


PLANT AND ANIMAL WHOLE GENOME SEQUENCING BEST PRACTICES



With Single Molecule, Real-Time (SMRT®) Sequencing and the Sequel™ System, you can easily and affordably generate high quality assemblies for even the most complex genomes. Users are regularly achieving plant and animal genome assemblies with megabase-size contig N50s and consensus accuracies >99.999%, resulting in the most complete genomes available today.


FROM GENOMIC DNA TO A COMPLETE GENOME



SAMPLE PREPARATION RECOMMENDATIONS

- Use recommend high-quality DNA input (>5 µg)
- Prepare long-insert SMRTbell™ library
 - 20 kb protocol¹ sufficient for most plant and animal genomes
 - >20 kb protocol² recommended for high-quality gDNA samples
 - >30 kb protocols^{3,4} available for genomes with long, complex repeats
- Size selection recommended
- Sequence to sufficient depth for accurate assembly; recommend:
 - 50-fold coverage for relatively homozygous diploids
 - >80-fold coverage for heterozygous or polyploid organisms
 - ~5 Gb generated per each SMRT Cell 1M on the Sequel System

SIZE SELECTION ENRICHES FOR LONGEST INSERTS




(A) Fragment Analyzer™ trace showing SMRTbell library sizes with no size-selection (black), a 4 kb selection (blue) or a 7 kb selection (red) performed with Blue Pippin™ (Sage Sciences). **(B)** Following SMRT Sequencing, the frequencies of reads spanning 7 kb regions increased by 38% for 4 kb size-selected libraries or 51% for 7 kb size-selected libraries⁵. Data provided by courtesy of Okinawa Institute of Advanced Sciences, Japan.

DATA ANALYSIS SOLUTIONS WITH SMRT® ANALYSIS AND PACBIO® DEVNET


- Exceptional *de novo* genome assemblies with megabase-size contig N50s
- High-quality consensus accuracies >99.999%
- Diploid- and polyploid-aware genome assembly; long reads provide phasing information for heterozygous genomes
- Open-source software
- Advanced data visualization and mining
- Assembly file in BAM and FASTA/Q formats allow for seamless input into downstream tools for further analysis

HIERARCHICAL GENOME ASSEMBLY PROCESS (HGAP)



HGAP utilizes all PacBio data using the longest reads for contiguity and all reads to generate high-quality *de novo* assemblies with high consensus accuracy (>QV50)⁶.

DIPLOID-AWARE GENOME ASSEMBLY PROCESS




In diploid-aware assembly, error-corrected reads are assembled using a string graph of the read overlaps, generating primary and alternative contigs that represent the alternative alleles, or structural variants (SVs), between the haplotypes⁷. Recommended for heterozygous, non-inbred genome assembly.

ACHIEVE HIGHLY CONTIGUOUS GENOME ASSEMBLIES


A Susie Reference sample



B Long-read assembly (Susie 3)



C Short-read assembly (gorGor3)



A new gorilla genome assembly was generated with SMRT Sequencing, resulting in a >150-fold improvement over previous assemblies. The rectangles in the treemaps in panels B and C are the largest contigs that cumulatively make up 300 Mb (~10%) of the assembly, showing the improved contiguity of the long-read assembly⁸.

KEY REFERENCES

1. Procedure and Checklist - 20 kb Template Preparation Using BluePippin™ Size-Selection System.
2. Procedure & Checklist >20 kb Template Preparation Using BluePippin™ Size-Selection System (15 - 20 kb Cutoff) for Sequel™ Systems.
3. Preparing >30 kb SMRTbell™ Libraries Using Needle Shearing and BluePippin™ Size-Selection System.
4. Procedure & Checklist - Preparing > 30 kb SMRTbell™ Libraries Using the Megaruptor® Shearing and BluePippin™ Size-Selection System.
5. Okinawa Institute of Advanced Sciences, Japan and Nippon Genetics (2014) Application Note: Long-fragment library preparation for PacBio RS P5-C3 Chemistry. Available at http://www.sagescience.com/wp-content/uploads/2014/10/application_note_2014_07_Eng.pdf
6. Chin, C.S. et al. (2013) Nonhybrid, finished microbial genome assemblies from long-read SMRT Sequencing data. *Nature Methods*. 10(6), 563-569.
7. Chin, C.S. et. al. (2016) Phased diploid genome assembly with single-molecule real-time sequencing. *Nature Methods*. 13(12), 1050-1054.
8. Gordon D. et al. (2016) Long-read sequence assembly of the gorilla genome. *Science*. 352 (6281), aae0344.

For Research Use Only. Not for use in diagnostic procedures. © Copyright 2017, Pacific Biosciences of California, Inc. All rights reserved. Information in this document is subject to change without notice. Pacific Biosciences assumes no responsibility for any errors or omissions in this document. Certain notices, terms, conditions and/or use restrictions may pertain to your use of Pacific Biosciences products and/or third party products. Please refer to the applicable Pacific Biosciences Terms and Conditions of Sale and to the applicable license terms at <http://www.pacb.com/legal-and-trademarks/terms-and-conditions-of-sale/>.

Pacific Biosciences, the Pacific Biosciences logo, PacBio, SMRT, SMRTbell, Iso-Seq, and Sequel are trademarks of Pacific Biosciences. BluePippin and SageELF are trademarks of Sage Science. NGS-go and NGSengine are trademarks of GenDx. Fragment Analyzer is a trademark of Advanced Analytical Technologies. All other trademarks are the sole property of their respective owners.



PN: BP102-011017