

## Introduction

Structural variants (genomic differences  $\geq 50$  base pairs) contribute to the evolution of organisms traits and human disease.




Figure 1. Common types of structural variation.

## Background

Most structural variants (SVs) are too small to detect with array comparative genomic hybridization but too large to reliably discover with short-read DNA sequencing. Recent studies in human genomes show that PacBio SMRT Sequencing sensitively detects structural variants<sup>1</sup>.




Figure 2. PacBio long reads have 5-fold increased sensitivity for structural variants compared to Illumina short reads.

While *de novo* assembly is the ideal method to identify variants in a genome, it requires high depth of coverage. Structural variant discovery using 10-fold coverage in humans analyzed with `pbsv` shows similar sensitivity for detecting variants.

Dataset	AK1 <sup>2</sup>	CHM <sup>3</sup>	NA12878 <sup>4</sup>
Analysis Method	<i>de novo</i> assembly	<i>de novo</i> assembly	<code>pbsv</code>
Fold Coverage	101-fold	41-fold	10-fold
Deletions (>50 bp)	7,358	6,111	8,209
Insertions (>50 bp)	10,077	9,638	11,350

Table 1. Structural variants in PacBio *de novo* human genome assemblies and low-coverage structural variants analysis.

## Sample Preparation




Figure 7. Recommendations for sample preparation and number of Sequel SMRT Cells for different genome sizes.

## Analysis: Map Reads, Chain Alignments, Call Variants

The analysis workflow to identify structural variants from low-coverage PacBio sequencing is: 1) map reads to the reference, 2) chain alignments, and 3) cluster indels to call variants.




Figure 3. NGM-LR is a read mapper designed for PacBio reads<sup>5</sup>. NGM-LR uses a convex gap penalty to model two sources of alignment gaps: biological indels and sequencing errors.




Figure 4. Comparison of read mapping spanning a deletion. (A) BWA, which uses a standard affine gap penalty, produces fragmented alignments at a deletion variant. (B) NGM-LR aligns the same PacBio reads with sharp boundaries at the deletion.




Figure 5. Chaining split alignments. Large gaps split NGM-LR alignments into primary and supplementary segments. Chaining connects collinear segments across large gaps. (A) A large deletion splits alignments of a read into two disjoint segments. (B) Chained alignments directly include a biological deletion, which simplifies visualization and variant calling.




Figure 6. Variant calling. To call structural variants from low-coverage sequencing, identify large deletion or insertion events in chained alignments, cluster nearby events that have similar length and sequence, and summarize into a call. Visualized in IGV v3.0 beta.

## Application of PBSV to *Drosophila*

PacBio RS II data<sup>6</sup> for *D. yakuba* mapped to v1.05 of Dyak genome with `pbsv align`. Variants called with `pbsv call` to generate BED and VCF formatted outputs.




Figure 8. Length Histograms of Insertions and Deletions in female of Cameroon strain CY21B3 (SRR1200825). Repeat Masker analysis identifies abundant 600-800bp SVs as Helitron repeats DNARep1\_Dyak and DNARep1\_DM.

## More X Insertions, Fewer Coding SVs

	INSERTIONS		DELETIONS	
	Observed	Expected	Observed	Expected
X	<b>434</b>	369	393	374
Autosomes	<b>1572</b>	<b>1638</b>	1638	1657
Coding	<b>81</b>	<b>374</b>	<b>52</b>	<b>380</b>
Noncoding	<b>1949</b>	<b>1656</b>	<b>2009</b>	<b>1681</b>

Table 2. Counts of SVs on the X vs Autosomes and Coding vs Noncoding sequence. `pbsv` BED file and Dyak v1.05 GFF file analyzed with `bedtools intersect`<sup>7</sup>. 2X<sup>2</sup> table values in **bold italic**;  $P < 0.05$ ,  $\chi^2$  test with Bonferonni correction. Expectations based on length of major chromosome arms or CDS versus non-CDS region lengths.

## Conclusion

- `pbsv align` uses the NGM-LR read mapper and alignment chaining to accurately map PacBio reads to a reference.
- `pbsv call` produces standard VCF and BED formats for custom downstream analysis.
- With yield of Sequel SMRT Cells 1M of 5-8 Gb, it is cost effective to screen populations for SVs using the `pbsv` workflow.
- Applications to humans and *Drosophila* identify thousands of insertions that cannot be detected with short-read technologies.

## References

- Huddleston J, et al. (2016). [Discovery and genotyping of structural variation from long-read haploid genome sequence data](#). *Genome Research*. doi:10.1101/gr.214007.116.
- Seo JS, et al. (2016). [De novo assembly and phasing of a Korean human genome](#). *Nature*. 538(7624), 243-247.
- Chaisson MJ, et al. (2015). [Resolving the complexity of the human genome using single-molecule sequencing](#). *Nature*. 517(7536), 608-611.
- Wenger AM, et al. (2016) <http://www.pacb.com/blog/identifying-structural-variants-na12878-low-fold-coverage-sequencing-pacbio-sequel-system/>
- Rescheneder P, et al. (2017). <https://github.com/phires/ngmlr>
- Rogers RL, et al. (2014). [Landscape of standing variation for tandem duplications in \*Drosophila yakuba\* and \*Drosophila simulans\*](#). *Mol Biol Evol*. 2014 31(7):1750-66.
- Quinlan AR and Hall IM. (2010). [BEDTools: a flexible suite of utilities for comparing genomic features](#). *Bioinformatics* 26(6), 841-842.

### Acknowledgements

The authors would like to thank Yuan Li, Luke Hickey, and Paul Peluso at PacBio.

### Conflict of Interest Statement

S.B.K. and A.M.W. are employees and shareholders of Pacific Biosciences, a company commercializing DNA sequencing technologies.