



# Technical overview – Kinnex kits for single-cell RNA, full-length RNA and 16S rRNA sequencing

Sequel II and IIe systems ICS v11.0

Revio system ICS v13.3

SMRT Link v25.3

PN 103-343-700 Rev 02 | December 2025

# Technical overview

## Kinnex kits for single-cell RNA, full-length RNA and 16S rRNA sequencing

1. Kinnex product introduction
2. Kinnex product configuration overview
3. Kinnex library preparation & sequencing workflow overview
4. Kinnex single-cell RNA workflow key highlights
5. Kinnex full-length RNA workflow key highlights
6. Kinnex 16S rRNA workflow key highlights
7. Technical documentation & applications support resources
8. APPENDIX 1 – Concatenating cDNA or gDNA amplicons using Kinnex kits to increase throughput
9. APPENDIX 2 – PacBio compatible Kinnex library preparation workflows
  - i. Preparing Kinnex single-cell libraries with Parse Evercode Whole Transcriptome kits
  - ii. Preparing Kinnex single-cell libraries with Twist exome enrichment kits
  - iii. Targeted enrichment of Kinnex full-length RNA libraries with IDT xGen hybridization and wash kit v3
10. APPENDIX 3 – SMRT Link run design procedure for Sequel IIe system



# Kinnex product introduction

# Kinnex kits for single-cell RNA, full-length RNA, and 16S rRNA sequencing

Kinnex kits offer scalable, cost-effective RNA sequencing solutions




## Kinnex single-cell RNA kit (103-072-200)

Single-cell RNA sequencing

4-plex using Kinnex barcoded adapters [16-fold concatenation]

100-120M reads – Revio system  
50-60M reads – Vega system  
30-40M reads – Sequel IIe system



Identify cell type-specific isoform expression




## Kinnex full-length RNA kit (103-072-100)

Full-length RNA sequencing

12-plex with the Iso-Seq express 2.0 kit<sup>1</sup> or 48-plex when combined with all four Kinnex barcoded adapters [8-fold concatenation]

50-60M reads – Revio system  
20-30M reads – Vega system  
16-20M reads – Sequel IIe system



Full-length transcripts from 1–10 kb for bulk RNA samples




## Kinnex 16S rRNA kit (103-072-000)

Full-length 16S rRNA for species identification


384-plex 16S amplicons or 1,536-plex when combined with all four Kinnex barcoded adapters [12-fold concatenation]

70-80M reads – Revio system  
30-40M reads – Vega system  
20-25M reads – Sequel IIe system




Taxonomic profiling for human, animal, and environmental samples

# Kinnex single-cell RNA kit for single-cell isoform sequencing




- 10x Chromium Single Cell 3' kit (v3.1) and 5' kit (v2)<sup>1</sup>
- 15–75 ng cDNA input
- 3,000 to 20,000 target cell recovery
- 2-day Kinnex library preparation using **Kinnex single-cell RNA kit**
- Barcoded Kinnex adapters support up to 4-plex multiplexing
- SMRT Link Run Design support for 'Kinnex single-cell RNA' application type option with auto-analysis (read segmentation + single-cell isoform analysis)
- SMRT Link single-cell Iso-Seq isoform-classification software to identify novel genes and isoforms
- Output compatible with tertiary single-cell analysis tools (e.g., *Seurat*, *Scanpy*, *Kana*)

# Kinnex full-length RNA kit for high-accuracy, full-length isoform sequencing



- Input 300 ng total RNA, RIN  $\geq 7$
- Generate up to 12-plex barcoded cDNA using **Iso-Seq express 2.0 kit (103-071-500)**
- 2-day Kinnex library preparation using **Kinnex full-length RNA kit (103-072-000)<sup>1</sup>**
- SMRT Link Run Design support for 'Kinnex full-length RNA' application type with auto-analysis (read segmentation + isoform analysis)
- SMRT Link Iso-Seq isoform-classification software to identify novel genes and isoforms with abundance information

# Kinnex 16S rRNA kit for full-length 16S sequencing



 **16S amplicon generation**

 **Kinnex library preparation**

 **PacBio sequencing**

 **Read segmentation + amplicon demultiplexing**


- Official protocol to generate barcoded 16S amplicons compatible with Kinnex 16S rRNA kit
- Protocol supports up to 384-plex multiplexing
- 2-day Kinnex library preparation using **Kinnex 16S rRNA kit (103-072-100)**
- SMRT Link Run Design support for 'Kinnex 16S rRNA' application type option with auto-analysis (read segmentation only)
- Demultiplex 16S amplicon barcodes in SMRT Link to generate per-sample read BAM files
- Analyze per-sample BAM files using GitHub tools or other custom 16S analysis pipeline



# Kinnex product configuration overview





# Kinnex product design overview and supported applications

Kinnex kits utilize the MAS-Seq concatenation method<sup>1</sup> to increase throughput on PacBio long-read systems



# Kinnex bundle kit product features

Kinnex bundle kit products are compatible with Sequel II/Ile, Vega and Revio systems

	Kinnex protocol	Kinnex bundle kit product	Concatenation factor	Sample multiplexing capacity			Throughput per SMRT Cell	SMRT Link support
				Amplicon multiplexing	Library multiplexing	Total multiplexing capacity		
Kinnex single-cell	Preparing Kinnex libraries using Kinnex single-cell RNA kit ( <a href="#">103-254-300</a> ) 	<b>Kinnex single-cell RNA kit</b> (103-072-200) 12 rxn 	16-fold	None	 4-plex	<b>4-plex</b>	<b>30 – 40 M</b> (Sequel Ile system) <b>50 – 60 M</b> (Vega system) <b>100 – 120 M</b> (Revio system)	<b>Full</b> Demux BC Read seg. Sc-Iso-Seq
Kinnex FL RNA	Preparing Kinnex libraries using the Kinnex full-length RNA kit ( <a href="#">103-238-700</a> ) 	<b>Kinnex full-length RNA kit</b> (103-072-000) 12 rxn 	8-fold	12-plex (Using Iso-Seq express 2.0 kit <sup>1</sup> ) 	 4-plex	<b>48-plex</b>	<b>16 – 20 M</b> (Sequel Ile system) <b>20 – 30 M</b> (Vega system) <b>50 – 60 M</b> (Revio system)	<b>Full</b> Demux BC Read seg. Iso-Seq
Kinnex 16S rRNA	Preparing Kinnex libraries from 16S rRNA amplicons ( <a href="#">103-238-800</a> ) 	<b>Kinnex 16S rRNA kit</b> (103-072-100) 12 rxn 	12-fold	384-plex (Using customer-supplied 16S PCR primers) 	 4-plex	<b>1,536-plex</b>	<b>20 – 25 M</b> (Sequel Ile system) <b>30 – 40 M</b> (Vega system) <b>70 – 80 M</b> (Revio system)	<b>Partial</b> Demux BC Read seg.

# Kinnex bundle kit product components

Kinnex bundle kit products are compatible with Sequel II/IIe, Vega and Revio systems

	Kinnex protocol	Kinnex bundle kit product	Kinnex bundle kit included components				Not included in Kinnex bundle kit but required for application-specific library prep workflow
Kinnex single-cell	<p>Preparing Kinnex libraries using Kinnex single-cell RNA kit (<a href="#">103-254-300</a>)</p> 	<p><b>Kinnex single-cell RNA kit</b> (103-072-200) 12 rxn</p> 	 Kinnex capture beads kit 103-076-000 (12 rxn)	 Kinnex single-cell concatenation kit 103-242-000 (12 rxn)	 Elution buffer 101-633-500	 SMRTbell cleanup beads 102-158-300	<p><b>10x Genomics Chromium Single Cell 3' or 5' reagent kits</b></p>
Kinnex FL RNA	<p>Preparing Kinnex libraries using the Kinnex full-length RNA kit (<a href="#">103-238-700</a>)</p> 	<p><b>Kinnex full-length RNA kit</b> (103-072-000) 12 rxn</p> 	 Kinnex PCR 8-fold kit 103-071-600 (12 rxn)	 Kinnex concatenation kit 103-071-800 (12 rxn)	 Elution buffer 101-633-500	 SMRTbell cleanup beads 102-158-300	 <p><b>Iso-Seq express 2.0 kit<sup>1</sup></b> 103-071-500 (24 rxn)</p>
Kinnex 16S rRNA	<p>Preparing Kinnex libraries from 16S rRNA amplicons (<a href="#">103-238-800</a>)</p> 	<p><b>Kinnex 16S rRNA kit</b> (103-072-100) 12 rxn</p> 	 Kinnex PCR 12-fold kit 103-071-700 (12 rxn)	 Kinnex concatenation kit 103-071-800 (12 rxn)	 Elution buffer 101-633-500	 SMRTbell cleanup beads 102-158-300	<p><b>16S gene amplification PCR primers (customer-supplied)<sup>2</sup></b></p>


<sup>1</sup> **Note:** Iso-Seq express 2.0 kit (103-071-500) is **not included** in Kinnex full-length RNA bundle kit and must be purchased separately from PacBio.

<sup>2</sup> Refer to *Procedure & checklist – Preparing Kinnex libraries from 16S rRNA amplicons* ([103-238-800](#)) for recommended 16S gene-specific forward and reverse PCR primer sequences to order.



# **Kinnex library preparation & sequencing workflow overview**

# Kinnex library preparation workflow overview<sup>1</sup>



# Kinnex library insert generation key workflow steps

## Single-cell cDNA sample preparation



Single-cell cDNA synthesis, amplification & cleanup



10x Forward (F) PCR primer

10x Reverse (R) PCR primer



Full-length single-cell 3' cDNA product (CBC = 10x cell barcode)

Min. cDNA input for TSO artifact removal step = 15 ng

OR



Full-length single-cell 5' cDNA product (CBC = 10x cell barcode)

Min. cDNA input for TSO artifact removal step = 15 ng

TSO artifact removal<sup>1</sup>



OR



Cleanup

Single-cell cDNA sample ready for Kinnex PCR step



## Full-length RNA sample preparation



Total RNA isolation & QC



Poly(A)+ mRNA template



Minimum total RNA input for cDNA synthesis = 300 ng

cDNA synthesis, amplification & cleanup

Can barcode up to 12 samples during cDNA amplification step using PacBio barcoded Fwd Iso-Seq primers

Forward (F) PCR primer (Barcoded Iso-Seq primer bcXX<sup>2</sup>)  
 BC TSO

Reverse (R) PCR primer (Iso-Seq cDNA amplification primer)



Barcoded full-length cDNA product

Pooling of barcoded cDNA samples & cleanup

If multiplexing, perform equal-mass pooling of barcoded cDNA samples


Pooled, full-length cDNA sample ready for Kinnex PCR step



## 16S RNA sample preparation



Bacterial metagenomic DNA isolation & QC



Recommended gDNA input for full-length 16S PCR amplification = 1-2 ng

PCR amplification of full-length 16S genes

Can barcode up to 384 samples during 16S gene amplification step using combinatorial dual-indexing with barcoded Fwd + Rev 16S gene-specific primers

Forward (F) PCR primer (Kinnex16S\_Fwd\_XX<sup>3</sup>)  
 Fwd BC

Reverse (R) PCR primer (Kinnex16S\_Rev\_XX<sup>3</sup>)  
 Rev BC



Dual-barcoded 16S amplicon product

Pooling of barcoded 16S PCR amplicons & cleanup

Perform equal-volume pooling of barcoded 16S samples

Pooled, full-length 16S amplicon sample ready for Kinnex PCR step



<sup>1</sup> Video tutorial is available: See *TSO artifact removal demo for PacBio MAS-Seq for 10x Single Cell 3' kit* [ Click to view the online video: [Link](#) ]


<sup>2</sup> Can use 12 different barcoded Iso-Seq Forward PCR primers to multiplex up to 12 samples per library prep Rx

<sup>3</sup> Can use 12 different 16S barcoded Forward PCR primers + 32 different 16S barcoded Reverse PCR primers to multiplex up to 384 samples per library prep Rx

# Application-specific Kinnex library insert generation workflow comparison

MAS-Seq for 10x single-cell 3'

Preparing MAS-Seq libraries using MAS-Seq for 10x Single Cell 3' kit (102-678-600)





Kinnex single-cell

Preparing Kinnex libraries using Kinnex single-cell RNA kit (103-254-300)




Kinnex FL RNA

Preparing Kinnex libraries using the Kinnex full-length RNA kit (103-238-700)



Kinnex 16S rRNA

Preparing Kinnex libraries from 16S rRNA amplicons (103-238-800)





# Kinnex PCR key workflow steps

Kinnex PCR step incorporates programmable segmentation adapter sequences into library insert samples

Application	Concatenation factor
Kinnex full-length RNA	8-fold
Kinnex 16S rRNA	12-fold
Kinnex single-cell RNA	16-fold

Set up 8/12/16 parallel PCR reactions per sample to generate amplified library inserts containing programmable sequences at both ends



Example Kinnex PCR setup for Kinnex single-cell RNA (16-fold concatenation factor)

Example: Amplified DNA products from Kinnex PCR 3 contain flanking segmentation adapter sequences that are **complementary** to the ends of DNA products from Kinnex PCR 2 & Kinnex PCR 4


Pool 8/12/16 Kinnex PCR products together and perform SMRTbell bead cleanup

Kinnex PCR-amplified library insert sample ready for Kinnex array formation step

# Application-specific Kinnex PCR workflow comparison

MAS-Seq for 10x single-cell 3'


Preparing MAS-Seq libraries using MAS-Seq for 10x Single Cell 3' kit (102-678-600)



~2 hrs

Kinnex single-cell

Preparing Kinnex libraries using Kinnex single-cell RNA kit (103-254-300)



~2 hrs

Kinnex FL RNA


Preparing Kinnex libraries using the Kinnex full-length RNA kit (103-238-700)



~2 hrs

Kinnex 16S rRNA

Preparing Kinnex libraries from 16S rRNA amplicons (103-238-800)



~1.5 hrs

Step	Temperature	Duration	Cycles
Initial Denaturation	98°C	3 min	1
Denaturation	98°C	20 s	
Annealing	68°C	30 s	9
Extension	72°C	4 min	
Final Extension	72°C	5 min	1
Hold	4°C	Hold	

Step	Temperature	Duration	Cycles
Initial Denaturation	98°C	3 min	1
Denaturation	98°C	20 s	
Annealing	68°C	30 s	9
Extension	72°C	4 min	
Final Extension	72°C	5 min	1
Hold	4°C	Hold	


Step	Temperature	Duration	Cycles
Initial Denaturation	98°C	3 min	1
Denaturation	98°C	20 s	
Annealing	68°C	30 s	9
Extension	72°C	4 min	
Final Extension	72°C	5 min	1
Hold	4°C	Hold	

Step	Temperature	Duration	Cycles
Initial Denaturation	98°C	3 min	1
Denaturation	98°C	20 s	
Annealing	68°C	30 s	9
Extension	72°C	90 s	
Final Extension	72°C	5 min	1
Hold	4°C	Hold	

Kinnex PCR thermal cycler program

# Kinnex array formation key workflow steps

Kinnex array formation step involves assembly of Kinnex PCR products (“segments”) containing a library insert + programmable “sticky ends” into a linear array



# Application-specific Kinnex array formation workflow comparison

MAS-Seq for 10x single-cell 3'


Preparing MAS-Seq libraries using MAS-Seq for 10x Single Cell 3' kit (102-678-600)



~4.5 hrs  
Final library

Kinnex single-cell

Preparing Kinnex libraries using Kinnex single-cell RNA kit (103-254-300)



~4.5 hrs  
Final library

Kinnex FL RNA

Preparing Kinnex libraries using the Kinnex full-length RNA kit (103-238-700)



~3 hrs  
Final library

Kinnex 16S rRNA

Preparing Kinnex libraries from 16S rRNA amplicons (103-238-800)



~3 hrs  
Final library


**Kinnex digestion + ligation performed as a single reaction** instead of two separate sequential reactions and SMRTbell bead cleanup is not required afterwards  
→ Kinnex full-length RNA & 16S rRNA protocols are more **streamlined**

<sup>1</sup> Kinnex digestion / ligation reaction does not require addition of ligation additive.

# Kinnex SMRTbell adapter design

Kinnex terminal adapters incorporate barcode sequences to enable up to 4-plex sample multiplexing at the library level


- Kinnex adapters contain **barcode sequences<sup>1</sup>** to enable (optional) sample multiplexing at the SMRTbell library level (**up to 4-plex**)
  - Forward and reverse Kinnex adapter pairs are pre-mixed in Kinnex concatenation kits
  - Kinnex concatenation kits contain a total of **4 barcoded Kinnex adapter mixes (bcM0001-bcM0004)** to enable multiplexing of up to **4 samples per SMRT Cell**



# Kinnex SMRTbell adapter design (cont.)

Kinnex terminal adapters use a new design that enables improved HiFi sequencing performance

- Kinnex adapters enable longer polymerase read length performance and improved sample loading efficiency
  - Improved overall HiFi data yield




- Kinnex adapter design requires a **different sequencing primer (Kinnex sequencing primer<sup>1</sup>)** for annealing reaction



# Application-specific Kinnex library prep input & expected step yield comparison


MAS-Seq for 10x single-cell 3'

Preparing MAS-Seq libraries using MAS-Seq for 10x Single Cell 3' kit (102-678-600)




Kinnex single-cell

Preparing Kinnex libraries using Kinnex single-cell RNA kit (103-254-300)




Kinnex FL RNA

Preparing Kinnex libraries using the Kinnex full-length RNA kit (103-238-700)





Kinnex 16S rRNA

Preparing Kinnex libraries from 16S rRNA amplicons (103-238-800)




# Example Femto Pulse DNA sizing QC results for final Kinnex SMRTbell libraries


**Kinnex single-cell RNA**  
(16-fold target concatenation)



**Kinnex full-length RNA**  
(8-fold target concatenation)



**Kinnex 16S rRNA**  
(12-fold target concatenation)



Kinnex library type	Kinnex library size range
Kinnex single-cell RNA	~10 – 17 kb
Kinnex full-length RNA	~12 – 16 kb
Kinnex 16S rRNA	~17 – 18 kb




# Kinnex sequencing preparation workflow overview

# Sample Setup & Run Design recommendations for Kinnex libraries

Follow SMRT Link run design instructions to prepare Kinnex libraries for sequencing

Select desired Kinnex application type in SMRT Link run design page



After specifying your application type, SMRT Link auto-fills selected Run Design parameter fields with default recommended values

## Kinnex single-cell RNA<sup>1</sup>

### Whole genome sequencing

- Human WGS
- Microbial assembly
- Other WGS

### RNA sequencing

- Iso-Seq method
- MAS-Seq single cell
- Kinnex single-cell RNA**
- Kinnex full-length RNA

### Viral sequencing

- Adeno-associated virus

### Metagenomics

- Kinnex 16S rRNA
- Full-length 16S rRNA sequencing

## Kinnex full-length RNA<sup>2</sup>

### Whole genome sequencing

- Human WGS
- Microbial assembly
- Other WGS

### RNA sequencing

- Iso-Seq method
- MAS-Seq single cell
- Kinnex single-cell RNA
- Kinnex full-length RNA**

### Viral sequencing

- Adeno-associated virus

### Metagenomics

- Kinnex 16S rRNA
- Full-length 16S rRNA sequencing

## Kinnex 16S rRNA<sup>2</sup>

### Whole genome sequencing

- Human WGS
- Microbial assembly
- Other WGS

### RNA sequencing

- Iso-Seq method
- MAS-Seq single cell
- Kinnex single-cell RNA
- Kinnex full-length RNA

### Viral sequencing

- Adeno-associated virus

### Metagenomics

- Kinnex 16S rRNA**
- Full-length 16S rRNA sequencing

SMRT Link module	Key setup parameters	Sequel II system recommended settings	Vega system recommended settings	Revio system recommended settings
Sample setup	Library type	Kinnex	Kinnex	Kinnex
	Primer	Kinnex sequencing primer	Kinnex sequencing primer	Kinnex sequencing primer
	Binding/Polymerase kit	Sequel II binding kit 3.2	Vega polymerase kit	Revio SPRQ polymerase kit
	Concentration on plate	40 – 60 pM	130 – 160 pM	130 – 160 pM
Runs → Run design	Use Adaptive Loading	YES	N/A	YES
	Library Type <sup>1</sup>	Kinnex	Kinnex	Kinnex
	Movie Acquisition Time	30 hrs	24 hrs	24 hrs
	Add Analysis	YES	YES	YES

# SMRT Link Run Design procedure for **Revio and Vega systems**

## Sample and run information

### Kinnex single-cell RNA

### Kinnex full-length RNA

### Kinnex 16S rRNA

▼ Plate 1, Well A01: Kinnex single-cell RNA library demo

▼ Plate 1, Well A01: Kinnex full-length RNA library demo

▼ Plate 1, Well A01: Kinnex 16S rRNA library demo

Application <small>Required</small>	Kinnex single-cell RNA	Application <small>Required</small>	Kinnex full-length RNA	Application <small>Required</small>	Kinnex 16S rRNA
Plate Well <small>Required</small>	Plate 1, Well A01	Plate Well <small>Required</small>	Plate 1, Well A01	Plate Well <small>Required</small>	Plate 1, Well A01
Well Name <small>Required</small>		Well Name <small>Required</small>		Well Name <small>Required</small>	
Well Comment		Well Comment		Well Comment	
Library Type <small>Required</small>	Kinnex	Library Type <small>Required</small>	Kinnex	Library Type <small>Required</small>	Kinnex
Insert Size (bp) <small>Required</small>	15000	Insert Size (bp) <small>Required</small>	16000	Insert Size (bp) <small>Required</small>	17000
Library Concentration (pM) <small>Required</small>		Library Concentration (pM) <small>Required</small>		Library Concentration (pM) <small>Required</small>	
Movie Acquisition Time (hours)	24	Movie Acquisition Time (hours)	24	Movie Acquisition Time (hours)	24

Select application type to autofill Library Type & Movie Acquisition Time recommended settings

Default = Kinnex library type (instead of Standard or AAV)<sup>1</sup>

### Standard SMRTbell library type containing standard SPK3 barcoded terminal adapters



Forward and reverse standard terminal adapters have the same structure

### Kinnex SMRTbell library type containing Kinnex barcoded terminal adapters



Forward and reverse Kinnex terminal adapters have different structures

**Library Type** field determines which adapter finding algorithm is used during post-primary analysis<sup>1</sup>

<sup>1</sup> **Note:** When sequencing a Kinnex library sample, if 'Standard' library type is mistakenly selected instead of 'Kinnex' then a higher missing adapter rate (> 95%) and a slight degradation in barcode demultiplexing performance (~93-96% barcoded HiFi read yield) will be observed.

# SMRT Link Run Design procedure for **Revio and Vega systems (cont.)**

## Sample indexing (barcoding) information

### Kinnex single-cell RNA

### Kinnex full-length RNA

### Kinnex 16S rRNA

**Samples**

Sample is indexed  YES  NO

**Indexes** Required: MAS adapter indexes

**Biosample names** Required: Interactively From a File

**Samples**

Sample is indexed  YES  NO

**Indexes** Required: MAS adapter indexes

**Biosample names** Required: Interactively From a File

**Samples**


Sample is indexed  YES  NO

**Indexes** Required: MAS adapter indexes

**Biosample names** Required: Interactively From a File

Default Indexes FASTA = MAS adapter indexes

Example complete Kinnex library molecule containing barcoded Kinnex terminal MAS adapters<sup>1</sup> at both ends



#### Barcode Selector and Sample Name Editor

Available Barcodes	Included Barcodes				
<input type="checkbox"/> bcM0001--bcM0001 <input type="checkbox"/> bcM0002--bcM0002 <input type="checkbox"/> bcM0003--bcM0003 <input type="checkbox"/> bcM0004--bcM0004	<table border="1"> <thead> <tr> <th>Barcode ID</th> <th>Bio Sample ID</th> </tr> </thead> <tbody> <tr> <td> </td> <td> </td> </tr> </tbody> </table>	Barcode ID	Bio Sample ID		
Barcode ID	Bio Sample ID				

#### Barcode Selector and Sample Name Editor

Available Barcodes	Included Barcodes						
<input type="checkbox"/> bcM0003--bcM0003 <input type="checkbox"/> bcM0004--bcM0004	<table border="1"> <thead> <tr> <th>Barcode ID</th> <th>Bio Sample ID</th> </tr> </thead> <tbody> <tr> <td>bcM0001--bcM0001</td> <td>Kinnex adapter-barcoded library 1</td> </tr> <tr> <td>bcM0002--bcM0002</td> <td>Kinnex adapter-barcoded library 2</td> </tr> </tbody> </table>	Barcode ID	Bio Sample ID	bcM0001--bcM0001	Kinnex adapter-barcoded library 1	bcM0002--bcM0002	Kinnex adapter-barcoded library 2
Barcode ID	Bio Sample ID						
bcM0001--bcM0001	Kinnex adapter-barcoded library 1						
bcM0002--bcM0002	Kinnex adapter-barcoded library 2						

SMRT Link



#### MAS SMRTbell barcoded adapter indexes

```
>bcM0001
ACAGTC
>bcM0002
ATGACG
>bcM0003
CACGTG
>bcM0004
CATCGC
```

<sup>1</sup> Four barcoded terminal Kinnex adapters (Kinnex adapter bcM0001-bcM0004) are available for Kinnex array formation step. Kinnex adapter barcode sequences can be downloaded from [SMRT Link](#) Data Management module.

# SMRT Link Run Design procedure for **Revio and Vega systems (cont.)**

## Data options

### Kinnex single-cell RNA

### Kinnex full-length RNA

### Kinnex 16S rRNA

Default = **NO** for  
**Include Base Kinetics**

▼ Data Options

Include Base Kinetics  YES  NO

Consensus Mode  MOLECULE  STRAND

Assign Data To Project

Default Consensus  
Mode = **MOLECULE**<sup>1</sup>

▼ Data Options

Include Base Kinetics  YES  NO

Consensus Mode  MOLECULE  STRAND

Assign Data To Project

Can leave **Include Base Kinetics** and **Consensus Mode**  
fields at their **default settings** for Kinnex library samples

▼ Data Options

Include Base Kinetics  YES  NO

Consensus Mode  MOLECULE  STRAND

Assign Data To Project

# SMRT Link Run Design analysis options

## Kinnex single-cell RNA

## Kinnex full-length RNA

## Kinnex 16S rRNA

Analysis Options **Default = YES for Add Analysis**

Analysis Options

Analysis Options

Add Analysis  YES  NO

Analysis Name Required: Kinnex\_Single-Cell\_RNA\_Demo\_Analysis\_Job\_Name

Select Analysis Workflow Required: Read Segmentation and Single-Cell Iso-Seq

Segmentation Adapter Set: MAS-Seq Adapter v1 (MAS16)

Primer Set Required: 10x Chromium single cell 3' cDNA primers

Reference Set Required: Human Genome hg38, with Gencode v39 annotations

Kit Type  10X 3' KIT  10X 5' KIT

Advanced Parameters

Add Analysis  YES  NO

Analysis Name Required: Kinnex\_Full-Length\_RNA\_Demo\_Analysis\_Job\_Name

Select Analysis Workflow Required: Read Segmentation and Iso-Seq

Segmentation Adapter Set: MAS-Seq Adapter v3 (MAS8)

Primer Set Required: Iso-Seq v2 Barcoded cDNA Primers

Reference Set: Human Genome hg38, with Gencode v39 annotations

Cluster of Barcoded Samples

Advanced Parameters

Add Analysis  YES  NO

Analysis Name Required: Kinnex\_Full-Length\_16S\_rRNA\_Demo\_Analysis\_Job\_Name

Select Analysis Workflow Required: Read Segmentation

Segmentation Adapter Set: MAS-Seq Adapter v2 (MAS12)

Advanced Parameters


Analysis Workflow is automatically filled in

Split arrayed HiFi reads at adapter positions, generating segmented reads (S-reads)


### Read Segmentation and Single-Cell Iso-Seq

### Read Segmentation and Iso-Seq

### Read Segmentation



Isoform-classification analysis to identify novel genes & isoforms



Isoform-classification analysis to identify novel genes & isoforms with abundance information (raw counts & normalized counts per million)

# SMRT Link Run Design analysis options (cont.)

## Kinnex single-cell RNA

**Analysis Options**

Add Analysis  YES  NO

Analysis Name Required: Kinnex\_Single-Cell\_RNA\_Demo\_Analysis\_Job\_Name

Select Analysis Workflow Required: Read Segmentation and Single-Cell Iso-Seq

Segmentation Adapter Set: MAS-Seq Adapter v1 (MAS16)

Primer Set Required: 10x Chromium single cell 3' cDNA primers

Reference Set Required: Human Genome hg38, with Gencode v39 annotations

Kit Type  10X 3' KIT  10X 5' KIT

Advanced Parameters

## Kinnex full-length RNA

**Analysis Options**

Add Analysis  YES  NO

Analysis Name Required: Kinnex\_Full-Length\_RNA\_Demo\_Analysis\_Job\_Name

Select Analysis Workflow Required: Read Segmentation and Iso-Seq

Segmentation Adapter Set: MAS-Seq Adapter v3 (MAS8)

Primer Set Required: Iso-Seq v2 Barcoded cDNA Primers

Reference Set: Human Genome hg38, with Gencode v39 annotations

Cluster of Barcoded Samples  Pool reads and cluster together

Advanced Parameters

## Kinnex 16S rRNA

**Analysis Options**

Add Analysis  YES  NO


Analysis Name Required: Kinnex\_Full-Length\_16S\_rRNA\_Demo\_Analysis\_Job\_Name

Select Analysis Workflow Required: Read Segmentation

Segmentation Adapter Set: MAS-Seq Adapter v2 (MAS12)

Advanced Parameters

Specify Segmentation Adapter Set that corresponds to the Kinnex library concatenation method used



# SMRT Link Run Design analysis options (cont.)

## Kinnex single-cell RNA

## Kinnex full-length RNA

**Analysis Options**

Add Analysis  YES  NO

Analysis Name Required: Kinnex\_Single-Cell\_RNA\_Demo\_Analysis\_Job\_Name

Select Analysis Workflow Required: Read Segmentation and Single-Cell Iso-Seq

Segmentation Adapter Set: MAS-Seq Adapter v1 (MAS16)

Primer Set Required: 10x Chromium single cell 3' cDNA primers

Reference Set Required: Human Genome hg38, with Gencode v39 annotations

Kit Type  10X 3' KIT  10X 5' KIT

Advanced Parameters

**Analysis Options**

Add Analysis  YES  NO

Analysis Name Required: Kinnex\_Full-Length\_RNA\_Demo\_Analysis\_Job\_Name

Analysis Workflow Required: Read Segmentation and Iso-Seq

Segmentation Adapter Set: MAS-Seq Adapter v3 (MAS8)

Primer Set Required: Iso-Seq v2 Barcoded cDNA Primers

Reference Set: Human Genome hg38, with Gencode v39 annotations

Cluster of Barcoded Samples  Pool reads and cluster together

Advanced Parameters

Specify Primer Set used for cDNA amplification

Specify Reference Set. Default sets are:

- Human Genome hg38, with Gencode v39 annotations
- Mouse Genome mm39, with Gencode vM28 annotations

Specify reference genome & annotation sets to align high quality isoforms to, and to collapse isoforms mapped to the same genomic loci.

10x Forward (F) PCR primer (pink)      10x Reverse (R) PCR primer (blue)

5' [CBC][UMI] TTTTTTTTTT LIBRARY INSERT CCC TSO 3'  
 3' [CBC][UMI] AAAAAAAAAA LIBRARY INSERT GGG TSO 5'

For Kinnex single-cell 3' RNA analysis, select '10x Chromium single cell 3' cDNA primers'

OR

5' [CBC][UMI][TSO]GGG LIBRARY INSERT AAAAAAAAAA 3'  
 3' [CBC][UMI][TSO]CCC LIBRARY INSERT TTTTTTTTTT 5'

For Kinnex single-cell 5' RNA analysis, select '10x Chromium single cell 5' cDNA primers'

Forward (F) PCR primer (Barcoded Iso-Seq primer bcXX<sup>2</sup>) (pink)      Reverse (R) PCR primer (Iso-Seq cDNA amplification primer) (blue)

5' BC TSO GGG LIBRARY INSERT AAAAAAAAAA 3'  
 3' BC TSO CCC LIBRARY INSERT TTTTTTTTTT 5'

For Kinnex full-length RNA analysis, select 'Iso-Seq v2 Barcoded cDNA primers'

Specify primer sequence file in FASTA format to identify cDNA primers for removal (include the 5' and 3' cDNA primers)



# SMRT Link Run Design analysis options (cont.)

## Kinnex single-cell RNA

**Analysis Options**

Add Analysis  YES  NO

Analysis Name Required: Kinnex\_Single-Cell\_RNA\_Demo\_Analysis\_Job\_Name

Select Analysis Workflow Required: Read Segmentation and Single-Cell Iso-Seq

Segmentation Adapter Set: MAS-Seq Adapter v1 (MAS16)

Primer Set Required: 10x Chromium single cell 3' cDNA

Reference Set Required: Human Genome hg38, with Gencode v39 annotations

Kit Type  10x 3' KIT  10x 5' KIT

Advanced Parameters

Specify Kit Type used for single-cell cDNA generation

Specification of Kit Type (10x 3' Kit or 10x 5' Kit) determines which set of 10x barcode sequences to use, and also affects UMI and single-cell barcode design settings<sup>1</sup>

5' [CBC][UMI]TTTTTTTTT LIBRARY INSERT CCC TSO 3'  
 3' [CBC][UMI]AAAAAAAAA LIBRARY INSERT GGG TSO 5'

For Kinnex single-cell 3' RNA analysis, select '10x 3' Kit'

OR

5' [CBC][UMI][TSO]GGG LIBRARY INSERT AAAAAAAAAA 3'  
 3' [CBC][UMI][TSO]CCC LIBRARY INSERT TTTTTTTTTT 5'

For Kinnex single-cell 5' RNA analysis, select '10x 5' Kit'

## Kinnex full-length RNA

**Analysis Options**

Add Analysis  YES  NO

Analysis Name Required: [blank]

Select Analysis Workflow Required: [blank]

Segmentation Adapter Set: MAS-Seq Adapter v1 (MAS8)

Primer Set Required: 10x 5' Barcoded cDNA Primers

Reference Set Required: Human Genome hg38, with Gencode v39 annotations

Cluster of Barcoded Samples  Pool reads and cluster together

Advanced Parameters

Specify how to perform read clustering for barcoded samples

Pool reads and cluster together



-- select --

Cluster reads separately


Pool reads and cluster together

Specification of Cluster of Barcoded Samples setting determines whether all FLNC reads will be pooled for clustering (Does not apply to non-barcoded samples.)


## Clustering options for Kinnex full-length RNA samples



### Analyze pooled samples



### Analyze by sample (default)



# SMRT Link Read Segmentation and Iso-Seq analysis video demonstration

Video demonstration of SMRT Link Read Segmentation and Iso-Seq application workflow for analysis of Kinnex full-length RNA samples

**Analysis Application** Required

Read Segmentation and Iso-Seq

Import Analysis Settings | Export

**Associated Inputs**

**Segmentation Adapter Set**

MAS-Seq Adapter v3 (MAS8)

**Primer Set** Required

Iso-Seq v2 Barcoded cDNA Primers

**Reference Set**

Human Genome hg38, with Gencode v39 annotations

**Cluster of Barcoded Samples**

Pool reads and cluster together

Advanced Parameters

**Analysis Name**

test

**Analysis Datasets**

Displaying rows 1 to 1 out of 1

ID	Name
21...	3230211_KPoS_64007_...

**Demo video for Read Segmentation and Iso-Seq workflow (SMRT Link v13.0+)**

- Workflow supports full-length isoform analysis for data generated on PacBio long-read systems using **Kinnex full-length RNA kit**
- End-to-end workflow begins with HiFi reads and **outputs full-length isoform classifications with supporting read count information**



# SMRT Link Read Segmentation and Demultiplex Barcodes video demonstration

Video demonstration of SMRT Link Read Segmentation and Demultiplex Barcodes workflow for analysis of Kinnex 16S rRNA samples

**Data Utility** Required

Read Segmentation

Import Analysis Settings Export

**Associated Inputs**

Segmentation Adapter Set

MAS-Seq Adapter v2 (MAS12)

Advanced Parameters

**Analysis Name**

test-ReadSeg

**Analysis Datasets**

Displaying rows 1 to 1 out of 1

ID	Name
21...	20231020-4_84028_13...

[Demo video](#) for Read Segmentation and Demultiplex Barcodes workflow (SMRT Link v13.0+)

- Demo video for analyzing Kinnex 16S rRNA data generated using Kinnex 16S rRNA kit in SMRT Link v13.0 and up






# **Kinnex single-cell RNA library preparation & sequencing workflow key highlights**

# Kinnex single-cell RNA library preparation procedure description

Procedure & checklist – Preparing Kinnex libraries using Kinnex single-cell RNA kit ([103-254-300](https://www.pacb.com/documentation/103-254-300)) describes the workflow for constructing single-cell Kinnex libraries from 10x Chromium 3' or 5' single cell cDNA using the **Kinnex single-cell RNA kit<sup>1</sup>** for library preparation and sequencing on PacBio long-read systems.

Overview	
Samples per kit	12
Workflow time	3 days for up to 12 samples
Number of SMRT® Cells per Kinnex library preparation	>8 SMRT Cells for Revio® using SPRQ™ chemistry >2 SMRT Cells for the Vega® or Revio® (non-SPRQ) systems >4 SMRT Cells for Sequel® II/IIE®

cDNA input	
Quantity	>15 ng 10x cDNA per sample cDNA concentration should be >1ng/μL with up to 15 μL in volume. See <a href="#">step 2.1</a> for 10x cDNA input requirement.
Average segment lengths	500–1,100 bp
Average 16-segment array lengths	10–17 kb



Preparing Kinnex™ libraries using Kinnex single-cell RNA kit PacBio

Procedure & checklist

[Before you begin](#)

This procedure describes the workflow for constructing single-cell Kinnex libraries from 10x Chromium 3' or 5' cDNA using the Kinnex single-cell RNA kit for library prep and sequencing on PacBio® Sequel® II, Sequel IIE, Vega™, and Revio® systems.

This kit is intended for use with single-cell cDNA generated using the 10x Chromium Next GEM Single Cell 3' kit (v3.1 or v4) or 10x Chromium Next GEM Single Cell 5' kit (v2, v3). For a full list of compatibility with 10x products, please consult the [10x website](#).

Overview	
Samples per kit	12
Workflow time	3 days for up to 12 samples
Number of SMRT® Cells per Kinnex library preparation	>8 SMRT Cells for Revio® using SPRQ™ chemistry >2 SMRT Cells for the Vega® or Revio® (non-SPRQ) systems >4 SMRT Cells for Sequel® II/IIE®

cDNA input	
Quantity	>15 ng 10x cDNA per sample cDNA concentration should be >1ng/μL with up to 15 μL in volume. See <a href="#">step 2.1</a> for 10x cDNA input requirement.
Average segment lengths	500–1,100 bp
Average 16-segment array lengths	10–17 kb

© 2025 PacBio. All rights reserved. Research use only. Not for use in diagnostic procedures.  
103-254-300 REV08 OCT2025 PacBio

- Kinnex single-cell RNA library prep protocol uses **Kinnex single-cell RNA kit**
- **Do not use** SMRTbell prep kit 3.0 with this protocol

PacBio [Documentation](https://www.pacb.com/documentation/103-254-300) (103-254-300)

# Comparison of *Kinnex single-cell RNA kit vs. MAS-Seq for 10x Single Cell 3' kit*



	Kinnex single-cell RNA kit (103-072-200)	MAS-Seq for 10x Single Cell 3' kit (102-659-600)
# reactions	12	8
Concatenation		16-fold
Compatibility	10x Chromium Single Cell 3' kit (v3.1) <b>10x Chromium Single Cell 5' kit (v2)</b>	10x Chromium Single Cell 3' kit (v3.1)
Sample multiplexing support	<b>4-plex using Kinnex barcoded adapters</b>	None
Workflow time		2 days
SMRT Link support	Yes (since SMRT Link v13.1)	Yes (since SMRT Link v11.1)
Sequencing primer	<b>Kinnex sequencing primer</b>	Standard sequencing primer
Sequencing OPLC	<b>40-60 pM (Sequel II/IIe system)</b> <b>130-160 pM (Vega/Revio system)</b>	85 pM (Sequel II/IIe system) 225 pM (Revio system)
Run time		Sequel II/IIe system: 30hr movie with adaptive loading Revio/Vega system: 24hr movie with adaptive loading
S-read yield		~30-40 M reads (Sequel II/IIe system) 50-60 M reads (Vega system) ~100-120 M reads (Revio system)


# Kinnex single-cell RNA experimental design considerations

## Kinnex single-cell RNA application use case recommendations for PacBio systems

	Sequel II and Ile systems	Vega system	Revio system + SPRQ chemistry
Experimental goal	Characterize alternative splicing in single cells / cell types		
Transcript reads per SMRT Cell	30-40 M reads	50-60 M reads	100-120 M reads
Sample multiplexing <sup>1</sup>	Not recommended	Not recommended	Up to 2 samples per Revio SMRT Cell (2-plex)
Cell input into 10x Chromium single cell 3' or 5' cDNA generation workflow	3,000 – 10,000 cells for running a single (non-multiplexed) sample on one Sequel II SMRT Cell 8M	3,000 – 10,000 cells for running a single (non-multiplexed) sample on one Vega SMRT Cell	<10,000 cells per sample if multiplexing 2 samples per Revio SMRT Cell (2-plex) 10,000 – 20,000 cells per sample if running a single (non-multiplexed) sample on 1–2 Revio SMRT Cells
Expected coverage	Obtain 3,000 – 10,000 unique reads/single cell	Obtain 3,000 – 10,000 unique reads/single cell	Obtain up to ~20,000 unique reads/single cell
Kinnex library prep protocol	Procedure & checklist – Preparing Kinnex libraries using Kinnex single-cell RNA kit (103-254-300)		
Single-cell cDNA input into Kinnex library prep workflow	15-75 ng of cDNA from (1) 10x Chromium Single Cell Universal 3' or 5' kit; Single Cell Next GEM (v3.1 3' or v2 5'); (2) 10x Visium HD Spatial Gene Expression kit; (3) Parse Evercode WT or WT mini kit; or (4) ArgenTag Single-Cell RNA Library Kit		
SMRT Link data analysis workflows	Read Segmentation and Single-cell Iso-Seq Analysis		
Community data analysis tools	Annotation & quantification: SQANTI3 / Differential analysis: TappAS / Fusion calling: pbfusion / Visualization: SWAN		


# Kinnex single-cell RNA method overview

## Single-cell cDNA sample preparation\*



Min. cDNA input for TSO artifact removal step = 15 ng

## Kinnex library prep, sequencing & analysis



**Preparing Kinnex™ libraries using Kinnex single-cell RNA kit**

**Before you begin**

This procedure describes the workflow for constructing single-cell Kinnex libraries from 10x Chromium 3' or 5' cDNA using the Kinnex single-cell kit for library prep and sequencing on PacBio® Sequel® II, Sequel IIe, Vega™, and Revio™ systems.

This kit is intended for use with amplified cDNA generated using the 10x Chromium Next Gen Single Cell 3' v2 (10x) or 10x Chromium Next Gen Single Cell 3' v2 (10x) for a full set of compatibility with the products, please consult the 10x website.

Parameter	Value
Samples per kit	12
Workflow time	2 days for up to 12 samples
Number of cDNA* Cells per Kinnex library preparation	10,000 cells for the 10x Chromium Next Gen Single Cell 3' v2 (10x) or 10x Chromium Next Gen Single Cell 3' v2 (10x) for a full set of compatibility with the products, please consult the 10x website.
Other notes	10x (10x cDNA per sample) cDNA concentration should be 100-1,000 ng/μL with up to 10% in volume, see also 10x website for 10x cDNA input requirement.
Quantity	500 - 1,000 ng
Average segment length	100 - 1,000 nt
Average fragment length	100 - 1,000 nt

**Procedure & checklist (103-254-300)**

\* Refer to [10x Genomics Support](https://www.10xgenomics.com/support/single-cell/next-gen/single-cell-3) website to download 10x Chromium single cell 3'/5' reagent kit user guides.

<sup>2</sup> Kinnex adapter barcode sequences can be downloaded from [SMRT Link](https://www.pacb.com/kinnex-data-management/) Data Management module.

# Kinnex single-cell RNA library preparation procedural notes

## I. Library insert generation



### TSO PCR

- Set up TSO PCR reactions **ON ICE** using the **CORRECT** TSO PCR primer set (3' 10x cDNA → 3' TSO PCR primer; 5' 10x cDNA → 5' TSO PCR primer)
  - PCR polymerase 3'→5' exonuclease activity negatively impacts amplification yield if prepared at room temp.

### TSO removal [ [Video demonstration](#) ]

- **IMPORTANT:** For bead capture steps, allow enough time for beads to magnetize as binding buffer is highly viscous; pipette mix with care and **avoid generating bubbles**
  - Wide-bore pipette tips are recommended to help minimize foaming (specifically when resuspending Kinnex capture beads)
- **Fully resuspend beads** during all wash steps to remove artifact cDNA effectively
- **Fully resuspend beads** before DNA quantification using Qubit dsDNA HS assay

### SMRTbell bead cleanup

- Prior to Kinnex array formation, perform a **1.5X** SMRTbell bead cleanup
  - If the cDNA contains smaller fragments <200 bp, it is recommended to increase the SMRTbell cleanup bead ratio to 1.8 – 2.0X

## III. Kinnex array formation



### Kinnex digestion

- Recommended input amount to proceed with Kinnex array formation is **5 µg** of Kinnex PCR amplicons
  - Proceeding with <3 µg is not recommended since lower input amounts may lead to insufficient final library yields to enable optimal sequencing results

## II. Kinnex PCR




### Kinnex PCR

- Set up Kinnex PCR reactions **ON ICE**
  - PCR polymerase 3'→5' exonuclease activity negatively impacts amplification yield if prepared at room temp.
- Perform Kinnex PCR using **25 ng** of cDNA after TSO removal step

# Example Kinnex single-cell RNA library preparation QC results


Kinnex single-cell 3' RNA library prepared with human cDNA

## Kinnex PCR DNA sizing QC (Single-cell 3' cDNA)



Example Bioanalyzer DNA sizing QC analysis results for Kinnex PCR products generated for a 10x Chromium single-cell 3' cDNA samples prepared from a human cell line (HG002).

## Final Kinnex single-cell RNA library QC



Example Femto Pulse DNA sizing QC analysis results for final Kinnex single-cell RNA library.

Final Kinnex library yield is typically sufficient to load the following number of SMRT Cells:

- >8 SMRT Cells for Revio using SPRQ chemistry
- >2 SMRT Cells for the Vega or Revio (non-SPRQ) systems
- >4 SMRT Cells for Sequel II/IIe systems


10x single cell 3' cDNA input	15 ng
cDNA input for Kinnex array formation	5000 ng
Post-nuclease treatment & final library cleanup yield (%) <sup>1</sup>	1100 ng (22.0%)

<sup>1</sup> Post-nuclease treatment & final cleanup yields typically ranged from ~10% to ~25% when using single-cell 3' cDNA samples for Kinnex single-cell RNA library construction.

# Example Kinnex single-cell RNA library preparation QC results (cont.)


Kinnex single-cell 5' RNA library prepared with human cDNA

## Kinnex PCR DNA sizing QC (Single-cell 5' cDNA)



Example Bioanalyzer DNA sizing QC analysis results for Kinnex PCR products generated for a 10x Chromium single-cell 5' cDNA samples prepared from a human cell line (HG002).

## Final Kinnex single-cell RNA library QC



Example Femto Pulse DNA sizing QC analysis results for final Kinnex single-cell RNA library.

Final Kinnex library yield is typically sufficient to load the following number of SMRT Cells:

- >8 SMRT Cells for Revio using SPRQ chemistry
- >2 SMRT Cells for the Vega or Revio (non-SPRQ) systems
- >4 SMRT Cells for Sequel II/IIe systems


10x single cell 5' cDNA input	15 ng
cDNA input for Kinnex array formation	5000 ng
Post-nuclease treatment & final library cleanup yield (%) <sup>1</sup>	1008 ng (20.2%)

<sup>1</sup> Post-nuclease treatment & final cleanup yields typically ranged from ~10% to ~25% when using single-cell 5' cDNA samples for Kinnex single-cell RNA library construction.


# Example sequencing performance for Kinnex single-cell RNA libraries prepared with human cDNA [ **Revio system + SPRQ chemistry<sup>1</sup>** ]

Kinnex single-cell RNA 5' library for PBMC single-cell cDNA sample from 10x Chromium Single Cell Universal 5' kit


## Raw Data Report



## HiFi Read Length



## Read Segmentation Metrics



- For Revio Kinnex single-cell RNA samples**
- Yield of HiFi reads is ideally >6 Million
  - Yield of segmented reads is ideally >100 Million
  - Mean S-read length will depend on size distribution of input 10x single cell cDNA (ideally >600 bp)
  - % of reads with full arrays is ideally >80%
  - Mean array size is ideally ≥15 segments

Raw Base Yield	1,539 Gb
Mean Polymerase Read Length	90.2 kb
P0	32%
P1	68%
P2	0%

Example sequencing metrics for a human Kinnex single-cell RNA 5' library sample run on a Revio system with Revio SPRQ polymerase kit / 140 pM on-plate loading concentration (OPLC) / 24-hrs movie time.

HiFi Reads	9.4 M
HiFi Base Yield	144.4 Gb
Mean HiFi Read Length	15.3 kb
Median HiFi Read Quality	Q29
HiFi Read Mean # of Passes	9

For human Kinnex single-cell RNA libraries, per-Revio SMRT Cell HiFi read counts were typically >6 Million depending on the final library insert size and P1 loading performance.


Input HiFi Reads	9,408,444
Segmented reads (S-reads)	144,801,149
Mean length of S-reads	971 bp
Percent of reads with full arrays	92.59%
Mean array size (concentration factor)	15.39

For Kinnex single-cell RNA libraries, per-Revio SMRT Cell segmentation read counts were typically >100 Million.


# Example sequencing performance for Kinnex single-cell RNA libraries prepared with human cDNA [ **Revio system + SPRQ chemistry<sup>1</sup>** ] (cont.)

Kinnex single-cell RNA 3' library for PBMC single-cell cDNA sample from 10x Chromium Single Cell Universal 3' kit


## Raw Data Report



## HiFi Read Length



## Read Segmentation Metrics



Raw Base Yield	1,570 Gb
Mean Polymerase Read Length	86.5 kb
P0	27%
P1	72%
P2	0%

Example sequencing metrics for a human Kinnex single-cell RNA 3' library sample run on a Revio system with Revio SPRQ polymerase kit / 140 pM on-plate loading concentration (OPLC) / 24-hrs movie time.

HiFi Reads	9.6 M
HiFi Base Yield	148.1 Gb
Mean HiFi Read Length	15.4 kb
Median HiFi Read Quality	Q29
HiFi Read Mean # of Passes	8

For human Kinnex single-cell RNA libraries, per-Revio SMRT Cell HiFi read counts were typically >6 Million depending on the final library insert size and P1 loading performance.


Input HiFi Reads	9,605,565
Segmented reads (S-reads)	149,894,918
Mean length of S-reads	963 bp
Percent of reads with full arrays	92.24%
Mean array size (concentration factor)	15.61

For Kinnex single-cell RNA libraries, per-Revio SMRT Cell segmentation read counts were typically >100 Million.


# Example sequencing performance for Kinnex single-cell RNA libraries prepared with human cDNA [ Vega system<sup>1</sup> ]

Kinnex single-cell RNA 3' library for PBMC single-cell cDNA sample from 10x Chromium Single Cell Universal 3' kit


## Raw Data Report



## HiFi Read Length



## Read Segmentation Metrics



**For Vega Kinnex single-cell RNA samples**

- Yield of HiFi reads is ideally >4 Million
- Yield of segmented reads is ideally >50 Million
- Mean S-read length will depend on size distribution of input 10x single cell cDNA (ideally >600 bp)
- % of reads with full arrays is ideally >80%
- Mean array size is ideally ≥15 segments

Mean Polymerase Read Length	94.1 kb
Loading Level	71%

Example sequencing metrics for a human Kinnex single-cell RNA 3' library sample run on a Vega system with Vega polymerase kit / 150 pM on-plate loading concentration (OPLC) / 24-hrs movie time.

HiFi Reads	4.4 M
HiFi Base Yield	64.9 Gb
Mean HiFi Read Length	14.8 kb
Median HiFi Read Quality	Q29
HiFi Read Mean # of Passes	9

For human Kinnex single-cell RNA libraries, per-Vega SMRT Cell HiFi read counts were typically >4 Million depending on the final library insert size and sample loading performance.

Input HiFi Reads	4,372,029
Segmented reads (S-reads)	68,570,669
Mean length of S-reads	924 bp
Percent of reads with full arrays	93.78%
Mean array size (concentration factor)	15.68


For Kinnex single-cell RNA libraries, per-Vega SMRT Cell segmentation read counts were typically >50 Million.



# **Kinnex full-length RNA library preparation & sequencing workflow key highlights**

# Kinnex full-length RNA library preparation procedure description

Procedure & checklist – Preparing Kinnex libraries using the Kinnex full-length RNA kit (103-238-700) describes the workflow for constructing Kinnex libraries from total RNA samples using the **Kinnex full-length RNA kit** for sequencing on PacBio long-read systems



PacBio [Documentation](#)  
(103-238-700)



Overview	
Samples	1–24
Workflow time	1.5 days (for up to 24 samples)
Number of SMRT® Cells per Kinnex library Prep	>8 SMRT Cells for Revio using SPRQ™ chemistry >2 SMRT Cells for the Vega or Revio (non-SPRQ) systems >4 SMRT Cells for Sequel II/Ie systems
RNA input	
Quality/size distribution	RIN (RNA integrity number) ≥7.0
Quantity	300 ng per library (minimum concentration 43 ng/μL per library)

- Kinnex full-length RNA library prep protocol uses **Kinnex full-length RNA kit** and **Iso-Seq express 2.0 kit**<sup>1</sup>
  - **Do not use** SMRTbell prep kit 3.0 with this protocol

## Kinnex full-length RNA library barcoding options

Kinnex full-length RNA library preparation procedure supports **up to 48-plex** sample multiplexing through combined use of:

- **12 different barcoded cDNA amplification PCR primers**
- **4 different barcoded Kinnex terminal SMRTbell adapters**



# Comparison of *Iso-Seq express 2.0 kit* vs. *Iso-Seq express oligo kit*



	Iso-Seq express 2.0 kit (103-071-500)	Iso-Seq express oligo kit (101-737-500)
Description	Includes Iso-Seq Express template switching oligo, barcoded cDNA PCR Primers, and other reagents needed for performing 1 <sup>st</sup> -strand cDNA synthesis and PCR amplification of cDNA products	Includes Iso-Seq express template switching oligo and non-barcoded cDNA PCR Primer to be used in conjunction with third-party reagents for performing 1 <sup>st</sup> -strand cDNA synthesis and PCR amplification of cDNA products
# reactions		24
Storage	Refer to product insert for storage instructions	-70°C to -80°C
Compatible SMRTbell library types	Standard (non-concatenated) library Kinnex (concatenated) library	Standard (non-concatenated) library
Sample multiplexing support	Includes 12 barcoded Iso-Seq PCR primers (bc01 – bc12) for up to 12-plex sample multiplexing	Requires additional purchase of barcoded PCR primers for cDNA amplification from a third-party vendor
Kit contents	Iso-Seq RT buffer Iso-Seq RT primer mix Iso-Seq RT enzyme mix Iso-Seq cDNA PCR mix	Iso-Seq template switch oligo Iso-Seq cDNA amplification primer Iso-Seq primers (bc01 – bc12)
SMRT Link support for barcoded cDNA primer Iso-Seq libraries	YES (select 'Iso-Seq v2 Barcoded cDNA Primers' for primer set to use for Iso-Seq analysis application)	YES (select 'Iso-Seq 12 Barcoded cDNA Primers' or 'Custom cDNA Primers' for primer set to use for Iso-Seq analysis application)


# Kinnex full-length RNA experimental design considerations

## Kinnex full-length RNA application use case recommendations for PacBio systems


Example application	Human genetics disease studies	Biopharma for identifying highly expressed targets	Plant & animal whole genome annotation
Experimental goal	Isoform discovery and quantification of moderate-to-rare transcripts	Isoform discovery of high expressed transcripts	Comprehensive transcript annotation in a species
Example study design	Disease vs. normal tissues with multiple replicates	Disease cohort with >20+ samples	Plant or animal with multiple tissue types
Target depth of coverage per sample	10 M reads per sample	5 M reads per sample	5 M reads per sample
Sample multiplexing <sup>1</sup>	<b>Sequel II/IIe system:</b> Up to 2 samples per SMRT Cell 8M (2-plex)	<b>Sequel II/IIe system:</b> Up to 3 samples per SMRT Cell 8M (3-plex)	
	<b>Vega system:</b> Up to 3 samples per Vega SMRT Cell (3-plex)	<b>Vega system:</b> Up to 6 samples per Vega SMRT Cell (6-plex)	
	<b>Revio system + SPRQ:</b> Up to 6 samples per Revio SMRT Cell (6-plex)	<b>Revio system + SPRQ:</b> Up to 12 samples per Revio SMRT Cell (12-plex)	
Expected data throughput (per SMRT Cell)	<b>Sequel II/IIe system:</b> Up to 20 M reads divided by <i>N</i> samples		
	<b>Vega system:</b> Up to 30 M reads divided by <i>N</i> samples		
	<b>Revio system + SPRQ:</b> Up to 60 M reads divided by <i>N</i> samples		
Kinnex library prep protocol	Procedure & checklist – Preparing Kinnex libraries using the Kinnex full-length RNA kit ( <a href="#">103-238-700</a> )		
Total RNA input into Kinnex library prep workflow	300 ng total RNA (RIN ≥7) for 1 <sup>st</sup> -strand cDNA synthesis		
SMRT Link data analysis workflows	Read Segmentation and Iso-Seq analysis application with option to “pool reads and cluster together” to get a master isoform classification file with per-sample full-length read counts		
Community data analysis tools	Annotation & quantification: PIGEON, SQANTI3 / Differential analysis: TappAS / Fusion calling: pbfusion / Visualization: SWAN		

# Kinnex full-length RNA method overview

## Full-length cDNA synthesis & amplification



## Kinnex library prep, sequencing & analysis



<sup>1</sup> Twelve barcoded Iso-Seq primers (Iso-Seq primer bc01–12) are available for cDNA amplification step.

<sup>2</sup> Kinnex adapter barcode sequences can be downloaded from [SMRT Link Data Management module](#).

# Kinnex full-length RNA library preparation procedural notes

## I. Library insert generation



### cDNA amplification

- **Set up on ice** and add PCR reaction to thermal cycler **after the lid has preheated to 105°C** to avoid digestion of primers by polymerase exonuclease activity
- Barcoded primers are used during cDNA amplification
  - 12 barcoded forward primers (bc01-bc12) available for use in combination with Iso-Seq cDNA amplification primer

### SMRTbell bead cleanup

- Previous non-Kinnex full-length RNA (bulk monomer) Iso-Seq protocol (*Procedure & checklist – Preparing Iso-Seq libraries using SMRTbell prep kit 3.0* [102-396-000]) offered 3 options for post-cDNA amplification SMRTbell bead cleanup: 0.82X, 0.86X, and 0.95X → For simplification, Kinnex full-length RNA protocol now only specifies to use **0.9X**

## II. Kinnex PCR



### Kinnex PCR

- Can transfer entire volume of primers to PCR tubes for ease of use with multi-channel pipettes (8 primer mix tubes)
- **Set up on ice** and add PCR reaction to thermal cycler **after lid has preheated to 105°C** to avoid digestion of primers by polymerase exonuclease activity

### Pooling of 8 Kinnex PCR products + SMRTbell bead cleanup

- **Pool exactly 23 µL** from each Kinnex PCR reaction for a total combined volume of **184 µL**
- Add exactly **193 µL** of SMRTbell cleanup beads (**1.05X**)
- Kinnex PCR mix significantly increases stringency of SMRTbell clean up beads, so **accurate pipetting is critical**

## III. Kinnex array formation



### Kinnex digestion

- Recommended input amount to proceed with Kinnex array formation is **5 µg** of Kinnex PCR amplicons
  - Proceeding with <3 µg is not recommended since lower input amounts may lead to insufficient final library yields to enable optimal sequencing results


### Final SMRTbell bead cleanup

- Perform **1.0X** SMRTbell bead cleanup on final library

# Example Kinnex full-length RNA library preparation QC results


Kinnex full-length RNA library prepared with human UHRR total RNA

## Amplified full-length cDNA QC



Example Bioanalyzer DNA sizing QC analysis results for amplified full-length cDNA generated from a universal human RNA reference (UHRR) total RNA sample.

## Final Kinnex full-length RNA library QC



Example Femto Pulse DNA sizing QC analysis results for final Kinnex full-length RNA library.

Final Kinnex library yield is typically sufficient to load the following number of SMRT Cells:

- >8 SMRT Cells for Revio using SPRQ chemistry
- >2 SMRT Cells for the Vega or Revio (non-SPRQ) systems
- >4 SMRT Cells for Sequel II/IIe systems

Total RNA input for cDNA synthesis	300 ng
------------------------------------	--------

cDNA input for Kinnex array formation	5900 ng
---------------------------------------	---------


Post-nuclease treatment & final library cleanup yield (%) <sup>1</sup>	1460 ng (24.7%)
--	-----------------

<sup>1</sup> Post-nuclease treatment & final cleanup yields typically ranged from ~10% to ~25% when using UHRR total RNA samples for Kinnex full-length RNA library construction.


# Example sequencing performance for Kinnex full-length RNA libraries prepared with human cDNA [ Revio system + SPRQ chemistry<sup>1</sup> ]

Kinnex full-length RNA library for universal human reference RNA (UHRR) sample


## Raw Data Report



## HiFi Read Length



## Read Segmentation Metrics



**For Revio Kinnex full-length RNA samples**

- Yield of HiFi reads is ideally >6 Million
- Yield of segmented reads is ideally >50 Million
- Mean S-read length will depend on size distribution of cDNA sample (ideally >1,500 bp)
- % of reads with full arrays is ideally >80%
- Mean array size is ideally >7.0 segments

Raw Base Yield	1,442 Gb
Mean Polymerase Read Length	71.9 kb
P0	20%
P1	80%
P2	1%

Example sequencing metrics for a Universal Human Reference RNA (UHRR) Kinnex full-length RNA library sample run on a Revio system with Revio SPRQ polymerase kit / 160 pM on-plate loading concentration (OPLC) / 24-hrs movie time.

HiFi Reads	10.4 M
HiFi Base Yield	146.6 Gb
Mean HiFi Read Length	14.0 kb
Median HiFi Read Quality	Q32
HiFi Read Mean # of Passes	7

For UHRR Kinnex full-length RNA libraries, per-Revio SMRT Cell HiFi read counts were typically >6 Million depending on the final library insert size and P1 loading performance.


Input HiFi Reads	10,404,037
Segmented reads (S-reads)	81,012,128
Mean length of S-reads	1,758 bp
Percent of reads with full arrays	94.55%
Mean array size (concentration factor)	7.79

For UHRR Kinnex libraries, per-Revio SMRT Cell segmentation read counts were typically >50 Million.


# Example sequencing performance for Kinnex full-length RNA libraries prepared with human cDNA [ Vega system<sup>1</sup> ]

Kinnex full-length RNA library for universal human reference RNA (UHRR) sample


## Raw Data Report



## HiFi Read Length



## Read Segmentation Metrics



**For Vega Kinnex full-length RNA samples**

- Yield of HiFi reads is ideally >2.5 Million
- Yield of segmented reads is ideally >20 Million
- Mean S-read length will depend on size distribution of cDNA sample (ideally >1,500 bp)
- % of reads with full arrays is ideally >80%
- Mean array size is ideally >7.0 segments

Mean Polymerase Read Length	95.3 kb
Loading level	61%

Example sequencing metrics for a Universal Human Reference RNA (UHRR) Kinnex full-length RNA library sample run on a Vega system with Vega polymerase kit / 130 pM on-plate loading concentration (OPLC) / 24-hrs movie time.

HiFi Reads	4.6 M
HiFi Base Yield	65.0 Gb
Mean HiFi Read Length	14.0 kb
Median HiFi Read Quality	Q35
HiFi Read Mean # of Passes	11

For UHRR Kinnex full-length RNA libraries, per-Vega SMRT Cell HiFi read counts were typically >2.5 Million depending on the final library insert size and sample loading performance.

Input HiFi Reads	4,642,616
Segmented reads (S-reads)	35,917,435
Mean length of S-reads	1,759 bp
Percent of reads with full arrays	94.14%
Mean array size (concentration factor)	7.74

For UHRR Kinnex libraries, per-Vega SMRT Cell segmentation read counts were typically >20 Million.



# **Kinnex 16S rRNA library preparation & sequencing workflow key highlights**

# Kinnex 16S rRNA library preparation procedure description

Procedure & checklist – Preparing Kinnex libraries from 16S rRNA amplicons (103-238-800) describes the workflow for constructing Kinnex libraries from full-length 16S amplicons using the **Kinnex 16S rRNA kit\*** for sequencing on PacBio long read systems



**Kinnex 16S rRNA kit**  
(103-072-100)

**Concatenation factor**  
12-fold

**Kinnex PCR 12-fold kit**  
103-071-700 (12 rxn)

**Kinnex concatenation kit**  
103-071-800 (12 rxn)

**SMRTbell cleanup beads**  
102-158-300

**Elution buffer**  
101-633-500

PacBio [Documentation](#)  
(103-238-00)



Overview	
Samples	1–384
Workflow time	1.5 days for up to 96 samples
	>8 SMRT Cells for Revio using SPRQ™ chemistry
Number of SMRT Cells per library Prep	>2 SMRT Cells for the Vega or Revio (non-SPRQ) systems
	>10 SMRT Cells for Sequel II/IIe
Metagenomic DNA mass	Total genomic DNA of 1–2 ng per metagenomic sample

- Kinnex full-length RNA library prep protocol uses **Kinnex 16S rRNA kit**
- **Do not use** SMRTbell prep kit 3.0 with this protocol

## Kinnex full-length 16s rRNA library barcoding options

Kinnex 16S rRNA library preparation procedure supports up to **1,536-plex** sample multiplexing through combined use of:

- 12 different 16S barcoded Forward PCR primers<sup>1</sup>
- 32 different 16S barcoded Reverse PCR primers<sup>1</sup>
- 4 different barcoded Kinnex terminal SMRTbell adapters



# Kinnex 16S rRNA experimental design considerations

## Kinnex 16S rRNA application use case recommendations for PacBio systems


	Sequel II and IIe systems	Vega system	Revio system + SPRQ chemistry			
Experimental goal	Determine the microbial diversity (phylogeny and taxonomy) of bacteria in a metagenomic sample					
Sample multiplexing <sup>1</sup>	Up to 384 samples per SMRT Cell 8M (384-plex)	Up to 1,152 samples per Vega SMRT Cell (1,152-plex)	Up to 1,536 samples per Revio SMRT Cell (1,536-plex)			
Expected coverage per sample <sup>2</sup>	96-plex	208 K	96-plex	416 K	96-plex	729 K
	192-plex	104 K	192-plex	208 K	192-plex	364 K
	384-plex	52 K	384-plex	104 K	384-plex	182 K
	768-plex	26 K	768-plex	52 K	768-plex	91 K
	1,152-plex	17 K	1,152-plex	34 K	1,152-plex	60 K
	1,536-plex	13 K	1,536-plex	26 K	1,536-plex	45 K
Kinnex library prep protocol	Procedure & checklist – Preparing Kinnex libraries from 16S rRNA amplicons (103-238-800)					
Metagenomic DNA input amount input into 16S gene amplification	1-2 ng of input gDNA per metagenomic sample					
16S amplicon DNA input into Kinnex library prep workflow	35 ng of purified pooled 16S amplicon DNA					
SMRT Link data analysis workflows	Read Segmentation					
Community data analysis tools	pb-16S-nf					

<sup>1</sup> Kinnex concatenation kit (103-071-800) can support up to 1,536-plex sample multiplexing through the combined use of 12 different 16S barcoded Forward PCR primers + 32 different 16S barcoded Reverse PCR primers and 4 different barcoded Kinnex terminal SMRTbell adapters during Kinnex 16s rRNA library construction.

<sup>2</sup> With proper full array formation and adequate sequencing, one SMRT Cell on the Sequel II, IIe, Vega and Revio systems are expected to achieve 20–25 million, 30-40 million and 70–80 million 16S sequences, respectively. For most 16S analysis applications, typically aim for ~30-50 K reads/sample.

# Kinnex 16S rRNA method overview

## Full-length 16S gene amplification




Bacterial gDNA fragments isolated from metagenomic samples

Recommended gDNA input for full-length 16S PCR amplification = 1-2 ng

PCR amplification of full-length 16S genes

Can barcode up to 384 samples during 16S gene amplification step using PacBio barcoded Fwd and Rev primers<sup>1</sup>



Sample plate layout for 96-plex PCR design using 12 different 16S Barcoded Forward Primers and 8 different 16S Barcoded Reverse Primers.

To multiplex 384 samples, use 12 barcoded Fwd primers + 32 barcoded Rev primers and set up four 96-well PCR reaction plates.

Forward (F) PCR primer (Kinnex16S\_Fwd\_XX<sup>1</sup>)

Reverse (R) PCR primer (Kinnex16S\_Rev\_XX<sup>1</sup>)



Dual-barcoded full-length 16S amplicon product

Contains Kinnex PCR handle 1

Equal-volume pooling & cleanup



Pooled barcoded full-length 16S amplicons for input into Kinnex library construction

Contains Kinnex PCR handle 2

## Kinnex library prep, sequencing & analysis




Amplified 16S products containing a different PacBio barcode at each end (Fwd BC ≠ Rev BC)

Kinnex PCR (12 parallel PCR Rx to append Kinnex segmentation adapters)



Pooling & Kinnex array formation

Barcoded Kinnex terminal adapter (bcM0001-bcM0004)<sup>2</sup>




Can also barcode samples during Kinnex array formation using up to 4 different barcoded Kinnex terminal SMRTbell adapters<sup>2</sup>

DNA damage repair & nuclease Tx




Anneal Kinnex sequencing primer / Bind polymerase / Cleanup (ABC)



Purified polymerase-bound 12-segment Kinnex library

Long-read sequencing & data analysis

Sequel II/IIe, Vega or Revo system



Third-party metagenomic secondary analysis tools

<sup>1</sup> 12 different 16S barcoded Forward PCR primers + 32 different 16S barcoded Reverse PCR primers are available for 16S gene amplification step to multiplex up to 384 samples.

<sup>2</sup> Kinnex adapter barcode sequences can be downloaded from [SMRT Link](#) Data Management module.

# Kinnex 16S rRNA library preparation procedural notes

## I. Library insert generation



### PCR amplification of full-length 16S genes

- Customer supplies Kapa PCR mix (HiFi HotStart ReadyMix) and oligos
- Up to **384-plex** can be done at this point using **combinatorial indexing**
- **Set up on ice** and add PCR reaction to thermal cycler **after the lid has preheated to 105°C** to avoid digestion of primers by polymerase exonuclease activity

## II. Kinnex PCR



### Kinnex PCR

- Can transfer entire volume of primers to PCR tubes for ease of use with multi-channel pipettes (12 primer mix tubes)
- **Set up on ice** and add PCR reaction to thermal cycler **after lid has preheated to 105°C** to avoid digestion of primers by polymerase exonuclease activity

### SMRTbell bead cleanup

- **Pool exactly 23 µL** from each Kinnex PCR reaction for a total combined volume of **276 µL**
- **Add exactly 304 µL** of SMRTbell clean up beads (**1.1X**)
- Kinnex PCR mix significantly increases stringency of SMRTbell clean up beads, so **accurate pipetting is critical**

## III. Kinnex array formation



### Kinnex digestion

- Recommended input amount to proceed with Kinnex array formation is **5 µg** of Kinnex PCR amplicons
  - Proceeding with <3 µg is not recommended since lower input amounts may lead to insufficient final library yields to enable optimal sequencing results


### Final SMRTbell bead cleanup

- Perform **1.0X** SMRTbell bead cleanup on final library

# Example Kinnex 16S rRNA library preparation QC results


Kinnex full-length 16S RNA library prepared from mock microbial community genomic DNA

## 16S amplicon DNA QC



Example Bioanalyzer DNA sizing QC analysis results for pooled 16S amplicon DNA samples generated from mock microbial community genomic DNA (ATCC MSA-1003 20 Strain Staggered Mix).

## Final Kinnex 16S rRNA library QC



Example Femto Pulse DNA sizing QC analysis results for final Kinnex 16S rRNA library.

### Example sample multiplexing design for 1,536-plex Kinnex 16S rRNA library

384-plex 16S PCR x 4-plex Kinnex adapter barcoding → Total sample multiplex level = **1,536-plex**

 <b>ZymoBIOMICS Fecal Reference</b> 96 barcoded samples	 <b>ZymoBIOMICS Gut Microbiome Standard</b> 96 barcoded samples	 <b>ATCC MSA-1002 20 Strain Even Mix</b> 96 barcoded samples	 <b>ATCC MSA-1003 20 Strain Staggered Mix</b> 96 barcoded samples	<b>X</b> 4 Kinnex barcoded adapters
--	--	---	---	-------------------------------------


gDNA input for 16S PCR	1.1 ng
16S amplicon DNA input Kinnex PCR products for Kinnex array formation	6000 ng
Post-nuclease treatment & final library cleanup yield (%) <sup>1</sup>	1080 ng (18%)

<sup>1</sup> Post-nuclease treatment & final cleanup yields typically ranged from ~10% to ~20% when using mock microbial community genomic DNA for Kinnex full-length 16S rRNA library construction.


# Example sequencing performance for Kinnex 16S rRNA libraries prepared from mock microbial community genomic DNA [ Revio system + SPRQ<sup>1</sup> ]

Kinnex full-length 16S library for mock microbial community gDNA sample mixture (384-plex)


## Raw Data Report



## HiFi Read Length



## Read Segmentation Metrics



**For Revio system Kinnex 16S rRNA samples**

- Yield of HiFi reads is ideally >6 Million
- Yield of segmented reads is ideally >70 Million
- Mean S-read length should be consistent with expected 16S amplicon length (~1.5 kb)
- % of reads with full arrays is ideally >80%
- Mean array size is ideally >11.0 segments

Raw Base Yield	1,400 Gb
Mean Polymerase Read Length	80.5
P0	30%
P1	69%
P2	3%

Example sequencing metrics for a Kinnex 16S rRNA library sample run on a Revio system with Revio SPRQ polymerase kit / 130 pM on-plate loading concentration (OPLC) / 24-hrs movie time.

HiFi Reads	7.7 M
HiFi Base Yield	144.5 Gb
Mean HiFi Read Length	18.7 kb
Median HiFi Read Quality	Q34
HiFi Read Mean # of Passes	7

For Kinnex 16S rRNA libraries, per-Revio SMRT Cell HiFi read counts were typically >6 Million depending on the final library insert size and P1 loading performance.


Input HiFi Reads	7,713,616
Segmented reads (S-reads)	91,255,042
Mean length of S-reads	1,558 bp
Percent of reads with full arrays	96.16%
Mean array size (concentration factor)	11.83

For Kinnex 16S rRNA libraries, per-Revio SMRT Cell segmentation read counts were typically >70 Million.


# Example sequencing performance for Kinnex 16S rRNA libraries prepared from mock microbial community genomic DNA [ Vega system<sup>1</sup> ]

Kinnex full-length 16S library for mock microbial community gDNA sample mixture (384-plex)


## Raw Data Report



## HiFi Read Length



## Read Segmentation Metrics



Mean Polymerase Read Length	64.1
Loading level	67%

Example sequencing metrics for a Kinnex 16S rRNA library sample run on a Vega system with Vega polymerase kit / 130 pM on-plate loading concentration (OPLC) / 24-hrs movie time.

HiFi Reads	3.5 M
HiFi Base Yield	65.0 Gb
Mean HiFi Read Length	18.8 kb
Median HiFi Read Quality	Q31
HiFi Read Mean # of Passes	6

For Kinnex 16S rRNA libraries, per-Vega SMRT Cell HiFi read counts were typically >2.5 Million depending on the final library insert size and sample loading performance.

Input HiFi Reads	3,464,764
Segmented reads (S-reads)	41,088,694
Mean length of S-reads	1,558 bp
Percent of reads with full arrays	95.86%
Mean array size (concentration factor)	11.86

For Kinnex 16S rRNA libraries, per-Vega SMRT Cell segmentation read counts were typically >30 Million.

# Improving sequencing performance of “difficult” 16S samples

Performing AMPure PB bead size-selection on Kinnex full-length 16S rRNA libraries can help improve sample loading of challenging metagenomic samples

Sample Name	<i>P1</i> % [ Revio system ]	Gb Yield	Mean Length	Mean QV
16S_collaborator_SOP	26	46	18,813 bp	Q29
<b>16S_collaborator_3.1X AMPure</b>	<b>80</b>	<b>87</b>	<b>18,851 bp</b>	<b>Q28</b>

Some bacterial 16S samples may have carry-over contaminants present leading to low sample (*P1*) loading on PacBio long-read systems

→ Using AMPure size-selection (3.1X 35% AMPure PB beads) can help mitigate this issue



# **Technical documentation & applications support resources**

# Technical resources for Kinnex library preparation, sequencing & data analysis

## Single-cell cDNA sample preparation literature & other resources

- 10x Genomics Chromium Next GEM Single Cell 3' v3.1 (Single Index) How-to Video [ [Link](#) ]
- 10x Genomics Chromium Single Cell 3' Reagent Kits User Guide – v3.1 ([CG000204](#))
- 10x Genomics Chromium Single Cell 5' Reagent Kits User Guide – v2 Chemistry Dual Index ([CG000331](#))

## Kinnex library preparation literature & other resources

- Application note – Kinnex 16S rRNA kit for full-length 16S sequencing ([102-326-601](#))
- Application note – Kinnex full-length RNA kit for isoform sequencing ([102-326-591](#))
- Application note – Kinnex single-cell RNA kit for single-cell isoform sequencing ([102-326-549](#))
- Brochure – Scalable, cost-effective RNA sequencing with PacBio Kinnex kits ([102-326-597](#))
- Procedure & checklist – Preparing Kinnex libraries using Kinnex single-cell RNA kit ([103-254-300](#))
- Procedure & checklist – Preparing Kinnex libraries using Kinnex full-length RNA kit ([103-238-700](#))
- Procedure & checklist – Preparing Kinnex libraries from 16s rRNA amplicons ([103-238-800](#))
- Technical overview – Kinnex kits for single-cell RNA, full-length RNA and 16S rRNA sequencing ([103-343-700](#))
- Technical overview – Kinnex library preparation for full-length 16S rRNA gene sequencing ([103-344-800](#))
- Technical overview – Kinnex library preparation using Kinnex full-length RNA kit ([103-344-700](#))
- Technical overview – Kinnex library preparation using Kinnex single-cell RNA kit ([103-344-600](#))
- Video tutorial – PacBio Kinnex single-cell RNA TSO artifact removal demo for Kinnex single-cell RNA kit [ [Link](#) ]
- Video tutorial – SMRT Link Sample Setup and Run Design setup procedure for Kinnex kits [ [Link](#) ]
- Whitepaper – Bulk and single-cell isoform sequencing for human disease research ([102-326-576](#))

# Technical resources for Kinnex library preparation, sequencing & data analysis (cont.)

## Data analysis resources

- Application note – Bioinformatics tools for full length isoform sequencing ([102-326-593](#))
- SMRT Link Kinnex full-length RNA troubleshooting guide ([103-552-100](#))
- SMRT Link Kinnex single-cell troubleshooting guide ([103-516-100](#))
- SMRT Link MAS-Seq troubleshooting guide ([102-994-400](#))
- SMRT Link software installation guide [ [Link](#) ]
- SMRT Link user guide [ [Link](#) ]
- SMRT Tools reference guide [ [Link](#) ]
- Video tutorial – Analyzing Kinnex 16S rRNA data in SMRT Link [ [Link](#) ]
- Video tutorial – Read Segmentation and Iso-Seq workflow in SMRT Link [ [Link](#) ]

## Publications and posters

- Wissel, D. et al. (2025) A systematic benchmark of high-accuracy PacBio long-read RNA sequencing for transcript-level quantification. BioRxiv [ [Link](#) ]
- Deng, E. et al. (2025) Systematic evaluation of single-cell RNA-seq analyses performance based on long-read sequencing platforms. J Advanced Res [ [Link](#) ]
- Al'Khafaji, A.M. et al. (2024) High-throughput RNA isoform sequencing using programmable cDNA concatenation. Nature biotechnology. [ [Link](#) ]
- Schertzer, M.D. et al. (2023) Cas13d-mediated isoform-specific RNA knockdown with a unified computational and experimental toolbox. Nature comm [ [Link](#) ]

## Webinars

- PacBio video (2024) – Kinnex explained – how concatenating smaller amplicons increases throughput for PacBio HiFi sequencing [ [Link](#) ]
- PacBio PRISM webinar (2024) – Let's stick together – exploring PacBio Kinnex kits [ [Link](#) ]
- PacBio video (2024) – Increasing the throughput of full-length 16S and RNA-Seq with Kinnex kits [ [Link](#) ]
- PacBio Iso-Seq social club webinar (2022) – TappAS for isoform differential expression analysis [ [Link](#) ]
- PacBio Iso-Seq social club webinar (2022) – Single-cell Iso-Seq applications in cancer and neurological disorders [ [Link](#) ]

# Technical resources for Kinnex library preparation, sequencing & data analysis (cont.)

## Example PacBio data sets

Application	Dataset	Data type	PacBio system
Kinnex single-cell RNA sequencing	Homo sapiens - PBMC 10x Chromium Single Cell 5' and 3' libraries [ <a href="#">Link</a> ]	HiFi long read	Sequel II, Revio and Vega systems
	Homo sapiens - HG002 (10x 5') [ <a href="#">Link</a> ]	HiFi long read	Revio system
Kinnex full-length RNA sequencing	Homo sapiens – Universal human reference RNA (UHRR) [ <a href="#">Link</a> ]	HiFi long read	Vega system
	Homo sapiens – Universal human reference RNA (UHRR) [ <a href="#">Link</a> ]	HiFi long read	Revio system – SPRQ chemistry
	Homo sapiens – Universal human reference RNA (UHRR) [ <a href="#">Link</a> ]	HiFi long read	Sequel II & Revio systems
	Homo sapiens – HG002 [ <a href="#">Link</a> ]	HiFi long read	Revio system
	Homo sapiens – Heart [ <a href="#">Link</a> ]	HiFi long read	Revio system
	Homo sapiens – Cerebellum [ <a href="#">Link</a> ]	HiFi long read	Revio system
	Homo sapiens – Brain [ <a href="#">Link</a> ]	HiFi long read	Revio system

# Technical resources for Kinnex library preparation, sequencing & data analysis (cont.)

## Example PacBio data sets

Application	Dataset	Data type	PacBio system
Kinnex 16S rRNA sequencing	ZymoBIOMICS Fecal Reference with TruMatrix Technology (human) [ <a href="#">Link</a> ]	HiFi long read	Sequel II & Revio systems
	Mixture: ZymoBIOMICS Gut Microbiome Standard, ZymoBIOMICS Fecal Reference with TruMatrix™ Technology, ATCC 20 Strain Even Mix Genomic Material, ATCC 20 Strain Staggered Mix Genomic Material [ <a href="#">Link</a> ]	HiFi long read	Sequel II & Revio systems



# **APPENDIX 1: Concatenating cDNA or gDNA amplicons using Kinnex kits to increase throughput**

# Concatenating cDNA or gDNA amplicons using Kinnex kits to increase throughput

Kinnex concatenation is a general method that can increase sequencing throughput for smaller amplicons<sup>1</sup>

## Benefits of Kinnex concatenation

- Increased throughput on PacBio long-read sequencers
- Retained HiFi accuracy despite throughput increase
- No change to secondary analysis – once reads are deconcatenated into S-reads, the S-reads represent the original, pre-concatenated amplicon and can be analyzed with established pipelines.

## When is Kinnex concatenation appropriate?

- The balance between the amplicon size and the concatenation factor, as well as additional Kinnex library generation cost, needs to be taken into consideration.
- HiFi sequencing produces optimal yield for inserts between 15–20 kb; therefore, the throughput advantage plateaus for larger amplicon sizes exceeding 3 kb.
- You can consider concatenating amplicons using Kinnex kits if:
  - The amplicons have an **average size between 200 bp – 3 kb**
  - The amplicons have **molecular ends that are either directly compatible or can be re-amplified to establish Kinnex compatibility**

Recommended Kinnex kit based on average amplicon sizes.

Average amplicon size	Example	Recommended Kinnex kit	Expected Kinnex library size
600 – 1000 bp	10x single-cell cDNA	Kinnex single-cell RNA kit (16-fold)	10 – 16 kb
1 – 2 kb	Full-length 16S	Kinnex 16S rRNA kit (12-fold)	~19 kb
2 – 3 kb	Bulk cDNA	Kinnex full-length RNA kit (8-fold)	15 – 20 kb
>3 kb	Not recommended for Kinnex concatenation		

Technical note  
**CONCATENATING AMPLICONS USING PACBIO KINNEX KITS TO INCREASE THROUGHPUT**

**Overview**  
This technical note describes the principles and advantages of concatenating amplicons using PacBio® Kinnex™ kits.

**Note:** This technical note is intended as a guide for best practices and to report typical example results (see Appendix). PacBio does not guarantee the success of concatenating amplicons that are not officially supported by the Kinnex kits, which currently include:

- Full-length cDNA generated using the Iso-Seq® express 2.0 kit with the Kinnex full-length RNA kit
- Full-length 16S sequences generated according to official protocol with the Kinnex 16S rRNA kit
- Single-cell cDNA generated using compatible single-cell platforms with the Kinnex single-cell RNA kit


How Kinnex works

**What is Kinnex?**  
The Kinnex kits are based on the method developed by Al'Khafaji et al. (2023). Originally termed multiplexed array isomorph sequencing (MAS-ISO-seq, or MAS-Seq), this method was developed to concatenate cDNAs into longer fragments suitable for long-read sequencing. It takes advantage of the disparity between optimal DNA fragment lengths for HiFi sequencing, 15–20 kb, and smaller sizes of transcript cDNA (1–10 kb with typical estimated average length of 2 kb for a human transcriptome). The original MAS-ISO-Seq method concatenated single-cell cDNA from the 10x platform, which has an average read length of 600–800 bp, and used a 15-fold concatenation array to increase throughput. PacBio commercialized the MAS-Seq method and adjusted the concatenation factors according to different average amplicon sizes (Table 1).

**Technical note – Concatenating cDNA or gDNA amplicons using Kinnex kits to increase throughput**  
[102-326-636](https://www.pacb.com/technical-note/concatenating-amplicons-using-kinnex-kits-to-increase-throughput)

# Kinnex concatenation workflow overview

Follow Kinnex library prep protocol documentation for specific details on concatenating cDNA or gDNA amplicons using Kinnex kits



**Kinnex concatenation workflow.** Amplicons must be generated or amplified to have Kinnex-compatible molecular ends before continuing to the Kinnex PCR and array formation step. Kinnex libraries should be sequenced with the appropriate sequencing chemistry and run configurations. Once de-concatenated using Read Segmentation in SMRT Link, the individual amplicons can be analyzed using amplicon-specific workflows.

## Kinnex concatenation procedural notes<sup>1</sup>

- To establish Kinnex compatibility, **amplicons are required to have Kinnex-compatible molecular ends** (see next section)
  - Once these are generated, choose the appropriate Kinnex kit based on the recommended concatenation factor listed in Table on previous slide and proceed with Kinnex PCR
- The **Kinnex PCR** steps consist of parallel PCR reactions per sample [i.e., 8, 12, or 16 reactions based on the Kinnex kit chosen] using premixed Kinnex primer pairs
  - The resulting PCRs generate amplified DNA products containing programmable sequences at both ends
- In the **Kinnex array formation** step, library inserts containing programmable ends are assembled to generate a linear array.
  - Further, the addition of barcoded Kinnex terminal adapters result in the formation of complete, full-length array SMRTbell templates along with partial arrays.
  - Subsequent nuclease treatment removes partial arrays to retain only full array SMRTbell templates for achieving optimal sequencing yield.
- **Note:** Kinnex terminal adapters are different from standard SMRTbell adapters and hence require the **Kinnex sequencing primer** (103-179-000) during the “Annealing, Binding, and Cleanup (ABC)” step for optimal sequencing results.
- Once HiFi reads are generated, **Read Segmentation** will produce the segmented reads (S-reads) that represent the original unconcatenated amplicons, which can be used for further analysis.

# Establishing Kinnex-compatible molecular ends

Amplicons must be generated or amplified to have Kinnex-compatible molecular ends before proceeding with Kinnex library prep procedure

- To be compatible with the Kinnex workflow, amplicons must be generated with sequence-defined ends as depicted in the underlined portion in the figure below
- Optional barcodes (such as sample indices, UMIs and single cell barcodes) should be placed internally between the Kinnex handles and the amplicon-specific primers.
- Kinnex handles may be present already in certain amplicons, such as 10x Single Cell Gene Expression libraries or the Kinnex 16S amplicons, or can be added by PCR amplification, such as for Parse Evercode single-cell libraries

Kinnex FWD primer CTACACGACGCTCTCCGATCT - [optional barcodes] - [amplicon specific FWD primers]

Kinnex REV primer AAGCAGTGGTATCAACGCAGAG - [optional barcodes] - [amplicon specific REV primers]

## Example 1. Kinnex 16S forward and reverse primer sequences

Kinnex 16S FWD01 CTACACGACGCTCTCCGATCT - GATCGAGTCA - AGRGTTYGATYMTGGCTCAG

Kinnex 16S REV13 AAGCAGTGGTATCAACGCAGAG - TCATCGACGT - RGYTACCTTGTACGACTT

## Example 2. Iso-Seq express 2.0 forward and reverse primer sequences

IsoSeqX bc01 FWD CTACACGACGCTCTCCGATCT - ACTACAC - GCAATGAAGTCGCAGGGTTGGG

IsoSeqX REV AAGCAGTGGTATCAACGCAGAGTAC

**Schematic for Kinnex-compatible primers.** Kinnex handles (5' to 3') are shown in black underline and must be present at the ends of the amplicons to be compatible with Kinnex concatenation. Optional barcodes can be included internally. Amplicon-specific primers (and optional internal barcode sequences) must be designed to avoid strong secondary structures in the context of Kinnex handles.

## Resources for Kinnex library prep<sup>1</sup>

### For using Kinnex full-length RNA kit (PN: 103-072-000) for 8-fold concatenation:

- Procedure & checklist – Preparing Kinnex libraries using the Kinnex full-length RNA kit ([103-238-700](#))
- Technical overview – Kinnex library preparation using Kinnex full-length RNA kit ([103-344-700](#))

### For using Kinnex 16S rRNA kit (PN: 103-072-100) for 12-fold concatenation:


- Procedure & checklist – Preparing Kinnex libraries using 16S rRNA amplicons ([103-238-800](#))
- Technical overview – Kinnex library preparation using Kinnex 16S rRNA kit ([103-344-800](#))

### For using Kinnex single-cell RNA kit (PN: 103-072-200) for 16-fold concatenation:

- Procedure & checklist – Preparing Kinnex libraries using Kinnex single-cell RNA kit ([103-254-300](#))
- Technical overview – Kinnex library preparation using Kinnex single-cell RNA kit ([103-344-600](#))



# **APPENDIX 2: PacBio compatible Kinnex library preparation workflows**






# Preparing Kinnex single-cell libraries with Parse Evercode Whole Transcriptome kits

# Preparing Kinnex single-cell libraries with Parse Evercode Whole Transcriptome kits: Getting started

Application-specific educational literature

Application-specific protocol documentation

Application-specific technical overviews



PacBio literature website [Link]

Application-specific brochures, informational guides and other product literature containing best practices recommendations for library preparation and data analysis workflows.

Technical note - Preparing Kinnex single-cell libraries with Parse Evercode Whole Transcriptome kits (102-326-610) [PacBio]

Technical documentation describing a modified workflow for preparing PacBio Kinnex single-cell libraries with single-cell or single-nucleus cDNA generated using Parse Biosciences Evercode Whole Transcriptome products.

Procedure & checklist - Procedure & checklist - Preparing Kinnex libraries using Kinnex single-cell RNA kit (103-254-300) [PacBio]

Technical documentation containing PacBio SMRTbell library construction details.

Technical overview: Kinnex library preparation using Kinnex single-cell RNA kit (103-344-600)

Technical overview presentations describe sample preparation details for constructing HiFi libraries for specific applications. Example sequencing performance data for a given application are also summarized.

# Technical reference for preparing Kinnex single-cell libraries with Parse Evercode Whole Transcriptome kits

*Technical note – Preparing Kinnex single-cell libraries with Parse Evercode Whole Transcriptome kits ([102-326-610](#)) describes a modified workflow for preparing PacBio Kinnex single-cell libraries with single-cell or single-nucleus cDNA generated using Parse Biosciences Evercode Whole Transcriptome products*

Technical note  
**PREPARING KINNEX™ SINGLE-CELL LIBRARIES WITH PARSE BIOSCIENCES EVERCODE™ WT KITS**

**Overview**  
This technical note describes the experimental conditions of a modified workflow for preparing PacBio® Kinnex single-cell libraries with single-cell or single-nucleus cDNA generated using Parse Biosciences Evercode Whole Transcriptome products, the WT and WT Mini kits. This workflow has not been tested using the Evercode WT Mega kit.


**Required materials and equipment**

- Refer to all materials described in the Kinnex single-cell protocol.<sup>1</sup>
- Single-cell or single-nucleus cDNA from the Parse Biosciences Evercode WT or Evercode WT Mini kits.
  - Note that the cDNA needs to be generated as described in the Parse protocols<sup>2</sup> for either the Evercode WT or Evercode WT Mini up through the end of those protocols' section 2. Then use either the 15 ng (from the WT protocol) or 60–75 ng (from the WT Mini protocol) of the generated cDNA as input into the workflow at right (see Figure 1).
- Parse Biosciences-specific amplification primers (Appendix 1).

**Library preparation**  
Library preparation will follow the PacBio Kinnex single-cell protocol<sup>1</sup> with modifications outlined below. Workflow modifications begin with single-cell or single-nucleus cDNA produced by the Parse Biosciences Evercode WT and WT Mini kits.

- Generate amplified cDNA per Parse Biosciences protocols. Complete the Evercode WT or Evercode WT Mini protocol up through section 2.
- Carry out QC as described in the Parse Biosciences user manual (see the figure at the end of section 2. Expected post-amplification sublibrary cDNA size distribution). While cDNA size ranges can vary, there should not be an excess of shorter (<500 bp) material.





**PacBio Technical note**  
*Preparing Kinnex single-cell libraries with Parse Evercode Whole Transcriptome kits ([102-326-610](#))*  
[ Select to use either single-cell or single-nucleus cDNA from Parse Biosciences Evercode WT or Evercode WT Mini kits ]


**PacBio Procedure & checklist**  
*Preparing Kinnex libraries using Kinnex single-cell RNA kit ([103-254-300](#))*

## Overview

- Parse whole transcriptome kits use split-pool combinatorial barcoding to enable scale-up of single cell projects to millions of cells or nuclei<sup>1</sup>
- This Kinnex library prep workflow has been tested with single-cell or single-nucleus cDNA generated using Parse Biosciences Evercode WT and WT Mini kit whole transcriptome products
- This workflow has not been tested using the Evercode WT Mega kit

# Workflow overview for preparing Kinnex single-cell libraries with Parse Evercode Whole Transcriptome kits

Follow Parse Bioscience whole transcriptome kit procedure<sup>1</sup> to generate single-cell cDNA samples suitable for Kinnex library prep and HiFi sequencing on PacBio long-read systems




**Parse single-cell cDNA generation procedure overview<sup>1</sup>**

Technical note – Preparing Kinnex single-cell libraries with Parse Evercode Whole Transcriptome kits (102-326-610)

**Technical note**  
PREPARING KINNEX® SINGLE-CELL LIBRARIES WITH PARSE BIOSCIENCES EVERCODE™ WT KITS

## Parse Evercode single-cell cDNA generation workflow<sup>1</sup> (Step 1)



**Parse single-cell cDNA generation procedure key steps**


- Single-cell cDNA material first needs to be generated as described in the Parse protocols<sup>1</sup> for either the **Evercode WT** or **Evercode WT Mini kit** up through the end of those protocols' **Section 2**.

**Evercode WT User Guides**  
[UMWT3300](#) or [UMWT3100](#)

- Use either **15 ng** (from the WT protocol) or **60–75 ng** (from the WT Mini protocol) of the generated cDNA as input into the **PacBio Kinnex single-cell library prep protocol** ([103-254-300](#)) starting at **Step 2: TSO PCR (with custom primers<sup>2</sup>)**.

# Workflow overview for preparing Kinnex single-cell libraries with Parse Evercode Whole Transcriptome kits (cont.)

Parse single-cell cDNA samples may be used for Kinnex library construction by following *PacBio Procedure & checklist – Preparing Kinnex libraries using Kinnex single-cell RNA kit* ([103-254-300](#))



**PacBio Kinnex library construction procedure<sup>1</sup>**


Procedure & checklist – Preparing Kinnex libraries using Kinnex single-cell RNA kit ([103-254-300](#))



- Kinnex SC library prep procedure key steps**
- Follow **PacBio Kinnex single cell protocol** ([103-254-300](#)) (**Steps 2 – 7**) with modifications below:
    - Generate amplified single-cell cDNA per Parse Biosciences protocols. **Complete the Evercode WT or Evercode WT Mini protocol up through Section 2.**
    - Carry out cDNA quantification and sizing QC as described in the Parse Biosciences user manual.
    - Begin **Kinnex single-cell protocol from step 2 (TSO PCR)** with the following modification: Follow step 2.1 and use **5 µL of Parse Biosciences-specific amplification primers** as described in Appendix 1 **instead of** Kinnex 5' or 3' capture primer mixes.<sup>1</sup>
    - After running TSO PCR program, proceed with the subsequent Kinnex workflow from cleanup with 1.5X SMRTbell cleanup beads in the TSO PCR step through the rest of the protocol.<sup>2</sup>

**1. Parse Evercode single-cell cDNA generation**

→ See **PacBio Technical note**  
 Preparing Kinnex single-cell libraries with Parse Evercode Whole Transcriptome kits ([102-326-610](#))



<sup>1</sup> This replacement is critical to adapt the Parse library termini to function with the Kinnex concatenation workflow.  
<sup>2</sup> Subsequent QC and library quantification can be carried out as described in the Kinnex single-cell protocol (step 7.21 in [103-254-300](#)).

# TSO PCR amplification primer oligo recommendations for preparing Kinnex single-cell libraries with Parse Evercode Whole Transcriptome kits


Order custom TSO PCR amplification primers as described below from your preferred oligo vendor

## Parse Biosciences-specific primers required for TSO PCR amplification step

- Order custom primers from your preferred oligo vendor
- An example of IDT code is shown at right
- Primers should be HPLC-purified
- Order these primers and reconstitute in TE buffer (10 mM Tris HCl pH 8, 0.1 mM EDTA).
- Use 10  $\mu$ M solutions of each of these primers in place of the yellow-capped Kinnex 3' capture primer mix or the red-capped Kinnex 5' capture primer mix in the Kinnex single-cell RNA kit (103-072-200)

### Parse Kinnex forward capture primer


/5Biosg/AA/ideoxyU/CTACACGACGCTCTT  
CCGATCTCAGACGTGTGCTCTTCCGATC



### Parse Kinnex reverse capture primer

AAGCAGTGGTATCAACGCAGAG

5' TSO 3'  
Parse Kinnex capture primer Rev




TSO PCR step in Kinnex single-cell RNA procedure uses a modified PCR primer (Parse Kinnex capture primer Fwd) to incorporate a biotin tag into desired cDNA products with the correct structure.<sup>1</sup>

# SMRT Link informatics and secondary analysis recommendations for Kinnex single-cell libraries generated with Parse Evercode Whole Transcriptome kits

Read Segmentation can be run in SMRT Link (or skera in the command line) to generate segmented reads

- Note: The SMRT Link single-cell Iso-Seq workflow does not currently support analysis of Parse Biosciences data.
- Instead, Read Segmentation can be run in SMRT Link (or skera in the command line) to generate segmented reads (S-reads).
- The S-reads will contain the original single-cell or single-nucleus cDNA-containing barcode information. Refer to Rebboah et al. (2021) (section Preprocessing of LR-Split-seq data)<sup>1</sup> and Parse Biosciences long-read data processing guide<sup>2</sup> (login required) for secondary analysis recommendations using command line and community tools.

## SMRT Link read segmentation workflow






# Preparing Kinnex single-cell libraries with Twist exome enrichment kits

# Preparing Kinnex single-cell libraries with Twist exome enrichment kits: Getting started

Application-specific educational literature


Application-specific protocol documentation

Application-specific technical overviews




PacBio literature website [Link]

Application-specific brochures, informational guides and other product literature containing best practices recommendations for library preparation and data analysis workflows.




Technical note – Preparing Kinnex single-cell libraries with Twist exome enrichment kits (102-326-633) [PacBio]

Technical documentation containing a modified version of Twist Bioscience hybrid capture target enrichment protocol for generating enriched cDNA libraries for PacBio sequencing.



Procedure & checklist – Preparing Kinnex libraries using Kinnex single-cell RNA kit (103-254-300) [PacBio]

Technical documentation containing PacBio SMRTbell library construction details.



Technical overview: Kinnex library preparation using Kinnex single-cell RNA kit (103-344-600)

Technical overview presentations describe sample preparation details for constructing HiFi libraries for specific applications. Example sequencing performance data for a given application are also summarized.

# Technical reference for preparing Kinnex single-cell libraries with Twist exome enrichment kits

Technical note – *Preparing Kinnex single-cell libraries with Twist exome enrichment kits* ([102-326-633](#)) describes a modified workflow for preparing PacBio Kinnex single-cell libraries with Twist exome enrichment kits

**Technical note**  
**PREPARING KINNEX™ SINGLE-CELL LIBRARIES WITH TWIST EXOME ENRICHMENT KITS**

**Overview**

This technical note describes the experimental conditions of a modified workflow preparing PacBio® Kinnex single-cell libraries with Twist exome enrichment kits. The purpose of using exome enrichment is to enrich for spliced mRNA in single-nuclei libraries<sup>1</sup>. The following workflow has been tested on *10x Single Cell 3' Gene Expression* (v3.1) libraries, but is also compatible with single-nuclei libraries generated with the following:

- *10x Single Cell Gene Expression 3'* (v3.1, v4)
- *10x Single Cell Immune Profiling 5'* (v2, v3)
- Parse Evercode WT and WT mini kits (see the Parse technical note)

The workflow has been tested with the following Twist Exome Enrichment Panels:

- Twist Human Comprehensive Exome Panel
- Twist RNA Exome Panel


**NOTE:** This is a single reaction workflow. This protocol does not support pooled indexed libraries such as those generated using Twist's UDI Adapter Systems.

**Required materials**  
 Please refer to the Kinnex single-cell protocol and respective single-cell provider manual for materials required.

Single-cell cDNA generation	
Single-nuclei cDNA library	Refer to single-cell vendor (10x, Parse)
10x cDNA primer*	10x 2000089
PacBio Kinnex library generation	
Kinnex single-cell RNA kit	Refer to Kinnex single-cell protocol
Custom 10x or Parse blockers†	Refer to Appendix 1
Twist exome enrichment protocol	
Twist Comprehensive Exome Panel	Twist 102031: 2 rxn Twist 102032: 12 rxn Twist 102033: 96 rxn
OR	
Twist RNA Exome Panel	Twist 107143: 2 rxn Twist 107144: 12 rxn Twist 107146: 96 rxn
Twist Universal Blockers	Twist 100856: 2 rxn Twist 100578: 12 rxn Twist 100767: 96 rxn
Twist Standard Hyb and Wash Kit v2 with Amp Mix	Twist 105559: 2 rxn Twist 105560: 12 rxn Twist 105561: 96 rxn

\* The 10x cDNA primer is included in 10x Single Cell Gene Expression kits. If stock kit primers are not available, please refer to Appendix 1 for custom ordering.  
 † Custom blockers are required to replace the Twist Universal Blockers to enable Kinnex compatibility. Please refer to Appendix 2 for custom blockers.

**PacBio**



**PacBio Procedure & checklist**  
 Preparing Kinnex libraries using Kinnex single-cell RNA kit ([103-254-300](#))

**PacBio Technical note**  
 Preparing Kinnex single-cell libraries with Twist exome enrichment kits ([102-326-633](#))  
 [ Select to use either Twist Human Comprehensive Exome Panel or Twist RNA Exome Panel ]

**PacBio Procedure & checklist**  
 Preparing Kinnex libraries using Kinnex single-cell RNA kit ([103-254-300](#))


## Overview

- The purpose of using Twist exome enrichment is to enrich for spliced mRNA in single-nuclei libraries<sup>1</sup>
- This workflow has been tested on *10x Single Cell 3' Gene Expression* (v3.1) libraries, but is also compatible with single-nuclei libraries generated with the following:
  - *10x Single Cell Gene Expression 3'* (v3.1, v4)
  - *10x Single Cell Immune Profiling 5'* (v2, v3)
  - Parse Biosciences Evercode WT and WT mini whole transcriptome kits (see Parse Technical note [102-326-610](#))
- The workflow has been tested with the following Twist Exome Enrichment Panels:
  - Twist Human Comprehensive Exome Panel
  - Twist RNA Exome Panel
- **NOTE:** This is a **single reaction** workflow. This protocol **does not** support pooled indexed libraries such as those generated using Twist's UDI Adapter Systems



# Workflow overview for preparing Kinnex single-cell libraries with Twist exome enrichment kits (cont.)

Follow PacBio Kinnex library prep protocol and modified Twist target enrichment procedure<sup>1</sup> to generate single-cell hybrid capture libraries suitable for HiFi sequencing on PacBio long-read systems



**Modified Twist hybrid capture procedure overview<sup>1</sup>**

Technical note – Preparing Kinnex single-cell libraries with Twist exome enrichment kits (102-326-633)


**Technical note**  
PREPARING KINNEX™ SINGLE-CELL LIBRARIES WITH TWIST EXOME ENRICHMENT KITS

**Overview**

This technical note provides the experimental conditions for modified workflow preparing PacBio Kinnex single-cell libraries with Twist exome enrichment kits. The procedure of single-cell library preparation is based on the modified Twist exome enrichment kit. The following procedure will generate a single-cell hybrid capture library suitable for HiFi sequencing on PacBio long-read systems.

**Required materials**

Material	Quantity
Input cDNA	100 ng
TSO PCR mix	100 μL
TSO artifact removal kit	100 μL
PCR to increase cDNA mass	100 μL
Twist exome capture & on-bead PCR	100 μL
Kinnex PCR	100 μL
Kinnex array formation	100 μL
DNA damage repair	100 μL
Nuclease treatment and cleanup	100 μL



**Twist exome capture & on-bead PCR procedure key steps<sup>1</sup>**

- **Step 5.1: Twist exome capture setup**
  - Concentrate amplified cDNA sample (generated previously in Step 4) with 1.5X SMRTbell cleanup beads and resuspend in 12 μL of hybridization blocker solution.
- **Step 5.2: Twist probe hybridization**
  - Prepare probe solution and incubate the hybridization reaction at 70°C for 16 hours in a thermal cycler with the lid at 85°C.
- **Step 5.3: Bind hybridization capture reaction to Streptavidin beads**
  - Prepare streptavidin beads and bind targets by adding hybridization reaction solution to the beads. Incubate for 5 minutes at 68°C. Agitation is not required. Do not vortex.
- **Step 5.4: On-bead PCR amplification of Twist enriched library**
  - Prepare PCR mix and keep sample on ice until the thermal cycler lid has heated to 105°C. After PCR amplification is complete, perform cleanup with 1.5X SMRTbell cleanup beads and evaluate DNA quantification QC using Qubit assay.
- If sample quantity is acceptable (**25 ng is required**), proceed to the **Kinnex PCR step of the Kinnex single-cell protocol (103-254-300)** and complete the remainder of the library construction.

# PCR amplification primer and hybridization blocker oligo recommendations for preparing Kinnex single-cell libraries with Twist exome enrichment kits

For optimal library prep yields and on-target hybrid capture performance, order custom primers as described below from your preferred oligo vendor

## 10x cDNA primer<sup>1</sup> required for PCR amplification step to increase cDNA mass [ Step 4 ]

- Order custom primers from your preferred oligo vendor.
- Standard desalt primers are acceptable.
- For each FWD and REV primer, add 10 µL of 10 µM to the PCR reaction to achieve a final concentration of 1 µM.

### 10x cDNA primer FWD

CTACACGACGCTCTTCCGATCT

### 10x cDNA primer REV

AAGCAGTGGTATCAACGCAGAG

## Custom 10x or Parse blocker oligos<sup>2</sup> required for preparing samples for Twist hybridization capture reaction [ Step 5.1 ]

- Order custom primers from your preferred oligo vendor.
- An example of IDT code is shown here.
- Blocker oligo primers should be HPLC-purified.
- Dilute blocker oligos to 10 µM and add 2 µL of each FWD and REV blocker.

### 10x Custom Blocker FWD

CTACACGACGCTCTTCCGATCT/3SpC3/

### 10x Custom Blocker REV

AAGCAGTGGTATCAACGCAGAG/3SpC3/

### Parse Custom Blocker FWD

CTACACGACGCTCTTCCGATCTCAGACGTGTGCTCTTCCGATC/3SpC3/

### Parse Custom Blocker REV

AAGCAGTGGTATCAACGCAGAG/3SpC3/

# Targeted enrichment of Kinnex full-length RNA libraries with IDT xGen hybridization and wash kit v3

# Targeted enrichment of Kinnex full-length RNA libraries with IDT xGen hybridization and wash kit v3: Getting started

Application-specific educational literature

Application-specific protocol documentation

Application-specific technical overviews

The screenshot shows the PacBio literature website interface. It features a search bar with the keyword 'Kinnex' entered. Below the search bar, there are several filter categories: 'Filter by Asset Type' (Technical note, Application and product notes, Brochure, Whitepapers), 'Filter by Application' (RNA sequencing, Sample + library preparation, Aggregomics, Biodiversity, Microbial sequencing methods, Plant + animal biology, Targeted sequencing), and 'Filter by Product' (Kinnex). A list of articles is displayed, including 'Application note – Kinnex single-cell RNA kit for single-cell isoform sequencing', 'Technical note – Preparing Kinnex single-cell libraries with Twist exome enrichment kits', 'Brochure – Kinnex', 'Technical note – Concatenating cDNA or gDNA amplicons using Kinnex kits to increase throughput', 'Technical note – Targeted enrichment of Kinnex full-length RNA libraries with IDT xGen hybridization and wash kit', 'Technical note – Preparing Kinnex single-cell libraries with Parse Exorcise Whole Transcriptome kits', 'Whitepaper – Transforming plant and animal research with full-length and single-cell sequencing', and 'Application note – Kinnex 16S rRNA kit for full-length 16S sequencing'.

The technical note document is titled 'Targeted enrichment of Kinnex full-length RNA libraries with IDT xGen hybridization and wash kit v3'. It includes an overview, a list of steps for generating and enriching libraries, a section on required materials and equipment, and a section on library preparation. The steps are: 1. Generate full-length cDNA using the Iso-Seq express 2.0 kit. 2. Enrich for targeted genes using the IDT xGen Hyb Wash v3 protocol. 3. Continue through the Kinnex PCR and pooling step (step 4) in the Kinnex full-length RNA protocol. 4. The targeted Kinnex library is compatible for sequencing on Sequel II/IIIe, Revio\*, and Vega\* systems. The required materials and equipment section lists: a comprehensive list of materials and equipment (refer to the PacBio Kinnex full-length RNA protocol and IDT xGen protocol), and for cDNA generation, use the Iso-Seq express 2.0 kit which is compatible with the IDT targeted enrichment and Kinnex full-length RNA library construction. The input to Iso-Seq express 2.0 is 300 ng of total RNA, ideally with RIN > 7. Lower-quality RNA may result in shorter cDNA, while total RNA lower than 100 ng may result in the incomplete capture of the full diversity of the targeted genes. The library preparation section states: Library preparation will follow the Kinnex full-length RNA protocol with modifications outlined here.

The procedure & checklist document is titled 'Preparing Kinnex™ libraries using the Kinnex full-length RNA kit'. It includes an overview, a 'Before you begin' section, and a table of RNA input requirements. The overview states: This procedure describes the workflow for constructing Kinnex full-length RNA libraries from total RNA samples for sequencing on PacBio® Sequel® II, Sequel IIIe, Vega\*, and Revio\* systems. The 'Before you begin' section includes a table with the following information: Samples (1-24), Workflow time (1-5 days for up to 24 samples), Number of SMART® Cells per Kinnex library Prep (x4 SMART Cells for Revio using SPRQ\* chemistry, x2 SMART Cells for the Vega or Revio (non-SPRQ) systems, x4 SMART Cells for Sequel II/IIIe systems). The RNA input table includes: Quality/size distribution (RIN (RNA integrity number) > 7.0), Quantity (300 ng per library (minimum concentration 43 ng/µL per library)).

The technical overview presentation slide is titled 'Technical overview – Kinnex library preparation using Kinnex single-cell RNA kit'. It includes the following information: Sequel II and IIe systems ICS v11.0, Revio system ICS v13.3, SMRT Link v25.3, and the contact information: PN 103-344-600 Rev 02 | December 2025.

The example sequencing performance data is presented in a table format. It includes a 'Raw Data Report' showing a heatmap of read counts, a 'HiFi Read Length' histogram, and a 'Read Segmentation Metrics' table. The 'Raw Data Report' shows a heatmap of read counts with a color scale from 0 to 1000. The 'HiFi Read Length' histogram shows a distribution of read lengths with a mean of 14.8 kb. The 'Read Segmentation Metrics' table includes the following data: Raw Data Yield (1.4 TB), HiFi Reads (16.8 B), Segmented Reads (12.8 million), Mean HiFi Read Length (14.8 kb), Percent of reads with Top Priority (94.0%), Mean HiFi Read Quality (Q32), Mean HiFi Reads per Transcript (7), Percent of reads with Top Priority (94.0%), and Mean HiFi Read Length (14.8 kb).

PacBio literature website [\[Link\]](#)

Application-specific brochures, informational guides and other product literature containing best practices recommendations for library preparation and data analysis workflows.

Technical note – Targeted enrichment of Kinnex full-length RNA libraries with IDT xGen hybridization and wash kit v3 [\(102-326-632\)](#) [\[PacBio\]](#)

Technical documentation describing the procedure for targeted enrichment using IDT xGen Hyb and Wash kit v3 for preparing Kinnex full-length RNA libraries for PacBio sequencing.

Procedure & checklist – Procedure & checklist – Preparing Kinnex libraries using Kinnex full-length RNA kit [\(103-238-700\)](#) [\[PacBio\]](#)

Technical documentation containing PacBio SMRTbell library construction details.

Technical overview: Kinnex library preparation using Kinnex full-length RNA kit [\(103-344-700\)](#)

Technical overview presentations describe sample preparation details for constructing HiFi libraries for specific applications. Example sequencing performance data for a given application are also summarized.

# Technical reference for targeted enrichment of Kinnex full-length RNA libraries with IDT xGen hybridization and wash kit v3

*Technical note – Targeted enrichment of Kinnex full-length RNA libraries with IDT xGen hybridization and wash kit (102-326-632)* describes the procedure for targeted enrichment using IDT xGen Hyb and Wash kit v3 for Kinnex full-length RNA libraries

Technical note

## TARGETED ENRICHMENT OF KINNEX FULL-LENGTH RNA LIBRARIES WITH IDT xGEN HYBRIDIZATION AND WASH KIT V3

**Overview**

For RNA sequencing, targeted enrichment is considered a favorable alternative to whole transcriptome sequencing because it can capture more isoform diversity of targeted genes, requiring less sequencing depth and enabling higher sample multiplexing. The PacBio® Iso-Seq® method for full-length RNA sequencing has been used to capture single genes<sup>1</sup>, highly similar paralogs or pseudogenes<sup>2</sup>, and low-abundance long non-coding RNAs for annotation purposes<sup>3</sup>. Targeted enrichment may be the right choice when the goal of sequencing is to:

- Identify and quantify transcripts from a known set of genes
- Characterize alternative 5' starts and 3' ends
- Detect lowly expressed transcripts
- Maximize sample multiplexing capacity
- Reduce bioinformatics analysis time


This technical note describes the procedure for targeted enrichment using IDT xGen Hyb and Wash kit v3 for Kinnex® full-length RNA libraries.

**Required materials and equipment**

- For a comprehensive list of materials and equipment, please refer to the PacBio Kinnex full-length RNA protocol<sup>4</sup> and the IDT xGen protocol<sup>5</sup>.
- For cDNA generation, use the Iso-Seq express 2.0 kit, which is compatible with the IDT targeted enrichment and Kinnex full-length RNA library construction. The input to Iso-Seq express 2.0 is 300 ng of total RNA, ideally with RIN ≥ 7. Lower-quality RNA may result in shorter cDNA, while total RNA lower than 100 ng may result in the incomplete capture of the full diversity of the targeted genes.

**Library preparation**


Library preparation will follow the Kinnex full-length RNA protocol with modifications outlined here.



**Figure 1.** Library workflow for generating targeted Kinnex full-length RNA libraries using the IDT xGen Hyb Wash v3 kit.

**SMRT® Link informatics and analysis recommendations**

Targeted Kinnex full-length RNA data can be analyzed the same way as whole transcriptome Kinnex datasets using the SMRT Link Read Segmentation and Iso-Seq workflow or through the command line. Note, however, that these workflows do not have specific analyses for targeted gene lists and will output gene and isoform information for all detected transcripts.

**PacBio Procedure & checklist**  
 Preparing Kinnex libraries using the Kinnex full-length RNA kit (103-238-700)

**PacBio Technical note**  
 Targeted enrichment of Kinnex full-length RNA libraries with IDT xGen hybridization and wash kit v3 (102-326-632)<sup>1</sup>


**PacBio Procedure & checklist**  
 Preparing Kinnex libraries using the Kinnex full-length RNA kit (103-238-700)

## Overview

- For RNA sequencing, **targeted enrichment** may be the right choice when the goal of sequencing is to:
  - Identify and quantify transcripts from a known set of genes
  - Characterize alternative 5' starts and 3' ends
  - Detect lowly expressed transcripts
  - Maximize sample multiplexing capacity
  - Reduce bioinformatics analysis time
- Technical note – Targeted enrichment of Kinnex full-length RNA libraries with IDT xGen hybridization and wash kit v3 (102-326-632)** describes the procedure for targeted enrichment using IDT xGen Hyb and Wash kit v3 for Kinnex full-length RNA libraries

# Workflow overview for targeted enrichment of Kinnex full-length RNA libraries with IDT xGen hybridization and wash kit v3

Follow PacBio Kinnex library prep protocol and IDT target enrichment procedure to generate full-length RNA hybrid capture libraries suitable for HiFi sequencing on PacBio long-read systems



**PacBio Kinnex library construction procedure<sup>1</sup>**

Procedure & checklist – Preparing Kinnex libraries using the Kinnex full-length RNA kit ([103-238-700](https://www.pacb.com/knowledge-base/preparing-kinnex-libraries-using-the-kinnex-full-length-rna-kit/))

Preparing Kinnex<sup>®</sup> libraries using Kinnex single-cell RNA kit

Procedure & checklist

Before you begin

This procedure describes the workflow for constructing single-cell RNA libraries from 10<sup>3</sup> Chromium V3 V2 2Dx using the Kinnex full-length RNA kit for targeted enrichment of PacBio Kinnex<sup>®</sup> libraries. It is intended for use with single-cell cDNA generated using the 10<sup>3</sup> Chromium V3 V2 2Dx Single Cell 2.0 Kit or the Chromium Next Single Cell V3 V2. For a full list of compatibility with this product, please consult the datasheet.

Parameter	Value
Sequencing kit	10 <sup>3</sup> Chromium V3 V2 2Dx
Sequencing flow	10 <sup>3</sup> Chromium V3 V2 2Dx
Number of single-cell runs per lane	48,000 Cells (to Kinnex <sup>®</sup> full-length RNA library)
Library preparation	10 <sup>3</sup> Chromium V3 V2 2Dx Single Cell 2.0 Kit or the Chromium Next Single Cell V3 V2


Library prep

Parameter	Value
Library prep	10 <sup>3</sup> Chromium V3 V2 2Dx
Library prep	10 <sup>3</sup> Chromium V3 V2 2Dx
Library prep	10 <sup>3</sup> Chromium V3 V2 2Dx
Library prep	10 <sup>3</sup> Chromium V3 V2 2Dx

A full list of compatibility with this product, please consult the datasheet.


**Kinnex full-length RNA library prep procedure key steps**

- For cDNA generation, use the Iso-Seq express 2.0 kit with 300 ng of total RNA input, ideally with RIN ≥ 7.<sup>1</sup>
- Library preparation will follow the **Kinnex full-length RNA protocol** ([103-238-700](https://www.pacb.com/knowledge-base/preparing-kinnex-libraries-using-the-kinnex-full-length-rna-kit/)) with modifications outlined here:
  - Generate full-length cDNA using the Iso-Seq express 2.0 kit by following **Steps 1-3** in the Kinnex full-length RNA protocol.
  - Enrich for targeted genes using the **IDT xGen Hyb Wash v3 protocol** (RUO24-3081\_001; also see **IDT xGen Technical note** ([102-326-632](https://www.idtdna.com/help-center/entry-detail/102-326-632))).
  - Continue through the **Kinnex PCR step** in the Kinnex full-length RNA protocol (page 13 Step 4) and complete the remainder of the library construction workflow



# Workflow overview for targeted enrichment of Kinnex full-length RNA libraries with IDT xGen hybridization and wash kit v3 (cont.)

Follow PacBio Kinnex library prep protocol and IDT target enrichment procedure to generate full-length RNA hybrid capture libraries suitable for HiFi sequencing on PacBio long-read systems



**IDT xGen hybrid capture procedure overview<sup>1</sup>**

Technical note – Targeted enrichment of Kinnex full-length RNA libraries with IDT xGen hybridization and wash kit v3 (102-326-632)


**Technical note**  
TARGETED ENRICHMENT OF KINNEK FULL-LENGTH RNA LIBRARIES WITH IDT xGEN HYBRIDIZATION AND WASH KIT V3

**Overview**

For this sequencing, targeted enrichment is considered a library preparation step. The IDT xGen Hybridization and Wash Kit v3 (EAP) for PacBio Targeted Kinnex protocol (RUO24-3081\_001) is used to enrich full-length RNA libraries for targeted enrichment. The IDT xGen Hybridization and Wash Kit v3 (EAP) for PacBio Targeted Kinnex protocol (RUO24-3081\_001) is used to enrich full-length RNA libraries for targeted enrichment. The IDT xGen Hybridization and Wash Kit v3 (EAP) for PacBio Targeted Kinnex protocol (RUO24-3081\_001) is used to enrich full-length RNA libraries for targeted enrichment.

**IDT xGEN hybrid capture procedure key steps<sup>1</sup>**

- Add Blocker Master Mix to each cDNA sample (500 ng<sup>2</sup>) and dry down
- Add Hybridization Master Mix to each sample well and incubate for desired hybridization time.<sup>3</sup>
- Prepare hybridization wash buffer and capture beads within 2 hours before use
- Add (streptavidin) capture beads to sample (30 min)
- Perform 3 rounds of washing with wash buffers
- Perform post-capture PCR reaction + cleanup
- Perform DNA quantification QC using Qubit assay and DNA sizing QC using Agilent TapeStation or other equivalent system



<sup>1</sup> For hybridization capture workflow details, refer to the IDT xGen Hybridization and Wash Kit v3 (EAP) for PacBio Targeted Kinnex protocol (RUO24-3081\_001).  
<sup>2</sup> 100 ng to 6 µg cDNA inputs can be used for short hybridizations (1-2 hrs). Generally, we recommend not to exceed 2.5 µg total input for 1-hr hybridization time for very large panels (>30 Mb).  
<sup>3</sup> This procedure supports short hybridization times of only 1 hr as well as the standard 4- or 16-hr incubations. We recommend users empirically test the shorter hybridization time for their workflows.


# SMRT Link informatics and secondary analysis recommendations for targeted Kinnex full-length RNA data

Use SMRT Link Read Segmentation and Iso-Seq application to perform isoform-classification analysis to identify novel genes & isoforms with abundance information

- Targeted Kinnex full-length RNA data can be analyzed the same way as whole transcriptome Kinnex datasets using the [SMRT Link Read Segmentation and Iso-Seq workflow](#) or through the command line.
- Note, however, that these workflows do not have specific analyses for targeted gene lists and will output gene and isoform information for all detected transcripts.



## SMRT Link read segmentation and Iso-Seq workflow<sup>1</sup>





# APPENDIX 3: SMRT Link run design procedure for Sequel IIe system

# SMRT Link Run Design procedure for Sequel I/II system

## Sample information and run information

- Select desired **Kinnex application** from the **Application** field drop-down menu
- The following fields are **auto-populated** with default recommended values and high-lighted in **green**:

- SMRTbell Adapter Design**
  - SMRTbell Kinnex Prep Kit
- Binding Kit**
  - Sequel II Binding Kit 3.2
- Sequencing Kit**
  - Sequel II Sequencing Plate 2.0 (4 rxn or 1 rxn)
- DNA Control Complex**
  - Sequel II DNA Internal Control Complex 3.2
- Movie Time per SMRT Cell**
  - 30 hrs
- Pre-Extension Time**
  - 2 hrs

**SMRTbell Adapter Design** field determines which adapter finding algorithm is used during post-primary analysis<sup>1</sup>

Default SMRTbell adapter design for Kinnex samples is **SMRTbell Kinnex Prep Kit**

Recommended OPLC for Sequel II/IIe Kinnex library samples is **40 – 60 pM**

Recommended movie time = **30 hrs**

Select desired Kinnex application type from drop-down menu

**SAMPLE 1: Kinnex full-length RNA library demo , A01, 30 hour movie, 16000 bp insert** [Copy] [Delete]

Import from Sample Setup [Select Sample]

**Application** (Required): Kinnex full-length RNA

**Well Sample Name** (Required): Kinnex full-length RNA library demo

**Bio Sample Name** (Required):

**Sample Comment**:

**Sample Well**: A01

**SMRTbell Adapter Design** (Required): SMRTbell® Kinnex Prep Kit

**Binding Kit** (Required): Sequel® II Binding Kit 3.2

**Sequencing Kit** (Required): Sequel® II Sequencing Plate 2.0 (4 rxn)

**DNA Control Complex**: Sequel® II DNA Internal Control Complex 3.2

**Insert Size (bp)** (Required): 16000

**Recommended Concentration on Plate (pM)**: 40 – 60 pM

**On-Plate Loading Concentration (pM)** (Required): 45

**Movie Time per SMRT Cell (hours)**: 30

**Use Pre-Extension**:  YES  NO

**Pre-Extension Time (hours)**: 2

CCS Analysis will be performed on-instrument to produce HiFi .bam files.

Example sample information entered into a Sequel I/II system run design worksheet for a Kinnex full-length RNA library sample.

# SMRT Link Run Design procedure for Sequel Ii system (cont.)

## Advanced options

- For all Kinnex library samples, leave the following **Advanced Options** fields at their **default settings**
  - Use Adaptive Loading**  
→ YES
  - Loading Target (P1 + P2)**  
→ 0.85
  - Maximum Loading Time**  
→ 2 hours
  - CCS Analysis Output - Include Low Quality Reads**  
→ NO
  - CCS Analysis Output - Include Kinetics Information**  
→ NO
  - Pre-Extension Time**  
→ 2 hrs
- If desired, specify to use an alternative project folder for the **Add Data to Project** field

Advanced Options

Use Adaptive Loading  YES  NO

Loading Target (P1 + P2) 0.85

Maximum Loading Time (hours) 2

CCS Analysis Output - Include Low Quality Reads  YES  NO

CCS Analysis Output - Include Kinetics Information  YES  NO

Add Data to Project

Example default Advanced Options settings entered into a Sequel Ii system run design worksheet for a Kinnex full-length RNA library sample.

# SMRT Link Run Design procedure for Sequel Ii system (cont.)

## Barcoded sample options

- For Kinnex library samples, can leave most **Barcoded Sample Options** fields at their **default settings**

Specify Bio Sample Names, either interactively or by downloading a CSV file (**Interactively** or **From a file**)

If desired, specify to perform barcode demultiplexing on-instrument or in SMRT Link (default = **On-instrument** for Sequel Ii system)

Can leave most of these fields at their default values

Barcoded Sample Options

Sample Is Barcoded  YES  NO

Barcode Set Required MAS SMRTbell barcoded adapters (v2)

Same Barcodes on Both Ends of Sequence ?  YES  NO

Assign Bio Sample Names to Barcodes ? Required

Demultiplex Barcodes  ON INSTRUMENT  IN SMRT LINK  DO NOT GENERATE

Example default Barcoded Sample Options settings entered into a Sequel Ii system run design worksheet for a Kinnex full-length RNA library sample.



[www.pacb.com](http://www.pacb.com)

Revision history (description)	Version	Date
Updated example Kinnex sequencing performance data for PacBio long-read systems (Revio system + SPRQ chemistry and Vega system). Also added information about PacBio compatible Kinnex library preparation workflows using third-party kit products.	02	December 2025

Research use only. Not for use in diagnostic procedures. © 2025 Pacific Biosciences of California, Inc. (“PacBio”). All rights reserved. Information in this document is subject to change without notice. PacBio assumes no responsibility for any errors or omissions in this document. Certain notices, terms, conditions and/or use restrictions may pertain to your use of PacBio products and/or third-party products. Refer to the applicable PacBio terms and conditions of sale and to the applicable license terms at [pacb.com/license](http://pacb.com/license). Pacific Biosciences, the PacBio logo, PacBio, Circulomics, Omniome, SMRT, SMRTbell, Iso-Seq, Sequel, Nanobind, SBB, Revio, Onso, Apton, Kinnex, PureTarget, SPRQ, and Vega are trademarks of PacBio.