



Technical overview – Kinnex library preparation using Kinnex full-length RNA kit

Sequel II and IIe systems ICS v11.0

Revio system ICS v13.3

SMRT Link v25.3



PN 103-344-700 Rev 02 | December 2025

Technical overview

Kinnex library preparation using Kinnex full-length RNA kit

1. Kinnex full-length RNA method overview
2. Kinnex full-length RNA library preparation workflow details
3. Kinnex full-length RNA sequencing preparation workflow details
4. Kinnex full-length RNA example sequencing performance data
5. Kinnex full-length RNA data analysis workflow overview
6. Technical documentation & applications support resources
7. APPENDIX 1 – Concatenating cDNA or gDNA amplicons using Kinnex kits to increase throughput
8. APPENDIX 2 – PacBio compatible Kinnex library preparation workflows
 - i. Targeted enrichment of Kinnex full-length RNA libraries with IDT xGen hybridization and wash kit v3
9. APPENDIX 3 – SMRT Link run design procedure for Sequel IIe system

Kinnex library preparation using Kinnex full-length RNA kit: Getting started



PacBio

Application note

Kinnex full-length RNA kit for isoform sequencing

Introduction

Alternative splicing (AS) in eukaryotic species generates functional diversity by expressing different combinations of exons in the same gene. Accurate characterization of full-length transcript isoforms generated by AS is critical for biological and disease studies. Bulk RNA-Seq using short reads cannot fully resolve isoform structures, as the complex nature of AS prohibits unambiguous transcript assembly with even the most sophisticated computational tools ([Spatz et al., 2019](#)). Long-read RNA-Seq using PacBio® technology (the Iso-Seq® method) eliminates the need for transcript assembly by sequencing full-length cDNAs and enables new discoveries across many applications (Figure 1).

PacBio RNA sequencing has been shown to offer the longest read lengths and highest quality reads that are required for isoform discovery and isoform quantification ([Zhang-Pajonk et al., 2018](#)). The Kinnex™ full-length RNA kit takes full-length cDNA generated by the Iso-Seq express 2.0 kit as input and outputs a sequencing-ready library that results in an 8-fold throughput increase compared to typical Iso-Seq libraries. Combined with the Read Segmentation & Iso-Seq analysis in SMRT™ Link software, PacBio offers cost-effective isoform sequencing that does not require orthogonal sequencing methods. SMRT Link software produces an isoform classification report with abundance information that can be used by [Isoform analysis tools](#).

© 2024 PacBio. All rights reserved. Research use only. Not for use in diagnostic procedures.
102-326-591 REV03 NOV2024

PacBio

Application note – Kinnex full-length RNA kit for isoform sequencing (102-326-591)

Summary overview of application-specific library preparation and data analysis workflow recommendations.

Preparing Kinnex™ libraries using the Kinnex full-length RNA kit

PacBio

Procedure & checklist

Before you begin

This procedure describes the workflow for constructing Kinnex® full-length RNA libraries from total RNA samples for sequencing on PacBio® Sequel II®, Sequel Ie, Vega®*, and Revio®* systems.

| Overview | |
|--|---|
| Samples | 1–24 |
| Workflow time | 1.5 days (for up to 24 samples) |
| Number of SMART® Cells per Kinnex library Prep | >8 SMART Cells for Revio using SPRQ™ chemistry |
| | >2 SMART Cells for the Vega or Revio (non-SPRQ) systems |
| | >4 SMART Cells for Sequel II/Ie systems |
| RNA input | |
| Quality/size distribution | RIN (RNA integrity number) ≥7.0 |
| Quantity | 300 ng per library (minimum concentration 43 ng/μL per library) |

* Not supported by PacBio

© 2025 PacBio. All rights reserved. Research use only. Not for use in diagnostic procedures.
109-238-700 REV07 JUN2025

Procedure & checklist – Preparing Kinnex libraries using the Kinnex full-length RNA kit (103-238-700)

Technical documentation containing application-specific library preparation protocol details.

[illegible]

Technical Overview – Kinnex library preparation using Kinnex full-length RNA kit (103-344-700)

Technical overview presentations describe sample preparation details for constructing Kinnex HiFi libraries for specific applications. Example sequencing performance data for a given application are also summarized.

cDNA synthesis & amplification

300 ng input total RNA per sample
RIN (RNA integrity number) ≥ 7.0
Sample multiplexing options available

Kinnex library preparation (Kinnex PCR 8-fold kit + Kinnex concatenation kit)

Use amplified cDNA to generate Kinnex library containing 8-segment array

Sample multiplexing up to 48-plex supported

SMRT sequencing

(Sequel Ile, Vega & Revio systems)

Perform ABC¹ and sequence Kinnex libraries on PacBio long-read systems

Data analysis (SMRT Link)

Use SMRT Link Read Segmentation data utility to split arrayed transcript HiFi reads


Use SMRT Link Iso-Seq analysis application to identify novel genes and isoforms with abundance information 3



Kinnex full-length RNA method overview

Kinnex full-length RNA method overview


Use Kinnex full-length RNA kit to perform high-accuracy, full-length isoform sequencing with PacBio long-read systems




- Input 300 ng total RNA, RIN ≥ 7
- Generate up to 12-plex barcoded cDNA using **Iso-Seq express 2.0 kit (103-071-500)**
- 2-day Kinnex library preparation using **Kinnex full-length RNA kit (103-072-000)**¹
- SMRT Link Run Design support for 'Kinnex full-length RNA' application type with auto-analysis (read segmentation + isoform analysis)
- SMRT Link Iso-Seq isoform-classification software to identify novel genes and isoforms with abundance information

Kinnex full-length RNA method overview (cont.)

Full-length cDNA synthesis & amplification



Kinnex library prep, sequencing & analysis



¹ Twelve barcoded Iso-Seq primers (Iso-Seq primer bc01–12) are available for cDNA amplification step.

² Kinnex adapter barcode sequences can be downloaded from [SMRT Link](#) Data Management module.

Kinnex full-length RNA library preparation procedure description

Procedure & checklist – Preparing Kinnex libraries using the Kinnex full-length RNA kit ([103-238-700](#)) describes the workflow for constructing Kinnex libraries from total RNA samples using the **Iso-Seq express 2.0 kit** and **Kinnex full-length RNA kit** for sequencing on PacBio long-read systems

| Overview | |
|---|--|
| Samples | 1–24 |
| Workflow time | 1.5 days (for up to 24 samples) |
| Number of SMRT® Cells per Kinnex library Prep | >8 SMRT Cells for Revio using SPRQ™ chemistry >2 SMRT Cells for the Vega or Revio (non-SPRQ) systems >4 SMRT Cells for Sequel II/Ile systems |

| RNA input | |
|---------------------------|---|
| Quality/size distribution | RIN (RNA integrity number) ≥7.0 |
| Quantity | 300 ng per library (minimum concentration 43 ng/μL per library) |



Preparing Kinnex™ libraries using the Kinnex full-length RNA kit

Procedure & checklist

Before you begin

This procedure describes the workflow for constructing Kinnex full-length RNA libraries from total RNA samples for sequencing on PacBio® Sequel® II, Sequel IIe, Vega™, and Revio® systems.

| Overview | |
|---|--|
| Samples | 1–24 |
| Workflow time | 1.5 days (for up to 24 samples) |
| Number of SMRT® Cells per Kinnex library Prep | >8 SMRT Cells for Revio using SPRQ™ chemistry >2 SMRT Cells for the Vega or Revio (non-SPRQ) systems >4 SMRT Cells for Sequel II/Ile systems |
| RNA input | |
| Quality/size distribution | RIN (RNA integrity number) ≥7.0 |
| Quantity | 300 ng per library (minimum concentration 43 ng/μL per library) |

© 2025 PacBio. All rights reserved. Research use only. Not for use in diagnostic procedures.
103-238-700 REV07 JUN2025

PacBio

PacBio [Documentation](#) ([103-238-700](#))








- Kinnex full-length RNA library prep protocol uses **Kinnex full-length RNA kit** and **Iso-Seq express 2.0 kit**
→ **Do not use** SMRTbell prep kit 3.0 with this protocol

Kinnex full-length RNA kit bundle and Iso-Seq express 2.0 kit components

Kinnex full-length RNA kit bundle and Iso-Seq express 2.0 kit provide full support for Kinnex library prep workflow





Iso-Seq express 2.0 kit (103-071-500)

Includes Iso-Seq Express template switching oligo, barcoded cDNA PCR Primers, and other reagents needed for performing 1st-strand cDNA synthesis and PCR amplification of cDNA products generated from input total RNA.

| Iso-Seq express 2.0 kit components | | |
|------------------------------------|---|--|
| Component | | Description |
| 1 |  | Iso-Seq RT buffer <ul style="list-style-type: none">For 1st-strand cDNA synthesis |
| 2 |  | Iso-Seq RT primer mix <ul style="list-style-type: none">For 1st-strand cDNA synthesis |
| 3 |  | Iso-Seq RT enzyme mix <ul style="list-style-type: none">For 1st-strand cDNA synthesis |
| 4 |  | Iso-Seq template switch oligo <ul style="list-style-type: none">For 1st-strand cDNA synthesis |
| 5 |  | Iso-Seq cDNA PCR mix <ul style="list-style-type: none">Enzyme nucleotide mix for cDNA amplification |
| 6 |  | Iso-Seq cDNA amplification primer <ul style="list-style-type: none">Reverse primer for cDNA amplification |
| 7 |  | Iso-Seq primers (bc01 – bc12) <ul style="list-style-type: none">Barcoded forward primers (bc01 – bc04) for cDNA amplification |

Kinnex full-length RNA kit bundle (103-072-000)

Includes Kinnex PCR kit, Kinnex concatenation and ancillary DNA cleanup reagents needed for incorporation of Kinnex segmentation adapters and Kinnex array formation for generating Kinnex full-length RNA libraries.


| Kinnex full-length RNA kit components | | |
|---------------------------------------|---|---|
| Component | | Description |
| 1 |  | Kinnex PCR 8-fold kit (12 rxn) <ul style="list-style-type: none">Contains reagents for Kinnex PCR to incorporate segmentation adapters |
| 2 |  | Kinnex concatenation kit (12 rxn) <ul style="list-style-type: none">Contains reagents for Kinnex array formation and SMRTbell template constructionIncludes barcoded Kinnex adapter mixes (bcM0001 – bcM0004) |
| 3 |  | SMRTbell cleanup beads <ul style="list-style-type: none">For DNA cleanup |
| 4 |  | Elution buffer <ul style="list-style-type: none">For DNA cleanup |

Kinnex full-length RNA library barcoding options for sample multiplexing

Kinnex full-length RNA library preparation procedure supports up to 48-plex sample multiplexing

Kinnex full-length RNA library preparation procedure supports **up to 48-plex** sample multiplexing through combined use of:

- 12 different barcoded cDNA amplification PCR primers (bc01 – bc12)
- 4 different barcoded Kinnex terminal SMRTbell adapters (bc01 – bc04)



Multiplexing best practices

Sample multiplexing can be achieved with one of the three following methods:


1. **Barcoded cDNA primers** using Iso-Seq primers bc01–12 in step 3 of the protocol. To multiplex, use the Iso-Seq cDNA amplification primer in combination with Iso-Seq primers bc01–12 to amplify samples. After SMRTbell cleanup, Iso-Seq samples can be pooled and brought through a single Kinnex PCR reaction. Each barcoded primer is sufficient for 2 reactions, with the Iso-Seq kit supporting a total of 24 reactions.
2. **Barcoded adapters** using Kinnex adapters bc01–04. In this case, use barcoded adapters at step 5 “Kinnex array
3. formation” in the workflow.
4. **A combination of the above two approaches** to achieve 48-plex.

Note: If not performing multiplexing, the same Iso-Seq primer barcodes and Kinnex adapter barcodes are still used, but without pooling.

Kinnex full-length RNA experimental design considerations

Kinnex full-length RNA application use case recommendations for PacBio systems

| Example application | Human genetics disease studies | Biopharma for identifying highly expressed targets | Plant & animal whole genome annotation |
|---|---|---|--|
| Experimental goal | Isoform discovery and quantification of moderate-to-rare transcripts | Isoform discovery of high expressed transcripts | Comprehensive transcript annotation in a species |
| Example study design | Disease vs. normal tissues with multiple replicates | Disease cohort with >20+ samples | Plant or animal with multiple tissue types |
| Target depth of coverage per sample | 10 M reads per sample | 5 M reads per sample | 5 M reads per sample |
| Sample multiplexing ¹ | Sequel II/Ile system: Up to 2 samples per SMRT Cell 8M (2-plex) | Sequel II/Ile system: Up to 3 samples per SMRT Cell 8M (3-plex) | |
| | Vega system: Up to 3 samples per Vega SMRT Cell (3-plex) | Vega system: Up to 6 samples per Vega SMRT Cell (6-plex) | |
| | Revio system + SPRQ: Up to 6 samples per Revio SMRT Cell (6-plex) | Revio system + SPRQ: Up to 12 samples per Revio SMRT Cell (12-plex) | |
| Expected data throughput (per SMRT Cell) | Sequel II/Ile system: Up to 20 M reads divided by <i>N</i> samples | | |
| | Vega system: Up to 30 M reads divided by <i>N</i> samples | | |
| | Revio system + SPRQ: Up to 60 M reads divided by <i>N</i> samples | | |
| Kinnex library prep protocol | Procedure & checklist – Preparing Kinnex libraries using the Kinnex full-length RNA kit (103-238-700) | | |
| Total RNA input into Kinnex library prep workflow | 300 ng total RNA (RIN ≥7) for 1 st -strand cDNA synthesis | | |
| SMRT Link data analysis workflows | Read Segmentation and Iso-Seq analysis application with option to “pool reads and cluster together” to get a master isoform classification file with per-sample full-length read counts | | |
| Community data analysis tools | Annotation & quantification: PIGEON, SQANTI3 / Differential analysis: TappAS / Fusion calling: pbfusion / Visualization: SWAN | | |



¹ Kinnex concatenation kit (103-071-800) can support up to 48-plex sample multiplexing through the combined use of 12 different barcoded cDNA amplification primers and 4 different barcoded Kinnex terminal SMRTbell adapters during Kinnex full-length RNA library construction.



Kinnex full-length RNA library preparation workflow details

Procedure & checklist – Preparing Kinnex libraries using the Kinnex full-length RNA kit (103-238-700)

Procedure & checklist [103-238-700](#) describes the workflow for constructing Kinnex libraries from total RNA samples using the **Iso-Seq express 2.0 kit** and **Kinnex full-length RNA kit** for sequencing on PacBio long-read systems

Procedure & checklist contents

1. Total RNA input QC recommendations and general best practices for reagent & sample handling.
2. Enzymatic workflow steps for cDNA synthesis and amplification.
3. Enzymatic workflow steps for construction of 8-segment Kinnex arrays from amplified cDNA.
4. Enzymatic workflow steps for DNA damage repair & nuclease treatment of Kinnex libraries.
5. Workflow steps for final cleanup of Kinnex SMRTbell libraries using SMRTbell cleanup beads.
6. Sample setup ABC¹ workflow steps to prepare Kinnex SMRTbell libraries for sequencing on Revio (+SPRQ) and Vega systems.

Preparing Kinnex™ libraries using the Kinnex full-length RNA kit

Procedure & checklist

Before you begin

This procedure describes the workflow for constructing Kinnex full-length RNA libraries from total RNA samples for sequencing on PacBio® Sequel® II, Sequel IIe, Vega™, and Revio® systems.

| Overview | |
|---|--|
| Samples | 1–24 |
| Workflow time | 1.5 days (for up to 24 samples) |
| Number of SMRT® Cells per Kinnex library Prep | >8 SMRT Cells for Revio using SPRQ™ chemistry >2 SMRT Cells for the Vega or Revio (non-SPRQ) systems >4 SMRT Cells for Sequel II/Ile systems |

| RNA input | |
|---------------------------|---|
| Quality/size distribution | RIN (RNA integrity number) ≥7.0 |
| Quantity | 300 ng per library (minimum concentration 43 ng/μL per library) |


© 2025 PacBio. All rights reserved. Research use only. Not for use in diagnostic procedures.
103-238-700 REV07 JUN2025

PacBio

PacBio [Documentation](#) ([103-238-700](#))

Kinnex full-length RNA library construction workflow overview

Procedure & checklist – Preparing Kinnex libraries using the Kinnex full-length RNA kit (103-238-700)



General best practices recommendations for preparing Kinnex full-length RNA libraries

Reagent and sample handling

Kinnex library prep reagents

- Room temperature is defined as any temperature in the range of 18–23°C for this protocol.
- Take care to accurately pipette SMRTbell cleanup beads because small changes in volume can significantly alter the size distribution of your sample.
- Equilibrate the SMRTbell cleanup beads at room temperature for 30 mins prior to use.
- In cDNA amplification and Kinnex PCR, **keep sample(s) on ice until thermal cycler lid has reached 105°C** to avoid digestion of primers by polymerase exonuclease activity.
- Ensure that the DNA damage repair mix is stored at -20°C to avoid poor library performance.**
- This workflow takes ~1.5 days to complete.
 - If a stop is necessary, refer to the workflow for **safe stopping points**.

| Temperature-sensitive reagents | | |
|-------------------------------------|------------|---|
| Iso-Seq express 2.0 kit 103-071-500 | | |
| | Tube color | Reagent |
| cDNA synthesis | Purple | Iso-Seq RT buffer 103-103-900 |
| | Orange | Iso-Seq RT primer mix 103-104-000 |
| | Yellow | Iso-Seq RT enzyme mix 103-104-100 |
| | Red | Iso-Seq cDNA PCR mix 103-104-200 |
| | Blue | Iso-Seq Express TSO 2.0 103-104-300 |
| | Green | Iso-Seq cDNA amplification primer 103-104-400 |
| | | Iso-Seq primer bc01 103-104-500 |
| | | Iso-Seq primer bc02 103-104-600 |
| | | Iso-Seq primer bc03 103-104-700 |
| | | Iso-Seq primer bc04 103-104-800 |
| | | Iso-Seq primer bc05 103-104-900 |
| | White | Iso-Seq primer bc06 103-105-000 |
| | | Iso-Seq primer bc07 103-105-100 |
| | | Iso-Seq primer bc08 103-105-200 |
| | | Iso-Seq primer bc09 103-105-300 |
| | | Iso-Seq primer bc10 103-105-400 |
| | | Iso-Seq primer bc11 103-105-500 |
| | | Iso-Seq primer bc12 103-105-600 |

| Temperature-sensitive reagents | |
|-----------------------------------|----------------------------------|
| Kinnex PCR 8-fold kit 103-071-600 | |
| Tube color | Reagent |
| Green | Kinnex PCR mix 103-107-700 |
| Orange | Kinnex primer mix A 103-107-800 |
| | Kinnex primer mix B 103-107-900 |
| | Kinnex primer mix C 103-108-000 |
| | Kinnex primer mix D 103-108-100 |
| | Kinnex primer mix E 103-108-200 |
| | Kinnex primer mix F 103-108-300 |
| | Kinnex primer mix G 103-108-400 |
| | Kinnex primer mix HQ 103-108-500 |

| Temperature-sensitive reagents | |
|--------------------------------------|--|
| Kinnex concatenation kit 103-071-800 | |
| Tube color | Reagent |
| Red | Kinnex enzyme 103-110-400 |
| Yellow | Kinnex ligase 103-110-500 |
| White | Kinnex array and repair buffer 103-110-300 |
| Green | DNA repair mix 103-110-000 |
| Light Purple | Nuclease buffer 103-110-200 |
| Light Green | Nuclease mix 103-110-100 |
| Blue | Kinnex adapter bc01 mix 103-109-600 |
| | Kinnex adapter bc02 mix 103-109-700 |
| | Kinnex adapter bc03 mix 103-109-800 |
| | Kinnex adapter bc04 mix 103-109-900 |

General best practices recommendations for preparing Kinnex full-length RNA libraries

Reagent and sample handling

Sequencing prep reagents (for sample setup ABC)


- Once thawed, place reaction buffers and sequencing primer on-ice prior to making master mix. The Loading buffer should be left at room-temperature.
- Note:** The Loading buffer is **light sensitive** and should be protected from light when not in use.
- Keep the following reagents on a cold block or ice:**
 - Sequencing polymerase
 - Sequencing control
- Bring the following reagents up to room temperature 30 minutes prior to use:**
 - Loading buffer
 - SMRTbell cleanup beads

Thaw the following reagents and room temperature:

| Component | Tube color |
|--------------------------|-------------|
| Annealing buffer | Light blue |
| Kinnex sequencing primer | Light green |
| Polymerase buffer | Yellow |
| Loading buffer | Green |
| Dilution buffer | Blue |

Input RNA quality control

Input RNA quality control is highly recommended before proceeding to the MAS-Seq library prep workflow



Left: Bioanalyzer electropherogram detailing the regions that are indicative of RNA quality. **Right:** Sample electro-pherograms corresponding to different RNA Integrity Number (RIN) scores. Samples range from intact (RIN 10), to degraded (RIN 2). Images from Agilent Application Note: RNA Integrity Number (RIN) – Standardization of RNA Quality Control ([5989-1165EN](#))

- 300 ng of total RNA per library (minimum concentration 43 ng/μL per library) is required for this procedure
- Sample QC of input total RNA samples should be assessed by measuring RNA Integrity Number (RIN) using a Bioanalyzer 2100 instrument (Agilent Technology) with RNA 6000 Nano kit¹
 - RIN ≥7.0 (ideally ≥8.0) is sufficient for Kinnex full-length RNA protocol
 - Samples with RIN <7.0 can be processed, but risk of significant underperformance or even failure is greatly increased
- RIN score (1 to 10) is related to ratio of the area under 28s and 18s fragment peaks and also takes into account signal intensity above baseline in the Inter-Region and Fast Region since this is where degradation products appear
- Higher RIN numbers are correlated with better overall sample quality and lower degradation



Step 1
procedural summary



DNA
sizing QC


Total RNA RIN ≥ 7

Proceed to **Step 2**




cDNA synthesis + cleanup

In this step, total RNA samples are converted to 1st-strand cDNA




Minimum total RNA input for cDNA synthesis = 300 ng



Step 2 procedural summary


2.1 Thawing reagents for first-strand cDNA synthesis



2.2 Primer annealing for first-strand cDNA synthesis



2.3 Reverse transcription and template switching




2.4 1.3X SMRTbell bead cleanup



Proceed to **Step 3**

cDNA synthesis + cleanup (cont.)

Procedural notes



2.1 Thawing reagents for first-strand synthesis

| ✓ | Step | Instructions |
|---|-------|---|
| | 2.1.1 | Quick-spin the Iso-Seq RT enzyme mix in the microcentrifuge to collect liquid, then place on ice. |

Thaw the following components at room temperature, briefly vortex to mix, then quick-spin to collect liquid and place on ice.

| Tube color | Reagent |
|------------|-----------------------|
| Orange | Iso-Seq RT primer mix |

- After thawing specified reagents **on ice**, perform a quick spin to collect liquid, then **place on ice**

2.2 Primer annealing for first-strand synthesis

| ✓ | Step | Instructions | | | | | | |
|-------|--------------------|--|---|------------|--------|-------|--------------------|-------|
| | | For each RNA sample to be processed, prepare reagent mix 1 on ice by adding the following components to each tube in the PCR strip tube. | | | | | | |
| | | <table><tr><th>✓</th><th>Components</th><th>Volume</th></tr><tr><td>2.2.1</td><td>Total RNA (300 ng)</td><td><7 µL</td></tr></table> | ✓ | Components | Volume | 2.2.1 | Total RNA (300 ng) | <7 µL |
| ✓ | Components | Volume | | | | | | |
| 2.2.1 | Total RNA (300 ng) | <7 µL | | | | | | |

- Set up primer annealing reaction mix **on ice** and then transfer to thermal cycler for incubation

2.3 Reverse transcription and template switching

| ✓ | Step | Instructions | | | | | | | | |
|-------|---|---|------------|------------|--------|--|--------|---|------|--|
| | | For each RNA sample, prepare reagent mix 2 on ice by adding the following components in the order and volume listed below. Adjust component volumes for the number of samples being prepared, plus 10% overage. | | | | | | | | |
| 2.3.1 | <table><tr><th>✓</th><th>Tube color</th><th>Components</th><th>Volume</th></tr><tr><td></td><td>Purple</td><td>Iso-Seq RT buffer (vortex briefly before use)</td><td>5 µL</td></tr></table> | ✓ | Tube color | Components | Volume | | Purple | Iso-Seq RT buffer (vortex briefly before use) | 5 µL | |
| ✓ | Tube color | Components | Volume | | | | | | | |
| | Purple | Iso-Seq RT buffer (vortex briefly before use) | 5 µL | | | | | | | |

- Set up reverse transcription reaction mix **on ice** and then transfer to thermal cycler for incubation



2.4 1.3X SMRTbell bead cleanup

| ✓ | Step | Instructions |
|---|-------|--|
| | 2.4.1 | For each sample, add 29 µL of elution buffer to the 21 µL reverse transcription and template switching reaction (Section 2.3) for a total volume of 50 µL. |
| | 2.4.2 | Add 65 µL (1.3X v/v) of resuspended, room-temperature SMRTbell cleanup beads. |
| | 2.4.3 | Mix beads by pipetting 10 times or until evenly distributed |

- Perform **1.3X** SMRTbell bead cleanup

cDNA amplification + cleanup


In this step, first-strand cDNA products are PCR-amplified and barcoded using barcoded Iso-Seq primers



¹ Twelve barcoded Iso-Seq primers (Iso-Seq primer bc01–12) are available for cDNA amplification step.
² If targeting longer cDNA, 86 µL of SMRTbell cleanup beads (0.86x v/v) can be used.

cDNA amplification + cleanup (cont.)

Procedural notes



3.2 Cleanup of amplified cDNA using 0.9X SMRTbell Cleanup beads


| ✓ | Step | Instructions |
|---|--------|---|
| | 3.2.1 | Add 55 µL of elution buffer to a new strip tube. Transfer 45 µL of PCR-amplified cDNA from Section 3.1 into the strip tube containing elution buffer for a final volume of 100 µL. Add 90 µL (0.9x v/v) of resuspended, room-temperature SMRTbell cleanup beads. The correct ratio of beads to sample is critical at this step. If targeting longer cDNA, 86 µL of SMRTbell cleanup beads (0.86x v/v) can be used. |
| | 3.2.2 | Mix beads by pipetting 10 times or until evenly distributed. |
| | 3.2.14 | Slowly pipette off the cleared supernatant without disturbing the beads. Transfer 24 µL of the supernatant to a new strip tube. Discard the old strip tube with beads. |
| | | Recommended: Measure concentration and size distribution of each cDNA sample. |
| | 3.2.15 | <ul style="list-style-type: none">Take a 1 µL aliquot from each strip tube. Dilute each aliquot with 4 µL of elution buffer.Measure DNA concentration with a Qubit fluorometer using the 1x dsDNA HS kit.Dilute 1:4 dilution further to 1.5 ng/µL based on the Qubit reading if needed.Run 1 µL on an Agilent Bioanalyzer using a High Sensitivity DNA kit. |
| | 3.2.16 | The expected recovery after cDNA amplification SMRTbell clean-up is >80 ng. A minimum of 55 ng of total cDNA is recommended to proceed with Kinnex PCR (Step 4). If less than 55 ng but more than 25 ng is recovered, proceed with Kinnex PCR but expect lower yields. Do not proceed with less than 25 ng. |

3.3 Pooling barcoded cDNA (skip if not multiplexing)

| ✓ | Step | Instructions |
|---|-------|---|
| | 3.3.1 | Using the concentration reading from the Qubit fluorometer, pool an equal mass of each barcoded cDNA sample for a total mass of 55 ng. Store any remaining purified, amplified barcoded cDNA at 4°C for future use. |
| | 3.3.2 | Quick-spin the tube strip in a microcentrifuge to collect liquid. |
| | 3.3.3 | Proceed to next step of the protocol. |

- Perform **0.9X** SMRTbell bead cleanup¹

- Perform DNA concentration QC using Qubit ds DNA HS assay and DNA sizing QC using Bioanalyzer




Example Bioanalyzer DNA sizing QC analysis results for amplified full-length cDNA generated from a universal human RNA reference (UHRR) total RNA sample.

- If multiplexing, perform **equal-mass pooling** of each barcoded cDNA sample after 0.9X SMRTbell bead cleanup


- Expected yield of purified cDNA product is **>80 ng**
- Minimum cDNA amount needed to proceed with Kinnex PCR = **55 ng** (do not proceed with <25 ng)

Kinnex PCR

In this step, incorporate programmable Kinnex segmentation adapter sequences into amplified cDNA products




Set up 8 parallel Kinnex PCR reactions per sample with premixed Kinnex primers to generate amplified cDNA products containing programmable sequences at both ends




Kinnex PCR (cont.)

Procedural summary




Step 4 procedural summary

4.1 Kinnex PCR



4.2 Pooling of 8 Kinnex PCR products and 1.05X SMRTbell cleanup




Proceed to **Step 5**

¹ Kinnex PCR 8-fold kit (103-071-600) contains eight different Kinnex primer mixes (Kinnex primer mix A/B/C/D/E/F/G/HQ).

Kinnex PCR (cont.)

Procedural notes



4.1 Kinnex PCR

✓

Step

Instructions

4.1.1

Thaw primers. The entire volume of primers can be transferred to an 8-tube strip for ease of use with a multi-channel pipette.

| 8X concatenation | Tube color | PN |
|------------------------|------------|-------------|
| 1 Kinnex primer mix A | Orange | 103-107-800 |
| 2 Kinnex primer mix B | | 103-107-900 |
| 3 Kinnex primer mix C | | 103-108-000 |
| 4 Kinnex primer mix D | | 103-108-100 |
| 5 Kinnex primer mix E | | 103-108-200 |
| 6 Kinnex primer mix F | | 103-108-300 |
| 7 Kinnex primer mix G | | 103-108-400 |
| 8 Kinnex primer mix HQ | | 103-108-500 |



4.1.2

Briefly vortex to mix, then quick-spin to collect liquid and place the primer mixes on ice.

4.1.3

Thaw the following components, briefly vortex to mix, then quick-spin to collect liquid and place on ice. Add the components **on ice** in a LoBind tube.

| Master mix components | Volume for 8X concatenation* |
|--|------------------------------|
| PCR-grade water | 88-X µL |
| Kinnex PCR mix | 110 µL |
| 55 ng of amplified cDNA from Step 3.2.16 | X µL |
| Total volume | 198 µL |

X = 55 (ng)/purified pooled DNA concentration from [step 3.2.16](#) (single-plex) or [step 3.3.3](#) (multiplex)
*10% overage included

4.1.4.

Aliquot **22.5 µL** of Master Mix 1 into each of the 8 PCR tubes (for 8X concatenation) **on ice**.

4.1.5

Add **2.5 µL** of Kinnex primer mix into each of 8 PCR tubes from step 4.1.4.

4.1.6

Thoroughly mix by pipetting up and down 10 times. Quick-spin to collect all liquid.

Set up the thermal cycler program listed below with the lid set to 105°C. **Keep sample(s) on ice until the lid is heated to 105°C.**

The duration of PCR is approximately 1 hour.

4.1.7

| Step | Temperature | Duration | Cycles |
|----------------------|-------------|----------|--------|
| Initial denaturation | 98°C | 3 min | 1 |
| Denaturation | 98°C | 20 s | |

• Can transfer entire volume of primers to PCR tubes for ease of use with multi-channel pipettes (8 primer mix tubes)

• Set up Kinnex PCR reactions **ON ICE**


• PCR polymerase 3'→5' exonuclease activity negatively impacts amplification yield if prepared at room temp.

• **Critical step! Correct setup of all 8 Kinnex PCR reactions is required** – any missing/incorrect MAS primer pairs will result in no/low SMRTbell yield

• Set up on ice and add PCR reaction to thermal cycler after lid has preheated to 105°C to avoid digestion of primers by polymerase exonuclease activity

Kinnex PCR (cont.)

Procedural notes



4.2 Pooling of 8 Kinnex PCR products and 1.05X SMRTbell cleanup

| ✓ | Step | Instructions |
|---|--------|---|
| | 4.2.1 | Add exactly 23 µL from each of the 8 PCR reactions into a 1.5 mL tube for a total volume of 184 µL . An equal volume of each PCR product is necessary for efficient array assembly. |
| | 4.2.2 | Add 193 µL (1.05X v/v) of resuspended, room-temperature SMRTbell cleanup beads to a tube of pooled Kinnex PCR amplicon. The correct ratio of beads to pooled sample is critical at this step. |
| | 4.2.3 | Pipette-mix the beads until evenly distributed. |
| | 4.2.4 | Quick-spin the tube in a microcentrifuge to collect liquid. |
| | 4.2.9 | Repeat the previous step. |
| | | Remove residual 80% ethanol: |
| | 4.2.10 | <ul style="list-style-type: none">Remove the tube from the magnetic separation rack.Quick-spin the tube in a microcentrifuge.Place the tube back in the magnetic separation rack until the beads separate fully from the solution.Remove residual 80% ethanol and discard. |
| | 4.2.11 | Remove the tube from the magnetic rack. Immediately add 40 µL of elution buffer to the tube and resuspend the beads by pipetting 10 times or until evenly distributed. |
| | 4.2.12 | Quick-spin the tube in a microcentrifuge to collect liquid. |
| | 4.2.13 | Incubate at room temperature for 5 minutes to elute DNA. |
| | 4.2.14 | Place tube in a magnetic separation rack until beads separate fully from the solution. |
| | 4.2.15 | Slowly aspirate the cleared eluate without disturbing the beads. Transfer eluate to a new tube strip . Discard old tube with beads. |
| | 4.2.16 | Make a 1:10 dilution of the sample in elution buffer and measure DNA concentration with a Qubit fluorometer using the 1x dsDNA HS kit. Typical yield is 6–18 µg. |

SAFE STOPPING POINT -- Store at 4°C or -20°C for long-term storage

- Pool **exactly 23 µL** from each Kinnex PCR reaction in a clean 1.5 mL DNA LoBind tube¹ for a total combined volume of 184 µL

- Add **exactly 193 µL** of SMRTbell cleanup beads (1.05X)
- Kinnex PCR mix significantly increases stringency of SMRTbell clean up beads, **so accurate pipetting is critical**



- Perform DNA concentration QC to verify there is sufficient yield of Kinnex PCR products (min. 4 µg) to proceed to Kinnex array formation step

Kinnex array formation


In this step, assemble cDNA transcripts (“segments”) containing programmable ends into a linear array

- 1 Input RNA quality control
- ↓ ○ Safe stop
- 2 cDNA synthesis + cleanup
- ↓
- 3 cDNA amplification + cleanup
- ↓ ○ Safe stop
- 4 Kinnex PCR + cleanup
- ↓
- 5 Kinnex array formation + cleanup
- ↓ ○ Safe stop
- 6 Nuclease treatment + cleanup


Treat pooled Kinnex PCR products with Kinnex enzyme to create single-stranded overhangs to enable subsequent directional assembly of cDNA transcripts into a linear, ordered array



Kinnex enzyme produces a single nucleotide gap at the location of a uracil (dU) residue within segmentation adapter sequences and enables generation of single-stranded overhangs



Barcoded Kinnex terminal adapters¹ are ligated to specific overhang sequences at array ends



Add Kinnex ligase to enable directional assembly of cDNA transcript segments containing “sticky ends” into a linear, ordered array




Kinnex library prep workflow supports incorporation of PacBio barcodes at the array formation step to enable up to 4-plex sample multiplexing at the SMRTbell library level

¹ Note: Four barcoded terminal Kinnex adapters (Kinnex adapter bcM0001-M0004) are available for Kinnex array formation step. Note: Kinnex concatenation workflow is not compatible with standard SMRTbell adapters from SMRTbell prep kit 3.0 and is also not compatible with SMRTbell barcoded adapter plate 3.0.

Kinnex array formation (cont.)

Kinnex terminal adapters incorporate barcode sequences to enable up to 4-plex sample multiplexing at the library level


- Kinnex adapters contain **barcode sequences¹** to enable (optional) sample multiplexing at the SMRTbell library level (**up to 4-plex**)
 - Forward and reverse Kinnex adapter pairs are pre-mixed in Kinnex concatenation kits
 - Kinnex concatenation kits contain a total of **4 barcoded Kinnex adapter mixes (bcM0001-bcM0004)** to enable multiplexing of up to **4 samples per SMRT Cell**



Kinnex array formation (cont.)

Kinnex terminal adapters use a new design that enables improved HiFi sequencing performance

- Kinnex adapters enable longer polymerase read length performance and improved sample loading efficiency
 - Improved overall HiFi data yield




- Kinnex adapter design requires a different sequencing primer (Kinnex sequencing primer¹) for annealing reaction




Kinnex array formation (cont.)

Procedural summary



Step 5 procedural summary

5.1 Kinnex array formation




5.2 1X SMRTbell cleanup



¹ Kinnex full-length RNA kit bundle (103-072-000) contains four different barcoded Kinnex adapter mixes (bcM0001–bcM0004) → Select one barcode to use per Kinnex library preparation reaction.

Kinnex array formation (cont.)

Procedural notes



5.1 Kinnex array formation

| ✓ | Step | Instructions | | | | | | | | | |
|------------|-------------------|--|------------|-------------|----------|-------|-------------------|--------|---|-----|------|
| | 5.1.0 | Place the DNA damage repair mix on ice at all times and immediately return the DNA damage repair mix back to the freezer (-20°C) after use. Improper storage and handling of the DNA damage repair mix may result in poor library performance and should not be used for subsequent reactions. | | | | | | | | | |
| | 5.1.1 | <p>In a 0.2 mL PCR tube, add 5 µg of sample from Step 4.2.15, in 39 µL of volume (128 ng/µL). Dilute with elution buffer going into this step if the sample is too concentrated.</p> <p>Add 2 µL of Kinnex adapter bc01–04 (select one barcode per library preparation).</p> <p>Note: if combining multiple barcoded Kinnex libraries for sequencing, make sure each library uses one of the 4 different Kinnex barcoded adapters.</p> | | | | | | | | | |
| | | <table><tr><th>Tube color</th><th>Components</th><th>Volume</th></tr><tr><td>Blue</td><td>Kinnex adapter bc</td><td>2 µL</td></tr></table> | Tube color | Components | Volume | Blue | Kinnex adapter bc | 2 µL | | | |
| Tube color | Components | Volume | | | | | | | | | |
| Blue | Kinnex adapter bc | 2 µL | | | | | | | | | |
| | 5.1.2 | <p>Add the following components in the listed order.</p> <p>If processing multiple samples, make a master mix with 10% overage. Pipette mix master mix.</p> | | | | | | | | | |
| | | <p>Add 17 µL of master mix to the PCR tube containing sample and Kinnex adapter. Pipette-mix and run the Kinnex primer digestion/ligation program with the lid set to 55°C.</p> | | | | | | | | | |
| | 5.1.3 | <table><tr><th>Step</th><th>Temperature</th><th>Duration</th></tr><tr><td>1</td><td>45°C</td><td>60 min</td></tr><tr><td>2</td><td>4°C</td><td>Hold</td></tr></table> | Step | Temperature | Duration | 1 | 45°C | 60 min | 2 | 4°C | Hold |
| Step | Temperature | Duration | | | | | | | | | |
| 1 | 45°C | 60 min | | | | | | | | | |
| 2 | 4°C | Hold | | | | | | | | | |
| | 5.1.4 | <p>After running the Kinnex primer digestion/ligation program, add 2 µL of DNA repair mix directly to the Kinnex primer digestion/ligation sample.</p> <table><tr><th>Tube color</th><th>Components</th><th>Volume</th></tr><tr><td>Green</td><td>DNA repair mix</td><td>2 µL</td></tr></table> | Tube color | Components | Volume | Green | DNA repair mix | 2 µL | | | |
| Tube color | Components | Volume | | | | | | | | | |
| Green | DNA repair mix | 2 µL | | | | | | | | | |
| | 5.1.5 | <p>Thoroughly mix by pipetting up and down 10 times and then quick spin to collect all liquid.</p> <p>Run the DNA Damage Repair Program with the lid set to >55°C.</p> | | | | | | | | | |
| | 5.1.6 | <table><tr><th>Step</th><th>Temperature</th><th>Duration</th></tr></table> | Step | Temperature | Duration | | | | | | |
| Step | Temperature | Duration | | | | | | | | | |

• **IMPORTANT!** Place DNA damage repair mix on ice at all times and immediately return the DNA damage repair mix back to the freezer (-20°C) after use

- Improper storage and handling of the DNA damage repair mix may result in poor library performance and should not be used for subsequent reactions

• Recommended input amount to proceed with Kinnex array formation is **5 µg** of Kinnex PCR amplicons (from Step 4)

- Proceeding with <3 µg is **not recommended** since lower input amounts may lead to insufficient final library yields to enable optimal sequencing results

• **IMPORTANT:** If combining multiple barcoded Kinnex libraries for sequencing, make sure each library uses one of the 4 different Kinnex barcoded adapters

• Run Kinnex primer digestion/ligation program with heated lid set to 55°C.



• Perform **DNA Damage Repair** step to repair nicked / damaged DNA sites within newly formed Kinnex array products

• After DNA Damage Repair step, perform **1X SMRTbell bead cleanup** at room temp. (**Step 5.2**)

- Add 1X v/v (60 µL) of resuspended, room temperature SMRTbell cleanup beads to each sample.

Nuclease treatment & cleanup


Perform nuclease treatment and final SMRTbell bead cleanup to remove incomplete SMRTbell template arrays



• After nuclease treatment, most remaining SMRTbell templates are complete (full-length) 8-segment transcript array molecules capped with Kinnex adapters


Nuclease treatment & cleanup (cont.)

Procedural summary




Step 6 procedural summary

5.1 Nuclease treatment



5.2 1X SMRTbell cleanup




Proceed to **ABC¹**

¹ Follow **Step 7** instructions in **Procedure & checklist** to prepare the Kinnex SMRTbell library (25 μ L) for sequencing on Revio (+SPRQ) or Vega system.

Nuclease treatment & cleanup (cont.)

Procedural notes




6.1 Nuclease treatment

| ✓ | Step | Instructions | | | | | | | | | | | | | | | |
|---------------------|-----------------|---|---------------------|-------------|----------|------|-----------|--------|--------------|-----------------|------|-------------|--------------|------|--------------|--|-------|
| | | Add the following components listed below to a new microcentrifuge tube. Adjust component volumes for the number of samples being prepared, plus 10% overage. Pipette mix master mix. For individual preps, add components directly to each sample from the previous step in the order and volume listed below. | | | | | | | | | | | | | | | |
| 6.1.1 | | <table><tr><th colspan="3">Nuclease master mix</th></tr><tr><th>Tube</th><th>Component</th><th>Volume</th></tr><tr><td>Light purple</td><td>Nuclease buffer</td><td>5 µL</td></tr><tr><td>Light green</td><td>Nuclease mix</td><td>5 µL</td></tr><tr><td colspan="2">Total volume</td><td>10 µL</td></tr></table> | Nuclease master mix | | | Tube | Component | Volume | Light purple | Nuclease buffer | 5 µL | Light green | Nuclease mix | 5 µL | Total volume | | 10 µL |
| Nuclease master mix | | | | | | | | | | | | | | | | | |
| Tube | Component | Volume | | | | | | | | | | | | | | | |
| Light purple | Nuclease buffer | 5 µL | | | | | | | | | | | | | | | |
| Light green | Nuclease mix | 5 µL | | | | | | | | | | | | | | | |
| Total volume | | 10 µL | | | | | | | | | | | | | | | |
| 6.1.2 | | Add 10 µL of Nuclease Master mix to each sample. The total volume should be 50 µL. | | | | | | | | | | | | | | | |
| 6.1.3 | | Thoroughly mix by pipetting up and down 10 times and then quick spin to collect all liquid. | | | | | | | | | | | | | | | |
| 6.1.4 | | Run the nuclease treatment program with the lid set to >47°C. <table><tr><th>Step</th><th>Temperature</th><th>Duration</th></tr><tr><td>1</td><td>37°C</td><td>15 min</td></tr><tr><td>2</td><td>4°C</td><td>Hold</td></tr></table> | Step | Temperature | Duration | 1 | 37°C | 15 min | 2 | 4°C | Hold | | | | | | |
| Step | Temperature | Duration | | | | | | | | | | | | | | | |
| 1 | 37°C | 15 min | | | | | | | | | | | | | | | |
| 2 | 4°C | Hold | | | | | | | | | | | | | | | |

- Perform nuclease treatment for 15 min

Nuclease treatment & cleanup (cont.)

Procedural notes



6.2 Final cleanup with 1X SMRTbell cleanup beads

| ✓ | Step | Instructions |
|---|--------|--|
| | 6.2.1 | Add 50 μ L SMRTbell cleanup beads to each sample from the previous step. Pipette-mix the beads until evenly distributed. |
| | 6.2.2 | Quick-spin the tube strip in a microcentrifuge to collect all liquid. |
| | 6.2.3 | Incubate at room temperature for 10 minutes to allow DNA to bind the beads. |
| | 6.2.4 | Place the tube strip in a magnetic separation rack until the beads separate fully from the solution. |
| | 6.2.5 | Slowly remove the cleared supernatant without disturbing the beads and discard the supernatant. |
| | 6.2.6 | Slowly dispense 200 μ L, or enough to cover the beads, of freshly prepared 80% ethanol into each tube. After 30 seconds, remove the 80% ethanol and discard. |
| | 6.2.7 | Repeat the previous step. |
| | | Remove residual 80% ethanol: |
| | 6.2.8 | <ul style="list-style-type: none"> Remove the tube strip from the magnetic separation rack. Quick-spin the tube strip in a microcentrifuge. Place the tube strip back in a magnetic separation rack until beads separate fully from |
| | 6.2.14 | Take a 1 μ L aliquot from each tube. Make a 1:5 dilution of the sample in elution buffer and measure DNA concentration with a Qubit fluorometer using the 1x dsDNA HS kit. Calculate the total mass. Expect 10-25% recovery of the starting Kinnex-PCR product. |
| | | Recommended: Further dilute each aliquot to 250 pg/ μ L with the Femto Pulse dilution buffer. Measure the final SMRTbell library size distribution with a Femto Pulse system. |
| | 6.2.15 | If required, normalize 25 μ L of library to <60 ng/ μ L. Using a concentration above 60 ng/ μ L into ABC will result in poor sequencing performance. |
| | 6.2.16 | Proceed to Section 7 to prepare library for sequencing with Vega or Revio (+SPRQ chemistry). or For sequencing on Revio (non-SPRQ) or Sequel II/e, proceed to SMRT Link Sample Setup for preparing samples. |
| | 6.2.17 | Store SMRTbell libraries at 4°C if sequencing within two weeks. Store long-term at -20°C. |

- Perform **DNA concentration QC** on final purified Kinnex RNA library using a Qubit dsDNA HS assay
 - Typical final SMRTbell library yield from 5 μ g of input DNA into Kinnex array formation is **~10 – 25%** – a much higher observed yield might suggest incomplete digestion of partial SMRTbell templates
 - Troubleshooting tip:** If SMRTbell library yield is higher than expected and sample loading is lower than expected, consider repeating nuclease treatment step




- Perform **DNA sizing QC** on final purified Kinnex full-length RNA library using a Femto Pulse system (expected final library insert size is **~12 – 16 kb**)

- Kinnex full-length RNA final SMRTbell library concentration must be **≤60 ng/ μ L** to proceed with sample setup (ABC¹) [Step 7]

→ Using a concentration above 60 ng/ μ L will result in lower loading during sequencing

Kinnex full-length RNA library prep inputs & expected step yields

Final Kinnex library yield is typically sufficient to load ≥ 2 SMRT Cells




Kinnex full-length RNA sequencing preparation workflow details

Sample Setup & Run Design recommendations for Kinnex libraries

Follow SMRT Link run design instructions to prepare Kinnex libraries for sequencing

Select desired Kinnex application type in SMRT Link run design page



Runs

After specifying your application type, SMRT Link auto-fills selected Run Design parameter fields with default recommended values

Kinnex single-cell RNA¹

Whole genome sequencing

- Human WGS
- Microbial assembly
- Other WGS

RNA sequencing

- Iso-Seq method
- MAS-Seq single cell
- Kinnex single-cell RNA**
- Kinnex full-length RNA

Viral sequencing

- Adeno-associated virus

Metagenomics

- Kinnex 16S rRNA
- Full-length 16S rRNA sequencing

Kinnex full-length RNA²

Whole genome sequencing

- Human WGS
- Microbial assembly
- Other WGS

RNA sequencing

- Iso-Seq method
- MAS-Seq single cell
- Kinnex single-cell RNA
- Kinnex full-length RNA**

Viral sequencing

- Adeno-associated virus

Metagenomics

- Kinnex 16S rRNA
- Full-length 16S rRNA sequencing

Kinnex 16S rRNA²

Whole genome sequencing

- Human WGS
- Microbial assembly
- Other WGS

RNA sequencing

- Iso-Seq method
- MAS-Seq single cell
- Kinnex single-cell RNA
- Kinnex full-length RNA

Viral sequencing

- Adeno-associated virus

Metagenomics

- Kinnex 16S rRNA**
- Full-length 16S rRNA sequencing

| SMRT Link module | Key setup parameters | Sequel II system recommended settings | Vega system recommended settings | Revio system recommended settings |
|-------------------|---------------------------|---------------------------------------|----------------------------------|-----------------------------------|
| Sample setup | Library type | Kinnex | | |
| | Primer | Kinnex sequencing primer | | |
| | Binding/Polymerase kit | Sequel II binding kit 3.2 | Vega polymerase kit | Revio SPRQ polymerase kit |
| | Concentration on plate | 40 – 60 pM | 130 – 160 pM | 130 – 160 pM |
| Runs → Run design | Use Adaptive Loading | YES | | |
| | Library Type ¹ | Kinnex | | |
| | Movie Acquisition Time | 30 hrs | 24 hrs | 24 hrs |
| | Add Analysis | YES | | |

SMRT Link Run Design procedure for Revio and Vega systems

Sample and run information

Kinnex single-cell RNA

▼ Plate 1, Well A01: Kinnex single-cell RNA library demo

| | |
|---|------------------------|
| Application <small>Required</small> | Kinnex single-cell RNA |
| Plate Well <small>Required</small> | Plate 1, Well A01 |
| Well Name <small>Required</small> | |
| Well Comment | |
| Library Type <small>Required</small> | Kinnex |
| Insert Size (bp) <small>Required</small> | 15000 |
| Library Concentration (pM) <small>Required</small> | |
| Movie Acquisition Time (hours) | 24 |

Kinnex full-length RNA

▼ Plate 1, Well A01: Kinnex full-length RNA library demo

| | |
|---|------------------------|
| Application <small>Required</small> | Kinnex full-length RNA |
| Plate Well <small>Required</small> | Plate 1, Well A01 |
| Well Name <small>Required</small> | |
| Well Comment | |
| Library Type <small>Required</small> | Kinnex |
| Insert Size (bp) <small>Required</small> | 16000 |
| Library Concentration (pM) <small>Required</small> | |
| Movie Acquisition Time (hours) | 24 |

Kinnex 16S rRNA

▼ Plate 1, Well A01: Kinnex 16S rRNA library demo

| | |
|---|-------------------|
| Application <small>Required</small> | Kinnex 16S rRNA |
| Plate Well <small>Required</small> | Plate 1, Well A01 |
| Well Name <small>Required</small> | |
| Well Comment | |
| Library Type <small>Required</small> | Kinnex |
| Insert Size (bp) <small>Required</small> | 17000 |
| Library Concentration (pM) <small>Required</small> | |
| Movie Acquisition Time (hours) | 24 |

Select application type to autofill Library Type & Movie Acquisition Time recommended settings

Default = Kinnex library type (instead of Standard or AAV)¹



Forward and reverse standard terminal adapters have the same structure



Library Type field determines which adapter finding algorithm is used during post-primary analysis¹

Forward and reverse Kinnex terminal adapters have different structures

¹ **Note:** When sequencing a Kinnex library sample, if 'Standard' library type is mistakenly selected instead of 'Kinnex' then a higher missing adapter rate (> 95%) and a slight degradation in barcode demultiplexing performance (~93-96% barcoded HiFi read yield) will be observed.

SMRT Link Run Design procedure for Revio and Vega systems (cont.)

Sample indexing (barcoding) information

Kinnex single-cell RNA

Samples

Sample is indexed ☒ YES ☐ NO

Indexes Required

MAS adapter indexes

Biosample names Required

Interactively

From a File

Default Indexes FASTA = MAS adapter indexes

Kinnex full-length RNA

Samples

Sample is indexed ☒ YES ☐ NO

Indexes Required

MAS adapter indexes

Biosample names Required

Interactively

From a File

Kinnex 16S rRNA

Samples

Sample is indexed ☒ YES ☐ NO

Indexes Required


MAS adapter indexes

Biosample names Required

Interactively

From a File

Example complete Kinnex library molecule containing barcoded Kinnex terminal MAS adapters¹ at both ends



Example interactive biosample name specification

Barcode Selector and Sample Name Editor

Available Barcodes

Filter...

☐

Barcode ID

☐

bcM0001--bcM0001

☐

bcM0002--bcM0002

☐

bcM0003--bcM0003

☐

bcM0004--bcM0004

Included Barcodes

Filter...

☐

Barcode ID

Bio Sample ID

Barcode Selector and Sample Name Editor

Available Barcodes

Filter...

☐

Barcode ID

☐

bcM0003--bcM0003

☐

bcM0004--bcM0004

Included Barcodes

Filter...

☐

Barcode ID

Bio Sample ID

☐

bcM0001--bcM0001

Kinnex adapter-barcoded library 1

☐

bcM0002--bcM0002

Kinnex adapter-barcoded library 2

SMRT Link



MAS SMRTbell barcoded adapter indexes

```
>bcM0001
ACAGTC
>bcM0002
ATGACG
>bcM0003
CACGTG
>bcM0004
CATCGC
```

¹ Four barcoded terminal Kinnex adapters (Kinnex adapter bcM0001-bcM0004) are available for Kinnex array formation step. Kinnex adapter barcode sequences can be downloaded from [SMRT Link](#) Data Management module.

SMRT Link Run Design procedure for Revio and Vega systems (cont.)

Data options

| Kinnex single-cell RNA | Kinnex full-length RNA | Kinnex 16S rRNA |
|--|--|--|
| <div><div>▼ Data Options</div><div><div>Include Base Kinetics</div><div><input type="radio"/> YES <input checked="" type="radio"/> NO</div></div><div><div>Consensus Mode</div><div><input checked="" type="radio"/> MOLECULE <input type="radio"/> STRAND</div></div><div><div>Assign Data To Project ?</div><div>General Project</div></div></div> | <div><div>▼ Data Options</div><div><div>Include Base Kinetics</div><div><input type="radio"/> YES <input checked="" type="radio"/> NO</div></div><div><div>Consensus Mode</div><div><input checked="" type="radio"/> MOLECULE <input type="radio"/> STRAND</div></div><div><div>Assign Data To Project ?</div><div>General Project</div></div></div> | <div><div>▼ Data Options</div><div><div>Include Base Kinetics</div><div><input type="radio"/> YES <input checked="" type="radio"/> NO</div></div><div><div>Consensus Mode</div><div><input checked="" type="radio"/> MOLECULE <input type="radio"/> STRAND</div></div><div><div>Assign Data To Project ?</div><div>General Project</div></div></div> |

Default = NO for Include Base Kinetics

Default Consensus Mode = MOLECULE¹

Can leave Include Base Kinetics and Consensus Mode fields at their default settings for Kinnex library samples

SMRT Link Run Design analysis options

Kinnex single-cell RNA

Analysis Options

Default = YES for Add Analysis

Add Analysis ☒ YES ☐ NO

Analysis Name Required Kinnex_Single-Cell_RNA_Demo_Analysis_Job_Name

Select Analysis Workflow Required Read Segmentation and Single-Cell Iso-Seq

Segmentation Adapter Set MAS-Seq Adapter v1 (MAS16)

Primer Set Required 10x Chromium single cell 3' cDNA primers

Reference Set Required Human Genome hg38, with Gencode v39 annotations

Kit Type ☒ 10X 3' KIT ☐ 10X 5' KIT

Advanced Parameters

Kinnex full-length RNA

Analysis Options

Add Analysis ☒ YES ☐ NO

Analysis Name Required Kinnex_Full-Length_RNA_Demo_Analysis_Job_Name

Select Analysis Workflow Required Read Segmentation and Iso-Seq

Segmentation Adapter Set MAS-Seq Adapter v3 (MAS8)

Primer Set Required Iso-Seq v2 Barcoded cDNA Primers

Reference Set Human Genome hg38, with Gencode v39 annotations

Cluster of Barcoded Samples Pool reads and cluster together

Advanced Parameters

Kinnex 16S rRNA

Analysis Options

Add Analysis ☒ YES ☐ NO

Analysis Name Required Kinnex_Full-Length_16S_rRNA_Demo_Analysis_Job_Name


Select Analysis Workflow Required Read Segmentation

Segmentation Adapter Set MAS-Seq Adapter v2 (MAS12)


Advanced Parameters

Split arrayed HiFi reads at adapter positions, generating segmented reads (S-reads)


Read Segmentation and Single-Cell Iso-Seq



Read Segmentation and Iso-Seq



Read Segmentation



SMRT Link Run Design analysis options (cont.)

Kinnex single-cell RNA

▼ Analysis Options

Add Analysis ☒ YES ☐ NO

Analysis Name Required

Kinnex_Single-Cell_RNA_Demo_Analysis_Job_Name

Select Analysis Workflow Required

Read Segmentation and Single-Cell Iso-Seq

Segmentation Adapter Set

MAS-Seq Adapter v1 (MAS16)

Primer Set Required

10x Chromium single cell 3' cDNA primers

Reference Set Required

Human Genome hg38, with Gencode v39 annotations

Kit Type ⓘ

☒ 10X 3' KIT ☐ 10X 5' KIT

Advanced Parameters

Kinnex full-length RNA

▼ Analysis Options

Add Analysis ☒ YES ☐ NO

Analysis Name Required

Kinnex_Full-Length_RNA_Demo_Analysis_Job_Name

Select Analysis Workflow Required

Read Segmentation and Iso-Seq

Segmentation Adapter Set

MAS-Seq Adapter v3 (MAS8)

Primer Set Required

Iso-Seq v2 Barcoded cDNA Primers

Reference Set

Human Genome hg38, with Gencode v39 annotations

Cluster of Barcoded Samples ⓘ

Pool reads and cluster together

Advanced Parameters

Kinnex 16S rRNA

▼ Analysis Options

Add Analysis ☒ YES ☐ NO

Analysis Name Required

Kinnex_Full-Length_16S_rRNA_Demo_Analysis_Job_Name

Select Analysis Workflow Required


Read Segmentation

Segmentation Adapter Set

MAS-Seq Adapter v2 (MAS12)

Advanced Parameters

Specify Segmentation Adapter Set that corresponds to the Kinnex library concatenation method used



SMRT Link Run Design analysis options (cont.)

Kinnex single-cell RNA

Analysis Options

Add Analysis

☒ YES ☐ NO

Analysis Name

Required

Kinnex_Single-Cell_RNA_Demo_Analysis_Job_Name

Select Analysis Workflow

Required

Read Segmentation and Single-Cell Iso-Seq

Segmentation Adapter Set

MAS-Seq Adapter v1 (MAS16)

Primer Set

Required

10x Chromium single cell 3' cDNA primers

Reference Set

Required

Human Genome hg38, with Gencode v39 annotations

Kit Type

☒ 10X 3' KIT ☐ 10X 5' KIT

Advanced Parameters

Specify Primer Set used for cDNA amplification

Kinnex full-length RNA

Analysis Options

Add Analysis

☒ YES ☐ NO

Analysis Name

Required

Kinnex_Full-Length_RNA_Demo_Analysis_Job_Name

Analysis Workflow

Required

Read Segmentation and Iso-Seq

Segmentation Adapter Set

MAS-Seq Adapter v3 (MAS8)

Primer Set

Required

Iso-Seq v2 Barcoded cDNA Primers

Reference Set

Human Genome hg38, with Gencode v39 annotations

Cluster of Barcoded Samples

☒ Pool reads and cluster together

Advanced Parameters

Specify Reference Set. Default sets are:

- Human Genome hg38, with Gencode v39 annotations
- Mouse Genome mm39, with Gencode vM28 annotations

Specify reference genome & annotation sets to align high quality isoforms to, and to collapse isoforms mapped to the same genomic loci.

10x Forward (F) PCR primer

10x Reverse (R) PCR primer

5'

[CBC][UMI]TTTTTTTT

LIBRARY INSERT

CCC

TSO

3'

3'

[CBC][UMI]AAAAAAAA

LIBRARY INSERT

GGG

TSO

5'

For Kinnex single-cell 3' RNA analysis, select '10x Chromium single cell 3' cDNA primers'

OR

5'

[CBC][UMI][TSO]GGG

LIBRARY INSERT

AAAAAAAA

3'

3'

[CBC][UMI][TSO]CCC

LIBRARY INSERT

TTTTTTTT

5'

For Kinnex single-cell 5' RNA analysis, select '10x Chromium single cell 5' cDNA primers'

Forward (F) PCR primer (Barcoded Iso-Seq primer bcXX²)

Reverse (R) PCR primer (Iso-Seq cDNA amplification primer)

5'

BC

TSO

GGG

LIBRARY INSERT

AAAAAAAA

3'

3'

BC

TSO

CCC


LIBRARY INSERT

TTTTTTTT

5'

For Kinnex full-length RNA analysis, select 'Iso-Seq v2 Barcoded cDNA primers'

Specify primer sequence file in FASTA format to identify cDNA primers for removal (include the 5' and 3' cDNA primers)



SMRT Link Run Design analysis options (cont.)

Kinnex single-cell RNA

▼

Analysis Options

Add Analysis

☒ YES ☐ NO

Analysis Name

Required

Kinnex_Single-Cell_RNA_Demo_Analysis_Job_Name

Select Analysis Workflow

Required

Read Segmentation and Single-Cell Iso-Seq

Segmentation Adapter Set

MAS-Seq Adapter v1 (MAS16)

Primer Set

Required

10x Chromium single cell 3' cDNA

Reference Set

Required

Human Genome hg38, with Gencode v39 annotations

Kit Type

☒ 10X 3' KIT ☐ 10X 5' KIT

Advanced Parameters

Specify Kit Type used for single-cell cDNA generation

Specification of Kit Type (10x 3' Kit or 10x 5' Kit) determines which set of 10x barcode sequences to use, and also affects UMI and single-cell barcode design settings¹

5' [CBC][UMI] TTTTTTTTTT LIBRARY INSERT CCC TSO 3'
3' [CBC][UMI] AAAAAAAAAA LIBRARY INSERT GGG TSO 5'

For Kinnex single-cell 3' RNA analysis, select '10x 3' Kit'

OR

5' [CBC][UMI][TSO]GGG LIBRARY INSERT AAAAAAAAAA 3'
3' [CBC][UMI][TSO]CCC LIBRARY INSERT TTTTTTTTTT 5'

For Kinnex single-cell 5' RNA analysis, select '10x 5' Kit'

Kinnex full-length RNA

▼

Analysis Options

Add Analysis

☒ YES ☐ NO

Analysis Name

Required

Select Analysis Workflow

Required

Segmentation Adapter Set

MAS-Seq Adapter v1 (MAS8)

Primer Set

Required

V2 Barcoded cDNA Primers

Reference Set

Required

Human Genome hg38, with Gencode v39 annotations

Cluster of Barcoded Samples

☒ Pool reads and cluster together

Advanced Parameters

Specify how to perform read clustering for barcoded samples

Pool reads and cluster together

▼


-- select --

Cluster reads separately

Pool reads and cluster together

Specification of Cluster of Barcoded Samples setting determines whether all FLNC reads will be pooled for clustering (Does not apply to non-barcoded samples.)

Clustering options for Kinnex full-length RNA samples



SMRT Link Read Segmentation and Iso-Seq analysis video demonstration

Video demonstration of SMRT Link Read Segmentation and Iso-Seq application workflow for analysis of Kinnex full-length RNA samples

Analysis Application Required

Read Segmentation and Iso-Seq

Import Analysis Settings Export

Associated Inputs

Segmentation Adapter Set

MAS-Seq Adapter v3 (MAS8)

Primer Set Required

Iso-Seq v2 Barcoded cDNA Primers

Reference Set

Human Genome hg38, with Gencode v39 annotations

Cluster of Barcoded Samples

Pool reads and cluster together

Advanced Parameters

Analysis Name

test


Analysis Datasets

Displaying rows 1 to 1 out of 1

| ID | Name |
|-------|------------------------|
| 21... | 3230211_KPoS_64007_... |

[Demo video](#) for Read Segmentation and Iso-Seq workflow (SMRT Link v13.0+)

- Workflow supports full-length isoform analysis for data generated on PacBio long-read systems using Kinnex full-length RNA kit
- End-to-end workflow begins with HiFi reads and outputs full-length isoform classifications with supporting read count information






Kinnex full-length RNA example sequencing performance data

Example Kinnex full-length RNA library preparation QC results


Kinnex full-length RNA library prepared with human UHRR total RNA

Amplified full-length cDNA QC



Example Bioanalyzer DNA sizing QC analysis results for amplified full-length cDNA generated from a universal human RNA reference (UHRR) total RNA sample.

Final Kinnex full-length RNA library QC



Example Femto Pulse DNA sizing QC analysis results for final Kinnex full-length RNA library.

Final Kinnex library yield is typically sufficient to load the following number of SMRT Cells:

- >8 SMRT Cells for Revio using SPRQ chemistry
- >2 SMRT Cells for the Vega or Revio (non-SPRQ) systems
- >4 SMRT Cells for Sequel II/Ile systems


| | |
|--|-----------------|
| Total RNA input for cDNA synthesis | 300 ng |
| cDNA input for Kinnex array formation | 5900 ng |
| Post-nuclease treatment & final library cleanup yield (%) ¹ | 1460 ng (24.7%) |

¹ Post-nuclease treatment & final cleanup yields typically ranged from ~10% to ~25% when using UHRR total RNA samples for Kinnex full-length RNA library construction.

Example sequencing performance for Kinnex full-length RNA libraries prepared with human cDNA [Revio system + SPRQ chemistry¹]

Kinnex full-length RNA library for universal human reference RNA (UHRR) sample


Raw Data Report



| | |
|-----------------------------|----------|
| Raw Base Yield | 1,442 Gb |
| Mean Polymerase Read Length | 71.9 kb |
| P0 | 20% |
| P1 | 80% |
| P2 | 1% |

Example sequencing metrics for a Universal Human Reference RNA (UHRR) Kinnex full-length RNA library sample run on a Revio system with Revio SPRQ polymerase kit / 160 pM on-plate loading concentration (OPLC) / 24-hrs movie time.


HiFi Read Length



| | |
|----------------------------|----------|
| HiFi Reads | 10.4 M |
| HiFi Base Yield | 146.6 Gb |
| Mean HiFi Read Length | 14.0 kb |
| Median HiFi Read Quality | Q32 |
| HiFi Read Mean # of Passes | 7 |

For UHRR Kinnex full-length RNA libraries, per-Revio SMRT Cell HiFi read counts were typically >6 Million depending on the final library insert size and P1 loading performance.

Read Segmentation Metrics



| | |
|--|------------|
| Input HiFi Reads | 10,404,037 |
| Segmented reads (S-reads) | 81,012,128 |
| Mean length of S-reads | 1,758 bp |
| Percent of reads with full arrays | 94.55% |
| Mean array size (concentration factor) | 7.79 |


For UHRR Kinnex libraries, per-Revio SMRT Cell segmentation read counts were typically >50 Million.

¹ HiFi read lengths, reads/data per SMRT Cell and other sequencing performance results can vary depending on DNA sample quality, insert size, sample loading performance & movie time. Note: Shorter library insert sizes (<15 kb), lower DNA quality samples, and suboptimal sample loading performance may result in reduced HiFi data yields per SMRT Cell.

Example sequencing performance for Kinnex full-length RNA libraries prepared with human cDNA [Vega system¹]

Kinnex full-length RNA library for universal human reference RNA (UHRR) sample


Raw Data Report



| | |
|-----------------------------|---------|
| Mean Polymerase Read Length | 95.3 kb |
| Loading level | 61% |

Example sequencing metrics for a Universal Human Reference RNA (UHRR) Kinnex full-length RNA library sample run on a Vega system with Vega polymerase kit / 130 pM on-plate loading concentration (OPLC) / 24-hrs movie time.


HiFi Read Length



| | |
|----------------------------|---------|
| HiFi Reads | 4.6 M |
| HiFi Base Yield | 65.0 Gb |
| Mean HiFi Read Length | 14.0 kb |
| Median HiFi Read Quality | Q35 |
| HiFi Read Mean # of Passes | 11 |

For UHRR Kinnex full-length RNA libraries, per-Vega SMRT Cell HiFi read counts were typically >2.5 Million depending on the final library insert size and sample loading performance.

Read Segmentation Metrics



| | |
|--|------------|
| Input HiFi Reads | 4,642,616 |
| Segmented reads (S-reads) | 35,917,435 |
| Mean length of S-reads | 1,759 bp |
| Percent of reads with full arrays | 94.14% |
| Mean array size (concentration factor) | 7.74 |

For UHRR Kinnex libraries, per-Vega SMRT Cell segmentation read counts were typically >20 Million.

¹ HiFi read lengths, reads/data per SMRT Cell and other sequencing performance results can vary depending on DNA sample quality, insert size, sample loading performance & movie time. Note: Shorter library insert sizes (<15 kb), lower DNA quality samples, and suboptimal sample loading performance may result in reduced HiFi data yields per SMRT Cell.




Kinnex full-length RNA data analysis workflow overview

Kinnex full-length RNA bioinformatics workflow overview

SMRT Link Read segmentation and Iso-Seq workflow processes HiFi reads generated from Kinnex full-length RNA libraries to produce classified isoforms with read counts that are compatible with tertiary analysis tools

SMRT Link read segmentation and Iso-Seq workflow¹



Read segmentation

- HiFi reads are segmented into individual segmented reads (**S-reads**) that represent the original cDNA sequences

Primer extraction

- Primers and polyA tails are removed, but also used to orient the read into 5' → 3' orientation

Isoform clustering

- FLNC reads are clustered by their sequencing similarity to produce isoform consensus sequences
- This step is the last step of Iso-Seq analysis if no genome is provided

Mapping

- If a genome is provided, isoform consensus sequences from the previous step are mapped and further collapsed by their exonic structures to produce isoforms as GFF files for visualization

Transcript classification²

- If an annotation (e.g., Gencode) is provided, isoforms are classified against it using pigeon (the PacBio implementation of SQANTI3) to identify known and novel genes/isoforms
- The Iso-Seq workflow can jointly analyze pooled sample reads to produce a unified isoform annotation with per-sample read counts, both raw and normalized as counts per million (CPM)

SMRT Link Read Segmentation and Iso-Seq analysis video demonstration

Video demonstration of SMRT Link Read Segmentation and Iso-Seq application workflow for analysis of Kinnex full-length RNA samples

Analysis Application Required

Read Segmentation and Iso-Seq

Import Analysis Settings Export

Associated Inputs

Segmentation Adapter Set

MAS-Seq Adapter v3 (MAS8)

Primer Set Required

Iso-Seq v2 Barcoded cDNA Primers

Reference Set

Human Genome hg38, with Gencode v39 annotations

Cluster of Barcoded Samples

Pool reads and cluster together

Advanced Parameters

Analysis Name

test


Analysis Datasets

Displaying rows 1 to 1 out of 1

| ID | Name |
|-------|------------------------|
| 21... | 3230211_KPoS_64007_... |

[Demo video](#) for Read Segmentation and Iso-Seq workflow (SMRT Link v13.0+)

- Workflow supports full-length isoform analysis for data generated on PacBio Sequel II/Ile and Revio systems using Kinnex full-length RNA kit
- End-to-end workflow begins with HiFi reads and outputs full-length isoform classifications with supporting read count information



Kinnex full-length RNA bioinformatics workflow recommendations

SMRT Link Read Segmentation and Iso-Seq workflow common considerations and recommendations for analysis of Kinnex full-length RNA data¹

Analysis recommendations for Iso-Seq data based on reference genome and annotation availability¹

- With SMRT Link v13, the **Read segmentation and Iso-Seq workflow** analysis application supports **human and mouse** reference genomes and annotations to produce classified isoforms with read counts.
- If working with **other organisms**, see table below for analysis recommendations

| Available reference or annotation | Analysis workflow recommendation |
|-------------------------------------|--|
| Human or mouse | <ul style="list-style-type: none">• Use the Iso-Seq workflow with preloaded human / mouse annotation to get mapped, unique isoforms with classifications and read count information (FASTA, GFF, TXT). |
| Model organism with good annotation | <ul style="list-style-type: none">• Run Iso-Seq workflow with uploaded reference genome to get mapped, unique isoforms (FASTA, GFF)• Generate pigeon-compliant annotation and use the command line for isoform classification with read count information (TXT) |
| Non-model organism with genome | <ul style="list-style-type: none">• Run Iso-Seq workflow with uploaded reference genome to get mapped, unique isoforms (FASTA, GFF) |
| No genome | <ul style="list-style-type: none">• Run Iso-Seq workflow without reference genome to get unique isoforms (FASTA) |

Kinnex full-length RNA bioinformatics workflow recommendations (cont.)

SMRT Link Read Segmentation and Iso-Seq workflow common considerations and recommendations for analysis of Kinnex full-length RNA data¹

Sequencing depth recommendations for Iso-Seq data based on experimental goals and study design

| Example application | Human genetics disease studies | Biopharma for identifying highly expressed targets | Plant & animal whole genome annotation |
|-------------------------------------|---|--|--|
| Experimental goal | Isoform discovery and quantification of moderate-to-rare transcripts | Isoform discovery of high expressed transcripts | Comprehensive transcript annotation in a species |
| Example study design | Disease vs. normal tissues with multiple replicates | Disease cohort with >20+ samples | Plant or animal with multiple tissue types |
| Target depth of coverage per sample | 10 M reads per sample | 5 M reads per sample | ≤5 M reads per tissue (of same species) |
| Sample multiplexing ¹ | Sequel II/Ile system: Up to 2 samples per SMRT Cell 8M (2-plex) | Sequel II/Ile system: Up to 3 samples per SMRT Cell 8M (3-plex) | |
| | Vega system: Up to 3 samples per Vega SMRT Cell (3-plex) | Vega system: Up to 6 samples per Vega SMRT Cell (6-plex) | |
| | Revio system + SPRQ: Up to 6 samples per Revio SMRT Cell (6-plex) | Revio system + SPRQ: Up to 12 samples per Revio SMRT Cell (12-plex) | |
| SMRT Link data analysis workflows | Read Segmentation and Iso-Seq analysis application with option to “pool reads and cluster together” to get a master isoform classification file with per-sample full-length read counts | | |

SMRT Link Read Segmentation and Iso-Seq analysis application setup

Specify **Read Segmentation and Iso-Seq** analysis application in SMRT Link¹

Read Segmentation and Iso-Seq analysis application processes HiFi reads generated from a Kinnex full-length RNA library to produce classified isoforms with read counts that are compatible with tertiary analysis tools

- Accepts **HiFi reads** (BAM format) as input
 - HiFi reads are reads generated with CCS analysis whose quality value is equal to or greater than 20
- HiFi reads should be generated using the Kinnex full-length RNA library preparation protocol ([103-238-700](#))
 - If the library is a regular (non-Kinnex) Iso-Seq monomer library without MAS-Seq concatenation, use the SMRT Link Iso-Seq Analysis workflow instead

Note on barcoded libraries

The **Read Segmentation and Iso-Seq** workflow will only process barcoded libraries at the cDNA level (such as using **Iso-Seq v2 Barcoded cDNA Primers** as part of the MAS-Seq for bulk Iso-Seq kit)

→ Demultiplexing of **barcoded adapters** (also part of the Kinnex full-length RNA kit) should first be performed by running the **Demultiplex Barcodes** data utility workflow in SMRT Link.

| ID | Name |
|-------|------------------------------------|
| 38820 | Kinnex_full_length_RNA_UHRR_sample |

SMRT Link Read Segmentation and Iso-Seq analysis application setup (cont.)

Specify **Read Segmentation** and **Single-Cell Iso-Seq** analysis application required associated inputs¹

PacBio SMRT Analysis

SMRT Analysis / Create New Analysis

1. Select Data 2. Select Analysis

Analysis Application Required

Read Segmentation and Iso-Seq

Import Analysis Settings Export

Associated Inputs

1 Segmentation Adapter Set

MAS-Seq Adapter v3 (MAS8)

2 Primer Set Required

Iso-Seq v2 Barcoded cDNA Primers

3 Reference Set

Human Genome hg38, with Gencode v39 annotatic

Cluster of Barcoded Samples

Cluster reads separately

Advanced Parameters

1. Segmentation Adapter Set (Default = MAS-Seq Adapter v3 (MAS8))

- Specify a FASTA file, provided by PacBio, containing segmentation adapters. If you need a custom segmentation adapter set, click Advanced Parameters and use a custom FASTA file formatted as described in the SMRT Link User Guide [documentation](#).

2. Primer Set (Required) (Default = Iso-Seq v2 Barcoded cDNA Primers)

- Specify a primer sequence file in FASTA format to identify cDNA primers for removal. The primer sequence includes the 5' and 3' cDNA primers
- Primer IDs must be specified using the suffix `_5p` to indicate 5' cDNA primers and the suffix `_3p` to indicate 3' cDNA primers. The 3' cDNA primer should not include the `Ts` and is written in reverse complement. (See the [SMRT Link User Guide](#) for example Iso-Seq v2 Barcoded cDNA Primer IDs and sequences)
- Each primer sequence must be unique

3. Reference Set (Required)

- Specify one of two default reference genome and annotation sets to align high quality isoforms to, and to collapse isoforms mapped to the same genomic loci. The default sets are `Human_hg38_Gencode_v39` and `Mouse_mm39_Gencode_vM28` annotations
- Alternatively, choose other reference genomes (but not with annotations) that were custom-uploaded to SMRT Link
- The Reference Set can be left blank. If blank, the workflow will stop after the isoform clustering step (`isoseq cluster`)

SMRT Link Read Segmentation and Iso-Seq analysis application setup (cont.)

Specify **Read Segmentation** and **Single-Cell Iso-Seq** analysis application required associated inputs¹

PacBio SMRT Analysis ▾

SMRT Analysis / Create New Analysis

1. Select Data

2. Select Analysis

Analysis Application Required

Read Segmentation and Iso-Seq ▾

⬆ Import Analysis Settings

⬇ Export ⓘ

Associated Inputs

Segmentation Adapter Set

MAS-Seq Adapter v3 (MAS8) ⌵

Primer Set Required

Iso-Seq v2 Barcoded cDNA Primers ⌵

Reference Set

Human Genome hg38, with Gencode v39 annotatic ⌵


Cluster of Barcoded Samples ⓘ

Cluster reads separately ▾

Advanced Parameters

4. Cluster of Barcoded Samples

- This option specifies barcoded samples that were barcoded at the cDNA level, where the (barcoded) cDNA primers are specified in the Primer Set option. This option does **not** address libraries that were barcoded using barcoded adapters
- Specify whether all FLNC reads will be pooled for clustering, then demultiplexed based on pooled result. **Note:** This setting does **not** apply to non-barcoded samples
- Specify **Pool reads and cluster together** if barcoded samples are from the **same** species, but different tissues, or samples of the same genes but different individuals. The samples are clustered with **all** barcodes pooled
- Specify **Cluster reads separately** if barcoded samples are from different species. The samples are clustered **separately** by barcode
- In either case, the samples on the results page are automatically named BioSample_1 through BioSample_N



Example SMRT Link Read Segmentation data utility processing results¹ for Kinnex full-length RNA library prepared with human UHRR sample

SMRT Link Read Segmentation data utility job report – Summary Metrics and Segmentation Statistics

Summary Metrics


| Value | Analysis Metric |
|------------|--|
| 10,404,037 | Reads |
| 81,012,128 | Segmented reads (S-reads) |
| 1,758.004 | Mean length of S-reads |
| 94.55 % | Percent of reads with full arrays |
| 7.79 | Mean array size (concatenation factor) |

Example Revio system + SPRQ data shown.

For Revio system with optimal sample *P1* loading:

- Yield of HiFi reads is ideally >6 Million
- Yield of segmented reads is ideally >50 Million
- Mean S-read length will depend on size distribution of cDNA sample (ideally >1,500 bp)
- % of reads with full arrays is ideally >80%
- Mean array size is ideally >7.0 segments

Segmentation Statistics



- Plot should ideally show a majority (>80%) of reads with array lengths = 8 segments

- **Reads:** Number of input arrayed HiFi reads
- **Segmented reads (S-reads):** Number of generated S-reads
- **Mean length of S-reads:** Mean read length of generated S-reads
- **Percent of reads with full arrays:** Percentage of input HiFi reads containing all adapter sequences in the order listed in the segmentation adapter FASTA file
- **Mean array size:** Mean number of fragments (or S-reads) found in input reads


Histogram distribution of the number of S-reads per HiFi read. (Example Revio system + SPRQ data shown.)

¹ HiFi read lengths, reads/data per SMRT Cell and other sequencing performance results can vary depending on DNA sample quality, insert size, *P1* loading performance & movie time. Note: Refer to *SMRT Link v13.1 Kinnex single-cell troubleshooting guide* (103-516-100) for example performance metrics typically achievable with Kinnex libraries under optimal *P1* loading conditions. For Revio system, we recommend aiming for ~60 – 80% *P1* loading for Kinnex full-length RNA libraries.

Example SMRT Link Read Segmentation data utility processing results¹ for Kinnex full-length RNA library prepared with human UHRR sample (cont.)


SMRT Link Read Segmentation data utility job report – Length of Reads and S-read Length Distribution

Length of Reads



Histogram distribution of the number of HiFi reads by read length, in base pairs. (Example Revio system + SPRQ data shown.)

S-read Length Distribution



Histogram distribution of the number of S-reads by HiFi read length, in base pairs. (Example Revio system + SPRQ data shown.)

¹ HiFi read lengths, reads/data per SMRT Cell and other sequencing performance results can vary depending on DNA sample quality, insert size, P1 loading performance & movie time. Note: Refer to *SMRT Link v13.1 Kinnex single-cell troubleshooting guide* ([103-516-100](#)) for example performance metrics typically achievable with Kinnex libraries under optimal P1 loading conditions. For Revio system, we recommend aiming for ~60 – 80% P1 loading for Kinnex full-length RNA libraries.

Example SMRT Link Iso-Seq analysis results for Kinnex full-length RNA library prepared with human UHRR sample

SMRT Link Iso-Seq analysis job report – Read Classification statistics

Summary Metrics

| Value | Analysis Metric |
|------------|--|
| 81,012,128 | Reads |
| 78,815,334 | Reads with 5' and 3' Primers |
| 78,699,872 | Non-Concatamer Reads with 5' and 3' Primers |
| 78,639,000 | Non-Concatamer Reads with 5' and 3' Primers and Poly-A Tail (FLNC Reads) |
| 1,643 | Mean Length of FLNC Reads |
| 6 | Unique Primers |
| 13,135,889 | Mean Reads per Primer |
| 13,959,398 | Max. Reads per Primer |
| 12,391,746 | Min. Reads per Primer |
| 2,196,794 | Reads without Primers |
| 96.85% | Percent Bases in Reads with Primers |
| 97.28% | Percent Reads with Primers |

- **Reads:** Total number of CCS reads
- **Reads with 5' and 3' Primers:** Number of CCS reads with 5' and 3' cDNA primers detected
- **Non-Concatemer Reads with 5' and 3' Primers:** Number of nonconcatemer CCS reads with 5' and 3' primers detected
- **Non-Concatemer Reads with 5' and 3' Primers and Poly-A Tail (FLNC Reads):** Number of non-concatemer CCS reads with 5' and 3' primers and polyA tails detected. This is usually the number for full-length, nonconcatemer (FLNC) reads, unless polyA tails are not present in the sample
- **Mean Length of FLNC Reads:** Mean length of the non-concatemer CCS reads with 5' and 3' primers and polyA tails detected
- **Unique Primers:** Number of unique primers in the sequence
- **Mean Reads per Primer:** Mean number of CCS reads per primer
- **Max. Reads per Primer:** Maximum number of CCS reads per primer
- **Min. Reads per Primer:** Minimum number of CCS reads per primer
- **Reads without Primers:** Number of CCS reads without a primer
- **Percent Bases in Reads with Primers:** Percentage of bases in CCS reads in the sequence data that contain primers
- **Percent Reads with Primers:** Percentage of CCS reads in the sequence data that contain primers

Example SMRT Link Iso-Seq analysis results for Kinnex full-length RNA library prepared with human UHRR sample (cont.)

SMRT Link Iso-Seq analysis job report – Read Classification statistics

Primer Data




| Bio Sample Name ↑ | Primer Name ⌵ | CCS Reads ⌵ | Mean Primer Quality ⌵ | Reads with 5' and 3' Primers ⌵ | Non-Concatamer Reads with 5' and 3' Primers ⌵ |
|--------------------------------|-----------------------------|--------------------------|------------------------------------|---|--|
| BioSample_1 | IsoSeqX_bc01_5p--IsoSeqX_3p | 13,800,150 | 99.7 | 13,800,150 | 13,787,052 |
| BioSample_2 | IsoSeqX_bc02_5p--IsoSeqX_3p | 13,959,398 | 99.7 | 13,959,398 | 13,937,187 |
| BioSample_3 | IsoSeqX_bc03_5p--IsoSeqX_3p | 12,728,468 | 99.7 | 12,728,468 | 12,707,044 |
| BioSample_4 | IsoSeqX_bc04_5p--IsoSeqX_3p | 12,404,920 | 99.7 | 12,404,920 | 12,386,009 |
| BioSample_5 | IsoSeqX_bc05_5p--IsoSeqX_3p | 12,391,746 | 99.7 | 12,391,746 | 12,372,742 |
| BioSample_6 | IsoSeqX_bc06_5p--IsoSeqX_3p | 13,530,652 | 99.7 | 13,530,652 | 13,509,838 |
| UNASSIGNED | No Primer | 2,196,794 | 0.0 | 0 | 0 |

- **Bio Sample Name:** Name of the biological sample associated with the primer
- **Primer Name:** A string containing the pair of primer indices associated with this biological sample
- **CCS Reads:** Number of CCS reads associated with the primer
- **Mean Primer Quality:** Mean primer quality associated with the primer
- **Reads with 5' and 3' Primers:** Number of CCS reads with 5' and 3' cDNA primers detected
- **Non-Concatamer Reads with 5' and 3' Primers:** Number of non-concatemer CCS reads with 5' and 3' primers detected
- **Non-Concatamer Reads with 5' and 3' Primers and Poly-A Tail:** Number of non-concatemer CCS reads with 5' and 3' primers and polyA tails detected. This is usually the number for full-length, non-concatemer (FLNC) reads, unless polyA tails are not present in the sample.

Example SMRT Link Iso-Seq analysis results for Kinnex full-length RNA library prepared with human UHRR sample (cont.)

SMRT Link Iso-Seq analysis job report – Read Classification statistics

Primer Read Statistics



Number Of Reads Per Primer: Maps the number of reads per primer, sorted by primer ranking


Primer Frequency Distribution: Maps the number of samples with primers by the number of reads with primers

Mean Read Length Distribution: Maps the read mean length against the number of samples with primers

Example SMRT Link Iso-Seq analysis results for Kinnex full-length RNA library prepared with human UHRR sample (cont.)


SMRT Link Iso-Seq analysis job report – Read Classification statistics

Primer Quality Scores




Primer Quality Score Distribution: Histogram of primer scores

Length of Full-length Non-Concatemer Reads



Length of Full-Length Non-Concatemer Reads: Per-sample histograms of the read length distribution of non-concatemer CCS reads with 5' and 3' primers and polyA tails detected



Example SMRT Link Iso-Seq analysis results for Kinnex full-length RNA library prepared with human UHRR sample (cont.)



SMRT Link Iso-Seq analysis job report – Transcript Clustering statistics

Summary Metrics

| Sample Name ^{IT} | Number of High-Quality Isoforms ^{IT} |
|---------------------------|---|
| BioSample_1 | 1,531,163 |
| BioSample_2 | 1,549,686 |
| BioSample_3 | 1,447,494 |
| BioSample_4 | 1,449,777 |
| BioSample_5 | 1,442,229 |
| BioSample_6 | 1,521,973 |
| All Samples | 2,930,150 |

- **Sample Name:** Sample name for which the following metrics apply
- **Number of High-Quality Isoforms:** Number of consensus isoforms that have an estimated accuracy **above** the specified threshold

Length of Consensus Isoforms



- **Length of Consensus Isoforms:** Per-sample histograms of the consensus isoform lengths and the distribution of isoforms exceeding a read length cutoff. Also includes a single histogram plot for all samples.

Example SMRT Link Iso-Seq analysis results for Kinnex full-length RNA library prepared with human UHRR sample (cont.)

SMRT Link Iso-Seq analysis job report – Transcript Mapping and Classification statistics

Summary Metrics (All samples)

| Sample Name | Genes | Genes, filtered | Known Genes, filtered | Isoforms | Isoforms, filtered | Known isoforms, filtered | Novel isoforms > 1TPM, filtered |
|-------------|---------|-----------------|-----------------------|----------|--------------------|--------------------------|---------------------------------|
| BioSample_1 | 403,607 | 25,116 | 21,839 | 829,922 | 236,783 | 64,884 | 8,731 |
| BioSample_2 | 405,157 | 25,172 | 21,978 | 840,868 | 240,486 | 65,292 | 8,933 |
| BioSample_3 | 375,913 | 24,601 | 21,629 | 778,487 | 223,422 | 64,022 | 7,682 |
| BioSample_4 | 369,570 | 24,459 | 21,607 | 781,424 | 231,051 | 63,712 | 9,386 |
| BioSample_5 | 367,898 | 24,382 | 21,574 | 778,483 | 230,618 | 63,540 | 9,473 |
| BioSample_6 | 393,111 | 24,862 | 21,827 | 827,709 | 241,390 | 64,789 | 9,234 |

- **Sample Name:** Sample name for which the following metrics apply
- **Total unique genes:** The total number of unique genes across all cells.
- **Total unique genes, filtered:** The total number of unique genes, after filtering out reads based on the SQANTI transcript filtering criteria.
- **Total unique isoforms:** The total number of unique isoforms across all cells
- **Total unique isoforms, filtered:** The total number of unique isoforms across all cells, after filtering out reads based on the SQANTI transcript filtering criteria.

Example SMRT Link Iso-Seq analysis results for Kinnex full-length RNA library prepared with human UHRR sample (cont.)

SMRT Link Iso-Seq analysis job report – Transcript Mapping and Classification statistics

Transcript Classification, filtered (All samples)¹

| Category IT | Count IT | CAGE Detected IT | CAGE Detected, (%) IT | polyA Detected IT | polyA Detected, (%) IT |
|----------------|----------|------------------|-----------------------|-------------------|------------------------|
| FSM | 202364 | 91165 | 45.05% | 108319 | 53.52% |
| ISM | 314609 | 40144 | 12.75% | 191810 | 60.96% |
| NIC | 183274 | 114701 | 62.58% | 96830 | 52.83% |
| NNC | 152653 | 92299 | 60.46% | 82649 | 54.14% |
| Antisense | 3074 | 599 | 19.48% | 1890 | 61.48% |
| Fusion | 5204 | 2937 | 56.43% | 2921 | 56.12% |
| More junctions | 135 | 76 | 56.29% | 82 | 60.74% |
| Genic intron | 0 | 0 | 0.00% | 0 | 0.00% |
| Genic genomic | 2059 | 802 | 38.95% | 1177 | 57.16% |
| Intergenic | 6606 | 509 | 7.70% | 5297 | 80.18% |

- **Category:** Transcript classification² assigned by the classification and filtering tool `pigeon`, based on the [SQANTI3](#) software
- **Count:** The number of transcripts, after filtering out reads based on the SQANTI filtering criteria, in a specific classification
- **CAGE Detected:** The number of transcripts where the transcription start site falls within 50 bp of an annotated CAGE (Cap Analysis of Gene Expression) peak site
- **CAGE Detected, (%):** The percentage of transcripts where the transcription start site falls within 50 bp of an annotated CAGE peak site
- **polyA Motif Detected:** The number of transcripts where a known polyA motif is detected upstream of the transcription end site
- **polyA Motif Detected, (%):** The percentage of transcripts where a known polyA motif is detected upstream of the transcription end site

Example SMRT Link Iso-Seq analysis results for Kinnex full-length RNA library prepared with human UHRR sample (cont.)


SMRT Link Iso-Seq analysis job report – Transcript Mapping and Classification statistics

Transcript Classification Plots

Isoform distributions across structural categories:


- Distribution of the % of isoforms by structural categories

Example Revio system + SPRQ data shown.



Filter out reads based on the SQANTI3 transcript filtering criteria


Transcript Classification Plots, Filtered




Structural categories by isoform lengths:

- Histogram display of the number of isoforms by their length in kb and their structural category

Example Revio system + SPRQ data shown.





Filter out reads based on the SQANTI3 transcript filtering criteria



Example SMRT Link Iso-Seq analysis results for Kinnex full-length RNA library prepared with human UHRR sample (cont.)

SMRT Link Iso-Seq analysis job report – Transcript Mapping and Classification statistics

Gene Saturation



Gene Saturation, all genes, filtered: Saturation plot showing the level of gene saturation for **all genes**, after filtering out reads based on the SQANTI transcript filtering criteria

Gene Saturation, known genes only, filtered: Saturation plot showing the level of gene saturation, for unique **known genes only** (genes annotated in the reference annotation) per cell, after filtering out reads based on the SQANTI transcript filtering criteria

Example SMRT Link Iso-Seq analysis results for Kinnex full-length RNA library prepared with human UHRR sample (cont.)

File Downloads tab

Edit Output File Name Prefix Example:analysis-Bio Sample 4-2110

| File ↑ | Size ↓ | Type ↓ |
|---|--------|------------------|
| Non-passing reads, unaligned | 4 GB | bam |
| Report read_segmentation | 3 KB | JsonReport |
| SMRT Link Log | 13 KB | log |
| Segmented Reads, passing, unaligned | 65 GB | bam |
| UHRR_verification_DL_bc4 (demux1) (Bio Sample 4) Segmented Reads | 20 KB | ConsensusReadSet |
| Unique mapped transcripts, GFF (All Samples) | 895 MB | gff |
| Unique mapped transcripts, classification TXT (All Samples) | 728 MB | txt |
| Unique mapped transcripts, filtered, GFF (All Samples) | 599 MB | gff |
| Unique mapped transcripts, filtered, classification TXT (All Samples) | 307 MB | txt |
| Unique mapped transcripts, filtered, junctions TXT (All Samples) | 624 MB | txt |
| Unique mapped transcripts, junctions TXT (All Samples) | 771 MB | txt |

Refer to [SMRT Link user guide](#) for descriptions of downloadable output files

- These files are useful for **visualizing** isoform structures in Integrative Genomics Viewer (IGV) / UCSC genome browser and enable understanding of why an isoform is novel/known, etc.
 - GFF file containing unique mapped transcripts after filtering
 - Text file containing unique mapped transcript classifications against annotations, after filtering
 - Text file containing information about unique mapped transcript junctions, after filtering

Files shown in the File Downloads tab are available on the analysis results page. Additional files are also available on the SMRT Link server in the analysis output directory.



Technical documentation & applications support resources

Technical resources for Kinnex full-length RNA library preparation, sequencing & data analysis

RNA sample preparation resources

- Note:** The products below have not been tested or validated by PacBio but are listed here as examples of third-party kits used by other PacBio customers for isolating total RNA for PacBio RNA sequencing (Iso-Seq) applications

| Kit type | Product name |
|-----------------------------|---|
| mRNA isolation | Ambion Poly(A) Purist MAG Kit [Link] |
| Total RNA isolation | Qiagen RNeasy Plus Kits [Link] |
| | Sigma Spectrum Plant Total RNA Kit [Link] |
| | iNtRON Easy Spin Total RNA [Link] |
| | TRIzol Reagent can be used to isolate total RNA from tissues or cells, including lipid-rich and difficult samples [Link] |
| RNA stabilization & storage | RNALater is an aqueous, nontoxic tissue storage reagent that rapidly permeates tissues to stabilize and protect cellular RNA [Link] |

Technical resources for Kinnex full-length RNA library preparation, sequencing & data analysis

RNA sample preparation resources (cont.)

Some important considerations to bear in mind when isolating total RNA for Iso-Seq analysis include the following:

- RNA sample has not been exposed to high temperatures (e.g.: $>65^{\circ}\text{C}$ for 1 hour can cause a detectable decrease in sequence quality) or pH extremes (<6 or >9).
- RNA sample has an OD260/OD280 ratio ~ 2.0 .
- RNA sample has an OD260/OD230 ratio ≥ 2.0 .
- RNA sample has a RIN number ≥ 7.0 (ideally recommend ≥ 8.0).
- RNA sample has not been exposed to intercalating fluorescent dyes or ultraviolet radiation. SYBR dyes are not RNA damaging, but do avoid ethidium bromide.
- RNA sample does not contain denaturants (e.g., guanidinium salts or phenol) or detergents (e.g., SDS or Triton-X100).
- RNA sample does not contain carryover contamination from the original organism / tissue (e.g., heme, humic acid, polyphenols, etc.).
- Only use RNase-free water supplied in the reagent kit or other suppliers.
- Make aliquots of the RNA sample and TSO to avoid excessive freeze-thaw cycles.
- Thaw RNA samples and TSO on ice before use – DO NOT leave on the benchtop.
- Avoid excessive pipetting and vortexing when working with RNA.
- **Note:** RNA samples should only be shipped with dry ice.

Technical resources for Kinnex full-length RNA library preparation, sequencing & data analysis (cont.)

Kinnex full-length RNA library preparation literature & other resources

- Application note – Kinnex full-length RNA kit for isoform sequencing ([102-326-591](#))
- Procedure & checklist – Preparing Kinnex libraries using Kinnex full-length RNA kit ([103-238-700](#))
- Technical overview – Kinnex kits for single-cell RNA, full-length RNA and 16S rRNA sequencing ([103-343-700](#))
- Technical overview – Kinnex library preparation using Kinnex full-length RNA kit ([103-344-700](#))
- Video tutorial – SMRT Link Sample Setup and Run Design setup procedure for Kinnex kits [[Link](#)]
- Whitepaper – Bulk and single-cell isoform sequencing for human disease research ([102-326-576](#))

Data analysis resources

- Application note – Bioinformatics tools for full length isoform sequencing ([102-326-593](#))
- SMRT Link Kinnex full-length RNA troubleshooting guide ([103-552-100](#))
- SMRT Link Kinnex single-cell troubleshooting guide ([103-516-100](#))
- SMRT Link MAS-Seq troubleshooting guide ([102-994-400](#))
- SMRT Link software installation guide [[Link](#)]
- SMRT Link user guide [[Link](#)]
- SMRT Tools reference guide [[Link](#)]
- Video tutorial – Read Segmentation and Iso-Seq workflow in SMRT Link [[Link](#)]

Publications

- Wissel, D. et al. (2025) A systematic benchmark of high-accuracy PacBio long-read RNA sequencing for transcript-level quantification. BioRxiv [[Link](#)]
- Al'Khafaji, A.M. et al. (2024) High-throughput RNA isoform sequencing using programmable cDNA concatenation. Nature biotechnology [[Link](#)]
- Pardo-Palacios, F. et al. (2024) Systematic assessment of long-read RNA-seq methods for transcript identification and quantification. Nature Biotech [[Link](#)]
- Pardo-Palacios, F., et al., (2024) SQANTI3: curation of long-read transcriptomes for accurate identification of known and novel isoforms. Nature Methods [[Link](#)]
- Schertzer, M.D. et al. (2023) Cas13d-mediated isoform-specific RNA knockdown with a unified computational and experimental toolbox. Nature comm [[Link](#)]

Technical resources for Kinnex full-length RNA library preparation, sequencing & data analysis (cont.)

Webinars

- PacBio webinar (2025) – Detection and quantification of transcript isoforms using high-depth Kinnex full-length RNA sequencing [[Link](#)]
- PacBio video (2024) – Kinnex explained – how concatenating smaller amplicons increases throughput for PacBio HiFi sequencing [[Link](#)]
- PacBio PRISM webinar (2024) – Let’s stick together – exploring PacBio Kinnex kits [[Link](#)]
- PacBio Iso-Seq social club webinar (2022) – TappAS for isoform differential expression analysis [[Link](#)]
- PacBio Iso-Seq social club webinar (2022) – Single-cell Iso-Seq applications in cancer and neurological disorders [[Link](#)]

Example PacBio data sets

| Application | Dataset | Data type | PacBio system |
|-----------------------------------|--|----------------|-------------------------------|
| Kinnex full-length RNA sequencing | Homo sapiens – Universal human reference RNA (UHRR) [Link] | HiFi long read | Vega system |
| | Homo sapiens – Universal human reference RNA (UHRR) [Link] | HiFi long read | Revio system – SPRQ chemistry |
| | Homo sapiens – Universal human reference RNA (UHRR) [Link] | HiFi long read | Sequel II & Revio systems |
| | Homo sapiens – HG002 [Link] | HiFi long read | Revio system |
| | Homo sapiens – Heart [Link] | HiFi long read | Revio system |
| | Homo sapiens – Cerebellum [Link] | HiFi long read | Revio system |
| | Homo sapiens – Brain [Link] | HiFi long read | Revio system |



APPENDIX 1: Concatenating cDNA or gDNA amplicons using Kinnex kits to increase throughput

Concatenating cDNA or gDNA amplicons using Kinnex kits to increase throughput

Kinnex concatenation is a general method that can increase sequencing throughput for smaller amplicons¹

Benefits of Kinnex concatenation

- Increased throughput on PacBio long-read sequencers
- Retained HiFi accuracy despite throughput increase
- No change to secondary analysis – once reads are deconcatenated into S-reads, the S-reads represent the original, pre-concatenated amplicon and can be analyzed with established pipelines.

When is Kinnex concatenation appropriate?

- The balance between the amplicon size and the concatenation factor, as well as additional Kinnex library generation cost, needs to be taken into consideration.
- HiFi sequencing produces optimal yield for inserts between 15–20 kb; therefore, the throughput advantage plateaus for larger amplicon sizes exceeding 3 kb.
- You can consider concatenating amplicons using Kinnex kits if:
 - The amplicons have an **average size between 200 bp – 3 kb**
 - The amplicons have **molecular ends that are either directly compatible or can be re-amplified to establish Kinnex compatibility**

Recommended Kinnex kit based on average amplicon sizes.

| Average amplicon size | Example | Recommended Kinnex kit | Expected Kinnex library size |
|-----------------------|--|--------------------------------------|------------------------------|
| 600 – 1000 bp | 10x single-cell cDNA | Kinnex single-cell RNA kit (16-fold) | 10 – 16 kb |
| 1 – 2 kb | Full-length 16S | Kinnex 16S rRNA kit (12-fold) | ~19 kb |
| 2 – 3 kb | Bulk cDNA | Kinnex full-length RNA kit (8-fold) | 15 – 20 kb |
| >3 kb | Not recommended for Kinnex concatenation | | |

Technical note

CONCATENATING AMPLICONS USING PACBIO KINNEX KITS TO INCREASE THROUGHPUT

Overview

This technical note describes the principles and advantages of concatenating amplicons using PacBio® Kinnex™ kits.

Note: This technical note is intended as a guide for best practices and to report typical example results (see Appendix). PacBio does not guarantee the success of concatenating amplicons that are not officially supported by the Kinnex kits, which currently include:

- Full-length cDNA generated using the Iso-Seq® express 2.0 kit with the Kinnex full-length RNA kit
- Full-length 16S sequences generated according to official protocol with the Kinnex 16S rRNA kit
- Single-cell cDNA generated using compatible single-cell platforms with the Kinnex single-cell RNA kit

How Kinnex works

Figure 1. How Kinnex works. Watch the Kinnex animation to learn more.

What is Kinnex?

The Kinnex kits are based on the method developed by Al'Khafaji et al. (2023). Originally termed multiplexed array isomorph sequencing (MAS-ISO-seq, or MAS-Seq), this method was developed to concatenate cDNAs into longer fragments suitable for long-read sequencing. It takes advantage of the disparity between optimal DNA fragment lengths for HiFi sequencing, 15–20 kb, and smaller sizes of transcript cDNA (1–10 kb with typical estimated average length of 2 kb for a human transcriptome). The original MAS-ISO-Seq method concatenated single-cell cDNA from the 10x platform, which has an average read length of 600–800 bp, and used a 15-fold concatenation array to increase throughput. PacBio commercialized the MAS-Seq method and adjusted the concatenation factors according to different average amplicon sizes (Table 1).

| Average amplicon size | Example | Recommended kit | Expected library insert size |
|-----------------------|--|--------------------------------------|------------------------------|
| 600–1000 bp | 10x single-cell cDNA | Kinnex single-cell RNA kit (16-fold) | 10–16 kb |
| 1–2 kb | Full-length 16S rRNA | Kinnex 16S rRNA kit (12-fold) | ~19 kb |
| 2–3 kb | Bulk cDNA | Kinnex full-length RNA kit (8-fold) | 15–20 kb |
| >3 kb | Not recommended for Kinnex concatenation | | |


Table 1. Recommended Kinnex kit based on average amplicon sizes.

PacBio

Technical note – Concatenating cDNA or gDNA amplicons using Kinnex kits to increase throughput ([102-326-636](#))

Kinnex concatenation workflow overview

Follow Kinnex library prep protocol documentation for specific details on concatenating cDNA or gDNA amplicons using Kinnex kits



Kinnex concatenation workflow. Amplicons must be generated or amplified to have Kinnex-compatible molecular ends before continuing to the Kinnex PCR and array formation step. Kinnex libraries should be sequenced with the appropriate sequencing chemistry and run configurations. Once de-concatenated using Read Segmentation in SMRT Link, the individual amplicons can be analyzed using amplicon-specific workflows.

Kinnex concatenation procedural notes¹

- To establish Kinnex compatibility, **amplicons are required to have Kinnex-compatible molecular ends** (see next section)
 - Once these are generated, choose the appropriate Kinnex kit based on the recommended concatenation factor listed in Table on previous slide and proceed with Kinnex PCR
- The **Kinnex PCR** steps consist of parallel PCR reactions per sample [i.e., 8, 12, or 16 reactions based on the Kinnex kit chosen] using premixed Kinnex primer pairs
 - The resulting PCRs generate amplified DNA products containing programmable sequences at both ends
- In the **Kinnex array formation** step, library inserts containing programmable ends are assembled to generate a linear array.
 - Further, the addition of barcoded Kinnex terminal adapters result in the formation of complete, full-length array SMRTbell templates along with partial arrays.
 - Subsequent nuclease treatment removes partial arrays to retain only full array SMRTbell templates for achieving optimal sequencing yield.
- **Note:** Kinnex terminal adapters are different from standard SMRTbell adapters and hence require the **Kinnex sequencing primer** (103-179-000) during the “Annealing, Binding, and Cleanup (ABC)” step for optimal sequencing results.
- Once HiFi reads are generated, **Read Segmentation** will produce the segmented reads (S-reads) that represent the original unconcatenated amplicons, which can be used for further analysis.

Establishing Kinnex-compatible molecular ends

Amplicons must be generated or amplified to have Kinnex-compatible molecular ends before proceeding with Kinnex library prep procedure

- To be compatible with the Kinnex workflow, amplicons must be generated with sequence-defined ends as depicted in the underlined portion in the figure below
- Optional barcodes (such as sample indices, UMIs and single cell barcodes) should be placed internally between the Kinnex handles and the amplicon-specific primers.
- Kinnex handles may be present already in certain amplicons, such as 10x Single Cell Gene Expression libraries or the Kinnex 16S amplicons, or can be added by PCR amplification, such as for Parse Evercode single-cell libraries

Kinnex FWD primer CTACACGACGCTCTTCCGATCT - [optional barcodes] - [amplicon specific FWD primers]

Kinnex REV primer AAGCAGTGGTATCAACGCAGAG - [optional barcodes] - [amplicon specific REV primers]

Example 1. Kinnex 16S forward and reverse primer sequences

Kinnex 16S FWD01 CTACACGACGCTCTTCCGATCT - GATCGAGTCA - AGRGTTYGATYMTGGCTCAG

Kinnex 16S REV13 AAGCAGTGGTATCAACGCAGAG - TCATCGACGT - RGYTACCTTGTACGACTT

Example 2. Iso-Seq express 2.0 forward and reverse primer sequences

IsoSeqX bc01 FWD CTACACGACGCTCTTCCGATCT - ACTACAC - GCAATGAAGTCGCAGGGTTGGG

IsoSeqX REV AAGCAGTGGTATCAACGCAGAGTAC

Schematic for Kinnex-compatible primers. Kinnex handles (5' to 3') are shown in black underline and must be present at the ends of the amplicons to be compatible with Kinnex concatenation. Optional barcodes can be included internally. Amplicon-specific primers (and optional internal barcode sequences) must be designed to avoid strong secondary structures in the context of Kinnex handles.

Resources for Kinnex library prep¹

For using Kinnex full-length RNA kit (PN: 103-072-000) for 8-fold concatenation:

- Procedure & checklist – Preparing Kinnex libraries using the Kinnex full-length RNA kit ([103-238-700](#))
- Technical overview – Kinnex library preparation using Kinnex full-length RNA kit ([103-344-700](#))

For using Kinnex 16S rRNA kit (PN: 103-072-100) for 12-fold concatenation:

- Procedure & checklist – Preparing Kinnex libraries using 16S rRNA amplicons ([103-238-800](#))
- Technical overview – Kinnex library preparation using Kinnex 16S rRNA kit ([103-344-800](#))

For using Kinnex single-cell RNA kit (PN: 103-072-200) for 16-fold concatenation:

- Procedure & checklist – Preparing Kinnex libraries using Kinnex single-cell RNA kit ([103-254-300](#))
- Technical overview – Kinnex library preparation using Kinnex single-cell RNA kit ([103-344-600](#))




APPENDIX 2: PacBio compatible Kinnex library preparation workflows

Targeted enrichment of Kinnex full-length RNA libraries with IDT xGen hybridization and wash kit v3

Targeted enrichment of Kinnex full-length RNA libraries with IDT xGen hybridization and wash kit v3: Getting started

Application-specific educational literature




PacBio literature website [\[Link\]](#)

Application-specific brochures, informational guides and other product literature containing best practices recommendations for library preparation and data analysis workflows.




Application-specific protocol documentation



Technical note – Targeted enrichment of Kinnex full-length RNA libraries with IDT xGen hybridization and wash kit v3 [\(102-326-632\)](#) [\[PacBio\]](#)

Technical documentation describing the procedure for targeted enrichment using IDT xGen Hyb and Wash kit v3 for preparing Kinnex full-length RNA libraries for PacBio sequencing.

Application-specific technical overviews



Technical overview: Kinnex library preparation using Kinnex full-length RNA kit [\(103-344-700\)](#)

Technical overview presentations describe sample preparation details for constructing HiFi libraries for specific applications. Example sequencing performance data for a given application are also summarized.

Technical reference for targeted enrichment of Kinnex full-length RNA libraries with IDT xGen hybridization and wash kit v3

Technical note – Targeted enrichment of Kinnex full-length RNA libraries with IDT xGen hybridization and wash kit ([102-326-632](#)) describes the procedure for targeted enrichment using IDT xGen Hyb and Wash kit v3 for Kinnex full-length RNA libraries

Technical note

TARGETED ENRICHMENT OF KINNEX FULL-LENGTH RNA LIBRARIES WITH IDT xGEN HYBRIDIZATION AND WASH KIT V3

Overview

For RNA sequencing, targeted enrichment is considered a favorable alternative to whole transcriptome sequencing because it can capture more isoform diversity of targeted genes, requiring less sequencing depth and enabling higher sample multiplexing. The PacBio® Iso-Seq® method for full-length RNA sequencing has been used to capture single genes¹, highly similar paralogs or pseudogenes², and low-abundance long non-coding RNAs for annotation purposes³. Targeted enrichment may be the right choice when the goal of sequencing is to:

- Identify and quantify transcripts from a known set of genes
- Characterize alternative 5' starts and 3' ends
- Detect lowly expressed transcripts
- Maximize sample multiplexing capacity
- Reduce bioinformatics analysis time

This technical note describes the procedure for targeted enrichment using IDT xGen Hyb and Wash kit v3 for Kinnex™ full-length RNA libraries.

Required materials and equipment

- For a comprehensive list of materials and equipment, please refer to the PacBio Kinnex full-length RNA protocol⁴ and the IDT xGen protocol⁵.
- For cDNA generation, use the Iso-Seq express 2.0 kit, which is compatible with the IDT targeted enrichment and Kinnex full-length RNA library construction. The input to Iso-Seq express 2.0 is 300 ng of total RNA, ideally with RIN ≥ 7. Lower-quality RNA may result in shorter cDNA, while total RNA lower than 100 ng may result in the incomplete capture of the full diversity of the targeted genes.

Library preparation

Library preparation will follow the Kinnex full-length RNA protocol with modifications outlined here.

1. Generate full-length cDNA using the Iso-Seq express 2.0 kit.

2. Enrich for targeted genes using the IDT xGen Hyb Wash v3 protocol.

3. Continue through the Kinnex PCR and pooling step (step 4) in the Kinnex full-length RNA protocol.

4. The targeted Kinnex library is compatible for sequencing on Sequel® II/IIe, Revio®, and Vega™ systems.

Full-length cDNA

Generate up to 12-plex cDNA using the Iso-Seq express 2.0

Target enrichment

Enrich for targeted genes using IDT customer probes and then xGen Hyb Wash v3 kit

Kinnex library

Generate Kinnex full-length RNA library using the Kinnex full-length RNA kit

Sequencing

Sequence on Sequel IIe, Vega, or Revio systems

Analyze

Analyze with the Read Segmentation and Iso-Seq workflow in SMRT® Link

Figure 1. Library workflow for generating targeted Kinnex full-length RNA libraries using the IDT xGen Hyb Wash v3 kit.

SMRT® Link informatics and analysis recommendations

Targeted Kinnex full-length RNA data can be analyzed the same way as whole transcriptome Kinnex datasets using the SMRT Link Read Segmentation and Iso-Seq workflow or through the command line. Note, however, that these workflows do not have specific analyses for targeted gene lists and will output gene and isoform information for all detected transcripts.

IDT

PacBio

```
graph TD; 1[1 Input RNA quality control] --> 2[2 cDNA synthesis + cleanup]; 2 --> 3[3 cDNA amplification + cleanup]; 3 --> 4[4 Target enrichment with IDT xGen kit]; 4 --> 5[5 Kinnex PCR]; 5 --> 6[6 Kinnex array formation + cleanup]; 6 --> 7[7 Nuclease treatment + cleanup]; 7 --> 8[8 Annealing, binding and cleanup ABC];
```

PacBio Procedure & checklist

Preparing Kinnex libraries using the Kinnex full-length RNA kit ([103-238-700](#))

PacBio Technical note

Targeted enrichment of Kinnex full-length RNA libraries with IDT xGen hybridization and wash kit v3 ([102-326-632](#))¹

PacBio Procedure & checklist

Preparing Kinnex libraries using the Kinnex full-length RNA kit ([103-238-700](#))

Overview


- For RNA sequencing, **targeted enrichment** may be the right choice when the goal of sequencing is to:
 - Identify and quantify transcripts from a known set of genes
 - Characterize alternative 5' starts and 3' ends
 - Detect lowly expressed transcripts
 - Maximize sample multiplexing capacity
 - Reduce bioinformatics analysis time
- **Technical note – Targeted enrichment of Kinnex full-length RNA libraries with IDT xGen hybridization and wash kit v3** ([102-326-632](#)) describes the procedure for targeted enrichment using IDT xGen Hyb and Wash kit v3 for Kinnex full-length RNA libraries

PacBio ¹ For step-by-step IDT hybridization capture workflow details, refer to the IDT xGen Hybridization and Wash Kit v3 (EAP) for PacBio Targeted Kinnex protocol ([RUO24-3081_001](#)).

82


Workflow overview for targeted enrichment of Kinnex full-length RNA libraries with IDT xGen hybridization and wash kit v3

Follow PacBio Kinnex library prep protocol and IDT target enrichment procedure to generate full-length RNA hybrid capture libraries suitable for HiFi sequencing on PacBio long-read systems



Workflow overview for targeted enrichment of Kinnex full-length RNA libraries with IDT xGen hybridization and wash kit v3 (cont.)

Follow PacBio Kinnex library prep protocol and IDT target enrichment procedure to generate full-length RNA hybrid capture libraries suitable for HiFi sequencing on PacBio long-read systems




IDT xGen hybrid capture procedure overview¹

Technical note – Targeted enrichment of Kinnex full-length RNA libraries with IDT xGen hybridization and wash kit v3 (102-326-632)



IDT xGEN hybrid capture procedure key steps¹

- Add Blocker Master Mix to each cDNA sample (500 ng²) and dry down
- Add Hybridization Master Mix to each sample well and incubate for desired hybridization time.³
- Prepare hybridization wash buffer and capture beads within 2 hours before use
- Add (streptavidin) capture beads to sample (30 min)
- Perform 3 rounds of washing with wash buffers
- Perform post-capture PCR reaction + cleanup
- Perform DNA quantification QC using Qubit assay and DNA sizing QC using Agilent TapeStation or other equivalent system



¹ For hybridization capture workflow details, refer to the IDT xGen Hybridization and Wash Kit v3 (EAP) for PacBio Targeted Kinnex protocol (RUO24-3081_001).




² 100 ng to 6 µg cDNA inputs can be used for short hybridizations (1-2 hrs). Generally, we recommend not to exceed 2.5 µg total input for 1-hr hybridization time for very large panels (>30 Mb).

³ This procedure supports short hybridization times of only 1 hr as well as the standard 4- or 16-hr incubations. We recommend users empirically test the shorter hybridization time for their workflows.

SMRT Link informatics and secondary analysis recommendations for targeted Kinnex full-length RNA data

Use SMRT Link Read Segmentation and Iso-Seq application to perform isoform-classification analysis to identify novel genes & isoforms with abundance information

- Targeted Kinnex full-length RNA data can be analyzed the same way as whole transcriptome Kinnex datasets using the [SMRT Link Read Segmentation and Iso-Seq workflow](#) or through the command line.
- Note, however, that these workflows do not have specific analyses for targeted gene lists and will output gene and isoform information for all detected transcripts.





APPENDIX 3: SMRT Link run design procedure for Sequel IIe system

SMRT Link Run Design procedure for Sequel IIe system

Sample information and run information

- Select desired **Kinnex application** from the **Application** field drop-down menu
- The following fields are **auto-populated** with default recommended values and high-lighted in **green**:

☐ **SMRTbell Adapter Design**

→ SMRTbell Kinnex Prep Kit

☐ **Binding Kit**

→ Sequel II Binding Kit 3.2

☐ **Sequencing Kit**

→ Sequel II Sequencing Plate 2.0 (4 rxn or 1 rxn)

☐ **DNA Control Complex**

→ Sequel II DNA Internal Control Complex 3.2

☐ **Movie Time per SMRT Cell**

→ 30 hrs

☐ **Pre-Extension Time**

→ 2 hrs

SMRTbell Adapter Design field determines which adapter finding algorithm is used during post-primary analysis¹

Default SMRTbell adapter design for Kinnex samples is **SMRTbell Kinnex Prep Kit**

Recommended OPLC for Sequel II/Ile Kinnex library samples is **40 – 60 pM**

Recommended movie time = **30 hrs**

Select desired Kinnex application type from drop-down menu

SAMPLE 1: Kinnex full-length RNA library demo , A01, 30 hour movie, 16000 Copy Delete

Import from Sample Setup Select Sample

Application Required Kinnex full-length RNA

Well Sample Name Required Kinnex full-length RNA library demo

Bio Sample Name

Sample Comment

Sample Well A01

SMRTbell Adapter Design Required SMRTbell® Kinnex Prep Kit

Binding Kit Required Sequel® II Binding Kit 3.2

Sequencing Kit Required Sequel® II Sequencing Plate 2.0 (4 rxn)

DNA Control Complex Sequel® II DNA Internal Control Complex 3.2

Insert Size (bp) Required 16000

Recommended Concentration on Plate (pM) 40 – 60 pM

On-Plate Loading Concentration (pM) Required 45

Movie Time per SMRT Cell (hours) 30

Use Pre-Extension ☒ YES ☐ NO

Pre-Extension Time (hours) 2

CCS Analysis will be performed on-instrument to produce HiFi .bam files.

Example sample information entered into a Sequel IIe system run design worksheet for a Kinnex full-length RNA library sample.

SMRT Link Run Design procedure for Sequel IIe system (cont.)

Advanced options

- For all Kinnex library samples, leave the following **Advanced Options** fields at their **default settings**
 - ☐ **Use Adaptive Loading**
→ YES
 - ☐ **Loading Target (P1 + P2)**
→ 0.85
 - ☐ **Maximum Loading Time**
→ 2 hours
 - ☐ **CCS Analysis Output - Include Low Quality Reads**
→ NO
 - ☐ **CCS Analysis Output - Include Kinetics Information**
→ NO
 - ☐ **Pre-Extension Time**
→ 2 hrs
- If desired, specify to use an alternative project folder for the **Add Data to Project** field

Advanced Options

Use Adaptive Loading ☒ YES ☐ NO

Loading Target (P1 + P2) 0.85

Maximum Loading Time (hours) 2

CCS Analysis Output - Include Low Quality Reads ☐ YES ☒ NO

CCS Analysis Output - Include Kinetics Information ☐ YES ☒ NO

Add Data to Project General Project

Example default Advanced Options settings entered into a Sequel IIe system run design worksheet for a Kinnex full-length RNA library sample.

SMRT Link Run Design procedure for Sequel IIe system (cont.)

Barcoded sample options

- For Kinnex library samples, can leave most **Barcoded Sample Options** fields at their **default settings**

Specify Bio Sample Names, either interactively or by downloading a CSV file (**Interactively** or **From a file**)

If desired, specify to perform barcode demultiplexing on-instrument or in SMRT Link (default = **On-instrument** for Sequel IIe system)

Can leave most of these fields at their default values

Barcoded Sample Options

Sample Is Barcoded ☒ YES ☐ NO

Barcode Set Required MAS SMRTbell barcoded adapters (v2)

Same Barcodes on Both Ends of Sequence ☒ YES ☐ NO

Assign Bio Sample Names to Barcodes Required ☒ Interactively ☐ From a File

Demultiplex Barcodes ☒ ON INSTRUMENT ☐ IN SMRT LINK ☐ DO NOT GENERATE

Example default Barcoded Sample Options settings entered into a Sequel IIe system run design worksheet for a Kinnex full-length RNA library sample.



www.pacb.com

| Revision history (description) | Version | Date |
|---|---------|---------------|
| Added procedural summary illustrations to clarify Kinnex library prep protocol steps, updated example Kinnex sequencing performance data for PacBio long-read systems (Revio system + SPRQ chemistry and Vega system) and incorporated information about PacBio compatible RNA sequencing workflows using third-party kit products. | 02 | December 2025 |

Research use only. Not for use in diagnostic procedures. © 2025 Pacific Biosciences of California, Inc. ("PacBio"). All rights reserved. Information in this document is subject to change without notice. PacBio assumes no responsibility for any errors or omissions in this document. Certain notices, terms, conditions and/or use restrictions may pertain to your use of PacBio products and/or third-party products. Refer to the applicable PacBio terms and conditions of sale and to the applicable license terms at pacb.com/license. Pacific Biosciences, the PacBio logo, PacBio, Circulomics, Omniome, SMRT, SMRTbell, Iso-Seq, Sequel, Nanobind, SBB, Revio, Onso, Apton, Kinnex, PureTarget, SPRQ, and Vega are trademarks of PacBio.