



Technical overview – Kinnex library preparation using Kinnex single-cell RNA kit

Sequel II and IIe systems ICS v11.0

Revio system ICS v13.3

SMRT Link v25.3

PN 103-344-600 Rev 02 | December 2025

Technical overview

Kinnex library preparation using Kinnex single-cell RNA kit

1. Kinnex single-cell RNA method overview
2. Kinnex single-cell RNA library preparation workflow details
3. Kinnex single-cell RNA sequencing preparation workflow details
4. Kinnex single-cell RNA example sequencing performance data
5. Kinnex single-cell RNA data analysis workflow overview
6. Technical documentation & applications support resources
7. APPENDIX 1 – Concatenating cDNA or gDNA amplicons using Kinnex kits to increase throughput
8. APPENDIX 2 – PacBio compatible Kinnex library preparation workflows
 - i. Preparing Kinnex single-cell libraries with Parse Evercode Whole Transcriptome kits
 - ii. Preparing Kinnex single-cell libraries with Twist exome enrichment kits
9. APPENDIX 3 – SMRT Link run design procedure for Sequel IIe system

Kinnex library preparation using Kinnex single-cell RNA kit: Getting started



PacBio

Application note

Kinnex single-cell RNA kit for single-cell isoform sequencing

Introduction

Understanding cell heterogeneity at the isoform level is critical for both basic and disease research. Short reads can only capture gene-level information, while other long-read technologies lack the accuracy for accurate unique molecular identifiers (UMI) and cell barcode (CBC) identification. PacBio[®] HiFi reads sequence full-length RNA isoforms along with single-cell barcode and UMI information, revealing extraordinary insight into single-cell biology.

The Kinnex[™] single-cell RNA kit takes as input single-cell cDNA and outputs a sequencing-ready library that results in a 16-fold throughput increase compared to regular single-cell Iso-Seq[®] libraries. Combined with isoform-aware single-cell analysis SMRT[™] Link software, PacBio offers cost-effective single-cell isoform sequencing that does not require orthogonal sequencing methods. The SMRT Link software supports bioinformatics analysis to produce an isoform-level single-cell data matrix compatible with tertiary analysis software.

© 2025 PacBio. All rights reserved. Research use only. Not for use in diagnostic procedures.
102-326-549 REV08 JULY 2025

PacBio

Application Note - Kinnex single-cell RNA for single-cell isoform sequencing (102-326-549)

Summary overview of application-specific library preparation and data analysis workflow recommendations.

PacBio

Preparing Kinnex[™] libraries using Kinnex single-cell RNA kit

Procedure & checklist

Before you begin

This procedure describes the workflow for constructing single-cell Kinnex libraries from 10x Chromium 3' or 5' cDNA using the Kinnex single-cell RNA kit for library prep and sequencing on PacBio[®] Sequel II, Sequel IIe, Vega[™], and Revio[™] systems.

This kit is intended for use with single-cell cDNA generated using the 10x Chromium Next GEM Single Cell 3' kit (v3.1 or v4) or 10x Chromium Next GEM Single Cell 5' kit (v2, v3). For a full list of compatibility with 10x products, please consult the [10x website](#).

Overview	
Samples per kit	12
Workflow time	3 days for up to 12 samples
Number of SMRT [™] Cells per Kinnex library preparation	>8 SMRT Cells for Revio [®] using SPRQ [™] chemistry >2 SMRT Cells for the Vega [®] or Revio [®] (non-SPRQ) systems >4 SMRT Cells for Sequel [®] II/III [®]

cDNA input	
Quantity	>15 ng 10x cDNA per sample cDNA concentration should be >1ng/μL with up to 15 μL in volume. See step 2.1 for 10x cDNA input requirement.
Average segment lengths	500–1,100 bp
Average 16-segment array lengths	10–17 kb

© 2025 PacBio. All rights reserved. Research use only. Not for use in diagnostic procedures.
102-254-300 REV08 OCT2025

PacBio

Procedure & checklist – Preparing Kinnex libraries using Kinnex single-cell RNA kit (102-254-300)

Technical documentation containing application-specific library preparation protocol details.

PacBio

Technical overview – Kinnex library preparation using Kinnex single-cell RNA kit

Sequel II and IIe systems ICS v11.0
Revio system ICS v13.3
SMRT Link v25.3

PN 103-344-600 Rev 02 | December 2025

Example sequencing performance for Kinnex single-cell RNA libraries prepared with human cDNA [Revio system + SPRQ chemistry¹]

Kinnex single-cell RNA 5' library for PBMC single-cell cDNA sample from 10x Chromium Single Cell Universal 5' kit

Raw Data Report		HiFi Read Length		Read Segmentation Metrics	
Raw Data Yield	1.53B (G)	HiFi Reads	9.4 M	Input HiFi Reads	9,400,000
Mean Polymerase Read Length	902.2 bp	HiFi Reads Yield	144.0 G/G	Segmentation Yield (16 Segments)	144.0 G/G
PI	95%	Mean HiFi Read Length	10,330	Mean Length of Segments	971 bp
P2	95%	Median HiFi Read Quality	Q29	Percent of reads with full entry	92.55%

© 2025 PacBio. All rights reserved. Research use only. Not for use in diagnostic procedures.

PacBio

Technical Overview: Kinnex library preparation using Kinnex single-cell RNA kit (103-344-600)

Technical overview presentations describe sample preparation details for constructing Kinnex HiFi libraries for specific applications. Example sequencing performance data for a given application are also summarized.

Single-cell cDNA synthesis

Use 10x Genomics Chromium or other system to perform single-cell cDNA synthesis (3,000 – 20,000 cells input)

Kinnex library preparation (Kinnex single-cell RNA kit)

Use ≥15 ng of single-cell cDNA input to generate Kinnex library containing 16-segment array

SMRT sequencing (Sequel IIe, Vega & Revio systems)

Perform ABC¹ and sequence Kinnex libraries on PacBio long-read systems

Data analysis (SMRT Link)


Use Read Segmentation data utility to split arrayed transcript HiFi reads
Use Single-cell Iso-Seq analysis application to identify novel genes and isoforms



Kinnex single-cell RNA method overview

Kinnex single-cell RNA method overview


Use Kinnex single-cell RNA kit to perform high-accuracy, single-cell isoform sequencing with PacBio long-read systems



- 10x Chromium Single Cell 3' kit (v3.1) and 5' kit (v2)¹
- 15–75 ng cDNA input
- 3,000 to 20,000 target cell recovery
- 2-day Kinnex library preparation using **Kinnex single-cell RNA kit**
- Barcoded Kinnex adapters support up to 4-plex multiplexing
- SMRT Link Run Design support for 'Kinnex single-cell RNA' application type option with auto-analysis (read segmentation + single-cell isoform analysis)
- SMRT Link single-cell Iso-Seq isoform-classification software to identify novel genes and isoforms
- Output compatible with tertiary single-cell analysis tools (e.g., *Seurat*, *Scanpy*, *Kana*)


Kinnex single-cell RNA method overview

Single-cell cDNA sample preparation¹



¹ Refer to [10x Genomics Support](https://www.10xgenomics.com/support) website to download 10x Chromium user guides and other documentation.

Kinnex library prep, sequencing & analysis



PacBio

Preparing MAS-Seq libraries using MAS-Seq for 10x Single Cell 3' Kit

Procedure & checklist (103-254-300)

Before you begin

This procedure describes the workflow for constructing single-cell MAS-Seq libraries using the 10x Chromium 3' Kit for library prep and sequencing on PacBio systems.


Reads are generated for use with single-cell cDNA generated using the 10x Chromium Next GEM Single Cell 3' Kit (10x) as described in the 10x Genomics Single Cell 3' Kit User Guide (https://www.10xgenomics.com/tech-papers/10x-chromium-single-cell-3-kit-user-guide).

This kit may be modified to work with the 10x Chromium Next GEM Single Cell 5' Kit and the Chromium Single Cell 5' Kit (https://www.10xgenomics.com/tech-papers/10x-chromium-single-cell-5-kit-user-guide).

Parameter	Value
Library type	10x Single Cell 3' Kit
Sequencing mode	10x Single Cell 3' Kit
Read length	10x Single Cell 3' Kit
Read depth	10x Single Cell 3' Kit
Read quality	10x Single Cell 3' Kit
Average read length	10x Single Cell 3' Kit
Average read depth	10x Single Cell 3' Kit

² ABC = Anneal sequencing primer / Bind polymerase / Complex cleanup

Single-cell cDNA sample preparation¹




10x & PacBio key protocol steps



Procedure & checklist (103-254-300)

Kinnex library prep, sequencing & analysis



¹ Refer to [10x Genomics Support](#) website to download 10x Chromium user guides and other documentation.

² Kinnex adapter barcode sequences can be downloaded from [SMRT Link](#) Data Management module.


Kinnex single-cell RNA library preparation procedure description

Procedure & checklist – Preparing Kinnex libraries using the Kinnex single-cell RNA kit (103-254-300) describes the workflow for constructing single-cell Kinnex libraries from 10x Chromium 3' or 5' cDNA using the Kinnex single-cell RNA kit (103-072-200) for library prep and sequencing on PacBio long-read systems

Overview	
Samples per kit	12
Workflow time	3 days for up to 12 samples
Number of SMRT® Cells per Kinnex library preparation	>8 SMRT Cells for Revio® using SPRQ™ chemistry >2 SMRT Cells for the Vega® or Revio® (non-SPRQ) systems >4 SMRT Cells for Sequel® II/Ile®
cDNA input	
Quantity	>15 ng 10x cDNA per sample cDNA concentration should be >1ng/μL with up to 15 μL in volume. See step 2.1 for 10x cDNA input requirement.
Average segment lengths	500–1,100 bp
Average 16-segment array lengths	10–17 kb



Kinnex single-cell RNA kit
103-072-200 (12 rxn)



Kinnex single-cell RNA library template (~10–17 kb)
Contains 16 concatenated full-length cDNA segments

For use with single-cell cDNA generated with 10x Chromium Next GEM Single Cell 3' kit v3.1 or 10x Chromium Next GEM Single Cell 5' kit v2, standard throughput¹

Kinnex single-cell RNA library preparation procedure supports up to **4-plex sample multiplexing** through use of 4 different barcoded Kinnex SMRTbell adapters²

PacBio

Preparing Kinnex™ libraries using Kinnex single-cell RNA kit

Procedure & checklist

Before you begin

This procedure describes the workflow for constructing single-cell Kinnex libraries from 10x Chromium 3' or 5' cDNA using the Kinnex single-cell RNA kit (103-072-200) for library prep and sequencing on PacBio® Sequel® II, Sequel Ile, and Revio™ systems.

This kit is intended for use with single-cell cDNA generated using the 10x Chromium Next GEM Single Cell 3' kit v3.1 or 10x Chromium Next GEM Single Cell 5' kit v2, standard throughput. It has not been tested for use on low throughput (LT) or high throughput (HT) kits which are currently unsupported.

Overview	
Samples per kit	12
Workflow time	3 days for up to 12 samples

cDNA input	
Quantity	>15 ng per library cDNA concentration should be >1ng/μL with up to 15 μL in volume. See step 2.1 for 10x cDNA input requirement.
Average segment lengths	500–1,000 bp
Average 16-segment array lengths	10–15 kb

© 2024 PacBio. All rights reserved. Research use only. Not for use in diagnostic procedures.
103-254-300 REV 02 MAR2024

PacBio

PacBio [Documentation \(103-254-300\)](#)

- Kinnex full-length RNA library prep protocol uses **Kinnex single-cell RNA kit**
→ **Do not use** SMRTbell prep kit 3.0 with this protocol

¹ This kit is intended for use with single-cell cDNA generated using the 10x Chromium Next GEM Single Cell 3' kit (v3.1 or v4) or 10x Chromium Next GEM Single Cell 5' kit (v2, v3). For a full list of compatibility with 10x products, please consult the [10x website](#).

² Kinnex adapter barcode sequences can be downloaded from [SMRT Link](#) Data Management module.

Kinnex single-cell RNA kit bundle components

Kinnex single-cell RNA kit bundle provides full support for Kinnex library prep workflow

Kinnex single-cell RNA kit (103-072-200)

Includes Kinnex PCR kit, Kinnex concatenation and ancillary DNA cleanup reagents needed for incorporation of Kinnex segmentation adapters and Kinnex array formation for generating Kinnex single-cell RNA libraries from input 10x Chromium Single Cell 5' and 3' cDNA.

Kinnex single-cell RNA kit components

Component	Description
1	 Kinnex capture beads kit (12 rxn) <ul style="list-style-type: none">Contains reagents for removing template-switch oligo (TSO) artifacts from single-cell cDNA
2	 Kinnex single cell concatenation kit (12 rxn) <ul style="list-style-type: none">Contains reagents for Kinnex array formation and SMRTbell template constructionIncludes barcoded Kinnex adapter mixes (bcM0001 – bcM0004)Also contains Kinnex capture primer oligos for TSO artifact removal
3	 SMRTbell cleanup beads <ul style="list-style-type: none">For DNA cleanup
4	 Elution buffer <ul style="list-style-type: none">For DNA cleanup

Kinnex single-cell RNA experimental design considerations

Kinnex single-cell RNA application use case recommendations for PacBio systems

	Sequel II and Ile systems	Vega system	Revio system + SPRQ chemistry
Experimental goal	Characterize alternative splicing in single cells / cell types		
Transcript reads per SMRT Cell	30-40 M reads	50-60 M reads	100-120 M reads
Sample multiplexing ¹	Not recommended	Not recommended	Up to 2 samples per Revio SMRT Cell (2-plex)
Cell input into 10x Chromium single cell 3' or 5' cDNA generation workflow	3,000 – 10,000 cells for running a single (non-multiplexed) sample on one Sequel II SMRT Cell 8M	3,000 – 10,000 cells for running a single (non-multiplexed) sample on one Vega SMRT Cell	<10,000 cells per sample if multiplexing 2 samples per Revio SMRT Cell (2-plex) 10,000 – 20,000 cells per sample if running a single (non-multiplexed) sample on 1–2 Revio SMRT Cells
Expected coverage	Obtain 3,000 – 10,000 unique reads/single cell	Obtain 3,000 – 10,000 unique reads/single cell	Obtain up to ~20,000 unique reads/single cell
Kinnex library prep protocol	Procedure & checklist – Preparing Kinnex libraries using Kinnex single-cell RNA kit (103-254-300)		
Single-cell cDNA input into Kinnex library prep workflow	15-75 ng of cDNA from (1) 10x Chromium Single Cell Universal 3' or 5' kit; Single Cell Next GEM (v3.1 3' or v2 5'); (2) 10x Visium HD Spatial Gene Expression kit; (3) Parse Evercode WT or WT mini kit; or (4) ArgenTag Single-Cell RNA Library Kit		
SMRT Link data analysis workflows	Read Segmentation and Single-cell Iso-Seq Analysis		
Community data analysis tools	Annotation & quantification: SQANTI3 / Differential analysis: TappAS / Fusion calling: pbfusion / Visualization: SWAN		



Kinnex single-cell RNA library preparation workflow details

Procedure & checklist – Preparing Kinnex libraries using the Kinnex single-cell RNA kit (103-254-300)

Procedure & checklist [103-254-300](#) describes the workflow for constructing single-cell Kinnex libraries from 10x Chromium 3' or 5' cDNA using the **Kinnex single-cell RNA kit** for library prep and sequencing on PacBio long-read systems

Procedure & checklist contents

1. **General best practices** for reagent & sample handling and **10x single cell cDNA input recommendations**.
2. Enzymatic workflow steps for **removal of template-switch oligo (TSO) artifacts** from input 10x single cell cDNA samples.
3. Enzymatic workflow steps for **construction of 16-segment Kinnex arrays** from 10x single cell cDNA.
4. Enzymatic workflow steps for **DNA damage repair & nuclease treatment** of Kinnex single-cell RNA SMRTbell libraries.
5. Workflow steps for **final cleanup of Kinnex single-cell RNA SMRTbell libraries** using SMRTbell cleanup beads.
6. **Sample setup ABC¹** workflow steps to prepare Kinnex SMRTbell libraries for sequencing on Revio (+SPRQ) and Vega systems.

Preparing Kinnex™ libraries using Kinnex single-cell RNA kit



Procedure & checklist

Before you begin

This procedure describes the workflow for constructing single-cell Kinnex libraries from 10x Chromium 3' or 5' cDNA using the Kinnex single-cell RNA kit for library prep and sequencing on PacBio® Sequel® II, Sequel IIe, Vega™, and Revio® systems.

This kit is intended for use with single-cell cDNA generated using the 10x Chromium Next GEM Single Cell 3' kit (v3.1 or v4) or 10x Chromium Next GEM Single Cell 5' kit (v2, v3). For a full list of compatibility with 10x products, please consult the [10x website](#).

Overview	
Samples per kit	12
Workflow time	3 days for up to 12 samples
Number of SMRT® Cells per Kinnex library preparation	>8 SMRT Cells for Revio® using SPRQ™ chemistry >2 SMRT Cells for the Vega® or Revio® (non-SPRQ) systems >4 SMRT Cells for Sequel® II/IIe®

cDNA input	
Quantity	>15 ng 10x cDNA per sample cDNA concentration should be >1ng/μL with up to 15 μL in volume. See step 2.1 for 10x cDNA input requirement.
Average segment lengths	500–1,100 bp
Average 16-segment array lengths	10–17 kb


© 2025 PacBio. All rights reserved. Research use only. Not for use in diagnostic procedures.
103-254-300 REV08 OCT2025



PacBio [Documentation](#) (103-254-300)

Kinnex single-cell RNA library construction workflow overview


Procedure & checklist – Preparing Kinnex libraries using Kinnex single-cell RNA kit (103-254-300)



General best practices recommendations for preparing Kinnex single-cell RNA libraries

cDNA Input

- Use an optimal input range of **3,000 – 10,000 cells*** for the 10x Chromium single cell 3' cDNA generation workflow
- Follow the best practices in the **10x Chromium user guide**.
- Input cDNA quality control is highly recommended before proceeding to the Kinnex single-cell RNA workflow



Refer to [10x Genomics Support](#) website to download 10x Chromium user guides and other documentation.



DNA sizing and quantitation QC

- Perform DNA concentration measurements with a **Qubit fluorometer** using the Qubit 1X dsDNA High Sensitivity (HS) Assay Kit.
- Perform DNA sizing measurements with a **Bioanalyzer system** using the High Sensitivity DNA Kit (**for input cDNA QC**) or with a **Femto Pulse system** using the Genomic DNA 165 kb Kit (**for final SMRTbell library QC**)



Qubit 4 fluorometer and 1X ds DNA High Sensitivity Assay Kit
(Thermo Fisher Scientific)



Bioanalyzer 2100 System
and **High Sensitivity DNA Kit**
(Agilent Technologies)



Femto Pulse System
and **Genomic DNA 165 Kit**
(Agilent Technologies)

General best practices recommendations for preparing Kinnex single-cell RNA libraries (cont.)

Reagent and sample handling

Kinnex library prep reagents

- Room temperature is defined as any temperature in the range of 18–23°C for this protocol.
- Thaw the repair buffer, nuclease buffer, Kinnex ligase buffer, Kinnex ligation additive, Kinnex adapter mix, and primers at room temperature.
- Briefly vortex reagent buffers and Kinnex adapters prior to use. Enzyme mixes do not require vortexing.
- Quick-spin all reagents in a microcentrifuge to collect liquid at tube bottom prior to use.
- Keep all temperature-sensitive reagents on ice.
- **Ensure that the DNA damage repair mix is stored at -20°C to avoid poor library performance.**
- **Note:** The Loading buffer is **light sensitive** and should be protected from light when not in use.
- Bring SMRTbell cleanup beads and Qubit 1X dsDNA HS reagents to room temperature for 30 minutes prior to use.
- Pipette-mix all bead binding and elution steps until beads are distributed evenly in solution.
- Wide-bore pipette tips help to minimize foaming specifically when resuspending Kinnex capture beads.
- Mix all library prep reactions by pipetting up and down 10 times or until fully resuspended.
- Samples can be stored at 4°C at all safe stopping points listed in the protocol.
- 1.5X SMRTbell cleanup is recommended before Kinnex array formation. If the cDNA contains smaller fragments <200 bp, it is recommended to increase the SMRTbell cleanup ratio to 2X.

Temperature-sensitive reagents:

Temperature-sensitive reagents			
Step used	Tube color	Reagent	
TSO PCR and Kinnex PCR	Green	Kinnex single-cell PCR mix 103-244-500	
	Yellow	Kinnex 3' capture primer mix 103-182-400	
	Red	Kinnex 5' capture primer mix 103-182-200	
	Orange		Kinnex primers premix (A-PQ)
			103-107-800 A
			103-107-900 B
			103-108-000 C
			103-108-100 D
			103-108-200 E
			103-108-300 F
			103-108-400 G
			103-153-000 H
			103-153-100 I
			103-153-200 J
			103-153-300 K
		103-153-400 L	
		103-153-500 M	
	103-153-600 N		
	103-153-700 O		
	103-153-800 PQ		
Kinnex array formation	Light green	Kinnex single-cell enzyme 103-243-800	
	Yellow	Kinnex single-cell ligase 103-244-000	
	White	Kinnex single-cell ligase buffer 103-244-100	
	Red	Kinnex single-cell ligation additive 103-244-400	
	Blue	Kinnex adapter mix bc01 103-109-600 bc02 103-109-700 bc03 103-109-800 bc04 103-109-900	
DNA damage repair Nuclease treatment	Green	DNA repair mix 103-110-000	
	Purple	Repair buffer 102-244-300	
	Light green	Nuclease mix 103-110-100	
	Light purple	Nuclease buffer 103-110-200	

General best practices recommendations for preparing Kinnex single-cell RNA libraries (cont.)

Reagent and sample handling

Sequencing prep reagents (for sample setup ABC)

- Once thawed, place reaction buffers and sequencing primer on-ice prior to making master mix. The Loading buffer should be left at room-temperature.
- **Note:** The Loading buffer is **light sensitive** and should be protected from light when not in use.
- **Keep the following reagents on a cold block or ice:**
 - Sequencing polymerase
 - Sequencing control
- **Bring the following reagents up to room temperature 30 minutes prior to use:**
 - Loading buffer
 - SMRTbell cleanup beads

Thaw the following reagents and room temperature:

Component	Tube color
Annealing buffer	Light blue
Kinnex sequencing primer	Light green
Polymerase buffer	Yellow
Loading buffer	Green
Dilution buffer	Blue

Input cDNA quality control

Input cDNA quality control is highly recommended before proceeding to the Kinnex single-cell RNA library prep workflow



- We recommend using an optimal range of **3,000 – 10,000 cells input** into 10x Chromium 3' or 5' single cell workflow¹
- Protocol requires a **minimum of 15 ng** of 10x Chromium single cell cDNA (**maximum of 75 ng** per library)
 - If your cDNA sample amounts are between **16 – 59 ng**, then **normalize** all samples to 15 ng
 - If your cDNA sample amounts are **>75 ng**, then **normalize** all samples to 75 ng
 - If your cDNA sample amounts are between **60 – 75 ng**, **normalization is not required**.
- Evaluate the size distribution of each input cDNA sample to determine whether it is suitable for the protocol (average cDNA fragment size should be between **500 – 1,500 bp**)
 - 10x single cell cDNA samples measured with a Bioanalyzer system typically show a peak at **~1 – 1.8 kb**

Step 1
procedural summary



DNA
sizing QC



Single-cell cDNA size
ideally >500 bp

Proceed to **Step 2**

✓	Step	Instructions
	1.1	Bring the Qubit 1X dsDNA HS working solution and standards to room temperature.
	1.2	Pulse vortex or pipette mix each sample to homogenize the DNA in solution.
	1.3	Quick spin each sample to collect liquid.
	1.4	Take a 1 µL aliquot from each sample.
	1.5	Measure DNA concentration with a Qubit fluorometer using the 1X dsDNA HS kit.
	1.6	Dilute each sample to 1.0-1.5 ng/µL in elution buffer or water, based on the Qubit reading.
	1.7	Measure DNA size distribution with a Bioanalyzer system using the High Sensitivity DNA Kit.
	1.8	Proceed to the next step of the protocol if sample quality is acceptable.

Input cDNA quality control (cont.)

Example Bioanalyzer DNA sizing QC results for single cell 3' cDNA prepared with the 10x Chromium system





Bioanalyzer DNA sizing QC analysis results for a 10x Chromium single cell 3' cDNA sample prepared from a human GM12878 cell line.

TSO PCR & TSO artifact removal


Perform Steps 2 & 3 to remove template switch oligo (TSO) priming artifacts generated during 10x cDNA synthesis

TSO priming artifacts can occur if the TSO acts as a nonspecific primer on poly(A)+ mRNA




- Up to ~50% of cDNA products from the 10x Chromium single cell cDNA preparation workflow may contain a TSO priming artifact instead of the correct structure

Example synthesis of 10x 3' cDNA products containing correct structure



Amplified (full-length) 10x cDNA products with correct structure (TSO sequence at one end)


Example synthesis of 10x 3' cDNA products containing TSO priming artifact




Amplified (non-full length) 10x cDNA products with TSO priming artifacts (TSO sequences at both ends)

TSO PCR & TSO artifact removal (cont.)


Role of Kinnex capture primers and Kinnex capture beads



TSO PCR step in Kinnex single-cell RNA procedure uses a modified PCR primer (**Kinnex capture primer Fwd**) to incorporate a biotin tag into desired cDNA products with the correct structure




TSO artifact removal step uses streptavidin-coated **Kinnex capture beads** to pull down and enrich for biotin-tagged cDNA products with the correct structure



TSO PCR

Procedural summary



Step 2 procedural summary

2.1 – 2.4 TSO PCR

- Normalize cDNA sample input to **15 ng** if it is **between 15 ng and 59 ng** using elution buffer
- Normalize cDNA sample input to **75 ng** if it is **higher than 75 ng** using elution buffer
- For cDNA amounts **between 60–75ng**, **proceed without normalizing**



TSO PCR program for 15 ng cDNA input

Step	Time	Temperature	Cycles
1	3 min	98°C	1
2	20 sec	98°C	
3	30 sec	65°C	5
4	4 min	72°C	
5	5 min	72°C	1
6	Hold	4°C	1

TSO PCR program for 60-75 ng cDNA input

Step	Time	Temperature	Cycles
1	3 min	98°C	1
2	20 sec	98°C	
3	30 sec	65°C	3
4	4 min	72°C	
5	5 min	72°C	1
6	Hold	4°C	1

2.5 – 2.20 1.5X SMRTbell cleanup



Proceed to the Step 3 of the protocol if sample quantity is acceptable (between 100 ng – 600 ng). If the yield is higher than 600 ng, proceed to TSO artifact removal step with 600 ng of amplified biotinylated DNA.

Proceed to **Step 3**

TSO PCR

Procedural notes



2. TSO PCR

Step	Instructions																				
2.1	<p>Normalize cDNA sample input to 15 ng if it is between 15 ng and 59 ng using elution buffer. Normalize cDNA sample input to 75 ng if it is higher than 75 ng using elution buffer.</p> <p>For cDNA amounts between 60–75ng, proceed without normalizing. Select either the Kinnex 3' or 5' capture primer mix depending on the 10x Genomics kit used. Set up the following PCR reaction on ice (RM1).</p> <p>Reaction Mix 1 (RM1):</p> <table border="1"> <thead> <tr> <th>Tube color</th> <th>Component</th> <th>Volume</th> </tr> </thead> <tbody> <tr> <td></td> <td>Nuclease-free water</td> <td>Up to 50 µL</td> </tr> <tr> <td>Green</td> <td>Kinnex single-cell PCR mix</td> <td>25 µL</td> </tr> <tr> <td>Red</td> <td>Kinnex 5' capture primer mix</td> <td rowspan="2">10 µL</td> </tr> <tr> <td>Yellow</td> <td>Kinnex 3' capture primer mix</td> </tr> <tr> <td></td> <td>10x 5' or 3' cDNA library (1–5 ng/µL)</td> <td>Up to 15 µL</td> </tr> <tr> <td></td> <td>Total volume</td> <td>50 µL</td> </tr> </tbody> </table>	Tube color	Component	Volume		Nuclease-free water	Up to 50 µL	Green	Kinnex single-cell PCR mix	25 µL	Red	Kinnex 5' capture primer mix	10 µL	Yellow	Kinnex 3' capture primer mix		10x 5' or 3' cDNA library (1–5 ng/µL)	Up to 15 µL		Total volume	50 µL
Tube color	Component	Volume																			
	Nuclease-free water	Up to 50 µL																			
Green	Kinnex single-cell PCR mix	25 µL																			
Red	Kinnex 5' capture primer mix	10 µL																			
Yellow	Kinnex 3' capture primer mix																				
	10x 5' or 3' cDNA library (1–5 ng/µL)	Up to 15 µL																			
	Total volume	50 µL																			
2.2	Pipette-mix RM1.																				
2.3	Quick spin RM1 in a microcentrifuge to collect liquid.																				
2.4	Select the TSO PCR program based on cDNA input. Keep sample on ice until thermal cycler lid has heated to 105°C.																				

TSO PCR program (15 ng input)

Heated lid set at 105°C

Step	Time	Temperature	Cycles
1	3 min	98°C	1
2	20 sec	98°C	1
3	30 sec	65°C	5
4	4 min	72°C	1
5	5 min	72°C	1
6	Hold	4°C	1

Or TSO PCR program (60-75 ng input)

Heated lid set at 105°C

Step	Time	Temperature	Cycles
1	3 min	98°C	1
2	20 sec	98°C	1
3	30 sec	65°C	3
4	4 min	72°C	1
5	5 min	72°C	1
6	Hold	4°C	1

• If needed, **normalize** cDNA sample input amounts to 15 ng or 75 ng

• **IMPORTANT!** Select either the Kinnex 3' or 5' capture primer mix depending on the 10x Genomics kit used

• Set up TSO PCR reactions **ON ICE**

• PCR polymerase 3'→5' exonuclease activity negatively impacts amplification yield if prepared at room temp.

• For **lower** cDNA sample inputs (15 ng), use a **higher** number of PCR cycles (5)

After completing **TSO PCR** step, perform cleanup with **1.5X SMRTbell cleanup beads** and proceed to **TSO artifact removal** (Step 3) if sample quantity is acceptable (**between 100 ng – 600 ng**)

• If the yield is higher than 600 ng, proceed to the next step with 600 ng of amplified biotinylated DNA.

TSO artifact removal

Procedural summary



¹ Scale up the volume of beads if processing more than 4 samples (with 10% overage). If preparing more than 40 μ L of beads, use a 1.5 mL LoBind tube instead of PCR tube.

² Scale up the volume of Kinnex capture binding buffer accordingly if preparing more than 40 μ L of beads.


Step 3 procedural summary

3.1 – 3.21 TSO artifact removal

3.1 – 3.4 Resuspend Kinnex capture beads by vortexing and transfer 10 μ L of beads per sample to a PCR tube




3.5 – 3.7 Resuspend Kinnex capture beads in Kinnex capture binding buffer (Perform resuspension 2 times)



3.8 – 3.10 Add purified TSO PCR sample (containing biotinylated DNA-fragments) from Step 2 to resuspended beads and incubate on rotator




3.11 – 3.14 Resuspend Kinnex capture beads/DNA-complex in Kinnex capture washing buffer (Perform wash 2 times)



TSO artifact removal

Procedural summary



¹ Scale up the volume of beads if processing more than 4 samples (with 10% overage). If preparing more than 40 μL of beads, use a 1.5 mL LoBind tube instead of PCR tube.

² Scale up the volume of Kinnex capture binding buffer accordingly if preparing more than 40 μL of beads.


Step 3 procedural summary (cont.)

3.1 – 3.21 TSO artifact removal (cont.)

3.15 – 3.16 Resuspend Kinnex capture beads/DNA-complex in nuclease-free water



3.17 – 3.20 Resuspend Kinnex capture beads/DNA-complex in elution buffer, add 2 μL Kinnex single-cell enzyme and run TSO artifact removal program



3.21 Place sample tube on magnet for 1 minute and move supernatant containing library to a fresh tube



3.22 – 3.38 1.5X SMRTbell cleanup



Proceed to **Step 4**

TSO artifact removal

Procedural notes



3. TSO artifact removal

Step	Instructions									
3.1	Bring Kinnex capture beads kit to room temperature. Resuspend the beads by vortexing.									
3.2	Transfer 10 μ L resuspended Kinnex capture beads per sample to a PCR tube. Scale up the amount of beads if processing more than 4 samples (with 10% overage). If preparing more than 40 μ L of beads, use a 1.5 mL LoBind tube instead of PCR tube.									
3.3	Place the tube on the magnet until the beads separate fully from the solution.									
3.4	Carefully remove and discard the supernatant while the tube remains on the magnet. Avoid touching the bead pellet with the pipette tip.									
3.5	<ul style="list-style-type: none">Remove the tube from the magnet.Add 40 μL Kinnex bead binding buffer along the inside wall of the tube where the beads are collected and gently resuspend by pipetting using wide bore tips. DO NOT VORTEX. <p>Note: the solution may be viscous. Highly recommend using wide bore tips to avoid foaming. When excess bubbles are present, lower cDNA recovery is expected.</p> <ul style="list-style-type: none">Quick-spin the tube in a microcentrifuge if needed. <p>Note: Scale up the volume of Kinnex capture binding buffer accordingly if preparing more than 40 μL of beads.</p>									
3.6	Place the tube on the magnet until the beads separate fully from the solution and remove the supernatant.									
3.18	Add 2 μ L Kinnex enzyme to the sample with capture beads to cleave the captured DNA products from Kinnex capture beads.									
3.19	Pipette-mix each sample and a very quick spin in a microcentrifuge to collect liquid.									
Run the TSO artifact removal program .										
TSO artifact removal program										
Heated lid set at $\geq 47^{\circ}\text{C}$										
3.20	<table border="1"><thead><tr><th>Step</th><th>Time</th><th>Temperature</th></tr></thead><tbody><tr><td>1</td><td>30 min</td><td>37°C</td></tr><tr><td>2</td><td>Hold</td><td>4°C</td></tr></tbody></table>	Step	Time	Temperature	1	30 min	37°C	2	Hold	4°C
Step	Time	Temperature								
1	30 min	37°C								
2	Hold	4°C								
3.21	Place the tube on the magnet for 1 minute and move the supernatant containing the library to a fresh tube.									

- Bring Kinnex capture beads to **room temperature** and resuspend by vortexing

- Critical step! For all Kinnex capture bead handling steps:** Pipette mix with care and avoid generating bubbles by using **wide bore tips** for mixing (**do not vortex**)
 - When excess bubbles are present, lower cDNA recovery is expected


- Add Kinnex enzyme to **cleave** captured cDNA products from Kinnex capture beads

- Keep the **supernatant** after treatment with Kinnex enzyme and placement on the magnet

After completing **TSO artifact removal** step, perform cleanup with **1.5X SMRTbell cleanup beads** and proceed to **Kinnex PCR** (Step 4) if sample quantity is acceptable (**minimum 25 ng**)


TSO PCR & TSO artifact removal (cont.)

TSO artifact video demonstration




Kinnex PCR

In this step, incorporate programmable segmentation adapter sequences into amplified cDNA products



Amplified cDNA products from PCR 3 contain terminal segmentation adapter sequences that are **complementary** to the ends of cDNA products from PCR 2 & PCR 4

Set up 16 parallel PCR reactions/sample with premixed Kinnex primers to generate amplified cDNA products containing programmable sequences at both ends.




Kinnex PCR (cont.)

Step 4 procedural summary

Procedural summary




3.1 Prepare Kinnex primers premix & 3.2 Kinnex PCR



Kinnex PCR Rx tube	12X concatenation	Tube color
1 [25 μ L Rx volume]	Kinnex primer mix A	Orange
2 [25 μ L Rx volume]	Kinnex primer mix B	
3 [25 μ L Rx volume]	Kinnex primer mix C	
4 [25 μ L Rx volume]	Kinnex primer mix D	
5 [25 μ L Rx volume]	Kinnex primer mix E	
6 [25 μ L Rx volume]	Kinnex primer mix F	
7 [25 μ L Rx volume]	Kinnex primer mix G	
...	...	
13 [25 μ L Rx volume]	Kinnex primer mix M	
14 [25 μ L Rx volume]	Kinnex primer mix N	
15 [25 μ L Rx volume]	Kinnex primer mix O	
16 [25 μ L Rx volume]	Kinnex primer mix PQ	



3.2 Pooling of 16 Kinnex PCR products and 1.1X SMRTbell cleanup



Proceed to **Step 5**

Kinnex PCR (cont.)

Procedural notes



4. Kinnex PCR

Step	Instructions										
	Set up the following PCR reaction mix per sample on ice. Reaction Mix 2 (RM2):										
4.1	<table border="1"> <thead> <tr> <th>Master mix components</th> <th>Volume for 16X concatenation*</th> </tr> </thead> <tbody> <tr> <td>PCR-grade water</td> <td>176-X μL</td> </tr> <tr> <td>Kinnex single-cell PCR mix</td> <td>220 μL</td> </tr> <tr> <td>25 ng of amplified cDNA from Step 3.38</td> <td>X μL</td> </tr> <tr> <td>Total volume</td> <td>396 μL</td> </tr> </tbody> </table> <p>X = 25 (ng) / purified cDNA concentrations from Step 3.38, if the cDNA yield from step 3.38 is greater than 25 ng, proceed with only 25 ng. *10% overage included</p>	Master mix components	Volume for 16X concatenation*	PCR-grade water	176-X μ L	Kinnex single-cell PCR mix	220 μ L	25 ng of amplified cDNA from Step 3.38	X μ L	Total volume	396 μL
Master mix components	Volume for 16X concatenation*										
PCR-grade water	176-X μ L										
Kinnex single-cell PCR mix	220 μ L										
25 ng of amplified cDNA from Step 3.38	X μ L										
Total volume	396 μL										
4.2	Quick-spin RM2 in a microcentrifuge to collect liquid.										
4.3	Add 22.5 μ L of RM2 to a new PCR tube on ice. Repeat this step to prepare a total of 16 tubes per sample (each containing 22.5 μ L of RM2). Add 2.5 μ L of Kinnex primers premix into each of 16 PCR tubes on ice according to the table below.										
4.4	<table border="1"> <thead> <tr> <th>PCR tube</th> <th>Kinnex primers premix</th> </tr> </thead> <tbody> <tr> <td>1</td> <td>Kinnex primers premix A</td> </tr> <tr> <td>2</td> <td>Kinnex primers premix B</td> </tr> <tr> <td>16</td> <td>Kinnex primers premix PQ</td> </tr> </tbody> </table>	PCR tube	Kinnex primers premix	1	Kinnex primers premix A	2	Kinnex primers premix B	16	Kinnex primers premix PQ		
PCR tube	Kinnex primers premix										
1	Kinnex primers premix A										
2	Kinnex primers premix B										
16	Kinnex primers premix PQ										
4.5	Pipette-mix each sample. The total volume of each tube should be 25.0 μ L.										
4.6	Quick-spin the PCR tubes in a microcentrifuge to collect liquid. Run the Kinnex PCR program with heated lid set to 105°C. Keep sample on ice until thermal cycler lid has heated to 105°C.										

- Set up Kinnex PCR reactions **ON ICE**
- PCR polymerase 3'→5' exonuclease activity negatively impacts amplification yield if prepared at room temp.

- **CRITICAL STEP!** Correct setup of all 16 Kinnex PCR reactions is required – any missing/incorrect Kinnex primer pairs will result in no/low SMRTbell yield

Kinnex PCR program

Heated lid set at 105°C


Step	Time	Temperature	Cycles
1	3 min	98°C	1
2	20 sec	98°C	
3	30 sec	68°C	9
4	4 min	72°C	
5	5 min	72°C	1
6	Hold	4°C	1

- After **Kinnex PCR** step, pool entire volume of all 16 reactions into a clean 1.5 mL LoBind tube¹ and perform cleanup with **1.5X SMRTbell cleanup beads**


¹ DNA concentration QC measurements for individual PCR reactions prior to pooling are not required.

Kinnex array formation

In this step, assemble cDNA transcripts (“segments”) containing programmable ends into a linear array



Treat pooled Kinnex PCR products with Kinnex enzyme to create single-stranded overhangs to enable subsequent directional assembly of cDNA transcripts into a linear, ordered array



Barcoded Kinnex terminal adapters¹ are ligated to specific overhang sequences at array ends

Kinnex enzyme produces a single nucleotide gap at the location of a uracil (dU) residue within segmentation adapter sequences and enables generation of single-stranded overhangs


Add Kinnex ligase to enable directional assembly of cDNA transcript segments containing “sticky ends” into a linear, ordered array

Kinnex library prep workflow supports incorporation of PacBio barcodes at the array formation step to enable up to 4-plex sample multiplexing at the SMRTbell library level

Kinnex array formation (cont.)

Kinnex terminal adapters incorporate barcode sequences to enable up to 4-plex sample multiplexing at the library level


- Kinnex adapters contain **barcode sequences¹** to enable (optional) sample multiplexing at the SMRTbell library level (up to 4-plex)
 - Forward and reverse Kinnex adapter pairs are pre-mixed in Kinnex concatenation kits
 - Kinnex concatenation kits contain a total of 4 **barcoded Kinnex adapter mixes (bcM0001-bcM0004)** to enable multiplexing of up to 4 samples per SMRT Cell



Kinnex array formation (cont.)

Kinnex terminal adapters use a new design that enables improved HiFi sequencing performance

- Kinnex adapters enable longer polymerase read length performance and improved sample loading efficiency
 - Improved overall HiFi data yield




- Kinnex adapter design requires a **different sequencing primer (Kinnex sequencing primer¹)** for annealing reaction




Kinnex array formation (cont.)

Procedural summary



Step 5 procedural summary

5.1 – 5.12 Kinnex array formation




5.13 – 5.27 1.2X SMRTbell cleanup



Proceed to **Step 6**

Kinnex array formation (cont.)

Procedural notes



5. Kinnex array formation

Step	Instructions												
5.1	In a 0.2 mL PCR tube, add 5 µg of sample from Step 4.22, in 47 µL of volume (106 ng/µl). Dilute with elution buffer going into this step if sample is too concentrated.												
5.2	Add 10 µL of Kinnex enzyme to create single-stranded extensions on PCR-amplified cDNA fragments to enable subsequent directional assembly of 16 PCR products.												
5.3	Pipette-mix each sample.												
5.4	Run the Kinnex primer digestion program .												
5.5	Add 3 µL of Kinnex adapter barcode 01–04 (use a single barcode per sample) and 20 µL of Kinnex ligation additive to each sample for a total volume of 80 µL.												
	<table border="1"><thead><tr><th>Tube color</th><th>Component</th><th>Volume</th></tr></thead><tbody><tr><td>Blue</td><td>Kinnex barcode adapter</td><td>3 µL</td></tr><tr><td>Red</td><td>Kinnex ligation additive</td><td>20 µL</td></tr><tr><td colspan="2">Total volume</td><td>23 µL</td></tr></tbody></table>	Tube color	Component	Volume	Blue	Kinnex barcode adapter	3 µL	Red	Kinnex ligation additive	20 µL	Total volume		23 µL
Tube color	Component	Volume											
Blue	Kinnex barcode adapter	3 µL											
Red	Kinnex ligation additive	20 µL											
Total volume		23 µL											
5.6	Pipette-mix each sample. Add the following components in the order and volume listed below to a new microcentrifuge tube. Adjust component volumes for the number of samples being prepared, plus 10% overage. For individual preps, add components directly to each sample in the order and volume listed below. Reaction Mix 3 (RM3):												
5.7	<table border="1"><thead><tr><th>Tube color</th><th>Component</th><th>Volume</th></tr></thead><tbody><tr><td>White</td><td>Kinnex single-cell ligase buffer</td><td>10 µL</td></tr><tr><td>Yellow</td><td>Kinnex single-cell ligase</td><td>10 µL</td></tr><tr><td colspan="2">Total volume</td><td>20 µL</td></tr></tbody></table>	Tube color	Component	Volume	White	Kinnex single-cell ligase buffer	10 µL	Yellow	Kinnex single-cell ligase	10 µL	Total volume		20 µL
Tube color	Component	Volume											
White	Kinnex single-cell ligase buffer	10 µL											
Yellow	Kinnex single-cell ligase	10 µL											
Total volume		20 µL											
5.12	Run the Kinnex array ligation program .												

- Recommended input amount to proceed with Kinnex array formation is **5 µg** of Kinnex PCR amplicons (from Step 4)
- Proceeding with <3 µg is **not recommended** since lower input amounts may lead to insufficient final library yields to enable optimal sequencing results

Kinnex primer digestion program

Heated lid set at 47°C

Step	Time	Temperature
1	30 min	37°C
2	Hold	4°C

- IMPORTANT:** Always add Kinnex adapters and Kinnex ligation additive to the sample **BEFORE** adding Kinnex single-cell ligase
- If combining multiple barcoded Kinnex libraries for sequencing, make sure each library uses one of the 4 different Kinnex barcoded adapters

Kinnex array ligation program



Heated lid set at 52°C

Step	Time	Temperature
1	60 min	42°C
2	Hold	4°C

- After [Kinnex array formation](#) step, perform cleanup with 1.2X SMRTbell cleanup beads and proceed to [DNA damage repair](#) (Step 6)

DNA damage repair & nuclease treatment

Perform DNA damage repair to repair nicked/damaged DNA sites and perform nuclease treatment to remove incomplete SMRTbell template arrays



• After nuclease treatment, most remaining SMRTbell templates are complete (full-length) 16-segment transcript array molecules capped with Kinnex adapters

DNA damage repair

Procedural summary



Step 6 procedural summary

6.1 – 6.8 DNA damage repair




6.9 – 6.22 1.2X SMRTbell cleanup



Proceed to **Step 7**

DNA damage repair

Procedural notes



6. DNA damage repair

✓ Step	Instructions																
6.1	In a 0.2 mL PCR tube, add 42 μ L of sample from Step 5.26																
	Prepare this reaction mix fresh before use. Place the DNA damage repair mix on ice at all times, and immediately return the DNA damage repair mix back to the freezer (-20°C) after use. Improper storage and handling of the DNA damage repair mix may result in poor library performance and should not be used for subsequent reactions.																
	Add the following components in the order and volume listed below to a new microcentrifuge tube. Adjust component volumes for the number of samples being prepared, plus 10% overage. For individual preps, add components directly to each sample in the order and volume listed below.																
6.2	Reaction Mix 4 (RM4):																
	<table border="1"><thead><tr><th>✓</th><th>Tube color</th><th>Component</th><th>Volume</th></tr></thead><tbody><tr><td></td><td>Purple</td><td>Repair buffer</td><td>6 μL</td></tr><tr><td></td><td>Green</td><td>DNA repair mix</td><td>2 μL</td></tr><tr><td></td><td></td><td>Total volume</td><td>8 μL</td></tr></tbody></table>	✓	Tube color	Component	Volume		Purple	Repair buffer	6 μ L		Green	DNA repair mix	2 μ L			Total volume	8 μL
✓	Tube color	Component	Volume														
	Purple	Repair buffer	6 μ L														
	Green	DNA repair mix	2 μ L														
		Total volume	8 μL														
6.3	Pipette-mix RM4.																
6.4	Quick-spin RM4 in a microcentrifuge to collect liquid.																
6.5	Add 8 μ L of RM4 to each sample. Total volume should equal 50 μ L.																
6.6	Pipette-mix each sample.																
6.7	Quick-spin the strip tube in a microcentrifuge to collect liquid.																
6.8	Run the DNA repair program .																

- **DO NOT** proceed with more than 5 μ g of Kinnex array input per individual DDR reaction (in a 50 μ L reaction volume)
- Adding more than 5 μ g per reaction can result in incomplete digestion of partial (non-full length) SMRTbell template arrays, resulting in low sample loading performance during sequencing

DNA damage repair program


Heated lid set at 47°C

Step	Time	Temperature
1	30 min	37°C
2	Hold	4°C

- After **DNA damage repair** step, perform cleanup with 1.2X SMRTbell cleanup beads and proceed to **Nuclease treatment** (Step 7)


Nuclease treatment and cleanup

Procedural summary



Step 7 procedural summary

7.1 – 7.7 Nuclease treatment



7.8 – 7.24 1.2X SMRTbell cleanup



Proceed to **ABC¹**

Nuclease treatment and cleanup (cont.)

Procedural notes



7.1 – 7.7 Nuclease treatment

Step	Instructions												
7.1	Add the following components in the order and volume listed below to a new microcentrifuge tube. Adjust component volumes for the number of samples being prepared, plus 10% overage. For individual preps, add components directly to each sample from the previous step in the order and volume listed below. Reaction Mix 5 (RM5):												
	<table border="1"><thead><tr><th>Tube color</th><th>Component</th><th>Volume</th></tr></thead><tbody><tr><td>Light purple</td><td>Nuclease buffer 103-110-200</td><td>5 µL</td></tr><tr><td>Light green</td><td>Nuclease mix 103-110-100</td><td>5 µL</td></tr><tr><td colspan="2">Total volume</td><td>10 µL</td></tr></tbody></table>	Tube color	Component	Volume	Light purple	Nuclease buffer 103-110-200	5 µL	Light green	Nuclease mix 103-110-100	5 µL	Total volume		10 µL
Tube color	Component	Volume											
Light purple	Nuclease buffer 103-110-200	5 µL											
Light green	Nuclease mix 103-110-100	5 µL											
Total volume		10 µL											
7.2	Pipette-mix RM5 .												
7.3	Quick-spin RM5 in a microcentrifuge to collect liquid.												
7.4	Add 10 µL of RM5 to each sample. Total volume should equal 50 µL.												
7.5	Pipette-mix each sample.												
7.6	Quick-spin the strip tube in a microcentrifuge to collect liquid.												
	Run the nuclease treatment program .												

• 1-hr nuclease treatment

Nuclease treatment program


Heated lid set at 47°C

Step	Time	Temperature
1	60 min	37°C
2	Hold	4°C

• After **Nuclease treatment** step, perform final cleanup with **1.2X SMRTbell cleanup beads** (Step 7.8)

Nuclease treatment and cleanup (cont.)

Procedural notes



7.8 – 7.24 Final cleanup with 1.2X SMRTbell cleanup beads

Step	Instructions
7.8	Add 60 μL (1.2X v/v) SMRTbell cleanup beads to each sample from the previous step. Pipette-mix the beads until evenly distributed.
7.9	Quick-spin the tube strip in a microcentrifuge to collect all liquid.
7.10	Incubate at room temperature for 10 minutes to allow DNA to bind beads.
7.11	Place the tube strip in a magnetic separation rack until beads separate fully from the solution.
7.12	Slowly remove the cleared supernatant without disturbing the beads. It is recommended to save the supernatant in another tube strip in case of poor DNA recovery.
7.13	Slowly dispense 200 μL , or enough to cover the beads, of freshly prepared 80% ethanol into each tube. After 30 seconds, remove the 80% ethanol and discard.
7.14	Repeat the previous step.
Remove residual 80% ethanol:	
	<ul style="list-style-type: none"> Remove tube strip from the magnetic separation rack. Quick spin tube strip in a microcentrifuge. Place tube strip back in a magnetic separation rack until beads separate fully from the solution.
7.15	Remove residual 80% ethanol and discard.



7.21	Take a 1 μL aliquot from each tube. Measure DNA concentration with a Qubit fluorometer using the 1x dsDNA HS kit. Calculate the total mass. Yield typically ranges from 500–1200 ng. Recommended: Further dilute each aliquot to 250 $\text{pg}/\mu\text{L}$ with Femto Pulse dilution buffer. Measure final SMRTbell library size distribution with a Femto Pulse system.
7.22	Normalize 25 μL of library to 20 $\text{ng}/\mu\text{L}$ prior to ABC. A concentration of less than 20 $\text{ng}/\mu\text{L}$ is acceptable if sample quantity is limited. Using a concentration above 20 $\text{ng}/\mu\text{L}$ will result in lower loading during sequencing.
7.23	Proceed to Section 8 to prepare library for sequencing on Vega or Revio (+SPRQ chemistry). Proceed to SMRT Link Sample Setup for preparing samples for sequencing on Revio (non-SPRQ) or Sequel II/e.
7.24	Store SMRTbell libraries at 4°C if sequencing within the week. Long-term storage should be at -20°C. Minimize freeze-thaw cycles when handling SMRTbell libraries.

SAFE STOPPING POINT - Store at 4°C



- Perform **DNA concentration QC** on final purified Kinnex RNA library using a Qubit dsDNA HS assay
 - Typical final SMRTbell library yield from 5 μg of input DNA into Kinnex array formation is **~10 – 25%** – a much higher observed yield might suggest incomplete digestion of partial SMRTbell templates
 - Troubleshooting tip:** If SMRTbell library yield is higher than expected and sample loading is lower than expected, consider repeating the nuclease treatment step






- Perform **DNA sizing QC** on final purified Kinnex single-cell RNA library using a Femto Pulse system (expected final library insert size is **~10 – 15 kb**)

- Normalize 25 μL of library to 20 $\text{ng}/\mu\text{L}$ prior to ABC¹. A concentration of less than 20 $\text{ng}/\mu\text{L}$ is acceptable if sample quantity is limited.
 - Using a concentration above 20 $\text{ng}/\mu\text{L}$ will result in lower loading during sequencing



Nuclease treatment and cleanup (cont.)

Example Femto Pulse DNA sizing QC results for Kinnex single-cell RNA library before nuclease treatment and after nuclease treatment & final cleanup



Kinnex single-cell RNA library prep inputs & expected step yields

Final Kinnex library yield is typically sufficient to load ≥ 2 SMRT Cells



Example Femto pulse DNA sizing QC analysis results for final Kinnex single-cell library generated for a 10x Chromium single-cell 3' cDNA sample prepared from a human cell line (HG002).

Example Kinnex single-cell RNA library prep yields

10x single cell 3' cDNA input	15 ng
cDNA input for Kinnex array formation	5000 ng
Post-nuclease treatment & final library cleanup yield (%) ¹	1100 ng (22.0%)

¹ Post-nuclease treatment & final cleanup yields typically ranged from ~10% to ~25% when using single-cell 3' cDNA samples for Kinnex single-cell RNA library construction.

Final Kinnex library yield is typically sufficient to load the following number of SMRT Cells:


- >8 SMRT Cells for Revio using SPRQ chemistry
- >2 SMRT Cells for the Vega or Revio (non-SPRQ) systems
- >4 SMRT Cells for Sequel II/IIe systems

Kinnex single-cell RNA sequencing preparation workflow details

Sample Setup & Run Design recommendations for Kinnex libraries

Follow SMRT Link run design instructions to prepare Kinnex libraries for sequencing

Select desired Kinnex application type in SMRT Link run design page



Runs

After specifying your application type, SMRT Link auto-fills selected Run Design parameter fields with default recommended values

Kinnex single-cell RNA¹

Whole genome sequencing

- Human WGS
- Microbial assembly
- Other WGS

RNA sequencing

- Iso-Seq method
- MAS-Seq single cell
- Kinnex single-cell RNA**
- Kinnex full-length RNA

Viral sequencing

- Adeno-associated virus

Metagenomics

- Kinnex 16S rRNA
- Full-length 16S rRNA sequencing

Kinnex full-length RNA²

Whole genome sequencing

- Human WGS
- Microbial assembly
- Other WGS

RNA sequencing

- Iso-Seq method
- MAS-Seq single cell
- Kinnex single-cell RNA
- Kinnex full-length RNA**

Viral sequencing

- Adeno-associated virus

Metagenomics

- Kinnex 16S rRNA
- Full-length 16S rRNA sequencing

Kinnex 16S rRNA²

Whole genome sequencing

- Human WGS
- Microbial assembly
- Other WGS

RNA sequencing

- Iso-Seq method
- MAS-Seq single cell
- Kinnex single-cell RNA
- Kinnex full-length RNA

Viral sequencing

- Adeno-associated virus

Metagenomics

- Kinnex 16S rRNA**
- Full-length 16S rRNA sequencing

SMRT Link module	Key setup parameters	Sequel II system recommended settings	Vega system recommended settings	Revio system recommended settings
Sample setup	Library type	Kinnex	Kinnex	Kinnex
	Primer	Kinnex sequencing primer	Kinnex sequencing primer	Kinnex sequencing primer
	Binding/Polymerase kit	Sequel II binding kit 3.2	Vega polymerase kit	Revio SPRQ polymerase kit
	Concentration on plate	40 – 60 pM	130 – 160 pM	130 – 160 pM
Runs → Run design	Use Adaptive Loading	YES	N/A	YES
	Library Type ¹	Kinnex	Kinnex	Kinnex
	Movie Acquisition Time	30 hrs	24 hrs	24 hrs
	Add Analysis	YES	YES	YES

SMRT Link Run Design procedure for **Revio and Vega systems**

Sample and run information

Kinnex single-cell RNA

Kinnex full-length RNA

Kinnex 16S rRNA

▼ Plate 1, Well A01: Kinnex single-cell RNA library demo

▼ Plate 1, Well A01: Kinnex full-length RNA library demo

▼ Plate 1, Well A01: Kinnex 16S rRNA library demo

Application <small>Required</small>	Kinnex single-cell RNA	Application <small>Required</small>	Kinnex full-length RNA	Application <small>Required</small>	Kinnex 16S rRNA
Plate Well <small>Required</small>	Plate 1, Well A01	Plate Well <small>Required</small>	Plate 1, Well A01	Plate Well <small>Required</small>	Plate 1, Well A01
Well Name <small>Required</small>		Well Name <small>Required</small>		Well Name <small>Required</small>	
Well Comment		Well Comment		Well Comment	
Library Type <small>Required</small>	Kinnex	Library Type <small>Required</small>	Kinnex	Library Type <small>Required</small>	Kinnex
Insert Size (bp) <small>Required</small>	15000	Insert Size (bp) <small>Required</small>	16000	Insert Size (bp) <small>Required</small>	17000
Library Concentration (pM) <small>Required</small>		Library Concentration (pM) <small>Required</small>		Library Concentration (pM) <small>Required</small>	
Movie Acquisition Time (hours)	24	Movie Acquisition Time (hours)	24	Movie Acquisition Time (hours)	24

Select application type to autofill Library Type & Movie Acquisition Time recommended settings

Default = Kinnex library type (instead of Standard or AAV)¹

Standard SMRTbell library type containing standard SPK3 barcoded terminal adapters



Forward and reverse standard terminal adapters have the same structure

Kinnex SMRTbell library type containing Kinnex barcoded terminal adapters



Library Type field determines which adapter finding algorithm is used during post-primary analysis¹

Forward and reverse Kinnex terminal adapters have different structures

¹ **Note:** When sequencing a Kinnex library sample, if 'Standard' library type is mistakenly selected instead of 'Kinnex' then a higher missing adapter rate (> 95%) and a slight degradation in barcode demultiplexing performance (~93-96% barcoded HiFi read yield) will be observed.

SMRT Link Run Design procedure for Revio and Vega systems (cont.)

Sample indexing (barcoding) information

Kinnex single-cell RNA

Kinnex full-length RNA

Kinnex 16S rRNA

Samples

Sample is indexed YES NO

Indexes Required: MAS adapter indexes

Biosample names Required: Interactively From a File

Samples

Sample is indexed YES NO

Indexes Required: MAS adapter indexes

Biosample names Required: Interactively From a File

Samples


Sample is indexed YES NO

Indexes Required: MAS adapter indexes

Biosample names Required: Interactively From a File

Default Indexes FASTA = MAS adapter indexes

Example complete Kinnex library molecule containing barcoded Kinnex terminal MAS adapters¹ at both ends



Example interactive biosample name specification

Barcode Selector and Sample Name Editor

Barcode ID	Bio Sample ID
<input type="checkbox"/> bcM0001--bcM0001	
<input type="checkbox"/> bcM0002--bcM0002	
<input type="checkbox"/> bcM0003--bcM0003	
<input type="checkbox"/> bcM0004--bcM0004	

Barcode Selector and Sample Name Editor

Barcode ID	Bio Sample ID
<input type="checkbox"/> bcM0003--bcM0003	
<input type="checkbox"/> bcM0004--bcM0004	
<input checked="" type="checkbox"/> bcM0001--bcM0001	Kinnex adapter-barcoded library 1
<input checked="" type="checkbox"/> bcM0002--bcM0002	Kinnex adapter-barcoded library 2

SMRT Link



Data Management

MAS SMRTbell barcoded adapter indexes

```
>bcM0001
ACAGTC
>bcM0002
ATGACG
>bcM0003
CACGTG
>bcM0004
CATCGC
```

¹ Four barcoded terminal Kinnex adapters (Kinnex adapter bcM0001-bcM0004) are available for Kinnex array formation step. Kinnex adapter barcode sequences can be downloaded from [SMRT Link](#) Data Management module.

SMRT Link Run Design procedure for **Revio and Vega systems (cont.)**

Data options

Kinnex single-cell RNA

Kinnex full-length RNA

Kinnex 16S rRNA

Default = **NO** for
Include Base Kinetics

▼ Data Options

Include Base Kinetics YES NO

Consensus Mode MOLECULE STRAND

Assign Data To Project

Default Consensus
Mode = **MOLECULE**¹

▼ Data Options

Include Base Kinetics YES NO

Consensus Mode MOLECULE STRAND

Assign Data To Project

Can leave **Include Base Kinetics** and **Consensus Mode**
fields at their **default settings** for Kinnex library samples

▼ Data Options

Include Base Kinetics YES NO

Consensus Mode MOLECULE STRAND

Assign Data To Project

SMRT Link Run Design analysis options

Kinnex single-cell RNA

Kinnex full-length RNA

Kinnex 16S rRNA

Analysis Options **Default = YES for Add Analysis**

Analysis Options

Analysis Options

Add Analysis YES NO

Analysis Name Required

Select Analysis Workflow Required

Segmentation Adapter Set

Primer Set Required

Reference Set Required

Kit Type 10X 3' KIT 10X 5' KIT

[Advanced Parameters](#)

Add Analysis YES NO

Analysis Name Required

Select Analysis Workflow Required

Segmentation Adapter Set

Primer Set Required

Reference Set

Cluster of Barcoded Samples

[Advanced Parameters](#)

Add Analysis YES NO

Analysis Name Required

Select Analysis Workflow Required

Segmentation Adapter Set

[Advanced Parameters](#)


Analysis Workflow is automatically filled in

Split arrayed HiFi reads at adapter positions, generating segmented reads (S-reads)


Read Segmentation and Single-Cell Iso-Seq

Read Segmentation and Iso-Seq

Read Segmentation



Isoform-classification analysis to identify novel genes & isoforms



Isoform-classification analysis to identify novel genes & isoforms with abundance information (raw counts & normalized counts per million)

SMRT Link Run Design analysis options (cont.)

Kinnex single-cell RNA

Analysis Options

Add Analysis YES NO

Analysis Name Required: Kinnex_Single-Cell_RNA_Demo_Analysis_Job_Name

Select Analysis Workflow Required: Read Segmentation and Single-Cell Iso-Seq

Segmentation Adapter Set: MAS-Seq Adapter v1 (MAS16)

Primer Set Required: 10x Chromium single cell 3' cDNA primers

Reference Set Required: Human Genome hg38, with Gencode v39 annotations

Kit Type 10X 3' KIT 10X 5' KIT

Advanced Parameters

Kinnex full-length RNA

Analysis Options

Add Analysis YES NO

Analysis Name Required: Kinnex_Full-Length_RNA_Demo_Analysis_Job_Name

Select Analysis Workflow Required: Read Segmentation and Iso-Seq

Segmentation Adapter Set: MAS-Seq Adapter v3 (MAS8)

Primer Set Required: Iso-Seq v2 Barcoded cDNA Primers

Reference Set: Human Genome hg38, with Gencode v39 annotations

Cluster of Barcoded Samples Pool reads and cluster together

Advanced Parameters

Kinnex 16S rRNA

Analysis Options

Add Analysis YES NO


Analysis Name Required: Kinnex_Full-Length_16S_rRNA_Demo_Analysis_Job_Name

Select Analysis Workflow Required: Read Segmentation

Segmentation Adapter Set: MAS-Seq Adapter v2 (MAS12)

Advanced Parameters

Specify Segmentation Adapter Set that corresponds to the Kinnex library concatenation method used



SMRT Link Run Design analysis options (cont.)

Kinnex single-cell RNA

Kinnex full-length RNA

Analysis Options

Add Analysis YES NO

Analysis Name Required: Kinnex_Single-Cell_RNA_Demo_Analysis_Job_Name

Select Analysis Workflow Required: Read Segmentation and Single-Cell Iso-Seq

Segmentation Adapter Set: MAS-Seq Adapter v1 (MAS16)

Primer Set Required: 10x Chromium single cell 3' cDNA primers

Reference Set Required: Human Genome hg38, with Gencode v39 annotations

Kit Type 10X 3' KIT 10X 5' KIT

Advanced Parameters

Analysis Options

Add Analysis YES NO

Analysis Name Required: Kinnex_Full-Length_RNA_Demo_Analysis_Job_Name

Analysis Workflow Required: Read Segmentation and Iso-Seq

Segmentation Adapter Set: MAS-Seq Adapter v3 (MAS8)

Primer Set Required: Iso-Seq v2 Barcoded cDNA Primers

Reference Set: Human Genome hg38, with Gencode v39 annotations

Cluster of Barcoded Samples Pool reads and cluster together

Advanced Parameters

Specify Primer Set used for cDNA amplification

Specify Reference Set. Default sets are:

- Human Genome hg38, with Gencode v39 annotations
- Mouse Genome mm39, with Gencode vM28 annotations

Specify reference genome & annotation sets to align high quality isoforms to, and to collapse isoforms mapped to the same genomic loci.



For Kinnex single-cell 3' RNA analysis, select '10x Chromium single cell 3' cDNA primers'

OR




For Kinnex single-cell 5' RNA analysis, select '10x Chromium single cell 5' cDNA primers'



For Kinnex full-length RNA analysis, select 'Iso-Seq v2 Barcoded cDNA primers'

Specify primer sequence file in FASTA format to identify cDNA primers for removal (include the 5' and 3' cDNA primers)



SMRT Link Run Design analysis options (cont.)

Kinnex single-cell RNA

Analysis Options

Add Analysis YES NO

Analysis Name Required: Kinnex_Single-Cell_RNA_Demo_Analysis_Job_Name

Select Analysis Workflow Required: Read Segmentation and Single-Cell Iso-Seq

Segmentation Adapter Set: MAS-Seq Adapter v1 (MAS16)

Primer Set Required: 10x Chromium single cell 3' cDNA

Reference Set Required: Human Genome hg38, with Gencode v39 annotations

Kit Type 10x 3' KIT 10x 5' KIT

Advanced Parameters

Specify Kit Type used for single-cell cDNA generation

Specification of Kit Type (10x 3' Kit or 10x 5' Kit) **determines which set of 10x barcode sequences to use**, and also affects **UMI and single-cell barcode design settings**¹

5' [CBC][UMI]TTTTTTTTT LIBRARY INSERT CCC TSO 3'
 3' [CBC][UMI]AAAAAAAAA LIBRARY INSERT GGG TSO 5'

For Kinnex single-cell 3' RNA analysis, select '10x 3' Kit'

OR

5' [CBC][UMI][TSO]GGG LIBRARY INSERT AAAAAAAAAA 3'
 3' [CBC][UMI][TSO]CCC LIBRARY INSERT TTTTTTTTTT 5'

For Kinnex single-cell 5' RNA analysis, select '10x 5' Kit'

Kinnex full-length RNA

Analysis Options

Add Analysis YES NO

Analysis Name Required: [blank]

Select Analysis Workflow Required: [blank]

Segmentation Adapter Set: MAS-Seq Adapter v1 (MAS8)

Primer Set Required: 10x 2 Barcoded cDNA Primers

Reference Set Required: Human Genome hg38, with Gencode v39 annotations

Cluster of Barcoded Samples Pool reads and cluster together

Advanced Parameters

Specify how to perform read clustering for barcoded samples

Pool reads and cluster together







-- select --

Cluster reads separately

Pool reads and cluster together

Specification of **Cluster of Barcoded Samples** setting **determines whether all FLNC reads will be pooled for clustering** (Does not apply to non-barcoded samples.)

Clustering options for Kinnex full-length RNA samples






Kinnex single-cell RNA example sequencing performance data

Example Kinnex single-cell RNA library preparation QC results


Kinnex single-cell 3' RNA library prepared with human cDNA

Kinnex PCR DNA sizing QC (Single-cell 3' cDNA)



Example Bioanalyzer DNA sizing QC analysis results for Kinnex PCR products generated for a 10x Chromium single-cell 3' cDNA samples prepared from a human cell line (HG002).

Final Kinnex single-cell RNA library QC



Example Femto Pulse DNA sizing QC analysis results for final Kinnex single-cell RNA library.

Final Kinnex library yield is typically sufficient to load the following number of SMRT Cells:

- >8 SMRT Cells for Revio using SPRQ chemistry
- >2 SMRT Cells for the Vega or Revio (non-SPRQ) systems
- >4 SMRT Cells for Sequel II/IIe systems


10x single cell 3' cDNA input	15 ng
cDNA input for Kinnex array formation	5000 ng
Post-nuclease treatment & final library cleanup yield (%) ¹	1100 ng (22.0%)

¹ Post-nuclease treatment & final cleanup yields typically ranged from ~10% to ~25% when using single-cell 3' cDNA samples for Kinnex single-cell RNA library construction.

Example Kinnex single-cell RNA library preparation QC results (cont.)


Kinnex single-cell 5' RNA library prepared with human cDNA

Kinnex PCR DNA sizing QC (Single-cell 5' cDNA)



Example Bioanalyzer DNA sizing QC analysis results for Kinnex PCR products generated for a 10x Chromium single-cell 5' cDNA samples prepared from a human cell line (HG002).

Final Kinnex single-cell RNA library QC



Example Femto Pulse DNA sizing QC analysis results for final Kinnex single-cell RNA library.

Final Kinnex library yield is typically sufficient to load the following number of SMRT Cells:

- >8 SMRT Cells for Revio using SPRQ chemistry
- >2 SMRT Cells for the Vega or Revio (non-SPRQ) systems
- >4 SMRT Cells for Sequel II/IIe systems


10x single cell 5' cDNA input	15 ng
cDNA input for Kinnex array formation	5000 ng
Post-nuclease treatment & final library cleanup yield (%) ¹	1008 ng (20.2%)

¹ Post-nuclease treatment & final cleanup yields typically ranged from ~10% to ~25% when using single-cell 5' cDNA samples for Kinnex single-cell RNA library construction.


Example sequencing performance for Kinnex single-cell RNA libraries prepared with human cDNA [**Revio system + SPRQ chemistry¹**]

Kinnex single-cell RNA 5' library for PBMC single-cell cDNA sample from 10x Chromium Single Cell Universal 5' kit

Raw Data Report




HiFi Read Length



- For Revio Kinnex single-cell RNA samples**
- Yield of HiFi reads is ideally >6 Million
 - Yield of segmented reads is ideally >100 Million
 - Mean S-read length will depend on size distribution of input 10x single cell cDNA (ideally >600 bp)
 - % of reads with full arrays is ideally >80%
 - Mean array size is ideally ≥15 segments

Read Segmentation Metrics



Raw Base Yield	1,539 Gb
Mean Polymerase Read Length	90.2 kb
P0	32%
P1	68%
P2	0%

Example sequencing metrics for a human Kinnex single-cell RNA 5' library sample run on a Revio system with Revio SPRQ polymerase kit / 140 pM on-plate loading concentration (OPLC) / 24-hrs movie time.

HiFi Reads	9.4 M
HiFi Base Yield	144.4 Gb
Mean HiFi Read Length	15.3 kb
Median HiFi Read Quality	Q29
HiFi Read Mean # of Passes	9

For human Kinnex single-cell RNA libraries, per-Revio SMRT Cell HiFi read counts were typically >6 Million depending on the final library insert size and P1 loading performance.


Input HiFi Reads	9,408,444
Segmented reads (S-reads)	144,801,149
Mean length of S-reads	971 bp
Percent of reads with full arrays	92.59%
Mean array size (concentration factor)	15.39

For Kinnex single-cell RNA libraries, per-Revio SMRT Cell segmentation read counts were typically >100 Million.


Example sequencing performance for Kinnex single-cell RNA libraries prepared with human cDNA [**Revio system + SPRQ chemistry¹**] (cont.)

Kinnex single-cell RNA 3' library for PBMC single-cell cDNA sample from 10x Chromium Single Cell Universal 3' kit


Raw Data Report



HiFi Read Length



Read Segmentation Metrics



For Revio Kinnex single-cell RNA samples

- Yield of HiFi reads is ideally >6 Million
- Yield of segmented reads is ideally >100 Million
- Mean S-read length will depend on size distribution of input 10x single cell cDNA (ideally >600 bp)
- % of reads with full arrays is ideally >80%
- Mean array size is ideally ≥15 segments

Raw Base Yield	1,570 Gb
Mean Polymerase Read Length	86.5 kb
P0	27%
P1	72%
P2	0%

Example sequencing metrics for a human Kinnex single-cell RNA 3' library sample run on a Revio system with Revio SPRQ polymerase kit / 140 pM on-plate loading concentration (OPLC) / 24-hrs movie time.

HiFi Reads	9.6 M
HiFi Base Yield	148.1 Gb
Mean HiFi Read Length	15.4 kb
Median HiFi Read Quality	Q29
HiFi Read Mean # of Passes	8

For human Kinnex single-cell RNA libraries, per-Revio SMRT Cell HiFi read counts were typically >6 Million depending on the final library insert size and P1 loading performance.


Input HiFi Reads	9,605,565
Segmented reads (S-reads)	149,894,918
Mean length of S-reads	963 bp
Percent of reads with full arrays	92.24%
Mean array size (concentration factor)	15.61

For Kinnex single-cell RNA libraries, per-Revio SMRT Cell segmentation read counts were typically >100 Million.


Example sequencing performance for Kinnex single-cell RNA libraries prepared with human cDNA [Vega system¹]

Kinnex single-cell RNA 3' library for PBMC single-cell cDNA sample from 10x Chromium Single Cell Universal 3' kit


Raw Data Report



HiFi Read Length



Read Segmentation Metrics



- For Vega Kinnex single-cell RNA samples**
- Yield of HiFi reads is ideally >4 Million
 - Yield of segmented reads is ideally >50 Million
 - Mean S-read length will depend on size distribution of input 10x single cell cDNA (ideally >600 bp)
 - % of reads with full arrays is ideally >80%
 - Mean array size is ideally ≥15 segments

Mean Polymerase Read Length	94.1 kb
Loading Level	71%

Example sequencing metrics for a human Kinnex single-cell RNA 3' library sample run on a Vega system with Vega polymerase kit / 150 pM on-plate loading concentration (OPLC) / 24-hrs movie time.

HiFi Reads	4.4 M
HiFi Base Yield	64.9 Gb
Mean HiFi Read Length	14.8 kb
Median HiFi Read Quality	Q29
HiFi Read Mean # of Passes	9

For human Kinnex single-cell RNA libraries, per-Vega SMRT Cell HiFi read counts were typically >4 Million depending on the final library insert size and sample loading performance.

Input HiFi Reads	4,372,029
Segmented reads (S-reads)	68,570,669
Mean length of S-reads	924 bp
Percent of reads with full arrays	93.78%
Mean array size (concentration factor)	15.68

For Kinnex single-cell RNA libraries, per-Vega SMRT Cell segmentation read counts were typically >50 Million.




Kinnex single-cell RNA data analysis workflow overview

Kinnex single-cell RNA bioinformatics workflow overview

SMRT Link Read Segmentation and Single-cell Iso-Seq workflow processes HiFi reads generated from Kinnex single-cell RNA libraries to produce classified isoforms with read counts that are compatible with tertiary analysis tools

SMRT Link read segmentation and single-cell Iso-Seq workflow¹



Read segmentation

- HiFi reads are segmented into individual segmented reads (**S-reads**) that represent the original cDNA sequences

Primer/UMI/barcode extraction

- Primers and polyA tails are removed, but also used to orient the read into 5' – 3' orientation
- Single-cell barcode and UMI information are extracted

Barcode correction & UMI deduplication

- Cell barcodes are corrected given an expected barcode list.
- Real cells — cell barcodes that represent encapsulated single cells (as opposed to ambient RNA) are also identified at this step.
- Reads are then deduplicated based on cell barcodes and UMIs.

Mapping and transcript classification²

- Deduplicated reads are mapped to the reference genome and classified against a transcript annotation (e.g., GENCODE).
- Finally, a gene- and isoform-level single-cell matrix is output for tertiary analysis.

SMRT Link Read Segmentation and Single-Cell Iso-Seq analysis application setup

Specify **Read Segmentation and Single-Cell Iso-Seq** analysis application type in SMRT Link¹

Enables automated analysis and functional characterization of full-length transcript isoforms with additional single-cell information, including single-cell barcodes & unique molecular identifiers (UMIs)

Analysis Application Required
Read Segmentation and Single-Cell Iso-Seq

Analysis Name
SMRT Analysis Demo - Creating a New Analysis

Analysis Datasets
Displaying rows 1 to 1 out of 1

ID #	Name #
59241	3pHG2_VERF_DLJZ_bc01

- Accepts **HiFi reads** (BAM format) as input
- HiFi reads are reads generated with CCS analysis whose quality value is equal to or greater than 20.

SMRT Link Read Segmentation and Single-Cell Iso-Seq analysis application setup (cont.)

Specify **Read Segmentation and Single-Cell Iso-Seq** analysis application required associated inputs¹

PacBio SMRT Analysis ▾

SMRT Analysis / Create New Analysis

1. Select Data 2. Select Analysis

Analysis Application Required

Read Segmentation and Single-Cell Iso-Seq

↑ Import Analysis Settings ↓ Export ⓘ

Associated Inputs

1 Segmentation Adapter Set
MAS-Seq Adapter v1 (MAS16)

2 Primer Set Required
10x Chromium single cell 3' cDNA primers

3 Reference Set Required
Human Genome hg38, with Gencode v39 annotatic

4 Kit Type ⓘ
 10X 3' KIT 10X 5' KIT

Advanced Parameters

1. Segmentation Adapter Set (Required)

- Specify a FASTA file, provided by PacBio, containing segmentation adapters. If you need a custom segmentation adapter set, click Advanced Parameters and use a custom FASTA file formatted as described in the SMRT Link User Guide [documentation](#).

2. Primer Set (Required)

- Specify a primer sequence file in FASTA format to identify cDNA primers for removal. The primer sequence includes the 5' and 3' cDNA primers.
- Primer IDs must be specified using the suffix `_5p` to indicate 5' cDNA primers and the suffix `_3p` to indicate 3' cDNA primers. The 3' cDNA primer should not include the Ts and is written in reverse complement.
- Each primer sequence must be unique.

3. Reference Set (Required)

- Specify one of two default reference genome and annotation sets to align high quality isoforms to, and to collapse isoforms mapped to the same genomic loci. The default sets are `Human_hg38_Gencode_v39` and `Mouse_mm39_Gencode_vM28`.

4. Segmentation Adapter Set (Required)

- Specify the 10x 3' Kit, or 10x 5' Kit. This determines which set of 10x primers and barcode sequences to use, and also affects the UMI and single-cell barcode design settings.

Example SMRT Link Read Segmentation data utility processing results¹ for Kinnex single-cell RNA libraries prepared with PBMC single cell cDNA

SMRT Link Read Segmentation data utility job report – Summary Metrics and Segmentation Statistics

Summary Metrics


Value	Analysis Metric
9,709,186	Reads
152,811,229	Segmented reads (S-reads)
928	Mean length of S-reads
94.89 %	Percent of reads with full arrays
15.74	Mean array size (concatenation factor)

Example Revio system (+SPRQ) data shown.

For Revio Kinnex single-cell RNA libraries under optimal sample loading conditions:

- Yield of HiFi reads is ideally >6 Million
- Yield of segmented reads is ideally >100 Million
- Mean S-read length will depend on size distribution of input 10x single cell cDNA (ideally >600 bp)
- % of reads with full arrays is ideally >80%
- Mean array size is ideally ≥ 15 segments

Segmentation Statistics



- Plot should ideally show a majority (>80%) of reads with array lengths = 16 segments

- **Reads:** Number of input arrayed HiFi reads
- **Segmented reads (S-reads):** Number of generated S-reads
- **Mean length of S-reads:** Mean read length of generated S-reads
- **Percent of reads with full arrays:** Percentage of input HiFi reads containing all adapter sequences in the order listed in the segmentation adapter FASTA file
- **Mean array size:** Mean number of fragments (or S-reads) found in input reads


Histogram distribution of the number of S-reads per HiFi read. (Example Revio system + SPRQ data shown.)

Example SMRT Link Read Segmentation data utility processing results for Kinnex single-cell RNA libraries prepared with PBMC single cell cDNA (cont.)

SMRT Link Read Segmentation data utility job report – Length of Reads and S-read Length


Length of Reads

- Mode should be **concordant** with the expected Kinnex single-cell RNA library insert size, with a majority of HiFi reads containing 16-segment arrays



Histogram distribution of the number of HiFi reads by read length, in base pairs. (Example Revio system + SPRQ data shown.)

S-read Length



Histogram distribution of the number of S-reads by HiFi read length, in base pairs. (Example Revio system + SPRQ data shown.)

Example SMRT Link Single-Cell Iso-Seq Analysis results for Kinnex single-cell RNA libraries prepared with PBMC single cell cDNA

SMRT Link Single-Cell Iso-Seq Analysis job report – Read Statistics

Summary Metrics

Value	Analysis Metric
152,811,229	Reads
SEGMENT	Read Type
151,565,748	Reads with 5' and 3' Primers with extracted UMIs and Barcodes
147,805,910	Non-Concatamer Reads with 5' and 3' Primers and Poly-A Tail (FLNC reads)
144,807,139	FLNC Reads with Valid Barcodes
147,466,710	FLNC Reads with Valid Barcodes, corrected
110,681,049	Reads after Barcode Correction and UMI Deduplication

- **Reads:** Total number of input reads for analysis.
- **Read Type:** Type of input reads - CCS, SEGMENT, or mixed if there are multiple input data sets with mixed data types.
- **Reads with 5' and 3' Primers with extracted UMIs and Barcodes:** The number of reads with 5' and 3' cDNA primers detected, and UMI/cell barcode information extracted. Also known as full-length tagged reads (FLT Reads).
- **Non-Concatamer Reads with 5' and 3' Primers and Poly-A Tail (FLNC Reads):** The number of non-concatemer reads with 5' and 3' primers and polyA tails detected after UMI/cell barcode information has been extracted.
- **FLNC Reads with Valid Barcodes:** Number of full-length non-concatemer reads that include valid single-cell barcodes.
- **FLNC Reads with Valid Barcodes, corrected:** Number of full-length non-concatemer reads that include valid single-cell barcodes, after barcode correction.
- **Reads after Barcode Correction and UMI Deduplication:** Number of deduplicated reads, after barcode correction.

Example Revio system + SPRQ data shown.

Example SMRT Link Single-Cell Iso-Seq Analysis results for Kinnex single-cell RNA libraries prepared with PBMC single cell cDNA (cont.)

SMRT Link Single-Cell Iso-Seq Analysis job report – Cell Statistics

Summary Metrics


Value	Analysis Metric
10,324	Estimated Number of Cells
90.65%	Reads in Cells
12,978	Mean Reads per Cell
9,194	Median UMIs per Cell

- There is no “correct” number of cells – this metric depends on what was specified in the 10x Chromium single cell workflow as the intended target cell recovery

Example Revio system + SPRQ data shown.

- **Estimated Number of Cells:** The estimated number of cells.
- **Reads in Cells:** The percentage of reads in cells.
- **Mean Reads per Cell:** The mean number of reads per cell.
- **Median UMIs per Cell:** The median number of unique molecular identifiers (UMIs) per cell.

Barcode Rank Plot



Displays the distribution of UMI counts and which barcodes were inferred to be associated with cells. The X-axis denotes barcodes ranked in decreasing order by UMI counts mapped to each barcode, and the Y-axis denotes the UMI count for the N-th ranked barcode. (Example Revio system + SPRQ data shown.)

Example SMRT Link Single-Cell Iso-Seq Analysis results for Kinnex single-cell RNA libraries prepared with PBMC single cell cDNA (cont.)

SMRT Link Single-Cell Iso-Seq Analysis job report – Transcript Statistics

Summary Metrics

Value	Analysis Metric
120,736,189	FLNC Reads Mapped Confidently to Genome
79,264,395	FLNC Reads Mapped Confidently to Transcriptome
1,902,966	Total Unique Genes
98,700	Total Unique Genes, filtered
37,474	Total Unique Genes, known genes only
28,955	Total Unique Genes, filtered, known genes only
4,399,569	Total Unique Transcripts
999,042	Total Unique Transcripts, filtered
103,849	Total Unique Transcripts, known transcripts only
93,952	Total Unique Transcripts, filtered, known transcripts only

Example Revio system + SPRQ data shown.

- **FLNC reads mapped confidently to genome:** The number of FLNC reads mapped to the reference genome. This number is calculated first based on the number of deduplicated reads mapped to the genome, then expanded to account for duplicate FLNC reads for each unique molecule.
- **FLNC reads mapped confidently to transcriptome:** The number of FLNC reads mapped to the reference genome in which the read is later associated with a transcript that is classified as one of the following: FSM, ISM, NIC, or NNC.
- **Total unique genes:** The total number of unique genes across all cells.
- **Total unique genes, filtered:** The total number of unique genes, after filtering out reads based on the SQANTI transcript filtering criteria.
- **Total unique genes, known genes only:** The total number of unique genes across all cells in which the gene is annotated in the reference annotation.
- **Total unique genes, filtered, known genes only:** The total number of unique genes (genes annotated in the reference annotation) across all cells, after filtering out reads based on the SQANTI transcript filtering criteria.
- **Total unique transcripts:** The total number of unique transcripts across all cells.
- **Total unique transcripts, filtered:** The total number of unique transcripts across all cells, after filtering out reads based on the SQANTI transcript filtering criteria.
- **Total unique transcripts, known transcripts only:** The total number of unique transcripts across all cells in which the gene the transcript belongs to is annotated in the reference annotation.
- **Total unique transcripts, filtered, known transcripts only:** The total number of unique transcripts across all cells, after filtering out reads based on the SQANTI transcript filtering criteria. Only transcripts associated with known genes (genes annotated in the reference annotation) are included.

Example SMRT Link Single-Cell Iso-Seq Analysis results for Kinnex single-cell RNA libraries prepared with PBMC single cell cDNA (cont.)

SMRT Link Single-Cell Iso-Seq Analysis job report – Transcript Statistics

Transcript Summary

Value	Analysis Metric
3,369	Median Genes per Cell
2,476	Median Genes per Cell, known genes only
4,522	Median Transcripts per Cell
2,145	Median Transcripts per Cell, known transcripts only
1,902,966	Total Unique Genes
37,474	Total Unique Genes, known genes only
4,399,569	Total Unique Transcripts
103,849	Total Unique Transcripts, known transcripts only



Filter out reads based on the SQANTI3 transcript filtering criteria¹

Example Revio system + SPRQ data shown.

Transcript Summary, Filtered

Value	Analysis Metric
1,822	Median Genes per Cell
1,787	Median Genes per Cell, known genes only
2,389	Median Transcripts per Cell
1,839	Median Transcripts per Cell, known transcripts only
98,700	Total Unique Genes
28,955	Total Unique Genes, known genes only
999,042	Total Unique Transcripts
93,952	Total Unique Transcripts, known transcripts only

- **Median genes per cell:** The median number of genes per cell.
- **Median genes per Cell, known genes only:** The median number of unique, known genes (genes annotated in the reference annotation) per input cell.
- **Median transcripts per cell:** The median number of transcripts per cell.
- **Median transcripts per cell, known transcripts only:** The median number of transcripts per cell. Only transcripts associated with known genes are included.

- **Total unique genes:** The total number of unique genes across all cells.
- **Total unique genes, known genes only:** The total number of unique, known genes (genes annotated in the reference annotation) across all cells.
- **Total unique transcripts:** The total number of unique transcripts across all cells.
- **Total unique transcripts, known transcripts only:** The total number of unique transcripts across all cells. Only transcripts associated with known genes are included.

Example SMRT Link Single-Cell Iso-Seq Analysis results for Kinnex single-cell RNA libraries prepared with PBMC single cell cDNA (cont.)

SMRT Link Single-Cell Iso-Seq Analysis job report – Transcript Statistics

Transcript Classification, filtered (All samples)

Category IT	Count IT	CAGE Detected IT	CAGE Detected, (%) IT	polyA Detected IT	polyA Detected, (%) IT
FSM	167573	103997	62.06%	58617	34.97%
ISM	314765	158530	50.36%	111958	35.56%
NIC	273808	190384	69.53%	109290	39.91%
NNC	514400	345302	67.12%	232608	45.21%
Antisense	17604	2071	11.76%	9938	56.45%
Fusion	16945	10215	60.28%	8019	47.32%
More junctions	282	158	56.02%	127	45.03%
Genic intron	0	0	0.00%	0	0.00%
Genic genomic	18361	9601	52.29%	7788	42.41%
Intergenic	63690	2479	3.89%	43701	68.61%


Example Revio system + SPRQ data shown.¹

- **Category:** Transcript classification² assigned by the classification and filtering tool `pigeon`, based on the [SQANTI3](#) software
- **Count:** The number of transcripts, after filtering out reads based on the SQANTI filtering criteria, in a specific classification
- **CAGE Detected:** The number of transcripts where the transcription start site falls within 50 bp of an annotated CAGE (Cap Analysis of Gene Expression) peak site
- **CAGE Detected, (%):** The percentage of transcripts where the transcription start site falls within 50 bp of an annotated CAGE peak site
- **polyA Motif Detected:** The number of transcripts where a known polyA motif is detected upstream of the transcription end site
- **polyA Motif Detected, (%):** The percentage of transcripts where a known polyA motif is detected upstream of the transcription end site


Example SMRT Link Single-Cell Iso-Seq Analysis results for Kinnex single-cell RNA libraries prepared with PBMC single cell cDNA (cont.)

SMRT Link Single-Cell Iso-Seq Analysis job report – Transcript Statistics

Transcript Classification Plots



Transcript Classification Plots, Filtered



Filter out reads based on the SQANTI3 transcript filtering criteria

Isoform distributions across structural categories:


- Distribution of the % of isoforms by structural categories.

Example Revio system + SPRQ data shown.


Structural categories by isoform lengths:

- Histogram display of the number of isoforms by their length in kb and their structural category.

Example Revio system + SPRQ data shown.




Filter out reads based on the SQANTI3 transcript filtering criteria



Example SMRT Link Single-Cell Iso-Seq Analysis results for Kinnex single-cell RNA libraries prepared with PBMC single cell cDNA (cont.)

SMRT Link Single-Cell Iso-Seq Analysis job report – Transcript Statistics

Gene Saturation



Gene Saturation, all genes, filtered: Saturation plot showing the level of gene saturation for **all genes**, after filtering out reads based on the SQANTI transcript filtering criteria

Gene Saturation, known genes only, filtered: Saturation plot showing the level of gene saturation, for unique **known genes only** (genes annotated in the reference annotation) per cell, after filtering out reads based on the SQANTI transcript filtering criteria

Example SMRT Link Single-Cell Iso-Seq Analysis results for Kinnex single-cell RNA libraries prepared with PBMC single cell cDNA (cont.)

File downloads tab

Edit Output File Name Prefix Example:analysis-Bio Sample 3-78360

File ↑	Size ↓	Type ↓
Non-passing reads, unaligned	4 GB	bam
Report read_segmentation	3 KB	JsonReport
SMRT Link Log	14 KB	log
Segmented Reads, passing, unaligned	57 GB	bam
Single-cell isoform and gene matrix, tar-gzipped	1 GB	tgz
Unique mapped transcripts, GFF	1 GB	gff
Unique mapped transcripts, classification TXT	788 MB	txt
Unique mapped transcripts, filtered, GFF	481 MB	gff
Unique mapped transcripts, filtered, classification TXT	246 MB	txt
Unique mapped transcripts, filtered, junctions TXT	451 MB	txt
Unique mapped transcripts, junctions TXT	785 MB	txt

Refer to [SMRT Link user guide](#) for descriptions of downloadable output files

- **Key output file!**
- Gzipped file containing Seurat-compatible isoform and gene matrix files

- These files are useful for **visualizing** isoform structures in Integrative Genomics Viewer (IGV) / UCSC genome browser and enable understanding of why an isoform is novel/known, etc.
 - GFF file containing unique mapped transcripts after filtering
 - Text file containing unique mapped transcript classifications against annotations, after filtering
 - Text file containing information about unique mapped transcript junctions, after filtering

Files shown in the File Downloads tab are available on the analysis results page. Additional files are also available on the SMRT Link server in the analysis output directory.



Technical documentation & applications support resources

Technical resources for Kinnex single-cell RNA library preparation, sequencing & data analysis

Single-cell cDNA sample preparation literature & other resources

- 10x Genomics Chromium Next GEM Single Cell 3' v3.1 (Single Index) How-to Video [[Link](#)]
- 10x Genomics Chromium Single Cell 3' Reagent Kits User Guide – v3.1 ([CG000204](#))
- 10x Genomics Chromium Single Cell 5' Reagent Kits User Guide – v2 Chemistry Dual Index ([CG000331](#))

Kinnex single-cell RNA library preparation literature & other resources

- Application note – Kinnex single-cell RNA for single-cell isoform sequencing ([102-326-549](#))
- Procedure & checklist – Preparing Kinnex libraries using Kinnex single-cell RNA kit ([102-254-300](#))
- Technical overview – Kinnex kits for single-cell RNA, full-length RNA and 16S rRNA sequencing ([103-343-700](#))
- Technical overview – Kinnex library preparation using Kinnex single-cell RNA kit ([103-344-600](#))
- Video tutorial – PacBio Kinnex single-cell RNA TSO artifact removal demo for Kinnex single-cell RNA kit [[Link](#)]
- Video tutorial – SMRT Link Sample Setup and Run Design setup procedure for Kinnex kits [[Link](#)]

Data analysis resources

- Application note – Bioinformatics tools for full length isoform sequencing ([102-326-593](#))
- SMRT Link Kinnex full-length RNA troubleshooting guide ([103-552-100](#))
- SMRT Link Kinnex single-cell troubleshooting guide ([103-516-100](#))
- SMRT Link MAS-Seq troubleshooting guide ([102-994-400](#))
- SMRT Link software installation guide [[Link](#)]
- SMRT Link user guide [[Link](#)]
- SMRT Tools reference guide [[Link](#)]

Technical resources for Kinnex single-cell RNA library preparation, sequencing & data analysis (cont.)

Publications

- Wissel, D. et al. (2025) A systematic benchmark of high-accuracy PacBio long-read RNA sequencing for transcript-level quantification. BioRxiv [[Link](#)]
- Al'Khafaji, A.M. et al. (2024) High-throughput RNA isoform sequencing using programmable cDNA concatenation. Nature biotechnology [[Link](#)]
- Deng, E. et al. (2025) Systematic evaluation of single-cell RNA-seq analyses performance based on long-read sequencing platforms. J Advanced Res [[Link](#)]

Webinars

- PacBio webinar (2025) - Detecting shared mis-splicing across blood cancers using long-read single-cell RNA sequencing [[Link](#)]
- PacBio video (2024) – Kinnex explained – how concatenating smaller amplicons increases throughput for PacBio HiFi sequencing [[Link](#)]
- PacBio PRISM webinar (2024) – Let's stick together – exploring PacBio Kinnex kits [[Link](#)]
- PacBio Iso-Seq social club webinar (2022) – TappAS for isoform differential expression analysis [[Link](#)]
- PacBio Iso-Seq social club webinar (2022) – Single-cell Iso-Seq applications in cancer and neurological disorders [[Link](#)]

Example PacBio data sets

Application	Dataset	Data type	PacBio system
Kinnex single-cell RNA sequencing	Homo sapiens - PBMC 10x Chromium Single Cell 5' and 3' libraries [Link]	HiFi long read	Sequel II & Revio systems
	Homo sapiens - HG002 (10x 5') [Link]	HiFi long read	Revio system



APPENDIX 1: Concatenating cDNA or gDNA amplicons using Kinnex kits to increase throughput

Concatenating cDNA or gDNA amplicons using Kinnex kits to increase throughput

Kinnex concatenation is a general method that can increase sequencing throughput for smaller amplicons¹

Benefits of Kinnex concatenation

- Increased throughput on PacBio long-read sequencers
- Retained HiFi accuracy despite throughput increase
- No change to secondary analysis – once reads are deconcatenated into S-reads, the S-reads represent the original, pre-concatenated amplicon and can be analyzed with established pipelines.

When is Kinnex concatenation appropriate?

- The balance between the amplicon size and the concatenation factor, as well as additional Kinnex library generation cost, needs to be taken into consideration.
- HiFi sequencing produces optimal yield for inserts between 15–20 kb; therefore, the throughput advantage plateaus for larger amplicon sizes exceeding 3 kb.
- You can consider concatenating amplicons using Kinnex kits if:
 - The amplicons have an **average size between 200 bp – 3 kb**
 - The amplicons have **molecular ends that are either directly compatible or can be re-amplified to establish Kinnex compatibility**

Recommended Kinnex kit based on average amplicon sizes.

Average amplicon size	Example	Recommended Kinnex kit	Expected Kinnex library size
600 – 1000 bp	10x single-cell cDNA	Kinnex single-cell RNA kit (16-fold)	10 – 16 kb
1 – 2 kb	Full-length 16S	Kinnex 16S rRNA kit (12-fold)	~19 kb
2 – 3 kb	Bulk cDNA	Kinnex full-length RNA kit (8-fold)	15 – 20 kb
>3 kb	Not recommended for Kinnex concatenation		

Technical note
CONCATENATING AMPLICONS USING PACBIO KINNEX KITS TO INCREASE THROUGHPUT

Overview
This technical note describes the principles and advantages of concatenating amplicons using PacBio® Kinnex™ kits.

Note: This technical note is intended as a guide for best practices and to report typical example results (see Appendix). PacBio does not guarantee the success of concatenating amplicons that are not officially supported by the Kinnex kits, which currently include:

- Full-length cDNA generated using the Iso-Seq® express 2.0 kit with the Kinnex full-length RNA kit
- Full-length 16S sequences generated according to official protocol with the Kinnex 16S rRNA kit
- Single-cell cDNA generated using compatible single-cell platforms with the Kinnex single-cell RNA kit


How Kinnex works

What is Kinnex?
The Kinnex kits are based on the method developed by Al'Khafaji et al. (2023). Originally termed multiplexed array isomorph sequencing (MAS-ISO-seq, or MAS-Seq), this method was developed to concatenate cDNAs into longer fragments suitable for long-read sequencing. It takes advantage of the disparity between optimal DNA fragment lengths for HiFi sequencing, 15–20 kb, and smaller sizes of transcript cDNA (1–10 kb with typical estimated average length of 2 kb for a human transcriptome). The original MAS-ISO-Seq method concatenated single-cell cDNA from the 10x platform, which has an average read length of 600–800 bp, and used a 15-fold concatenation array to increase throughput. PacBio commercialized the MAS-Seq method and adjusted the concatenation factors according to different average amplicon sizes (Table 1).

Technical note – Concatenating cDNA or gDNA amplicons using Kinnex kits to increase throughput
[\(102-326-636\)](https://www.pacb.com/technical-note/concatenating-amplicons-using-kinnex-kits-to-increase-throughput/)

Kinnex concatenation workflow overview

Follow Kinnex library prep protocol documentation for specific details on concatenating cDNA or gDNA amplicons using Kinnex kits



Kinnex concatenation workflow. Amplicons must be generated or amplified to have Kinnex-compatible molecular ends before continuing to the Kinnex PCR and array formation step. Kinnex libraries should be sequenced with the appropriate sequencing chemistry and run configurations. Once de-concatenated using Read Segmentation in SMRT Link, the individual amplicons can be analyzed using amplicon-specific workflows.

Kinnex concatenation procedural notes¹

- To establish Kinnex compatibility, **amplicons are required to have Kinnex-compatible molecular ends** (see next section)
 - Once these are generated, choose the appropriate Kinnex kit based on the recommended concatenation factor listed in Table on previous slide and proceed with Kinnex PCR
- The **Kinnex PCR** steps consist of parallel PCR reactions per sample [i.e., 8, 12, or 16 reactions based on the Kinnex kit chosen] using premixed Kinnex primer pairs
 - The resulting PCRs generate amplified DNA products containing programmable sequences at both ends
- In the **Kinnex array formation** step, library inserts containing programmable ends are assembled to generate a linear array.
 - Further, the addition of barcoded Kinnex terminal adapters result in the formation of complete, full-length array SMRTbell templates along with partial arrays.
 - Subsequent nuclease treatment removes partial arrays to retain only full array SMRTbell templates for achieving optimal sequencing yield.
- **Note:** Kinnex terminal adapters are different from standard SMRTbell adapters and hence require the **Kinnex sequencing primer** (103-179-000) during the “Annealing, Binding, and Cleanup (ABC)” step for optimal sequencing results.
- Once HiFi reads are generated, **Read Segmentation** will produce the segmented reads (S-reads) that represent the original unconcatenated amplicons, which can be used for further analysis.

Establishing Kinnex-compatible molecular ends

Amplicons must be generated or amplified to have Kinnex-compatible molecular ends before proceeding with Kinnex library prep procedure

- To be compatible with the Kinnex workflow, amplicons must be generated with sequence-defined ends as depicted in the underlined portion in the figure below
- Optional barcodes (such as sample indices, UMIs and single cell barcodes) should be placed internally between the Kinnex handles and the amplicon-specific primers.
- Kinnex handles may be present already in certain amplicons, such as 10x Single Cell Gene Expression libraries or the Kinnex 16S amplicons, or can be added by PCR amplification, such as for Parse Evercode single-cell libraries

Kinnex FWD primer CTACACGACGCTCTTCCGATCT - [optional barcodes] - [amplicon specific FWD primers]
Kinnex REV primer AAGCAGTGGTATCAACGCAGAG - [optional barcodes] - [amplicon specific REV primers]

Example 1. Kinnex 16S forward and reverse primer sequences

Kinnex 16S FWD01 CTACACGACGCTCTTCCGATCT - GATCGAGTCA - AGRGTTYGATYMTGGCTCAG
Kinnex 16S REV13 AAGCAGTGGTATCAACGCAGAG - TCATCGACGT - RGYTACCTTGTACGACTT

Example 2. Iso-Seq express 2.0 forward and reverse primer sequences

IsoSeqX bc01 FWD CTACACGACGCTCTTCCGATCT - ACTACAC - GCAATGAAGTCGCAGGGTTGGG
IsoSeqX REV AAGCAGTGGTATCAACGCAGAGTAC

Schematic for Kinnex-compatible primers. Kinnex handles (5' to 3') are shown in black underline and must be present at the ends of the amplicons to be compatible with Kinnex concatenation. Optional barcodes can be included internally. Amplicon-specific primers (and optional internal barcode sequences) must be designed to avoid strong secondary structures in the context of Kinnex handles.

Resources for Kinnex library prep¹

For using Kinnex full-length RNA kit (PN: 103-072-000) for 8-fold concatenation:

- Procedure & checklist – Preparing Kinnex libraries using the Kinnex full-length RNA kit ([103-238-700](#))
- Technical overview – Kinnex library preparation using Kinnex full-length RNA kit ([103-344-700](#))

For using Kinnex 16S rRNA kit (PN: 103-072-100) for 12-fold concatenation:


- Procedure & checklist – Preparing Kinnex libraries using 16S rRNA amplicons ([103-238-800](#))
- Technical overview – Kinnex library preparation using Kinnex 16S rRNA kit ([103-344-800](#))

For using Kinnex single-cell RNA kit (PN: 103-072-200) for 16-fold concatenation:

- Procedure & checklist – Preparing Kinnex libraries using Kinnex single-cell RNA kit ([103-254-300](#))
- Technical overview – Kinnex library preparation using Kinnex single-cell RNA kit ([103-344-600](#))



APPENDIX 2: PacBio compatible Kinnex library preparation workflows







Preparing Kinnex single-cell libraries with Parse Evercode Whole Transcriptome kits

Preparing Kinnex single-cell libraries with Parse Evercode Whole Transcriptome kits: Getting started

Application-specific educational literature

Application-specific protocol documentation

Application-specific technical overviews



PacBio literature website [Link]

Application-specific brochures, informational guides and other product literature containing best practices recommendations for library preparation and data analysis workflows.

Technical note – Preparing Kinnex single-cell libraries with Parse Evercode Whole Transcriptome kits (102-326-610) [PacBio]

Technical documentation describing a modified workflow for preparing PacBio Kinnex single-cell libraries with single-cell or single-nucleus cDNA generated using Parse Biosciences Evercode Whole Transcriptome products.

Procedure & checklist – Procedure & checklist – Preparing Kinnex libraries using Kinnex single-cell RNA kit (103-254-300) [PacBio]

Technical documentation containing PacBio SMRTbell library construction details.

Technical overview: Kinnex library preparation using Kinnex single-cell RNA kit (103-344-600)

Technical overview presentations describe sample preparation details for constructing HiFi libraries for specific applications. Example sequencing performance data for a given application are also summarized.

Technical reference for preparing Kinnex single-cell libraries with Parse Evercode Whole Transcriptome kits

Technical note – Preparing Kinnex single-cell libraries with Parse Evercode Whole Transcriptome kits ([102-326-610](#)) describes a modified workflow for preparing PacBio Kinnex single-cell libraries with single-cell or single-nucleus cDNA generated using Parse Biosciences Evercode Whole Transcriptome products

Technical note
PREPARING KINNEX™ SINGLE-CELL LIBRARIES WITH PARSE BIOSCIENCES EVERCODE™ WT KITS

Overview
This technical note describes the experimental conditions of a modified workflow for preparing PacBio® Kinnex single-cell libraries with single-cell or single-nucleus cDNA generated using Parse Biosciences Evercode Whole Transcriptome products, the WT and WT Mini kits. This workflow has not been tested using the Evercode WT Mega kit.


Required materials and equipment

- Refer to all materials described in the Kinnex single-cell protocol.¹
- Single-cell or single-nucleus cDNA from the Parse Biosciences Evercode WT or Evercode WT Mini kits.
 - Note that the cDNA needs to be generated as described in the Parse protocols² for either the Evercode WT or Evercode WT Mini up through the end of those protocols' section 2. Then use either the 15 ng (from the WT protocol) or 60–75 ng (from the WT Mini protocol) of the generated cDNA as input into the workflow at right (see Figure 1).
- Parse Biosciences-specific amplification primers (Appendix 1).

Library preparation
Library preparation will follow the PacBio Kinnex single-cell protocol¹ with modifications outlined below. Workflow modifications begin with single-cell or single-nucleus cDNA produced by the Parse Biosciences Evercode WT and WT Mini kits.

- Generate amplified cDNA per Parse Biosciences protocols. Complete the Evercode WT or Evercode WT Mini protocol up through section 2.
- Carry out QC as described in the Parse Biosciences user manual (see the figure at the end of section 2. Expected post-amplification sublibrary cDNA size distribution). While cDNA size ranges can vary, there should not be an excess of shorter (<500 bp) material.





PacBio Technical note
Preparing Kinnex single-cell libraries with Parse Evercode Whole Transcriptome kits ([102-326-610](#))
[Select to use either single-cell or single-nucleus cDNA from Parse Biosciences Evercode WT or Evercode WT Mini kits]


PacBio Procedure & checklist
Preparing Kinnex libraries using Kinnex single-cell RNA kit ([103-254-300](#))

Overview

- Parse whole transcriptome kits use split-pool combinatorial barcoding to enable scale-up of single cell projects to millions of cells or nuclei¹
- This Kinnex library prep workflow has been tested with single-cell or single-nucleus cDNA generated using Parse Biosciences Evercode WT and WT Mini kit whole transcriptome products
- This workflow has not been tested using the Evercode WT Mega kit

Workflow overview for preparing Kinnex single-cell libraries with Parse Evercode Whole Transcriptome kits

Follow Parse Bioscience whole transcriptome kit procedure¹ to generate single-cell cDNA samples suitable for Kinnex library prep and HiFi sequencing on PacBio long-read systems



Parse single-cell cDNA generation procedure overview¹

Technical note – Preparing Kinnex single-cell libraries with Parse Evercode Whole Transcriptome kits (102-326-610)


Technical note
PREPARING KINNEX® SINGLE-CELL LIBRARIES WITH PARSE BIOSCIENCES EVERCODE™ WT KITS

Parse single-cell cDNA generation procedure key steps

- Single-cell cDNA material first needs to be generated as described in the Parse protocols¹ for either the **Evercode WT** or **Evercode WT Mini kit** up through the end of those protocols' **Section 2**.

Evercode WT User Guides
[UMWT3300](#) or [UMWT3100](#)


- Use either **15 ng** (from the WT protocol) or **60–75 ng** (from the WT Mini protocol) of the generated cDNA as input into the **PacBio Kinnex single-cell library prep protocol** ([103-254-300](#)) starting at **Step 2: TSO PCR (with custom primers²)**.



¹ For single-cell cDNA generation procedure details, refer to **Evercode WT User Guides** ([UMWT3300](#) or [UMWT3100](#)) available from the [Parse Biosciences website](#).
² Perform **TSO PCR (Step 2)** with Parse Biosciences-specific amplification primers (see Appendix 1 in **Parse Technical note** [102-326-610](#)).

Workflow overview for preparing Kinnex single-cell libraries with Parse Evercode Whole Transcriptome kits (cont.)

Parse single-cell cDNA samples may be used for Kinnex library construction by following PacBio *Procedure & checklist – Preparing Kinnex libraries using Kinnex single-cell RNA kit* ([103-254-300](#))



PacBio Kinnex library construction procedure¹

Procedure & checklist – Preparing Kinnex libraries using Kinnex single-cell RNA kit ([103-254-300](#))



Kinnex SC library prep procedure key steps

- Follow **PacBio Kinnex single cell protocol** ([103-254-300](#)) (**Steps 2 – 7**) with modifications below:
 - Generate amplified single-cell cDNA per Parse Biosciences protocols. **Complete the Evercode WT or Evercode WT Mini protocol up through Section 2.**
 - Carry out cDNA quantification and sizing QC as described in the Parse Biosciences user manual.
 - Begin **Kinnex single-cell protocol from step 2 (TSO PCR)** with the following modification: Follow step 2.1 and use **5 µL of Parse Biosciences-specific amplification primers** as described in Appendix 1 **instead of** Kinnex 5' or 3' capture primer mixes.¹
 - After running TSO PCR program, proceed with the subsequent Kinnex workflow from cleanup with 1.5X SMRTbell cleanup beads in the TSO PCR step through the rest of the protocol.²

1. Parse Evercode single-cell cDNA generation

→ See **PacBio Technical note**
Preparing Kinnex single-cell libraries with Parse Evercode Whole Transcriptome kits ([102-326-610](#))

Kinnex single-cell library construction workflow (Steps 2 – 7)

¹ This replacement is critical to adapt the Parse library termini to function with the Kinnex concatenation workflow.
² Subsequent QC and library quantification can be carried out as described in the Kinnex single-cell protocol (step 7.21 in [103-254-300](#)).

TSO PCR amplification primer oligo recommendations for preparing Kinnex single-cell libraries with Parse Evercode Whole Transcriptome kits


Order custom TSO PCR amplification primers as described below from your preferred oligo vendor

Parse Biosciences-specific primers required for TSO PCR amplification step

- Order custom primers from your preferred oligo vendor
- An example of IDT code is shown at right
- Primers should be HPLC-purified
- Order these primers and reconstitute in TE buffer (10 mM Tris HCl pH 8, 0.1 mM EDTA).
- Use 10 μ M solutions of each of these primers in place of the yellow-capped Kinnex 3' capture primer mix or the red-capped Kinnex 5' capture primer mix in the Kinnex single-cell RNA kit (103-072-200)

Parse Kinnex forward capture primer


/5Biosg/AA/ideoxyU/CTACACGACGCTCTT
CCGATCTCAGACGTGTGCTCTTCCGATC



Parse Kinnex reverse capture primer

AAGCAGTGGTATCAACGCAGAG

5' TSO 3'
Parse Kinnex capture primer Rev




TSO PCR step in Kinnex single-cell RNA procedure uses a modified PCR primer (Parse Kinnex capture primer Fwd) to incorporate a biotin tag into desired cDNA products with the correct structure.¹

SMRT Link informatics and secondary analysis recommendations for Kinnex single-cell libraries generated with Parse Evercode Whole Transcriptome kits

Read Segmentation can be run in SMRT Link (or skera in the command line) to generate segmented reads

- Note: The SMRT Link single-cell Iso-Seq workflow does not currently support analysis of Parse Biosciences data.
- Instead, Read Segmentation can be run in SMRT Link (or skera in the command line) to generate segmented reads (S-reads).
- The S-reads will contain the original single-cell or single-nucleus cDNA-containing barcode information. Refer to Rebboah et al. (2021) (section Preprocessing of LR-Split-seq data)¹ and Parse Biosciences long-read data processing guide² (login required) for secondary analysis recommendations using command line and community tools.

SMRT Link read segmentation workflow






Preparing Kinnex single-cell libraries with Twist exome enrichment kits

Preparing Kinnex single-cell libraries with Twist exome enrichment kits: Getting started

Application-specific educational literature


Application-specific protocol documentation

Application-specific technical overviews



[PacBio literature website \[Link\]](#)

Application-specific brochures, informational guides and other product literature containing best practices recommendations for library preparation and data analysis workflows.




[Technical note – Preparing Kinnex single-cell libraries with Twist exome enrichment kits \(102-326-633\) \[PacBio\]](#)

Technical documentation containing a modified version of Twist Bioscience hybrid capture target enrichment protocol for generating enriched cDNA libraries for PacBio sequencing.



[Procedure & checklist – Procedure & checklist – Preparing Kinnex libraries using Kinnex single-cell RNA kit \(103-254-300\) \[PacBio\]](#)

Technical documentation containing PacBio SMRTbell library construction details.



[Technical overview: Kinnex library preparation using Kinnex single-cell RNA kit \(103-344-600\)](#)

Technical overview presentations describe sample preparation details for constructing HiFi libraries for specific applications. Example sequencing performance data for a given application are also summarized.

Technical reference for preparing Kinnex single-cell libraries with Twist exome enrichment kits

Technical note – *Preparing Kinnex single-cell libraries with Twist exome enrichment kits* ([102-326-633](#)) describes a modified workflow for preparing PacBio Kinnex single-cell libraries with Twist exome enrichment kits

Technical note
PREPARING KINNEX™ SINGLE-CELL LIBRARIES WITH TWIST EXOME ENRICHMENT KITS

Overview

This technical note describes the experimental conditions of a modified workflow preparing PacBio® Kinnex single-cell libraries with Twist exome enrichment kits. The purpose of using exome enrichment is to enrich for spliced mRNA in single-nuclei libraries¹. The following workflow has been tested on *10x Single Cell 3' Gene Expression* (v3.1) libraries, but is also compatible with single-nuclei libraries generated with the following:

- *10x Single Cell Gene Expression 3'* (v3.1, v4)
- *10x Single Cell Immune Profiling 5'* (v2, v3)
- Parse Evercode WT and WT mini kits (see the Parse technical note)

The workflow has been tested with the following Twist Exome Enrichment Panels:

- Twist Human Comprehensive Exome Panel
- Twist RNA Exome Panel


NOTE: This is a single reaction workflow. This protocol does not support pooled indexed libraries such as those generated using Twist's UDI Adapter Systems.

Required materials
 Please refer to the Kinnex single-cell protocol and respective single-cell provider manual for materials required.

Single-cell cDNA generation	
Single-nuclei cDNA library	Refer to single-cell vendor (10x, Parse)
10x cDNA primer*	10x 2000089
PacBio Kinnex library generation	
Kinnex single-cell RNA kit	Refer to Kinnex single-cell protocol
Custom 10x or Parse blockers†	Refer to Appendix 1
Twist exome enrichment protocol	
Twist Comprehensive Exome Panel	Twist 102031: 2 rxn Twist 102032: 12 rxn Twist 102033: 96 rxn
OR	
Twist RNA Exome Panel	Twist 107143: 2 rxn Twist 107144: 12 rxn Twist 107146: 96 rxn
Twist Universal Blockers	Twist 100856: 2 rxn Twist 100578: 12 rxn Twist 100767: 96 rxn
Twist Standard Hyb and Wash Kit v2 with Amp Mix	Twist 105559: 2 rxn Twist 105560: 12 rxn Twist 105561: 96 rxn

* The 10x cDNA primer is included in 10x Single Cell Gene Expression kits. If stock kit primers are not available, please refer to Appendix 1 for custom ordering.
 † Custom blockers are required to replace the Twist Universal Blockers to enable Kinnex compatibility. Please refer to Appendix 2 for custom blockers.

PacBio



PacBio Procedure & checklist
 Preparing Kinnex libraries using Kinnex single-cell RNA kit ([103-254-300](#))

PacBio Technical note
 Preparing Kinnex single-cell libraries with Twist exome enrichment kits ([102-326-633](#))
 [Select to use either Twist Human Comprehensive Exome Panel or Twist RNA Exome Panel]


PacBio Procedure & checklist
 Preparing Kinnex libraries using Kinnex single-cell RNA kit ([103-254-300](#))

Overview

- The purpose of using Twist exome enrichment is to enrich for spliced mRNA in single-nuclei libraries¹
- This workflow has been tested on *10x Single Cell 3' Gene Expression* (v3.1) libraries, but is also compatible with single-nuclei libraries generated with the following:
 - *10x Single Cell Gene Expression 3'* (v3.1, v4)
 - *10x Single Cell Immune Profiling 5'* (v2, v3)
 - Parse Biosciences Evercode WT and WT mini whole transcriptome kits (see Parse Technical note [102-326-610](#))
- The workflow has been tested with the following Twist Exome Enrichment Panels:
 - *Twist Human Comprehensive Exome Panel*
 - *Twist RNA Exome Panel*
- **NOTE:** This is a **single reaction** workflow. This protocol **does not** support pooled indexed libraries such as those generated using Twist's UDI Adapter Systems

Workflow overview for preparing Kinnex single-cell libraries with Twist exome enrichment kits

Follow PacBio Kinnex library prep protocol and modified Twist target enrichment procedure to generate single-cell hybrid capture libraries suitable for HiFi sequencing on PacBio long-read systems



PacBio Kinnex library construction procedure¹

Procedure & checklist – Preparing Kinnex libraries using Kinnex single-cell RNA kit ([103-254-300](#))

Preparing Kinnex[®] libraries using Kinnex single-cell RNA kit

Procedure & checklist

Before you begin

This procedure describes the workflow for constructing single-cell Kinnex libraries from 10x Chromium 3' v2.2 cDNA using the Kinnex single-cell RNA kit to target enrichment on the Kinnex Single-Cell Sequencing Kit (103-254-300).

Things to consider for use with single-cell cDNA generated using the 10x Chromium Next GEM Single Cell v2.2 kit or the Chromium Next GEM Single Cell v2.0 kit. For a full list of compatibility with this product, please consult the 10x website.

Parameter	Value
Sequencing kit	10x Genomics Next GEM Single Cell v2.2
Sequencing flow	10x Genomics Next GEM Single Cell v2.2
Number of single-cell cDNA per run	48,000 (cDNA per run) (10x Genomics)
Library construction kit	10x Genomics Next GEM Single Cell v2.2

Library prep

1. Using the 10x Genomics cDNA concentration and bead-to-cDNA ratio, add 10 µL of volume. See [10x Genomics](#) for details.

2. Add 10x Genomics Next GEM Single Cell v2.2

3. Add 10x Genomics Next GEM Single Cell v2.2

4. Add 10x Genomics Next GEM Single Cell v2.2

5. Add 10x Genomics Next GEM Single Cell v2.2

6. Add 10x Genomics Next GEM Single Cell v2.2

7. Add 10x Genomics Next GEM Single Cell v2.2

8. Add 10x Genomics Next GEM Single Cell v2.2

9. Add 10x Genomics Next GEM Single Cell v2.2

10. Add 10x Genomics Next GEM Single Cell v2.2

11. Add 10x Genomics Next GEM Single Cell v2.2

12. Add 10x Genomics Next GEM Single Cell v2.2

13. Add 10x Genomics Next GEM Single Cell v2.2

14. Add 10x Genomics Next GEM Single Cell v2.2

15. Add 10x Genomics Next GEM Single Cell v2.2

16. Add 10x Genomics Next GEM Single Cell v2.2

17. Add 10x Genomics Next GEM Single Cell v2.2

18. Add 10x Genomics Next GEM Single Cell v2.2

19. Add 10x Genomics Next GEM Single Cell v2.2

20. Add 10x Genomics Next GEM Single Cell v2.2

21. Add 10x Genomics Next GEM Single Cell v2.2

22. Add 10x Genomics Next GEM Single Cell v2.2

23. Add 10x Genomics Next GEM Single Cell v2.2

24. Add 10x Genomics Next GEM Single Cell v2.2

25. Add 10x Genomics Next GEM Single Cell v2.2

26. Add 10x Genomics Next GEM Single Cell v2.2

27. Add 10x Genomics Next GEM Single Cell v2.2

28. Add 10x Genomics Next GEM Single Cell v2.2

29. Add 10x Genomics Next GEM Single Cell v2.2

30. Add 10x Genomics Next GEM Single Cell v2.2

31. Add 10x Genomics Next GEM Single Cell v2.2

32. Add 10x Genomics Next GEM Single Cell v2.2

33. Add 10x Genomics Next GEM Single Cell v2.2

34. Add 10x Genomics Next GEM Single Cell v2.2

35. Add 10x Genomics Next GEM Single Cell v2.2

36. Add 10x Genomics Next GEM Single Cell v2.2

37. Add 10x Genomics Next GEM Single Cell v2.2

38. Add 10x Genomics Next GEM Single Cell v2.2

39. Add 10x Genomics Next GEM Single Cell v2.2

40. Add 10x Genomics Next GEM Single Cell v2.2

41. Add 10x Genomics Next GEM Single Cell v2.2

42. Add 10x Genomics Next GEM Single Cell v2.2

43. Add 10x Genomics Next GEM Single Cell v2.2

44. Add 10x Genomics Next GEM Single Cell v2.2

45. Add 10x Genomics Next GEM Single Cell v2.2

46. Add 10x Genomics Next GEM Single Cell v2.2

47. Add 10x Genomics Next GEM Single Cell v2.2

48. Add 10x Genomics Next GEM Single Cell v2.2

49. Add 10x Genomics Next GEM Single Cell v2.2

50. Add 10x Genomics Next GEM Single Cell v2.2

51. Add 10x Genomics Next GEM Single Cell v2.2

52. Add 10x Genomics Next GEM Single Cell v2.2

53. Add 10x Genomics Next GEM Single Cell v2.2

54. Add 10x Genomics Next GEM Single Cell v2.2

55. Add 10x Genomics Next GEM Single Cell v2.2

56. Add 10x Genomics Next GEM Single Cell v2.2

57. Add 10x Genomics Next GEM Single Cell v2.2

58. Add 10x Genomics Next GEM Single Cell v2.2

59. Add 10x Genomics Next GEM Single Cell v2.2

60. Add 10x Genomics Next GEM Single Cell v2.2

61. Add 10x Genomics Next GEM Single Cell v2.2

62. Add 10x Genomics Next GEM Single Cell v2.2

63. Add 10x Genomics Next GEM Single Cell v2.2

64. Add 10x Genomics Next GEM Single Cell v2.2

65. Add 10x Genomics Next GEM Single Cell v2.2

66. Add 10x Genomics Next GEM Single Cell v2.2

67. Add 10x Genomics Next GEM Single Cell v2.2

68. Add 10x Genomics Next GEM Single Cell v2.2

69. Add 10x Genomics Next GEM Single Cell v2.2

70. Add 10x Genomics Next GEM Single Cell v2.2

71. Add 10x Genomics Next GEM Single Cell v2.2

72. Add 10x Genomics Next GEM Single Cell v2.2

73. Add 10x Genomics Next GEM Single Cell v2.2

74. Add 10x Genomics Next GEM Single Cell v2.2

75. Add 10x Genomics Next GEM Single Cell v2.2

76. Add 10x Genomics Next GEM Single Cell v2.2

77. Add 10x Genomics Next GEM Single Cell v2.2

78. Add 10x Genomics Next GEM Single Cell v2.2

79. Add 10x Genomics Next GEM Single Cell v2.2

80. Add 10x Genomics Next GEM Single Cell v2.2

81. Add 10x Genomics Next GEM Single Cell v2.2

82. Add 10x Genomics Next GEM Single Cell v2.2

83. Add 10x Genomics Next GEM Single Cell v2.2

84. Add 10x Genomics Next GEM Single Cell v2.2

85. Add 10x Genomics Next GEM Single Cell v2.2

86. Add 10x Genomics Next GEM Single Cell v2.2

87. Add 10x Genomics Next GEM Single Cell v2.2

88. Add 10x Genomics Next GEM Single Cell v2.2

89. Add 10x Genomics Next GEM Single Cell v2.2

90. Add 10x Genomics Next GEM Single Cell v2.2

91. Add 10x Genomics Next GEM Single Cell v2.2

92. Add 10x Genomics Next GEM Single Cell v2.2

93. Add 10x Genomics Next GEM Single Cell v2.2

94. Add 10x Genomics Next GEM Single Cell v2.2

95. Add 10x Genomics Next GEM Single Cell v2.2

96. Add 10x Genomics Next GEM Single Cell v2.2

97. Add 10x Genomics Next GEM Single Cell v2.2


98. Add 10x Genomics Next GEM Single Cell v2.2

99. Add 10x Genomics Next GEM Single Cell v2.2

100. Add 10x Genomics Next GEM Single Cell v2.2

Kinnex SC library prep procedure key steps

- Follow Kinnex protocol instructions until the end of the TSO artifact removal step
 - If using compatible 10x products, follow the **Kinnex single-cell library protocol** ([103-254-300](#)) through the end of **Step 3** (TSO artifact removal)¹
- For **Step 4** (PCR to increase cDNA mass), this additional PCR step is required to generate sufficient cDNA mass for Twist exome enrichment (500 ng)
- After performing **Twist exome capture and on-bead PCR** (**Step 5**), proceed to the Kinnex PCR step of the **Kinnex single-cell library protocol** ([103-254-300](#)) and complete the remainder of the library construction workflow



PCR amplification primer and hybridization blocker oligo recommendations for preparing Kinnex single-cell libraries with Twist exome enrichment kits

For optimal library prep yields and on-target hybrid capture performance, order custom primers as described below from your preferred oligo vendor

10x cDNA primer¹ required for PCR amplification step to increase cDNA mass [Step 4]

- Order custom primers from your preferred oligo vendor.
- Standard desalt primers are acceptable.
- For each FWD and REV primer, add 10 µL of 10 µM to the PCR reaction to achieve a final concentration of 1 µM.

10x cDNA primer FWD

CTACACGACGCTCTTCCGATCT

10x cDNA primer REV

AAGCAGTGGTATCAACGCAGAG

Custom 10x or Parse blocker oligos² required for preparing samples for Twist hybridization capture reaction [Step 5.1]

- Order custom primers from your preferred oligo vendor.
- An example of IDT code is shown here.
- Blocker oligo primers should be HPLC-purified.
- Dilute blocker oligos to 10 µM and add 2 µL of each FWD and REV blocker.

10x Custom Blocker FWD

CTACACGACGCTCTTCCGATCT/3SpC3/

10x Custom Blocker REV

AAGCAGTGGTATCAACGCAGAG/3SpC3/

Parse Custom Blocker FWD

CTACACGACGCTCTTCCGATCTCAGACGTGTGCTCTTCCGATC/3SpC3/

Parse Custom Blocker REV

AAGCAGTGGTATCAACGCAGAG/3SpC3/



APPENDIX 3: SMRT Link run design procedure for Sequel IIe system

SMRT Link Run Design procedure for Sequel IIe system

Sample information and run information

- Select desired **Kinnex application** from the **Application** field drop-down menu
- The following fields are **auto-populated** with default recommended values and high-lighted in **green**:

- SMRTbell Adapter Design**

→ SMRTbell Kinnex Prep Kit

- Binding Kit**

→ Sequel II Binding Kit 3.2

- Sequencing Kit**

→ Sequel II Sequencing Plate 2.0 (4 rxn or 1 rxn)

- DNA Control Complex**

→ Sequel II DNA Internal Control Complex 3.2

- Movie Time per SMRT Cell**

→ 30 hrs

- Pre-Extension Time**

→ 2 hrs

SMRTbell Adapter Design field determines which adapter finding algorithm is used during post-primary analysis¹

Default SMRTbell adapter design for Kinnex samples is **SMRTbell Kinnex Prep Kit**

Recommended OPLC for Sequel II/IIe Kinnex library samples is **40 – 60 pM**

Recommended movie time = **30 hrs**

Select desired Kinnex application type from drop-down menu

SAMPLE 1: Kinnex full-length RNA library demo , A01, 30 hour movie, 16000 bp insert

Import from Sample Setup

Application Required: Kinnex full-length RNA

Well Sample Name Required: Kinnex full-length RNA library demo

Bio Sample Name

Sample Comment

Sample Well: A01

SMRTbell Adapter Design Required: SMRTbell® Kinnex Prep Kit

Binding Kit Required: Sequel® II Binding Kit 3.2

Sequencing Kit Required: Sequel® II Sequencing Plate 2.0 (4 rxn)

DNA Control Complex: Sequel® II DNA Internal Control Complex 3.2

Insert Size (bp) Required: 16000

Recommended Concentration on Plate (pM): 40 – 60 pM

On-Plate Loading Concentration (pM) Required: 45

Movie Time per SMRT Cell (hours): 30

Use Pre-Extension: YES NO

Pre-Extension Time (hours): 2

CCS Analysis will be performed on-instrument to produce HiFi .bam files.

Example sample information entered into a Sequel IIe system run design worksheet for a Kinnex full-length RNA library sample.

SMRT Link Run Design procedure for Sequel IIe system (cont.)

Advanced options

- For all Kinnex library samples, leave the following **Advanced Options** fields at their **default settings**
 - Use Adaptive Loading**
→ YES
 - Loading Target (P1 + P2)**
→ 0.85
 - Maximum Loading Time**
→ 2 hours
 - CCS Analysis Output - Include Low Quality Reads**
→ NO
 - CCS Analysis Output - Include Kinetics Information**
→ NO
 - Pre-Extension Time**
→ 2 hrs
- If desired, specify to use an alternative project folder for the **Add Data to Project** field

Advanced Options

Use Adaptive Loading YES NO

Loading Target (P1 + P2) 0.85

Maximum Loading Time (hours) 2

CCS Analysis Output - Include Low Quality Reads YES NO

CCS Analysis Output - Include Kinetics Information YES NO

Add Data to Project

Example default Advanced Options settings entered into a Sequel IIe system run design worksheet for a Kinnex full-length RNA library sample.

SMRT Link Run Design procedure for Sequel Ii system (cont.)

Barcoded sample options

- For Kinnex library samples, can leave most **Barcoded Sample Options** fields at their **default settings**

Specify Bio Sample Names, either interactively or by downloading a CSV file (**Interactively** or **From a file**)

If desired, specify to perform barcode demultiplexing on-instrument or in SMRT Link (default = **On-instrument** for Sequel Ii system)

Can leave most of these fields at their default values

Barcoded Sample Options

Sample Is Barcoded YES NO

Barcode Set Required MAS SMRTbell barcoded adapters (v2)

Same Barcodes on Both Ends of Sequence ? YES NO

Assign Bio Sample Names to Barcodes ? Required

Demultiplex Barcodes ON INSTRUMENT IN SMRT LINK DO NOT GENERATE

Example default Barcoded Sample Options settings entered into a Sequel Ii system run design worksheet for a Kinnex full-length RNA library sample.



www.pacb.com

Revision history (description)	Version	Date
Added procedural summary illustrations to clarify Kinnex library prep protocol steps, updated example Kinnex sequencing performance data for PacBio long-read systems (Revio system + SPRQ chemistry and Vega system) and incorporated information about PacBio compatible RNA sequencing workflows using third-party kit products.	02	December 2025

Research use only. Not for use in diagnostic procedures. © 2025 Pacific Biosciences of California, Inc. (“PacBio”). All rights reserved. Information in this document is subject to change without notice. PacBio assumes no responsibility for any errors or omissions in this document. Certain notices, terms, conditions and/or use restrictions may pertain to your use of PacBio products and/or third-party products. Refer to the applicable PacBio terms and conditions of sale and to the applicable license terms at pacb.com/license. Pacific Biosciences, the PacBio logo, PacBio, Circulomics, Omniome, SMRT, SMRTbell, Iso-Seq, Sequel, Nanobind, SBB, Revio, Onso, Apton, Kinnex, PureTarget, SPRQ, and Vega are trademarks of PacBio.