



Technical overview: Adeno-associated virus (AAV) library preparation using SMRTbell prep kit 3.0

Sequel II and IIe systems ICS v11.0 / SMRT Link v11.0

PN 102-390-400 Version 01 (April 2022)

Adeno-associated virus (AAV) library preparation using SMRTbell prep kit 3.0

Technical overview

1. PacBio sequencing for quality control of gene therapy methods
2. AAV library sample preparation workflow overview
3. AAV library sample preparation workflow details
4. AAV library sequencing preparation workflow overview
5. AAV data analysis workflow overview
6. AAV library example performance data
7. Technical documentation & applications support resources



AAV Sequencing: How to get started

Application-specific educational literature

Application-specific Procedure & checklist


Application-specific technical overviews

Library construction, sequencing & analysis




PacBio literature website ([Link](#))

Lists application-specific brochures, informational guides and and other product literature containing best practices recommendations for sample preparation and data analysis workflows.




Procedure & Checklist – Preparing multiplexed AAV SMRTbell libraries using SMRTbell prep kit 3.0 ([102-126-400](#))

Technical documentation containing sample library construction and sequencing preparation protocol details



Technical Overview: AAV Library preparation using SMRTbell prep kit 3.0 ([102-390-400](#))

Technical overview presentations describe sample preparation details for constructing HiFi libraries for specific applications. Example sequencing performance data for a given application are also summarized.






PacBio sequencing for quality control of gene therapy methods

PacBio sequencing for quality control of gene therapy methods

PacBio is the only sequencing technology able to support a broad range of gene therapy use cases



High accuracy and long read lengths enables a broad range of gene therapy methods



Discovery & Engineering



Vector QC



Editing Efficiency



Safety

Adeno-associated virus (AAV) for gene therapy


Recombinant adeno-associated viral (rAAV) vectors are promising tools for gene therapy

AAV structure is comprised of a non-enveloped capsid protein shell surrounding and protecting a small, single-stranded DNA genome of ~4.8 kb


Why Use AAV as a Vector

- Ability to generate recombinant AAV (rAAV)
- Form episomal concatemers in the cell nucleus
- Reduced host genome integration
- Proven in various therapeutic applications
- Low immunogenicity
- Strong vector persistence
- Long-term transferred gene (transgene) expression
- Non-toxic

A



B



Adeno-Associated Virus (AAV) structure. **A.** Cartoon illustration showing simulated AAV size and 3-D structure. **B.** Electron microscope image of purified AAV vector particles.

Recombinant AAV vector genome design

Single-stranded AAV (ssAAV) and self-complementary AAV (scAAV)

- Wild-type AAV genome consists of the viral *rep* and *cap* genes between two **inverted terminal repeats (ITRs)**
- Recombinant AAV vector is produced by co-transfection of
 1. **Vector plasmid** containing an ITR-flanked transgene cassette;
 2. **Packaging plasmid** that encodes the *rep* and *cap* genes of a specific AAV serotype; and
 3. **Helper plasmid** that supplies the essential adenovirus helper genes


The vector plasmid may encode for:

- **Single-stranded DNA (ssAAV)** containing ~4.5 kb of novel transgene sequence; or
- **Self-complementary DNA (scAAV)** containing up to ~2.2 kb of novel transgene sequence in duplex form

Wild-type AAV



Recombinant AAV




Recombinant AAV Vector QC Example


PacBio HiFi reads allow for complete sequencing of viral vector particles. This allows for unambiguous detection of undesired plasmid or host gDNA that has been encapsidated

Alignments of heterogeneous vector populations to the pcis-plasmid reference

Vector QC is central to understanding that the vectors are packaging the **correct** DNA





Undesirable reads containing plasmid backbone sequences



Characterization of chimeric reads that map to the human genome

Undesirable reads containing human host sequences





AAV library sample preparation workflow overview

AAV sample preparation procedure description

Procedure & checklist – Preparing multiplexed adeno-associated virus (AAV) libraries using SMRTbell prep kit 3.0 ([102-126-400](#)) describes a library preparation procedure for sequencing multiplexed AAV samples on the Sequel II and IIe Systems



SMRTbell Prep Kit 3.0
(102-182-700)



SMRTbell Barcoded Adapter Plate 3.0
(102-009-200)

Procedure & checklist highlights

- Compatible with single-stranded AAV (**ssAAV**) and self-complementary AAV (**scAAV**)
- AAV DNA sample extraction is performed using third-party methods
- Multiplexing is performed using **SMRTbell barcoded adapter plate 3.0** (102-009-200)
- A total of **1 µg of pooled AAV DNA** is required for SMRTbell library preparation. The input DNA requirement per individual AAV sample depends on the multiplex level and can range from **250 ng per sample for a 4-plex**, to **42 ng per sample for a 24-plex** preparation.

Preparing multiplexed AAV SMRTbell libraries using SMRTbell prep kit 3.0

PacBio

Procedure & checklist

Before you begin

This protocol describes how to prepare multiplexed AAV libraries for sequencing on the Sequel®, Sequel II and Sequel IIe Systems. Multiplexing achieved using SMRTbell Barcoded adapters.

Once barcoded, samples can be pooled and purified for sequencing on the Sequel, Sequel II and Sequel IIe Systems. A total of 1 µg of pooled AAV DNA is required for SMRTbell library preparation. The input requirements per individual AAV sample depends on the multiplex level and can range from 250 ng for a 4-plex, to 42 ng per sample for a 24-plex preparation.

Preparation of AAV DNA for SMRTbell library preparation

AAV DNA may have several structures or formats: single-stranded DNA (ss), double-stranded DNA with some non-complementary regions (ds), or single-stranded self-complementary (sc) DNA, with an inverted repeat in the middle. All structures generally have inverted terminal repeats (ITR) at the ends.

SMRTbell template formation requires double-stranded DNA with blunt ends onto which hairpin SMRTbell adapters can be ligated, following A-tailing. Creating appropriate ds input DNA requires extra steps if the sample contains a substantial ss region without ds blunt ends, as shown on the left in Figure 1, below.

If the sample contains self-complementary molecules that have only one site accessible for hairpin adapter ligation, as shown on the right in Figure 1, it is important to collect and analyze data appropriately. Non-complementary ds molecules also require additional analysis.




Figure 1: Structures of ss (single-stranded) and sc (self-complementary) DNA molecules. The ssDNA structure at left requires a second strand (as shown in yellow) for SMRTbell formation.

© 2022 PacBio. All rights reserved. Research use only. Not for use in diagnostic procedures.
PN 102-126-400 EA V1 APR2022

PacBio

PacBio [Documentation](#) ([102-126-400](#))

APPLICATIONS
VIRAL SEQUENCING
AAV Sequencing



AAV sample preparation & sequencing workflow overview

Workflow summary for constructing SMRTbell libraries suitable for sequencing on the Sequel and Sequel II/IIe Systems for multiplexed AAV applications




AAV DNA extraction & QC

- Perform AAV DNA sample extraction using third-party methods (Samples can be comprised of a mixture of ssAAV and scAAV species)
- ≥ 1 μg of pooled AAV DNA is required per SMRT Cell 8M
- Input DNA requirement per AAV sample can range from 250 ng for a 4-plex, to 42 ng per sample for a 24-plex.



Sequencing

- Follow SMRT Link Sample Setup instructions for primer annealing, polymerase binding, complex cleanup and sample loading




HiFi Read

PacBio HiFi reads
achieve 99.9% accuracy



Data analysis

- Perform CCS analysis and de-multiplex barcodes on-instrument (Sequel IIe system only) or in SMRT Link
- Align demultiplexed AAV data against a reference and visualize results using command line tools available in PacBio GitHub




SMRTbell library construction

- *Procedure & checklist – Preparing multiplexed adeno-associated virus (AAV) libraries using SMRTbell prep kit 3.0 (102-126-400)*
- Multiplex AAV samples using SMRTbell Barcoded Adapter Plate 3.0 (102-009-200).
- Pool barcoded AAV samples and purify final multiplexed SMRTbell library using SMRTbell cleanup beads



AAV DNA extraction for SMRTbell library preparation

- Refer to **Procedure & checklist – Preparing multiplexed adeno-associated virus (AAV) libraries using SMRTbell prep kit 3.0** ([102-126-400](#)) for recommendations on third-party AAV DNA extraction methods
- AAV DNA may have several structures or formats:
 - Single-stranded (ss) DNA;
 - Double-stranded (ds) DNA with some non-complementary regions; or
 - Single-stranded self-complementary (sc) DNA, with an inverted repeat in the middle
- All structures generally have inverted terminal repeats (ITR) at the ends (some may have mutated ITRs or other short hairpin sequences present)
- Note:** SMRTbell template synthesis requires **dsDNA** with blunt ends onto which hairpin SMRTbell adapters can be ligated
 - Creating appropriate dsDNA input for SMRTbell library construction requires an extra thermal annealing step if the sample contains substantial ss regions without ds blunt ends, as shown on the left in the figure for ssAAV samples



Structures of single-stranded (ss) (left) and self-complementary (sc) (right) AAV DNA molecules. The ssDNA structure at left requires a second strand (as shown in yellow) for SMRTbell library construction.

AAV SMRTbell library construction workflow overview

DNA repair & A-tailing


- After extracting AAV DNA, perform DNA damage repair, polish fragment ends and carry out A-tailing reaction using SMRTbell prep kit 3.0 (102-182-700)

Barcoded adapter ligation & cleanup

- Use SMRTbell barcoded adapter plate (102-009-200) to prepare multiplexed AAV samples

Nuclease treatment & cleanup

- Perform nuclease treatment and cleanup to remove incomplete SMRTbell templates/un-ligated adapters and obtain a purified AAV library sample ready for pooling & sequencing



Structure of AAV SMRTbell library template molecules

Comparison of standard SMRTbell library template structure *versus* ssAAV SMRTbell templates and scAAV SMRTbell templates


Standard SMRTbell template

- Linear dsDNA insert molecules are ligated to SMRTbell adapters **at each end** during library construction




Single-stranded AAV (ssAAV) SMRTbell template

- ssAAV molecules can be converted to dsDNA structures by performing intermolecular annealing or 2nd-strand synthesis
- Forward (+) and reverse (–) sequences are on separate strands, bracketed by T-shaped inverted terminal repeat (ITR) sequences and ligated to SMRTbell adapters **at each end** during library construction



Self-complementary AAV (scAAV) SMRTbell template

- scAAV molecules feature a linear self-complementary sequence containing a mutant ITR (mITR) (or short hairpin DNA, shDNA) in the middle of the molecule
- The two ITRs at the open end of the scAAV vector genome are ligated to a **single** SMRTbell adapter





AAV library sample preparation workflow details

Procedure & checklist – Preparing multiplexed AAV SMRTbell libraries using SMRTbell prep kit 3.0 (102-126-400)

Procedure & checklist [102-126-400](#) describes a method for constructing SMRTbell libraries using SMRTbell prep kit 3.0 (SPK 3.0) and SMRTbell barcoded adapter plate 3.0 that are suitable for generating HiFi reads on the Sequel II and IIe Systems for AAV sequencing applications

Procedure & checklist contents

1. AAV DNA input requirements per SMRT Cell 8M.
2. Recommendations for AAV DNA extraction and QC.
3. Recommendations for thermal annealing of ssAAV vectors.
4. Enzymatic workflow steps for preparation of multiplexed AAV SMRTbell libraries using SMRTbell prep kit 3.0 (102-182-700) and SMRTbell barcoded adapter plate 3.0 (102-009-200).
5. Guidance for pooling barcoded AAV SMRTbell libraries for multiplexed sequencing on a single SMRT Cell 8M.

Preparing multiplexed AAV SMRTbell libraries using SMRTbell prep kit 3.0

PacBio

Procedure & checklist

Before you begin

This protocol describes how to prepare multiplexed AAV libraries for sequencing on the Sequel®, Sequel II and Sequel IIe Systems. Multiplexing achieved using SMRTbell Barcoded adapters.

Once barcoded, samples can be pooled and purified for sequencing on the Sequel, Sequel II and Sequel IIe Systems. A total of 1 µg of pooled AAV DNA is required for SMRTbell library preparation. The input requirements per individual AAV sample depends on the multiplex level and can range from 250 ng for a 4-plex, to 42 ng per sample for a 24-plex preparation.

Preparation of AAV DNA for SMRTbell library preparation

AAV DNA may have several structures or formats: single-stranded DNA (ss), double-stranded DNA with some non-complementary regions (ds), or single-stranded self-complementary (sc) DNA, with an inverted repeat in the middle. All structures generally have inverted terminal repeats (ITR) at the ends.

SMRTbell template formation requires double-stranded DNA with blunt ends onto which hairpin SMRTbell adapters can be ligated, following A-tailing. Creating appropriate ds input DNA requires extra steps if the sample contains a substantial ss region without ds blunt ends, as shown on the left in Figure 1, below.

If the sample contains self-complementary molecules that have only one site accessible for hairpin adapter ligation, as shown on the right in Figure 1, it is important to collect and analyze data appropriately. Non-complementary ds molecules also require additional analysis.




Figure 1: Structures of ss (single-stranded) and sc (self-complementary) DNA molecules. The ssDNA structure at left requires a second strand (as shown in yellow) for SMRTbell formation.



© 2022 PacBio. All rights reserved. Research use only. Not for use in diagnostic procedures.
PN 102-126-400 EA V1 APR2022

PacBio

PacBio [Documentation](#) (102-126-400)

AAV SMRTbell construction workflow overview

**** Note:** Extracted AAV DNA samples typically contain a **mixture** of both ssAAV DNA and scAAV DNA



Workflow Step

DNase I treatment

AAV DNA extraction

Thermal annealing for
ssAAV vectors

Normalize DNA input
& cleanup

Repair & A-tailing

Adapter ligation & cleanup

Nuclease treatment &
cleanup

Pool libraries & concentrate

DNase I treatment

DNase I treatment removes any non-encapsulated DNA from the vector

- We recommended using DNase I from NEB (M0303S) to treat extracted AAV DNA samples
- The capsid shell is resistant to DNase I treatment, therefore, DNase I will not degrade the encapsulated DNA.
- Add 20 units of DNase I to 7×10^{11} vector genomes (in 50 μL containing 1.5×10^{13} genome copies/mL)

Reagent	Volume (μL)
AAV vector	50
Nuclease-free Water	120
Reaction Buffer	20
DNase I (RNase-free)	10
<hr/>	
Total Volume	200



Workflow Step

DNase I treatment



AAV DNA extraction



Thermal annealing for ssAAV vectors



Normalize DNA input & cleanup



Repair & A-tailing



Adapter ligation & cleanup



Nuclease treatment & cleanup




Pool libraries & concentrate



rAAV DNA extraction

Extract DNA from DNase I treated sc/ss AAV vector by using the PureLink Viral RNA/DNA Mini Kit following the manufacturer's instructions

- **PureLink Viral RNA/DNA Mini kit** from Invitrogen (12280050) is recommended for AAV DNA extraction



- PureLink kit allows efficient lysis of viral particles at elevated temperatures and selective binding of viral nucleic acids to the silica matrix under highly denaturing condition (~1.5 hours)
- PureLink kit **removes fragments <200 bp**
 - If retention of DNA fragments <200 bp is desired, use an alternative extraction method such as phenol / chloroform / isoamyl alcohol (see Tran *et al.*)
- Perform DNA quantification QC on the extracted AAV DNA using a Qubit DNA assay kit and DNA sizing QC using a Bioanalyzer tool
- **Note:** Most extracted AAV DNA preparations typically contain a **mixture** of ssAAV and scAAV DNA



Example Bioanalyzer sizing QC analysis for an extracted scAAV DNA sample. The major peak at ~2 kb corresponds to the expected size range (1.8 kb) of this scAAV DNA sample.

Workflow Step

DNase I treatment



AAV DNA extraction



Thermal annealing for
ssAAV vectors



Normalize DNA input
& cleanup



Repair & A-tailing



Adapter ligation & cleanup



Nuclease treatment &
cleanup





Pool libraries & concentrate

Thermal annealing for ssAAV vectors

Perform thermal annealing to anneal + / – strands in ssAAV DNA to convert ssAAV vectors into a double-stranded DNA form suitable for SMRTbell library construction

- Thermal annealing step is only required for ssAAV vectors
 - Skip this step if working with scAAV vectors.
- To perform thermal annealing:
 - Add 1 vol. of ssAAV DNA sample to 9 vol. of 1X annealing buffer stock = 25 mM NaCl, 10 mM Tris-HCl [pH 8.5], 0.5 mM EDTA [pH 8]
 - Heat the diluted ssAAV DNA sample (prepared above) at 95°C for 5 min and then cool to 25°C (1 min for every –1°C) on a thermocycler
- Perform DNA quantification QC on the thermally annealed ssAAV DNA using a Qubit DNA assay kit and DNA sizing QC using a Bioanalyzer tool

ssAAV



Thermal annealing



Double-stranded DNA structure
(may contain non-complementary strands)

if your sample contains substantial single-stranded DNA regions without double-stranded blunt ends as shown above for ssAAV, perform a thermal annealing step to create dsDNA suitable for SMRTbell library construction.

Workflow Step

DNase I treatment



AAV DNA extraction



Thermal annealing for
ssAAV vectors



Normalize DNA input
& cleanup



Repair & A-tailing



Adapter ligation & cleanup



Nuclease treatment &
cleanup



Pool libraries & concentrate

Normalize DNA input & cleanup

Normalize the amount of AAV DNA for each sample to go into SMRTbell library construction to generate more even sequencing coverage across multiplexed samples

- Use a **combined total of $\geq 1 \mu\text{g}$** of AAV DNA per SMRT Cell 8M
 - This is the **combined total mass** of samples that will be multiplexed
- **Per sample input amounts should be $1 \mu\text{g}$ / number of samples**
 - Use an equal mass of input DNA per AAV sample for multiplexing
 - Recommended per sample input AAV DNA amounts for different multiplex levels are shown in the table below

Multiplex level	DNA Input
4-plex	250 ng per sample
8-plex	125 ng per sample
12-plex	83 ng per sample
16-plex	63 ng per sample
24-plex	42 ng per sample

- Perform clean-up using 1.3X SMRTbell cleanup beads to concentrate the sample in an elution volume of $47 \mu\text{L}$

Workflow Step

DNase I treatment



AAV DNA extraction



Thermal annealing for
ssAAV vectors



Normalize DNA input
& cleanup



Repair & A-tailing



Adapter ligation & cleanup



Nuclease treatment &
cleanup



Pool libraries & concentrate

Repair & A-tailing

Perform DNA damage repair, end repair and A-tailing reactions in a single tube

- For **multiplexed** preps:
 - Prepare a reagent master mix by adding the required components in the order and volume listed in the table below in a new tube
 - Adjust the buffer and enzyme component volumes for the total number of samples being prepared, plus 10% overage
- For **non-multiplexed (individual)** preps:
 - Add components directly to the sample from the previous step at the specified volumes and skip the preparation of a reagent master mix

✓	Tube	Reagent	Volume
	Previous	Contents from previous step	46 µL
	Purple	Repair buffer	8 µL
	Blue	End repair mix	4 µL
	Green	DNA repair mix	2 µL
Total volume			60 µL



SMRTbell Prep Kit 3.0
(102-182-700)

- Incubate at 37°C for 30 min followed by 65°C for 5 min

Workflow Step

DNase I treatment



AAV DNA extraction



Thermal annealing for
ssAAV vectors



Normalize DNA input
& cleanup



Repair & A-tailing



Adapter ligation & cleanup



Nuclease treatment &
cleanup



Pool libraries & concentrate

Adapter ligation & cleanup

Ligate either barcoded or non-barcoded SMRTbell adapters to repaired and A-tailed AAV dsDNA to create SMRTbell templates

- For **multiplexed** preps:

- Use **barcoded adapters** included with SMRTbell barcoded adapter plate 3.0 (102-009-200)

SMRTbell Barcoded Adapter Plate 3.0 (102-009-200) contains 96 ready-to-use barcoded SMRTbell adapters



- For **non-multiplexed (individual)** preps:

- Use the **standard (non-barcoded) adapter** included with SMRTbell prep kit 3.0 (102-182-700)

SMRTbell Prep Kit 3.0 (102-182-700) contains the standard (non-barcoded) SMRTbell adapter



- Add 4 μ L of adapter per sample, followed by the components listed in the table at right
- Incubate at 20°C for 30 min
- Perform clean-up using 1.3X SMRTbell cleanup beads

Reaction Mix 2 (RM2)

✓	Tube	Component	Volume per Sample
	Yellow	Ligation mix	30 μ L
	Red	Ligation enhancer	1 μ L
Total volume			31 μ L

Workflow Step

DNase I treatment



AAV DNA extraction



Thermal annealing for
ssAAV vectors



Normalize DNA input
& cleanup



Repair & A-tailing



Adapter ligation & cleanup



Nuclease treatment &
cleanup



Pool libraries & concentrate

Nuclease treatment & cleanup

Perform a nuclease treatment step to reduce levels of undesirable reaction products that are not completely intact, circular SMRTbell templates

- For **multiplexed** preps:
 - Prepare a reagent master mix by adding the required components in the order and volume listed in the table below in a new tube
 - Adjust the buffer and enzyme component volumes for the total number of samples being prepared, plus 10% overage
- For **non-multiplexed (individual)** preps:
 - Add components directly to the sample from the previous step at the specified volumes and skip the preparation of a reagent master mix

Reaction Mix 3 (RM3)			
✓	Tube	Component	Volume per sample
	Light purple	Nuclease Buffer	5 µL
	Light green	Nuclease Mix	5 µL
Total Volume			10 µL



SMRTbell Prep Kit 3.0
(102-182-700)

- Incubate at 37°C for 15 min
- Perform clean-up using 1.3X SMRTbell cleanup beads
- OPTIONAL: Take 1 µL and measure DNA concentration with the Qubit dsDNA HS assay kit to check for variable sample loss prior to pooling barcoded samples

Workflow Step

DNase I treatment



AAV DNA extraction



Thermal annealing for
ssAAV vectors



Normalize DNA input
& cleanup



Repair & A-tailing



Adapter ligation & cleanup



Nuclease treatment &
cleanup




Pool libraries & concentrate

Pool barcoded AAV samples & concentrate

Pool barcoded AAV SMRTbell library samples for sequencing on a single SMRT Cell

- Pooling libraries that have **similar average insert sizes** will provide more even sequence coverage between samples
 - Increasing differences in average molecular weights will lead to increasing differences in coverage across libraries because of differences in molarity for SMRT Cell loading
 - For sensitive experiments where balanced coverage is critical, consider pooling libraries in equal molar amounts
- 1.3X cleanup using SMRTbell cleanup beads to concentrate
- Perform DNA quantification QC on the final pooled AAV library using a Qubit DNA assay kit and DNA sizing QC using a Bioanalyzer tool




Example Bioanalyzer sizing QC analysis for a pooled 8-plex scAAV library sample. The major peak at ~2.6 kb corresponds to the expected size range (2.1 kb) of this scAAV DNA sample.



AAV library sequencing preparation workflow overview

Sample Setup & Run Design recommendations for AAV libraries

In SMRT Link Sample Setup & Run Design, select 'Viral Sequencing' / 'Adeno-Associated Virus' for application type



Binding kit 3.1 & cleanup beads (102-333-400) is recommended for preparing AAV samples for sequencing.

- We recommend using **Sequel II binding kit 3.1 & cleanup beads** (102-333-400) to perform ABC (anneal primer / bind polymerase / clean up complex) with AAV samples
- Refer to **Quick reference card – Loading and pre-extension time recommendations for the Sequel II and Ile systems** ([101-769-100](#)) for updates to ABC workflow for specific applications

Sequel II binding kit 3.1 & cleanup beads (102-333-400) includes the following components:

- Sequencing primer 3.1
- Sequel II polymerase 2.1
- SMRTbell cleanup beads for complex cleanup
- DNA internal control 3.1 (defined 2 kb template bound to Polymerase 2.1)
- Supports ≥ 24 binding reactions, and up to 4 SMRT Cells 8M per binding reaction (96 cells total), depending on use case, sample size and concentration

AAV Sample Setup guidance

Use SMRT Link Sample Setup High-Throughput (HT) mode and follow instructions to perform ABC (anneal primer / bind polymerase / clean up complex) using recommended settings for AAV samples

PacBio

Sample Setup

smark (Lab Tech)

?

Sample Setup

New Calculation

Edit/Print

Import

Export

Delete

Version:

Classic

High-Throughput

Sample Setup HT for Sequel II and Sequel IIe

Displaying rows 1 to 1 out of 1

AAV

<input type="checkbox"/>	Name	Date Created ↓	Number of Samples	Comment	Created By	Locked
<input type="checkbox"/>	Example AAV Sample Setup HT m...	2022-04-12, 01:31:07 PM	2	This batch includes Pooled_AAV_Sample_01 Pooled_AAV_Sample_02	smark	false

- Sample Setup High-Throughput mode provides a simplified, streamlined workflow to efficiently process either one sample or multiple samples with similar library properties (such as mean insert size and DNA concentration) in parallel
- You can also export the calculated values to a CSV file for laboratory automation

Note: We recommend using an on-plate loading concentration (OPLC) of 70 – 150 pM for AAV samples

< Sample Group >

Actions

Copy Remove Lock Automate

Name

Example AAV Sample Setup

Comment

This batch includes Pooled_AAV_Sample_01 Pooled_AAV_Sample_02

Application

Adeno-Associated Virus

Binding Kit

Sequel II Binding Kit 3.1

Number of Samples

2

samples

SMRT Cells per Sample

1

cells

Available Volume per Sample

10

uL

Insert Size

3500

bp

Sample Concentration

12

ng/uL

Cleanup Anticipated Yield

60

%

Recommended Concentration on Plate

40-150

pM

Specify Concentration on Plate

90

pM

Minimum Pipetting Volume

1

uL

Warnings

Example Sample Setup HT mode worksheet for a batch comprised of two AAV samples.

AAV Run Design guidance

Follow SMRT Link Run Design instructions to set up a sequencing run using recommended settings for AAV samples

- **IMPORTANT:** Select **Adeno-Associated Virus** from the Application field drop-down menu in SMRT Link Run Design
- **Modified adapter calling + heteroduplex detection (HD)** are automatically enabled by default when AAV application type is selected
- Specify **YES** for '**Detect and Resolve Heteroduplex Reads**' regardless of whether you sample contains ssAAV, scAAV, or a mixture of ssAAV and scAAV DNA

IMPORTANT: AAV application type specification in Run Design is **required** for proper analysis of AAV data. If inadvertently not specified, the AAV sample must be re-run.

IMPORTANT: Heteroduplex detection is **required** for proper analysis of AAV data. If inadvertently not specified in Run Design, the AAV sample must be **re-run**.

The screenshot displays the PacBio Run Design interface. The 'Run Information' section on the left includes fields for System Type (SEQUEL II selected), Run Name (Example_AAV_Run_Design), Run Comments, Experiment Name, Experiment ID, and Estimated Run Duration. The 'Sample Information' section on the right shows a dropdown for 'SAMPLE 1: Pooled_AAV_Sample_01, A01, 24 hour movie, 3500 bp insert'. Below this, the 'Import from Sample Setup (Classic only)' section has a 'Select Sample' button. The 'Application' field is set to 'Adeno-Associated Virus'. Other fields include Well Sample Name (Pooled_AAV_Sample_01), Bio Sample Name, Sample Comment, Sample Well (A01), Template Prep Kit (SMRTbell® Prep Kit 3.0), Binding Kit (Sequel® II Binding Kit 3.1), Sequencing Kit (Sequel® II Sequencing Plate 2.0 (4 rxn)), DNA Control Complex (Sequel® II DNA Internal Control Complex 3.1), Insert Size (bp) (3500), Recommended Concentration on Plate (pM) (40-150 pM), On-Plate Loading Concentration (pM) (90), Movie Time per SMRT Cell (hours) (24), Use Pre-Extension (YES selected), Pre-Extension Time (hours) (2), and Detect and Resolve Heteroduplex Reads (YES selected). A note at the bottom states: 'CCS Analysis will be performed on-instrument to produce HiFi .bam files.'

Example sample information entered into Run Design for sequencing an AAV sample.

AAV Run Design guidance (cont.)

OPTIONAL: Run Design setup procedure for automated barcode demultiplexing of pooled AAV library samples

1. Sample is Barcoded: YES
2. Barcode Set: Select 'SMRTbell barcoded adapter plate 3.0 (bc2001-bc2096)'
3. Same Barcodes on Both Ends of Sequence: YES
4. Assign a Biological Sample Name to each barcoded sample using one of two ways: From a (CSV) File or Interactively
5. Specify if barcode demultiplexing is to be performed on-instrument (Sequel IIe system only) or in SMRT Link. (Optionally specify Do Not Generate.)

The screenshot shows the PacBio Run Design web interface. The 'Sample Information' section is active, and the 'Advanced Options' tab is selected. Within this tab, the 'Barcoded Sample Options' are highlighted with a red box and numbered 1 through 5:

- 1. Sample Is Barcoded: YES (selected)
- 2. Barcode Set: SMRTbell Barcoded Adapter Plate 3.0 (Required)
- 3. Same Barcodes on Both Ends of Sequence: YES (selected)
- 4. Assign Bio Sample Names to Barcodes: Interactively (selected)
- 5. Demultiplex Barcodes: ON INSTRUMENT (selected)

Below the highlighted options, there are 'Copy' and 'Delete' buttons, and an 'Add Sample' button at the bottom. The left sidebar shows 'Run Information' with fields for Run Name, Run Comments, Experiment Name, and Experiment ID, and a list of 'Run Reagents / Consumables'.

Example barcoding information entered into Run Design for sequencing a pooled AAV sample.




AAV data analysis workflow overview

AAV data analysis workflow overview


After sequencing AAV samples using the AAV sequencing mode in Run Design, perform AAV data analysis QC using command line tools available in GitHub

Workflow Step




AAV application type specification in Run Design is IMPORTANT!

AAV sequencing mode in Run Design applies a modified adapter calling algorithm during post-primary analysis and enables automated heteroduplex (HD) read detection and resolution during CCS analysis




Modified adapter calling is automatically performed for AAV samples containing scAAV DNA (or a mixture of scAAV and ssAAV*)

AAV sequencing mode in Run Design automatically applies a modified adapter calling algorithm to enable correct annotation of full-pass subreads for library molecules containing only one SMRTbell adapter



Heteroduplex detection and splitting is automatically performed for AAV samples containing ssAAV DNA (or a mixture of ssAAV and scAAV)


AAV sequencing mode in Run Design automatically applies a heteroduplex detection (HD) and splitting algorithm during CCS analysis of library molecules containing non-complementary regions



AAV data analysis QC recommendations

Map AAV HiFi reads to AAV reference genome using pbmm2 command line tool available in GitHub

Workflow Step




Map HiFi Reads to AAV Genome

- Example **pbmm2** command:
 - `pbmm2 -preset isoseq -sort -j <threads> <aav_genome.fa> <reads.bam> <mapped.sorted.bam>`
 - We use the isoseq preset to capture large deletion events (denoted by N in CIGAR)

AAV data analysis QC recommendations (cont.)

Summarize and visualize AAV alignment data using `summarize_AAV_alignment.py` and `plotAAVreport.R` command line tools available in GitHub

Workflow Step



Summarize & plot alignment results

- After alignment, `summarize_AAV_alignment.py` is used to generate 3 csv files that give you information at varying levels of detail
 - `<prefix>.summary.csv`: mapping summary
 - Read name/length, if the read is mapped, if the alignment is supplementary, mapping coords, and mapping identity
 - `<prefix>.per_read.csv`: AAV type assignments
 - Read name/length, if that read has a primary/supp alignment, the assigned AAV type (sc, ss, unknown), and subtype (full length or partial)
 - `<prefix>.nonmatch_stat.csv`: breakdown of mistakes found in each read
 - Read name, error position, type of error (I, D, X), length of error
- Plot the alignment results using `plotAAVreport.R`
 - `<prefix>._AAV_report.pdf` report
 - Distribution of mapped reference lengths, mapped identity to reference, fraction of reads that map to the reference start/end
 - Distribution of non-matches (deletions, insertions, substitutions) by reference position
 - Distribution of read lengths by assigned AAV type (ssAAV, scAAV, unknown)




AAV library example performance data

Example library QC and sequencing performance for AAV libraries


scAAV DNA extraction QC metrics

scAAV-CBA-eGFP Plasmid Map



scAAV target region: 662:2739 (2077 bp) (scAAV-CBA-eGFP, Vector Biolabs)

scAAV DNA Extraction QC




Example Bioanalyzer sizing QC analysis for an extracted scAAV DNA sample. The major peak at ~2.3 kb corresponds to the expected size range (~2.1kb) of this scAAV species (scAAV-CBA-eGFP, Vector Biolabs)

Example extracted DNA yield for an scAAV DNA sample extracted with the PureLink Viral RNA/DNA Mini kit from Invitrogen (12280050)

scAAV viral sample input	scAAV extracted DNA
1.5 x 10 ¹² genome copies [1.5 x 10 ¹³ GC/mL (100µl)]	1000 ng


Example library QC and sequencing performance for AAV libraries (cont.)

SMRTbell library QC and primary sequencing metrics





Extracted scAAV DNA input	1000 ng
Extracted scAAV DNA cleanup (%)	658 ng (65%)
Post-ligation & cleanup (%)	598 ng (59%)
Post-nuclease treatment & cleanup (%)	208 ng (20%)

Cleanup steps were performed using 1.3X SMRTbell cleanup beads.



Raw Base Yield	251.2 Gb
Mean Polymerase Read Length	46.0 kb
P0	30.6%
P1	68.1%
P2	1.3%

Example sequencing metrics for a scAAV sample run with Polymerase 2.1 / 90 pM on-plate concentration / 24-h movie time / 2-h Pre-Extension Time.




Mean HiFi Read Length	3,785 bp
Median HiFi Read Quality	Q32
HiFi Read Mean # of Passes	15

Total number of HiFi reads generated per SMRT Cell 8M is typically ≥2.5 Million for AAV samples achieve >50% P1 loading.

Example library QC and sequencing performance for AAV libraries (cont.)

Barcode demultiplexing metrics for a pooled (8-plex) scAAV sample


CCS Analysis Report



Demultiplexing Report


Barcode	Reads	Avg. Barcode Quality
bc1001	276,940	97.0
bc1002	264,996	97.0
bc1003	202,163	98.0
bc1008	188,205	97.0
bc1009	248,151	98.0
bc1010	257,343	98.0
bc1011	303,674	97.0
bc1012	250,317	98.0
None	854	0.0

Coverage Report for Demultiplexed scAAV Samples




Plot of demultiplexed AAV HiFi reads aligned to the correct position of the reference transfer plasmid. scAAV target region: 662:2739 (2077 bp) (scAAV-CBA-eGFP, Vector Biolabs)

Target Region Coverage per Demultiplexed Barcode




Example plotAAVreport.R report output

Summarize and visualize AAV alignment data using `summarize_AAV_alignment.py` and `plotAAVreport.R` command line tools available in GitHub




Distribution of mapped reference lengths and fraction of reads that map to the reference start/end for a scAAV control sample (scAAV-CBA-eGFP, Vector Biolabs).



Distribution of non-matches (deletions, insertions, substitutions) by reference position for a scAAV sample (scAAV-CBA-eGFP, Vector Biolabs).

Example plotAAVreport.R report output (cont.)

Summarize and visualize AAV alignment data using `summarize_AAV_alignment.py` and `plotAAVreport.R` command line tools available in GitHub



Distribution of mapped identities to reference and distribution of non-matches (deletions, insertions, substitutions) and their lengths by reference position for a scAAV sample (scAAV-CBA-eGFP, Vector Biolabs).

Err Type	Err Length	Count	Frequency (%)
D	1–10	33123	32.14
D	11–100	32	0.03
D	100–500	3	0.00
I	1–10	16532	16.04
I	11–100	66	0.06
I	100–500	1	0.00
I	>500	2	0.00
X	1–10	53289	51.71

Assigned Type	Assigned Subtype	Count	Frequency (%)
ssAAV	full	8	0.08
ssAAV	partial	8462	84.61
unknown	full	9	0.09
unknown	partial	1522	15.22

Frequency of non-matches (deletions, insertions, substitutions) by length [top] and assigned AAV type (ssAAV, scAAV or unknown) [bottom] for a ssAAV sample (pAV-CMV-GFP).



Technical documentation & applications support resources

Subtitle

Technical resources for AAV library preparation, sequencing & data analysis

AAV DNA extraction literature

- Tran, N.T. et al. (2020) AAV-genome population sequencing of vectors packaging CRISPR components reveals design-influenced heterogeneity. *Molecular Therapy Methods & Clinical Development*. 18:639 – 651. [[Link](#)]
- Guerin, K. et al. (2020) A novel next-generation sequencing and analysis platform to assess the identity of recombinant adeno-associated viral preparations from viral DNA extracts. *Human Gene Therapy*. 31:664 – 678 [[Link](#)]
- Lecompte et al. (2015) Advanced characterization of DNA molecules in rAAV vector preparations by single-stranded virus next-generation sequencing. *Molecular Therapy Nucleic Acids*. 4:E260. [[Link](#)]
- Gao, G. and Sena-Esteves, M. (2012). Introducing genes into mammalian cells: Viral vectors. *Molecular cloning: A laboratory manual*, Volume 2 New York: Cold Spring Harbor Laboratory Press. Pp. 1209 – 1313.

Sample preparation literature

- Procedure & checklist – Preparing multiplexed AAV SMRTbell libraries using SMRTbell prep kit 3.0 ([102-126-400](#))
- Quick reference card – Loading and pre-extension recommendations for the Sequel II and IIe systems ([101-769-100](#))
- Overview – Sequel systems application options and sequencing recommendations ([101-851-300](#))
- Technical note: Preparing DNA for PacBio HiFi sequencing – Extraction and quality control ([101-061920](#))
- Technical overview: AAV library preparation using SMRTbell prep kit 3.0 ([102-390-400](#))

Technical resources for AAV library preparation, sequencing & data analysis (cont.)

Data analysis resources

- Tutorial: Analyzing AAV data (GitHub): <https://github.com/Magdoll/AAV/wiki/Tutorial:-Analyzing-AAV-Data>
- SMRT Link v11.0 user guide ([102-278-200](#))
- SMRT Tools v11.0 reference guide ([102-278-500](#))
- Sequel II and IIe systems: Data files ([102-144-100](#))

Example PacBio data sets

AAV sequencing application	Dataset	Data type	PacBio system
scAAV (8-plex)	2021-scAAV/rh10-CBA-eGFP	HiFi Reads	Sequel II System
ssAAV (Single sample)	2021-pAV-CMV-GFP	HiFi Reads	Sequel II System

Posters

- Ranade, S. et al. (2018). High-throughput SMRT sequencing of clinically relevant targets. [[Link](#)]

Technical resources for AAV library preparation, sequencing & data analysis (cont.)

Publications

- Tran, N.T. et al. (2020) AAV-genome population sequencing of vectors packaging CRISPR components reveals design-influenced heterogeneity. *Molecular Therapy Methods & Clinical Development*. 18:639 – 651. [[Link](#)]
- Tai, P.W. et al. (2018) Adeno-associated virus genome population sequencing achieves full vector genome resolution and reveals human-vector chimeras. *Molecular Therapy Methods & Clinical Development*. 9:130–141. [[Link](#)]
- Paulk, N.K. et al. (2018) Bioengineered AAV Capsids with Combined High Human Liver Transduction In Vivo and Unique Humoral Seroreactivity. *Molecular Therapy*. 26:289 – 303. [[Link](#)]
- Xie, J. et al., 2017. Short DNA hairpins compromise recombinant adeno-associated virus genome homogeneity. *Molecular Therapy*. 25(6):1363–1374. [[Link](#)]
- Hüser, D. et al. (2014) Adeno-associated virus type 2 wild-type and vector-mediated genomic integration profiles of human diploid fibroblasts analyzed by third-generation PacBio DNA sequencing. *Journal of virology*, 88:11253–11263.[[Link](#)]

Webinars

- PacBio Journal Club (2021): AAV-genome population sequencing of vectors packaging CRISPR components reveals design-influenced heterogeneity [[Link](#)]
- Genewiz Cell and Gene Therapy Treatments Event (2021): Gene editing validation with single-molecule resolution using highly accurate HiFi reads [[Link](#)]
- PacBio Webinar (2019): Highly accurate SMRT sequencing for gene editing applications [[Link](#)]



www.pacb.com

Research use only. Not for use in diagnostic procedures. © 2022, PacBio. All rights reserved. Pacific Biosciences, the PacBio logo, PacBio, Circulomics, Omnione, SMRT, SMRTbell, Iso-Seq, Sequel, Nanobind, and SBB are trademarks of Pacific Biosciences of California, Inc. (PacBio). All other trademarks are the sole property of their respective owners.