





# COMPREHENSIVE HUMAN GENOMIC VARIANT DETECTION WITH HiFi LONG-READ SEQUENCING

PacBio® HiFi sequencing enables greater variant detection in more regions of the genome than alternative technologies, providing exceptional insight at any depth. Choosing the right level of sequencing is a critical first step for project planning, and HiFi sequencing delivers:


- Industry-leading small and structural variant calling performance
- Excellent small variant (SNV and indel) calling at 20x HiFi coverage
- Highly accurate structural variant detection with as little as 10x coverage

## The benefits of HiFi sequencing

 <b>Mappability</b> Access the entire genome	 <b>Phasing</b> Haplotype resolution	 <b>Exceptional accuracy</b> Call more variants with fewer errors	 <b>Methylation detection</b> Single-molecule 5mC insights
--	--	---	--

## The 20x HiFi human genome is best-in-class

Titration of a HG002 dataset sequenced on the Revio® system demonstrates industry-leading variant calling performance of 20x HiFi sequencing compared to 30x Illumina and 60x Oxford Nanopore Technologies (ONT) call sets. The power of a 20x HiFi genome is further supported by a recent independent study showing that 20x HiFi coverage is sufficient to call 96.2% of the clinically relevant germline variants identified at 30x.<sup>1</sup>



- 20x HiFi exceeds
- SNV accuracy of 99.8%
  - Indel accuracy of 60x ONT

Figure 1. Variant calling performance for single nucleotide variants (SNVs) and insertions and deletions (INDELs, < 50 bp) against the GIAB v4.2.1 benchmark<sup>2</sup>. HG002 PacBio HiFi coverage titration shown in magenta and purple (SNV-solid, INDEL-dashed) with each down-sampled depth marked as a point. Performance for the 30x Illumina DRAGEN 3.7.5<sup>3</sup> and 60x ONT R10.4.1<sup>4</sup> call sets is shown in orange and blue, respectively.

## Total benchmark errors by dataset




Figure 2. Total small variant (SNV and INDEL) calling errors for each technology. HiFi errors at 20x coverage are shown in magenta, with 30x Illumina DRAGEN 3.7.5 in orange, and 60x ONT R10.4.1 in blue. False negatives are represented as a lighter shade, while false positives are shown in the darker shade.

For structural variants, as little as 10x HiFi sequencing outperforms Illumina and ONT




Figure 3. Structural variant (SV) calling performance against the GIAB HG002 telomere-to-telomere (T2T) genome benchmark<sup>5</sup>. Magenta line represents HG002 HiFi performance (F1-score) with each down-sampled data point marked at 10x (93.2%), 12x, 14x, 16x, 18x, 20x (96.7%), 25x, 30x (97.2%), 35x, and 40x (97.4%). Performance for 60x ONT R10.4.1 data is shown in blue (F1-score, 93.1%) and 35x Illumina DRAGEN v4.2.4 data is shown in orange (F1-score, 55.5%).

The PacBio HG002 titrated HiFi dataset was generated with a library prepared using the automated HiFi prep kit 96 workflow<sup>6</sup> and sequenced on a single Revio SMRT<sup>®</sup> Cell using SPRQ<sup>™</sup> chemistry. Sequencing generated 146 Gb of mapped HiFi data and was down sampled to coverage levels ranging from 8x to 40x. HiFi data was processed using the PacBio WGS Variant Pipeline to generate small variant calls (SNVs and INDELS)<sup>7</sup>. Additional call sets, including 30x Illumina DRAGEN v3.7.5<sup>8</sup> and 60x ONT (Oxford Nanopore Technology) R10.4.1 SUP data<sup>4</sup>, were retrieved from publicly available sources. Small variant calls were benchmarked against the GIAB v4.2.1 reference set<sup>2</sup> using the hap.py framework<sup>9</sup>. Structural variant calls were generated from HiFi (Sawfish v0.12.4)<sup>9</sup>, 3.5x Illumina (DRAGEN v4.2.4)<sup>10</sup>, and 60x Nanopore R10.4.1 (Sniffles v2.0.7) data, and subsequently benchmarked with Truvari<sup>11</sup> against the T2T Q100 HG002 truth set<sup>12</sup>.

Details on the analysis including the Revio SPRQ chemistry aligned bam files can be found here:  
<https://github.com/PacificBiosciences/pb-benchmarks/tree/main/SPRQ-NOV-2024>

Visit [pacb.com/revio](https://pacb.com/revio) and [pacb.com/wgs](https://pacb.com/wgs) for more information on workflow and automation options, or connect with a PacBio scientist today.

### KEY REFERENCES

- Höps W., et al. (2024). HiFi long-read genomes for difficult-to-detect clinically relevant variants. *medRxiv*, 2024-09.
- Wagner, J., et al. (2022). Benchmarking challenging small variants with linked and long reads. *Cell Genomics*, 2(5).
- Rossi, M., et al. (2023). Unlocking the full potential of Illumina genomes: The journey to enhanced variant calling quality with DRAGEN informatics and high-quality sequencing. Illumina Genomics Research Hub, published June 29, 2023. ([https://human-pangenomics.s3.amazonaws.com/publications/PANGENOME\\_2022/DeepVariant/samples/HG002/HG002.30x\\_novaseq\\_pcrfree.dragen.vcf.gz](https://human-pangenomics.s3.amazonaws.com/publications/PANGENOME_2022/DeepVariant/samples/HG002/HG002.30x_novaseq_pcrfree.dragen.vcf.gz))
- Oxford Nanopore Open Data. <https://labs.epi2me.io/giab-2023.05/> Accessed September 2024.
- Telomere-to-telomere consortium HG002 Q100 project: <https://github.com/marbl/HG002>
- PacBio Guide & overview – Automated HiFi prep 96 and HiFi ABC for the Hamilton NGS STAR MOA system.
- PacBio HiFi Human WGS WDL pipeline: <https://github.com/PacificBiosciences/HiFi-human-WGS-WDL>
- Illumina Hap.py framework: <https://github.com/Illumina/hap.py>
- PacBio Sawfish: <https://github.com/PacificBiosciences/sawfish>
- Behera, S., et al. (2024). Comprehensive and accurate genome analysis at scale using DRAGEN accelerated algorithms. *bioRxiv*.
- English, A.C., et al. (2022). Truvari: refined structural variant comparison preserves allelic diversity. *Genome Biology*, 23(1), 271.
- HG002 T2T draft genome benchmark: [https://ftp-trace.ncbi.nlm.nih.gov/ReferenceSamples/giab/data/AshkenazimTrio/analysis/NIST\\_HG002\\_DraftBenchmark\\_defrabbV0.018-20240716/](https://ftp-trace.ncbi.nlm.nih.gov/ReferenceSamples/giab/data/AshkenazimTrio/analysis/NIST_HG002_DraftBenchmark_defrabbV0.018-20240716/)

Research use only. Not for use in diagnostic procedures. © 2024 Pacific Biosciences of California, Inc. ("PacBio"). All rights reserved. Information in this document is subject to change without notice. PacBio assumes no responsibility for any errors or omissions in this document. Certain notices, terms, conditions and/or use restrictions may pertain to your use of PacBio products and/or third-party products. Refer to the applicable PacBio terms and conditions of sale and to the applicable license terms at [pacb.com/license](https://pacb.com/license). Pacific Biosciences, the PacBio logo, PacBio, Circulomics, Omniome, SMRT, SMRTbell, Iso-Seq, Sequel, Nanobind, SBB, Revio, Onso, Apton, Kinnex, PureTarget, and SPRQ are trademarks of PacBio.