

## Application note

# Robust detection of somatic variants from tumor-normal samples with highly accurate long-read whole genome sequencing

## Introduction

Highly accurate long-read sequencing has enabled more complete germline variant detection and the completion of the human genome. In contrast to germline mutations, somatic variant detection has presented distinct challenges due to lower variant frequencies, yielding a lower signal-to-noise ratio and requiring higher sequencing depth. However, the exceptional accuracy and long read lengths of PacBio® HiFi sequencing on the Revio® and Vega™ systems is increasingly being applied for the robust detection of complex variants that were previously inaccessible with short reads or less accurate long reads (Vasan et al., 2019; Nattestad et al., 2018), now with the availability of the higher throughput needed to detect variants present at lower allele frequencies.

Short-read sequencing limits the ability to reconstruct important variation in cancer genomes, including complex structural variation and repetitive regions (Cortes-Ciriano et al., 2022). Inaccurate nanopore long-read sequencing faces challenges in the detection of small variants, such as single nucleotide variants (SNVs) (Olson et al., 2022). Paired tumor-normal WGS studies with HiFi **highly accurate long reads** can detect and phase a wide range of cancer-specific genetic variation, including SNVs, structural variants (SVs), deletions and insertions (indels), copy number variations (CNVs) and methylation, in a single assay. This Application note provides the workflow for the detection of somatic small variants, structural variants, and methylation for paired tumor-normal samples with HiFi whole genome sequencing (WGS) (Figure 1).

## Experimental results

HiFi WGS of the COLO829 melanoma cell line performed on the Revio system reveals genomic (Figure 1a) and epigenetic (Figure 1b) variation relative to the matched normal cell line.




Figure 1. IGV view of a) WGS and (b) methylation data for COLO829 tumor and matched normal COLO829-BL cell lines. Arrows indicate noted tumor variants.

Due to the high accuracy and precision of HiFi sequencing, paired-tumor WGS was able to reveal a 12 kb deletion in the *PTEN* tumor suppressor gene, along with an upstream SNV. Additionally, 5mC methylation calling reveals hypermethylation upstream of *KLLN*, a p53-regulated DNA replication inhibitor. Epigenetics plays a large role in the tumor progression through the regulation of oncogenes and tumor suppressors. Therefore, methylation calling is an important part of this workflow that can reveal valuable insights for cancer research (Hao et al., 2017).

In an analysis using the breast cancer cell line HCC1395, the accuracy of two tumor/normal sequencing depth ratios were compared for SNVs at varying allele frequencies (Figure 2). Both 30X tumor/30X normal and 60X tumor/30X normal coverage for HiFi WGS sequencing approach 0.99 F1 scores down to 0.2 variant allele frequency (VAF). At lower VAF, 60X/30X provided better variant detection, providing guidance for sequencing depth required as a function of VAF.




Figure 2. F1 scores for SNV detection in breast cancer cell line HCC1395 at 30X tumor/30X normal and 60X tumor/30X normal coverage.

Similarly, both sequencing depth ratios exhibit high recall (>83%) of a set of 62 structural variants from the [Valle-Inclan et al., 2022 truth dataset](#) using the Severus caller (Figure 3).




Figure 3. Recall of 62 structural variants improve with increasing COLO829 tumor purity.


## Tumor-normal WGS workflow overview

### Extract DNA from blood and/or tissue



It is recommended to use extraction protocols according to the [Nanobind® CBB kit](#) or [Nanobind PanDNA kit](#). The protocols have been validated for many tissue types including skeletal muscle, cerebellum, and kidney tumor and require 5–25 mg of sample (input material quantity may vary depending on tissue) and 200–400 µL of human whole blood.

### Library prep



It is recommended to construct single WGS libraries from extracted gDNA using the [SMRTbell® prep kit 3.0](#) and with 500 ng of DNA input per Revio SMRT® Cell with SPRQ™ chemistry (lower DNA input has been shown to yield high quality libraries with minimal decline in sequencing depth). High-quality genomic DNA will maximize sequencing coverage per SMRTbell library and reduce sequencing costs.

### HiFi sequencing



The resulting libraries can be prepared for sequencing on the [Revio](#) or [Vega](#) systems following instructions in the Sample Setup module of SMRT® Link. With 90% of bases at  $\geq Q30$ , both HiFi sequencing systems provide whole genome sequencing and on-board 5mC and 6mA methylation calling. One Revio SMRT Cell has an expected yield of 90 Gb, and can support 1,300 human genomes at 30x coverage per year.

### Analysis



The [HiFi somatic WDL](#) is a tumor-normal variant calling pipeline that consolidates variant callers for small and large variants into a single workflow. Together these tools identify and phase single nucleotide variants, small insertions and deletions, copy number variants, structural variants, and methylation.

## Conclusion

These results demonstrate the ability of the PacBio somatic tumor-normal WGS workflow to accurately detect biologically relevant variants, including methylation, across sample and sequencing conditions. This variant-calling performance is driven by the high accuracy of HiFi sequencing and is particularly important in cancers where allele frequencies are much lower. As such, this workflow represents the achievement of a more complete interrogation of cancer genomes, including the detection of both small and complex structural variants.

## Resources and references

### Resources

[Guide & overview](#) – Nanobind CBB kit

[Guide & overview](#) – Nanobind PanDNA kit

[Procedure & checklist](#) – Preparing whole genome and metagenome libraries using SMRTbell prep kit 3.0

[Dataset](#) – HCC1395 tumor and matched normal WGS on the Revio system

[Dataset](#) – COLO829 tumor and matched normal WGS on the Revio system

### References

- Cortés-Ciriano, I., Gulhan, D. C., Lee, J. J. K., Melloni, G. E., & Park, P. J. (2022). Computational analysis of cancer genome sequencing data. *Nature Reviews Genetics*, 23(5), 298-314.  
<https://doi.org/10.1038/s41576-021-00431-y>
- Hao, X., et al. (2017). DNA methylation markers for diagnosis and prognosis of common cancers. *Proceedings of the National Academy of Sciences*, 114(28), 7414-7419.  
<https://doi.org/10.1073/pnas.1703577114>
- Mandelker, D., & Ceyhan-Birsoy, O. (2020). Evolving significance of tumor-normal sequencing in cancer care. *Trends in Cancer*, 6(1), 31-39.  
<https://doi.org/10.1016/j.trecan.2019.11.006>
- Nattestad, M., et al. (2018). Complex rearrangements and oncogene amplifications revealed by long-read DNA and RNA sequencing of a breast cancer cell line. *Genome Research*, 28(8), 1126-1135.  
<https://doi.org/10.1101/gr.231100.117>
- Olson, N. D., et al., (2022). PrecisionFDA Truth Challenge V2: Calling variants from short and long reads in difficult-to-map regions. *Cell Genomics*, 2(5).  
<https://doi.org/10.1016/j.xgen.2022.100129>
- Valle-Inclan, J. E., et al. (2022). A multi-platform reference for somatic structural variation detection. *Cell Genomics*, 2(6).  
<https://doi.org/10.1016/j.xgen.2022.100139>
- Vasan, N., et al. (2019). Double PIK3CA mutations in cis increase oncogenicity and sensitivity to PI3Kα inhibitors. *Science*, 366(6466), 7.  
<https://doi.org/10.1126/science.aaw9032>

Research use only. Not for use in diagnostic procedures. © 2025 Pacific Biosciences of California, Inc. ("PacBio"). All rights reserved. Information in this document is subject to change without notice. PacBio assumes no responsibility for any errors or omissions in this document. Certain notices, terms, conditions and/or use restrictions may pertain to your use of PacBio products and/or third-party products. Refer to the applicable PacBio terms and conditions of sale and to the applicable license terms at [pacb.com/license](https://pacb.com/license). Pacific Biosciences, the PacBio logo, PacBio, Circulomics, Omniome, SMRT, SMRTbell, Iso-Seq, Sequel, Nanobind, SBB, Revio, Onso, Apton, Kinnex, PureTarget, SPRQ, and Vega are trademarks of PacBio.