



PACIFIC
BIOSCIENCES®



Back to the lab package

Contents

- Getting started!
- Useful guides and webinars
- Instrument start up and technical support
- New software and chemistry release
- Appendix



Getting started!

PACBIO SEQUENCING 101 SERIES



PacBio Sequencing - How it works
 PacBio sequencing instruments are powered by Single Molecule, Real-Time (SMRT) Sequencing technology. Learn how to go from DNA to discovery and explore the versatility of SMRT Sequencing across a broad range of applications, including whole genome sequencing, variant detection, RNA sequencing, to resolve complex populations, targeted sequencing, and epigenetics.
<https://www.youtube.com/watch?v= ID8JyAbwEo>

- PacBio's sequencing 101 series provides information on SMRT sequencing technology, genomics and much more!

Explore other posts in the sequencing 101 series

OVERVIEW – SEQUEL SYSTEMS APPLICATION OPTIONS AND SEQUENCING RECOMMENDATIONS

Document provides high-level application workflow guidance and links to protocols for preparing samples for sequencing on the Sequel Systems and analysis.



Application	Whole Genome Sequencing					RNA Sequencing		Metagenomics		Targeted Sequencing		
	De Novo Assembly - HiFi Reads	De Novo Assembly - for Low DNA Input	De Novo Assembly and Variant Detection - for Ultra-Low DNA Input	Microbial De Novo Assembly	Variant Detection	Structural Variation Detection	Iso-Seq Method	Single-Cell Iso-Seq	Full-length VES RNA Sequencing	Shotgun Metagenomic Profiling or Assembly	Amplicon Sequencing	No-Amp Targeted Sequencing
Experimental Design	<p>With 1 SMRT Cell 8M you can:</p> <ul style="list-style-type: none"> Produce reference quality assemblies for genomes up to 2 Gb Produce reference quality assemblies for genomes up to 1 Gb. Multiplex up to 2 small genomes on the Sequel II System. Produce reference quality assemblies for genomes up to 500 Mb. Sequence up to 48 microbes With 2 SMRT Cells 8M, Call CNVs, iCNVs, and SVs in a 3 Gb genome Call SVs for up to 2 samples with ~3 Gb genomes Characterize alternative splicing/transcriptome a genome with full length transcripts Characterize alternative splicing/transcriptome a genome with full length transcripts up to 2M full length reads with cell barcode and UMI information Multiplex up to 2 small genomes on the Sequel II System <p>Minimum Recommended Coverage: 10- to 15-fold HiFi read coverage per haplotype (for HiFi); >30-fold HiFi read coverage per diploid sample for de novo assembly; >15-fold HiFi read coverage for human variant detection; >30-fold Unique Molecular Coverage (UMC) per microbial genome; >15-fold HiFi read coverage of a human genome; 4- to 25-fold UMC per microbial genome; 10,000 unique reads/single cell for 300 cells; 1,000 unique reads/single cell for 300 cells; 5,000 unique reads/single cell for 300 cells</p> <p>Library Insert Size: 15-20 kb; ~15 kb; 10-12 kb; 10-15 kb; 15-20 kb; >15 kb; <2 kb to >3 kb; <2 kb to >3 kb</p>											
Procedure and Checklist Reference	<p>Prepare HiFi SMRTbell Libraries using SMRTbell Express Template Prep Kit 2.0; Prepare HiFi Libraries from Low DNA Input using SMRTbell Express Template Prep Kit 2.0; Prepare HiFi SMRTbell Libraries from Ultra-Low DNA Input using SMRTbell Express Template Prep Kit 2.0; Prepare Multiplexed, Microbial Libraries using SMRTbell Express Template Prep Kit 2.0; Prepare HiFi SMRTbell Libraries using SMRTbell Express Template Prep Kit 2.0 for Structural Variant Detection; Iso-Seq Express Template Preparation for Single and Sequel II Systems; Prepare Single-Cell Iso-Seq Libraries using SMRTbell Express Template Prep Kit 2.0; Amplify and Prepare Single-Cell SMRTbell Libraries</p>											
SMRTbell From Plate Preparation	<p>Minimum Input Amount: 15 µg; >400 ng per 1 Gb genome size (single sample); >300 ng per 600 Mb genome size (2-pile); 5-20 ng per 500 Mb genome size; 1 µg per microbe; 15 µg; 3 µg; 300 ng total RNA for 1st Strand cDNA Synthesis; >100 ng cDNA AFTER reamplification; 50</p> <p>Recommended PacBio Template Prep Kit: Express TPK 2.0 + SMRTbell Enzyme Channel Kit + Sequencing primer v2; Express TPK 2.0 + Sequencing primer v4; Express TPK 2.0 + Sequencing primer v4; Express TPK 2.0 + SMRTbell Enzyme Channel Kit + Sequencing primer v2; Express TPK 2.0 + SMRTbell Enzyme Channel Kit + Sequencing primer v2; Express TPK 2.0 + Sequencing primer v2; Express TPK 2.0 + Sequencing primer v4; Express TPK 2.0 + Sequencing primer v4</p> <p>Multiplexing/SMRT Cell: NA; NA; NA; Up to 48 microbes / SMRT Cell 8M; Up to 16 microbes / SMRT Cell 1M; NA; Up to 2 human samples/ SMRT Cell 8M; N/A SMRT Cell 1M; The protocol supports up to 12 barcodes available; Detects cell barcodes and UMIs; Up to 5M</p>											
Sequencing Preparation and Yield Performance	<p>Sequel II Binding Kit: 2.0; 2.0; 2.0; 2.0; 2.0; 2.0; 2.0 / 2.1* Sequel II Sequencing Plate: 2.0; 2.0; 2.0; 2.0; 2.0; 2.0; 2.0 Sequencing Mode: CCS; CCS; CCS; CLR; CCS; CLR; CCS Movie Collection Time: 30 h; 15 / 30 h; 30 h; 15 h; 30 h; 15 h; 24 h; 24 h Notes: *Use Binding Kit 2.1 as default and 2.0 for long transcripts</p>											
Sequencing Preparation and Yield Performance	<p>Sequel Binding kit: 3.0; 3.0; NA; 3.0; 3.0; 3.0; 3.0 Sequel Sequencing Plate: 3.0; 3.0; NA; 3.0; 3.0; 3.0; 3.0 Sequencing Mode: CCS; CCS; NA; CLR; CCS; CLR; CCS Movie Collection Time: 20 h; 20 h; NA; 10 h; 20 h; 10 h; 20 h; 20 h Notes: NA</p>											
Other Compatible (Contact PacBio for Details)	<p>SMRT Analysis GUI Applications: Mark PCR Duplicates, Trim cDNA Amplification Adapters, Structural Variant Calling; Demultiplex Barcodes followed by Microbial Assembly analysis; CCS with Mapping; Structural Variant Calling; Iso-Seq; Iso-Seq</p> <p>PacBio Connect Tools: Recommended: IPA, hifiasm, hicanu; Recommended: IPA, hifiasm, hicanu; Recommended: IPA, hifiasm, hicanu, Deep Variant; Recommended: Google Deep Variant; GitHub Tools</p>											

Application

De Novo Assembly - HiFi Reads **De Novo Assembly - for Low DNA Input**

Experimental Design

With 1 SMRT Cell 8M you can:

- Produce reference quality assemblies for genomes up to 2 Gb
- Produce reference quality assemblies for genomes up to 1 Gb. Multiplex up to 2 small genomes on the Sequel II System

Minimum Recommended Coverage

- 10- to 15-fold HiFi read coverage per haplotype
- 10- to 15-fold HiFi read coverage per haplotype

Library Insert Size

- 15 - 20 kb
- ~15 kb

Procedure and Checklist Reference

- [Preparing HiFi SMRTbell Libraries using SMRTbell Express Template Prep Kit 2.0](#)
- [Preparing HiFi Libraries from Low DNA Input using SMRTbell Express Template Prep Kit 2.0](#)


LIBRARY CONSTRUCTION

PRIMER ANNEALING AND POLY BINDING

Libraries can be stored @4C for 1 week (elution buffer, 101-633-500)

Libraries can be stored @-20C for ~6 months

Libraries can be stored @4C for 1 week



* For <3kb fragments, amplicons or transcripts use 2.1
For >3kb fragments, amplicons or transcripts use 2.0/2.2

PREPARATION FOR LOADING ON INSTRUMENT

	# Part Number	# Reactions	# Per Run	
SMRT Cell 8M Tray	101-389-001	4	1	Unused cells stored @RT
Sequel II Seq kit 2.0	101-820-200	4	1	Partially used kit can be stored @4C up to 100 hours
Sequel SMRT Cell Oil	100-621-300	5	1	Mitigates evaporation during runs on the SMRT Cell
Tube septa	100-667-700	10	1	Tube septa for closing reagent tubes on the system
Sample Plate Foil	100-667-400	100	1	Foil seals for covering Sample Plates
Sample Plate	000-448-888	50	1	96-well PCR plates for sample placement on the system
Sequel Mixing Plate	100-667-500	120	1	96-well deep well plates for reagent and sample mixing on the system
Sequel Pipet Tips v2	100-667-601	3840	10-12	Robotic pipet tips compatible for liquid handling on the system



Useful guides and webinars

DNA SAMPLE PREP AND RUN PERFORMANCE GUIDES AND APPLICATION BRIEFS

Sample Preparation Technical Notes

[Link](#)

Run Performance Evaluation Guide

Guide - Step-By-Step Run Performance Evaluation

Introduction

This guide provides information on how to troubleshoot sub-optimal performance of PacBio® SMRT® sequencing runs using the Internal Control and primary metrics immediately available upon run completion. It is intended to help customers understand and interpret the metrics most important to raw sequencing performance and downstream secondary analysis applications.

Step 1: Evaluate the Performance of the DNA Internal Control

What is the Sequel DNA Internal Control Complex and Why Use It?

The DNA Internal Control Complex is a previously prepared bound complex used as a spike-in sequencing control on Sequel® Systems. It is composed of a fixed insert of 1966 bps with ligated SMRTbell® adapters, annealed primer and a bound polymerase. This control complex is spiked into the bound sample at the end of SMRT® Link Sample Setup and is intended to be a known ideal sample for monitoring the sequencing performance of Sequel Systems. Poor or unexpected performance of this control could indicate potential issues with the instrument, sequencing reagents, or consumables.

The sequence of the control does not have identity to any known organism and therefore can be easily detected and separated from sample data (see Additional Information for sequence). Control read filtering occurs on the instrument and all detected control reads are sent to the scraps.bam file. A few low-quality control reads may pass through the filter to the subreads.bam file, but this leak rate should be rare (less than 1% of the total number of control reads from a given SMRT® Cell).

How to Prepare the DNA Internal Control Complex for Sequencing

The DNA Internal Control Complex is provided as a stock solution and requires a serial dilution prior to adding it directly to the sample. It is critical to appropriately follow the directions for preparing the control complex as indicated in SMRT Link Sample Setup or the application-specific Procedure & Checklist.

The Control Complex should be kept on ice when handling and, after use, promptly returned to storage at -20°C to maintain stability. The stock tube, and all subsequent dilutions, must be fully homogenized by gentle finger tapping and spinning down before use. Do not reuse preparations of the control from previous runs. Typical Control Complex read counts are shown in the "Sequel Systems Control Concordance: Too Low?" section (Table 1 for SMRT Cell 1M and Table 2 for SMRT Cell 8M).

How to use the DNA Internal Control Sequencing Performance Metrics for Troubleshooting

When a sequencing performance issue is observed for a sample, users should first evaluate the primary sequencing metrics of the DNA Internal Control. The Run QC table (Figure 1) in SMRT Link contains primary metrics for the control.

Control		Concordance	
Poly RT Mean (tsp)	Total Reads	Mean	Mode
55559	1420	0.85	0.87
49102	1730	0.84	0.87

Figure 1: DNA Internal Control sequencing metrics table reported in SMRT Link Run QC

Page 1 Part Number 101-993-600 Version 01 (September 2020)

[Link](#)

Best Practices Application Briefs

WHOLE GENOME SEQUENCING FOR DE NOVO ASSEMBLY BEST PRACTICES

Single Molecule, Real-Time (SMRT®) Sequencing on the Sequel® II System enables easy and affordable generation of high-quality de novo assemblies. With megabase size contig N50s, accuracies >99.99%, and phased haplotypes, you can do more biology – capturing undetected SNVs, fully intact genes, and regulatory elements embedded in complex regions.

CHOOSE THE SEQUENCING MODE THAT'S RIGHT FOR YOUR PROJECT

PacBio is the only sequencing technology to offer highly accurate long reads (HiFi reads).

HiFi Reads

- Highly accurate long reads with minimum accuracy of Q20 (99%)
- Small file sizes and fast analysis time
- Assemble up to a 2 Gb genome in a single SMRT Cell 8M for ~\$1,300*
- Run up to 200 samples (2 Gb) per year, per system!

Long Reads

- Longest reads, with half of data >50 kb, and maximum read lengths up to 175 kb
- Shorter sequencing run times
- Assemble up to a 3 Gb genome in a single SMRT Cell 8M for ~\$1,300*
- Run up to 400 samples (3 Gb) per year, per system!

*Read lengths, reads/data per SMRT Cell 8M, and other sequencing performance results vary based on sample quality/type and insert size. Prices, listed in USD, are approximate and may vary by region. Pricing includes library and sequencing reagents run on a Sequel II System and does not include instrument amortization or other reagents.

GENERATE COMPLETE AND ACCURATE GENOME ASSEMBLIES


HiFi reads and long reads both provide assemblies with high base pair accuracy and gene completeness. (A) Base pair accuracy as measured against a reference is >Q40, even outside of known gene regions. (B) BUSCO and larger, species-specific gene sets are >94% complete in both human and rice genome assemblies.

pacb.com/wgs


[Link](#)

FLEXIBLE MULTIPLEXING AND BARCODING SOLUTIONS

Barcoded M13 Primers




Barcoded Universal Primers



Barcoded Adapters

Adapter Ligation (SMRTbell Library Preparation)



Barcoded Primers



DNA SAMPLE PREPARATION ONLINE RESOURCE

Literature resource for sample collection and DNA extraction protocol references

The listing below is a collection of publications by the scientific community describing extraction protocols for high-molecular weight DNA followed by PacBio sequencing. When possible, the links point directly to the methods section (or supplementary information).

Plants

- [Animals](#)
- [Plants](#)
- [Fungi](#)
- [Protists](#)

If you have protocols you would like to share, or have questions about DNA extraction for PacBio sequencing, contact ExtractDNA@pacb.com.

Animals


- L. Invertebrates
 - L. Panova2016 – DNA extraction protocols for whole-genome sequencing in marine organisms
 - L. *microinvertebrates*
 - L. Laumer2020 – protocols for diverse meiofauna species, including *C. elegans* & *Castrella truncata* (SMRT-seq presentation)
 - L. *arthropods*
 - L. arachnids
 - L. Guerrero2019 – The Pacific Biosciences de novo assembled genome dataset from a parthenogenetic New Zealand wild population of the longhorned tick, *Haemaphysalis longicornis* Neumann, 1901
 - L. Liu2019 – DNA Methylation Patterns in the Social Spider, *Stegodyphus dumicola*
 - L. Burgess2018 – Draft genome assembly of the sheep scab mite, *Pteroptes ovis*
 - L. Randall2018 – The draft genome assembly of *Dermatophagoides pteronyssinus* supports identification of novel allergen isoforms in *Dermatophagoides* species
 - L. Cramaro2017 – Genome scaffolding and annotation for the pathogen vector *Isodetes ricinus* by ultra-long single molecule sequencing
 - L. crustaceans
 - L. Liu2019 – Crab muscle extraction & sequencing (Circularomics application note)
 - L. Baldwin-Brown2018 – A New Standard for Crustacean Genomes: The Highly Contiguous, Annotated Genome Assembly of the Clam Shrimp *Enlimnadia texana* Reveals HOX Gene Order and Identifies the Sex Chromosome
 - L. Gonçalves Athanasio2016 – Optimisation of DNA extraction from the crustacean *Daphnia*

Methods

Sample collection

A female yellowbelly pufferfish (Fig. 2), reared in the fish breeding centre of Fujian Normal University in Fuzhou City of Fujian Province was used for genome sequencing and assembly. Fresh white muscle, eye, skin, gonad, gut, liver, kidney, blood, gall bladder and air bladder tissues were collected and quickly frozen in liquid nitrogen for one hour. White muscle tissues were used for DNA sequencing for genome assembly, while all tissues were used for transcriptome sequencing.

Fig. 2



A picture of the yellowbelly pufferfish used in the genome sequencing and assembly.

DNA and RNA sequencing

Genomic DNA from white muscle tissue was extracted using the

at genomic DNA for long-read sequencing of single molecules

ee Genome With Single-Molecule Sequencing Uncovers a Recent Burst of Genes

sembly reveals the structure of the *Arabidopsis thaliana* Nd-1 genome and its gene set

sive, high-throughput plant genomic DNA extraction method suitable for

r consist

ng and

ect ru

olly pine mega-genome using long-read single-molecule sequencing

ice for the Fast-Growing Microalga *Picoclorum celeri*

Undaria pinnatifida. Chromosome-Level Assembly Using PacBio and Hi-C

ne suggests convergent functions of homeobox genes in algae and land plants (algae)

hile-genome sequencing in marine organisms (algae)

www.ExtractDNAforPacBio.com

Fungi

- L. Faure2019 – Long Read Genome Sequence of the Sugar Beet Rhizosphere Mycoparasite *Pythium oligandrum*
- L. Nagappan2018 – Improved nucleic acid extraction protocols for *Candida bombinis*, *G. miniaticinctum* and *G. torulosa*
- L. Schwessinger2017 – Extraction of high molecular weight DNA from fungal rust spores for long read sequencing
- L. Solomon2016 – Robust and effective methodologies for cryopreservation and DNA extraction from anaerobic gut fungi
- L. Sonnenberg2016 – A detailed analysis of the recombination landscape of the button mushroom *Agaricus bisporus* var. *bisporus*

PacBio does not assume responsibilities/guarantees for these external publications/protocols, but we are happy to help as best as we can to guide / connect. Please contact ExtractDNA@pacb.com for more discussions around your particular species & sequencing project!

DNA AND SAMPLE PREPARATION FOR PACBIO HIFI SEQUENCING – BEST PRACTICES AND FAQ



In this *webinar* Greg Young covers the points below:

- *Preparing DNA for PacBio HiFi sequencing and selecting the right HiFi library preparation protocol for your project*
- *Interpreting sequencing metrics and troubleshooting*
- *Sharing Frequently Asked Questions*

PACBIO DATA DEEPLDIVE: A CLOSER LOOK AT HIFI SEQUENCING




In this *webinar* Dan Browne covers the points below:

- Exploring the content of PacBio HiFi BAM files (base calls, quality scores and kinetic information)
- Providing examples of how HiFi reads can be filtered
- Going through recommendations for downstream analysis - for example, how HiFi reads are demultiplexed, assembled and mapped to reference sequences

SMRTLINK USER GUIDE

SMRT Link User Guide



[Link](#)

Sequel Ii System: Location of HiFi reads files

Sequel[®] Ii System: Location of HiFi Reads Files

Introduction

This document describes how to locate the `hifi_reads` files generated by SMRT[®] Link when you perform an on-instrument CCS analysis on the Sequel[®] Ii System.

Note: This document applies **only** to the Sequel Ii System.

HiFi Reads Generation

An on-instrument CCS analysis generates a `reads.bam` file and transfers it to the network server. The `reads.bam` file contains HiFi Reads and non-HiFi Reads, and should **not** be used unfiltered as input for tools that expect \geq QV 20. SMRT Link **automatically** launches an Export Reads analysis on the `reads.bam` to filter out the HiFi Reads, and generates the following HiFi data files by default:

- `<Movie_Name>.hifi_reads.fastq.gz` - FASTQ file containing HiFi Reads
- `<Movie_Name>.hifi_reads.fasta.gz` - FASTA file containing HiFi Reads
- `<Movie_Name>.hifi_reads.bam` - BAM file containing HiFi Reads

If **not** using SMRT Link for subsequent analysis, please use these three files as input with any third-party analysis tools.

Finding the `hifi_reads` Files Generated Using On-Instrument CCS

1. In Run QC, click the desired run, then click the sample name to view the CCS Data Set.

No.	Name	Movie Title (L)	Status	Total Da.
AD1	2M_Lambda_OCC	2	Comp...	16.63
BD1	2M_Lambda_OCC	2	Comp...	16.43
CD1	2M_Lambda_OCC	2	Comp...	20.71

2. Click **Analyses** in the left-side panel.

Page 1

[Link](#)




Instrument start up and technical support

INSTRUMENT SELF-TEST

A self-test is recommended for dormant instruments. The robot, optics, chiller and pipettor are tested. In addition, the instrument will perform a work deck clean up and start a 2 minute acquisition with the test cell (found on the workdeck) testing data transfer to data storage location. The self-test takes 45 minutes to complete.


1



The screenshot shows the instrument's main interface with a menu open. The menu items are: Home, Admin, Tools, Restart..., Shutdown..., and Self-Test... (highlighted). The background shows a table with columns for Cell, Well, Sample, and Bases (Ob). The 'Cell 1 Stats' section displays two graphs: 'Polymerase Read Length' and 'Longest Subread Length', both showing a sharp peak at the beginning of the x-axis (Length in bp). Below the graphs is a 'Productivity' table:

	P0	P1	P2
Productivity	46.2%	52.5%	1.3%
	(2,706,071)	(4,205,244)	(105,358)

2




Instrument Self-Test

This self-test takes approximately 45 minutes

CANCEL START

3



Instrument Self-Test

Running...

CANCEL START

Instrument Self-Test

Completed: **All Tests Passed**


- PASSED: The robot control board reports no errors
- PASSED: The instrument can report all stage control registers
- PASSED: Robot control board 24V supply can be verified
- PASSED: Robot control board 36V supply can be verified
- PASSED: Robot can home axes
- PASSED: Robot can move axis using a servo loop
- PASSED: Robot can move XY to limits
- PASSED: Robot can move with variable speeds without losing precision
- PASSED: Robot board test date and time are recorded
- PASSED: Instrument can communicate with the pipettor
- PASSED: Pipettor can move the plunger
- PASSED: Pipettor self-tests pass
- PASSED: Pipettor can pick up a tip
- PASSED: Verify pipettor leakage with no alarms
- PASSED: Illumination control board communication
- PASSED: Illumination control board reports no errors
- PASSED: Instrument can report all illumination control registers
- PASSED: Illumination control board 12V and 5V supplies are nominal
- PASSED: Illumination control board revision matches installed FRU assembly
- PASSED: Lightbrush AT motor moves without losing encoder counts
- PASSED: Lightbrush CT motor moves without losing encoder counts
- PASSED: Lightbrush Pitch motor moves without losing encoder counts
- PASSED: Lightbrush Yaw motor moves without losing encoder counts
- PASSED: Lightbrush Pattern Rotation has full range of motion
- PASSED: Illumination board test date and time are recorded
- PASSED: Instrument can perform a dry run acquisition and data transfer
- PASSED: Post-run illumination checks performed without alarms
- PASSED: Instrument can move the test cell from the stage to its home location

DONE

HOW TO CONTACT PACBIO FOR TECHNICAL SUPPORT

We have an experienced team of technical support scientists, engineers, and bioinformaticians to assist you with all of your support needs!

1. Log into the **PacBio Customer Portal** (<http://www.pacbioportal.com/>) from PacBio's corporate website and create a new case



- **Preferred route** for most customers because PacBio Tech Support and your local field team (FAS and FSE) are immediately notified via email of your issue
- Facilitates the exchange of troubleshooting and diagnostic data and maintains a record of all communications between our support team and yourself





2. Email: support@pacb.com


3. Tech Support Toll-free Telephone Number (Mon. – Fri., 7:00 a.m. to 5:00 p.m. PST)
+1.877.920.PACB [+1 (877) 920-7222]

HOW TO LOG INTO THE CUSTOMER PORTAL

1. Click on the '[Customer Login](#)' link at the bottom of any PacBio webpage
 - Or go to <https://mypacbio.force.com/customer/s/login/> directly
2. Enter your Customer Portal login credentials.
 - If you do not remember your password, click on the "[Forgot your password?](#)" link on the Customer Portal login webpage.
 - For future Users needing accounts, contact your Local FAS to request access or simply click on the "[Request Access](#)" link on the [PacBio Technical Support](#) webpage and the request will be routed to your local PacBio Support Team for activation within 24 hours.



YOUR PACBIO LOCAL ACCOUNT TEAM



Jonathan Seddon
Territory Account
Manager
UK & IRE



Carla Simoes
Sequencing Application
Sales Specialist
UK & IRE



Kirrenjit Kaur
Field Application Specialist
UK & IRE



Johnathan Smith,
Field Service Engineer
Northern Europe



Bas Tolhuis
BFX Field Application Specialist
EMEA



New software and chemistry release

NEW SMRT LINK V10.1 SOFTWARE RELEASE: KEY FEATURES & IMPROVEMENTS

- Supports complete **SMRT Link analysis workflow on Amazon Cloud** (SMRT Link Cloud)
- Features several **new and improved analysis applications** for HiFi de novo assembly, SV calling, multiplexed Iso-Seq analysis and SARS-CoV-2 full-viral genome sequencing
- Reduces HPC requirements to enable **lower-cost data analysis & storage configurations**
- Provides a **simplified user experience** for run setup and includes **usability improvements** to support high-throughput sequencing



NEW SEQUEL II BINDING KIT 2.2

Part Number	New Products: Name and Description
102-089-000	The Sequel II Binding Kit 2.2 and Internal Control 1.0 contains reagent supplies to bind prepared DNA template libraries to the Sequel II Polymerase 2.2 in preparation for HiFi sequencing. The result is a DNA polymerase template complex ready for use with Adaptive Loading on the Sequel II Systems. Sequel II Polymerase 2.2 should be used only with Sequel II Sequencing Kit 2.0. Reagent quantities support approximately 24 binding reactions, and up to 70 SMRT Cell 8M runs depending on use case, sample size and concentration. Note Sequel II polymerase 2.2 is not recommended for inserts <3 kb on the Sequel II and Sequel IIe Systems. This item also includes the Internal Control 1.0, used as a positive control for loading and sequencing.
101-932-600	The SMRTbell Enzyme Clean Up Kit 2.0 contains 5 enzymes in a single tube and an enzyme reaction buffer. The kit is designed to reduce the template molecules that are not intact to enhance HiFi sequencing. The kit contains sufficient reagents to perform 18 reactions.
102-067-400	Sequencing Primer v5 is a sequencing primer used to anneal to HiFi SMRTbell templates generated with the SMRTbell Express Template Prep Kit 2.0.
102-104-700	Sequel II HiFi Bundle - 18 2.0 contains a bundle of PacBio consumables required to perform the de novo assembly and/or variant detection applications using HiFi reads on the Sequel II and Sequel IIe Systems including template preparation and polymerase binding. The bundle specifically includes the following kits and their associated quantities: SMRTbell Express Template Prep Kit 2.0 (PN: 100-938-900 / Qty: 1), SMRTbell Enzyme Cleanup Kit 2.0 (PN: 101-932-600 / Qty: 1), Sequencing Primer V5 (PN:102-067-400 / Qty: 2), AMPure PB Beads (PN: 100-265-900 / Qty: 1), Sequel II Binding Kit 2.2 and Internal Control 1.0 (PN:102-089-000 / Qty: 1).
102-088-900	The HiFi Express Template Prep Kit 2.0 Bundle contains one each: SMRTbell Express Template Preparation Kit 2.0, SMRTbell Enzyme Clean Up Kit 2.0 and Sequencing Primer v5. Each kit contains sufficient reagents to perform 18 standard reactions.

New binding kit to support HiFi Whole Genome Sequencing applications (*de novo* assembly and variant detection). New chemistry further increases the accuracy of our HiFi reads.

Binding kit 2.2 features:

- Faster polymerase
- Improved version of Enzyme Cleanup kit
- New primer v5



Appendix

Application	Template prep	Sequencing primer	Binding kit
De Novo Assembly HiFi Reads	<ul style="list-style-type: none"> SMRTbell Express Template 2.0 (100-938-900) 18 standard libraries SMRTbell Enzyme Clean Up 2.0 (101-932-600) NEW 18 reactions 	<ul style="list-style-type: none"> V5 (102-067-400) NEW 10 reactions 	<ul style="list-style-type: none"> 2.2 (102-089-000) NEW 24 reactions
De Novo Assembly Low Input	<ul style="list-style-type: none"> SMRTbell Express Template 2.0 (100-938-900) 18 standard libraries SMRTbell Enzyme Clean Up (101-746-400) 18 reactions Barcoded Overhang Adapter 8A (101-628-400) 6 reactions each barcode Barcoded Overhang Adapter 8B (101-628-500) 6 reactions each barcode 	<ul style="list-style-type: none"> V4 (included in prep kit) 	<ul style="list-style-type: none"> 2.0 (101-842-900) 24 reactions
De Novo & Variant Detection Ultra Low Input	<ul style="list-style-type: none"> SMRTbell gDNA Sample Amplification (101-980-000) 18 standard libraries SMRTbell Express Template 2.0 (100-938-900) 18 standard libraries 	<ul style="list-style-type: none"> V4 (included in prep kit) 	<ul style="list-style-type: none"> 2.0 (101-842-900) 24 reactions
Microbial De Novo Assembly	<ul style="list-style-type: none"> SMRTbell Express Template 2.0 (100-938-900) 48 microbial genomes (12kb) Barcoded Overhang Adapter 8A (101-628-400) 6 reactions each barcode Barcoded Overhang Adapter 8B (101-628-500) 6 reactions each barcode 	<ul style="list-style-type: none"> V4 (included in prep kit) 	<ul style="list-style-type: none"> 2.0 (101-842-900) 24 reactions
Variant Detection	<ul style="list-style-type: none"> SMRTbell Express Template 2.0 (100-938-900) 18 standard libraries SMRTbell Enzyme Clean Up 2.0 (101-932-600) NEW 18 reactions 	<ul style="list-style-type: none"> V5 (102-067-400) NEW 10 reactions 	<ul style="list-style-type: none"> 2.2 (102-089-000) NEW 24 reactions

Application	Template prep	Sequencing primer	Binding kit
Structural Variant Detection	<ul style="list-style-type: none"> SMRTbell Express Template 2.0 (100-938-900) 18 standard libraries Barcoded Overhang Adapter 8A (101-628-400) 6 reactions each barcode Barcoded Overhang Adapter 8B (101-628-500) 6 reactions each barcode 	<ul style="list-style-type: none"> V2 (101-847-900) 10 reactions 	<ul style="list-style-type: none"> 2.0 (101-842-900) 24 reactions
Iso-Seq	<ul style="list-style-type: none"> Iso-Seq Express Oligo kit (101-737-500) 24 reactions SMRTbell Express Template 2.0 (100-938-900) 18 standard libraries 	<ul style="list-style-type: none"> V4 (included in prep kit) 	<ul style="list-style-type: none"> 2.0 (101-842-900) 24 reactions 2.1 (101-843-000)* 24 reactions
Full length 16S rRNA Seq	<ul style="list-style-type: none"> SMRTbell Express Template 2.0 (100-938-900) 96 samples 	<ul style="list-style-type: none"> V4 (included in prep kit) 	<ul style="list-style-type: none"> 2.1 (101-843-000) 24 reactions
Shotgun Metagenomics	<ul style="list-style-type: none"> SMRTbell Express Template 2.0 (100-938-900) 18 standard libraries SMRTbell Enzyme Clean Up (101-746-400) 18 reactions Barcoded Overhang Adapter 8A (101-628-400) 6 reactions each barcode Barcoded Overhang Adapter 8B (101-628-500) 6 reactions each barcode 	<ul style="list-style-type: none"> V2 (101-847-900) 10 reactions 	<ul style="list-style-type: none"> 2.0 (101-842-900) 24 reactions

* Use Binding kit 2.1 as default and Binding kit 2.0 for transcripts >3kb

	Template prep	Sequencing primer	Binding kit
Amplicon Sequencing	<ul style="list-style-type: none"> • SMRTbell Express Template 2.0 (100-938-900) 96 amplicon samples • Barcoded Overhang Adapter 8A (101-628-400) 6 reactions each barcode • Barcoded Overhang Adapter 8B (101-628-500) 6 reactions each barcode • Bar Univ F/R Primers 96 Plate (101-629-100) 96 barcodes 	<ul style="list-style-type: none"> • V4 (included in prep kit) 	<ul style="list-style-type: none"> • 2.0 (101-842-900) 24 reactions • 2.1 (101-843-000)** 24 reactions
No-Amp Targeted Sequencing	<ul style="list-style-type: none"> • No-Amp Accessory kit (101-788-900) 20 reactions • SMRTbell Enzyme Clean Up (101-746-400) 18 reactions 	<ul style="list-style-type: none"> • V4 (included in accessory kit) 	<ul style="list-style-type: none"> • 2.0 (101-842-900) 24 reactions

** Use Binding kit for 500bp-3kb and Binding kit 2.0 for >3kb

Other reagents	
Ampure PB, 5mL (100-265-900)	For library construction and cleanup of polymerase-bound complexes
Elution Buffer, 50mL (101-633-500)	For library prep and barcoding workflows requiring additional clean-up steps
MagBead Binding Buffer (101-046-400)	Used to purify DNA/polymerase complexes after binding of polymerase – functionally equivalent to Complex Dilution Buffer for the Internal Control and final loading dilution steps

Reagents for instrument loading	# Reactions	# Per run	
SMRT Cell 8M Tray (101-389-001)	4	1	Unused cells stored @RT
Sequel II Seq kit 2.0 (101-820-200)	4	1	Partially used kit can be stored @4C up to 100 hours
Sequel SMRT Cell Oil (100-621-300)	5	1	Mitigates evaporation during runs on the SMRT Cell
Tube septa (100-667-700)	10	1	Tube septa for closing reagent tubes on the system
Sample Plate Foil (100-667-400)	100	1	Foil seals for covering Sample Plates
Sample Plate (000-448-888)	50	1	96-well PCR plates for sample placement on the system
Sequel Mixing Plate (100-667-500)	120	1	96-well deep well plates for reagent and sample mixing on the system
Sequel Pipet Tips v2 (100-667-601)	3840	10-12	Robotic pipet tips compatible for liquid handling on the system



PACBIO®

www.pacb.com

For Research Use Only. Not for use in diagnostic procedures. © Copyright 2021 by Pacific Biosciences of California, Inc. All rights reserved. Pacific Biosciences, the Pacific Biosciences logo, PacBio, SMRT, SMRTbell, Iso-Seq, and Sequel are trademarks of Pacific Biosciences. Pacific Biosciences does not sell a kit for carrying out the overall No-Amp Targeted Sequencing method. Use of these No-Amp methods may require rights to third-party owned intellectual property. FEMTO Pulse and Fragment Analyzer are trademarks of Agilent Technologies Inc.

All other trademarks are the sole property of their respective owners.