



# Streamlined SMRTbell® library generation using addition-only, single tube strategy for all library types reduces time to results.

Jenny Ekholm<sup>1)</sup>, Jenny Gu<sup>1)</sup>, Jeff Burke<sup>2)</sup>, Michelle Kim<sup>2)</sup>, Justin Sibert<sup>2)</sup>, Kelvin Liu<sup>2)</sup>  
Pacific Biosciences, 1305 O'Brien Drive, Menlo Park, CA 94025, Circulomics Inc. 701 E Pratt St, Baltimore, MD 21202

## ABSTRACT

We have streamlined the SMRTbell library generation protocols with improved workflows to deliver seamless end-to-end solutions from sample to analysis. A key improvement to the development of a single-tube reaction strategy that shortened hands-on time needed to generate each SMRTbell library, reduced time-consuming AMPure purification steps, and minimized sample-handling induced gDNA damage to improve the integrity of long-insert SMRTbell templates for sequencing. The improved protocols support all large-insert genomic libraries, multiplexed microbial genomes, and amplicon sequencing. These advances enable completion of library preparation in less than a day (approximately 4 hours) and opens opportunities for automated library preparation for large-scale projects. Here we share data summarizing performance of the new SMRTbell Express Template Kit 2.0 representing our solutions for 10 kb and >50 kb large-insert genomic libraries, complete microbial genome assemblies, and high-throughput amplicon sequencing. The improved throughput of the Sequel System with read lengths up to 30 kb and high consensus accuracy (>99.999% accuracy) makes sequencing with high-quality results increasingly accessible to the community.

## EXPRESS 2.0 WORKFLOWS




Figure 1. Detailed Express 2.0 workflows for long-insert and multiplexed microbial libraries.

## RECOMMENDED gDNA INPUT REQUIREMENTS

Insert Size Target	gDNA Selection	gDNA Quality Required	Shearing Method	Range of Sheared and Concentrated Input gDNA (µg)	Recommended Sheared and Concentrated DNA Amount	Insert Size Target	gDNA Selection	gDNA Quality Required	Shearing Method	Recommended microbial gDNA Amount
>10 kb	No	Made >20 kb	g-TUBE	1.0 to 2.0 µg	1.0 µg	>12 kb	Optional	Preferably >20 kb	g-TUBE	1.0 µg
10-15 kb	Yes	Made >40 kb	g-TUBE or Neogenex®	2.0 to 5.0 µg	3.0 µg					
100 kb	Yes	Made >20 kb	25G needle or Neogenex®	3.0 to 5.0 µg	3.0 µg					

## EXPECTED SEQUENCING PERFORMANCE

Library Type with 3.8 sequencing chemistry	Polymersase Read Length (Average)	Polymersase Read Length (NSD)	Longest observed average (NSD)	Library Type with 3.8 sequencing chemistry	Polymersase Read Length (Average)	Polymersase Read Length (NSD)	Longest observed average (NSD)
>50 kb libraries with no size selection	up to 30 kb	up to 40 kb	up to 35 kb	>10 kb insert	up to 45 kb	up to 85 kb	up to 8 kb
10 kb libraries with no size selection	up to 24 kb	up to 30 kb	up to 25 kb	>12 kb insert	up to 45 kb	up to 85 kb	up to 11 kb

Sequencing performance, made using PacBio Cell and other reagents tested only based on sample quality and insert size.

## LARGE-INSERT gDNA QC AND SEQUENCING RESULTS




Table 3. Large-insert gDNA library sequencing yield for the different library target sizes.

## MULTIPLEXED MICROBIAL gDNA QC AND SEQUENCING RESULTS




Table 4. Detailed information on the 10 microbial genomes tested.



Figure 5. Sequencing results for pooled microbial genome libraries. We demonstrated an average concentration of 8.0 to 20.000 additional to your sequencing with optimal results.

## Multiplexed Microbial Genome Assembly Results

Barcode ID	Sample ID	Expected Genome Size (kb)	Assembled Genomes (Gen)	GC Content (%)	# Contigs (non-redundant)	Proportion of Non-Reference DNA	Finished Contigs (k)	Max. Contig Length (bp)	NS Contig Length (bp)	Sum of Contig Length (bp)	Pre-Assembled Yield (%)	Filtered Subread Coverage (x)	Assembly Notes
BC1001	E. coli K12 (Jugens)	4,833,240	1	50.8%	1	52	1	4,842,498	4,842,498	4,842,498	95.1%	84	Complete chromosomal assembly
BC1002	E. coli K12 (Circulomics)	4,833,240	1	50.8%	1	57	1	4,842,500	4,842,500	4,842,500	92.9%	55	Complete chromosomal assembly
BC1005	K. pneumoniae	5,781,501	1	57.0%	1	48	5	5,435,746	5,432,746	5,748,500	92.5%	50	Complete chromosomal assembly, and 140 kb and 85 kb plasmid assemblies. 118 kb plasmid captured in 2 contigs. Misassembled 2 kb plasmid.
BC1010	B. cereus	5,427,083	1	51.3%	1	35	2	5,408,315	5,408,315	5,423,588	92.7%	59	Complete chromosomal and 16 kb plasmid assembly.
BC1012	L. monocytogenes	3,032,208	1	37.9%	1	46	2	3,043,149	3,043,149	3,137,529	93.7%	68	Complete chromosomal assembly
BC1015	S. aureus	5,082,553	1	51.0%	1	42	1	4,913,454	4,913,454	4,913,454	92.1%	53	Complete chromosomal assembly. Missing eight expected plasmids.
BC1016	N. meningitidis	2,194,951	1	51.0%	1	50	1	2,213,947	2,213,947	2,213,947	92.2%	74	Complete chromosomal and 27 kb plasmid assemblies.
BC1018	S. aureus	2,806,345	1	52.3%	1	16	2	2,778,800	2,778,800	2,806,350	92.6%	62	Complete chromosomal and 10 kb plasmid assemblies.
BC1019	E. coli Strain W	5,005,347	1	50.8%	1	44	2	4,988,327	4,988,327	5,004,300	93.0%	65	Complete chromosomal and 103 kb plasmid assemblies. Misassembled 5 kb plasmid.
BC1022	S. aureus HRP107	2,981,405	1	51.3%	1	33	2	2,982,756	2,982,756	2,984,872	93.4%	62	Complete chromosomal and 24 kb plasmid assemblies.

Table 5. Multiplexed Microbial Genome Assembly Results. Summary of genome assembly results with respect to GC content and genomic complexity.

## CONCLUSION

The SMRTbell Express Template Prep Kit 2.0 product family provides streamlined, expedited workflows for improved ease of use with better recovery and yield. This enables the completion of library generation within a day for sequencing.

For large-insert gDNA library, SMRTbell library construction can be completed in 4 hours.

For multiplexed microbial genome assemblies, the key improvement is the change in shearing condition to generate ~12 kb insert libraries, and the much simplified workflow to complete library preparation from gDNA to SMRTbell library in ~ 8.5 hrs, depending on number of samples being processed. Customers will appreciate the reduced AMPure purification steps required by the protocol, and the single-tube addition strategy opens up opportunities to explore automated library generation.

Forthcoming this year in Summer 2019, we anticipate new solutions for our full-length transcript Iso-Seq application and amplicon sequencing to join the Express TPK 2.0 product family.

## REFERENCES

- Procedure A Checklist - Preparing gDNA Libraries Using the SMRTbell Express Template Preparation Kit 2.0
- Procedure A Checklist - Preparing Multiplexed Microbial Libraries Using SMRTbell Express Template Prep Kit 2.0
- Analysis Procedure - Multiplexed Microbial Assembly with SMRT Link v8.0.0 and Express Template Prep Kit 2.0
- Express Microbial Multiplexing Calculator