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KEY REFERENCES

1. Olson, N.D. et al. (2021) **precisionFDA Truth Challenge V2: Calling variants from short- and long-reads in difficult-to-map regions.** *bioRxiv preprint*.
2. Wenger, A. M. et al. (2019) **Accurate circular consensus long-read sequencing improves variant detection and assembly of a human genome.** *Nature Biotechnology*. 37, 1155–1162.
3. Farrow, E. et al. (2021) **Applications of Third Generation Sequencing in Unsolved Disease.** ACMG conference presentation.
4. **Procedure + checklist – Preparing HiFi SMRTbell libraries using SMRTbell express template prep kit 2.0.** PacBio documentation.
5. **Overview – Sequel systems application options and sequencing recommendations.** PacBio documentation.
6. **Technical note: Preparing DNA for PacBio HiFi sequencing – Extraction and quality control.** PacBio literature.
7. **pbsv – PacBio structural variant (SV) calling and analysis tools.** PacBio GitHub.
8. **Product brochure: SMRT Link – Explore and analyze your data with confidence.** PacBio literature.
9. Poplin, R. et al. (2018) **A universal SNP and small-indel variant caller using deep neural networks.** *Nature Biotechnology*. 36, 983–987.
10. Martin, M. et al. (2016) **WhatsHap: fast and accurate read-based phasing.** *bioRxiv preprint*.

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Information in this document is subject to change without notice. PacBio assumes no responsibility for any errors or omissions in this document. Certain notices, terms, conditions and/or use restrictions may pertain to your use of PacBio products and/or third party products. Refer to the applicable PacBio terms and conditions of sale and to the applicable license terms at <http://www.pacb.com/legal-and-trademarks/terms-and-conditions-of-sale/>. PacBio, the PacBio logo, PacBio, SMRT, SMRTbell, Iso-Seq, and Sequel are trademarks of PacBio. PacBio does not sell a kit for carrying out the overall No-Amp targeted sequencing method. Use of the No-Amp method may require rights to third-party owned intellectual property. SageELF is a trademark of Sage Science. Femto Pulse and Fragment Analyzer are trademarks of Agilent Technologies Inc. All other trademarks are the sole property of their respective owners.

© 2022 Pacific Biosciences of California, Inc. All rights reserved. For research use only.
102-193-604 V1 02DEC2021

PacBio