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### Long-read fusion gene detection tools

Several fusion gene detection tools have been developed specifically for long-read RNA data, including [CTAT-LR-fusion](#), [pbfusion](#), [PB\\_FLIP](#), and [JAFFAL](#).

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Study	Samples	Method	Key findings
Nattestad et al. (2018)	SKBR3 cell line	WGS, Iso-Seq	SVs, 2- and 3-hop fusion genes
Dondi et al. (2023)	Ovarian cancer	Single-cell Iso-Seq with Kinnex*	Cell-type specific tumor isoforms, patient-specific fusion gene
Qin et al. (2024)	SeraCare FusionStandard	Single-cell Iso-Seq with Kinnex*	CTAT-LR-fusion long-read fusion detection tool
Black et al. (2024)	Chronic lymphocytic leukemia (CLL) samples with BTK resistance	Single-cell Iso-Seq with Kinnex*	Improved single cell genotyping and clonal evolution tracing
Wedemeyer et al. (2024)	Patient-derived PIK3CA-altered fibroblast	Targeted single-cell Iso-Seq with Kinnex*	Improved single cell genotyping
Veiga et al. (2022)	Breast cancer cell lines and tumor samples	Iso-Seq	Tumor-associated novel isoforms
Kohli et al. (2017)	Prostate cancer	Targeted Iso-Seq	Co-expression of AR-V9 with AR-V7
Huang et al. (2022)	Gastric cancer	Iso-Seq	Novel isoforms with alternative promoters associated with survival outcome
Li et al. (2024)	Colorectal cancer	Single-cell Iso-Seq	Tumor-specific isoforms, neoepitopes
Wijeratne et al. (2024)	HBR_SIRV, pediatric glioma sample	Iso-Seq with Kinnex*	Differentially expressed novel isoforms

\*Includes use of the original or modified versions of the MAS-Iso-Seq method, the MAS-Seq for 10x Single Cell 3' kit, and all Kinnex RNA kits.

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© 2025 PacBio. All rights reserved. Research use only. Not for use in diagnostic procedures.

102-326-634 REV01 MAR2025

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- Al'Khafaji, A. M., et al. (2023). High-throughput RNA isoform sequencing using programmed cDNA concatenation. *Nature Biotechnology*, 1-5. <https://doi.org/10.1038/s41587-023-01815-7>
- Anczukow, O., et al. (2024). Steering research on mRNA splicing in cancer towards clinical translation. *Nature Reviews Cancer*, 24(12), 887-905. <https://doi.org/10.1038/s41568-024-00750-2>
- Black et al. (2024). Long-read single-cell RNA sequencing enables the study of cancer subclone-specific genotype and phenotype in chronic lymphocytic leukemia. *BioRxiv*. <https://doi.org/10.1101/2024.03.15.585298>
- Dondi, A., Lischetti, U., et al. (2023). Detection of isoforms and genomic alterations by high-throughput full-length single-cell RNA sequencing in ovarian cancer. *Nature Communications*, 14(1), 7780. <https://doi.org/10.1038/s41467-023-43387-9>
- Huang, K. K., et al. (2021). Long-read transcriptome sequencing reveals abundant promoter diversity in distinct molecular subtypes of gastric cancer. *Genome Biology*, 22, 1-24. <https://doi.org/10.1186/s13059-021-02261-x>
- Kohli, M., et al. (2017). Androgen receptor variant AR-V9 is coexpressed with AR-V7 in prostate cancer metastases and predicts abiraterone resistance. *Clinical Cancer Research*, 23(16), 4704-4715. <https://doi.org/10.1158/1078-0432.CCR-17-0017>
- Li, Z., Zhang, B., Chan, J. J., et al. (2024). An isoform-resolution transcriptomic atlas of colorectal cancer from long-read single-cell sequencing. *Cell Genomics*, 4(9). <https://doi.org/10.1016/j.xgen.2024.100641>
- Nattestad, M., et al. (2018). Complex rearrangements and oncogene amplifications revealed by long-read DNA and RNA sequencing of a breast cancer cell line. *Genome Research*, 28(8), 1126-1135. <https://doi.org/10.1101/gr.231100.117>
- Qin et al. (2024). CTAT-LR-fusion: accurate fusion transcript identification from long and short read isoform sequencing at bulk or single cell resolution. *BioRxiv*. <https://doi.org/10.1101/2024.02.24.581862>

Veiga, D. F., Nesta, et al. (2022). A comprehensive long-read isoform analysis platform and sequencing resource for breast cancer. *Science Advances*, 8(3), eabg6711. <https://doi.org/10.1126/sciadv.abg6711>

Wedemeyer, M. A., et al. (2024). Defining the transcriptome of PIK3CA-altered cells in a human capillary malformation using single cell long-read sequencing. *Scientific Reports*, 14(1), 25440. <https://doi.org/10.1038/s41598-024-72167-8>

Wijeratne, S., et al. (2024). Full-length isoform concatenation sequencing to resolve cancer transcriptome complexity. *BMC Genomics*, 25(1), 122. <https://doi.org/10.1186/s12864-024-10021-x>

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Application note – [Kinnex full-length RNA kit for isoform sequencing](#)

Application note – [Kinnex single-cell RNA kit for single-cell isoform sequencing](#)

Webinar – [The Full-Length Transcriptomic Atlas of Human Colorectal Cancer from Single-Cell Isoform Sequencing](#) (speaker: Zhongxiao Li)

Webinar – [The RNA isoform landscape of cancer](#) (speaker: Olga Anczukow)

Webinar – [Detecting cancer-related RNA dysregulation with long-read sequencing](#) (speaker: Arthur Dondi)

Webinar – [Understanding clonal evolution using game theory and single-cell long-read isoform analysis](#) (speaker: Nathan Salomonis)

Research use only. Not for use in diagnostic procedures. © 2025 Pacific Biosciences of California, Inc. ("PacBio"). All rights reserved. Information in this document is subject to change without notice. PacBio assumes no responsibility for any errors or omissions in this document. Certain notices, terms, conditions and/or use restrictions may pertain to your use of PacBio products and/or third-party products. Refer to the applicable PacBio terms and conditions of sale and to the applicable license terms at [pacb.com/license](https://pacb.com/license). Pacific Biosciences, the PacBio logo, PacBio, Circulomics, Omniome, SMRT, SMRTbell, Iso-Seq, Sequel, Nanobind, SBB, Revio, Onso, Apton, Kinnex, PureTarget, SPRQ, and Vega are trademarks of PacBio.