

# Procedure & Checklist - Preparing SMRTbell® Libraries using PacBio® Barcoded Universal Primers for Multiplexing Amplicons

## Before You Begin

This procedure describes a 2-step PCR process for generating up to 96 barcoded amplicons for multiplexed sequencing on the Sequel System and Sequel II System. The first-round of PCR requires internal primers that contain a combination of universal and target-specific sequences. The target-specific primers tailed with universal sequences are designed and supplied by users. In the second-round of PCR, barcodes are incorporated by using universal sequences tailed with 16 bp PacBio barcode sequences.

A ready-to-use reagent kit containing 96 Barcoded Universal Primers (BUP) in a plate format is commercially available from Pacific Biosciences. Using this kit, barcode sequences are incorporated after the second round of PCR so that the same barcode sequence is found on both ends of the PCR products (Figure 1). Barcoded samples are then pooled as one sample for SMRTbell library construction using the SMRTbell Express Template Prep Kit 2.0. To multiplex >96 samples using the BUP approach, contact PacBio Technical Support ([support@pacb.com](mailto:support@pacb.com)) or your local Field Applications Specialist.

Figure 1 below summarizes the amplification workflow.




Figure 1: Barcodes are incorporated into the PCR amplicon via a two-step PCR approach. The first step is to amplify the region of interest with target-specific primers tailed with Forward and Reverse universal sequences. The second step requires re-amplification of the PCR product with pairs of Forward and Reverse Barcoded Universal Primers from the Barcoded Univ. F/R Primers Plate-96v2 kit. Purple regions correspond to Universal sequences and Blue regions corresponds to 16-bp Barcode sequences. The final barcoded PCR amplicon product contains the same barcode sequence on both ends.

## Workflow




Figure 2: Workflow for preparing multiplexed BUP amplicon libraries using the SMRTbell Express Template Prep Kit 2.0

## Required Materials

Item	Where Used	Vendor	Part Number
Target-Specific F/R Primers tailed with Universal Sequences (Customer-supplied)	PCR Amplification (1st-Round)	Oligo Synthesis Company	N/A
Bar. Univ. F/R Primers Plate-96v2	PCR Amplification (2nd-Round)	PacBio	101-629-100
High Fidelity DNA Polymerase Master Mix (Two recommended mixes are listed below)			
PCR Master Mix Option 1: Phusion® High-Fidelity PCR Master Mix with HF Buffer	PCR Amplification	Thermo Scientific	F531S
PCR Master Mix Option 2: KAPA HiFi HotStart ReadyMix PCR Kit	PCR Amplification	Roche	KK2600, KK2601, or KK2602
SMRTbell® Express Template Prep Kit 2.0	Library Prep	PacBio	100-938-900
AMPure® PB Beads	Library Purification	PacBio	100-265-900
Sequel® System or Sequel® II System Binding and Internal Control Kit (Recommended kits are listed below)			
Sequel® II Bind Kit 2.1 and Int Ctrl 1.0	Sequencing short (<3 kb) amplicon samples on the Sequel II System	PacBio	101-843-000
Sequel® II Bind Kit 2.0 and Int Ctrl 1.0	Sequencing long (≥3 kb) amplicon samples on the Sequel II System	PacBio	101-842-900
Sequel® Binding and Internal Ctrl Kit 3.0	Sequencing amplicon samples on the Sequel System	PacBio	101-626-600
Sequel® System or Sequel® II System Sequencing Kit (Recommended kits are listed below)			
Sequel® II Sequencing Kit 2.0	Supports 4 sequencing reactions on the Sequel II System	PacBio	101-820-200
Sequel® Sequencing Kit 3.0 (8-rxn)	Supports 8 sequencing reactions on the Sequel System	PacBio	101-597-800
Sequel® Sequencing Kit 3.0 (4-rxn)	Supports 4 sequencing reactions on the Sequel System	PacBio	101-597-900
Sequel® System or Sequel® II System SMRT® Cells			
SMRT® Cell 8M Tray	Sequencing on the Sequel II System	PacBio	101-389-001
SMRT® Cell 1M v3 Tray	Sequencing on the Sequel System (Max. 10-hour movie collection time)	PacBio	101-531-000
SMRT® Cell 1M v3 LR Tray	Sequencing on the Sequel System (Max. 20-hour movie collection time)	PacBio	101-531-001
8- or 12-Multichannel Pipettor	High Throughput Pipetting	Any	Vendor-specific
Qubit™ 4 Fluorometer	DNA concentration measurement	ThermoFisher	Q33238
Qubit™ 1x dsDNA HS Assay Kit	DNA concentration measurement	ThermoFisher	Q33230

Table 1. List of Required Materials and Equipment.

## Designing Target-Specific Primers Tailed with Universal Sequences

PacBio highly recommends complying with the following requirements when ordering oligos from your oligo synthesis provider:

1. The 5'-end of the primer **must** be blocked (e.g., 5AmMC6) to prevent amplicons carried over from the first round of PCR from forming SMRTbell templates during library construction. Failure to block the 5'-end reduces yield of barcoded reads.
2. Oligos are desalted.
3. Use the following primer format:

Primer Type	Universal Sequence	Target-Specific Primer	Primer to Order
Forward_Internal_PCR_Primer	/5AmMC6/ gcagtcgaacatgtagctgactcagggtcac	FOR_EXAMPLE1	/5AmMC6/ gcagtcgaacatgtagctgactcagggtcac FOR_EXAMPLE1
Reverse_Internal_PCR_Primer	/5AmMC6/ tggatcactgtgcaagcatcacatcgtag	REV_EXAMPLE2	/5AmMC6/ tggatcactgtgcaagcatcacatcgtag REV_EXAMPLE2

Table 2. Recommended Primer Format for Ordering Oligos

## Best Practices for Generating High-Quality PCR Products

Clean, target-specific PCR products are extremely important for obtaining high-quality sequence data. Non-specific products can represent a substantial percentage of the sequencing reads if they are not removed. To minimize the presence of non-specific products, consider the following recommendations for generating high-quality amplicons suitable for SMRTbell library preparation and sequencing.

1. KAPA HiFi HotStart 2x ReadyMix PCR Reagent. Thaw on ice and mix well before use. Note: All KAPA HiFi HotStart reagents and reactions must be set up and kept on ice until PCR; the high proofreading activity of the enzyme will rapidly degrade primers at room temperature.
2. Begin with high-quality nucleic acids and work in a clean environment.
  - a. If extracted nucleic acids must be stored, freeze at high concentrations in appropriately-buffered solutions.
  - b. To minimize possible contamination and degradation caused by multiple freeze/thaw cycles, sub-aliquot DNA into smaller volumes for storage. For DNA samples, DNASTable<sup>®</sup> Plus from Biomatrica may be used to help preserve extracted DNA.
  - c. Set up PCR reactions in an environment free from sources of non-specific primer and template contaminants; ideally a laminar flow hood, using dedicated pre-PCR pipettor, tips and reagents.
3. Use PCR reagents and conditions for generating target-specific, full-length amplicons.
  - a. Use the highest fidelity polymerase compatible with your PCR amplification system.
  - b. Use desalted or HPLC-purified oligos; damaged bases at the ends of the amplicons cannot be repaired by DNA Damage Repair enzymes.
  - c. Optimize PCR conditions to minimize total time spent at high (>65°C) temperatures, particularly during denaturation.
  - d. PCR extension time should be long enough to ensure complete extension, taking into consideration the polymerase used and target amplicon size. For mixed samples with similar targets, it is important to complete extension at every step to avoid generating chimeric products in subsequent steps. As a general guideline, use extension times of one minute per 1000 base pairs (e.g. 3 minutes for a 3 kb product).
4. Use the lowest number of cycles required for obtaining adequate yields (ng) of PCR products to proceed with SMRTbell library construction. Avoid over-amplification.
5. If non-specific products are present, optimize PCR conditions or perform AMPure PB Bead-based size selection to enrich for PCR amplicons with the desired target size (see recommendations below).

## Recommendations for AMPure PB Bead Purifications

It is highly recommended to purify amplicons before SMRTbell library construction to remove PCR reagents, buffers, primer dimers and short non-specific PCR products. Depending on the amplicon size, the required concentration (volumetric ratio) of AMPure PB beads for purification varies. Use the table below to determine the appropriate concentration of AMPure PB beads to use for your sample.

Amplicon Size	AMPure PB Bead Ratio
250-500 bp	1.8X
500-1000 bp	1.0X
1-3 kb	0.6X
3-10 kb	0.45X
15 kb	0.45X

Table 3. AMPure PB Bead Recommendations Based on Amplicon Size

## DNA Input Requirements for Library Construction

When planning your amplification experiments, always consider the total input DNA amount required for library construction. The table below summarizes the total mass of pooled PCR products required to go into the DNA Damage Repair step for different amplicon sizes.

Amplicon Size	Input DNA Amounts per Pool (ng)
250-500 bp	250 – 500*
500-1000 bp	250 – 500*
1-3 kb	500 – 1000
3-10 kb	1000 – 2000
15 kb	1500 – 3000

\*If input is greater than 500 ng, use 2x Overhang Adapter in the ligation reaction.

Table 4. Input DNA Recommendations Based on Amplicon Size

Note that the input requirements shown in Table 4 refer to the total pooled mass of DNA, not the mass of individual members of the pool. Samples should be present at equimolar concentrations after pooling (see Best Practices for Equimolar Pooling below).

If necessary, replicate PCR reactions should be set up to obtain the required amount of DNA product. This approach also minimizes PCR sampling bias for samples containing heterogeneous templates.

It is also important not to exceed the recommended DNA input amounts shown in Table 4 when working with highly multiplexed samples (>96-plex). For high-multiplexing projects, it may be desirable to pool amplicons, purify with AMPure PB beads, measure the concentration, and then aliquot the appropriate amount for library construction.

## Best Practices for Equimolar Pooling

For studies targeting a single consensus sequence per sample, amplicons may be multiplexed to leverage the throughput capacity of a single SMRT Cell. However, pooling is generally recommended for amplicons of similar sizes (i.e., within +/-15% of the mean size).

1. Ideally, amplicons should be AMPure PB bead purified prior to pooling. (See step 4 on Page 9.)
2. To obtain equal representation of each amplicon in the data, it is important to pool samples in equimolar amounts. To do this, purify with AMPure PB and quantify each amplicon using the Agilent® Bioanalyzer System, Agilent TapeStation or Agilent Fragment Analyzer™ system.
3. Remove non-specific PCR products (contaminating bands) prior to pooling. The presence of non-specific products in the pool will impact sequencing data yield.
  - a. If amplicons contain secondary bands that are <1.5 kb, it may be possible to remove them from >3 kb amplicons using AMPure PB bead purification using an appropriate concentration of beads.
  - b. If the contaminating bands are close in size to or are larger than the desired amplicon, or are greater than 1.5 kb (i.e, they cannot be removed by AMPure PB purification), size selection using an automated size selection tool or other gel-based method may be necessary.
  - c. If removal of contaminating bands is not possible, we recommend re-optimization of the amplification reaction using more stringent PCR conditions.
  - d. Always determine the concentration of the amplicon **target band** or **peak only** and use this value to calculate the mass or volume of the amplicon sample to be used during pooling.
  - e. If presence of contaminating bands is determined to be acceptable (i.e., their presence has minimum impact on the sequencing yield of the desired target), you may choose to include the amplicon in the sample pool. In such cases, however, it may be necessary to increase the relative input amounts of such amplicons (containing non-specific products) during pooling in order to achieve adequate sequencing data yields for each amplicon in the sample pool.
4. For higher multiplexing (e.g. ≥96-plex), purifying and quantifying each individual PCR product may be difficult or impractical. A QC method that may work well is to load samples on an agarose gel to view the PCR products prior to pooling. This QC method may work well if PCR conditions are fully optimized to generate clean specific PCR products consistently.
  - a. PCR products that show the same band intensity on a gel may be pooled by volume or mass. To do this, include control fragments of known concentrations when loading samples to perform agarose gel electrophoresis. The pooled samples must meet the minimum DNA input requirements for SMRTbell library construction listed in Table 4.
  - b. For samples that show weak signals on a gel, increase the volume or mass used during pooling.

## First-Round PCR Using Target Specific Primers Tailed with a Universal Sequence

The procedure below was optimized using high fidelity Phusion Hot Start II for PCR. KAPA HiFi HotStart ReadyMix may also be used. PCR optimizations are highly recommended before proceeding with the second-round PCR step.

1. Prepare a master mix by including 25% overage to account for pipetting errors.

**Note:** Use the master mix set-up below for experiments requiring amplifications of the same target region (same locus) on different samples. For experiments requiring **different** target regions (different loci), individual PCR reactions using locus-specific universal sequence-tailed primers must be performed by following the 1X reaction in the table below.

Component	Stock Concentration	Final Concentration	1X (µL)	N Plex (µL)
Nuclease-Free Water			9.5	N x 9.5 x 1.25
Phusion HSII High Fidelity PCR Master Mix	2X	1X	12.0	N x 12 x 1.25
Universal sequence-tailed Primer 1 <sup>a</sup> and Primer 2 <sup>a</sup>	2 µM each	0.2 µM	2.5	N x 2.5 x 1.25
Total Volume			24.0	

a. 5' end must be blocked (e.g., 5AmMC6)

2. Mix gently by tapping the tube. Quick spin the tube.
3. Transfer 24 µL of the master mix to appropriate PCR tubes.
4. Add 1.0 µL DNA (1 ng/µL) to the PCR mix so that the final reaction volume is 25.0 µL. Optimization is highly recommended to determine the appropriate amount of DNA template to use.
5. Mix gently by tapping the tubes. Quick spin the tubes.
6. Perform amplification using the following cycling parameters. Success of the amplification depends on the specific sample and primers used; therefore, optimizations are highly recommended.

Step	Temperature	Time
1	98°C	30 seconds <sup>a</sup>
2	98°C <sup>b</sup>	15 seconds
3	N <sup>c</sup>	15 seconds
4	72°C	X min <sup>d</sup>
5	Repeat steps 2 to 4 (20 cycles total) <sup>e</sup>	-
6	72°C	7 minutes
7	4°C	Hold

a. Initial denaturation is dependent on template complexity (30 seconds – 3 minutes).

b. Temperature is dependent on the polymerase.

c. N in step 3 is defined by the user depending on Primers 1 and 2.

d. X in step 4 is defined by the user depending on the insert size. 1 min should be allowed for every 1 kb of the target product.

e. The number of cycles in step 5 can be increased to 20-25 cycles for visual inspection on an agarose gel.

7. After amplification, perform visual inspection of the PCR products on agarose gel.
8. Proceed to the second-round of PCR.

## Second-Round PCR Using Barcoded Universal Primers

The procedure below was optimized using high fidelity Phusion Hot Start II for PCR. Other high-fidelity polymerases may be used. Optimizations are highly recommended before proceeding with SMRTbell library construction.

1. Prepare a master mix by including 25% overage to account for pipetting errors. The BUP primers are not added in the master mix. They are added to each tube in step 4 below.

Component	Stock Concentration	Final Concentration	1X	N Plex <sup>a</sup>
Nuclease-free Water			10.75	N x 10.75 x 1.25
Phusion HSII High Fidelity PCR Master Mix	2X	1X	10.75	N x 10.75 x 1.25
Total Volume			21.5	

a. Include 25% overage in the Master Mix to account for pipetting errors.

2. Mix gently by tapping the tube. Quick spin the tube.
3. Transfer 21.5 µL of the master mix to appropriate PCR tubes.
4. Add 2.5 µL of the desired Barcoded Universal Primer pair from the Barcoded Univ. F/R Primers Plate-96v2 to each tube so that the final reaction volume is 24.0 µL.
5. Add 1.0 µL of Round 1 PCR product (1 ng/µL) to the PCR mix. Optimization is highly recommended to determine the appropriate amount of round 1 PCR product to use.
6. Mix gently by tapping the tubes. Quick spin the tubes.
7. Perform amplification using the following cycling parameters. Success of the amplification depends on the specific sample and primers; therefore, optimizations are highly recommended.

Step	Temperature	Time
1	98°C	30 seconds
2	98°C	15 seconds
3	64°C <sup>a</sup>	15 seconds
4	72°C	X min <sup>b</sup>
5	Repeat steps 2 to 4 (20 cycles total) <sup>c</sup>	-
6	72°C	7 minutes
7	4°C	Hold

a. Recommended for Barcoded Universal Primers.

b. Defined by the user depending on the insert size.

c. The number of cycles can be increased to 20-25 cycles for visual inspection on an agarose gel.

8. After amplification, perform visual inspection of the PCR products on an agarose gel.
9. Proceed to pooling.

## Sample Pooling

STEP	✓	Sample Pooling	Notes
1		Refer to Table 4 on page 5 for minimum and maximum DNA input amounts of the pooled sample required for library construction. The total DNA input requirements depend on the average amplicon insert size.	
2		Pool samples in equimolar amounts. See “Best Practices for Equimolar Pooling” for recommendations for pooling samples.	



## Concentrate Pooled PCR Amplicons Using AMPure PB Beads

Refer to Table 3 on page 5 for recommendations for AMPure purification per amplicon size. If the total volume of the pooled sample is <100 µL, bring the volume to a minimum volume of 100 µL with 1X Elution Buffer.

STEP	✓	Concentrate DNA	Notes
1		<p>Add the appropriate volume of AMPure PB beads to the pooled amplicons. The required volume depends on amplicon size. Refer to Table 3.</p> <p>Note that the beads must be brought to room temperature before use and all AMPure PB bead purification steps should be performed at room temperature.</p> <p>Before using, mix the bead reagent well until the solution appears homogenous. Pipette the reagent slowly since the bead mixture is viscous and precise volumes are critical to the purification process.</p>	
2		Pipette mix 15 times. It is important to mix well.	
3		Quickly spin down the tube (for 1 second) to collect the beads.	
4		Incubate the mix on bench top for 5 minutes at room temperature.	
5		Spin down the tube (for 1 second) to collect beads.	
6		Place the tube in a magnetic bead rack until the beads collect to the side of the tube and the solution appears clear. The actual time required to collect the beads to the side depends on the volume of beads added.	
7		<p>With the tube still on the magnetic bead rack, slowly pipette off the cleared supernatant and save in another tube. Avoid disturbing the beads.</p> <p>If the DNA is not recovered at the end of this procedure, you can add equal volumes of AMPure PB beads to the saved supernatant and repeat the AMPure PB bead purification steps to recover the DNA.</p>	
8		<p>Wash the beads with freshly prepared 80% ethanol.</p> <p>Note that 80% ethanol is hygroscopic and should be prepared FRESH to achieve optimal results. Also, 80% ethanol should be stored in a tightly capped polypropylene tube for no more than 3 days.</p> <ul style="list-style-type: none"> <li>– Do not remove the tube from the magnetic rack.</li> <li>– Use a sufficient volume of 80% ethanol to fill the tube (1.5 mL for 1.5 mL tube or 2 mL for 2 mL tube). Slowly dispense the 80% ethanol against the side of the tube opposite the beads.</li> <li>– Do not disturb the beads.</li> <li>– After 30 seconds, pipette and discard the 80% ethanol.</li> </ul>	
9		Repeat <a href="#">step 8</a> .	
10		<p>Remove residual 80% ethanol.</p> <ul style="list-style-type: none"> <li>– Remove the tube from the magnetic bead rack and spin. Both the beads and any residual 80% ethanol will be at the bottom of the tube.</li> <li>– Place the tube back on the magnetic bead rack and allow the beads to separate.</li> <li>– Pipette off any remaining 80% ethanol.</li> </ul>	
11		Check for any remaining droplets in the tube. If droplets are present, repeat <a href="#">step 10</a> .	
12		The volume to use for elution is <b>48 µL</b> (the volume to go into “DNA Damage Repair”) plus 1 µL for QC – see <a href="#">step 14</a> below).	




13		<p>Add the Elution Buffer volume to your beads. Pipette mix 15 times. It is important to mix well.</p> <ul style="list-style-type: none"> <li>- Place at 37°C for 15 minutes to elute the DNA from the beads.</li> <li>- Spin the tube down, then place the tube back on the magnetic bead rack.</li> <li>- Let beads separate fully. Then without disturbing the beads, transfer the supernatant to a new 1.5 ml Lo-Bind tube.</li> <li>- Discard the beads.</li> </ul>	
14		<p>Verify your DNA amount and concentration using a Qubit quantitation platform.</p> <ul style="list-style-type: none"> <li>- Measure the DNA concentration using a Qubit fluorometer.</li> <li>- Using 1 µL of the eluted sample, make a 1:10 dilution in EB.</li> <li>- Use 1 µL of this 1:10 dilution to measure the DNA concentration using the Qubit dsDNA HS Assay kit according to the manufacturer's recommendations.</li> </ul>	
15		<p>Actual recovery per µL and total available sample material: _____</p>	

## SMRTbell Library Construction

The amount of input DNA required in this step depends on the size of the amplicons. See Table 4 for input DNA requirements.

### DNA Damage Repair


1. Prepare the following reaction.

Reagent (Reaction Mix 1)	Tube Cap Color	Volume	✓	Notes
DNA Prep Buffer		7.0 µL		
Pooled and Purified PCR Product		47.0 µL		
NAD		1.0 µL		
DNA Damage Repair Mix v2		2.0 µL		
Total Volume		57.0 µL		

2. Pipette mix 10 times. It is important to mix well.
3. Spin down the contents of the tube with a quick spin in a microfuge.
4. Incubate at 37°C for 30 minutes, then return the reaction to 4°C. Proceed to the next step.

### End-Repair/A-tailing





1. Prepare the following reaction.

Reagent (Reaction Mix 2)	Tube Cap Color	Volume	✓	Notes
Reaction Mix 1		57.0 µL		
End Prep Mix		3.0 µL		
Total Volume		60.0 µL		

2. Pipette mix 10 times. It is important to mix well.
3. Spin down the contents of the tube with a quick spin in a microfuge.
4. Incubate at 20°C for 30 minutes.
5. Incubate at 65°C for 30 minutes, then return the reaction to 4°C. Proceed to the next step.

### Adapter Ligation

1. Prepare the following reaction, adding the components below in the order listed.

Reagent (Reaction Mix 3)	Tube Cap Color	Volume	✓	Notes
Reaction Mix 2		60.0 µL		
Overhang Adapter v3		5.0 µL		
Ligation Mix		30.0 µL		
Ligation Additive		1.0 µL		
Ligation Enhancer		1.0 µL		
Total Volume		97.0 µL		

1. Pipette mix 10 times. It is important to mix well.
2. Spin down the contents of the tube with a quick spin in a microfuge.
3. Incubate at 20°C for 60 minutes, then return the reaction to 4°C. Proceed to the next step.

## Purification of SMRTbell Templates

Bring to a volume of 100  $\mu\text{L}$  with 1X EB before AMPure PB bead purification. For AMPure PB bead concentration guidelines, refer to Table 3.

STEP	✓	Purify DNA	Notes
1		<p>Add the appropriate volume of AMPure PB beads to the pooled amplicons. The required volume depends on amplicon size. Refer to Table 3.</p> <p>Note that the beads must be brought to room temperature before use and all AMPure PB bead purification steps should be performed at room temperature.</p> <p>Before using, mix the bead reagent well until the solution appears homogenous. Pipette the reagent slowly since the bead mixture is viscous and precise volumes are critical to the purification process.</p>	
2		Mix the bead/DNA solution thoroughly by pipette mixing 15 times. It is important to mix well.	
3		Quickly spin down the tube (for 1 second) to collect the beads.	
4		Incubate samples on bench top for 5 minutes at room temperature.	
5		Spin down the tube (for 1 second) to collect beads.	
6		Place the tube in a magnetic bead rack to collect the beads to the side of the tube.	
7		Slowly pipette off cleared supernatant and save (in another tube). Avoid disturbing the beads.	
8		<p>Wash the beads with freshly prepared 80% ethanol.</p> <p>Note that 80% ethanol is hygroscopic and should be prepared FRESH to achieve optimal results. Also, 80% ethanol should be stored in a tightly capped polypropylene tube for no more than 3 days.</p> <ul style="list-style-type: none"> <li>– Do not remove the tube from the magnetic rack.</li> <li>– Use a sufficient volume of 80% ethanol to fill the tube (1.5 mL for a 1.5 mL DNA LoBind tube).</li> <li>– Slowly dispense the 80% ethanol against the side of the tube opposite the beads.</li> <li>– Do not disturb the beads.</li> <li>– After 30 seconds, pipette and discard the 80% ethanol.</li> </ul>	
9		Repeat <a href="#">step 8</a> .	
10		<p>Remove residual 80% ethanol.</p> <ul style="list-style-type: none"> <li>– Remove the tube from the magnetic bead rack and spin. Both the beads and any residual 80% ethanol will be at the bottom of the tube.</li> <li>– Place the tube back on the magnetic bead rack and allow beads to separate.</li> <li>– Pipette off any remaining 80% ethanol.</li> </ul>	
11		Check for any remaining droplets in the tube. If droplets are present, repeat step 10.	
12		<p>Immediately add <b>20 <math>\mu\text{L}</math></b> of Elution Buffer volume to your beads. Pipette mix 15 times. It is important to mix well.</p> <ul style="list-style-type: none"> <li>– Elute the DNA by letting the mix incubate at 37 °C for 15 minutes. This is important to maximize recovery of high molecular weight DNA.</li> <li>– Spin the tube down, then place the tube back on the magnetic bead rack.</li> <li>– Let beads separate fully. Then without disturbing the beads, transfer supernatant to a new 1.5 ml Lo-Bind tube.</li> <li>– Discard the beads.</li> </ul>	
13		<p>Verify your DNA amount and concentration using a Qubit quantitation platform.</p> <ul style="list-style-type: none"> <li>– Measure the DNA concentration using a Qubit fluorometer.</li> <li>– Using 1 <math>\mu\text{L}</math> of the eluted sample, make a 1:10 dilution in EB.</li> <li>– Use 1 <math>\mu\text{L}</math> of this 1:10 dilution to measure the DNA concentration using a Qubit fluorometer and the dsDNA HS Assay kit according to the manufacturer's recommendations.</li> </ul>	
14		<p>Actual recovered DNA SMRTbell concentration (ng/<math>\mu\text{L}</math>): _____</p> <p>Total recovered DNA SMRTbell amount (ng): _____</p>	

## Prepare for Sequencing

### Sequel System:

- For primer annealing, follow the instructions in SMRT Link Sample Setup v8.0 or higher. Use Sequencing Primer v4 for primer annealing.
- For polymerase binding, follow the instructions in SMRT Link Sample Setup v8.0 or higher. Use Sequel I Binding Kit 3.0 for polymerase binding
- For detailed recommendations for sequencing specific amplicon library insert size ranges, refer to the Quick Reference Card – Diffusion Loading and Pre-Extension Time Recommendations for the Sequel System [here](#).

### Sequel II System:

#### Amplicons ≥3 kb:

- For primer annealing, follow the instructions in SMRT Link Sample Setup v8.0 or higher. Use Sequencing Primer v4 for primer annealing.
- For polymerase binding, follow the instructions in SMRT Link Sample Setup v8.0 or higher. Use Sequel II Binding Kit 2.0 for polymerase binding.
- For complex cleanup, PacBio recommends using AMPure PB beads.
- For detailed recommendations for sequencing amplicon libraries >3 kb, refer to the Quick Reference Card – Diffusion Loading and Pre-Extension Time Recommendations for the Sequel II System [here](#).

#### Amplicons <3 kb:

- Sample setup recommendations for sequencing amplicon libraries <3 kb:

Sequencing Primer	Sequencing Primer v4
Binding Kit	Sequel II Binding Kit 2.1
Binding Time	1 Hour
Sample Complex Cleanup	AMPure PB Beads or ProNex® Beads*
Sequencing Kit	Sequel II Sequencing Plate 2.0
On-Plate Loading Conc (pM)	100 pM

\* [ProNex beads](#) can also be used for complex cleanup and may lead to improved library recovery efficiencies for short amplicon samples .

- Recommended Run parameters:

Movie Time (hours)	10
Pre-Extension Time	0.5 Hours
Immobilization Time	2 Hours (default)

- If using ProNex beads for complex cleanup, the instructions provided in the “Purification of Polymerase Bound SMRTbell® Complexes” section of the SMRT Link Sample Setup (v8.0 or higher) calculation worksheet should be modified by following the changes indicated below:

### ProNex Beads

- Equilibrate the ~~AMPure PB beads~~ and Complex Dilution Buffer to room temperature.
- Add Complex Dilution Buffer to binding reaction.

	Barcoded Sample	✓	
Volume of Binding Reaction (uL)	76.0 uL		
Volume of Complex Dilution Buffer (uL)	<del>177.0 uL</del>		Bring to 100 µl
Total Volume (uL)	<del>253.0 uL</del>		100 µl

- Measure the concentration with a Qubit® High Sensitivity kit. Enter concentration and volume of diluted bound complex into the

	Barcoded Sample	✓	
Volume of Diluted Bound Complex (uL)	<span style="background-color: yellow;">    </span> uL		
Bound Complex Concentration (ng/uL)	<span style="background-color: yellow;">    </span> ng/uL		

### ProNex Beads

- Add ~~AMPure PB beads~~ and gently pipette-mix. Pipette-mixing with wide orifice pipette tips is recommended.

	Barcoded Sample	✓	
Volume of <del>AMPure PB</del> beads (uL)	<del>204.0 uL</del>		2X volume or 200 µL

### ProNex Beads

- To bind the polymerase-bound complexes to ~~AMPure PB beads~~, incubate the mixture on the benchtop for 5 minutes. Note: Longe
- Place the tube in a magnetic bead rack until the beads collect to the side of the tube and the solution appears clear. Discard the s
- Immediately resuspend the beads in room temperature Complex Dilution Buffer and pipette-mix. Pipette-mixing with wide orifice

	Barcoded Sample	✓	
Volume of Complex Dilution Buffer (uL)	<del>127.0 uL</del>		30 µL or more

@ 37°C

- To elute the polymerase-bound complexes, incubate the sample on the benchtop for at least 15 minutes at ~~room temperature~~.
- Place the tube in a magnetic bead rack until the beads collect to the side of the tube and the solution appears clear.
- Transfer eluate to a new LoBind microcentrifuge tube. Place the tube (containing the now-purified polymerase-bound complexes)
- Measure the recovered DNA concentration with a Qubit® High Sensitivity kit. Enter concentration and volume of the purified co

## Sequencing Data Analysis

### Demultiplexing Analysis:

- Use the Demultiplex Barcodes analysis application in SMRT Link v8.0 (or higher) to demultiplex barcoded samples.
- Support FASTA file containing 96 PacBio barcode sequences can be downloaded [here](#).
- For details about how to use the Demultiplex Barcodes analysis application, refer to the [SMRT Link User Guide](#).

Revision History (Description)	Version	Date
BETA DRAFT	BETA	May 2019
Initial Release	1	June 2019
Updated to support both Sequel and Sequel II Systems	2	April 2020

For Research Use Only. Not for use in diagnostic procedures. © Copyright 2019 - 2020, Pacific Biosciences of California, Inc. All rights reserved. Information in this document is subject to change without notice. Pacific Biosciences assumes no responsibility for any errors or omissions in this document. Certain notices, terms, conditions and/or use restrictions may pertain to your use of Pacific Biosciences products and/or third party products. Please refer to the applicable Pacific Biosciences Terms and Conditions of Sale and to the applicable license terms at <http://www.pacificbiosciences.com/licenses.html>. Pacific Biosciences, the Pacific Biosciences logo, PacBio, SMRT, SMRTbell, Iso-Seq and Sequel are trademarks of Pacific Biosciences. BluePippin and SageELF are trademarks of Sage Science, Inc. NGS-go and NGSengine are trademarks of GenDx. FEMTO Pulse and Fragment Analyzer are trademarks of Agilent Technologies. All other trademarks are the sole property of their respective owners.