


Extracting HMW DNA from mammalian brain tissue using Nanobind[®] kits

Procedure & checklist



User supplied equipment and reagent list

| Equipment | Model |
|--|--|
| Nanobind® tissue kit | PacBio® (102-302-100) |
| Magnetic tube rack | Thermo Fisher DynaMag-2 (12321D) |
| TissueRuptor II | Qiagen (9002755) |
| Surgical scalpel | Fisher Scientific (22-079-712) |
| ThermoMixer | Eppendorf (5382000023) |
| Platform rocker | Thermo Scientific (M48725Q) |
| Mini-centrifuge | Ohaus Mini-Centrifuge (FC5306) |
| 1.5 mL Protein LoBind microcentrifuge tubes* | Eppendorf (022431081) |
| 2.0 mL Protein LoBind microcentrifuge tubes* | Eppendorf (022431102) |
| 14 mL round bottom tubes | Fisher Scientific (14-956-3B) |
| Wide bore 200 µL pipette tips | USA Scientific (1011-8410) |
| Ethanol (96–100%) | |
| Isopropanol (100%) | |
| UV/Vis | Thermo Fisher Scientific NanoDrop 2000 |
| Fluorescent DNA quantitation | Thermo Qubit 3.0, dsDNA BR and RNA BR Assay Kits |

*Eppendorf Protein LoBind tubes (Eppendorf #022431081 and #022431102) are highly recommended for all extractions to reduce protein contamination from tube carryover. Protein LoBind tubes are more effective in reducing carryover contamination than DNA LoBind tubes or other tubes and will result in improved UV purity.

Kit storage

RNase A and Buffer CT should be stored at 4°C upon arrival.

Nanobind Disks and all other buffers should be stored at room temperature (15–30°C).

Safety precautions

Buffer BL3 and Buffer CW1 contain guanidine hydrochloride. Warning! Guanidine hydrochloride is harmful if swallowed or inhaled and causes skin and eye irritation. DO NOT mix with bleach or acidic solutions.

Product use

Nanobind tissue kits are intended for research use only.

Mammalian brain

This Procedure & checklist describes an example extraction of HMW DNA from mammalian brain tissue. This example provides sample specific details along with any modifications to the standard protocol.

For the standard Nanobind HMW tissue kit protocols, see either the [Procedure & checklist – Extracting HMW DNA from standard Dounce homogenizer tissue using Nanobind kits](#) or [Procedure & checklist – Extracting HMW DNA from animal tissue using TissueRuptor](#). Please use the [Guide & overview – Nanobind tissue kit](#) to determine which protocol is suitable for a given sample type and to find general information regarding the kit.

Sample notes

- The DNA content of brain samples will vary depending on the region of the brain biopsied.
- Many regions of the brain will result in high amounts of DNA per mg input because of regionally high nuclei density; however, some regions such as white matter are comparatively lower in nuclei density.

Protocol notes

- This Procedure & checklist uses the **Nanobind tissue kit** (102-302-100).
- This Procedure & checklist describes DNA extraction from 22 mg of mouse brain cortex using TissueRuptor homogenization for disruption.
- Dounce homogenization would also be fine for brain disruption.
- Brain tissue benefits from centrifugation speeds in steps 6 and 8 that are faster (6,000 x *g*) than in the standard tissue protocol.
- Brain tissue in this Procedure & checklist did not require mincing prior to disruption due to its loose, liquid-like consistency.
- If the brain tissue is sampled from a region with very high nuclei density, the supernatant in step 19 will have very concentrated DNA that forms a gel-like matrix – **this gel-like matrix contains all the DNA**. Transfer the entire gel-like matrix in step 19.

Protocol

1. Place a 14 mL round bottom tube on ice and chill the centrifuge to 4°C.
2. Transfer 22 mg mouse brain cortex to the 14 mL round bottom tube. Keep the tube on ice during the entire disruption process.
 - Our brain sample was frozen and, when it thawed, it had a liquid-like consistency and did not require mincing.
 - If your sample does not have a liquid-like consistency, finely mince the tissue to $\leq 1 \text{ mm}^3$ pieces using a scalpel.
3. Add 750 μL of cold Buffer CT.
 - Buffer CT should be placed on ice when removed from the refrigerator.
4. Submerge the TissueRuptor probe tip in the buffer and blend at max speed for 10 s.
 - Blending the mouse brain cortex generated a moderate amount of foam. All of the foam was transferred in the next step.
5. Transfer homogenate and all foam to a 2 mL Protein LoBind microcentrifuge tube.
6. Pellet homogenate by centrifuging at 6,000 x g and 4°C for 5 min. Discard supernatant.
7. Add 1 mL cold Buffer CT and pipette mix 10X with a wide bore P200 pipette to resuspend tissue.
8. Pellet homogenate by centrifuging at 6,000 x g and 4°C for 5 min. Discard supernatant.
9. Pulse vortex pellet 1s x 2 times (max setting) to dislodge pellet.
10. Add 20 μL of Proteinase K to the previous pellet.
11. Add 150 μL of Buffer CLE3 and pipette mix 10X with a wide bore P200 pipette.
12. Incubate on a ThermoMixer at 55°C and 900 rpm for 30 min.
13. Spin the tube on a mini-centrifuge for 2 s to remove liquid from the cap.
14. Add 20 μL of RNaseA.
15. Incubate on a ThermoMixer at 55°C and 900 rpm for 30 min.
16. Spin the tube on a mini-centrifuge for 2 s to remove liquid from the cap.
17. Add 60 μL of Buffer SB and pulse vortex for 1s x 5 times (max setting) to mix.
18. Centrifuge at 10,000 x g and RT (15–30°C) for 5 min.
19. Transfer up to 250 μL supernatant to a new 1.5 mL Protein LoBind microcentrifuge tube using a wide bore P200 pipette. (Discard the 2 mL Protein LoBind microcentrifuge tube containing the precipitated pellet.)
 - Typical supernatant volumes will be 225–250 μL .
 - A very small pellet was barely visible.
 - Due to the high DNA content of this brain sample, the supernatant contained a gel-like matrix that comprised the DNA. All of the matrix was transferred.
20. Add 50 μL of Buffer BL3 to the previous supernatant and inversion mix 10X.
 - The solution became cloudy but cleared up in step 23.

Quick tip

Thorough tissue disruption is key to efficient lysis. It is also important to keep the tissue cold during the entire disruption process.

Quick tip

The 2 mL tube is essential for efficient lysis because of its shape; the narrow taper of a 1.5 mL tube prevents proper mixing of the lysate during subsequent thermomixing.

Quick tip

If there are still visible, undigested tissue pieces after step 12, the incubation may be extended up to 2 h. However, if tissue is appropriately disrupted in steps 1-4, then 30 min should be sufficient.

Quick tip

The narrow taper of the 1.5 mL tube is essential for proper removal of wash buffer in steps 31 & 32 and for thorough recovery of eluate in step 35.

21. Spin the tube on a mini-centrifuge for 2 s to remove liquid from the cap.
22. Add Nanobind disk to lysate and add 300 μ L of isopropanol. Inversion mix 10X.
 - The Nanobind disk must be added before isopropanol.
 - A large, cloudy mass appeared upon addition of isopropanol and inversion mixing; this adhered to the Nanobind disk and became clear during the next step.
23. Mix on a platform rocker at 20 rpm for 15 min at RT.
24. Place tube rack on the magnetic base using the method described in the [Guide & overview – Nanobind tissue kit](#) Magnetic Rack Handling Procedure section.
25. Discard the supernatant with a pipette using the method described in the [Guide & overview – Nanobind tissue kit](#) Pipetting section, taking care to avoid pipetting the DNA or contacting the Nanobind disk.
26. Add 500 μ L of Buffer CW1, remove tube rack from magnetic base, inversion mix 4X, replace the tube rack on the magnetic base and discard the supernatant.
27. Repeat step 26.
28. Add 500 μ L of Buffer CW2, remove tube rack from magnetic base, inversion mix 4X, replace the tube rack on the magnetic base, and discard the supernatant.
29. Repeat step 28.
30. Pipette out any residual liquid from the tube cap.
31. Spin the tube on a mini-centrifuge for 2 s. With the tube rack already on the magnetic base and right-side-up, place tube on tube rack and remove residual liquid.
 - If the Nanobind disk is blocking the bottom of the tube, gently push it aside with the tip of the pipette towards the magnet.
32. Repeat step 31.
33. Add 75 μ L of Buffer EB directly onto the Nanobind disk and incubate at RT for 10 min.
 - The Nanobind disk does not need to be fully immersed in Buffer EB – it need only be wetted and sitting atop the liquid.
34. Collect DNA by transferring eluate to a new 1.5 mL microcentrifuge tube using a wide bore P200 pipette.
 - Either Protein LoBind or DNA LoBind tubes can be used in this step.
 - Avoid Axygen tubes as these have been shown to interfere with PacBio sequencing.
35. Spin the tube containing the Nanobind disk on a mini-centrifuge for 5 s. Use a standard P200 pipette to combine any additional liquid that comes off the disk with the previous eluate. Repeat if necessary.
 - For 22 mg of brain cortex, step 35 had to be performed 3 times to get all the DNA off of the disk.
 - Tissue types with high DNA content such as certain brain regions may result in a dense, clear gel that adheres strongly to the Nanobind disk. **This clear gel is DNA.** For these tissue types, this spin step is critical for recovering all the DNA. Repeat until all the clear gel has spun off of the Nanobind disk.
36. Pipette mix 5X with a standard P200 pipette to homogenize the eluate and disrupt any unsolubilized “jellies” that may be present.

Quick tip

The Nanobind disk only needs to be wetted in the elution step: **THE DISK DOES NOT NEED TO BE FULLY SUBMERGED IN BUFFER EB.**

Quick tip

This 5 s spin is **CRITICAL** for recovering the DNA. We do not recommend a 2nd elution.

- Limited pipette mixing will not noticeably reduce DNA size or sequencing read lengths but is important for accurate quantitation and consistent sequencing performance.
- Take care to disrupt any regions that feel more viscous than other regions.

37. Let eluate rest overnight at RT to allow DNA to solubilize.

- Visible “jellies” should disperse after resting.
- The extracted HMW DNA can be heterogeneous. This is normal and is one of the challenges of working with HMW DNA. The bigger the DNA, the more this will be apparent.

38. Following overnight rest, pipette mix 5X with a standard P200 pipette and perform triplicate NanoDrop measurements by sampling the top, middle, and bottom of the eluate.

- If the concentration %CV exceeds 30%, 5X pipette mix with a standard P200 pipette and allow DNA to rest at RT for 1 hour to overnight. Take care to disrupt any regions that feel more viscous than other regions. Remeasure with NanoDrop.
- Limited pipette mixing will not noticeably reduce DNA size or sequencing read lengths but is important for accurate quantitation and consistent sequencing performance.
- We routinely see A260/A280 in the range of 1.81–1.90 and A260/A230 in the range of 1.51–2.17 for mouse brain samples.

39. Use Qubit dsDNA BR assay to determine DNA concentration.

- We recommend making multiple measurements from the top, middle, and bottom of the eluate for an accurate DNA concentration reading.

40. Run pulsed field gel electrophoresis (PFGE) or Agilent Femto Pulse to size the HMW DNA.


Quick tip

The DNA will solubilize after resting at RT or by coaxing it into solution using standard P200 pipette mixing. For samples that need to be used immediately, we recommend needle shearing.

DNA extraction yield and purity


- Mouse brain samples result in moderate to very high extraction yields, depending on which portion of the brain is sampled.
- Recovery and purities obtained from mouse brain cortex were very good.

| Sample | Sample input | dsDNA yield | % RNA | A260/A280 | A260/A230 |
|--------------------|--------------|-------------|-------|-----------|-----------|
| Mouse brain cortex | 22 mg | 5.3 µg | 14 | 1.85 | 2.03 |



DNA size

- Size of DNA extracted from mouse brain cortex is 50–300 kb.



PFGE of DNA extracted from mouse brain cortex.

| Revision history (description) | Version | Date |
|---------------------------------------|----------------|---------------|
| Initial release | 01 | July 2022 |
| Minor updates throughout | 02 | December 2022 |

Research use only. Not for use in diagnostic procedures. © 2022 Pacific Biosciences of California, Inc. ("PacBio"). All rights reserved. Information in this document is subject to change without notice. PacBio assumes no responsibility for any errors or omissions in this document. Certain notices, terms, conditions and/or use restrictions may pertain to your use of PacBio products and/or third-party products. Refer to the applicable PacBio terms and conditions of sale and to the applicable license terms at pacb.com/license. Pacific Biosciences, the PacBio logo, PacBio, Circulomics, Omniome, SMRT, SMRTbell, Iso-Seq, Sequel, Nanobind, and SBB are trademarks of PacBio. All other trademarks are the sole property of their respective owners.