

# Generating PureTarget™ repeat expansion panel libraries

## Procedure & checklist

## Before you begin

---

This procedure describes the workflow for generating PureTarget repeat expansion libraries.

### Overview

Samples	8–24 (processed in batches of 8)
Library prep time	8 hours +/- 2 hours for up to 24 samples
Annealing binding clean up time	1 hour +/- 10 minutes

### DNA input


Quantity	2 µg in Buffer EB, TE buffer (pH 8), or nuclease-free water
DNA size distribution	50% ≥30 kb

### Sample multiplexing

Sequel® II systems	Up to 24 samples
Revo™ system	Up to 48 samples

## Workflow overview

---



## Required materials and equipment

DNA sizing	
Femto Pulse system	Agilent Technologies, Inc. M5330AA
Femto Pulse gDNA 165kb analysis kit	Agilent Technologies, Inc. FP-1002-0275
DNA quantitation	
Qubit fluorometer	ThermoFisher Scientific Q33238
Qubit 1X dsDNA HS assay kit	ThermoFisher Scientific Q33230
Target enrichment and library preparation	
PureTarget™ repeat expansion panel	PacBio® 103-390-400
PureTarget™ beads kit (store at 4C upon arrival)	PacBio® 103-234-800*
<ul style="list-style-type: none"> <li>• SMRT boost beads</li> <li>• SMRT boost buffer</li> </ul>	
PureTarget™ repeat expansion targeting kit (store at -20C upon arrival)	PacBio® 103-234-700*
<ul style="list-style-type: none"> <li>• Cas9 Buffer</li> <li>• Phosphatase</li> <li>• Cas9 Nuclease</li> <li>• dA tail buffer</li> <li>• dATP (10mM)</li> <li>• <i>Taq</i> DNA polymerase</li> <li>• gRNA mix (5 µM)</li> </ul>	
SMRTbell® prep kit 3.0 (store at -20C upon arrival)	
<ul style="list-style-type: none"> <li>• Repair buffer</li> <li>• Ligation mix</li> <li>• Ligation enhancer</li> <li>• Nuclease buffer</li> <li>• Nuclease mix</li> <li>• <i>End repair mix (not used in this protocol)</i></li> <li>• <i>DNA repair mix (not used in this protocol)</i></li> </ul>	PacBio® 102-141-700*
SMRTbell® Cleanup Beads	PacBio® 102-158-300*
Low TE buffer (pH 8.0)	PacBio® 102-178-400*
SMRTbell® adapter index plate 96A (for barcoding)	PacBio® 102-009-200
Revio™ polymerase kit (store at -20C upon arrival)	PacBio® 102-793-100
<ul style="list-style-type: none"> <li>• Annealing buffer</li> </ul>	

- Standard sequencing primer
- *Kinnex™ sequencing primer (not used in this protocol)*
- Polymerase buffer
- Sequencing polymerase
- Dilution buffer
- Sequencing control
- Loading buffer

Sequel® II binding kit 3.2 (store at -20C upon arrival)

- Sequel® II annealing buffer
- Sequel® II primer 3.2
- *Kinnex™ sequencing primer (not used in this protocol)*
- Sequel® II polymerase dilution buffer
- Sequel® II DNA polymerase 2.2
- Sequel® II ABC buffer
- Sequel® II DNA internal control complex 3.2
- Sequel® II loading buffer 3.2

PacBio® 102-194-100

200 Proof ethanol, molecular biology or ACS grade

Any major lab supplier (MLS)

Mini-tube rotator

Any MLS (e.g., Fisher Scientific 05-450-127)

8-channel pipettes

Any MLS

0.2 mL 8-tube strips

USA Scientific TempAssure 1402-4708

Microcentrifuge

Any MLS

Magnetic separation rack compatible with 0.2 mL 8-tube strips

Any MLS (e.g., V&P Scientific VP 772F4-1 8-strip or Permagen MSR812 24-strip)

Magnetic separation rack compatible with 1.5 mL tubes

Any MLS (e.g., Invitrogen DynaMag™-2 Magnet 12321D)

Thermocycler compatible with 0.2 mL tube strips

Any MLS

Nuclease-free (NF) water

Any MLS, molecular biology grade

1.5 mL DNA LoBind tubes

Eppendorf 022431021

\*Sold together as part of the PureTarget™ repeat expansion panel kit (103-390-400)

# General best practices

---

## DNA input

For optimal performance, this protocol requires high-quality, high molecular weight (HMW) human gDNA with at least 50% of the mass of DNA in molecules at  $\geq 30$  kb in length, or genome quality number (GQN) of  $\geq 5$  at 30 kb based on the Agilent Femto Pulse system.

We recommend starting with high-quality, HMW genomic DNA extracted with the Nanobind<sup>®</sup> CBB kit (PacBio 102-301-900) for DNA from human cells. For human whole blood we recommend using the Nanobind PanDNA kit (PacBio 103-260-000) and RBC lysis as the extraction method.

The recommended total mass of DNA is **2  $\mu$ g per sample** to ensure there are sufficient gene copies to load and maximize sequencing coverage. This protocol is suitable for 1–4  $\mu$ g per sample.

## Multiplexing samples

### Sequel II/ Sequel IIe systems

Up to 24 samples can be barcoded and sequenced per SMRT<sup>®</sup> Cell on the Sequel II and Sequel IIe systems. These samples should be processed in batches of 8, 16, or 24 samples.

### Revio system

Up to 48 samples can be barcoded and sequenced per SMRT Cell on the Revio system. These samples should be processed in batches of 8, 16, 24, 32, 40, or 48 samples.

For Revio runs with 8, 16, or 24 samples, follow this protocol as is.

For Revio runs with 32, 40, or 48 samples, process two workflows in parallel: 1) 24 samples following this protocol and 2) an additional 8, 16, or 24 samples following this protocol.

**Each PureTarget repeat expansion panel kit supports the preparation of 24 samples. Therefore, two library prep kits are required to generate a Revio run of 32, 40, or 48 samples.**

## Reagent and sample handling

Room temperature is defined as any temperature in the range of 18–23°C for this protocol.

Mix all reagents well prior to use. Vortex-mix all buffers prior to use. Do not vortex enzymes.

Quick-spin all reagents in a microcentrifuge to collect liquid at the bottom prior to use.

Keep all temperature-sensitive reagents on ice prior to use. Combine and mix master mixes involving temperature-sensitive reagents at room temperature and keep on ice until use.

## Safety precautions

Refer to the Safety Data Sheet (SDS) for information on reagent hazards and protocols for safe handling, use, storage, and disposal.

Temperature-sensitive reagents		
Step used	Tube	Reagent
Dephosphorylation	Blue	Phosphatase
Cas9 digestion	Green	Cas9 Nuclease
Cas9 digestion	Purple	gRNA mix
dA Tailing	Light Blue	<i>Taq</i> DNA Polymerase
dA Tailing	Yellow	dATP
Adapter ligation		SMRTbell adapter plate
Adapter ligation	Yellow	Ligation mix
Adapter ligation	Red	Ligation enhancer
Nuclease treatment	Light green	Nuclease mix

Bring SMRTbell cleanup beads, SMRT boost beads, SMRT boost buffer, and Qubit reagents to room temperature for 30 minutes prior to use.

Pipette-mix all reactions by pipetting up and down 10 times. Use full-volume pipette mixing of all reactions to ensure thorough mixing of all reaction components.

Samples can be stored at -20°C at all safe stopping points listed in the protocol.

# Thermocycler conditions

---

Program thermocycler(s) prior to beginning the protocol for the first time.

## 1. Dephosphorylation program

Set the lid temperature to 95°C. If the lid temperature is not programmable, it is acceptable to leave at 95–105°C.

Step	Time	Temperature
1	10 min	37°C
2	3 min	80°C
3	Hold	4°C

## 2. Cas9 gRNA complex formation program

For this and the rest of thermocycler programs, set the lid temperature to 75°C. If the lid temperature is not programmable, it is acceptable to leave at 95–105°C.

Step	Time	Temperature
1	10 min	37°C
2	Hold	4°C

## 3. Cas9 digestion program

Step	Time	Temperature
1	60 min	37°C
2	Hold	4°C

## 4. dA-tail program

Step	Time	Temperature
1	20 min	37°C
2	5 min	72°C
3	Hold	4°C

## 5. Adapter ligation program

Step	Time	Temperature
1	30 min	20°C
2	Hold	4°C

## 6. Nuclease treatment program

Step	Time	Temperature
1	60 min	37°C
2	Hold	4°C

## Workflow steps

### 1. Input DNA quality control and dilution

Before you begin, evaluate the quantity and size distribution of input DNA using Agilent Femto Pulse system to determine whether it is suitable for the protocol.

✓	Step	Instructions
	1.1	Bring the Qubit 1X dsDNA HS working solution and standards to <b>room temperature</b> . Pulse vortex and/or pipette-mix each sample 5 times to homogenize the DNA in solution.
	1.2	For viscous input DNA, it is important to homogenize the extracted DNA prior to start of the protocol. To homogenize the DNA, pulse-vortex 5 times and/or pipette-mix full sample volume 5 times, up and down with standard (not wide bore) tips. These steps will maintain HMW of your DNA but will improve accuracy of quantification and subsequent handling.
	1.3	Quick-spin each sample to collect liquid.
	1.4	Take a <b>1 µL</b> aliquot from each sample and dilute with <b>9 µL</b> of <b>elution buffer or water</b> .
	1.5	Measure DNA concentration with a Qubit fluorometer using the <b>1X dsDNA HS kit</b> .
	1.6	Dilute each aliquot to <b>250 pg/µL</b> in Femto Pulse dilution buffer based on the Qubit reading.
	1.7	Measure DNA size distribution with a Femto Pulse system using the <b>gDNA 165 kb analysis kit</b> .
	1.8	Aliquot or bring 1–4 µg DNA to a final volume of <b>67 µL per sample</b> with nuclease-free water and transfer to an 8 tube PCR strip.

**SAFE STOPPING POINT - Store at 4°C**



## 2. Dephosphorylation

This step enables dephosphorylation of genomic DNA 5' and 3' ends, which prevents subsequent adapter ligation to non-targeted genomic DNA ends.

✓	Step	Instructions																								
		Add the following components in the order and volume listed below to a new microcentrifuge tube. Adjust component volumes for the number of samples being prepared, plus 15% overage according to the table below.																								
		<b>Reaction Mix 1 (RM1)</b>																								
2.1	✓	<table border="1"> <thead> <tr> <th>Tube color</th> <th>Component</th> <th>Volume per sample</th> <th>Volume per 8-plex <small>With 15% overage</small></th> <th>Volume per 16-plex <small>With 15% overage</small></th> <th>Volume per 24-plex <small>With 15% overage</small></th> </tr> </thead> <tbody> <tr> <td>Red</td> <td>Cas9 buffer</td> <td>8 <math>\mu</math>L</td> <td>73.6 <math>\mu</math>L</td> <td>147.2 <math>\mu</math>L</td> <td>220.8 <math>\mu</math>L</td> </tr> <tr> <td>Blue</td> <td>Phosphatase</td> <td>5 <math>\mu</math>L</td> <td>46 <math>\mu</math>L</td> <td>92 <math>\mu</math>L</td> <td>138 <math>\mu</math>L</td> </tr> <tr> <td colspan="2"><b>Total volume</b></td> <td><b>13 <math>\mu</math>L</b></td> <td><b>119.6 <math>\mu</math>L</b></td> <td><b>239.2 <math>\mu</math>L</b></td> <td><b>358.8 <math>\mu</math>L</b></td> </tr> </tbody> </table>	Tube color	Component	Volume per sample	Volume per 8-plex <small>With 15% overage</small>	Volume per 16-plex <small>With 15% overage</small>	Volume per 24-plex <small>With 15% overage</small>	Red	Cas9 buffer	8 $\mu$ L	73.6 $\mu$ L	147.2 $\mu$ L	220.8 $\mu$ L	Blue	Phosphatase	5 $\mu$ L	46 $\mu$ L	92 $\mu$ L	138 $\mu$ L	<b>Total volume</b>		<b>13 <math>\mu</math>L</b>	<b>119.6 <math>\mu</math>L</b>	<b>239.2 <math>\mu</math>L</b>	<b>358.8 <math>\mu</math>L</b>
Tube color	Component	Volume per sample	Volume per 8-plex <small>With 15% overage</small>	Volume per 16-plex <small>With 15% overage</small>	Volume per 24-plex <small>With 15% overage</small>																					
Red	Cas9 buffer	8 $\mu$ L	73.6 $\mu$ L	147.2 $\mu$ L	220.8 $\mu$ L																					
Blue	Phosphatase	5 $\mu$ L	46 $\mu$ L	92 $\mu$ L	138 $\mu$ L																					
<b>Total volume</b>		<b>13 <math>\mu</math>L</b>	<b>119.6 <math>\mu</math>L</b>	<b>239.2 <math>\mu</math>L</b>	<b>358.8 <math>\mu</math>L</b>																					
2.2		Pipette-mix <b>RM1</b> up and down 10 times (do not vortex). Pipette-mix full volume.																								
2.3		Quick-spin <b>RM1</b> in a microcentrifuge to collect liquid.																								
2.4		Add <b>13 <math>\mu</math>L</b> of the <b>RM1</b> to each sample from step 1.8 for total reaction volume of <b>80 <math>\mu</math>L</b> .																								
2.5		Pipette-mix each sample up and down 10 times. Pipette-mix full volume.																								
2.6		Quick-spin in a microcentrifuge to collect liquid.																								
		Run the <b>dephosphorylation</b> thermocycler program. Set the lid temperature to 95°C. If the lid temperature is not programmable, it is acceptable to leave at 95–105°C.																								
2.7		<table border="1"> <thead> <tr> <th>Step</th> <th>Time</th> <th>Temperature</th> </tr> </thead> <tbody> <tr> <td>1</td> <td>10 min</td> <td>37°C</td> </tr> <tr> <td>2</td> <td>3 min</td> <td>80°C</td> </tr> <tr> <td>3</td> <td>Hold</td> <td>4°C</td> </tr> </tbody> </table>	Step	Time	Temperature	1	10 min	37°C	2	3 min	80°C	3	Hold	4°C												
Step	Time	Temperature																								
1	10 min	37°C																								
2	3 min	80°C																								
3	Hold	4°C																								
2.8		Proceed to the next step of the protocol.																								

### 3. Cas9 gRNA complex formation

This step enables Cas9-gRNA ribonucleoprotein complex formation required for subsequent targeting and digestion of genomic DNA in Step 4.

✓	Step	Instructions																																																								
		Add the following components in the order and volume listed below to a new PCR tube. Adjust component volumes for the number of samples being prepared, plus 25% overage according to the table below.																																																								
		<table border="1"> <thead> <tr> <th colspan="7">Reaction Mix 2 (RM2)</th> </tr> <tr> <th>✓</th> <th>Reagent</th> <th>Component</th> <th>Volume per sample</th> <th>Volume per 8-plex</th> <th>Volume per 16-plex</th> <th>Volume per 24-plex</th> </tr> </thead> <tbody> <tr> <td></td> <td></td> <td></td> <td></td> <td>With 25% overage</td> <td>With 25% overage</td> <td>With 25% overage</td> </tr> <tr> <td>3.1</td> <td></td> <td>NF water</td> <td>1.75 µL</td> <td>17.5 µL</td> <td>35 µL</td> <td>52.5 µL</td> </tr> <tr> <td></td> <td>Red</td> <td>Cas9 buffer</td> <td>0.5 µL</td> <td>5 µL</td> <td>10 µL</td> <td>15 µL</td> </tr> <tr> <td></td> <td>Green</td> <td>Cas9 Nuclease</td> <td>0.5 µL</td> <td>5 µL</td> <td>10 µL</td> <td>15 µL</td> </tr> <tr> <td></td> <td>Purple</td> <td>gRNA mix (5 µM)</td> <td>2 µL</td> <td>20 µL</td> <td>40 µL</td> <td>60 µL</td> </tr> <tr> <td></td> <td></td> <td><b>Total volume</b></td> <td><b>4.75 µL</b></td> <td><b>47.5 µL</b></td> <td><b>95 µL</b></td> <td><b>142.5 µL</b></td> </tr> </tbody> </table>	Reaction Mix 2 (RM2)							✓	Reagent	Component	Volume per sample	Volume per 8-plex	Volume per 16-plex	Volume per 24-plex					With 25% overage	With 25% overage	With 25% overage	3.1		NF water	1.75 µL	17.5 µL	35 µL	52.5 µL		Red	Cas9 buffer	0.5 µL	5 µL	10 µL	15 µL		Green	Cas9 Nuclease	0.5 µL	5 µL	10 µL	15 µL		Purple	gRNA mix (5 µM)	2 µL	20 µL	40 µL	60 µL			<b>Total volume</b>	<b>4.75 µL</b>	<b>47.5 µL</b>	<b>95 µL</b>	<b>142.5 µL</b>
Reaction Mix 2 (RM2)																																																										
✓	Reagent	Component	Volume per sample	Volume per 8-plex	Volume per 16-plex	Volume per 24-plex																																																				
				With 25% overage	With 25% overage	With 25% overage																																																				
3.1		NF water	1.75 µL	17.5 µL	35 µL	52.5 µL																																																				
	Red	Cas9 buffer	0.5 µL	5 µL	10 µL	15 µL																																																				
	Green	Cas9 Nuclease	0.5 µL	5 µL	10 µL	15 µL																																																				
	Purple	gRNA mix (5 µM)	2 µL	20 µL	40 µL	60 µL																																																				
		<b>Total volume</b>	<b>4.75 µL</b>	<b>47.5 µL</b>	<b>95 µL</b>	<b>142.5 µL</b>																																																				
3.2		Pipette-mix <b>RM2</b> up and down 10 times (do not vortex). Pipette-mix full volume.																																																								
3.3		Ensure <b>RM2</b> is fully mixed and quick-spin in a microcentrifuge to collect liquid. <b>Proceed to incubating the RM2 Master Mix using thermocycler program below. The full master mix volume for 8, 16, or 24 samples should be incubated in one PCR tube.</b>																																																								
		Run the <b>Cas9 gRNA complex formation</b> thermocycler program.																																																								
3.4		Set the lid temperature to 75°C. If the lid temperature is not programmable, it is acceptable to leave at 95–105°C.																																																								
		<table border="1"> <thead> <tr> <th>Step</th> <th>Time</th> <th>Temperature</th> </tr> </thead> <tbody> <tr> <td>1</td> <td>10 min</td> <td>37°C</td> </tr> <tr> <td>2</td> <td>Hold</td> <td>4°C</td> </tr> </tbody> </table>	Step	Time	Temperature	1	10 min	37°C	2	Hold	4°C																																															
Step	Time	Temperature																																																								
1	10 min	37°C																																																								
2	Hold	4°C																																																								
3.5		Proceed to the next step of the protocol.																																																								

## 4. Cas9 digestion and cleanup

This step enables digestion of double-stranded DNA at targeted regions using Master Mix RM2 prepared in Step 3.

✓	Step	Instructions																														
		Add the following components in the order and volume listed below to a new microcentrifuge tube. Adjust component volumes for the number of samples being prepared, plus 15% overage according to the table below. <b>RM2 is Cas9 gRNA complex from step 3 after thermocycler program.</b>																														
		<b>Reaction Mix 3 (RM3)</b>																														
4.1	✓	<table border="1"> <thead> <tr> <th>Reagent</th> <th>Component</th> <th>Volume per sample</th> <th>Volume per 8-plex</th> <th>Volume per 16-plex</th> <th>Volume per 24-plex</th> </tr> </thead> <tbody> <tr> <td></td> <td>NF water</td> <td>12.75 <math>\mu</math>L</td> <td>117.3 <math>\mu</math>L</td> <td>234.6 <math>\mu</math>L</td> <td>351.9 <math>\mu</math>L</td> </tr> <tr> <td>Red</td> <td>Cas9 buffer</td> <td>1.5 <math>\mu</math>L</td> <td>13.8 <math>\mu</math>L</td> <td>27.6 <math>\mu</math>L</td> <td>41.4 <math>\mu</math>L</td> </tr> <tr> <td></td> <td>RM2 (Cas9 gRNA complex from Step 3)</td> <td>4.75 <math>\mu</math>L</td> <td>43.7 <math>\mu</math>L</td> <td>87.4 <math>\mu</math>L</td> <td>131.1 <math>\mu</math>L</td> </tr> <tr> <td colspan="2"><b>Total volume</b></td> <td><b>19 <math>\mu</math>L</b></td> <td><b>174.8 <math>\mu</math>L</b></td> <td><b>349.6 <math>\mu</math>L</b></td> <td><b>524.4 <math>\mu</math>L</b></td> </tr> </tbody> </table>	Reagent	Component	Volume per sample	Volume per 8-plex	Volume per 16-plex	Volume per 24-plex		NF water	12.75 $\mu$ L	117.3 $\mu$ L	234.6 $\mu$ L	351.9 $\mu$ L	Red	Cas9 buffer	1.5 $\mu$ L	13.8 $\mu$ L	27.6 $\mu$ L	41.4 $\mu$ L		RM2 (Cas9 gRNA complex from Step 3)	4.75 $\mu$ L	43.7 $\mu$ L	87.4 $\mu$ L	131.1 $\mu$ L	<b>Total volume</b>		<b>19 <math>\mu</math>L</b>	<b>174.8 <math>\mu</math>L</b>	<b>349.6 <math>\mu</math>L</b>	<b>524.4 <math>\mu</math>L</b>
Reagent	Component	Volume per sample	Volume per 8-plex	Volume per 16-plex	Volume per 24-plex																											
	NF water	12.75 $\mu$ L	117.3 $\mu$ L	234.6 $\mu$ L	351.9 $\mu$ L																											
Red	Cas9 buffer	1.5 $\mu$ L	13.8 $\mu$ L	27.6 $\mu$ L	41.4 $\mu$ L																											
	RM2 (Cas9 gRNA complex from Step 3)	4.75 $\mu$ L	43.7 $\mu$ L	87.4 $\mu$ L	131.1 $\mu$ L																											
<b>Total volume</b>		<b>19 <math>\mu</math>L</b>	<b>174.8 <math>\mu</math>L</b>	<b>349.6 <math>\mu</math>L</b>	<b>524.4 <math>\mu</math>L</b>																											
4.2		Pipette-mix <b>RM3</b> up and down 10 times (do not vortex). Pipette-mix full volume.																														
4.3		Quick-spin in a microcentrifuge to collect liquid.																														
4.4		Add <b>19 <math>\mu</math>L</b> of the RM3 to each sample from Step 2 for total reaction volume of <b>99 <math>\mu</math>L</b> .																														
4.5		Pipette-mix each sample up and down 10 times. Pipette-mix full volume.																														
4.6		Spin down the tube strip for 15–30 seconds in a microcentrifuge to collect liquid and remove bubbles.																														
		Run the <b>Cas-9 digestion</b> thermocycler program.																														
4.7		Set the lid temperature to 75°C. If the lid temperature is not programmable, it is acceptable to leave at 95–105°C.																														
		<table border="1"> <thead> <tr> <th>Step</th> <th>Time</th> <th>Temperature</th> </tr> </thead> <tbody> <tr> <td>1</td> <td>60 min</td> <td>37°C</td> </tr> <tr> <td>2</td> <td>Hold</td> <td>4°C</td> </tr> </tbody> </table>	Step	Time	Temperature	1	60 min	37°C	2	Hold	4°C																					
Step	Time	Temperature																														
1	60 min	37°C																														
2	Hold	4°C																														
4.8		Quick-spin in a microcentrifuge to collect liquid.																														
		<b>Cleanup with 1X SMRTbell cleanup beads</b>																														
4.9		Add <b>100 <math>\mu</math>L</b> of resuspended, room-temperature SMRTbell cleanup beads to each sample.																														
4.10		Pipette-mix the beads slowly 8–10x until evenly distributed. <b>Bead clumping may occur and is not a concern, however, avoid over-pipetting as it may cause DNA/bead mixture to stick to the pipette tip.</b>																														
4.11		Quick-spin the tube strip in a microcentrifuge to collect all liquid from the sides of the tubes.																														

- 4.12 Leave at **room temperature** for **10 minutes** to allow DNA to bind beads.
- 4.13 Place the tube strip in a magnetic separation rack for 3–5 minutes until beads separate fully from the solution.
- 4.14 Slowly pipette off the cleared supernatant without disturbing the beads. Discard the supernatant.
- 4.15 Slowly dispense **200  $\mu$ L**, or enough to cover the beads, of **freshly prepared 80% ethanol** into each tube. After **30 seconds**, pipette off the 80% ethanol and discard.
- 4.16 Repeat the previous step.
- Remove residual 80% ethanol:
- Remove the tube strip from the magnetic separation rack.
  - Quick-spin the tube strip in a microcentrifuge.
- 4.17
- Place the tube strip back in a magnetic separation rack until beads separate fully from the solution.
  - Pipette off residual 80% ethanol and discard.
- 4.18 Remove the tube strip from the magnetic rack. **Immediately** add **41  $\mu$ L** of **low-TE buffer** to each tube and resuspend the beads.
- 4.17 Quick-spin the tube strip in a microcentrifuge.
- 4.18 Leave at **room temperature** for **5 minutes** to elute DNA.
- 4.19 Place the tube strip in a magnetic separation rack for 3–5 minutes until beads separate fully from the solution.
- 4.20 Slowly pipette off the cleared supernatant without disturbing the beads. Transfer supernatant to a **new tube strip**. Discard old tube strip with beads.
- Optional QC step: Take a **1  $\mu$ L** aliquot from each sample and dilute with **9  $\mu$ L** of **elution buffer or water**. Measure DNA concentration with a Qubit fluorometer using the **1X dsDNA HS kit**.
- 4.21
- Expect recovery of 50-100% per samples relative to starting mass.**
- 4.22 Proceed to the next step of the protocol.

**SAFE STOPPING POINT - Store at -20°C**

## 5. dA-tail

This step enables A-tailing of DNA 3' ends after Cas9-gRNA digestion at targeted regions.

✓	Step	Instructions																																																								
		Add the following components in the order and volume listed below to a new microcentrifuge tube. Adjust component volumes for the number of samples being prepared, plus 15% overage according to the table below.																																																								
		<table border="1"> <thead> <tr> <th colspan="7">Reaction Mix 4 (RM4)</th> </tr> <tr> <th>✓</th> <th>Reagent</th> <th>Component</th> <th>Volume per sample</th> <th>Volume per 8-plex</th> <th>Volume per 16-plex</th> <th>Volume per 24-plex</th> </tr> <tr> <td></td> <td></td> <td></td> <td></td> <td>With 15% overage</td> <td>With 15% overage</td> <td>With 15% overage</td> </tr> </thead> <tbody> <tr> <td></td> <td></td> <td>NF water</td> <td>2.5 <math>\mu</math>L</td> <td>23 <math>\mu</math>L</td> <td>46 <math>\mu</math>L</td> <td>69 <math>\mu</math>L</td> </tr> <tr> <td></td> <td>Orange</td> <td>dA tail buffer</td> <td>5 <math>\mu</math>L</td> <td>46 <math>\mu</math>L</td> <td>92 <math>\mu</math>L</td> <td>138 <math>\mu</math>L</td> </tr> <tr> <td></td> <td>Yellow</td> <td>dATP (10 mM)</td> <td>1.25 <math>\mu</math>L</td> <td>11.5 <math>\mu</math>L</td> <td>23 <math>\mu</math>L</td> <td>34.5 <math>\mu</math>L</td> </tr> <tr> <td></td> <td>Light Blue</td> <td><i>Taq</i> DNA polymerase</td> <td>1.25 <math>\mu</math>L</td> <td>11.5 <math>\mu</math>L</td> <td>23 <math>\mu</math>L</td> <td>34.5 <math>\mu</math>L</td> </tr> <tr> <td></td> <td></td> <td><b>Total volume</b></td> <td><b>10 <math>\mu</math>L</b></td> <td><b>92 <math>\mu</math>L</b></td> <td><b>184 <math>\mu</math>L</b></td> <td><b>276 <math>\mu</math>L</b></td> </tr> </tbody> </table>	Reaction Mix 4 (RM4)							✓	Reagent	Component	Volume per sample	Volume per 8-plex	Volume per 16-plex	Volume per 24-plex					With 15% overage	With 15% overage	With 15% overage			NF water	2.5 $\mu$ L	23 $\mu$ L	46 $\mu$ L	69 $\mu$ L		Orange	dA tail buffer	5 $\mu$ L	46 $\mu$ L	92 $\mu$ L	138 $\mu$ L		Yellow	dATP (10 mM)	1.25 $\mu$ L	11.5 $\mu$ L	23 $\mu$ L	34.5 $\mu$ L		Light Blue	<i>Taq</i> DNA polymerase	1.25 $\mu$ L	11.5 $\mu$ L	23 $\mu$ L	34.5 $\mu$ L			<b>Total volume</b>	<b>10 <math>\mu</math>L</b>	<b>92 <math>\mu</math>L</b>	<b>184 <math>\mu</math>L</b>	<b>276 <math>\mu</math>L</b>
Reaction Mix 4 (RM4)																																																										
✓	Reagent	Component	Volume per sample	Volume per 8-plex	Volume per 16-plex	Volume per 24-plex																																																				
				With 15% overage	With 15% overage	With 15% overage																																																				
		NF water	2.5 $\mu$ L	23 $\mu$ L	46 $\mu$ L	69 $\mu$ L																																																				
	Orange	dA tail buffer	5 $\mu$ L	46 $\mu$ L	92 $\mu$ L	138 $\mu$ L																																																				
	Yellow	dATP (10 mM)	1.25 $\mu$ L	11.5 $\mu$ L	23 $\mu$ L	34.5 $\mu$ L																																																				
	Light Blue	<i>Taq</i> DNA polymerase	1.25 $\mu$ L	11.5 $\mu$ L	23 $\mu$ L	34.5 $\mu$ L																																																				
		<b>Total volume</b>	<b>10 <math>\mu</math>L</b>	<b>92 <math>\mu</math>L</b>	<b>184 <math>\mu</math>L</b>	<b>276 <math>\mu</math>L</b>																																																				
5.1																																																										
5.2		Pipette-mix <b>RM4</b> up and down 10 times (do not vortex). Pipette-mix full volume.																																																								
5.3		Quick-spin <b>RM4</b> in a microcentrifuge to collect liquid.																																																								
5.4		Add <b>10 <math>\mu</math>L</b> of the <b>RM4</b> to <b>40 <math>\mu</math>L</b> of each sample for a total reaction volume of <b>50 <math>\mu</math>L</b> .																																																								
5.5		Pipette-mix each sample up and down 10 times. Pipette-mix full volume.																																																								
5.6		Quick-spin the tube strip in a microcentrifuge to collect liquid.																																																								
		Run the <b>dA-tail</b> thermocycler program.																																																								
		Set the lid temperature to 75°C. If the lid temperature is not programmable, it is acceptable to leave at 95–105°C.																																																								
5.7		<table border="1"> <thead> <tr> <th>Step</th> <th>Time</th> <th>Temperature</th> </tr> </thead> <tbody> <tr> <td>1</td> <td>20 min</td> <td>37°C</td> </tr> <tr> <td>2</td> <td>5 min</td> <td>72°C</td> </tr> <tr> <td>3</td> <td>Hold</td> <td>4°C</td> </tr> </tbody> </table>	Step	Time	Temperature	1	20 min	37°C	2	5 min	72°C	3	Hold	4°C																																												
Step	Time	Temperature																																																								
1	20 min	37°C																																																								
2	5 min	72°C																																																								
3	Hold	4°C																																																								
5.8		Proceed to the next step of the protocol.																																																								

## 6. Adapter ligation and cleanup

This step ligates the indexed SMRTbell adapter to the ends of each targeted DNA fragment.

✓	Step	Instructions
6.1		To a PCR strip with 50 $\mu\text{L}$ DNA from Step 5, add <b>4 <math>\mu\text{L}</math> SMRTbell barcoded adapter 3.0</b> to each sample. Tap-mix or pipette up and down 10 times (do not vortex).

Add the following components in the order and volume listed below to a new microcentrifuge tube. Adjust component volumes for the number of samples being prepared, plus 15% overage according to the table below.

Aspirate and dispense viscous ligation mix slowly to reduce liquid sticking to the exterior and interior wall of the pipette tip.

Reaction Mix 5 (RM5)							
6.2	✓	Reagent	Component	Volume per sample	Volume Per 8-plex With 15% overage	Volume Per 16-plex With 15% overage	Volume Per 24-plex With 15% overage
			NF water	2 $\mu\text{L}$	18.4 $\mu\text{L}$	36.8 $\mu\text{L}$	55.2 $\mu\text{L}$
		Purple	Repair buffer	8 $\mu\text{L}$	73.6 $\mu\text{L}$	147.2 $\mu\text{L}$	220.8 $\mu\text{L}$
		Yellow	Ligation mix	30 $\mu\text{L}$	276 $\mu\text{L}$	552 $\mu\text{L}$	828 $\mu\text{L}$
		Red	Ligation enhancer	1 $\mu\text{L}$	9.2 $\mu\text{L}$	18.4 $\mu\text{L}$	27.6 $\mu\text{L}$
			<b>Total volume</b>	<b>41 <math>\mu\text{L}</math></b>	<b>377.2 <math>\mu\text{L}</math></b>	<b>754.4 <math>\mu\text{L}</math></b>	<b>1131.6 <math>\mu\text{L}</math></b>

6.3	Pipette-mix <b>RM5</b> up and down 10 times (do not vortex). Pipette-mix full volume. Ensure slow mixing and dispensing of RM5 to prevent liquid adhering to inside wall of tip and volume loss.
-----	--

6.4 Quick-spin **RM5** in a microcentrifuge to collect liquid.

6.5	Add <b>41 <math>\mu\text{L}</math></b> of <b>RM5</b> to each sample from previous step for a total volume of <b>95 <math>\mu\text{L}</math></b> .
-----	---

6.6 Pipette-mix each sample up and down 10 times. Pipette-mix full volume.

6.7	Quick-spin the tube strip in a microcentrifuge to collect liquid.
-----	---

Run the **adapter ligation** thermocycler program.

Set the lid temperature to 75°C. If the lid temperature is not programmable, it is acceptable to leave at 95–105°C.

6.8	Step	Time	Temperature
	1	30 min	20°C
	2	Hold	4°C

### Cleanup with 1X SMRTbell cleanup beads

- 6.9 Add **95  $\mu$ L** of resuspended, room-temperature SMRTbell cleanup beads to each sample.
- 6.10 Pipette-mix the beads 10 times until evenly distributed.
- 6.11 Quick-spin the tube strip in a microcentrifuge to collect all liquid from the sides of the tubes.
- 6.12 Leave at **room temperature** for **10 minutes** to allow DNA to bind beads.
- 6.13 Place the tube strip in a magnetic separation rack for 3–5 minutes until beads separate fully from the solution.
- 6.14 Slowly pipette off the cleared supernatant without disturbing the beads. Discard the supernatant.
- 6.15 Slowly dispense **200  $\mu$ L**, or enough to cover the beads, of **freshly prepared 80% ethanol** into each tube. After **30 seconds**, pipette off the 80% ethanol and discard.
- 6.16 Repeat the previous step.
- Remove residual 80% ethanol:
- Remove the tube strip from the magnetic separation rack.
  - Quick-spin the tube strip in a microcentrifuge.
  - Place the tube strip back in a magnetic separation rack until beads separate fully from the solution.
  - Pipette off residual 80% ethanol and discard.
- 6.17
- 6.18 Remove the tube strip from the magnetic rack. **Immediately** add **41  $\mu$ L** of **elution buffer** to each tube and resuspend the beads.
- 6.19 Quick-spin the tube strip in a microcentrifuge.
- 6.20 Leave at **room temperature** for **5 minutes** to elute DNA.
- 6.21 Place the tube strip in a magnetic separation rack until beads separate fully from the solution.
- 6.22 Slowly pipette off the cleared supernatant without disturbing the beads. Transfer supernatant to a **new tube strip**. Discard the old tube strip with beads.
- 6.23 Optional QC step: Take a **1  $\mu$ L** aliquot from each sample and dilute with **9  $\mu$ L** of **elution buffer or water**. Measure DNA concentration with a Qubit fluorometer using the **1X dsDNA HS kit**.
- Expect recovery of 50-100% per samples relative to starting mass.**
- 6.23 Proceed to the next step of the protocol.

**SAFE STOPPING POINT - Store at -20°C**

## 7. Nuclease treatment

This step removes DNA fragments that have not formed SMRTbell templates. This step removes the vast majority of starting DNA fragments.

✓	Step	Instructions																																										
		Add the following components in the order and volume listed below to a new microcentrifuge tube. Adjust component volumes for the number of samples being prepared, plus 15% overage according to the table below.																																										
		<table border="1"> <thead> <tr> <th colspan="7">Reaction Mix 6 (RM6)</th> </tr> <tr> <th>✓</th> <th>Reagent</th> <th>Component</th> <th>Volume per sample</th> <th>Volume per 8-plex</th> <th>Volume per 16-plex</th> <th>Volume per 24-plex</th> </tr> <tr> <td></td> <td></td> <td></td> <td></td> <td>With 15% overage</td> <td>With 15% overage</td> <td>With 15% overage</td> </tr> </thead> <tbody> <tr> <td></td> <td>Light Purple</td> <td>Nuclease buffer</td> <td>5 <math>\mu</math>L</td> <td>46 <math>\mu</math>L</td> <td>92 <math>\mu</math>L</td> <td>138 <math>\mu</math>L</td> </tr> <tr> <td></td> <td>Light Green</td> <td>Nuclease mix</td> <td>5 <math>\mu</math>L</td> <td>46 <math>\mu</math>L</td> <td>92 <math>\mu</math>L</td> <td>138 <math>\mu</math>L</td> </tr> <tr> <td></td> <td></td> <td><b>Total volume</b></td> <td><b>10 <math>\mu</math>L</b></td> <td><b>92 <math>\mu</math>L</b></td> <td><b>184 <math>\mu</math>L</b></td> <td><b>276 <math>\mu</math>L</b></td> </tr> </tbody> </table>	Reaction Mix 6 (RM6)							✓	Reagent	Component	Volume per sample	Volume per 8-plex	Volume per 16-plex	Volume per 24-plex					With 15% overage	With 15% overage	With 15% overage		Light Purple	Nuclease buffer	5 $\mu$ L	46 $\mu$ L	92 $\mu$ L	138 $\mu$ L		Light Green	Nuclease mix	5 $\mu$ L	46 $\mu$ L	92 $\mu$ L	138 $\mu$ L			<b>Total volume</b>	<b>10 <math>\mu</math>L</b>	<b>92 <math>\mu</math>L</b>	<b>184 <math>\mu</math>L</b>	<b>276 <math>\mu</math>L</b>
Reaction Mix 6 (RM6)																																												
✓	Reagent	Component	Volume per sample	Volume per 8-plex	Volume per 16-plex	Volume per 24-plex																																						
				With 15% overage	With 15% overage	With 15% overage																																						
	Light Purple	Nuclease buffer	5 $\mu$ L	46 $\mu$ L	92 $\mu$ L	138 $\mu$ L																																						
	Light Green	Nuclease mix	5 $\mu$ L	46 $\mu$ L	92 $\mu$ L	138 $\mu$ L																																						
		<b>Total volume</b>	<b>10 <math>\mu</math>L</b>	<b>92 <math>\mu</math>L</b>	<b>184 <math>\mu</math>L</b>	<b>276 <math>\mu</math>L</b>																																						
	7.1																																											
	7.2	Pipette-mix <b>RM6</b> up and down 10 times (do not vortex). Pipette-mix full volume.																																										
	7.3	Quick-spin <b>RM6</b> in a microcentrifuge to collect liquid.																																										
	7.4	Add <b>10 <math>\mu</math>L</b> of <b>RM6</b> to <b>40 <math>\mu</math>L</b> of each sample. The total volume should equal <b>50 <math>\mu</math>L</b> .																																										
	7.5	Pipette-mix each sample up and down 10 times. Pipette-mix full volume.																																										
	7.6	Quick-spin the tube strip in a microcentrifuge to collect liquid.																																										
		Run the <b>nuclease treatment</b> thermocycler program.																																										
	7.7	Set the lid temperature to 75°C. If the lid temperature is not programmable, it is acceptable to leave at 95–105°C.																																										
		<table border="1"> <thead> <tr> <th>Step</th> <th>Time</th> <th>Temperature</th> </tr> </thead> <tbody> <tr> <td>1</td> <td>60 min</td> <td>37°C</td> </tr> <tr> <td>2</td> <td>Hold</td> <td>4°C</td> </tr> </tbody> </table>	Step	Time	Temperature	1	60 min	37°C	2	Hold	4°C																																	
Step	Time	Temperature																																										
1	60 min	37°C																																										
2	Hold	4°C																																										
	7.8	Proceed to the next step of the protocol (pooling and cleanup with 1x SMRTbell cleanup beads).																																										



## 8. Sample pooling

Pool SMRTbell templates **in units of 8** samples after nuclease treatment of individual samples

✓	Step	Instructions for pooling, bead binding, washing, and sample elution
	8.1	In a 1.5 mL DNA LoBind tube combine nuclease treated libraries from step 7.7 in groups of 8 ( <b>8 x 50 µL</b> ) for a final sample volume of <b>400 µL</b> .
	8.2	Add <b>400 µL</b> of SMRTbell cleanup beads to each pooled nuclease treated sample.
	8.3	Pipette-mix up and down 10 times until the beads are evenly distributed. Pipette-mix full volume.
	8.4	Quick-spin the tube strip in a microcentrifuge to collect all liquid from the sides of the tubes.
	8.5	Leave at <b>room temperature</b> for <b>10 minutes</b> to allow DNA to bind beads.
	8.6	Place the tube strip in a magnetic separation rack until beads separate fully from the solution.
	8.7	Slowly pipette off the cleared supernatant without disturbing the beads. Discard the supernatant.
	8.8	Slowly dispense <b>1 mL</b> of <b>freshly prepared 80% ethanol</b> into each tube. After <b>30 seconds</b> , pipette off the 80% ethanol and discard.
	8.9	Repeat the previous step.
		Remove residual 80% ethanol:
	8.10	<ul style="list-style-type: none"> <li>• Remove the tube strip from the magnetic separation rack.</li> <li>• Quick-spin the tube strip in a microcentrifuge.</li> <li>• Place the tube strip back in a magnetic separation rack until beads separate fully from the solution.</li> <li>• Pipette off residual 80% ethanol and discard.</li> </ul>
	8.11	Remove the tube strip from the magnetic rack. <b>Immediately</b> add <b>100 µL</b> of <b>elution buffer</b> to each tube and resuspend the beads by pipetting up and down 10 times.
	8.13	Leave at <b>room temperature</b> for <b>5 minutes</b> to elute DNA.
	8.14	Place the tube strip in a magnetic separation rack until beads separate fully from the solution.
	8.15	Slowly pipette off the cleared supernatant without disturbing the beads. Transfer supernatant to a <b>1.5 mL DNA LoBind tube</b> . Discard old tube with beads.
Proceed to the next step of the protocol.		

## 9. SMRT boost bead wash of SMRTbell templates

The SMRT boost bead wash will prepare the library for sequencing.

Wash pooled SMRTbell templates with SMRT boost beads.

Perform one SMRT boost bead cleanup per final pool of up to 24 samples. For >24 samples you need to perform two SMRT boost bead cleanups with 24 samples each.

Bring SMRT boost buffer and SMRT boost beads to room-temperature prior to use.

Vortex room-temperature SMRT boost buffer prior to use. Vortex room-temperature SMRT boost beads prior to use and spin down to collect.

✓	Step	Instructions for bead washing and sample collection			
		Prepare the beads by washing in SMRT boost buffer			
		<ul style="list-style-type: none"> <li>To a 1.5 mL DNA LoBind tube, first add 1 mL of <b>room-temperature SMRT boost buffer</b> and then 100 <math>\mu</math>L of resuspended, <b>room-temperature SMRT boost beads</b>.</li> <li>Pulse-vortex 5 times to mix.</li> <li>Quick-spin the tube and magnetically separate for 2 minutes until the supernatant is clear.</li> </ul>			
9.1		<ul style="list-style-type: none"> <li>Remove the supernatant and wash once more by adding <b>1 mL of SMRT boost buffer</b>.</li> <li>Pulse-vortex 5 times to mix.</li> <li>Quick-spin the tube and magnetically separate for 2 minutes until the supernatant is clear. Remove the supernatant.</li> <li>Resuspend the beads in <b>200 <math>\mu</math>L SMRT boost buffer</b> by pulse-vortexing 5 times to mix.</li> <li>Quick-spin to collect the beads.</li> </ul>			
		Prepare the SMRTbell templates for SMRT boost bead washing by mixing together one or more 8-plex pools (each at 100 $\mu$ L) from Step 8. Bring the total volume to <b>300 <math>\mu</math>L with Elution Buffer</b> if combining less than 3 pools of 8, as shown in the table below.			
		<b># of Samples in Batch</b>	<b>8</b>	<b>16</b>	<b>24</b>
9.2		8-plex pool 1	100 $\mu$ L	100 $\mu$ L	100 $\mu$ L
		8-plex pool 2	0 $\mu$ L	100 $\mu$ L	100 $\mu$ L
		8-plex pool 3	0 $\mu$ L	0 $\mu$ L	100 $\mu$ L
		Elution Buffer	200 $\mu$ L	100 $\mu$ L	0 $\mu$ L
		<b>Total sample volume</b>	<b>300 <math>\mu</math>L</b>	<b>300 <math>\mu</math>L</b>	<b>300 <math>\mu</math>L</b>
9.3		Add <b>300 <math>\mu</math>L of SMRTbell templates</b> to <b>200 <math>\mu</math>L SMRT boost beads</b> in wash buffer from step 9.1. Pulse-vortex 5 times to mix. Spin down to collect.			
9.4		Gently <b>rotate-mix for 30 minutes at RT</b> using a rotator at low speed (~10 rpm).			
9.5		Spin down to collect and magnetically separate for 2 minutes until the supernatant is clear.			
9.6		Aliquot <b>500 <math>\mu</math>L</b> of the <b>SMRTbell-containing supernatant</b> into a fresh 1.5 mL LoBind tube.			

### Cleanup with 1X SMRTbell cleanup beads

- 9.7 Add **500  $\mu$ L** of resuspended, room-temperature SMRTbell cleanup beads to **500  $\mu$ L SMRTbell-containing supernatant** from step 9.6.
- 9.8 Pipette-mix the beads 8–10 times until evenly distributed.
- 9.9 Quick-spin the tube in a microcentrifuge to collect all liquid from the sides of the tubes.
- 9.10 Leave at **room temperature** for **10 minutes** to allow DNA to bind beads.
- 9.11 Place the tube in a magnetic separation rack until beads separate fully from the solution.
- 9.12 Slowly pipette off the cleared supernatant without disturbing the beads. Discard the supernatant.
- 9.13 Slowly dispense **1 mL** of **freshly prepared 80% ethanol** into each tube. After **30 seconds**, pipette off the 80% ethanol and discard.
- 9.14 Repeat the previous step.
- Remove residual 80% ethanol:
- Remove the tube from the magnetic separation rack.
- 9.15
- Quick-spin the tube in a microcentrifuge.
  - Place the tube back in a magnetic separation rack until beads separate fully from the solution.
  - Pipette off residual 80% ethanol and discard.
- 9.16 Remove the tube from the magnetic rack. **Immediately** add **100  $\mu$ L** of **elution buffer** to each tube and resuspend the beads.
- 9.17 Leave at **room temperature** for **5 minutes** to elute DNA.
- 9.18 Place the tube in a magnetic separation rack until beads separate fully from the solution.
- 9.19 Slowly pipette off the cleared supernatant without disturbing the beads. Transfer supernatant to a **1.5 mL DNA LoBind tube**. Discard old tube with beads.

Proceed to the next step of the protocol.

**SAFE STOPPING POINT - Store at -20°C**

## 10. Concentrating samples for ABC

Concentrate up to 24 SMRTbell templates into 15  $\mu\text{L}$  volume for ABC.

✓	Step	Instructions for bead binding, washing, and sample elution
	10.1	Add <b>100 <math>\mu\text{L}</math></b> SMRTbell cleanup beads to a 1.5mL DNA LoBind tube with <b>100 <math>\mu\text{L}</math> of SMRTbell templates</b> .
	10.2	Pipette-mix <b>slowly</b> up and down 10 times until the beads are evenly distributed.
	10.3	Quick-spin the tube in a microcentrifuge to collect all liquid from the sides of the tubes.
	10.4	Leave at <b>room temperature</b> for <b>10 minutes</b> to allow DNA to bind beads.
	10.5	Place the tube in a magnetic separation rack until the beads separate fully from the solution.
	10.6	Slowly pipette off the cleared supernatant without disturbing the beads. Discard the supernatant.
	10.7	Slowly dispense <b>250 <math>\mu\text{L}</math></b> , or enough to cover the beads, of <b>freshly prepared 80% ethanol</b> into each tube. After <b>30 seconds</b> , pipette off the 80% ethanol and discard.
	10.8	Repeat the previous step.
		Remove residual 80% ethanol:
	10.9	<ul style="list-style-type: none"> <li>Remove the tube from the magnetic separation rack.</li> <li>Quick-spin the tube in a microcentrifuge.</li> <li>Place the tube back in a magnetic separation rack until beads separate fully from the solution.</li> <li>Pipette off residual 80% ethanol and discard.</li> </ul>
	10.10	Remove the tube from the magnetic rack. <b>Immediately</b> add <b>16 <math>\mu\text{L}</math> of elution buffer</b> to each tube and resuspend the beads.
	10.11	Leave at <b>room temperature</b> for <b>5 minutes</b> to elute DNA.
	10.12	Place the tube in a magnetic separation rack until beads separate fully from the solution.
	10.13	Slowly pipette off the cleared eluate (supernatant) without disturbing the beads. Transfer supernatant to a <b>1.5 mL DNA LoBind tube</b> . Discard old tube with beads.
		QC step: Take <b>1 <math>\mu\text{L}</math></b> of eluted DNA for quantification and measure DNA concentration with a Qubit fluorometer using the <b>1X dsDNA HS kit</b> . <b>Do NOT use Qubit dsDNA BR assay kit as the concentration may be too low to measure. Do NOT dilute sample 1:10 as concentration may be too low to measure.</b>
	10.14	<p><b>Expect recovery of 0.5% or less (range 0.02% - 0.5%) relative to input starting mass.</b> For example, starting with 2 <math>\mu\text{g}</math> input per sample (or 16 <math>\mu\text{g}</math> per pool of 8 samples), the final mass recovered at this step is expected to be less than 80 ng total (range 3 ng–80 ng). Note, however, it is possible to observe recoveries outside of this range and still get good sequencing yield.</p> <p>Proceed to the next step of the protocol.</p> <p><b>SAFE STOPPING POINT - Store at -20°C</b></p>

## 11a. Annealing, binding & cleanup (ABC) for Revio System

Use the entire volume of 15  $\mu\text{L}$  pooled SMRTbell templates per ABC preparation and loading.

$\leq 24$  pooled SMRTbell templates (15  $\mu\text{L}$  total) go into Preparation A; The rest of the SMRTbell templates (also pooled in 15  $\mu\text{L}$  total) go into Preparation B.

✓	Step	Instructions																																																
		<p><b>Note:</b> Always use these values for each pooled preparation. Do not adjust based on measured concentration, value, or plex level.</p> <p><b>Annealing Sequencing primer</b> Combine the following components in a new low-binding tube and pipette to mix.</p> <table border="1"> <thead> <tr> <th></th> <th>Preparation A</th> <th>Preparation B</th> </tr> </thead> <tbody> <tr> <td># of Samples in Preparation</td> <td>8 or 16 or 24</td> <td>8 or 16 or 24</td> </tr> <tr> <td>Sample (SMRTbell templates)</td> <td>15 <math>\mu\text{L}</math></td> <td>15 <math>\mu\text{L}</math></td> </tr> <tr> <td>Annealing buffer</td> <td>7.5 <math>\mu\text{L}</math></td> <td>7.5 <math>\mu\text{L}</math></td> </tr> <tr> <td>Sequencing primer</td> <td>7.5 <math>\mu\text{L}</math></td> <td>7.5 <math>\mu\text{L}</math></td> </tr> <tr> <td><b>Total Volume</b></td> <td><b>30 <math>\mu\text{L}</math></b></td> <td><b>30 <math>\mu\text{L}</math></b></td> </tr> </tbody> </table> <p>Incubate at room temperature for 15 minutes then proceed to the next step.</p> <p><b>Binding sequencing polymerase</b></p> <p><b>Dilute sequencing polymerase</b> Combine the following components in a single low-bind tube and pipette to mix. The prepared volume of diluted polymerase is sufficient to process all specified samples in a batch.</p> <table border="1"> <thead> <tr> <th></th> <th>Preparation A</th> <th>Preparation B</th> </tr> </thead> <tbody> <tr> <td># of Samples in preparation</td> <td>8 or 16 or 24</td> <td>8 or 16 or 24</td> </tr> <tr> <td>Polymerase stock</td> <td>3.5 <math>\mu\text{L}</math></td> <td>3.5 <math>\mu\text{L}</math></td> </tr> <tr> <td>Polymerase buffer</td> <td>26.5 <math>\mu\text{L}</math></td> <td>26.5 <math>\mu\text{L}</math></td> </tr> <tr> <td><b>Total Volume</b></td> <td><b>30 <math>\mu\text{L}</math></b></td> <td><b>30 <math>\mu\text{L}</math></b></td> </tr> </tbody> </table> <p>Diluted polymerase must be used immediately.</p> <p><b>Bind sequencing polymerase</b> Add annealed sample to diluted polymerase and finger tap or pipette to mix.</p> <table border="1"> <thead> <tr> <th></th> <th>Preparation A</th> <th>Preparation B</th> </tr> </thead> <tbody> <tr> <td># of Samples in preparation</td> <td>8 or 16 or 24</td> <td>8 or 16 or 24</td> </tr> <tr> <td>Annealed sample</td> <td>30 <math>\mu\text{L}</math></td> <td>30 <math>\mu\text{L}</math></td> </tr> <tr> <td>Diluted Polymerase</td> <td>30 <math>\mu\text{L}</math></td> <td>30 <math>\mu\text{L}</math></td> </tr> <tr> <td><b>Total Volume</b></td> <td><b>60 <math>\mu\text{L}</math></b></td> <td><b>60 <math>\mu\text{L}</math></b></td> </tr> </tbody> </table>		Preparation A	Preparation B	# of Samples in Preparation	8 or 16 or 24	8 or 16 or 24	Sample (SMRTbell templates)	15 $\mu\text{L}$	15 $\mu\text{L}$	Annealing buffer	7.5 $\mu\text{L}$	7.5 $\mu\text{L}$	Sequencing primer	7.5 $\mu\text{L}$	7.5 $\mu\text{L}$	<b>Total Volume</b>	<b>30 <math>\mu\text{L}</math></b>	<b>30 <math>\mu\text{L}</math></b>		Preparation A	Preparation B	# of Samples in preparation	8 or 16 or 24	8 or 16 or 24	Polymerase stock	3.5 $\mu\text{L}$	3.5 $\mu\text{L}$	Polymerase buffer	26.5 $\mu\text{L}$	26.5 $\mu\text{L}$	<b>Total Volume</b>	<b>30 <math>\mu\text{L}</math></b>	<b>30 <math>\mu\text{L}</math></b>		Preparation A	Preparation B	# of Samples in preparation	8 or 16 or 24	8 or 16 or 24	Annealed sample	30 $\mu\text{L}$	30 $\mu\text{L}$	Diluted Polymerase	30 $\mu\text{L}$	30 $\mu\text{L}$	<b>Total Volume</b>	<b>60 <math>\mu\text{L}</math></b>	<b>60 <math>\mu\text{L}</math></b>
	Preparation A	Preparation B																																																
# of Samples in Preparation	8 or 16 or 24	8 or 16 or 24																																																
Sample (SMRTbell templates)	15 $\mu\text{L}$	15 $\mu\text{L}$																																																
Annealing buffer	7.5 $\mu\text{L}$	7.5 $\mu\text{L}$																																																
Sequencing primer	7.5 $\mu\text{L}$	7.5 $\mu\text{L}$																																																
<b>Total Volume</b>	<b>30 <math>\mu\text{L}</math></b>	<b>30 <math>\mu\text{L}</math></b>																																																
	Preparation A	Preparation B																																																
# of Samples in preparation	8 or 16 or 24	8 or 16 or 24																																																
Polymerase stock	3.5 $\mu\text{L}$	3.5 $\mu\text{L}$																																																
Polymerase buffer	26.5 $\mu\text{L}$	26.5 $\mu\text{L}$																																																
<b>Total Volume</b>	<b>30 <math>\mu\text{L}</math></b>	<b>30 <math>\mu\text{L}</math></b>																																																
	Preparation A	Preparation B																																																
# of Samples in preparation	8 or 16 or 24	8 or 16 or 24																																																
Annealed sample	30 $\mu\text{L}$	30 $\mu\text{L}$																																																
Diluted Polymerase	30 $\mu\text{L}$	30 $\mu\text{L}$																																																
<b>Total Volume</b>	<b>60 <math>\mu\text{L}</math></b>	<b>60 <math>\mu\text{L}</math></b>																																																
11a.1																																																		
11a.2																																																		
11a.3																																																		

Incubate at room temperature for 15 minutes. The bound complex can be stored at 4°C for 4 weeks.

### Purification of Polymerase Bound SMRTbell Complexes

1. Equilibrate the SMRTbell cleanup beads and the loading buffer to room temperature at least 30 min before use.
2. Add the following buffer volumes to each sample in each batch, as indicated:

	Preparation A	Preparation B
<b># of Samples in preparation</b>	8 or 16 or 24	8 or 16 or 24
<b>Binding Reaction</b>	60 µL	60 µL
<b>Dilution Buffer</b>	40 µL	40 µL
<b>Total Volume</b>	<b>100 µL</b>	<b>100 µL</b>

3. Add the indicated volume of SMRTbell cleanup beads to each sample in each preparation and gently pipette-mix. Incubate on the benchtop for 10 minutes.

	Preparation A	Preparation B
<b># of Samples in preparation</b>	8 or 16 or 24	8 or 16 or 24
<b>Binding Reaction</b>	100 µL	100 µL
<b>SMRTbell cleanup Beads</b>	120 µL	120 µL
<b>Total Volume</b>	<b>220 µL</b>	<b>220 µL</b>

11a.4

4. Place each tube in a magnetic bead rack until the beads collect to the side of the tube and the solution appears clear. Discard the supernatant. **DO NOT** wash the collected bead pellet with ethanol.
5. Immediately resuspend the beads in the indicated volumes of room temperature Loading buffer and pipette-mix:

	Preparation A	Preparation B
<b># of Samples in preparation</b>	8 or 16 or 24	8 or 16 or 24
<b>Loading Buffer</b>	49 µL	49 µL

6. To elute the polymerase-bound complexes, incubate the samples on the benchtop for at least 10 minutes at room temperature.
7. Place each tube in a magnetic bead rack until the beads collect to the side of the tube and the solution appears clear.
8. Transfer eluates to new low-binding tube. Place on ice **and protect from light**.

### Internal control dilution

**Prepare only one control dilution, regardless of number of samples (up to 48 samples).**

11a.5

1st Dilution. Mix well by flicking the tube by hand and pulse-spin to collect contents and keep on ice.

Reagent	Internal Control
Dilution buffer	19 $\mu$ L
Sequencing control	1.0 $\mu$ L
<b>Total volume</b>	<b>20 <math>\mu</math>L</b>

2<sup>nd</sup> Dilution. Mix well by flicking the tube by hand and pulse-spin to collect contents and keep on ice.

11a.6

Reagent	Internal Control
Dilution buffer	19 $\mu$ L
Sequencing control	1.0 $\mu$ L
<b>Total volume</b>	<b>20 <math>\mu</math>L</b>

3<sup>rd</sup> Dilution. Mix well by flicking the tube by hand and pulse-spin to collect contents and keep on ice.

11a.7

Reagent	Internal Control
Dilution buffer	19 $\mu$ L
Sequencing control	1.0 $\mu$ L
<b>Total Volume</b>	<b>20 <math>\mu</math>L</b>

### Final Loading Dilution

Combine the following and **protect from light**:

#### If loading $\leq 24$ samples

11a.8

# of Samples in preparation	$\leq 24$ -plex
Prepared Sample (Preparation A)	48.5 $\mu$ L
Loading buffer	48.5 $\mu$ L
Diluted internal control (Dilution 3)	3 $\mu$ L
<b>Total volume</b>	<b>100 <math>\mu</math>L</b>

Load 100  $\mu$ L of sample per well and/or store at 4C for up to 24 hours before use.

#### If loading $> 24$ samples

# of Samples in preparation	$> 24$ -plex
Prepared sample (Preparation A)	48.5 $\mu$ L
Prepared sample (Preparation B)	48.5 $\mu$ L
Loading buffer	0 $\mu$ L
Diluted internal control (Dilution 3)	3 $\mu$ L
<b>Total volume</b>	<b>100 <math>\mu</math>L</b>

Load 100  $\mu$ L of sample per well and/or store at 4°C for up to 24 hours before use.

## 11b. Annealing, binding, + cleanup (ABC) for Sequel II systems

Use entire volume of 15  $\mu\text{L}$  pooled SMRTbell templates per ABC batch and loading.

✓	Step	Instructions										
		<p><b>Note:</b> Always use these values for each pooled batch. Do not adjust based on measured concentration, value, or plex level of batch.</p> <p><b>Annealing Sequencing primer</b></p> <p>Combine the following components in a new low-binding tube and pipette to mix.</p>										
11b.1		<table border="1"> <thead> <tr> <th># of Samples in Batch</th> <th>8 or 16 or 24</th> </tr> </thead> <tbody> <tr> <td>Sample (SMRTbell templates)</td> <td>15 <math>\mu\text{L}</math></td> </tr> <tr> <td>Annealing Buffer</td> <td>7.5 <math>\mu\text{L}</math></td> </tr> <tr> <td>Sequel II Primer 3.2</td> <td>7.5 <math>\mu\text{L}</math></td> </tr> <tr> <td><b>Total volume</b></td> <td><b>30 <math>\mu\text{L}</math></b></td> </tr> </tbody> </table> <p>Incubate at room temperature for 15 minutes then proceed to the next step.</p> <p><b>Binding Sequencing polymerase</b></p> <p><b>Dilute Sequencing polymerase</b></p> <p>Combine the following components in a single low-bind tube and pipette to mix. The prepared volume of diluted polymerase is sufficient to process all specified samples in a batch.</p>	# of Samples in Batch	8 or 16 or 24	Sample (SMRTbell templates)	15 $\mu\text{L}$	Annealing Buffer	7.5 $\mu\text{L}$	Sequel II Primer 3.2	7.5 $\mu\text{L}$	<b>Total volume</b>	<b>30 <math>\mu\text{L}</math></b>
# of Samples in Batch	8 or 16 or 24											
Sample (SMRTbell templates)	15 $\mu\text{L}$											
Annealing Buffer	7.5 $\mu\text{L}$											
Sequel II Primer 3.2	7.5 $\mu\text{L}$											
<b>Total volume</b>	<b>30 <math>\mu\text{L}</math></b>											
11b.2		<table border="1"> <thead> <tr> <th># of Samples in Batch</th> <th>8 or 16 or 24</th> </tr> </thead> <tbody> <tr> <td>Polymerase Stock</td> <td>3.5 <math>\mu\text{L}</math></td> </tr> <tr> <td>Sequel II polymerase dilution buffer</td> <td>26.5 <math>\mu\text{L}</math></td> </tr> <tr> <td><b>Total Volume</b></td> <td><b>30 <math>\mu\text{L}</math></b></td> </tr> </tbody> </table> <p>Diluted polymerase must be used immediately.</p> <p><b>Bind Sequencing polymerase</b></p> <p>Add annealed sample to diluted polymerase and finger tap or pipette to mix.</p>	# of Samples in Batch	8 or 16 or 24	Polymerase Stock	3.5 $\mu\text{L}$	Sequel II polymerase dilution buffer	26.5 $\mu\text{L}$	<b>Total Volume</b>	<b>30 <math>\mu\text{L}</math></b>		
# of Samples in Batch	8 or 16 or 24											
Polymerase Stock	3.5 $\mu\text{L}$											
Sequel II polymerase dilution buffer	26.5 $\mu\text{L}$											
<b>Total Volume</b>	<b>30 <math>\mu\text{L}</math></b>											
11b.3		<table border="1"> <thead> <tr> <th># of Samples in Batch</th> <th>8 or 16 or 24</th> </tr> </thead> <tbody> <tr> <td>Annealed sample</td> <td>30 <math>\mu\text{L}</math></td> </tr> <tr> <td>Diluted polymerase</td> <td>30 <math>\mu\text{L}</math></td> </tr> <tr> <td><b>Total volume</b></td> <td><b>60 <math>\mu\text{L}</math></b></td> </tr> </tbody> </table> <p>Incubate at room temperature for 15 minutes. Bound complex can be stored at 4°C for 4 weeks.</p>	# of Samples in Batch	8 or 16 or 24	Annealed sample	30 $\mu\text{L}$	Diluted polymerase	30 $\mu\text{L}$	<b>Total volume</b>	<b>60 <math>\mu\text{L}</math></b>		
# of Samples in Batch	8 or 16 or 24											
Annealed sample	30 $\mu\text{L}$											
Diluted polymerase	30 $\mu\text{L}$											
<b>Total volume</b>	<b>60 <math>\mu\text{L}</math></b>											



### Purification of polymerase bound SMRTbell complexes

1. Equilibrate the cleanup beads and Sequel II loading buffer 3.2 to room temperature at least 30 min before use.
2. Add the following buffer volumes to each sample in each batch, as indicated:

# of Samples in Batch	8 or 16 or 24
Binding reaction	60 $\mu$ L
ABC buffer	40 $\mu$ L
<b>Total Volume</b>	<b>100 <math>\mu</math>L</b>

3. Add the indicated volume of cleanup beads to each sample in each batch and gently pipette-mix. Incubate on the benchtop for 10 minutes.

# of Samples in Batch	8 or 16 or 24
Binding reaction	100 $\mu$ L
Cleanup beads	120 $\mu$ L
<b>Total Volume</b>	<b>220 <math>\mu</math>L</b>

11b.4

4. Place each tube in a magnetic bead rack until the beads collect to the side of the tube and the solution appears clear. Discard the supernatant. **DO NOT** wash the collected bead pellet with ethanol.
5. Immediately resuspend the beads in the indicated volumes of room temperature Sequel II loading buffer 3.2 and pipette-mix:

# of Samples in Batch	8 or 16 or 24
Sequel II loading buffer 3.2	50 $\mu$ L

6. To elute the polymerase-bound complexes, incubate the samples on the benchtop for at least 10 minutes at room temperature.
7. Place each tube in a magnetic bead rack until the beads collect to the side of the tube and the solution appears clear.
8. Transfer eluates to new low-binding tube. Place on ice **and protect from light**.

#### Internal control dilution

1st Dilution. Mix well by flicking the tube by hand and pulse-spin to collect contents and keep on ice.

11b.5

Reagent	Internal Control
ABC buffer	19 $\mu$ L
Sequel II DNA internal control complex 3.2	1.0 $\mu$ L
<b>Total volume</b>	<b>20 <math>\mu</math>L</b>

- 11b.6 2nd Dilution. Mix well by flicking the tube by hand and pulse-spin to collect contents and keep on ice.

Reagent	Internal Control
ABC buffer	19 $\mu$ L
Sequel II DNA internal control complex 3.2	1.0 $\mu$ L
<b>Total volume</b>	<b>20 <math>\mu</math>L</b>

3rd Dilution. Mix well by flicking the tube by hand and pulse-spin to collect contents and keep on ice.

11b.7

Reagent	Internal Control
ABC buffer	19 $\mu$ L
Sequel II DNA internal control complex 3.2	1.0 $\mu$ L
<b>Total volume</b>	<b>20 <math>\mu</math>L</b>

### Final loading dilution

Combine the following and **protect from light**:

11b.8

# of Samples in batch	Sample 1 (24-plex or less)
Prepared sample	50 $\mu$ L
Sequel II loading buffer 3.2	67 $\mu$ L
Diluted internal control (dilution 3)	3 $\mu$ L
<b>Total volume</b>	<b>120 <math>\mu</math>L</b>

Load 115  $\mu$ L of sample per well and/or store at 4°C for up to 24 hours before use.

**PROTOCOL COMPLETE**

Revision history (description)	Version	Date
Initial release	01	March 2024

Research use only. Not for use in diagnostic procedures. © 2024 Pacific Biosciences of California, Inc. ("PacBio"). All rights reserved. Information in this document is subject to change without notice. PacBio assumes no responsibility for any errors or omissions in this document. Certain notices, terms, conditions and/or use restrictions may pertain to your use of PacBio products and/or third-party products. Refer to the applicable PacBio terms and conditions of sale and to the applicable license terms at [pacb.com/license](https://pacb.com/license). Pacific Biosciences, the PacBio logo, PacBio, Circulomics, Omniome, SMRT, SMRTbell, Iso-Seq, Sequel, Nanobind, SBB, Revio, Onso, Apton, Kinnex, and PureTarget are trademarks of PacBio.

© 2024 PacBio. All rights reserved. Research use only. Not for use in diagnostic procedures.

103-329-400 REV01 MAR2024

