


Introduction

There are many sequencing-based approaches to understanding complex metagenomic communities spanning targeted amplification to whole-sample shotgun sequencing. While targeted approaches provide valuable data at low sequencing depth, they are limited by primer design and PCR amplification. Whole-sample shotgun experiments generally use short-read, second-generation sequencing, which results in data processing difficulties. For example, reads less than 1 Kb in length will likely not cover a complete gene or region of interest, and will require assembly. This not only introduces the possibility of incorrectly combining sequence from different community members, it requires a high depth of coverage. As such, rare community members may not be represented in the resulting assembly.

Circular-consensus, single molecule, real-time (SMRT®) Sequencing reads in the 1-2 kb range, with >99% accuracy can be efficiently generated for low amounts of input DNA. 10 ng of input DNA sequenced in 4 SMRT Cells would generate >100,000 such reads. While throughput is low compared to second-generation sequencing, the reads are a true random sampling of the underlying community, since SMRT Sequencing has been shown to have no sequence-context bias. Long read lengths mean that that it would be reasonable to expect a high number of the reads to include gene fragments useful for analysis.

Profiling Using Circular Consensus Reads



BEI Mock Staggered Community






Figure 1. Comparison of CCS-gene prediction and other mapping strategies. (a) Distribution of taxonomic hits to the Silva Database⁴ using CCS gene predictions (blue) and 16S sequences (red). Filtered at 0.0005% abundance. Pie charts on nodes show total of all downstream assignments. * Indicates known members of the mock community. * Incorrect assignment likely due to gene sharing between *Salmonella*, *Shigella* and *E. coli*. (b) Mapping predicted genes gives comparable abundance numbers to mapping of raw sequence reads.

Stool Sample from a 4 Day Old Healthy Infant




Clostridium difficile Patient, Fecal Microbiota Transplant (FMT)

(a) Class level comparison with published results³ (not same sample)



(b)



(c)




Figure 3. CCS profile analysis of pre- and post-fecal transplant samples. (a) Class level assignments compare well with previously published results from similar samples characterized using microarrays and 16S techniques. (b) Full species-level classification shows very distinct community profiles pre- and post-transplant. (c) Profiling using SEED⁵ functional classifications, here pre-FMT secretion system genes are seen that are not present post-FMT.

Conclusions

- Metagenomic CCS profiling offers distinct advantages over both 16S and shotgun assembly methods.
- While having a high tolerance for sample input problems such as low input quantities and fragmented DNA, CCS profiling allows species level taxonomic classification and functional studies.
- Using a mock community, we demonstrate abundance measurements comparable to 16S quantification.
- Using primary samples we show comparable results to both 16S and microarray data, while allowing finer grain species level classification and meaningful functional insight.

References

- Hyatt D. *et al.*, Gene and translation initiation site prediction in metagenomic sequences. *Bioinformatics*. 2012 Sep 1;28(17):2223-2230.
- Huson D.H. *et al.*, Integrative analysis of environmental sequences using MEGAN 4. *Genome Res*. 2011. 21:1552-1560.
- Shankar V. *et al.* Species and genus level resolution analysis of gut microbiota in *Clostridium difficile* patients following fecal microbiota transplantation. *Microbiome*. 2014 Apr 21;2:13.
- Quast C. *et al.* (2013) The SILVA ribosomal RNA gene database project: improved data processing and web-based tools. *Nucl. Acids Res*. 41 (D1): D590-D596.
- Overbeek R. *et al.* The subsystems approach to genome annotation and its use in the project to annotate 1000 genomes. *Nucleic Acids Res*. 2005 Oct 7;33(17):5691-702.

