

Note: This Instrument Control Software release applies **only** to Revio systems.

New features

- Added support for SMRT® Link v13.1.
- Updated `lima` demultiplexing parameters to remove `--shared-prefix`. The updated parameters match those used in SMRT Link, and provide a minor improvement in precision.
- Updated adapter-calling to support SMRTbell® adapter index plates 96B, 96C, and 96D.

Known issues

- At times the touchscreen may freeze and be unresponsive during initialization. The workaround is to hard-reset the instrument to restart initialization by **either** pressing the power button, **or** by unplugging the instrument. Contact PacBio Technical Support if the issue persists.

Note: Contact your FSE (Field Service Engineer) or PacBio Technical Support to perform the Instrument Control Software installation.

Research use only. Not for use in diagnostic procedures. © 2024 Pacific Biosciences of California, Inc. ("PacBio"). All rights reserved. Information in this document is subject to change without notice. PacBio assumes no responsibility for any errors or omissions in this document. Certain notices, terms, conditions and/or use restrictions may pertain to your use of PacBio products and/or third-party products. Refer to the applicable PacBio terms and conditions of sale and to the applicable license terms at pacb.com/license.

Pacific Biosciences, the PacBio logo, PacBio, Circulomics, Omniome, SMRT, SMRTbell, Iso-Seq, Sequel, Nanobind, SBB, Revio, Onso, Apton, Kinnex and PureTarget are trademarks of PacBio.

© 2024 PacBio. All rights reserved. For research use only. Not for use in diagnostic procedures.

PN 103-407-800 Version 01 (April 2024)