

Supported instruments and chemistries

SMRT Link Cloud v25.3 supports the following:

Revio® system	Instrument software v13.3, all chemistries
Vega™ system	Instrument software v1.1.0, all chemistries

The **Sequel® II and IIe system** is not supported by SMRT Link Cloud. A local installation of SMRT Link v13.1 is required. Please see <https://www.pacb.com/support/software-downloads/> for details.

Supported data transfer options

- Direct transfers to network storage over an encrypted connection provided by SSH.
- Direct file transfer to cloud storage for Amazon S3, Google Cloud Storage, and Microsoft Azure Blob Storage, and S3-compatible API, including Oracle Cloud Object Storage and Cloudflare R2.

Notice: SMRT Link v25.3 will be the last release with the unencrypted rsync daemon transfer for Revio. Once disabled, file transfers using this option will fail. PacBio strongly recommends transitioning to a supported transfer option.

SMRT Link Cloud upgrades are performed automatically for all users during the maintenance window.

New Features and Updates

Overall

- API access is now enabled, and API keys can be easily generated from the SMRT Link user interface. See the SMRT Link Cloud user guide at <https://www.pacb.com/support/software-downloads/>.
- New PureTarget™ carrier application for use with the new PureTarget carrier panel 1.0.

Sample Setup

- The improved Loading Calculator now displays samples in fixed order and allows viewing of all samples without having to resize the table.
- The Loading Calculator now includes the lock feature for each sample, enabling “Import from Sample Setup” during Run Design. In Run Design, the sequencing plate must be selected before a sample can be imported.

Fixed issues

- **Run Design:** Fixed an issue where when using "Import from Sample Setup" the Application and Bio Sample name were not always populated.
- **Run Design:** Fixed an issue, where when importing Run Design CSV, Bio Sample names with unallowed characters such as spaces would create a run design with empty Bio Sample names instead of producing an error as expected.

Known issues

- **Sample Setup:** For Loading Calculator samples, Sample Setup CSVs cannot be imported or exported.
- **Run Design:** By design, Plate 1 and Plate 2 must be the same part number for Revio run designs. When using both the 4 rxn Revio SPRQ™ sequencing plate and 1 rxn Revio SPRQ sequencing plate, separate run designs must be created for each product.
- **Run Design:** For the Kinnex™ single-cell RNA application, when assigning Bio Sample names interactively, the populated default for "bcM0001–bcM0001" is "Bio Sample 1" and erroneously contains a space. Spaces are not an allowed character for Bio Sample names. The default Bio Sample names must be edited to exclude spaces to save the run design.
- **Run Design:** When interactively assigning Bio Sample names to barcodes, selecting multiple cells and pressing DEL will delete the Run Design. Only one cell should be deleted at a time.

Research use only. Not for use in diagnostic procedures. © 2025 Pacific Biosciences of California, Inc. ("PacBio"). All rights reserved. Information in this document is subject to change without notice. PacBio assumes no responsibility for any errors or omissions in this document. Certain notices, terms, conditions and/or use restrictions may pertain to your use of PacBio products and/or third-party products. Refer to the applicable PacBio terms and conditions of sale and to the applicable license terms at pacb.com/license.

Pacific Biosciences, the PacBio logo, PacBio, Circulomics, Omniome, SMRT, SMRTbell, Iso-Seq, Sequel, Nanobind, SBB, Revio, Onso, Apton, Kinnex, PureTarget, SPRQ and Vega are trademarks of PacBio.

See <https://github.com/broadinstitute/cromwell/blob/develop/LICENSE.txt> for Cromwell redistribution information.