
This system version supports Sequel II Systems **only**.

- | | |
|---------------------|--|
| New Features | Sequel® Instrument Control Software (ICS) <ul style="list-style-type: none">• Support for Sequel II Sequencing Chemistries 2.0 and 2.1. |
| Improvements | Sequel® Instrument Control Software (ICS) <ul style="list-style-type: none">• Numerous improvements to reliability, robotics workflow, and error handling.• Decreased convergence time for Active Laser Power.• Pre-extension time indication added to the instrument UI. |
| Fixed Issues | Sequel® Instrument Control Software (ICS) <ul style="list-style-type: none">• After system startup or restart, the system software version menu now displays the correct PA version. |
| Known Issues | <ul style="list-style-type: none">• When performing a run using two Sequencing Plates, each plate must have a separate Sequel SMRT Cell Oil tube. |

For Research Use Only. Not for use in diagnostic procedures. © Copyright 2019, Pacific Biosciences of California, Inc. All rights reserved. Information in this document is subject to change without notice. Pacific Biosciences assumes no responsibility for any errors or omissions in this document. Certain notices, terms, conditions and/or use restrictions may pertain to your use of Pacific Biosciences products and/or third party products. Please refer to the applicable Pacific Biosciences Terms and Conditions of Sale and the applicable license terms at <http://www.pacb.com/legal-and-trademarks/product-license-and-use-restrictions/>.

Pacific Biosciences, the Pacific Biosciences logo, PacBio, SMRT, SMRTbell, Iso-Seq and Sequel are trademarks of Pacific Biosciences. BluePippin and SageELF are trademarks of Sage Science, Inc. NGS-go and NGSengine are trademarks of GenDx. FEMTO Pulse and Fragment Analyzer are trademarks of Agilent Technologies Inc. All other trademarks are the sole property of their respective owners.

P/N 101-774-400 Version 3 (October 2019)