

---

This system version supports Sequel Systems **only**.

- New Features**    **Sequel® Instrument Control Software (ICS)**
- Automatic detection and recovery of inserts with missing adapters.
- Improvements**    **Sequel® Instrument Control Software (ICS)**
- Numerous improvements to reliability, robotics workflow, and error handling.
  - Improved support for the Auto Analysis feature in SMRT® Link v8.0.
  - Decreased convergence time for Active Laser Power.
  - Pre-extension time indication added to the instrument UI.
- Fixed Issues**    **Sequel® Instrument Control Software (ICS)**
- The on-instrument expiration time for Sequencing Plate 3.0, 8-reaction version is no longer incorrectly set to 72 hours. The correct expiration time is 100 hours.
  - The virtual keyboard no longer occasionally covers input fields and selections when using the Admin dialog or when using the Quick Run functionality.
  - The instrument self-test no longer becomes intermittently unresponsive.
- Known Issues**    **Sequel® Instrument Control Software (ICS)**
- When performing a run using two Sequencing Plates, each plate must have a separate Sequel SMRT Cell Oil tube.

For Research Use Only. Not for use in diagnostic procedures. © Copyright 2015 - 2020, Pacific Biosciences of California, Inc. All rights reserved. Information in this document is subject to change without notice. Pacific Biosciences assumes no responsibility for any errors or omissions in this document. Certain notices, terms, conditions and/or use restrictions may pertain to your use of Pacific Biosciences products and/or third party products. Please refer to the applicable Pacific Biosciences Terms and Conditions of Sale and the applicable license terms at <http://www.pacb.com/legal-and-trademarks/product-license-and-use-restrictions/>.

Pacific Biosciences, the Pacific Biosciences logo, PacBio, SMRT, SMRTbell, Iso-Seq and Sequel are trademarks of Pacific Biosciences. BluePippin and SageELF are trademarks of Sage Science, Inc. NGS-go and NGSengine are trademarks of GenDx. FEMTO Pulse and Fragment Analyzer are trademarks of Agilent Technologies Inc. All other trademarks are the sole property of their respective owners.

P/N 100-715-100 Version 15 (January 2020)