Improved liquid biopsy assay performance using sequencing by binding (SBB)

Daniel Nasko, Phillip Pham, Stuti Joshi, Kristi Kim, Nairi Pezeshkian, Christopher Kingsley, Daniel Portik, Young Kim, Alex Sockell, and Jonas Korlach PacBio, 1305 O'Brien Drive, Menlo Park, CA 94025

Introduction	Decreased noise at variant sites	More accurate VAF estimation
Liquid biopsy is revolutionizing the field of early cancer detection research through non-invasive detection of tumor DNA in the blood. However, existing liquid biopsy assays are limited in their sensitivity for ctDNA detection at low variant allele frequencies (VAFs). Here we describe the application of the PacBio Onso short-read sequencing system to enable detection of	Increased sequencing accuracy achieved by SBB reduces noise due to errors at variant sites. In SBS data, all 3 alternative alleles have a similar number of base calls, obscuring the true variant. In SBB data, the lack of errors reveals the variant is a $T \rightarrow C$.	SBB also performs better at estimating the true VAF than SBS with only half the sequencing depth, with less variability in allele frequency estimates. $\frac{0.10\% \text{ VAF}}{\text{True positive VAF distributions}}$

ctDNA at low VAFs using the SeraCare Complete ctDNA Mutation Mix reference sample.

PacBie









Figure 5. Distribution of observed VAFs for true positive variants in the 0.1% VAF sample for SBB at 12,000X and SBS at 24,000X.

Improved recall at low VAFs

Improved accuracy with SBB results in improved recall for low frequency VAFs (0.05%, 0.1%) compared to SBS, even

Figure 1. (A) SBB results in improved empirical Q scores compared to SBS due to decoupling of the interrogation and incorporation steps during each sequencing cycle. **(B)** Empirical Q score for top 99% of bases with SBB vs. SBS.

Experimental design



Figure 3. IGV plots (left) and base counts (right) for SBS **(A)** and SBB **(B)**.

Reduction in false positives

Reduced false positive variant calls in the WT sample for SBB lowers the limit of detection and improves recall when specificity is set at 100%.



with half the sequencing depth.



Figure 6. Sensitivity for variant detection at VAFs 0.0% (WT), 0.05%, 0.1%, and 0.25% for SBB at 12,000X and SBS at 24,000X.

Conclusions

We observe superior sensitivity for ctDNA detection using SBB compared to SBS at low VAFs (0.05%, 0.10%) at comparable sequencing depth. Furthermore, SBB requires significantly less sequencing to achieve comparable sensitivity results to SBS. Taken together, our results demonstrate the potential of SBB to improve upon existing methods of liquid biopsy and better enable research on early cancer detection.

Figure 2. Overview of experimental workflow. Post-capture libraries were split into 2 aliquots for the comparison, one of which was converted using the Onso library conversion kit to add Onso-compatible adapter sequences.

Figure 4. Observed VAFs for WT sample (false positives, FP) and 0.1% VAF sample (true positives, TP). Dashed box indicates variants that would be excluded with this specificity cutoff for each technology.

Research use only. Not for use in diagnostic procedures. © 2023 Pacific Biosciences of California, Inc. ("PacBio"). All rights reserved. Information in this document is subject to change without notice. PacBio assumes no responsibility for any errors or omissions in this document. Certain notices, terms, conditions and/or use restrictions may pertain to your use of PacBio products and/or third-party products. Refer to the applicable PacBio terms and conditions of sale and to the applicable license terms at pacb.com/license. Pacific Biosciences, the PacBio logo, PacBio, Circulomics, Omniome, SMRT, SMRTbell, Iso-Seq, Sequel, Nanobind, SBB, Revio, and Onso are trademarks of PacBio.