

The background of the slide is a blurred image of a laboratory setting. In the foreground, a multi-well microplate is visible, with several wells containing a bright pink liquid. A pipette tip is positioned above one of the wells, with a single drop of the pink liquid about to fall. The overall scene is brightly lit, with a soft focus on the background, emphasizing the laboratory equipment and the color of the liquid.

PacBio

Technical overview – HiFi library preparation using HiFi prep kits for high-throughput sequencing on PacBio long-read systems

Sequel II and IIe systems ICS v11.0

Revio system ICS v13.1

SMRT Link v13.1


PN 103-424-600 Rev 01 | April 2024

HiFi library preparation using HiFi prep kits for high-throughput sequencing on PacBio long-read systems

Technical Overview

1. HiFi library prep products and workflow options for WGS applications
2. HiFi prep kit 96 workflow overview for WGS applications
3. HiFi prep kit 96 example sequencing performance data
4. HiFi plex prep kit 96 workflow overview for WGS and amplicon sequencing applications
5. HiFi plex prep kit 96 example sequencing performance data
6. Technical documentation & applications support resources

HiFi library preparation using HiFi prep kits for high-throughput sequencing on PacBio long-read systems: Getting started



PacBio | HI-FI PREP + PLEX KITS

Scale up + stay ahead with HiFi prep kits

PacBio HiFi prep kits offer scalable library preparation solutions for projects and genomes of all sizes. Maximize your sequencing on the Revio™ system by preparing up to 96 libraries at a time with flexible, automated workflows including DNA shearing and size selection. Now is the time to do more with long reads using the new family of HiFi prep kits.

HiFi prep kits

HiFi prep kit 96: 24–96 SMRTbell™ libraries using only 2 µg of DNA per sample. Ideal for human, plant, and animal WGS applications.

HiFi plex prep kit 96: 24–96 multiplexed SMRTbell libraries from 20–200 ng of gDNA per sample or 20–250 ng per amplicon. Ideal for microbial, metagenomic, targeted, or low coverage WGS applications.

Workflow steps

SRE HT	HiFi prep kit 96	HiFi plex prep kit 96
	Automated shearing	Automated shearing
	HiFi prep kit 96	HiFi plex prep kit 96
	Pooling	Pooling
Nucleic acid treatment	HiFi prep kit 96	HiFi plex prep kit 96
	HiFi prep kit 96	HiFi plex prep kit 96
Sequencing	HiFi prep kit 96	HiFi plex prep kit 96
	HiFi prep kit 96	HiFi plex prep kit 96

Sequencing performance on the Revio system

HiFi reads	HiFi prep kit 96	HiFi plex prep kit 96
HiFi reads yield	101.22 Gb	87.64 Gb
HiFi read length (NSO)	18,016 bp	14,033 bp
HiFi read quality (median)	Q34	Q40
Base quality ≥Q30	92.78%	94.57%

Learn more: pacb.com/wgs

PacBio

Brochure – HiFi prep kits (102-326-608)

Summary overview of application-specific DNA sample extraction, library preparation and/or data analysis workflow recommendations.

PacBio

Generating PureTarget™ repeat expansion panel libraries

Procedure & checklist

Before you begin

This procedure describes the workflow for generating PureTarget repeat expansion panel libraries.

Overview	
Samples	8–24 (processed in batches of 8)
Library prep time	8 hours + 1–2 hours for up to 24 samples
Annealing binding clean up time	1 hour +/- 10 minutes

DNA input	
Quantity	2 µg in Buffer EB, TE buffer (pH 8), or nuclease-free water
DNA size distribution	50% >30 kb

Sample multiplexing	
Sequel® II systems	Up to 24 samples
Revio™ system	Up to 48 samples

© 2024 PacBio. All rights reserved. Research use only. Not for use in diagnostic procedures. 103-329-490 REV01 MAR2024

PacBio

Procedure & checklist – Preparing whole genome libraries using the HiFi prep kit 96 (103-420-700)

Technical documentation containing application-specific library preparation protocol details.

PacBio

Technical overview – HiFi library preparation using HiFi prep kits for high-throughput sequencing on PacBio long-read systems

Sequel II and IIe systems ICS v11.0
Revio system ICS v13.1
SMRT Link v13.1

PN 103-424-600 Rev 01 | April 2024

HiFi sequencing performance of HiFi prep kit 96 libraries


Example HiFi read length and yield performance obtained for HiFi prep kit 96 libraries constructed from gDNA isolated from human and animal whole blood samples using Nanobind PanDNA kit

© 2024 PacBio. All rights reserved. Research use only. Not for use in diagnostic procedures. 103-329-490 REV01 MAR2024

PacBio

Technical Overview – Technical overview – HiFi library preparation using HiFi prep kits for high-throughput sequencing on PacBio long-read systems (103-424-600)

Technical overview presentations describe workflow details for constructing PacBio SMRTbell libraries for specific applications. Example sequencing performance data for a given application are also summarized.





HiFi library prep products and workflow options for WGS applications

HiFi library prep workflow recommendations for low-throughput whole genome sequencing applications



≤24 samples per month

Low-throughput HiFi WGS workflow

DNA sample extraction




DNA sample QC




DNA size selection




DNA shearing




DNA repair & A-tailing



Adapter ligation



Nuclease treatment



Library cleanup or size selection



Library construction

Library construction

Protocol documentation or reference



Guide & overview – Nanobind PanDNA kit ([103-394-800](#))



Refer to third-party user guide documentation



Procedure & checklist – Removing short DNA fragments with the Short Read Eliminator (SRE) kit ([102-982-300](#)) [For manual workflows]



Guide and overview – Short Read Eliminator (SRE), DNA shearing, and cleanup for the Hamilton Microlab Prep system ([103-424-100](#))

SPK 3.0 protocol reference for manual & automated workflows



Procedure & checklist – Preparing whole genome and metagenome sequencing libraries using SMRTbell prep kit 3.0 ([102-166-600](#))

Recommended equipment & consumables



Nanobind PanDNA kit (103-260-000 / 24 RXN)



Qubit 4 fluorometer system
Femto Pulse system



Optional: Short read eliminator (SRE) kit (102-208-300 / 24 RXN)



Hamilton MicroLab Prep (Pipette-based shearing)



SMRTbell prep kit 3.0 (102-182-700 / 24 RXN)



Animal/plant/human WGS:¹
3.1X (35%) AMPure PB SS (>4 kb)

Microbial/metagenomic WGS:
1X SMRTbell bead cleanup



¹ Can optionally perform gel cassette-based size selection to enrich for SMRTbell library inserts >10 kb.


HiFi library prep workflow recommendations for high-throughput whole genome sequencing applications



>24 samples per month

High-throughput HiFi WGS workflow

DNA sample extraction




DNA sample QC




DNA size selection




DNA shearing




DNA repair & A-tailing



Adapter ligation



Nuclease treatment



Library cleanup or size selection



Library construction

Library construction

Protocol documentation or reference



Guide & overview – Nanobind HT kits ([103-028-100](#))¹



Refer to third-party user guide documentation

HiFi prep kit 96 protocol references for manual & automated workflows



Procedure & checklist – Preparing whole genome libraries using the HiFi prep kit 96 ([103-420-700](#))



Guide & overview – Automated HiFi prep 96 and HiFi annealing, binding, and cleanup for the Hamilton NGS Microlab STAR system ([103-425-700](#))

HiFi plex prep kit 96 protocol references for manual & automated workflows



Procedure & checklist – Preparing multiplexed whole genome and amplicon libraries using the HiFi plex prep kit 96 ([103-418-800](#))



Guide & overview – Automated HiFi plex prep 96 for the Hamilton NGS Microlab STAR system ([103-425-800](#))

Recommended equipment & consumables



Nanobind HT CBB kit / Nanobind HT 1 mL whole blood kit



Microplate reader + Quant-iT dsDNA HS assay; Femto Pulse system



SRE HT kit (103-124-500 / 96 RXN) using Hamilton automation system



Pipette-based shearing using Hamilton automation system

HiFi prep kit 96 (103-381-200 / 96 RXN); **or**



Revio HiFi prep kit 96 (103-382-200 / 96 RXN) [Includes AMPure PB + Revio polymerase kit 96]; **or**

HiFi plex prep kit 96 (103-381-300 / 96 RXN) + SMRTbell adapter index plate 96A/B/C/D



Animal/plant/human WGS:
3.1X (35%) AMPure PB SS (>4 kb)








Microbial/metagenomic WGS:
1X SMRTbell bead cleanup

¹ Refer to PacBio [Documentation](#) website for specific Nanobind HT Procedures & checklists.


Overview of HiFi prep kit product bundles and components

PacBio HiFi prep kits offer scalable library prep solutions for projects and genomes of all sizes

Product	Components	Supported workflow steps	Description
SRE HT kit [96 rxn] 103-124-500 	<ul style="list-style-type: none"> • Buffer SRE • Buffer LTE 	<ul style="list-style-type: none"> <input type="checkbox"/> Genomic DNA size selection 	<ul style="list-style-type: none"> • 96 short read eliminator (SRE) reactions for performing rapid high-pass size selection of genomic DNA samples to remove DNA fragments <10 kb (with progressive deletion of fragments <25 kb)
HiFi prep kit 96 bundle [96 rxn] 103-381-200 	<ul style="list-style-type: none"> • HiFi prep kit 96 • SMRTbell cleanup beads-85mL • SMRTbell adapter index plate 96A • Elution buffer (50 mL) • Buffer LTE HT (50 mL) 	<ul style="list-style-type: none"> <input type="checkbox"/> Pipette-based DNA shearing <input type="checkbox"/> SMRTbell library construction 	<ul style="list-style-type: none"> • 24–96 SMRTbell libraries using only 2 µg of gDNA per sample • Ideal for human, plant, and animal WGS applications
Revio HiFi prep kit 96 bundle [96 rxn] 103-382-200 	<ul style="list-style-type: none"> • HiFi prep kit 96 • SMRTbell cleanup beads-85mL • SMRTbell adapter index plate 96A • Elution buffer (50 mL) • SRE HT kit (96 RXN) • AMPure PB (5 mL) • Revio polymerase kit 96¹ 	<ul style="list-style-type: none"> <input type="checkbox"/> gDNA size selection (SRE) <input type="checkbox"/> Pipette-based DNA shearing <input type="checkbox"/> SMRTbell library construction <input type="checkbox"/> SMRTbell library size selection (AMPure PB) <input type="checkbox"/> Sequencing preparation (ABC – Anneal sequencing primer / Bind polymerase / Complex cleanup) 	<ul style="list-style-type: none"> • 24–96 SMRTbell libraries using only 2 µg of gDNA per sample • Includes Revio polymerase kit 96 for performing sequencing primer annealing, polymerase binding and complex cleanup reactions to prepare samples for sequencing on a Revio system • Ideal for human, plant, and animal WGS applications
HiFi plex prep kit 96 bundle [96 rxn] 103-381-300 	<ul style="list-style-type: none"> • HiFi plex prep kit 96 • SMRTbell cleanup beads-52mL • Buffer LTE HT (50 mL) • Elution buffer (50 mL) 	<ul style="list-style-type: none"> <input type="checkbox"/> Pipette-based DNA shearing <input type="checkbox"/> Multiplexed SMRTbell library construction 	<ul style="list-style-type: none"> • 24–96 multiplexed SMRTbell libraries from 50–300 ng of gDNA per sample or 20–250 ng per amplicon sample • Ideal for microbial, metagenomic, targeted, or low coverage WGS (2-3X coverage per sample) applications
SMRTbell adapter index plate 96A/B/C/D [96 rxn] 102-009-200 (A) / 102-547-800 (B) / 102-547-900 (C) / 102-548-000 (D) 	<ul style="list-style-type: none"> • SMRTbell adapter index plate 96A/B/C/D 	<ul style="list-style-type: none"> <input type="checkbox"/> Indexed adapter ligation (during multiplexed SMRTbell library construction) 	<ul style="list-style-type: none"> • SMRTbell adapter index plate contains 96 indexed adapters in plate format (one sample per indexed adapter) • Plate 96A: bc2001–bc2096 / Plate 96B: bc2097–bc2192 / Plate 96C: bc2193–bc2288 / Plate 96D: bc2289–bc2384

Overview of HiFi prep kit bundle supported workflow steps



Revio HiFi prep kit 96 bundle supports library preparation and sequencing preparation (ABC) workflows



Sequencing preparation workflow overview for HiFi prep kit libraries

For HiFi prep kit 96 WGS libraries bound with Revio polymerase kit or Sequel II binding kit 3.2, follow SMRT Link Sample Setup calculator instructions for ABC and final loading dilution procedure

HiFi prep kit bundle	Polymerase kit ¹	Sample setup workflow & procedural reference	
HiFi prep kit 96 bundle 103-381-200	Revio polymerase kit (102-817-600) <i>or</i> Sequel II binding kit 3.2 (102-333-300)	Anneal sequencing primer, bind polymerase, complex cleanup (ABC)	<input type="checkbox"/> Follow SMRT Link v13.1+ Sample Setup
		Final loading dilution procedure	<input type="checkbox"/> Follow SMRT Link v13.1+ Sample Setup




¹ HiFi prep kit 96 bundle (103-381-200) does not include any polymerase kit component; Revio polymerase kit (102-817-600) or Sequel II binding kit 3.2 (102-333-300) must be purchased separately.

Sequencing preparation workflow overview for HiFi prep kit libraries (cont.)

For HiFi prep kit 96 WGS libraries bound with Revio polymerase kit 96¹, follow SMRT Link Revio polymerase kit 96 Loading Calculator instructions for final loading dilution procedure


HiFi prep kit bundle	Polymerase kit ¹	Sample setup workflow & procedural reference	
Revio HiFi prep kit 96 bundle 103-382-200	Revio polymerase kit 96 (103-253-600)	Anneal sequencing primer, bind polymerase, complex cleanup (ABC)	<input type="checkbox"/> Follow Procedure & checklist or Guide & overview
		Sample pooling (optional)	<input type="checkbox"/> Follow SMRT Link v13.1+ Pooling Calculator
		Final loading dilution procedure	<input type="checkbox"/> Follow SMRT Link v13.1+ Loading Calculator



Procedure & checklist – Preparing whole genome libraries using the HiFi prep kit 96 (103-420-700)




Guide & overview – Automated HiFi prep 96 and HiFi annealing binding and cleanup for the Hamilton NGS Microlab STAR system (103-425-700)



Sequencing preparation workflow overview for HiFi plex prep kit libraries

For HiFi plex prep kit 96 WGS or amplicon libraries bound with Revio polymerase kit or Sequel II binding kit 3.1/3.2, follow SMRT Link Sample Setup calculator instructions for ABC and final loading dilution procedure






HiFi prep kit bundle	Polymerase kit ²	Sample setup workflow & procedural reference	
HiFi plex prep kit 96 bundle ¹ 103-381-300	Revio polymerase kit (102-817-600) <i>or</i> Sequel II binding kit 3.1 (102-333-400) <i>or</i> Sequel II binding kit 3.2 (102-333-300)	Anneal sequencing primer, bind polymerase, complex cleanup (ABC)	<input type="checkbox"/> Follow SMRT Link v13.1+ Sample Setup
		Final loading dilution procedure	<input type="checkbox"/> Follow SMRT Link v13.1+ Sample Setup




The screenshot shows the PacBio Sample Setup web interface. The top navigation bar includes the PacBio logo, 'Sample Setup', and links for Notifications, Settings, Help, and the user's name (smark, Lab Tech). The main content area features a 'Sample Setup' header with a gear icon and a list of actions: '+ Add Calculation', 'Edit', 'Import', 'Export', and 'Delete'. Below this is a table with columns for 'Sample name', 'Date created', and 'Comment'. A dropdown menu is open, showing three options: 'Sequel II binding kit 2.1/2.2', 'Sequel II binding kit 3.1/3.2, Revio polymerase kit' (which is checked), and 'Revio polymerase kit 96'. The table shows one row with the sample name 'My Batch of Samples', a date of '2024-03-25, 05:15:58 PM', and a comment 'smark Revio polymerase kit'. A search bar and a note 'Displaying rows 1 to 13 out of 112' are also visible.

HiFi prep kit 96 bundle (103-381-200) configuration details

Contains HiFi library preparation reagents for processing up to 96 samples

Product	Components	Description
HiFi prep kit 96 bundle (103-381-200)	 HiFi prep kit 96 (103-122-600)	<ul style="list-style-type: none"> HiFi library preparation reagents for 96 samples. SMRTbell cleanup beads and low TE buffer included. Note: Size-selection reagents are sold separately.¹
	 SMRTbell cleanup beads-85mL (103-294-600)	<ul style="list-style-type: none"> Paramagnetic beads supplied at 85 mL to selectively bind DNA fragments. The beads are used for: 1) cleanup – removing excess primers, nucleotides, salts, and enzymes; and 2) buffer exchanges.
	 SMRTbell adapter index plate 96A (102-009-200)	<ul style="list-style-type: none"> SMRTbell adapter index plate 96A contains 96 indexed adapters in plate format (one sample per indexed adapter). It includes indexes bc2001–bc2096.
	 Elution buffer (50 mL) (101-633-500)	<ul style="list-style-type: none"> Elution buffer supplied at 50 mL is used with the SMRTbell library template preparation and barcoding associated kits for workflows requiring additional elution buffer and AMPure PB beads cleanup steps.
	 Buffer LTE HT (50 mL) (103-306-100)	<ul style="list-style-type: none"> Low TE buffer supplied at 50 mL to support 96 reactions. This buffer consists of 10 mM Tris-HCl (pH 9.0) and 0.1 mM EDTA.

Supported workflow steps










HiFi prep kit 96 bundle
(103-381-200)

¹ For high-throughput animal/plant/human WGS applications, we recommend performing size selection on input genomic DNA using SRE HT kit (103-124-500) to remove DNA fragments <10 kb and on HiFi SMRTbell libraries using AMPure PB beads (100-265-900) to remove residual library insert fragments <5 kb.

² HiFi prep kit 96 bundle includes Buffer LTE for diluting samples for DNA shearing step. For high-throughput animal/plant/human WGS applications, we recommend performing DNA shearing using a pipette-based shearing method with a third-party liquid handling system.

Revio HiFi prep kit 96 bundle (103-382-200) configuration details

Contains HiFi library preparation reagents for processing up to 96 samples and includes Revio polymerase kit

Product	Components	Description
Revio HiFi prep kit 96 bundle (103-382-200)	 HiFi prep kit 96 (103-122-600)	<ul style="list-style-type: none"> HiFi library preparation reagents for 96 samples. SMRTbell cleanup beads and low TE buffer included.
	 SMRTbell cleanup beads-85mL (103-294-600)	<ul style="list-style-type: none"> Paramagnetic beads supplied at 85 mL to selectively bind DNA fragments. The beads are used for: 1) cleanup – removing excess primers, nucleotides, salts, and enzymes; and 2) buffer exchanges.
	 SMRTbell adapter index plate 96A (102-009-200)	<ul style="list-style-type: none"> SMRTbell adapter index plate 96A contains 96 indexed adapters in plate format (one sample per indexed adapter). It includes indexes bc2001–bc2096.
	 Elution buffer (50 mL) (101-633-500)	<ul style="list-style-type: none"> Elution buffer supplied at 50 mL is used with the SMRTbell library template preparation and barcoding associated kits for workflows requiring additional elution buffer and AMPure PB beads cleanup steps.
	 SRE HT kit (103-124-500)	<ul style="list-style-type: none"> 96 short read eliminator (SRE) reactions for performing rapid high-pass size selection¹ of genomic DNA samples to remove DNA fragments <10 kb (with progressive deletion of fragments <25 kb).
	 AMPure PB (5 mL) (100-265-900)	<ul style="list-style-type: none"> AMPure PB contains specially formulated paramagnetic beads for the SMRT sequencing workflow.¹
	 Revio polymerase kit 96 (103-253-600)	<ul style="list-style-type: none"> Reagents for binding sequencing polymerase to SMRTbell library. Includes 96 reactions, each supporting one Revio SMRT Cell.

Supported workflow steps







Revio HiFi prep kit 96 bundle
(103-382-200)

¹ For high-throughput animal/plant/human WGS applications, we recommend performing size selection on input genomic DNA using SRE HT kit (103-124-500) to remove DNA fragments <10 kb and on HiFi SMRTbell libraries using AMPure PB beads (100-265-900) to remove residual library insert fragments <3 kb.

² Revio HiFi prep kit 96 bundle includes Buffer LTE for diluting samples for DNA shearing step. For high-throughput animal/plant/human WGS applications, we recommend performing DNA shearing using a pipette-based shearing method with a third-party liquid handling system.

HiFi plex prep kit 96 bundle (103-381-300) configuration details

Contains HiFi library preparation reagents for multiplexing 96 samples

Product	Components	Description
HiFi plex prep kit 96 bundle (103-381-300)	 HiFi plex prep kit 96 (103-122-800)	<ul style="list-style-type: none"> HiFi library preparation reagents for multiplexing 96 samples. SMRTbell® cleanup beads are included. Note: SMRTbell adapter indexes and SMRTbell library size-selection reagents are sold separately.
	 SMRTbell cleanup beads-52mL (103-294-500)	<ul style="list-style-type: none"> Paramagnetic beads supplied at 52 mL to selectively bind DNA fragments. The beads are used for: 1) cleanup – removing excess primers, nucleotides, salts, and enzymes; and 2) buffer exchanges.
	 Elution buffer (50 mL) (101-633-500)	<ul style="list-style-type: none"> Elution buffer supplied at 50 mL is used with the SMRTbell library template preparation and barcoding associated kits for workflows requiring additional elution buffer and AMPure PB beads cleanup steps.
	 Buffer LTE HT (50 mL) 103-306-100	<ul style="list-style-type: none"> Low TE buffer supplied at 50 mL to support 96 reactions. This buffer consists of 10 mM Tris-HCl (pH 9.0) and 0.1 mM EDTA.


Supported workflow steps



Note: SMRTbell adapter index plate must be purchased separately

Choose one or more of the following products to perform sample indexing:

- SMRTbell adapter index plate 96A (102-009-200) – Contains indexes bc2001–bc2096
- SMRTbell adapter index plate 96B (102-547-800) – Contains indexes bc2097–bc2192
- SMRTbell adapter index plate 96C (102-547-900) – Contains indexes bc2193–bc2288
- SMRTbell adapter index plate 96D (102-548-000) – Contains indexes bc2289–bc2384





HiFi prep kit 96 workflow overview for WGS applications

Procedure & checklist – Preparing whole genome libraries using the HiFi prep kit 96 (103-420-700)

Procedure & checklist [103-420-700](#) describes the steps for constructing whole genome sequencing (WGS) libraries from genomic DNA using HiFi prep kit 96 and describes the polymerase binding steps using Revio polymerase kit 96. This workflow is intended as a high-throughput library prep method and has been optimized for use with liquid handler automation¹.

Overview	
Applications	WGS of human, animal, or plant samples
Samples	24–96 using automation (1–96 when doing manual preps)
Minimum automated batch size	24
Maximum automated batch size	96

DNA input and fragment size recommendations	
gDNA per Revio SMRT® Cell	2 µg
DNA shearing	Automated pipette-tip shearing
Target fragment lengths	15–20 kb
Size selection	SRE on gDNA, and 3.1X (35% v/v) AMPure® PB on HiFi library
Average library recovery	15% of genomic DNA input (dependent on input gDNA quality)

Note: Procedure & checklist 103-420-700 includes instructions for HiFi prep kit 96 **SMRTbell library construction workflow** *and* **sequencing preparation (ABC²) workflow using Revio polymerase kit 96²**

PacBio

Preparing whole genome libraries using the HiFi prep kit 96

Procedure & checklist

Overview

This procedure describes the steps for constructing whole genome sequencing (WGS) libraries from genomic DNA using the HiFi prep kit 96 and describes the polymerase binding steps using the Revio™ polymerase kit 96. This workflow is intended as a high-throughput library prep method and has been optimized for use with liquid handler automation.

Overview		
Applications	WGS of human, animal, or plant samples	
Samples	24–96 using automation (1–96 when doing manual preps)	
Minimum automated batch size	24	
Maximum automated batch size	96	
Workflow time	Automation time	
	Hamilton NGS STAR	Hamilton Microlab Prep
SRE	3.5 hours for 96 samples	3 hours for 24 samples
Shearing	10 min for 24–96 samples	22 min for 24 samples
Library prep	6.5 hours for 96 samples (start from post shearing cleanup)	1.5 hours for 24 samples (post-shearing cleanup only)
Anneal, bind, cleanup (ABC)	2.5 hours for 96 samples	N/A
Average total time	13 hours	5 hours

DNA input and fragment size recommendations	
gDNA per Revio SMRT® Cell	2 µg
DNA shearing	Automated pipette-tip shearing
Target fragment lengths	15–20 kb
Size selection	SRE on gDNA, and 3.1X (35% v/v) AMPure® PB on HiFi library
Average library recovery	15% of genomic DNA input (dependent on input gDNA quality)

© 2024 PacBio. All rights reserved. Research use only. Not for use in diagnostic procedures.
103-420-700 REV01 MAR2024

PacBio

¹ If using the HiFi prep kit 96 or Revio HiFi prep kit 96 in an **automated** high-throughput library prep workflow with a Hamilton liquid handling system, also refer to **Guide & overview – Automated HiFi prep 96 and HiFi annealing, binding, and cleanup for the Hamilton NGS Microlab STAR system** ([103-425-700](#)) for specific details surrounding automation equipment setup recommendations.


² If using **Revio polymerase kit 96 (103-253-600)** for polymerase binding: For primer annealing, polymerase binding & complex cleanup (ABC) steps, follow sample setup instructions for HiFi prep kit 96 libraries in **Procedure & checklist – Preparing whole genome libraries using the HiFi prep kit 96** ([103-420-700](#)). To perform the final dilution step after completing ABC with Revio polymerase kit 96, follow the instructions provided in the **SMRT Link Loading Calculator** to prepare HiFi prep kit 96 samples for sequencing. **Note:** If using **Revio polymerase kit (102-817-600)** or **Sequel II binding kit 3.2 (102-333-300)** for polymerase binding, then follow sample setup ABC and final loading dilution instructions provided in **SMRT Link Sample Setup Calculator for Sequel II binding kit 3.1/3.2, Revio polymerase kit.**

Guide & overview – Automated HiFi prep 96 and HiFi annealing, binding, and cleanup for the Hamilton NGS Microlab STAR system (103-425-700)¹

Guide & overview 103-425-700 describes the automated workflow for constructing whole genome sequencing (WGS) libraries from genomic DNA using the Revio HiFi prep kit 96 and the Revio polymerase kit 96. The SRE HT kit, HiFi Prep 96 kit and the Revio Polymerase 96 kit are designed for a minimum of 24 and maximum of 96 samples per automated run.

Overview	
Applications	WGS of human, animal, or plant samples
Samples	24–96 using automation
Minimum automated batch size	24
Maximum automated batch size	96
Workflow time	Automation time
SRE	3.5 hours for 96 samples
Shearing	10 min for 24–96 samples
Library prep	6.5 hours for 96 samples (start from post-shearing cleanup)
Anneal, bind, cleanup (ABC)	2.5 hours for 96 samples
Average total time	13 hours

Note: Guide & overview 103-425-700 includes instructions for **automated** HiFi prep kit 96 **SMRTbell library construction workflow** **and** sequencing preparation (ABC²) workflow using Revio polymerase kit 96²




¹ If using the HiFi prep kit 96 or Revio HiFi prep kit 96 in a **manual** high-throughput library prep workflow, refer to **Procedure & checklist – Preparing whole genome libraries using the HiFi prep kit 96 (103-420-700)**.

² If using **Revio polymerase kit 96 (103-253-600)** for polymerase binding: For primer annealing, polymerase binding & complex cleanup (ABC) steps, follow sample setup instructions for HiFi prep kit 96 libraries in **Guide & overview – Automated HiFi prep 96 and HiFi annealing, binding, and cleanup for the Hamilton NGS Microlab STAR system (103-425-700)**. To perform the final dilution step after completing ABC with Revio polymerase kit 96, follow the instructions provided in the **SMRT Link Loading Calculator** to prepare HiFi prep kit 96 samples for sequencing. **Note:** If using **Revio polymerase kit (102-817-600)** or **Sequel II binding kit 3.2 (102-333-300)** for polymerase binding, then follow sample setup ABC and final loading dilution instructions provided in **SMRT Link Sample Setup Calculator for Sequel II binding kit 3.1/3.2, Revio polymerase kit**.

HiFi prep kit 96 workflow overview for human/animal/plant WGS applications

Automation of HiFi library prep workflow enables high-throughput processing of 24 to 96 samples in 2 days



General best practices recommendations for preparing WGS libraries using HiFi prep kit 96

Automation

- The HiFi prep kit 96 and workflow was **designed for use with NGS liquid handling automation systems**
- As a result, this protocol (***Preparing whole genome libraries using the HiFi prep kit 96*** ([103-425-700](#)) is intended to describe the SRE, shearing, library prep enzymatic reactions, and bead cleanups for automation method development, or in certain instances manual preparation
- **Note:** This protocol was developed using the **Hamilton NGS STAR MOA 96 system**
 - Refer to ***Guide & overview – Automated HiFi prep 96 and HiFi annealing, binding, and cleanup for the Hamilton NGS Microlab STAR system*** ([103-425-700](#))
- Because of differences between automation instruments, modifications not described herein may be needed to adapt the protocol to user-specific instrumentation
 - Please visit www.pacb.com or contact your local support team for a list of instruments with a PacBio qualified method



Hamilton NGS Star workstation

Genomic DNA sample extraction

- **PacBio Nanobind DNA extractions kits** are recommended to ensure sufficient mass and quality of high-molecular weight (HMW) DNA for this protocol

Nanobind HT CBB kit (102-762-700; 96 rxn)	Nanobind HT 1 mL blood kit (102-762-800; 96 rxn)	Nanobind PanDNA kit (103-260-000; 24 rxn)
<ul style="list-style-type: none">• For high-throughput HMW DNA extraction from up to 200 μL human/mammalian blood, non-mammalian animal blood¹, cultured cells, and bacteria• Expected HMW DNA yield: 3–15 μg for blood and cultured mammalian cells and 2–10 μg for bacteria	<ul style="list-style-type: none">• For high-throughput HMW DNA extraction from 1 mL human blood• Expected HMW DNA yield: 3–70 μg	<ul style="list-style-type: none">• For HMW DNA extraction from cells, blood, bacteria, tissues, insects, and plant nuclei• Expected HMW DNA yield: 3–26 μg



PacBio Nanobind kits

General best practices recommendations for preparing WGS libraries using HiFi prep kit 96

DNA sizing QC

- Agilent **Femto Pulse system**¹ is highly recommended for the accurate sizing of genomic DNA samples
- Femto Pulse system enables simple, rapid sizing QC of genomic DNA and SMRTbell libraries, and conserves sample by using femtogram ranges of input DNA
 - Resolves fragments 1,300 bp to 165 kb using gDNA 165 kb Analysis kit (can resolve 100 – 6,000 bp using Ultra Sensitivity NGS kit)
 - Requires <1 ng of sample DNA
 - Can analyze up to 12 samples in <1.5 hrs
 - Outputs quality metrics such as Genomic Quality Number (GQN)² to quickly score integrity of HMW gDNA



Femto Pulse system
(Agilent Technologies)

DNA quantification QC

- For high-throughput DNA quantification QC workflows, we recommend using the Quant-iT 1X dsDNA high sensitivity assay kit³ (Thermo Fisher Scientific) with the Varioskan LUX multimode microplate reader (Thermo Fisher Scientific) (or similar instrument)
 - Varioskan LUX multimode microplate reader enables rapid, specific and accurate determination of nucleic acid concentrations in a wide range
 - Quant-iT assay is well-adapted to high-throughput use (typically in 96-well or 384-well plates) and is highly selective for dsDNA over RNA (in the 0.2–100 ng range, the fluorescence signal is linear with the amount of DNA)
 - **Note:** We do not recommend quantification with UV-Vis Spectrophotometers (e.g. NanoDrop) that measure all nucleic acids in a sample. For example, measuring all nucleic acid will inflate the true concentration of gDNA in samples



Varioskan LUX multimode
microplate reader
(Thermo Fisher Scientific)

¹ See *Product Note – HiFi WGS sequencing with the Agilent Femto Pulse system* ([102-326-561](#)) for more details.

² See *Application Note – Quality Metrics for Nucleic Acids with the Agilent Fragment Analyzer and Femto Pulse Systems* (Agilent [5994-0521EN](#))

³ Alternatively, for lower-throughput applications DNA quantification QC may be performed without a microplate reader using the **Qubit 1X dsDNA high sensitivity assay kit** with a **Qubit fluorometer**.

Note: Do not use a Qubit Flex fluorometer when performing DNA quantification QC on polymerase-bound SMRTbell library in Loading Buffer 96 since concentration readings will not be accurate.

General best practices recommendations for preparing WGS libraries using HiFi prep kit 96

Recommended genomic DNA input amount and quality

- 70% or more of the DNA should be ≥ 10 kb for this protocol
 - This corresponds to a genome quality number (GQN) of 7.0 or higher at 10 kb
- Recommended DNA input amounts will vary by starting gDNA quality

DNA quality	90% >10 kb	80% >10 kb	70% >10 kb
gDNA input into SRE size selection step	2 – 3 μ g	3 – 4 μ g	4 – 5 μ g

Note: The maximum input gDNA mass tolerated by shearing and library enzymatic reactions in this HiFi prep kit 96 protocol is **3 μ g per reaction**
→ Perform parallel library prep reactions if using >3 μ g input gDNA

- Overall SMRTbell library construction yield is dependent on input gDNA quality and size
 - The recovery from input gDNA to completed SMRTbell library typically ranges between **10 – 25%**
- Starting with **2 μ g** of input gDNA (going into SRE size selection step) will, on average, provide enough library to load **1 Revio SMRT Cell**

Mean library insert size	Library mass needed to load one Revio SMRT Cell at 250 pM OPLC
15,000 bp	243 ng
18,000 bp	292 ng
20,000 bp	341 ng

Note: It is recommended to use at least 2 μ g of input gDNA even if only using one Revio SMRT Cell the following reasons:

- Ensure adequate SRE recovery
- Final mean library size is not known prior beginning protocol.
- Available library for re-sequencing or topping off coverage if necessary






- If targeting larger insert sizes or working with lower quality DNA, start with at least 3 μ g of input gDNA (going into SRE size selection step) to ensure adequate library for optimal SMRT Cell loading

General best practices recommendations for preparing WGS libraries using HiFi prep kit 96

Reagent handling

- Room temperature is defined as any temperature in the range of 18 – 25°C for this protocol

HiFi prep kit 96

Thaw these reagents at room temperature	Keep these temperature-sensitive reagents on ice	Bring these reagents to room temperature 30 minutes prior to use	Bring these reagents to room temperature 1.5 hours prior to use
<input type="checkbox"/> Repair buffer 96	 <input type="checkbox"/> End repair mix 96	<input type="checkbox"/> AMPure PB beads	<input type="checkbox"/> SMRTbell cleanup beads-85 mL ¹
<input type="checkbox"/> Nuclease buffer 96	 <input type="checkbox"/> DNA repair mix 96	<input type="checkbox"/> Elution buffer	
<input type="checkbox"/> SMRTbell adapter index plate	 <input type="checkbox"/> Ligation mix 96	<input type="checkbox"/> dsDNA quantification reagents	
	 <input type="checkbox"/> Ligation enhancer 96		
	 <input type="checkbox"/> Nuclease mix 96		









¹ Alternatively, can bring SMRTbell cleanup beads-85 mL to room temperature the night before if starting protocol in the morning.

- Once thawed, reaction buffers and adapter index plate may be stored on a cold block, at 4°C, or on-ice prior to making master mix or placing on the liquid handler work deck
- Briefly spin down all reagent tubes in a microcentrifuge to collect all liquid at the bottom
- Briefly vortex then spin down SMRTbell adapter index plate in a centrifuge with a plate adapter to collect all liquid at the bottom of the wells
- Shake/vortex SMRTbell cleanup beads and AMPure PB beads immediately before use

General best practices recommendations for preparing WGS libraries using HiFi prep kit 96

Reagent handling (cont.)

Revio polymerase kit 96

Thaw these reagents at room temperature		Keep these reagents on a cold block or on ice		Bring these reagents to room temperature 30 minutes prior to use		Bring these reagents to room temperature 1.5 hours prior to use	
	Annealing buffer 96		Sequencing polymerase 96		Loading buffer 96	SMRTbell cleanup beads-85 mL ¹	
	Standard sequencing primer 96		Sequencing control 96				
	Polymerase buffer 96						
	Loading buffer 96						
	Dilution buffer 96						

¹ Alternatively, can bring SMRTbell cleanup beads-85 mL to room temperature the night before if starting protocol in the morning.

- Once thawed, reaction buffers and sequencing primer may be stored on a cold block, at 4°C, or on-ice prior to making master mix or placing on the liquid handler work deck
- Loading buffer 96 should be left at room-temperature
- **Note:** Loading buffer 96 is light sensitive and should be protected from light when not in use

General best practices recommendations for preparing WGS libraries using HiFi prep kit 96

Multiplexing samples

- All libraries constructed for each DNA sample using this protocol will include a **SMRTbell adapter index ('barcode')**
- Use the SMRT Link (v13.1+) Sample Setup **Pooling Calculator tool** to help determine appropriate volumes to use for multiplexing libraries prepared with HiFi prep kit 96 and bound with Revo polymerase kit 96
- Prior to pooling HiFi libraries together please consider the following guidelines:
 - Only **pool samples with similar genome sizes to ensure balanced coverage**
 - Ensure that samples to be pooled have a **similar mean insert size and similar insert length size distribution**¹
 - Aim to pool samples in an equal molar concentration for best balanced coverage
- It is recommended to **pool HiFi libraries post-ABC** (i.e., **after** performing primer annealing, polymerase binding and complex cleanup) for the following reasons:
 - **Prevent an inhibitor in one sample from affecting the polymerase binding of all samples** in a pool
 - **Ability to quickly pool different libraries together on additional runs** to “top off” coverage (any un-pooled complexed library is available for future sequencing runs without having to re-do ABC)

Pooling Calculator

Number of samples to be multiplexed ↑ Import Export Print

Pooled library target volume (μL)

Concentration output units

Pooled library concentration (ng/μL)

Buffer volume (μL) to add to pooled sample: **8.66**

Sample name	Conc. (ng/μL)	Pooling volume (μL)
Sample 1	3.8	23.03
Sample 2	4.78	18.31


Close

SMRT Link Sample Setup Pooling Calculator tool can be used to calculate the required volumes of HiFi prep kit library samples needed for pooling when performing multiplexed sequencing on a single SMRT Cell.

General best practices recommendations for preparing WGS libraries using HiFi prep kit 96

Sequencing preparation (ABC) and polymerase-bound library storage

- **Note:** *Procedure & checklist* – *Preparing whole genome libraries using HiFi prep kit 96* ([103-420-700](#)) brings the final SMRTbell library through the primer annealing, polymerase binding, and complex cleanup (ABC) sample setup steps
 - Do not need to use SMRT Link Sample Setup to perform ABC calculations for HiFi prep kit 96 libraries.
 - Follow SMRT Link Revio polymerase kit 96 **Loading Calculator** tool to perform final loading procedure for HiFi prep kit 96 libraries
- **Sequencing polymerase is stable once bound to the HiFi prep kit library and can be stored at 4°C or frozen at -20°C.**
- Stored polymerase-bound library shows equivalent loading to freshly prepared bound libraries up to the recommendations listed below.
 - Recommended polymerase-bound storage:¹**
 - **Polymerase-bound library is stable at 4°C for 1 month**
 - **Frozen polymerase-bound library is stable for at least 6 months**
- Stored polymerase-bound library needs to be **protected from light** while stored.




HiFi prep kit 96 library prep inputs & expected step yields

Final HiFi library yield is typically sufficient to load ≥ 1 SMRT Cells

	Cleanup	DNA QC	Step input	Step output / Yield (%)
DNA sample extraction	N/A	N/A	200 μ L human or animal whole blood / 1 mL human whole blood / 1×10^6 cells cultured mammalian cells	Nanobind HT CBB kit : 3 – 15 μ g HMW DNA Nanobind HT 1 mL blood kit 3 – 70 μ g HMW DNA
DNA sample QC	N/A	Quant-iT/Qubit dsDNA HS assay Femto Pulse system	Nanobind HT-extracted HMW DNA	DNA quantification QC \rightarrow Aim for ≥ 3 μ g HMW DNA DNA sizing QC \rightarrow Ideally GQN(10 kb) ≥ 7.0
SRE 1 Short Read Eliminator on gDNA	N/A	N/A	40-100 ng/ μ L DNA in a total volume of 50 μ L Input DNA mass depends on expected recovery	Expect ~75% step recovery or 75% total recovery
Shear 2 Pipette DNA shearing	1X SMRTbell cleanup beads	Post-shearing QC with Quant-iT/Qubit (optional) & Femto Pulse	≤ 10 ng/ μ L DNA in 300 μ L (3 μ g total input DNA mass) into Step 2	Expect ~80% step recovery or 60% total recovery; target DNA shear size is ~15-20 kb
3 Post-shearing cleanup				
Library construction 4 DNA repair & A-tailing	N/A	N/A	49 μ L sheared DNA	60 μ L repaired & A-tailed DNA
5 Adapter ligation	1X SMRTbell cleanup beads	N/A	60 μ L post-repaired & A-tailed DNA into Step 5	Expect ~80% step recovery or 48% total recovery
6 Post-ligation cleanup				
7 Nuclease treatment	N/A	N/A	40 μ L of post-ligation cleanup sample	Expect ~40% step recovery or 19% total recovery
8 AMPure PB bead cleanup	3.1X (35%) AMPure PB	Quant-iT/Qubit dsDNA HS assay Femto Pulse system (optional)	50 μ L of nuclease-treated sample	Expect 80% step recovery or 15% total recovery (can range from ~10% - 25%)
ABC 9 Annealing, binding & cleanup	1X SMRTbell cleanup beads	Quant-iT/Qubit dsDNA assay ¹	DNA concentration must be less than 60 ng/ μ l to proceed with ABC	Expect ~70% step recovery relative to starting library input mass for ABC

DNA sample extraction

Perform automated high-throughput HMW DNA extraction using Nanobind HT kits¹



ThermoFisher
SCIENTIFIC



KingFisher Apex

HAMILTON



Microlab NIMBUS Presto

Automated solutions built on Nanobind technology

Nanobind HT kits use **magnetic disk processing** to automate lysis, binding, washing and elution steps and are compatible with instruments from Hamilton and Thermo Fisher

- Hamilton NIMBUS Presto is a walkaway solution with **automated** plate filling
- Thermo Fisher KingFisher instruments are **semi-automated** with manual plate filling and limited user interaction during the run

Available Nanobind HT kits



Nanobind HT CBB kit (102-762-700; 96 rxn)

- For up to 200 μ L human/mammalian blood, non-mammalian animal blood², cultured cells, and bacteria
- Expected HMW DNA yield: 3–15 μ g for blood and cultured mammalian cells and 2–10 μ g for bacteria

Nanobind HT 1 mL blood kit (102-762-800; 96 rxn)


- For 1 mL human blood
- Expected HMW DNA yield: 3–70 μ g

¹ See *Technical overview — Automated high-throughput HMW DNA extraction for PacBio long-read sequencing using Nanobind HT kits* ([103-401-700](#)).

² Note: For animals with nucleated red blood cells, only 2.5 – 20 μ L of blood is needed for HMW DNA extraction using Nanobind HT CBB kit.

DNA sample extraction (cont.)


Available Nanobind HT procedures



Automation system	Procedure & checklist	Sample type
KingFisher Apex	Extracting HMW DNA using the Nanobind HT CBB kit for mammalian cultured cells on KingFisher Apex system [102-996-100]	Mammalian cells
	Extracting HMW DNA using the Nanobind HT CBB kit for 200 μ L human whole blood on KingFisher Apex system [102-995-700]	Human whole blood
	Extracting HMW DNA using the Nanobind HT CBB kit for bacteria on the KingFisher Apex system [103-377-600]	Cultured bacteria
	Extracting HMW DNA using the Nanobind HT CBB kit for non-human mammalian blood (NHMB) on the KingFisher Apex system [103-397-300]	Non-human mammalian blood
	Extracting HMW DNA using the Nanobind HT CBB kit for nucleated red blood cells (nRBCs) on the KingFisher Apex system [103-377-800]	Non-mammalian blood (nucleated RBCs)
	Extracting HMW DNA using Nanobind HT 1 mL blood kit for human whole blood on KingFisher Apex system [102-995-300]	Human whole blood
Hamilton NIMBUS Presto	Extracting HMW DNA using the Nanobind HT CBB kit for mammalian cultured cells on Hamilton NIMBUS Presto system [102-996-400]	Mammalian cells
	Extracting HMW DNA using the Nanobind HT CBB kit for 200 μ L human whole blood on Hamilton NIMBUS Presto system [102-996-000]	Human whole blood
	Extracting HMW DNA using the Nanobind HT CBB kit for bacteria on Hamilton NIMBUS Presto system [103-397-400]	Cultured bacteria
	Extracting HMW DNA using the Nanobind HT CBB kit for non-human mammalian blood (NHMB) on the Hamilton NIMBUS Presto system [103-377-700]	Non-human mammalian blood
	Extracting HMW DNA using the Nanobind HT CBB kit for nucleated red blood cells (nRBCs) on the Hamilton NIMBUS Presto system [103-397-500]	Non-mammalian blood (nucleated RBCs)
	Extracting HMW DNA using the Nanobind HT 1 mL blood kit for human whole blood on Hamilton NIMBUS Presto system [102-995-600]	Human whole blood

DNA sample QC

Perform high-throughput DNA QC using a microplate reader¹ and a Femto Pulse system



DNA quantification QC



Varioskan LUX multimode microplate reader
(Thermo Fisher Scientific)

Varioskan LUX multimode microplate reader enables rapid, specific and accurate determination of nucleic acid concentrations in high-throughput workflows¹

- Perform DNA concentration measurements using **Quant-iT dsDNA assay²** (Thermo Fisher Scientific)
- Quant-iT assay is well-adapted to high-throughput use, typically in **96-well** or **384-well plates**
- Quant-iT assay is highly selective for dsDNA over RNA, and in the range of **0.2–100 ng**, the fluorescence signal is linear with the amount of DNA

¹ Alternatively, for lower-throughput applications DNA quantification QC may be performed without a microplate reader using the Qubit 1X dsDNA high sensitivity assay kit. **Note:** Do not use a Qubit Flex fluorometer when performing DNA quantification QC on polymerase-bound SMRTbell library in Loading buffer 96 (from Revo polymerase kit 96 PN 103-253-600) or Loading Buffer (from Revo polymerase kit PN 102-817-600) since concentration readings will not be accurate.

DNA sizing QC




Femto Pulse system
(Agilent Technologies)

Femto Pulse system enables simple, rapid sizing QC of genomic DNA and SMRTbell libraries, and conserves sample by using femtogram ranges of input DNA

- Use the **Femto Pulse gDNA 165 kb analysis kit** (FP-1002-0275)
- Dilute samples to **250 pg/uL**
- **70% or more of the DNA should be ≥10 kb** for this protocol. This corresponds to a genome quality number (GQN)³ of 7.0 or higher at 10 kb.

Short Read Eliminator on gDNA

Perform high-throughput size selection on input genomic DNA using SRE HT kit to remove <10 kb fragments



1. Short Read Eliminator


Step	Instructions for SRE on gDNA
1.1	Bring the DNA samples to a concentration between 40 and 100 ng/μL in a total volume of 50 μL using Buffer LTE in a 1.5 ml LoBind tube or hard-shell plate. The DNA input requirement for shearing is <3 μg so gDNA input into SRE will depend on expected recovery.
1.2	Add 50 μL of Buffer SRE to each sample. If working in a plate format, heat seal with foil. Vortex/shake to mix for 5 seconds at max speed.
1.3	Incubate the sample for 1 hour at 50°C in a heat block or thermal cycler. After incubation, if using a plate format, ensure that it is compatible with a 300 μL elution. If not, transfer to the appropriate deep well plate after incubation and seal with an adhesive seal.
1.4	Load plate or tube (with the hinge facing toward the outside of the rotor) into centrifuge.
1.5	Centrifuge a tube at 10,000 rcf for 30 minutes. Centrifuge a plate at >2250 rcf (max 3220 rcf) for 1 hour. <ul style="list-style-type: none"> If using a centrifuge with temperature control (i.e., cooling function), turn this function off by specifying a target temperature set point higher than ambient room temperature (e.g., 29°C or 30°C).
1.6	Carefully remove supernatant without disturbing the pellet. <ul style="list-style-type: none"> Leaving up to 10 μL is acceptable to be sure the pellet is not disturbed.
1.7	Add 300 μL of Buffer LTE to the tube and incubate at room temperature for 10 minutes.
1.8	After incubation, pipette-mix 20 times and vortex/shake the tube/sealed plate for 15s to ensure that the DNA is properly re-suspended and mixed.
1.9	Quantify the resuspension to measure DNA recovery. If the recovery is lower than 50%, repeat pipette-mixing 20 times and vortex/shake. If the recovery is greater than 50%, proceed to next step (DNA shearing).
1.10	Proceed to automated DNA shearing. It is recommended to proceed to DNA shearing within 2 weeks of performing SRE

SAFE STOPPING POINT - Store at 4°C

IMPORTANT! Use SRE on only genomic DNA. Attempting to use SRE on HiFi libraries (post-SMRTbell library construction) will result in poor recoveries and potential loss of the entire library.

- DNA samples should be between 40 to 100 ng/μL in a total volume of 50 μL using Buffer LTE
- DNA input requirement into shearing is <3 μg so gDNA input into SRE will depend on expected recovery

SRE workflow for HiFi prep kit 96 protocol¹




¹ **Note:** SRE size selection procedure described in high-throughput HiFi prep kit 96 protocols ([103-420-700](#) and [103-425-700](#)) uses slightly different starting sample concentration and volume ranges from the low-throughput SRE protocol described in *Procedure & checklist – Removing short DNA fragments with the Short Read Eliminator (SRE) kit* ([102-982-300](#)).

² Please refer to *Automated HiFi prep 96 and HiFi annealing, binding, and cleanup for the Hamilton NGS Microlab STAR system* ([103-425-700](#)) for details on required consumables.

Pipette DNA shearing

Perform automated high-throughput DNA shearing for WGS samples using Hamilton automation¹



2. Pipette DNA shearing

Step	Instructions for automated DNA shearing on Hamilton systems
2.1	Adjust DNA concentration to ≤ 10 ng/ μ L, if necessary (e.g. if more than 3 μ g of gDNA was recovered from SRE). Use Buffer LTE to dilute samples. Bring all samples up to 300 μ L in a 0.8 mL, 96 DeepWell plate (Thermo Fisher Scientific AB0859).
Parameters for shearing on the Microlab Prep, or Hamilton assay-ready workstations are listed below. These parameters should already be part of the installed method on the instrument.	
Parameter	Setting
DNA concentration	≤ 10 ng/ μ L
Volume of Buffer LTE	300 μ L
Number of mixes	300 cycles
Pipette mixing speed	500 μ L/se
Mix volume	83% volume
Liquid following, cLLD	On
Pipette tip	300 μ L CO-RE II tips (filtered, black, non-sterile)
2.2	
2.3	Place the plate on the appropriate work deck position and start the shearing procedure. Optional: measure DNA concentration with a Qubit fluorometer using the 1x dsDNA HS kit once the shearing procedure is complete.
2.4	Recommended: Further dilute each aliquot to 250 pg/ μ L with Femto Pulse dilution buffer. Measure the final SMRTbell library size distribution with a Femto Pulse system to ensure efficient shearing.
2.5	Proceed to the 1X SMRTbell cleanup bead procedure to concentrate samples for library preparation.

• For automated pipette-based shearing using Hamilton systems, DNA samples should be at **<10 ng/ μ L** in a total volume of **300 μ L** in a 0.8 mL, 96 DeepWell plate (use Buffer LTE to dilute samples)



Hamilton NGS STAR/STARlet/STAR V



Hamilton Microlab Prep

Pipette-based DNA shearing can be completed within ~10 min for up to 96 samples using Hamilton liquid handling instruments.


• **Note:** Post-shear cleanup step should be performed on the **same day** as the library construction steps to ensure there is sufficient volume of SMRTbell cleanup beads to complete the entire library prep workflow

IMPORTANT!

- A mean fragment size between **15 to 20 kb** is recommended for this protocol
- In addition, the distribution of fragment sizes should be narrow and generally **between 10 to 30 kb**
 - Fragments that are too short produce less yield per read, and fragments that are too long may result in lower read accuracy and are less likely to produce HiFi reads
- Deviating from the concentration and automation settings specified for this HiFi prep kit 96 workflow is not recommended and will result in under-sheared DNA

Post-shearing cleanup

Perform post-shearing cleanup using 1X SMRTbell cleanup beads



3. Post-shearing cleanup

Step	Instructions for SMRTbell cleanup bead step
3.1	Add 300 μL (1.0X) of resuspended, room-temperature SMRTbell cleanup beads to each sample. Note: If using less than 300 μL , add 1.0X (v/v) concentration of SMRTbell cleanup beads.
3.2	Pipette-mix the sample until the beads are evenly distributed.
3.3	Leave at room temperature for 10 minutes to allow DNA to bind to the beads.
3.4	Place samples on an appropriate magnet and allow beads to separate fully from the solution. The solution should be clear and beads pelleted to the magnet.
3.5	Slowly pipette off the supernatant without disturbing the beads. Discard the supernatant.
3.6	Slowly dispense 200 μL , or enough to cover the beads, of freshly prepared 80% ethanol to each sample. After 30 seconds , pipette off the 80% ethanol and discard.
3.7	Repeat the previous step.
3.8	Remove residual 80% ethanol: <ul style="list-style-type: none">Remove the samples from the magnetQuick-spin to collect liquid at the bottom of the tube or well.Place the tube or plate back in a magnetic separation rack until beads separate fully from the solution.Pipette off residual 80% ethanol and discard. Alternatively, air dry samples for 1 minute to allow residue ethanol to evaporate. Do not let the bead pellet completely dry out.
3.9	Remove samples from the magnet and immediately resuspend the beads with 49 μL of elution buffer .
3.10	Resuspend by pipetting mixing until beads are evenly distributed in solution.
3.11	Leave samples at room temperature for 5 minutes to elute DNA off beads.
3.12	Place samples back on the magnet and allow beads to separate fully from the solution. The solution should be clear and beads pelleted to the magnet before aspirating the supernatant.
3.13	Slowly pipette off the cleared supernatant without disturbing the beads. Transfer supernatant to a new tube strip or plate. Discard old tube strip or plate with beads.


SAFE STOPPING POINT - Store at 4°C

- **IMPORTANT!** Allow SMRTbell cleanup beads to come up to room temperature by bringing them out of 4°C storage at least 1.5 hrs prior to use
- If performing cleanup in the morning, beads may be left out at room temperature overnight

Note: Post-shear cleanup step should be performed on the **same day** as the library construction steps to ensure that there is a sufficient volume of SMRTbell cleanup beads to complete the entire library prep workflow¹

DNA repair & A-tailing

Repair sites of DNA damage and prepare sheared DNA for ligation to SMRTbell adapter



4. DNA repair & A-tailing

✓ Step Instructions for repair and A-tailing step

Prepare the appropriate volume of master mix with 15% overage using the per reaction volumes listed below.

Repair mix		
Tube	Component	Per rxn vol.
Purple	Repair buffer 96	8 µL
Blue	End repair mix 96	2 µL
Green	DNA repair mix 96	1 µL
Total volume		11 µL

4.1

4.2 Slowly pipette mix the **Repair mix** and quick-spin to collect liquid at the bottom of the tube. If bubbles form during mixing, pulse spin to remove.

4.3 Add **11 µL** of the **Repair mix** to each sample. Total reaction volume should be **60 µL**.

4.4 Pipette-mix the reactions and quick-spin to collect liquid at the bottom of the well/tube. If using a plate format, seal with a heated foil seal.

Run the **Repair and A-tailing** thermocycler program. Set lid temperature to 75°C if programmable.

Step	Time	Temperature
4.5 1	30 min	37°C
2	5 min	65°C
3	Hold	4°C

4.6 Proceed directly to the next step of the protocol.

- Prepare a reaction master mix by adding the required components in the order and volume listed to a new microcentrifuge tube¹
- Adjust component volumes for the number of samples being prepared (e.g., 24, 48, 72 or 96¹), plus 15% overage


- Run **Repair and A-tailing** on-deck thermal cycler program
- Set the lid temperature to 75°C

Adapter ligation & cleanup

Ligate SMRTbell adapter to the ends of each DNA fragment



Revio HiFi prep kit 96 includes SMRTbell adapter index plate 96A (102-009-200), which contains indexes bc2001–bc2096



5. Adapter ligation

✓ Step Instructions for adapter ligation

5.1 Add 4 μL of indexed adapter to each sample from the previous step. Any of the adapters from the four SMRTbell adapter index plates can be used (e.g. 96A, 96B, 96C, or 96D).

Component	Per rxn vol.
SMRTbell adapter index plate 96(A, B, C, or D)	4 μL

Prepare the appropriate volume of master mix with 10% overage using the per reaction volumes listed below.

Ligation mix			
Tube	Component	Volume	
Yellow	Ligation mix 96	20 μL	
Red	Ligation enhancer 96	1 μL	
Total volume		21 μL	

5.5 Pipette mix each reaction to ensure reagents are thoroughly mixed. Quick spin plate if necessary to collect liquid at the bottom of the well.

5.6 Run the **Adapter ligation** thermocycler program. Set lid temperature to 75°C if programmable.

Step	Time	Temperature
1	30 min	20°C
2	Hold	4°C

5.7 Proceed to post-ligation cleanup step

- Add 4 μL of indexed adapter to each sample from the previous step.
- Any of the adapters from the four SMRTbell adapter index plates can be used (e.g. 96A, 96B, 96C, or 96D)

- Prepare a reaction master mix by adding the required components in the order and volume listed to a new microcentrifuge tube¹
- Adjust component volumes for the number of samples being prepared (e.g., 24, 48, 72 or 96¹), plus 10% overage

- Run **Adapter ligation** on-deck thermal cycler program
- Set the lid temperature to 75°C

6. Post-ligation 1X SMRTbell bead cleanup

✓ Step Instructions for SMRTbell cleanup bead step


6.1 Add 85 μL (1.0X) of resuspended, room-temperature SMRTbell cleanup beads to each sample.

6.2 Pipette mix the sample until the beads are evenly distributed.

- **IMPORTANT!** Allow SMRTbell cleanup beads to come up to room temperature by bringing them out of 4°C storage at least 1.5 hours prior to beginning

Nuclease treatment

Remove unligated DNA fragments and leftover SMRTbell adapters from the sample



7. Nuclease treatment

7.1 Prepare the appropriate volume of master mix with 10% overage using the per reaction volumes listed below. **Prepare the master mix immediately before use to ensure optimal activity.**

Tube	Component	Volume
Light purple	Nuclease buffer 96	5 µL
Light green	Nuclease mix 96	5 µL
Total volume		10 µL

7.2 Pipette mix **Nuclease mix** and quick spin to collect liquid.

7.3 Add **10 µL** of **Nuclease mix** to each sample. Total volume should equal **50 µL**.

7.4 Pipette-mix each sample and quick-spin to collect liquid.

7.5 Run the **Nuclease treatment** thermocycler program. Set lid temperature to 75°C if programmable.

Step	Time	Temperature
1	15 min	37°C
2	Hold	4°C


7.6 Proceed to the next step of the protocol.

- Prepare a reaction master mix by adding the required components in the order and volume listed to a new microcentrifuge tube¹
- Adjust component volumes for the number of samples being prepared (e.g., 24, 48, 72 or 96¹), plus 10% overage

- Run **Nuclease treatment** on-deck thermal cycler program
- Set the lid temperature to 75°C

Diluted AMPure PB cleanup and size selection

AMPure PB bead size cleanup and selection step will clean the library and deplete DNA fragments <3 kb



8. Diluted AMPure PB cleanup and size selection

Step	Instructions for AMPure PB bead cleanup
8.1	Make a 35% v/v dilution of AMPure PB beads by adding 1.75 mL of resuspended AMPure PB beads to 3.25 mL of elution buffer. The 35% dilution can be stored at 4°C for 30 days. Note: The AMPure PB dilution may be scaled as appropriate for smaller/larger scale projects.
8.2	Add 3.1X v/v (155 µL) of resuspended, room-temperature 35% AMPure PB beads to each sample from the previous step.
8.3	Pipette-mix the beads until evenly distributed.
8.4	Leave at room temperature for 20 minutes to allow DNA to bind beads.
8.5	Place sample on an appropriate magnet and allow beads separate fully from the solution.
8.6	Slowly pipette off the cleared supernatant without disturbing the beads.
8.7	Slowly dispense 200 µL , or enough to cover the beads, of freshly prepared 80% ethanol into each sample. After 30 seconds , pipette off the 80% ethanol and discard.
8.8	Repeat the previous step.
8.10	Remove samples from the magnet and immediately add 25 µL of elution buffer to each sample.
8.11	Pipette-mix the beads until evenly distributed.
8.12	Leave at room temperature for 5 minutes to elute DNA of the beads.
8.13	Place samples on the magnet and allow the beads separate fully from the solution.
8.14	Slowly pipette off the cleared supernatant without disturbing the beads. Transfer supernatant to a new tube.
8.15	Optional: Take a 1 µL aliquot from each tube and dilute with 9 µL of elution buffer or water . Measure DNA concentration with a Qubit fluorometer using the 1x dsDNA HS kit. Calculate the total mass. The final overall recovery should be 10–25% as measured from gDNA input to completed SMRTbell library (includes SRE, shearing, and library prep). DNA concentration must be less than 60 ng/µL to proceed to ABC; however, libraries typically are at <40 ng/µL after the SMRTbell library preparation process. Optional: Further dilute each aliquot to 250 pg/µL with Femto Pulse dilution buffer. Measure final SMRTbell library size distribution with a Femto Pulse system.

SAFE STOPPING POINT - Store at 4°C

- Prepare a **35% (v/v) dilution of AMPure PB beads** using elution buffer
- 35% AMPure PB solution can be stored at **4°C for 30 days**
- **Note:** The AMPure PB dilution procedure may be scaled as appropriate for smaller-/larger-scale projects (each sample requires 155 µL of 35% AMPure PB beads)

- **Optional:** Perform **DNA concentration QC** on final purified HiFi prep kit 96 SMRTbell library using a Qubit dsDNA HS assay or Quant-iT ds DNA HS assay
 - Final overall recovery should be **10 – 25%** as measured from gDNA input to completed SMRTbell library (includes SRE, shearing, and library prep)
- **Note:** Final HiFi prep kit 96 SMRTbell library concentration must be **<60 ng/ µL** to proceed with annealing, binding & cleanup (ABC)¹
 - Using a concentration above 60 ng/µL will result in lower *P1* loading during sequencing
- **Optional:** Perform **DNA sizing QC** on final purified HiFi prep kit 96 SMRTbell library using a Femto Pulse system




Varioskan LUX multimode microplate reader (Thermo Fisher Scientific)



Femto Pulse system (Agilent Technologies)

Annealing, Binding, and Cleanup (ABC)

Perform primer annealing, polymerase binding and complex cleanup to prepare samples for sequencing



9.1. Anneal sequencing primer

Step	Instructions															
9.1	Prepare the appropriate volume of master mix with 10% overage using the per reaction volumes listed below.															
<table border="1"> <thead> <tr> <th colspan="3">Annealing mix</th> </tr> <tr> <th>Tube</th> <th>Component</th> <th>Volume</th> </tr> </thead> <tbody> <tr> <td>Light blue</td> <td>Annealing buffer 96</td> <td>12.5 μL</td> </tr> <tr> <td>Light green</td> <td>Standard sequencing primer 96</td> <td>12.5 μL</td> </tr> <tr> <td colspan="2">Total volume</td> <td>25 μL</td> </tr> </tbody> </table>		Annealing mix			Tube	Component	Volume	Light blue	Annealing buffer 96	12.5 μ L	Light green	Standard sequencing primer 96	12.5 μ L	Total volume		25 μ L
Annealing mix																
Tube	Component	Volume														
Light blue	Annealing buffer 96	12.5 μ L														
Light green	Standard sequencing primer 96	12.5 μ L														
Total volume		25 μ L														
9.2	Pipette-mix the Annealing mix and quick spin to collect liquid.															
9.3	Add 25 μ L of the Annealing mix to each library. Total volume should equal 50 μ L.															
9.4	Pipette-mix each sample and quick spin to collect liquid.															
9.5	Incubate at room temperature for 15 minutes .															
9.6	During primer incubation, prepare the polymerase dilution (see below) and store on ice.															

- For primer annealing reaction, each sample requires **25 μ L of Annealing mix solution**
- Scale up the appropriate volume of Annealing mix solution required to process all samples by adding 1 volume of **Standard sequencing primer 96 + 1 volume of Annealing buffer 96** to a new tube



Annealed SMRTbell library

9.7. Bind polymerase

Step	Instructions															
9.7	To prepare the polymerase, add the following components to a new microcentrifuge tube on ice. Adjust component volumes for the number of samples being prepared, plus 10% overage.															
<table border="1"> <thead> <tr> <th colspan="3">Polymerase Dilution</th> </tr> <tr> <th>Tube</th> <th>Component</th> <th>Volume</th> </tr> </thead> <tbody> <tr> <td>Yellow</td> <td>Polymerase buffer 96</td> <td>47 μL</td> </tr> <tr> <td>Purple</td> <td>Sequencing polymerase 96</td> <td>3 μL</td> </tr> <tr> <td colspan="2">Total volume</td> <td>50 μL</td> </tr> </tbody> </table>		Polymerase Dilution			Tube	Component	Volume	Yellow	Polymerase buffer 96	47 μ L	Purple	Sequencing polymerase 96	3 μ L	Total volume		50 μ L
Polymerase Dilution																
Tube	Component	Volume														
Yellow	Polymerase buffer 96	47 μ L														
Purple	Sequencing polymerase 96	3 μ L														
Total volume		50 μ L														
9.8	Pipette mix the polymerase dilution and quick-spin to collect liquid.															
9.9	Add 50 μ L of polymerase dilution to primer annealed sample. Total volume should equal 100 μ L.															
9.10	Pipette-mix each sample and quick-spin to collect liquid.															
9.11	Incubate at room temperature for 15 minutes .															
9.12	Proceed immediately to the next step of the protocol to remove excess polymerase.															


- For binding reaction, each annealed sample requires **50 μ L of Revio sequencing polymerase 96 working solution**
- Scale up the appropriate volume of Revio polymerase 96 working solution required to process all samples by diluting the **Revio polymerase 96 stock 16.7-fold in Polymerase buffer 96**



Polymerase-bound SMRTbell library

Annealing, Binding, and Cleanup (ABC) (cont.)

Perform primer annealing, polymerase binding and complex cleanup to prepare samples for sequencing



9.13. Clean up polymerase-bound complexes


Step	Instructions
9.13	Add 100 μL of resuspended, room-temperature SMRTbell cleanup beads to each sample
9.14	Pipette-mix the beads until evenly distributed and quick-spin if necessary to collect all liquid from the sides of the tube.
9.15	Leave at room temperature for 10 minutes to allow DNA to bind beads
9.16	Place sample on an appropriate magnet and allow beads to separate fully from the solution
9.17	Slowly pipette off the cleared supernatant without disturbing the beads. Discard the supernatant. DO NOT USE EtOH . Proceed immediately to the elution. It is important not to let the beads dry out.
9.18	Remove sample from the magnet and immediately add 50 μL of Loading Buffer 96 to each tube and resuspend the beads by pipette mixing.
9.19	Quick-spin the samples to collect any liquid from the sides of the tube.
9.20	Leave at room temperature for 5 minutes to elute DNA
9.21	Place sample on magnet and allow beads to separate fully from the solution.
9.22	Slowly pipette off the cleared supernatant without disturbing the beads. Transfer supernatant to a new tube . Discard the old tube with beads
9.23	Use 1 μL of sample to measure DNA concentration with a Qubit fluorometer using the 1x dsDNA HS kit. Important: The Qubit Flex instrument is not compatible with measuring polymerase-bound library in Loading Buffer 96. Concentration readings will not be accurate.
9.24	Proceed to the Loading Calculator in SMRT Link v13.1 or higher to calculate the final dilution for adding the sample to the Revio sequencing plate.

PROTOCOL COMPLETE

Storing polymerase-bound HiFi prep kit 96 libraries

- Polymerase-bound libraries can be stored at **4°C** for **1 month**, or at **-20°C** for **>6 months** prior to sequencing.
- Polymerase-bound libraries can withstand **>4 freeze-thaw cycles** without affecting sequencing performance.

• For complex cleanup reaction, each sample requires **100 μ L** of **SMRTbell cleanup beads**



• Each sample is eluted in **50 μ L** of **Loading Buffer 96**

• Perform **DNA concentration QC** on final purified polymerase-bound HiFi prep kit 96 SMRTbell library using a Qubit dsDNA HS assay or Quant-iT ds DNA HS assay

- Post-ABC recovery yield (relative to library input amount for primer annealing reaction) should be **~70%**

• **IMPORTANT!** The Qubit Flex instrument is not compatible with measuring polymerase-bound library in Loading Buffer 96 → Concentration readings will not be accurate.


• Prepare DNA internal control and perform final loading dilution procedure using **SMRT Link Loading Calculator** tool



SMRT Link Loading Calculator & Pooling Calculator for HiFi prep kit 96 libraries

Use **SMRT Link Revio polymerase kit 96 Loading Calculator** to calculate amount of polymerase-bound library to add to Revio sequencing plate and optionally use **SMRT Link Pooling Calculator**^{1,2} to calculate sample pooling volumes

Sample setup workflow for preparing Revio HiFi prep kit 96 libraries for sequencing



Procedure & checklist

Preparing whole genome libraries using the HiFi prep kit 96 (103-420-700)




ABC instructions for manual HiFi prep kit 96 workflows

Or

Guide & overview

Automated HiFi prep 96 and HiFi annealing binding and cleanup for the Hamilton NGS Microlab STAR system (103-425-700)




ABC instructions for automated HiFi prep kit 96 workflows



SMRT Link v13.1+

- Optional: Use **SMRT Link Pooling Calculator** to calculate required sample pooling volumes if performing multiplexed sequencing^{1,2}

- Use **SMRT Link Loading Calculator** to calculate final sample dilution when you are starting with **Revio polymerase kit 96-bound SMRTbell libraries** prepared using liquid handler automation or previously prepared manually



¹ If sequencing **HiFi plex prep kit 96 SMRTbell libraries**, you do not need to use this SMRT Link Pooling Calculator tool since samples are typically **already pooled** during SMRTbell library construction.
² If sequencing **HiFi prep kit 96 SMRTbell libraries**, this SMRT Link Pooling Calculator tool can optionally be used for pooling samples **either before** performing ABC **or after** performing ABC.

SMRT Link Pooling Calculator & Loading Calculator workflow overview

Pooling & Loading Calculator workflow overview for SMRTbell libraries bound with Revio polymerase kit 96

Use SMRT Link Revio polymerase kit 96 Loading Calculator & Pooling Calculator

Sample Setup

Sequel II binding kit 2.1/2.2
Sequel II binding kit 3.1/3.2, Revio polymerase kit
✓ Revio polymerase kit 96

Specify all required information in SMRT Link calculator tools

SMRT Link Pooling Calculator

Optional if using the Revio polymerase kit 96 **and** pooling multiple biological samples together in one well

Specify # of samples to be multiplexed, pooled library target volume, pooled library concentration and other parameters

Pooling Calculator

Number of samples to be multiplexed: 2

Pooled library target volume (µL): 50

Concentration output units: (ng/µL)

Pooled library concentration (ng/µL): 4

Buttons: Import, Export, Print

SMRT Link Loading Calculator

Specify the # of sample wells per Revio sequencing plate, followed by the concentration, average insert size & loading concentration of each sample

Revio sequencing plates

Plate 1: Number of wells: 4

Plate 2: Number of wells: 0

Plate 1

Well ID*	Sample name*	Conc. (ng/uL)*	Avg. insert size (bp)*	Loading conc. (pM)*	Comment
A01	Sample P1S1				
B01	Sample P1S2				
C01	Sample P1S3				
D01	Sample P1S4				

Follow provided instructions to process samples in the lab




Pool samples (OPTIONAL)

Prepare DNA internal control

Perform final sample dilution at target loading concentration

Load sample into Revio sequencing plate¹



SMRT Link Pooling Calculator procedure for Revio HiFi prep kit 96 libraries

Pooling Calculator procedure for performing equal-mass pooling of HiFi prep kit 96 libraries bound with Revio polymerase kit 96

1. Specify number of samples to be multiplexed

- Can specify a value between 2 and 384 samples

2. Specify pooled library target volume (μL)

- For post-ABC library samples, typically specify a pooled library target volume $\geq 50 \mu\text{L}$
- For pre-ABC library samples, typically specify a pooled library target volume $\geq 25 \mu\text{L}$

3. Specify concentration output units

- Specify $\text{ng}/\mu\text{L}$ to perform **equal-mass pooling**

4. Specify pooled library concentration

- Typically, specify a pooled library concentration value that is lower than the least concentrated sample

5. Specify sample concentrations

- Specify the concentration of each sample in $\text{ng}/\mu\text{L}$

6. Print or export instructions (optional)

- To print the calculation(s) and instructions, click the **Print** button.

Pooling Calculator

Number of samples to be multiplexed ? 3 1

Pooled library target volume (μL) 50 2

Concentration output units ($\text{ng}/\mu\text{L}$) \downarrow 3

Pooled library concentration ($\text{ng}/\mu\text{L}$) 5 4

Buffer volume (μL) to add to pooled sample: **3.48**

Import Export Print 6

Sample name	Conc. ($\text{ng}/\mu\text{L}$)	Pooling volume (μL)
Sample 1	6.2	13.44
Sample 2	4.8	17.36
Sample 3	5.3	15.72

5

Pooling Calculator
outputs required pooling
volume for each sample

SMRT Link Pooling Calculator procedure for Revio HiFi prep kit 96 libraries (cont.)

Pooling Calculator procedure for performing equal-molar pooling of HiFi prep kit 96 libraries bound with Revio polymerase kit 96

1. Specify number of samples to be multiplexed

- Can specify a value between 2 and 384 samples

2. Specify pooled library target volume (μL)

- For post-ABC library samples, typically specify a pooled library target volume ≥ 50 μL
- For pre-ABC library samples, typically specify a pooled library target volume ≥ 25 μL

3. Specify concentration output units

- Specify nM or pM to perform **equal-molar pooling**

4. Specify pooled library concentration

- Typically, specify a pooled library concentration value that is lower than the least concentrated sample

5. Specify sample concentrations

- Specify the concentration of each sample in nM or pM

6. Specify sample insert sizes

- Specify the insert size of each sample in bp

7. Print or export instructions (optional)

- To print the calculation(s) and instructions, click the [Print](#) button.

Pooling Calculator

Number of samples to be multiplexed 1

Pooled library target volume (μL) 2

Concentration output units (pM) 3

Pooled library concentration (pM) 4

Buffer volume (μL) to add to pooled sample: **12.87**

Library name	Conc. (ng/μL)	Insert size (bp)	Conc. (pM)	Pooling volume (μL)
Sample 1	6.2	15839	602	10.38
Sample 2	4.8	16025	461	13.56
Sample 3	5.3	17193	474	13.19

5

6

Pooling Calculator outputs required pooling volume for each sample

SMRT Link Loading Calculator procedure for Revio HiFi prep kit 96 libraries

Entering sample information

1. Specify number of samples to use

- Specify the number (1-4) of sample wells to use per Revio sequencing plate
- Note:** If you are using only one Revio sequencing plate, specify 0 for Plate 2

2. Enter information for first sample well

- Sample name
- Concentration (ng/μL)
- Average insert size (in base pairs)
- Loading concentration (in pM)
- Comments (optional)
- Note:** If using a partially-used Revio sequencing plate, can delete a Well ID by clicking on the 'x' button at right-hand side of table

3. Repeat Step 2 for additional sample wells

- Note:** All sample wells must be filled in for the instructions to display.

4. Print instructions (optional)

- To print the calculation(s) and instructions, click the [Print](#) button.

Sample Setup / Sample Calculation
Revio polymerase kit 96

Use this tool to calculate the final sample dilution for the Revio sequencing plate when you are starting with previously-prepared polymerase-bound SMRTbell libraries. Input the number of sample wells per Revio sequencing plate, followed by the concentration, average insert size, and loading concentration of each sample. The tool will return instructions for making the final dilution for each sample well.

Revo sequencing plates ⓘ

Plate 1 Number of wells 4

Plate 2 Number of wells 0

Plate 1

Well ID*	Sample name*	Conc. (ng/uL)*	Avg. insert size (bp)*	Loading conc. (pM)*	Comment
☐ A01	Sample P1S1	5.2	17362	225	HiFi prep kit 96 workflow
☐ B01	Sample P1S2	3.1	16934	225	HiFi prep kit 96 workflow
☐ C01	Sample P1S3	3.8	15904	225	HiFi prep kit 96 workflow
☐ D01	Sample P1S4	4.5	18003		

Pooling Calculator Print

1. Sequencing control dilution

2. Final loading dilution

3. Sample loading on sequencing plate

Loading Calculator outputs instructions for final loading dilution procedure

SMRT Link Loading Calculator procedure for Revio HiFi prep kit 96 libraries (cont.)

Calculator outputs and instructions

1. Sequencing control dilution

1. Sequencing control dilution

Perform three sequential dilution steps of **Sequencing control** using **Dilution buffer**. Use a new Lo-bind tube for each dilution step. The volume is sufficient for all samples in all batches.

Step	Instructions	Tube	Component	Volume
Perform a first dilution by adding Dilution buffer to Sequencing control.				
1.1		blue	Dilution buffer	19.0 uL
		red	Sequencing control	1.0 uL
			Total volume	20.0 uL
Pipette-mix the dilution. Quick-spin to collect contents. Keep on ice.				
Perform a second dilution by adding Dilution buffer to Dilution 1.				
1.3		blue	Dilution buffer	19.0 uL
			Dilution 1	1.0 uL
			Total volume	20.0 uL
Pipette-mix the dilution. Quick-spin to collect contents. Keep on ice.				
Perform a third dilution by adding Dilution buffer to Dilution 2.				
1.5		blue	Dilution buffer	19.0 uL
			Dilution 2	1.0 uL
			Total volume	20.0 uL
1.6	Pipette-mix the dilution. Quick-spin to collect contents.			
1.7	Discard Dilution 1 and Dilution 2.			
1.8	Keep Dilution 3 on ice.			

Three small-volume dilution steps are performed for better compatibility with automation systems



2. Final loading dilution

2. Final loading dilution

For each sample in each batch, add the following components in the order and volume listed to a new Lo-bind tube or tube strip.

Step	Tube	Plate:Well	Sample P1S1 1 : A01	Sample P1S2 1 : B01	Sample P1S3 1 : C01	Sample P1S4 1 : D01
2.1	green	Prepared Sample	48.8 uL	79.9 uL	61.2 uL	58.5 uL
		Loading buffer	48.2 uL	17.1 uL	35.5 uL	38.5 uL
		Diluted sequencing control (Dilution 3)	3.0 uL	3.0 uL	3.0 uL	3.0 uL
		Total volume	100.0 uL	100.0 uL	100.0 uL	100.0 uL
2.2	Pipette-mix each sample.					
2.3	Protect samples from light.					
2.4	Discard any unused Dilution 3.					

Add 3 µL (per Revio SMRT Cell) of serially diluted (1:8000) Revio sequencing control to final on-plate sample dilution

SAFE STOPPING POINT - Store protected from light at 4°C for up to 24 hours.

Final loading dilution total volume (per Revio SMRT Cell) is 100 µL

Final diluted samples should be used within 24 hrs




3. Sample loading on sequencing plate

3. Sample loading on sequencing plate

Step	Instructions
3.1	Load 95 uL of sample per well on a Revio sequencing plate.

Load 95 µL of final diluted sample per well on a Revio sequencing plate.





HiFi prep kit 96 example sequencing performance data

HiFi sequencing performance of HiFi prep kit 96 libraries

Example HiFi read length and yield performance obtained for HiFi prep kit 96 libraries constructed from human HG002 genomic DNA samples¹

Plate A 96 pool




Plate B 96 pool




Plate C 96 pool





	Plate A	Plate B	Plate C
SMRTbell adapter index plate	96A	96B	96C
HiFi yield	111 Gb	121 Gb	117 Gb
HiFi read length	15 kb	16 kb	15 kb
Average Revio SMRT Cells/2 µg gDNA	2.1 ± 0.4	1.9 ± 0.3	1.8 ± 0.3
Automation workflow time (gDNA → ABC)	12.4 hrs	11.6 hrs	12 hrs

¹ Revio system data shown. **Note:** HiFi read lengths, reads/data per SMRT Cell and other sequencing performance results can vary depending on DNA sample quality, insert size, *P1* loading performance & movie time. Shorter library insert sizes (<15 kb), lower DNA quality samples, and suboptimal *P1* loading performance may result in HiFi data yields <90 Gb per Revio SMRT Cell.


HiFi sequencing performance of HiFi prep kit 96 libraries

Example HiFi read length and yield performance obtained for HiFi prep kit 96 libraries constructed from gDNA isolated from human and animal whole blood samples using Nanobind PanDNA kit¹

Human and animal whole blood samples



Tuna whole blood sample





HiFi plex prep kit 96 workflow overview for WGS and amplicon sequencing applications

Procedure & checklist – Preparing multiplexed whole genome and amplicon libraries using the HiFi plex prep kit 96 (103-418-800)

Procedure & checklist 103-418-800 describes the workflow for constructing multiplexed whole genome sequencing (WGS) and amplicon libraries using the HiFi Plex Prep Kit 96 for sequencing on PacBio systems. Automation is highly recommended for the first half of the protocol up to and including pooling; however, the entire procedure may be performed manually. Please see instrument-specific protocols for automation details¹.

Overview		
Applications	<ul style="list-style-type: none"> • Microbial WGS • Metagenome shotgun sequencing • Amplicon sequencing • Low-pass WGS 	
Samples	24–96 per kit	
Minimum batch size supported	4 x 24	
Maximum batch size supported	96	
DNA input and fragment size recommendations		
	gDNA	Amplicons
Per sample input	50–300 ng	20–200 ng
DNA shearing	Automated pipette-tip shearing	N/A
Target fragment lengths	13–20 kb*	Any size >1 kb

*Smaller fragment lengths can be used when working with lower quality DNA samples (see gDNA quality recommendations).

Note: To prepare HiFi plex prep kit 96 samples for sequencing on PacBio long-read systems, follow ABC² workflow instructions provided in **SMRT Link Sample Setup** calculator tool

¹ If using the HiFi plex prep kit 96 in an **automated** high-throughput library prep workflow with a Hamilton liquid handling system, refer to **Guide & overview – Automated HiFi plex prep 96 for the Hamilton NGS Microlab STAR system** (103-425-800).

² If using **Revio polymerase kit (102-817-600)** or **Sequel II binding kit 3.2 (102-333-300)** for polymerase binding, follow sample setup ABC instructions provided in **SMRT Link Sample Setup calculator tool for Sequel II binding kit 3.1/3.2, Revio polymerase kit**.to perform primer annealing, polymerase binding & complex cleanup (ABC) steps and final loading dilution procedure to prepare HiFi plex prep kit 96 libraries for sequencing (do not use SMRT Link Loading Calculator).

PacBio

Preparing multiplexed whole genome and amplicon libraries using the HiFi plex prep kit 96

Procedure & checklist

Overview

This procedure describes the workflow for constructing multiplexed whole genome sequencing (WGS) and amplicon libraries using the HiFi Plex Prep Kit 96 for sequencing on PacBio® systems. Automation is highly recommended for the first half of the protocol, prior to pooling; however, the entire procedure may be performed manually. Please see instrument-specific protocols for automation details.

Overview		
Applications	<ul style="list-style-type: none"> • Microbial WGS • Metagenome shotgun sequencing • Amplicon sequencing • Low-pass WGS 	
Samples	24–96 per kit	
Minimum batch size supported	4 x 24	
Maximum batch size supported	96	
	Hamilton NGS STAR	Hamilton Microlab Prep
Shearing	10 min for 24–96 samples	22 min for 24 samples
Post-shearing cleanup	1 hour for 96 samples	1.5 hours for 24 samples (Post-shearing cleanup only)
Automated steps (from ER/DDR to pooling)	4 hours for 96 samples	N/A
Manual steps (from post-ligation cleanup to post-nuclease cleanup)	1.5 hours for 96 samples (4 pools of 24 samples)	N/A
Average total time	7 hours	~2 hours for shearing and cleanup only
	gDNA	Amplicons
Per sample input	50–300 ng	20–200 ng
DNA shearing	Automated pipette-tip shearing	N/A
Target fragment lengths	13–20 kb*	Any size >1 kb

*Smaller fragment lengths can be used when working with lower quality DNA samples (see gDNA quality recommendations).

© 2024 PacBio. All rights reserved. Research use only. Not for use in diagnostic procedures.
103-418-800 REV01 MAR2024


PacBio

Guide & overview – Automated HiFi plex prep 96 for the Hamilton NGS Microlab STAR system (103-425-800)¹

Guide & overview 103-425-800 describes the workflow for PacBio high-throughput semi-automated long read sample prep utilizing the HiFi Prep Plex 96 kit that is designed for a minimum of 24 and maximum of 96 samples per automated run.

Overview	
Applications	<ul style="list-style-type: none">• Microbial WGS• Metagenome shotgun sequencing• Amplicon sequencing• Long read low pass WGS
Samples	24 – 96 per kit
Minimum batch size supported	4 x 24
Maximum batch size supported	96
	Hamilton NGS STAR
Shearing	10 min for 24–96 samples
Post-shearing cleanup	1 hour for 96 samples
Automated steps (ER/DDR to pooling)	4 hours for 96 samples
Manual steps (Post-ligation cleanup to post-nuclease cleanup)	1.5 hours for 96 samples (4 pools of 24 samples)
Average total time	7 hours

Note: To prepare HiFi plex prep kit 96 samples for sequencing on PacBio long-read systems, follow ABC² workflow instructions provided in **SMRT Link Sample Setup** calculator tool



¹ If using the HiFi plex prep kit 96 in an **automated** high-throughput library prep workflow with a liquid handling system, refer to **Guide & overview – Automated HiFi plex prep 96 for the Hamilton NGS Microlab STAR system** (103-425-800).

² If using **Revio polymerase kit (102-817-600)** or **Sequel II binding kit 3.2 (102-333-300)** for polymerase binding, follow sample setup ABC instructions provided in **SMRT Link Sample Setup calculator tool for Sequel II binding kit 3.1/3.2, Revio polymerase kit**.to perform primer annealing, polymerase binding & complex cleanup (ABC) steps and final loading dilution procedure to prepare HiFi plex prep kit 96 libraries for sequencing (do not use SMRT Link Loading Calculator).

HiFi plex prep kit 96 workflow overview for multiplexed WGS & amplicon applications

Automation of HiFi library prep workflow enables high-throughput processing of 24 to 96 samples in 2 days

Key workflow steps	Automation	Time	Protocol documentation or reference	Recommended equipment & consumables
DNA sample extraction		2 – 2.5 hrs	Guide & overview – Nanobind HT kits (103-028-100) [UPDATED]	Nanobind HT CBB kit / Nanobind HT 1 mL whole blood kit
DNA sample QC		2 hrs	Refer to third-party user guide documentation	DNA concentration QC: • Varioskan LUX microplate reader DNA sizing QC: • Femto Pulse system
DAY 1				
1 Pipette DNA shearing		10 min	Protocol for manual & automated workflows	Pipette-based shearing with Hamilton Microlab Prep or NGS STAR
↓ ○ Safe stop				
2 Post-shearing cleanup ¹		1 hr	Procedure & checklist – Preparing multiplexed whole genome and amplicon libraries using the HiFi plex prep kit 96 (103-418-800)	
↓ ○ Safe stop				
3 DNA repair & A-tailing		0.5 hrs	Guide & overview for automated workflows	
↓				
4 Adapter ligation, termination & pooling		3.5 hrs	Guide & overview – Automated HiFi plex prep 96 for the Hamilton NGS Microlab STAR system (103-425-800)	
↓				
5 Pool cleanup		45 min	Supported automation systems Thermo Fisher KingFisher Apex Hamilton Microlab NIMBUS Presto Hamilton Microlab Prep Hamilton NGS Star	HiFi plex prep kit 96 (103-381-300 / 96 RXN) + SMRTbell adapter index plate 96A/B/C/D with Hamilton
↓				
6 Nuclease treatment		15 min		
↓ ○ Safe stop				
7 Final cleanup		45 min		1X SMRTbell bead cleanup

Library construction

General best practices recommendations for preparing multiplexed WGS and amplicon libraries using HiFi plex prep kit 96

Automation

- The HiFi plex prep kit 96 and workflow was **designed to work with NGS liquid handling automation up to the sample pooling step**
 - Once samples are pooled, the **remaining steps are prepared manually**
- **Note:** This protocol was developed using the **Hamilton NGS STAR MOA 96 system**
 - Refer to ***Guide & overview – Automated HiFi plex prep 96 for the Hamilton NGS Microlab STAR system*** (103-425-800)
- Because of differences between automation instruments, modifications not described herein may be needed to adapt the protocol to your specific instrumentation
 - Please visit www.pacb.com or contact your local support team for a list of instruments with a PacBio qualified method



Hamilton NGS Star workstation

General best practices recommendations for preparing multiplexed WGS and amplicon libraries using HiFi plex prep kit 96 (cont.)

DNA sample extraction

Genomic DNA extraction from cultured bacteria

PacBio Nanobind DNA extractions kits are recommended to ensure sufficient mass and quality of high-molecular weight (HMW) DNA for use in HiFi plex prep kit 96 protocol

Nanobind HT CBB kit (102-762-700; 96 rxn)	Nanobind PanDNA kit (103-260-000; 24 rxn)
<ul style="list-style-type: none">For high-throughput HMW DNA extraction from up to 200 µL human/mammalian blood, non-mammalian animal blood¹, cultured cells, and bacteriaExpected HMW DNA yield: 3–15 µg for blood and cultured mammalian cells and 2–10 µg for bacteria	<ul style="list-style-type: none">For HMW DNA extraction from cells, blood, bacteria, tissues, insects, and plant nucleiExpected HMW DNA yield: 3–26 µg

Genomic DNA extraction from metagenomic samples

Note: The products below have not been tested or validated by PacBio but are listed here as examples of third-party kits used by other PacBio customers for isolating genomic DNA for PacBio metagenomic sequencing applications

Sample type	Third-party product or kit
Fecal and soil	QIAGEN DNeasy PowerSoil Pro (PN 47014)
	QIAGEN PowerFecal Pro (PN 51804)
	QIAGEN DNeasy PowerClean Pro Cleanup Kit (PN 12997-50) <ul style="list-style-type: none">If needed, can be used after extracting DNA with PowerSoil or PowerFecal kits to further improve sequencing performance
Saliva	DNA Genotek Oragene OG 500 collection tubes (PN OG-500) <ul style="list-style-type: none">Recommended for collection of saliva samples

Amplicon DNA generation

- Note:** Using gel-extracted amplicon products may result in lower sequencing performance due to the damage inherently caused by intercalating dyes such as ethidium bromide and exposure to UV radiation.
- Sequencing amplicons stained with SYBR dyes from ThermoFisher Scientific is untested, and therefore cannot be recommended.
- If working with a gel-extracted product that has been stained with a dye, it is recommended to bring it through additional rounds of amplification to remove damage and/or dyes prior to library prep and sequencing.

General best practices recommendations for preparing multiplexed WGS and amplicon libraries using HiFi plex prep kit 96 (cont.)

DNA sizing QC

- Agilent **Femto Pulse system**¹ is highly recommended for the accurate sizing of genomic DNA samples
- Femto Pulse system enables simple, rapid sizing QC of genomic DNA and SMRTbell libraries, and conserves sample by using femtogram ranges of input DNA
 - Resolves fragments 1,300 bp to 165 kb using gDNA 165 kb Analysis kit (can resolve 100 – 6,000 bp using Ultra Sensitivity NGS kit)
 - Requires <1 ng of sample DNA
 - Can analyze up to 12 samples in <1.5 hrs
 - Outputs quality metrics such as Genomic Quality Number (GQN)² to quickly score integrity of HMW gDNA



Femto Pulse system
(Agilent Technologies)

DNA quantification QC

- For high-throughput DNA quantification QC workflows, we recommend using the Quant-iT 1X dsDNA high sensitivity assay kit³ (Thermo Fisher Scientific) with the Varioskan LUX multimode microplate reader (Thermo Fisher Scientific)
 - Varioskan LUX multimode microplate reader enables rapid, specific and accurate determination of nucleic acid concentrations in a wide range
 - Quant-iT assay is well-adapted to high-throughput use (typically in 96-well or 384-well plates) and is highly selective for dsDNA over RNA (in the 0.2–100 ng range, the fluorescence signal is linear with the amount of DNA)
 - **Note:** We do not recommend quantification with UV-Vis Spectrophotometers (e.g. NanoDrop) that measure all nucleic acids in a sample. For example, measuring all nucleic acid will inflate the true concentration of gDNA in samples



Varioskan LUX multimode
microplate reader
(Thermo Fisher Scientific)

¹ See *Product Note – HiFi WGS sequencing with the Agilent Femto Pulse system* ([102-326-561](#)) for more details.

² See *Application Note – Quality Metrics for Nucleic Acids with the Agilent Fragment Analyzer and Femto Pulse Systems* (Agilent [5994-0521EN](#))

³ Alternatively, for lower-throughput applications DNA quantification QC may be performed without a microplate reader using the **Qubit 1X dsDNA high sensitivity assay kit** with a **Qubit fluorometer**.

General best practices recommendations for preparing multiplexed WGS and amplicon libraries using HiFi plex prep kit 96

Recommended DNA input amount and quality

- Recommended DNA input amounts will vary by starting DNA sample type

Genomic DNA	Amplicon DNA
50–300 ng per sample	20–200 ng per sample

IMPORTANT!

- Do not exceed >300 ng of gDNA per sample going into the Repair and A-tailing steps.
→ Too much DNA may overwhelm enzymatic reactions and lead to poor library recovery.
- At least 24 samples need to be used when using lower input amounts (e.g., 50 ng).
- Samples must be pooled after ligation


- For genomic DNA samples, 70% or more of the DNA should be ≥ 10 kb for this protocol
 - This corresponds to a genome quality number (GQN) of 7.0 or higher at 10 kb
- Remove all RNA from genomic DNA samples prior to beginning
 - Residual RNA can inhibit sequencing polymerase binding and can therefore lead to low loading

General best practices recommendations for preparing multiplexed WGS and amplicon libraries using HiFi plex prep kit 96


DNA shearing

- This protocol recommends shearing genomic DNA using automated liquid handler systems to a size between **13–20 kb**
 - If shearing DNA samples using a **Hamilton Microlab Prep** liquid handling instrument, refer to **Technical note – High throughput DNA shearing using Hamilton Microlab Prep** (102-326-606) or **Guide and overview – Short Read Eliminator (SRE), DNA shearing, and cleanup for the Hamilton Microlab Prep system** (103-424-100) for details about third-party consumables requirements
 - If shearing DNA samples using a **Hamilton NGS Star MOA** liquid handling workstation, refer to **Guide & overview – Automated HiFi plex prep 96 for the Hamilton NGS Microlab STAR system** (103-425-800) for details about third-party consumables requirements
- **Microbial and metagenomic samples often have degraded DNA** where the majority of fragments are **already <15 kb** in length
 - To better balance the number of reads between samples it may be necessary to shear all samples to a mean size of **~10 kb**
- For shearing DNA samples to a target fragment size below **<15 kb** for PacBio HiFi sequencing, we recommend using the following equipment options:

SPEX SamplePrep 1600 MiniG homogenizer ¹	
Parameter	Settings
Sheared fragment size	7–10 kb
Speed	1500 RPM
Time	3 minutes
DNA input	300 ng – 3 µg
Volume	Up to 300 µL



MP Bio FastPrep 96 homogenizer ²	
Parameter	Settings
Sheared fragment size	7–10 kb
Speed	1800 RPM
Time	60 seconds
DNA input	300 ng – 3 µg
Volume	50 µL



¹ See **Technical note – HT plate shearing microbial inserts** ([102-326-575](#)) for description of SamplePrep 1600 MiniG system for PacBio workflows.






² See **Technical note – HT DNA shearing for HiFi whole genome sequencing from whole blood samples** ([102-326-579](#)) for description of FastPrep 96 system for PacBio workflows.

General best practices recommendations for preparing multiplexed WGS and amplicon libraries using HiFi plex prep kit 96

Reagent handling

- Room temperature is defined as any temperature in the range of 18 – 25°C for this protocol

HiFi plex prep kit 96

Thaw these reagents at room temperature	Keep these temperature-sensitive reagents on ice	Bring these reagents to room temperature 30 minutes prior to use	Bring these reagents to room temperature 1.5 hours prior to use
<input type="checkbox"/> Repair buffer M96	 <input type="checkbox"/> End repair mix M96	<input type="checkbox"/> 1X dsDNA HS assay reagents	<input type="checkbox"/> SMRTbell cleanup beads-85 mL ¹
<input type="checkbox"/> Nuclease buffer M96	 <input type="checkbox"/> DNA repair mix M96		
<input type="checkbox"/> SMRTbell adapter index plate	 <input type="checkbox"/> Ligation mix M96		
<input type="checkbox"/> Stop solution M96	 <input type="checkbox"/> Ligation enhancer M96		
	 <input type="checkbox"/> Nuclease mix M96		


¹ Alternatively, can bring SMRTbell cleanup beads-85 mL to room temperature the night before if starting protocol in the morning.

- Once thawed, reaction buffers and adapter index plate may be stored on ice
- Briefly spin down all reagent tubes to collect all liquid at bottom
- Briefly vortex then spin down SMRTbell adapter index plate in a centrifuge with a plate adapter to collect all liquid at the bottom of the wells
- Shake/vortex SMRTbell cleanup beads immediately before use
- Pipette-mix bead binding and elution steps until beads are distributed evenly
- Pipette-mix all SMRTbell prep reactions by pipetting up and down 10 times.
- Samples can be stored at 4°C at all safe stopping points listed in the protocol.

General best practices recommendations for preparing multiplexed WGS and amplicon libraries using HiFi plex prep kit 96

Multiplexing samples


- **Note:** The procedure requires any of the four available SMRTbell adapter index plates:
 - ❑ SMRTbell adapter index plate 96A (bc2001–bc2096)
 - ❑ SMRTbell adapter index plate 96B (bc2097–bc2192)
 - ❑ SMRTbell adapter index plate 96C (bc2193–bc2288)
 - ❑ SMRTbell adapter index plate 96D (bc2289–bc2384)
- To balance the number of reads per sample, consider the following:
 - Shear all gDNA samples to **similar mean fragment sizes and distributions**
 - **Normalize DNA input** across all samples
- To pool more than 96 samples per SMRT Cell, combine indexed libraries constructed using multiple HiFi plex prep kits following the final cleanup
 - Each prep kit needs to use a different SMRTbell adapter index plate
- A total of 384 samples can be pooled for sequencing on a single Revio SMRT Cell¹
 - A total of 1,536 (4 cells x 384 samples) can be sequenced in a single Revio system 4-cell run



General best practices recommendations for preparing multiplexed WGS and amplicon libraries using HiFi plex prep kit 96

Sequencing preparation (ABC) and polymerase-bound library storage

- **Note: Procedure & checklist** – *Preparing multiplexed whole genome and amplicon libraries using the HiFi plex prep kit 96* (103-418-800) does not include instructions for the primer annealing, polymerase binding, and complex cleanup (ABC) sample setup steps
 - Follow **SMRT Link Sample Setup Sequel II binding kit 3.1/3.2**, **Revio polymerase kit** calculator tool to perform ABC and final loading dilution procedure for HiFi plex prep kit 96 libraries
- Sequencing polymerase is stable once bound to the HiFi prep kit library and can be stored at 4°C or frozen at -20°C.
- Stored polymerase-bound library shows equivalent loading to freshly prepared bound libraries up to the recommendations listed below.
 - Recommended polymerase-bound storage:¹**
 - Polymerase-bound library is stable at 4°C for 1 month
 - Frozen polymerase-bound library is stable for at least 6 months
- Please note that the stored polymerase-bound library needs to be **protected from light** while stored.



HiFi plex prep kit 96 library prep inputs & expected step yields

Final HiFi library yield is typically sufficient to load ≥ 1 SMRT Cells


	Cleanup	DNA QC	Step input	Step output / Yield (%)
DNA sample extraction	N/A	N/A	200 μ L human or animal whole blood / 1 mL human whole blood / 1×10^6 cells cultured mammalian cells	Nanobind HT CBB kit : 3 – 15 μ g HMW DNA Nanobind HT 1 mL blood kit 3 – 70 μ g HMW DNA
DNA sample QC	N/A	Quant-iT/Qubit dsDNA HS assay Femto Pulse system	Nanobind HT-extracted HMW DNA	DNA quantification QC \rightarrow Aim for ≥ 3 μ g HMW DNA DNA sizing QC \rightarrow Ideally GQN(10 kb) ≥ 7.0
1 Pipette DNA shearing	N/A	Post-shearing QC with Quant-iT/Qubit (optional) & Femto Pulse	≤ 10 ng/ μ L DNA in 300 μ L (3 μ g total input DNA mass) into Step 2	Expect up to 80% step recovery (60% total recovery) Target DNA shear size is ~15-20 kb
2 Post-shearing cleanup	1X SMRTbell cleanup beads	N/A	≤ 300 μ L sheared DNA	4.5 μ L purified sheared DNA
3 DNA repair & A-tailing	N/A	N/A	4.5 μ L purified sheared DNA	30 μ L repaired & A-tailed DNA
4 Adapter ligation & pooling	N/A	N/A	30 μ L post-repaired & A-tailed DNA	960 μ L for 24-plex of pooled libraries (24 x 40 μ L per sample)
5 Pool cleanup	1X SMRTbell cleanup beads	N/A	960 μ L for 24-plex of pooled libraries	40 μ L of purified pooled libraries
6 Nuclease treatment	N/A	N/A	40 μ L of purified pooled libraries	50 μ L of nuclease-treated sample
7 Final cleanup	1X (for >3 kb) or 1.3X (for <3 kb) SMRTbell cleanup beads	Quant-iT/Qubit dsDNA HS assay Femto Pulse system	50 μ L of nuclease-treated sample	Final DNA concentration must be <60 ng/ μ L for >10 kb libraries, <20 ng/ μ L for 3-10 kb libraries and <10 ng/ μ L for <3 kb libraries to proceed with ABC ¹

DNA shearing

Library construction

DNA sample extraction

Perform automated high-throughput HMW DNA extraction using Nanobind HT kits¹



Available Nanobind HT kits



Automated solutions built on Nanobind technology

Nanobind HT kits use **magnetic disk processing** to automate lysis, binding, washing and elution steps and are compatible with instruments from Hamilton and Thermo Fisher

- Hamilton NIMBUS Presto is a walkaway solution with **automated** plate filling
- Thermo Fisher KingFisher instruments are **semi-automated** with manual plate filling and limited user interaction during the run

Nanobind HT CBB kit (102-762-700; 96 rxn)

- For up to 200 μ L human/mammalian blood, non-mammalian animal blood², cultured cells, and bacteria
- Expected HMW DNA yield: 3–15 μ g for blood and cultured mammalian cells and 2–10 μ g for bacteria

Nanobind HT 1 mL blood kit (102-762-800; 96 rxn)


- For 1 mL human blood
- Expected HMW DNA yield: 3–70 μ g

¹ See *Technical overview — Automated high-throughput HMW DNA extraction for PacBio long-read sequencing using Nanobind HT kits* ([103-401-700](#)).

² Note: For animals with nucleated red blood cells, only 2.5 – 20 μ L of blood is needed for HMW DNA extraction using Nanobind HT CBB kit.

DNA sample extraction (cont.)


Available Nanobind HT procedures for high-throughput HMW DNA extraction from bacteria



Automation system	Procedure & checklist	Sample type
KingFisher Apex	Extracting HMW DNA using the Nanobind HT CBB kit for mammalian cultured cells on KingFisher Apex system [102-996-100]	Mammalian cells
	Extracting HMW DNA using the Nanobind HT CBB kit for 200 μ L human whole blood on KingFisher Apex system [102-995-700]	Human whole blood
	Extracting HMW DNA using the Nanobind HT CBB kit for bacteria on the KingFisher Apex system [103-377-600]	Cultured bacteria
	Extracting HMW DNA using the Nanobind HT CBB kit for non-human mammalian blood (NHMB) on the KingFisher Apex system [103-397-300]	Non-human mammalian blood
	Extracting HMW DNA using the Nanobind HT CBB kit for nucleated red blood cells (nRBCs) on the KingFisher Apex system [103-377-800]	Non-mammalian blood (nucleated RBCs)
	Extracting HMW DNA using Nanobind HT 1 mL blood kit for human whole blood on KingFisher Apex system [102-995-300]	Human whole blood
Hamilton NIMBUS Presto	Extracting HMW DNA using the Nanobind HT CBB kit for mammalian cultured cells on Hamilton NIMBUS Presto system [102-996-400]	Mammalian cells
	Extracting HMW DNA using the Nanobind HT CBB kit for 200 μ L human whole blood on Hamilton NIMBUS Presto system [102-996-000]	Human whole blood
	Extracting HMW DNA using the Nanobind HT CBB kit for bacteria on Hamilton NIMBUS Presto system [103-397-400]	Cultured bacteria
	Extracting HMW DNA using the Nanobind HT CBB kit for non-human mammalian blood (NHMB) on the Hamilton NIMBUS Presto system [103-377-700]	Non-human mammalian blood
	Extracting HMW DNA using the Nanobind HT CBB kit for nucleated red blood cells (nRBCs) on the Hamilton NIMBUS Presto system [103-397-500]	Non-mammalian blood (nucleated RBCs)
	Extracting HMW DNA using the Nanobind HT 1 mL blood kit for human whole blood on Hamilton NIMBUS Presto system [102-995-600]	Human whole blood

DNA sample QC

Perform high-throughput DNA QC using a microplate reader¹ and a Femto Pulse system



DNA quantification QC



Varioskan LUX multimode microplate reader
(Thermo Fisher Scientific)

Varioskan LUX multimode microplate reader enables rapid, specific and accurate determination of nucleic acid concentrations in high-throughput workflows¹

- Perform DNA concentration measurements using **Quant-iT dsDNA assay²** (Thermo Fisher Scientific)
- Quant-iT assay is well-adapted to high-throughput use, typically in **96-well** or **384-well plates**
- Quant-iT assay is highly selective for dsDNA over RNA, and in the range of **0.2–100 ng**, the fluorescence signal is linear with the amount of DNA

¹ Alternatively, for lower-throughput applications DNA quantification QC may be performed without a microplate reader using the Qubit 1X dsDNA high sensitivity assay kit.

DNA sizing QC




Femto Pulse system
(Agilent Technologies)

Femto Pulse system enables simple, rapid sizing QC of genomic DNA and SMRTbell libraries, and conserves sample by using femtogram ranges of input DNA

- Use the **Femto Pulse gDNA 165 kb analysis kit** (FP-1002-0275)
- Dilute samples to **250 pg/uL**
- **70% or more of the DNA should be ≥10 kb** for optimal results. This corresponds to a genome quality number (GQN)³ of 7.0 or higher at 10 kb.

Pipette DNA shearing

Perform automated high-throughput shearing for HMW gDNA samples using Hamilton automation^{1,2}



1. Pipette DNA shearing

Step	Instructions
1.1	Use Buffer LTE HT to bring all samples up to 300 μL total volume in a 0.8 mL, 96 DeepWell plate (Thermo Fisher Scientific AB0859). The concentration of each sample must be <10 ng/μL.
Parameters for shearing on the Microlab Prep, or Hamilton assay-ready workstations. These parameters should already be part of the installed method on the instrument.	
Parameter	Setting
1.2 DNA concentration	≥ 10 ng/ μL
1.2 Volume of Buffer LTE	300 μL
1.2 Number of mixes	300 cycles
1.2 Pipette mixing speed	500 $\mu\text{L}/\text{sec}$
1.2 Liquid following	83% volume
1.2 Pipette tip	300 μL CO-RE II tips (filtered, black, non-sterile)
1.3	Place the plate on the appropriate work deck position and start the shearing procedure.
Optional: measure DNA concentration with a Qubit fluorometer using the 1x dsDNA HS kit once shearing procedure is complete.	
1.4	Recommended: Further dilute each aliquot to 250 pg/ μL with Femto Pulse dilution buffer. Measure the final SMRTbell library size distribution with a Femto Pulse system to ensure efficient shearing.
1.5	Proceed to the 1X SMRTbell cleanup bead procedure to concentrate samples for library preparation.

- For automated pipette-based shearing using Hamilton systems, DNA samples should be at **<10 ng/ μL** in a total volume of **300 μL** in a 0.8 mL, 96 DeepWell plate (use Buffer LTE to dilute samples)



Hamilton NGS
STAR/STARlet/STAR V



Hamilton
Microlab Prep

Pipette-based DNA shearing can be completed within ~10 min for up to 96 samples using Hamilton liquid handling instruments.

IMPORTANT!


- Please follow the exact shearing conditions outlined below for the automated DNA shearing method using Hamilton automation.
- Deviating from these conditions may result in undersheared or unsheared DNA.

¹ This section describes the procedure for shearing HMW gDNA samples with the Hamilton Microlab Prep or Hamilton assay ready workstations (NGS STAR MOA, STARlet, and STAR V). It may be possible to shear DNA using other NGS liquid handler systems – Please [contact](#) PacBio Technical Support for updated information on all qualified DNA shearing methods.

² **Note:** If using a SPEX SamplePrep 1600 MiniG homogenizer or MP Bio FastPrep 96 homogenizer to shear lower-quality gDNA samples to a target fragment size <15 kb, follow the alternative equipment settings provided in *Procedure & checklist – Preparing multiplexed whole genome and amplicon libraries using the HiFi plex prep kit 96* (103-418-800).

Post-shearing cleanup

Perform post-shearing cleanup using 1X SMRTbell cleanup beads



2. Post-shearing cleanup

Step	Instructions for SMRTbell cleanup bead step
2.1	Sheared genomic DNA: add 300 μL (1.0X) of resuspended, room-temperature SMRTbell cleanup beads to each sample. Amplicon DNA: if clean-up or concentration is required for amplicons, add 1.0X (amplicons >3 kb) or 1.3X (amplicons <3 kb) (v/v) concentration of SMRTbell cleanup beads.
2.2	Pipette-mix the sample until the beads are evenly distributed. If necessary, quick-spin the samples to collect liquid.
2.3	Leave at room temperature for 10 minutes to allow DNA to bind beads.
2.4	Place samples on a magnet and allow beads to separate fully from the solution. The solution should be clear and beads pelleted to the magnet.
2.5	Slowly pipette off the supernatant without disturbing the beads. Discard the supernatant.
2.6	Slowly dispense 200 μL , or enough to cover the beads, of freshly prepared 80% ethanol to each sample. After 30 seconds , pipette off the 80% ethanol and discard.
2.7	Repeat the previous step.
2.8	Remove residual 80% ethanol: <ul style="list-style-type: none">Remove the samples from the magnet and quick-spin to collect liquid.Place samples back on the magnet and wait until beads separate fully from the solution.Carefully pipette-off the residual 80% ethanol without disturbing the bead pellet and discard.
2.9	Remove samples from the magnet and immediately add 25.5 μL of elution buffer . Resuspend by pipette mixing until beads are evenly distributed in the solution. Quick-spin samples if necessary to collect liquid.
2.10	Leave samples at room temperature for 5 minutes to elute DNA off beads.
2.11	Place samples back on the magnet and allow beads to separate fully from the solution. The solution should be clear, and beads pelleted to the magnet before proceeding.
2.12	Slowly pipette off the cleared supernatant without disturbing the beads. Transfer supernatant to a new plate or tube strip. Discard old plate or tube with beads.


SAFE STOPPING POINT - Store at 4°C

• **IMPORTANT!** Allow SMRTbell cleanup beads to come up to room temperature by bringing them out of 4°C storage at least 1.5 hrs prior to beginning

Note: Post-shear cleanup step should be performed on the **same day** as the library construction steps to ensure that there is a sufficient volume of SMRTbell cleanup beads to complete the entire library prep workflow¹

DNA repair & A-tailing

Repair sites of DNA damage and prepare sheared DNA for ligation to SMRTbell adapter



3. DNA repair & A-tailing

3.1 Prepare the appropriate volume of master mix with 15% overage using the per reaction volumes listed below.

Repair mix		
Tube	Component	Volume
Purple	Repair buffer M96	4 μ L
Blue	End repair mix M96	1 μ L
Green	DNA repair mix M96	0.5 μ L
Total volume		5.5 μ L

- Prepare a reaction master mix by adding the required components in the order and volume listed to a new microcentrifuge tube¹
- Adjust component volumes for the number of samples being prepared (e.g., 24-plex, 48-plex, 72-plex or 96-plex¹), plus 15% overage

3.2 Slowly pipette-mix the **repair mix** and quick spin to collect liquid at the bottom of the tube. If bubbles form during mixing, pulse-spin to remove bubbles.

3.3 Add **5.5 μ L** of the repair mix to each sample. The total reaction volume should be **30 μ L**.

3.4 Pipette-mix the reactions and quick-spin to collect liquid at the bottom of the well.

3.5 Run the **Repair and A-tailing** thermocycler program. Set lid temperature to 75°C if programmable.


Step	Time	Temperature
1	30 min	37°C
2	5 min	65°C
3	Hold	4°C

- Run **Repair and A-tailing** on-deck thermal cycler program
- Set the lid temperature to 75°C

3.6 Proceed to the next step of the protocol.

Adapter ligation, termination & pooling

In this step, ligate SMRTbell indexed adapter to the ends of each DNA fragment, stop the ligation reaction and then pool indexed samples together



4. Adapter ligation, termination & pooling

✓ Step Instructions for adapter index ligation

Add **4 μL** of indexed adapter to each sample from the previous step.

4.1

Component	Volume
SMRTbell adapter index plate 96 (A, B, C, or D)	4 μL

4.2 Prepare the appropriate volume of master mix using the per reaction volumes listed below. Add 10% overage if preparing the reactions manually.

Tube	Component	Volume
Yellow	Ligation mix M96	10 μL
Red	Ligation enhancer M96	0.5 μL
Total volume		10.5 μL

4.3 Pipette-mix the ligation mix and quick-spin to collect liquid.

4.4 Add **10.5 μL** of ligation mix to each sample from previous step. The total volume should be **44.5 μL** .

4.5 Pipette-mix each sample thoroughly to ensure the mix is evenly distributed.

4.6 Quick-spin to collect liquid.

Run the **Adapter ligation** thermocycler program. Set the lid temperature to 75°C if programmable.

Step	Time	Temperature
4.7 1	30 min	20°C
2	Hold	4°C

Proceed to the next step to terminate the ligation reaction prior to pooling.

4.8 Terminate the ligation reaction by adding **15.5 μL** of **Stop Solution M96** to each sample. Pipette mix thoroughly.

4.9 Pool **40 μL** from up to 24 reactions into a **2 mL LoBind tube**. The total pooled volume for 24 samples should be **960 μL** . If processing 96 samples, use four 2 mL LoBind tubes.

4.10 Proceed to the next step of the protocol.

- Add 4 μL of indexed adapter to each sample from the previous step.
- Any of the adapters from the four SMRTbell adapter index plates can be used (e.g. 96A, 96B, 96C, or 96D)

- Prepare a reaction master mix by adding the required components in the order and volume listed to a new microcentrifuge tube¹
- Adjust component volumes for the number of samples being prepared (e.g., 24-plex, 48-plex, 72-plex or 96-plex¹), plus 10% overage


- Run **Adapter ligation** on-deck thermal cycler program
- Set the lid temperature to 75°C

- **Terminate** ligation reaction by adding 15.5 μL of Stop Solution M96 to each sample and mixing thoroughly

- **Pool** 40 μL from up to 24 reactions into a 2 mL LoBind tube (total pooled volume for 24 samples = 960 μL)
- If processing 96 samples, use four 2 mL LoBind tubes

Pool cleanup

Perform 1X SMRTbell cleanup bead purification of pooled libraries



5. SMRTbell cleanup bead purification of pooled libraries


✓ Step	Instructions
5.1	Add 960 µL (1.0X) of resuspended, room-temperature SMRTbell cleanup beads to each pool. <i>If using less than 960 µL, add 1.0X (v/v) concentration of SMRTbell cleanup beads.</i>
5.2	Pipette-mix or invert the sample until the beads are evenly distributed. Quick-spin the samples to collect liquid.
5.3	Leave at room temperature for 10 minutes to allow DNA to bind beads.
5.4	Place samples on a magnet and allow beads to separate fully from the solution. The solution should be clear and beads pelleted to the magnet.
5.5	Slowly pipette off the supernatant without disturbing the beads. Discard the supernatant.
5.6	Slowly dispense 200 µL, or enough to cover the beads, of freshly prepared 80% ethanol to each sample. After 30 seconds, pipette off the 80% ethanol and discard.
5.7	Repeat the previous step.
5.8	Remove residual 80% ethanol: <ul style="list-style-type: none">Remove the samples from the magnet and quick-spin to collect liquid.Place samples back on the magnet and wait until beads separate fully from the solution.Carefully pipette off the residual 80% ethanol without disturbing the bead pellet and discard.
5.9	Remove samples from the magnet and immediately add 40 µL of elution buffer . Resuspend by pipetting mixing until beads are evenly distributed in solution. Quick-spin samples if necessary to collect liquid.
5.10	Leave samples at room temperature for 5 minutes to elute DNA off beads.
5.11	Place samples back on the magnet and allow beads to separate fully from the solution. The solution should be clear, and beads pelleted to the magnet before proceeding.
5.12	Slowly pipette off the cleared supernatant without disturbing the beads. Transfer supernatant to a new tube strip. Discard the old tube with beads.
5.13	Proceed to the next step (nuclease treatment), or store samples at 4°C.

• **IMPORTANT!** Allow SMRTbell cleanup beads to come up to room temperature by bringing them out of 4°C storage at least 1.5 hours prior to beginning

SAFE STOPPING POINT - Store at 4°C

Nuclease treatment

Remove unligated DNA fragments and leftover SMRTbell adapters from the sample



6. Nuclease treatment

✓ Step Instructions for nuclease treatment

Add the following components in the order and volume listed below to each of the pools from the previous step.

Nuclease mix		
✓ Tube	Component	Volume
Light purple	Nuclease buffer M96	5 μ L
Light green	Nuclease mix M96	5 μ L
Total volume		10 μ L

6.1

6.2 Pipette-mix and spin down to collect liquid. The total volume should equal **50 μ L**.

Run the **Nuclease treatment** thermocycler program. Set lid temperature to 75°C if programmable.

Step	Time	Temperature
1	15 min	37°C
2	Hold	4°C


6.3

- For Nuclease treatment reaction, each sample requires 10 μ L of **Nuclease mix** solution
- Prepare the appropriate volume of Nuclease mix solution required to process all samples by adding 1 volume of **Nuclease buffer M96** and 1 volume of **Nuclease mix M96** to a new tube¹

- Run **Nuclease treatment** thermocycler program
- Set the lid temperature to 75°C (if lid temp. is not programmable, leave at 95–105°C)

Final cleanup

Perform SMRTbell cleanup bead purification of final libraries



7. SMRTbell cleanup bead purification of final libraries

Step	Instructions
7.1	For DNA >3 kb , add 50 µL (1.0X) of resuspended, room-temperature SMRTbell cleanup beads to each sample. For DNA <3kb , add 65 µL (1.3x) of resuspended, room-temperature SMRTbell cleanup beads to the sample.
7.2	Pipette-mix the sample until the beads are evenly distributed. Quick spin the samples to collect liquid.
7.3	Leave at room temperature for 10 minutes to allow DNA to bind beads.
7.4	Place samples on a magnet and allow beads to separate fully from the solution. The solution should be clear and beads pelleted to the magnet.
7.5	Slowly pipette off the supernatant without disturbing the beads. Discard the supernatant.
7.10	Leave samples at room temperature for 5 minutes to elute DNA off beads.
7.11	Place samples back on the magnet and allow beads to separate fully from the solution. The solution should be clear, and beads pelleted to the magnet before proceeding.
7.12	Slowly pipette off the cleared supernatant without disturbing the beads. Transfer supernatant to a new tube strip. Discard tube with beads.
7.13	Recommended: quality control checkpoint. Use 1 µL from each pool to measure the concentration and size distribution of the completed library.
7.14	Please note: elutions from the bead cleanup can be pooled together at this stage prior to proceeding to ABC in SMRT® Link Sample Setup. The input library concentration into ABC must be <60 ng/µL for >10 kb libraries, <20 ng/µL for 3–10 kb libraries, and <10 ng/µL for <3 kb libraries. Failure to dilute library prior to ABC may result in low loading.

PROTOCOL COMPLETE

Note: To prepare HiFi plex prep kit 96 samples for sequencing on PacBio long-read systems, follow annealing, binding & cleanup (ABC) workflow instructions provided in **SMRT Link Sample Setup** calculator tool

- **For DNA >3 kb**
→ Add 50 µL (1.0X) of resuspended, room-temperature SMRTbell cleanup beads to sample
- **For DNA <3kb**
→ Add 65 µL (1.3x) of resuspended, room-temperature SMRTbell cleanup beads to sample

- Perform **DNA concentration QC** on final purified HiFi plex prep kit 96 SMRTbell library using a Qubit dsDNA HS assay or Quant-iT ds DNA HS assay
- Perform **DNA sizing QC** on final purified HiFi plex prep kit 96 SMRTbell library using a Femto Pulse system



Varioskan LUX multimode microplate reader
(Thermo Fisher Scientific)



Femto Pulse system
(Agilent Technologies)

- **Note:** Final HiFi plex prep kit 96 SMRTbell library concentration must be **<60 ng/µL for >10 kb; <20 ng/µL for 3-10 kb and <10 ng/L for <3 kb** to proceed with ABC
→ Failure to dilute libraries to the recommended concentration range will result in lower *P1* loading during sequencing





HiFi plex prep kit 96 example sequencing performance data

HiFi sequencing performance of HiFi plex prep kit 96 libraries

Example HiFi sequencing performance obtained for HiFi plex prep kit 96 libraries constructed from *E. coli* gDNA sheared using SPEX SamplePrep 1600 MiniG homogenizer¹

Sample ²	Polymerase RL (bp)	Yield (Gb)	HiFi Mean Length (bp)	Mean QV	P1 %
96_plex E.coli_Cell 1 (250 pM)	81,529	67.8	7,477	Q45	61
96_plex E.coli_Cell 2 (250 pM)	96,977	73.1	7,810	Q46	59




¹ Revio system data shown. **Note:** HiFi read lengths, reads/data per SMRT Cell and other sequencing performance results can vary depending on DNA sample quality, insert size, P1 loading performance & movie time. Shorter library insert sizes (<15 kb), lower DNA quality samples, and suboptimal P1 loading performance may result in HiFi data yields <90 Gb per Revio SMRT Cell.

² Samples were sequenced on a Revio system using 250 pM OPLC with 24-hours movie time.

HiFi sequencing performance of HiFi plex prep kit 96 libraries (cont.)

Example HiFi sequencing performance obtained for HiFi plex prep kit 96 libraries constructed from microbial gDNA samples sheared using Hamilton NGS STAR workstation¹

Metric	Result
HiFi reads	7.3 M
HiFi reads yield	87.64 Gb
HiFi read length (mean)	11.99 kb
HiFi read length (N50)	14,033 kb
Median read quality	Q40
Base quality \geq Q30	94.57%
Barcodes detected	96
Barcoded HiFi reads %	99.81%



- HMW DNA isolated from 96 different microbial isolates and sheared on Hamilton NGS STAR
- All samples prepared with HiFi plex prep kit 96 with no size selection used
- Single 96-plex sample pool was processed through ABC workflow with Revio polymerase kit
- Sequenced on Revio system using 200 pM on-plate loading concentration with 24 hrs movie time



Technical documentation & applications support resources

Technical resources for high-throughput HiFi library preparation, sequencing & data analysis

DNA extraction literature & other resources

- Guide & overview – Nanobind HT kits for automated HMW DNA extraction ([103-028-100](#))
- Guide & overview – Nanobind PanDNA kit ([103-394-800](#))
- Technical note – High-throughput DNA extraction ([102-326-611](#))
- Technical overview – Automated high-throughput HMW DNA extraction for PacBio long-read sequencing using Nanobind HT kits ([103-401-700](#))

DNA shearing literature & other resources

- Guide & overview – Short Read Eliminator (SRE), DNA shearing, and cleanup for the Hamilton Microlab Prep system ([103-424-100](#))
- Technical note – High-throughput DNA shearing for HiFi whole genome sequencing from whole blood samples [MP Biomedicals FastPrep-96] ([102-326-579](#))
- Technical note – High-throughput DNA shearing for long-read microbial WGS [SPEX SamplePrep 1600 MiniG] ([102-326-575](#))
- Technical note – High throughput DNA shearing using Hamilton Microlab Prep ([102-326-606](#))

HiFi SMRTbell library preparation literature & other resources

- Brochure – HiFi prep kits ([102-326-608](#))
- Guide & overview – Automated HiFi plex prep 96 for the Hamilton NGS Microlab STAR system ([103-425-800](#))
- Guide & overview – Automated HiFi prep 96 and HiFi annealing, binding, and cleanup for the Hamilton NGS Microlab STAR system ([103-425-700](#))
- Procedure & checklist – Preparing multiplexed whole genome and amplicon libraries using the HiFi plex prep kit 96 ([103-418-800](#))
- Procedure & checklist – Preparing whole genome libraries using the HiFi prep kit 96 ([103-420-700](#))
- Procedure & checklist – Preparing whole genome and metagenome sequencing libraries using SMRTbell prep kit 3.0 ([102-166-600](#))
- Technical Overview – Technical overview – HiFi library preparation using HiFi prep kits for high-throughput sequencing on PacBio long-read systems ([103-424-600](#))

Technical resources for high-throughput HiFi library preparation, sequencing & data analysis (cont.)

Third-party automation instrumentation literature

- Application note – Automation of Long-Read Sequencing Library Preparation with PacBio SMRTbell prep kit 3.0 on Hamilton NGS STAR MOA (Hamilton [AN-2305-05](#))
- Technical note – Automated extraction of High Molecular Weight (HMW) DNA with PacBio Nanobind technology on the Hamilton NIMBUS Presto Assay Ready Workstation (Hamilton [AN-2205-05](#))
- Technical note – Automated Isolation of High Molecular Weight (HMW) DNA from Human Blood Samples with PacBio Nanobind Technology on the Hamilton NIMBUS Presto – Next Level Preparation of Extracts for Long-Read Sequencing (Hamilton [AN-2212-03](#))

HiFi data analysis literature & other resources

- SMRT Link software installation guide [[Link](#)]
- SMRT Link user guide [[Link](#)]
- SMRT Tools reference guide [[Link](#)]

HiFi WGS applications literature & other resources

- Application note – Consolidated analysis tools with the PacBio WGS Variant Pipeline ([102-326-588](#))
- Application note – Robust detection of somatic variants from tumor-normal samples with highly accurate long-read whole genome sequencing ([102-326-582](#))
- Overview – Human genomics: Unlock your next great discovery with HiFi sequencing ([102-326-536](#))

Webinars

- PacBio webinar (2024) – More samples, lower costs, less time: New PacBio HiFi prep kits + microbial WGS and antimicrobial resistance [[Link](#)]



www.pacb.com

Research use only. Not for use in diagnostic procedures. © 2024 Pacific Biosciences of California, Inc. ("PacBio"). All rights reserved. Information in this document is subject to change without notice. PacBio assumes no responsibility for any errors or omissions in this document. Certain notices, terms, conditions and/or use restrictions may pertain to your use of PacBio products and/or third-party products. Refer to the applicable PacBio terms and conditions of sale and to the applicable license terms at pacb.com/license. Pacific Biosciences, the PacBio logo, PacBio, Circulomics, Omniome, SMRT, SMRTbell, Iso-Seq, Sequel, Nanobind, SBB, Revio, Onso, Apton, Kinnex, and PureTarget are trademarks of PacBio.