PacBi

Technical overview – Nanobind HT kits for automated HMW DNA extraction


High-throughput HMW DNA extraction solutions for PacBio long-read sequencing

Nanobind HT kits for automated HMW DNA extraction

Technical Overview

- 1. Nanobind HT method overview
- 2. Nanobind HT workflow details
- 3. Nanobind HT example performance data
- 4. Technical documentation & applications support resources

Nanobind HT HMW DNA extraction for SMRT sequencing: Getting started



Nanobind HT HMW DNA extraction workflow overview



Available Nanobind HT HMW DNA extraction protocols

Select the appropriate Nanobind HT Procedure & checklist based on the type of automation system, sample type, desired sample throughput and available input sample amount

Automation			Input	Max #	Nanobind HT kit (96-RXN)	
system	Procedure & checklist	Sample type	Sample amount	samples per run	Kit	Part no.
KingFisher Duo	Nanobind HT HMW DNA extraction – Cultured cells – KingFisher Duo [102-996-200]	Mammalian cells	1x10 ⁶ cells*	12	Nanobind HT CBB kit	102-762-700
-	Nanobind HT HMW DNA extraction – 200 $_{\mu}L$ whole blood – KingFisher Duo [102-995-800]	Whole blood	200 μL	12		102-762-700
	Nanobind HT HMW DNA extraction – 1 mL whole blood – KingFisher Duo [102-995-400]	Whole blood	1 mL	6	Nanobind HT 1 mL whole blood kit	102-762-800
KingFisher Flex	Nanobind HT HMW DNA extraction – Cultured cells – KingFisher Flex [102-996-300]		1x10 ⁶ cells*	96	Nanobind HT CBB kit	102-762-700
	Nanobind HT HMW DNA extraction – 200 μL whole blood – KingFisher Flex [102-995-900]	Whole blood	200 μL	96	Nanobinu III CBB Kit	102 102 100
	Nanobind HT HMW DNA extraction – 1 mL whole blood – KingFisher Flex [102-995-500]	Whole blood	1 mL	24	Nanobind HT 1 mL whole blood kit	102-762-800
KingFisher Apex	Nanobind HT HMW DNA extraction – Cultured cells – KingFisher Apex [102-996-100]	Mammalian cells	1x10 ⁶ cells*	96	Nanobind HT CBB kit	102-762-700
	Nanobind HT HMW DNA extraction – 200 $_{\mu}L$ whole blood – KingFisher Apex [102-995-700]	Whole blood	200 μL	96		102-762-700
	Nanobind HT HMW DNA extraction – 1 mL whole blood – KingFisher Apex [102-995-300]	Whole blood	1 mL	24	Nanobind HT 1 mL whole blood kit	102-762-800
Hamilton NIMBUS	Nanobind HT HMW DNA extraction – Cultured cells – Hamilton NIMBUS Presto [102-996-400]	Mammalian cells	1x10 ⁶ cells*	96	Nanobind HT CBB kit	102-762-700
Presto	Nanobind HT HMW DNA extraction – 200 μ L whole blood – Hamilton NIMBUS Presto [102-996-000]	Whole blood	200 μL	96		102-762-700
	Nanobind HT HMW DNA extraction – 1 mL whole blood – Hamilton NIMBUS Presto [102-995-600]	Whole blood	1 mL	24	Nanobind HT 1 mL whole blood kit	102-762-800

 * Cell samples are typically resuspended in 50 μL of 1X PBS for Nanobind HT processing.

Pace See PacBio's <u>Documentation</u> website for the most up-to-date list of Nanobind HT HMW DNA sample preparation protocols.

Supported automation platforms for high-throughput HMW DNA extraction using Nanobind HT kits

Nanobind HT HMW DNA extraction kits are optimized for use with Thermo Fisher Scientific KingFisher sample purification systems and Hamilton Microlab NIMBUS automated liquid handlers










		o Fisher Sc Fisher Duo I			o Fisher So ngFisher F			o Fisher Sc ngFisher Ap		Microla	Hamilton ab NIMBUS	Presto
Instrument design		benchtop sa on system	mple	 Benchtop sample purification system 		ification	 Benchtop sample purification system with touchscreen 		 Hamilton liquid handler with integrated KingFisher Presto sample purification system 		Presto	
Nanobind HT throughput	anobind HT throughput • 6 or 12 samples / run • 24 o		• 24 or 96	samples / rui	n	• 24 or 96 samples / run		• 24 or 96 samples / run				
Nanobind HT workflow automation level			equires limited	 Semi-automated Manual plate filling + requires limited user interaction after run start 		 Semi-automated Manual plate filling + requires limited user interaction after run start 		 Fully automated Automated plate filling + fully walk- away after run start 				
Nanobind HT processing time	12 (1x10 ⁶ cells)	12 (200 μL)	<mark>6</mark> (1 mL)	96 (1x10 ⁶ cells)	96 (200 μL)	24 (1 mL)	96 (1x10 ⁶ cells)	<mark>96</mark> (200 μL)	24 (1 mL)	96 (1x10 ⁶ cells)	<mark>96</mark> (200 μL)	24 (1 mL)
Total time	1 hr 20 min	1.5 hrs	1 hr 55 min	1 hr 40 min	1 hr 50 min	2 hrs 15 min	1 hr 40 min	1 hr 50 min	2 hrs 15 min	2 hrs 25 min	2 hrs 35 min	2 hrs 30 min
Automation run time	1 hr 5 min	1 hr 15 min	1 hr 45 min	1 hr 5 min	1 hr 15 min	1 hr 45 min	1 hr 5 min	1 hr 15 min	1 hr 45 min	2 hrs	2 hrs 10 min	2 hrs 10 min
Hands-on time	15 min	15 min	10 min	35 min	35 min	30 min	35 min	35 min	30 min	25 min	25 min	20 min


Nanobind HT HMW DNA extraction workflow overview



* Refer to appropriate vendor websites to download user guides and other documentation for specific automation systems.


How Nanobind HMW DNA extraction kits work with KingFisher systems


Nanobind kits enable extraction of high-molecular weight DNA from common samples as well as more challenging samples such as animal tissue, insects, fungi, and plants



Rapid magnetic purification

- Rapid and simple bind, wash, and elute protocol
- Can perform manual processing using a magnetic separation rack
- Nanobind disks are automation compatible for high-throughput applications





Extracted HMW DNA bound to a Nanobind disk

8


How Nanobind HMW DNA extraction kits work with KingFisher systems (cont.)

Kingfisher sample purification systems employ inverse magnetic particle processing (MPP) technology to enable automated transfer and processing of magnetic Nanobind disks in a microplate format.

Kingfisher instruments




- Automate sample extraction by moving magnetic particles or objects (not liquids)
- Inverse MPP uses magnetic rods covered with a disposable tip comb and PCR or deepwell plates
- Function without any liquid dispensing or aspiration parts or devices
- Enable quicker reaction times and a more efficient wash process




How Nanobind HMW DNA extraction kits work with KingFisher systems (cont.)

Standard manual processing with Nanobind disks


Liquid is moved into and out of a tube using a pipette while the magnetic Nanobind disk remains immobilized inside the tube





Automated inverse magnetic particle processing (MPP) with Nanobind disks


Rather than moving liquids, **magnetic Nanobind disks are moved** from plate to plate using magnetic rods covered with a disposable plastic tip comb



Schematic illustration of KingFisher rotary platform modified *after* Ficarro et al., Anal Chem. 2009 June 1; 81(11): 10 4566–4575. doi:10.1021/ac9004452.

PacBi

PacBio scripts support high-throughput HMW DNA extraction using Nanobind HT kits with semi-automated and fully-automated workflows





Nanobind HT Procedure & checklist documentation overview

Twelve Procedures & checklists to support 3 high-throughput HMW DNA extraction workflows on 4 automation systems

Nanobind HT Procedure & checklist protocols describe procedures for automated HMW DNA extraction from mammalian cells and 200 μ L/1 mL of blood using a Kingfisher Duo/Flex/Apex system or Hamilton NIMBUS Presto system

Procedure & checklist contents

- 1. Equipment and reagent list
- 2. Sample input requirements
- 3. Automated HMW (50 kb 300+ Mb) DNA extraction protocol
- 4. QC procedure
- 5. Expected QC results
- 6. Troubleshooting FAQ

Nanobind HT Procedures & checklists have been validated using human whole blood samples and mammalian cells.



Extracting HMW DNA using the Nanobind[®] HT CBB kit for mammalian cultured cells on the KingFisher Apex system Procedure & checklist **PacBi**

This procedure describes the workflow for high-throughput automated extraction of HMW (50–300 kb) DNA from cultured mammalian cells using the Thermo Fisher KingFisher Apex robotic instrument. This procedure requires the Nanobind HT CBB kt (102-762-700) and is recommended for HiFi sequencing.

The Nanobind HT CBB kit has enough reagents for 96 extractions to be run in one of the following formats: 1 run x 96 samples, 2 runs x 48 samples, or 4 runs x 24 samples. We do not recommend running fewer than 24 samples per run as the kit is designed to accommodate dead volumes for a maximum of 4 runs (4 runs x 24 samples).

Required equipment and materials

Equipment/reagent	Manufacturer (part number)
Nanobind HT CBB kit	PacBio® (102-762-700)
KingFisher Apex System	Thermo Fisher Scientific (5400930, includes Apex 96 deep well magnet head)
KingFisher Apex 96 deep-well magnet head	Thermo Fisher Scientific (24079930)
KingFisher Apex 96 deep-well heating block	Thermo Fisher Scientific (24075920)
KingFisher 96 deep-well plates, barcoded	Thermo Fisher Scientific (95040450B)
KingFisher 96 deep-well tip combs, barcoded	Thermo Fisher Scientific (97002534B)
Ethanol (96–100%)	Any major lab supplier (MLS)
Isopropanol (100%)	Any MLS
UV/Vis	Thermo Fisher Scientific NanoDrop 2000
Fluorescent DNA Quantification	Thermo Qubit 3.0, dsDNA BR and RNA BR Assay Kits

Before you begin

Prior to starting


Buffer CW1 and CW2 are supplied as concentrates. This kit uses CW1 with a 60% final ethanol concentration. This kit uses CW2 with a 60% final ethanol concentration. Before using, add the appropriate amount of ethanol (96-100%) to Buffer CW1 and Buffer CW2 as indicated on the bottles.

See PacBio sample preparation Documentation

PacBie See PacBio's Documentation website for the most up-to-date list of Nanobind HT HMW DNA sample preparation protocols.

Equipment & materials required for Nanobind HT HMW DNA extraction

For each type of sample purification system, use the correct matching magnet heads, heating blocks, plates and other consumables as listed in the appropriate Nanobind HT procedure



Sample input requirements for Nanobind HT HMW DNA extraction

	Human whole blood	Mammalian cells						
Nanobind HT kit	Nanobind HT CBB kit (96 Rx) (102-762-700) Nanobind HT 1 mL whole blood kit (96 Rx) (102-762-800)	Nanobind HT CBB kit (96 Rx) (102-762-700)						
Storage conditions	 Use Potassium EDTA (K₂EDTA) anticoagulant Sodium heparin (NaHep) and citrate (NaCit) also performed well in limited testing Store at 4°C for ≤2 days to prevent sample degradation 	• N/A						
	 Cells: White blood cell (WBC) count should be ≥4 × 10⁶ cells/mL* to achieve ≥3 ug DNA yield for input into SMRTbell library prep 	 Use 1x10⁶ diploid human cells* or equivalent Determine cell counts using a hemocytometer or cell counter 						
Sample input amount	 Blood: Use either 200 μL or 1 mL of input human whole blood 200 μL: DNA yield is 1 – 12 μg depending on donor WBC conc. 1 mL: DNA yield is 3 – 70 μg depending on donor WBC conc. 	 For non-diploid/non-human cells: Scale cell input to contain 5–25 µg of DNA Warning: >25 µg inputs can cause Nanobind disks to be "dropped" in Lysis/Binding solution and/or cause well-to-well contamination 						
Sample freezing / thawing conditions	 Freeze blood samples as quickly as possible after being drawn Aliquot blood samples to avoid repeated freeze-thaws Thaw blood samples in water bath/dry block heater at 37 °C (15 min Mix by inverting tube >15 times immediately prior to use Note: Improperly thawed and mixed samples may result in inconsistent DNA yield and purity 	 Cell pellets should be frozen dry with as much liquid removed as possible. No cryoprotectant needed 						
No systematic difference observed in DNA QC or sequencing results between fresh and frozen samples								

Before you begin

End-to-end overall Nanobind HT workflow

Prepare sample, storage & reagent plates

Select appropriate automation script to run

Load prepared plates and other consumables onto work deck

Run selected Nanobind HT DNA extraction protocol

Store extracted DNA samples and perform QC

Semi-automated workflow (KingFisher Duo/Flex/Apex)





Fully-automated workflow (Hamilton NIMBUS Presto)



Prior to beginning Nanobind HT protocol

- Ensure that the proper automation script has been installed (see <u>Nanobind HT kit Guide & overview</u> "Programs" section)
 - To obtain KingFisher automation scripts, contact PacBio Technical Support at support@pacb.com
 - To obtain Hamilton NIMBUS Presto automation scripts, contact <u>Hamilton Technical Support</u> at <u>roboticservice@hamiltoncompany.com</u>
- Ensure that instrument is set up with the correct magnet head and heating block.

Before you begin (cont.)

End-to-end overall Nanobind HT workflow

Prepare sample, storage & reagent plates

Semi-automated workflow (KingFisher Duo/Flex/Apex)





Fully-automated workflow (Hamilton NIMBUS Presto)



Select appropriate automation script to run

Load prepared plates and other consumables onto work deck

Run selected Nanobind HT DNA extraction protocol

Store extracted DNA samples and perform QC

PacBi

Prior to beginning Nanobind HT protocol

For Hamilton NIMBUS Presto only:

- Prepare a *.xls sample worklist
- Column 1 = "Sample_ID"

•

- List one sample ID for each sample to be processed
- Column 2 = "Sample_Position"
 - List the corresponding well location for each sample

Sample worklist enables NIMBUS liquid handler to **automatically** prepare lysis reactions on instrument work deck using a loaded sample plate

	А	В	С	D
1	Sample_ID	Sample_Position		
2	Sample1	A1		
3	Sample2	B1		
4	Sample3	C1		
5	Sample4	D1		
6	Sample5	A2		
7	Sample6	B2		
8	Sample7	C2		
9	Sample8	D2		
10	Sample9	A3		
11	Sample10	B3		
12	Sample11	C3		
13	Sample12	D3		
14	Sample13	A4		
15	Sample14	B4		

Sample, storage & reagent plate preparations

NIMBUS Presto workflow requires upfront preparation of **fewer plates** compared to Duo/Flex/Apex workflows

End-to-end overall Nanobind HT workflow Prepare sample, storage & reagent plates Load prepared plates and other consumables onto work deck Run selected Nanobind HT DNA

Store extracted DNA samples

and perform QC

PacBi

Semi-automated workflow (KingFisher Duo/Flex/Apex)

Manually prepare deep-well plates with required quantities of Nanobind reagents & consumables

Example list of 7 DW plates to be prepared for Kingfisher Apex.

Plate	Plate name	Reagent	Notes
1	Lysis / binding	Sample / lysis / binding solution	Prepare only after all other plates have been prepared (see below)
2	Nanobind storage	1 Nanobind disk per well*	200 μL / cells \rightarrow 3 mm disk 1 mL \rightarrow 5 mm disk
3	Wash 1	Buffer CW1	200 μL sample / cells $\ \rightarrow$ 700 μL 1 mL sample $\ \rightarrow$ 2 mL
4	Wash 2	Buffer CW2	200 μL sample / cells \rightarrow 700 μL 1 mL sample \rightarrow 2 mL
5	Wash 3	Buffer CW3	200 μL sample / cells \rightarrow 700 μL 1 mL sample \rightarrow 2 mL
6	Elution	Buffer EB	200 μL sample / cells \rightarrow 100 μL 1 mL sample \rightarrow 200 μL
7	Tip comb storage	For 24- or 96-tip comb	Insert tip comb into a DW plate for loading onto instrument

Plate 1 (lysis / binding plate) preparation procedure :

- Add Proteinase K
- 2. Add whole blood sample
- 3. Add Buffer BL3
- 4. Add RNase A

Note: Sample and reagents MUST be added to each well in the order described in the Procedure

Example list of 3 DW plates to be prepared for NIMBUS Presto.

Fully-automated workflow

(Hamilton NIMBUS Presto)

Plate	Plate name	Reagent	Notes		
1	Sample plate	Samples for analysis	Prepare only after all other plates have been prepared		
2	Nanobind storage	1 Nanobind disk per well*	200 μL sample / cells \rightarrow 3 mm disk 1 mL sample \rightarrow 5 mm disk		
3	Tip comb storage	For 24- or 96-tip comb	Insert tip comb into a DW plate for loading onto instrument		

Plate 1 (sample plate) preparation procedure:


1. Add whole blood sample to appropriate well positions as specified in *.xls sample worklist.




1	A	В	С	D
1	Sample_ID	Sample_Position		
2	Sample1	A1		
3	Sample2	81		
4	Sample3	C1		
5	Sample4	D1		
6	Sample5	A2		
7	Sample6	B2		
8	Sample7	C2		
9	Sample8	D2		
10	Sample9	A3		
11	Sample10	B3		
12	Sample11	C3		
13	Sample12	D3		. الس ، الدال

* Nanobind disks do not need to be perfectly centered, but please ensure they are at the bottom of the well and not stuck to the side walls.

Automation script selection



Work deck preparation



Automation script execution

End-to-end overall Nanobind HT workflow

Prepare sample, storage & reagent plates

Select appropriate automation script to run

Load prepared plates and other consumables onto work deck

Run selected Nanobind HT DNA extraction protocol

Store extracted DNA samples and perform QC


PacBi

Semi-automated workflow (KingFisher Duo/Flex/Apex)

The run will start <u>automatically</u> when final plate is loaded and **Next** button is pressed.



After run starts, screen shows time until user action required, remaining run time, and current protocol step



User action required: When prompted, remove Lysis/Binding plate (Plate 1) from instrument and manually add IPA to each well. After adding IPA, return plate to instrument and resume program.



Note: Add isopropanol (IPA) gently against the **side** of the well into the Lysis / Binding solution.*

Fully-automated workflow (Hamilton NIMBUS Presto)

Close front cover of instrument to protect against environmental contamination and click **Ok** to start run



NIMBUS Presto laptop screen will request final verification that correct magnet head is installed prior to starting run



No user action required until run completes

* Adding IPA directly to the Lysis / Binding solution may affect extraction purity.

Sample storage and QC

End-to-end overall Nanobind HT workflow

Prepare sample, storage & reagent plates

Select appropriate automation script to run

Load prepared plates and other consumables onto work deck

Run selected Nanobind HT DNA extraction protocol

Store extracted DNA samples and perform QC

PacBi

Semi-automated workflow (KingFisher Duo/Flex/Apex)

Fully-automated workflow (Hamilton NIMBUS Presto)

When run ends, remove each plate from instrument as prompted.












Manually transfer eluates from Elution Plate to a new storage plate or storage tubes



Elution plate New storage plate tube to correct to corr

KingFisher Duo/Flex/Apex automation scripts are designed to leave Nanobind disks in elution plate after program completes

NIMBUS Presto automation scripts are designed to place Nanobind disks back into the disk storage plate, so the elution plate should not contain any disks inside

Pipette mix extracted DNA sample 10 times with a standard P200 pipette to homogenize and disrupt any unsolubilized "jellies" that may be present. Take care to disrupt any regions that feel more viscous than other regions.



Note: Limited pipette mixing will not noticeably reduce DNA size or sequencing read lengths but is important for accurate quantitation and consistent sequencing performance

Let eluate rest overnight at room temp. to allow DNA to solubilize. Following overnight rest, pipette mix 10 times with a standard P200 pipette. Visible "jellies" should disperse after resting.

Recommended QC procedure for HMW DNA samples extracted using Nanobind HT kits

Use recommended tools to perform DNA QC evaluation of high-molecular weight DNA



Qubit fluorometer

DNA concentration QC

Use Qubit dsDNA broad range (BR) assay kit for DNA concentration QC

 We recommend using a Qubit fluorometer (Thermo Fisher Scientific) with the dsDNA BR assay kit since it provides more reliable concentration measurements of (unsheared) HMW DNA

Optional: Determine RNA concentration using the Qubit RNA BR assay kit (Thermo Fisher Scientific)

• We recommend taking a single measurement to obtain an approximate RNA concentration reading



Nanodrop

DNA purity QC

Perform UV/VIS assay measurements with a Nanodrop to determine total nucleic acid conc. as well as purity (A260/A280, A260/230)

- Note: Inconsistent Nanodrop concentration or spectrophotometric readings may be obtained if the DNA is very heterogeneous or contains large amounts of unsolubilized "jellies"
- For information on how to disrupt particularly viscous regions, refer to the kit Guide & overview "Heterogeneity and viscosity" section or the "Troubleshooting FAQ" section in the Procedures & checklists





DNA sizing QC

Use a Femto Pulse system with Genomic DNA 165 kb Kit (Agilent Technologies) for HMW genomic DNA sizing QC

- We recommend diluting samples to 250 pg/µL (use wide-bore pipette for making dilutions and gentle mixing)
- · Avoid mixing with a standard pipette! This will shear dilute solutions of DNA.

Expected QC results for HMW DNA samples extracted using Nanobind HT kits

Example ranges typically observed for DNA concentration, DNA purity and DNA sizing QC metrics



Qubit fluorometer

DNA concentration QC

- 200 µL human whole blood typically yields ~3 12 µg of DNA depending on donor white blood cell (WBC) count
- 1 mL of human whole blood typically yields ~3 70 µg depending on donor WBC count
- 1x10⁶ GM24385 cells typically yields ~4 10 µg
- 1x10⁶ MCF-7 cells typically yields ~12 18 μg



Nanodrop

DNA purity QC

- 260/280 ratios should consistently be 1.8 2.0
- 260/230 ratio can vary from 1.3 2.2
- Samples with UV purities within expected ranges should sequence well
- UV purities outside of these ranges may indicate abnormalities in extraction process



Femto Pulse system

DNA sizing QC

 Mode of extracted human cell and whole blood DNA measured on Femto Pulse system (Agilent Technologies) is typically 100 kb+

Refer to Nanobind HT Procedure & checklist documentation for guidance on troubleshooting specific types of DNA extraction quality issues

Nanobind HT HMW DNA extraction example performance data




Nanobind HT HMW DNA extraction performance on KingFisher Apex and Hamilton NIMBUS Presto

Example HMW DNA extraction QC results

	Human whole blood	Mammalian cells
•	200 μL or 1 mL fresh whole blood (max. 2 days storage at 4°C) or frozen	 1 x 10⁶ diploid cells (fresh culture or frozen cell pellet)
•	For donor WBC count $\ge 4 \times 10^6$ cells/mL, expected HMW DNA yield is 3 – 12 µg for 200 µL blood sample and 3 – 70 µg for 1 mL blood sample	- Expected HMW DNA yield is 5 – 15 μ g

Sample	Robot	Input sample	DNA yield (µg)	A260/ A280 ratio	A260/ A230 ratio	DNA mean size (kb)	GQN at 10 kb
Mammalian cell (HG001)	Nimbus	10 ⁶ cells	6.4	1.92	2.04	167	9.4
Mammalian cell (HG001)	Apex	10 ⁶ cells	7.0	1.88	2.08	113	9.6
Whole human blood	Nimbus	200 µL	7.0	1.83	1.96	98	9.4
Whole human blood	Apex	200 µL	6.1	1.83	1.80	124	9.7
Whole human blood	Nimbus	1 mL	32.6	1.86	1.86	94	9.3
Whole human blood	Apex	1 mL	36.4	1.89	2.18	115	9.4



Sizing QC results exceed SPK3 WGS library preparation input DNA quality requirements

- Mean DNA fragment size >90 kb
- >93% of DNA fragments longer than 10 kb (GQN at 10 kb)

Input DNA quality requirements for Procedure & checklist – Preparing whole genome and metagenome sequencing libraries using SMRTbell prep kit 3.0 (102-166-600)

DNA Quality	Human, plant, and animal samples	Microbial and Metagenomic samples
DNA size distribution (Femto Pulse system)	50% ≥30 kb & 90% ≥10 kb	90% ≥7 kb


Nanobind HT HMW DNA extraction performance on KingFisher Apex and Hamilton NIMBUS Presto

Example HiFi sequencing performance results


Sample	Robot	Input sample	HiFi mean read length (bp)	HiFi yield (Gb)
Mammalian cell (HG001)	Nimbus	10 ⁶ cells	16,563	41.9
Mammalian cell (HG001)	Apex	10 ⁶ cells	16,544	35.4
Whole human blood	Nimbus	200 µL	18,067	37.3
Whole human blood	Apex	200 µL	13,026	34.8
Whole human blood	Nimbus	1 mL	17,163	36.8
Whole human blood	Apex	1 mL	15,863	36.1

SPK3 WGS libraries constructed with Nanobind HT extracted DNA show excellent HiFi sequencing performance

- HiFi yield >35 Gb
- Mean HiFi read length 13 18 kb
- Limited number of reads <8 kb



HiFi read length distribution





PacBi

Technical documentation & applications support resources



Technical resources for Nanobind HT HMW DNA extraction

PacBio Nanobind HT HMW DNA sample preparation literature

- Brochure Nanobind HT HMW DNA extraction (<u>102-326-565</u>)
- Guide & overview Nanobind HT kits for automated HMW DNA extraction (103-028-100)
- Overview Nanobind HT HMW DNA extraction-robotic procedures (<u>103-032-000</u>)
- Procedure & Checklist : Nanobind HT HMW DNA extraction 1 mL whole blood KingFisher Apex (102-995-300)
- Procedure & Checklist : Nanobind HT HMW DNA extraction 1 mL whole blood KingFisher Duo (102-995-400)
- Procedure & Checklist : Nanobind HT HMW DNA extraction 1 mL whole blood KingFisher Flex (102-995-500)
- Procedure & Checklist : Nanobind HT HMW DNA extraction 1 mL whole blood Hamilton NIMBUS Presto (102-995-600)
- Procedure & Checklist : Nanobind HT HMW DNA extraction 200 μL whole blood KingFisher Apex (102-995-700)
- Procedure & Checklist : Nanobind HT HMW DNA extraction 200 μL whole blood KingFisher Duo (102-995-800)
- Procedure & Checklist : Nanobind HT HMW DNA extraction 200 μL whole blood KingFisher Flex (102-995-900)
- Procedure & Checklist : Nanobind HT HMW DNA extraction 200 μL whole blood Hamilton NIMBUS Presto (102-996-000)
- Procedure & Checklist : Nanobind HT HMW DNA extraction Cultured cells KingFisher Apex (102-996-100)
- Procedure & Checklist : Nanobind HT HMW DNA extraction Cultured cells KingFisher Duo (<u>102-996-200</u>)
- Procedure & Checklist : Nanobind HT HMW DNA extraction Cultured cells KingFisher Flex (<u>102-996-300</u>)
- Procedure & Checklist : Nanobind HT HMW DNA extraction Cultured cells Hamilton NIMBUS Presto (102-996-400)
- Technical overview Nanobind HT kits for automated HMW DNA extraction (<u>103-020-800</u>)

PacBi

Technical resources for Nanobind HT HMW DNA extraction (cont.)

Third-party automation instrumentation literature

- Hamilton Technical Note Automated extraction of High Molecular Weight (HMW) DNA with PacBio Nanobind technology on the Hamilton NIMBUS Presto Assay Ready Workstation (<u>AN-2205-05</u>)
- Hamilton Technical Note Automated Isolation of High Molecular Weight (HMW) DNA from Human Blood Samples with PacBio Nanobind Technology on the Hamilton NIMBUS Presto - Next Level Preparation of Extracts for Long-Read Sequencing (<u>AN-2212-03</u>)

Posters

• PacBio poster (2022) – High-throughput workflow for human whole genome sequencing using PacBio HiFi [Link]

Nanobind HT consumable product literature

- Nanobind HT 1 mL blood kit (PN <u>102-762-800</u>)
 - Product insert: Nanobind HT 1 mL blood kit 4C
 - Product insert: Nanobind HT 1 mL blood kit RT
- Nanobind HT CBB kit (PN <u>102-762-700</u>)
 - Product insert Nanobind HT CBB kit 4C
 - Product insert Nanobind HT CBB kit RT
- PacBio sample prep consumables website [Link]

PacBio

www.pacb.com

Research use only. Not for use in diagnostic procedures. © 2023 Pacific Biosciences of California, Inc. ("PacBio"). All rights reserved. Information in this document is subject to change without notice. PacBio assumes no responsibility for any errors or omissions in this document. Certain notices, terms, conditions and/or use restrictions may pertain to your use of PacBio products and/or third-party products. Refer to the applicable PacBio terms and conditions of sale and to the applicable license terms at pacb.com/license. Pacific Biosciences, the PacBio logo, PacBio, Circulomics, Omniome, SMRT, SMRTbell, Iso-Seq, Sequel, Nanobind, and SBB are trademarks of PacBio. All other trademarks are the sole property of their respective owners.