

The PacBio logo is displayed in a bold, pink, sans-serif font. A single, large pink droplet is positioned at the end of the word 'Bio', appearing to fall from the top right corner of the frame. The background is a blurred laboratory setting with a multi-well plate containing pink liquid in the foreground and a pipette tip visible in the upper right.

## Technical overview

# PureTarget library preparation using PureTarget kit 2.0

Vega system ICS v1.1

Revio system ICS v13.3

SMRT Link v25.3


PN 103-742-600 Rev 01 | September 2025

# Technical overview

## PureTarget library preparation using PureTarget kit 2.0

1. PureTarget method & products overview
2. Library preparation workflow details
3. Sequencing preparation workflow details
4. Data analysis workflow overview
5. Example sequencing performance data
6. Technical documentation & applications support resources
7. Appendix

# PureTarget library preparation using PureTarget kit 2.0: Getting started



**Application-specific literature**

**PacBio**

Application note

### Comprehensive genotyping with the PureTarget repeat expansion panel and HiFi sequencing

**Introduction**

Tandem repeats are regions of the genome consisting of repetitive units of specific DNA sequences. These regions are hypermutable: they can increase in length across generations and may have variable lengths within somatic tissues of an individual (Paulson 2018). Many tandem repeats become pathogenic when they exceed a length threshold, which varies from gene to gene, resulting in mutations called repeat expansions (REs) (Ibáñez et al., 2022). Repeat expansions have been linked to dozens of diseases and cancer, most notably neurodegenerative disorders like Huntington's disease, Fragile-X disorder, spinocerebellar ataxia, and myotonic dystrophies (DePieris & Mandel, 2021). Disease severity and age of onset of these conditions are often associated with their repeat length (Ibáñez et al., 2022).

Though common, these regions are challenging to characterize and as such, the majority of patients with rare neurological diseases remain undiagnosed (Ibáñez et al., 2022). Repeat expansions have historically been

profiled with Southern blotting and PCR-based assays (Tartien, 2003), and more recently short-read bioinformatics methods (Dobchenko et al., 2020; Dashnow et al., 2022), but these methods are limited by throughput and read length, respectively. Given that many repeat expansions are longer than the typical length of short reads, comprehensive genotyping of REs requires a high-throughput long-read sequencing approach that can reliably manage the high structural variability of these regions. Powered by the exceptional accuracy of HiFi sequencing and the Tandem Repeat Genotyping Tool (TRGT), the PacBio® PureTarget™ repeat expansion panel offers more comprehensive genotyping for 20 of the most important repeat expansions for human health. This application note demonstrates the performance of the PureTarget repeat expansion panel, and presents the PureTarget panel as a scalable, more comprehensive solution for profiling repeat expansions, compared to legacy genotyping and next-generation sequencing methods.

© 2024 PacBio. All rights reserved. Research use only. Not for use in diagnostic procedures. 102-326-614 REV01 MAR2024 **PacBio**

**Application note – Comprehensive genotyping with the PureTarget repeat expansion panel and HiFi sequencing (102-326-614)**

**Application note – Consolidate challenging genes with PureTarget carrier screen panel (102-326-653)**

**Technical note – A practical guide to amplification-free PureTarget custom panels (102-326-652)**

**Application-specific protocol**

**PacBio**

### Generating PureTarget™ libraries with PureTarget kit 24 – manual protocol

Procedure & checklist

**Before you begin**

This procedure describes the workflow for manually generating PureTarget libraries.

Overview	
Samples	8–48 (processed in batches of 8)
Library prep time	8 hours +/- 2 hours for up to 24 samples
DNA input	
Quantity	1 – 4 µg in Elution buffer, TE buffer (pH 8.0, 0.1mM EDTA), or nuclease-free water
DNA size distribution	50% ≥30 kb on Femto Pulse system; DIN >7 on TapeStation system
Sample multiplexing	
Vega™ system	Up to 48 samples
Revo™ system with SPRQ	Up to 48 samples

© 2025 PacBio. All rights reserved. Research use only. Not for use in diagnostic procedures. PN 103-740-700 REV01 SEP2025 **PacBio**

**Procedure & checklist – Generating PureTarget libraries with PureTarget kit 24 – Manual protocol (103-740-700)**

**Procedure & checklist – Generating PureTarget libraries with PureTarget kit 96 – Automation protocol (103-740-800)**

**Guide & overview – Automated PureTarget kit 96 for the Hamilton NGS STAR MOA system (103-740-900)**

**Application-specific technical overview**

**PacBio**

### Technical overview PureTarget library preparation using PureTarget kit 2.0

Vega system ICS v1.1.3  
Revo system ICS v13.3  
SMRT Link v25.3

PN 103-742-600 Rev 01 | September 2025

**Example 96-plex PureTarget repeat expansion panel library sequencing performance (Revo system + SPRQ chemistry)**

96-plex PureTarget repeat expansion panel example data for human whole blood DNA samples\*

Raw Data Report	HiFi Read Length	Target Enrichment Summary Metrics
Mean Expansion Read Length: 81.3 kb	4.8 kb Mean HiFi Read Length	Total Reads: 11,612,262,440
95%: 100%	141 reads (0.3%) are longer than 4.0 kb	Total Reads: 2,586,700
90%: 100%		Median Read Length: 4,102 bp
85%: 100%		Median Read Quality: Q48
80%: 100%		Median Read Quality: Q48
75%: 100%		Target Regions: 36
70%: 100%		Control Regions: 95

© 2025 PacBio. All rights reserved. Research use only. Not for use in diagnostic procedures. PN 103-742-600 REV01 SEP2025 **PacBio**

**Technical Overview – PureTarget library preparation using PureTarget kit 2.0 (103-742-600)**

Technical overview presentations describe workflow details for constructing PacBio SMRTbell libraries for specific applications. Example sequencing performance data for a given application are also summarized.

**Library preparation, sequencing & analysis**

**gDNA sample extraction & QC (Nanobind kit)**

1 – 2 µg DNA per sample  
Genome quality number (30 kb) ≥5.0  
Use Nanobind PanDNA kit for blood extraction

**PureTarget library preparation (Manual or automated workflow)**

PureTarget repeat expansion / carrier / custom panel options available  
PureTarget kit 24 (manual)  
PureTarget kit 96 (automated)

**SMRT sequencing (Vega & Revo systems)**

Vega system: Up to 48-plex  
Revo system + SPRQ: Up to 96-plex

**Data analysis (SMRT Link + PacBio variant calling tools)**


SMRT Link support for coverage QC analysis and repeat expansion genotyping



# PureTarget sequencing method overview

# CRISPR-Cas9 technology overview

CRISPR-Cas9 system comprises a guide RNA (gRNA or sgRNA) and Cas9 nuclease, which together form a ribonucleoprotein (RNP) complex that can introduce a site-specific double-strand break in DNA<sup>1,2</sup>




PureTarget uses the **CRISPR-Cas9 system** to generate targeted **native DNA libraries** for HiFi sequencing. This amplification-free approach retains epigenetic signals and has no PCR artifacts or errors.

## Guide RNA configurations


### Two-component guide RNA (gRNA)

Annealed [crRNA : tracrRNA] complex



### Single guide RNA (sgRNA)

Synthetic crRNA-tracrRNA single-molecule hybrid construct



sgRNA is used in PureTarget library prep workflow

<sup>1</sup> Image modified from: <https://horizondiscovery.com/en/applications/gene-editing>

<sup>2</sup> CRISPR (= *clustered regularly interspaced short palindromic repeats*) is a family of DNA sequences found in the genomes of prokaryotic organisms such as bacteria that play a role in the anti-viral defense system of these organisms.

# CRISPR-Cas9 technology overview (cont.)

CRISPR-Cas9 system comprises a guide RNA (gRNA or sgRNA) and Cas9 nuclease, which together form a ribonucleoprotein (RNP) complex that can introduce a site-specific double-strand break in DNA<sup>1,2</sup>





Illustration of **Cas9 nuclease** (blue), programmed by the **tracrRNA** (violet) : **crRNA** (olive) complex (= **guide RNA**) cutting both strands of genomic DNA 5' of the protospacer-adjacent motif (**PAM**) (red).

## Guide RNA configurations


### Two-component guide RNA (gRNA)

Annealed [crRNA : tracrRNA] complex



### Single guide RNA (sgRNA)

Synthetic crRNA-tracrRNA single-molecule hybrid construct




sgRNA is used in PureTarget library prep workflow

<sup>1</sup> Image modified from: <https://horizondiscovery.com/en/applications/gene-editing>

<sup>2</sup> CRISPR (= clustered regularly interspaced short palindromic repeats) is a family of DNA sequences found in the genomes of prokaryotic organisms such as bacteria that play a role in the anti-viral defense system of these organisms.

# How CRISPR-Cas9 is used in PureTarget native DNA library prep workflow


For PureTarget HiFi sequencing, CRISPR-Cas9 system is used to selectively enrich for regions of interest in a genomic sample prior to long-read sequencing of native DNA molecules<sup>1</sup>



## Excising a target region with two flanking gRNA-Cas9 complexes

- Cas9 nuclease, in close association with appropriate guide RNA (gRNA) oligonucleotides, identifies and then **cleaves a specific recognition site on each side** of the target region
- Different regions of interest can be **simultaneously targeted** in a single CRISPR-Cas9 digestion reaction by using multiple sets of gRNA pairs

## Target enrichment with CRISPR-Cas9



<sup>1</sup> Tsai, Y. C., et al. (2022). Multiplex CRISPR/Cas9-Guided No-Amp targeted sequencing panel for spinocerebellar ataxia repeat expansions. In Genomic Structural Variants in Nervous System Disorders (pp. 95-120). New York, NY: Springer US.

<sup>2</sup> **Note:** PureTarget guide RNAs are designed to cut a 5-kb region in the human reference genome; however, some DNA fragment sizes that are sequenced may be much larger with repeat expansion alleles.

# PureTarget 2.0 consumables for manual and automated library prep workflows

PureTarget 2.0 panel products enable flexibility in content and scale



## Automated PureTarget library prep

## HiFi Sequencing

### PureTarget kit 96 bundle (103-708-000)<sup>1</sup>



PureTarget kit 2.0  
(103-632-900)



PureTarget cleanup kit  
(103-708-100)<sup>2</sup>



HiFi plex prep kit 96  
(103-122-800)

SMRTbell cleanup beads  
(102-158-300) [2 x 10 mL]

Elution buffer  
(101-633-500) [1 x 50 mL]

### Panel gRNA reagents



PureTarget repeat expansion panel 2.0  
(103-633-100); or

PureTarget carrier panel  
(103-633-200); or

PureTarget control panel  
(103-633-300)

### Indexed SMRTbell adapters



SMRT adapter index plate 96A/B/C/D  
(102-009-200)

### PacBio long-read systems



Revio system + SPRQ  
96-plex

*With automation, a single Revio system can process 96,000 samples per year with full system utilization using PureTarget kit 96*



## Manual PureTarget library prep

## HiFi Sequencing

### PureTarget kit 24 bundle (103-707-900)



PureTarget kit 2.0  
(103-632-900)



PureTarget cleanup kit  
(103-708-100)<sup>2</sup>



SMRTbell prep kit 3.0  
(102-141-700)

SMRTbell cleanup beads  
(102-158-300) [1 x 10 mL]

Elution buffer  
(101-633-500) [1 x 50 mL]

### Panel gRNA reagents



PureTarget repeat expansion panel 2.0  
(103-633-100); or

PureTarget carrier panel  
(103-633-200); or

PureTarget control panel  
(103-633-300)

### Indexed SMRTbell adapters



SMRT adapter index plate 96A/B/C/D  
(102-009-200)

### PacBio long-read systems



Revio system + SPRQ  
8- to 48-plex




Vega system  
8- to 48-plex

*Vega system can process ~10,000 samples a year using manual workflow at 48-plex*


# PureTarget 2.0 product portfolio & workflow overview

Updated product portfolio supports PureTarget repeat expansion panel 2.0, PureTarget carrier screening & PureTarget custom [ manual + automated ] workflows



# PureTarget 2.0 product portfolio & workflow overview (cont.)


Updated product portfolio supports PureTarget repeat expansion panel 2.0, PureTarget carrier screening & PureTarget custom [ manual + automated ] workflows



<sup>1</sup> **Note:** Human saliva samples are only supported for manual library prep workflows using PureTarget kit 24.  
<sup>2</sup> PureTarget kit 2.0 does not include any guide RNA mix.

# PureTarget 2.0 product portfolio & workflow overview (cont.)

Updated product portfolio supports PureTarget repeat expansion panel v2, PureTarget carrier screening & PureTarget custom [ manual + automated ] workflows



<sup>1</sup> **Note:** Human saliva samples are only supported for **manual** library prep workflows using PureTarget kit 24.

<sup>2</sup> PureTarget kit 2.0 does not include any guide RNA mix.

<sup>3</sup> PureTarget custom control panel contains a gRNA mix to target 3 positive control genes (1 X-linked + 2 autosomal); customers may additionally spike in custom gRNA mixes for custom panel designs.

# PureTarget manual library preparation procedure reference [ For up to 48-plex ]

Procedure & checklist – Generating PureTarget libraries with PureTarget kit 24 – Manual protocol ([103-740-700](#)) describes the manual workflow for generating up to 48-plex PureTarget libraries using the **PureTarget kit 24** reagent bundle for sequencing on PacBio long-read systems

Overview	
Samples	8–48 (processed in batches of 8)
Library prep time	8 hours +/- 2 hours for up to 24 samples
DNA input	
Quantity	1 – 4 µg in Elution buffer, TE buffer (pH 8.0, 0.1mM EDTA), or nuclease-free water
DNA size distribution	50% ≥30 kb on Femto Pulse system; DIN >7 on TapeStation system
Sample multiplexing	
Vega™ system	Up to 48 samples
Revio® system with SPRQ	Up to 48 samples



**PureTarget library template (~4 – 5 kb)**  
Expanded alleles can be ≥20 kb<sup>1</sup>

PacBio

Generating PureTarget™ libraries with PureTarget kit 24 – manual protocol

Procedure & checklist

---

**Before you begin**

This procedure describes the workflow for manually generating PureTarget libraries.

Overview	
Samples	8–48 (processed in batches of 8)
Library prep time	8 hours +/- 2 hours for up to 24 samples
DNA input	
Quantity	1 – 4 µg in Elution buffer, TE buffer (pH 8.0, 0.1mM EDTA), or nuclease-free water
DNA size distribution	50% ≥30 kb on Femto Pulse system; DIN >7 on TapeStation system
Sample multiplexing	
Vega™ system	Up to 48 samples
Revio® system with SPRQ	Up to 48 samples

© 2025 PacBio. All rights reserved. Research use only. Not for use in diagnostic procedures.  
PN 103-740-700 REV01 SEP2025

PacBio

**Note:** Procedure [103-740-700](#) includes instructions for manual **PureTarget library construction workflow** *and* sequencing preparation (**ABC<sup>2</sup>**) workflow


<sup>1</sup> For normal alleles, resulting sequences are ~4–5 kb in length but reads for expanded alleles may be longer. **Note:** Sequencing results for PureTarget samples with large expansions indicate that it is possible to span repeats up to 35 kb in length in a single HiFi read.

<sup>2</sup> For primer annealing, polymerase binding & complex cleanup (ABC) and final dilution steps, follow sample setup instructions for PureTarget libraries in **Procedure & checklist – Generating PureTarget libraries with PureTarget kit 24 – Manual protocol** ([103-740-700](#)) – Do not use SMRT Link Sample Setup.

# PureTarget automated library preparation procedure reference [ For 96-plex ]

Procedure & checklist – Generating PureTarget libraries with PureTarget kit 96 – Automation protocol ([103-740-800](#)) describes the automated workflow for generating 96-plex PureTarget libraries using the **PureTarget kit 96** reagent bundle<sup>1</sup> for sequencing on PacBio long-read systems

Overview	
Samples	96
Library prep time	16 hours with the Hamilton NGS STAR MOA system
DNA input	
Quantity	1 - 1.5 µg in Elution buffer, TE buffer (pH 8, 0.1mM EDTA), or nuclease-free water
DNA size distribution	50% ≥30 kb on Femto Pulse system; DIN >7 on TapeStation system
Sample multiplexing	
Revio® system with SPRQ™ chemistry	96 samples



**PureTarget library template (~4 – 5 kb)**  
Expanded alleles can be ≥20 kb<sup>2</sup>

PacBio

Generating PureTarget™ libraries with PureTarget kit 96 – automation protocol

Procedure & checklist

**Before you begin**

This procedure describes the workflow for generating PureTarget repeat expansion or carrier libraries.

Overview	
Samples	96
Library prep time	16 hours with the Hamilton NGS STAR MOA system
DNA input	
Quantity	1 - 1.5 µg in Elution buffer, TE buffer (pH 8, 0.1mM EDTA), or nuclease-free water
DNA size distribution	50% ≥30 kb on Femto Pulse system; DIN >7 on TapeStation system
Sample multiplexing	
Revio® system with SPRQ™ chemistry	96 samples

**Note:** Procedure [103-740-800](#) includes instructions for automated **PureTarget library construction workflow** *and* sequencing preparation (**ABC<sup>3</sup>**) workflow

<sup>1</sup> PureTarget kit 96 bundle (103-708-000) is intended to support automated PureTarget library preparation workflows and is not intended to support manual workflows.

<sup>2</sup> For normal alleles, resulting sequences are ~4–5 kb in length but reads for expanded alleles may be longer. **Note:** Sequencing results for PureTarget samples with large expansions indicate that it is possible to span repeats up to 35 kb in length in a single HiFi read.

<sup>3</sup> For primer annealing, polymerase binding & complex cleanup (ABC) and final dilution steps, follow sample setup instructions for PureTarget libraries in *Procedure & checklist – Generating PureTarget libraries with PureTarget kit 96 – Automation protocol* ([103-740-800](#)) – Do not use SMRT Link Sample Setup.

# PureTarget automated 96-plex library preparation guide & overview reference for Hamilton NGS STAR MOA system

Guide & overview – Automated PureTarget kit 96 for the Hamilton NGS STAR MOA system ([103-740-900](https://www.pacbio.com/resources/technical_documents/103-740-900)) describes the automated workflow for generating 96-plex PureTarget libraries using the **Hamilton NGS STAR MOA liquid handling system** for sequencing on PacBio long-read systems<sup>1</sup>

Overview	
Samples	96
Workflow step	Automation time
gDNA Repair & SMRTbell cleanup (Optional*)	2.5 hours
Dephosphorylation	1 hour
Cas9 digestion & SMRTbell cleanup	2 hours
dA-tail	0.5 hours
Adapter ligation, pool & SMRTbell cleanup	2 hours
Nuclease treatment 1, pool & SMRTbell cleanup	2 hours
Nuclease treatment 2 & SMRTbell cleanup	2 hours
PureTarget Library cleanup	1 hour
ABC	1 hour
Average total time	< 16 hours (including off-deck prep)
DNA input	
DNA input Quantity	1.0–1.5 µg in Elution buffer, TE buffer (pH 8, 0.1mM EDTA), or nuclease-free water
DNA size distribution	50% ≥30 kb on Femto Pulse system; DIN >7 on TapeStation system



Hamilton NGS STAR MOA system



## Automated PureTarget™ kit 96 for the Hamilton NGS STAR MOA system

[Guide & overview](#)

© 2025 PacBio. All rights reserved. Research use only. Not for use in diagnostic procedures.

**Note:** This Hamilton NGS STAR MOA system method supports automated PureTarget library preparation for a **single batch of 96 samples – batch sizes less than 96 are not supported**

# PureTarget 2.0 consumables for manual and automated library prep workflows

PureTarget 2.0 panel products enable flexibility in content and scale



## Automated PureTarget library prep

## Sequencing

### PureTarget kit 96 bundle (103-708-000)<sup>1</sup>

### Panel gRNA reagents

### Indexed SMRTbell adapters

### PacBio long-read systems



PureTarget kit 2.0  
(103-632-900)



PureTarget cleanup kit  
(103-708-100)



HiFi plex prep kit 96  
(103-122-800)

SMRTbell cleanup beads  
(102-158-300) [2 x 10 mL]

Elution buffer  
(101-633-500) [1 x 50 mL]



PureTarget repeat expansion panel 2.0  
(103-633-100); *or*

PureTarget carrier panel  
(103-633-200); *or*

PureTarget control panel  
(103-633-300)



SMRT adapter index plate 96A/B/C/D  
(102-009-200)



Revio system + SPRQ  
96-plex



## Manual PureTarget library prep

## Sequencing

### PureTarget kit 24 bundle (103-707-900)

### Panel gRNA reagents

### Indexed SMRTbell adapters

### PacBio long-read systems



PureTarget kit 2.0  
(103-632-900)



PureTarget cleanup kit  
(103-708-100)



SMRTbell prep kit 3.0  
(102-141-700)

SMRTbell cleanup beads  
(102-158-300) [1 x 10 mL]

Elution buffer  
(101-633-500) [1 x 50 mL]



PureTarget repeat expansion panel 2.0  
(103-633-100); *or*

PureTarget carrier panel  
(103-633-200); *or*

PureTarget control panel  
(103-633-300)



SMRT adapter index plate 96A/B/C/D  
(102-009-200)



Revio system + SPRQ  
8- to 48-plex



Vega system  
8- to 48-plex

# PureTarget kit 24 bundle<sup>1,2</sup> and PureTarget cleanup kit components

PureTarget kit 24 bundle includes PureTarget kit 2.0 and SPK 3.0 to support manual PureTarget library prep

## PureTarget kit 2.0 (103-632-900)

- Reagents for PureTarget library preparation
- **Note:** This kit **does not** include guide RNA reagents<sup>3</sup>
- Supports 24 samples for manual PureTarget library prep workflows when bundled with SMRTbell prep kit 3.0

### PureTarget kit 2.0 components

Component	Description
1	 <b>Cas9 buffer</b> <ul style="list-style-type: none"> <li>• For Cas 9 digestion of gDNA</li> </ul>
2	 <b>Phosphatase</b> <ul style="list-style-type: none"> <li>• For blocking gDNA fragment ends</li> </ul>
3	 <b>Cas9 nuclease</b> <ul style="list-style-type: none"> <li>• For Cas 9 digestion of gDNA</li> </ul>
4	 <b>PureTarget nuclease mix</b> <ul style="list-style-type: none"> <li>• For nuclease Tx of SMRTbell libraries</li> </ul>
5	 <b>dA tail buffer</b> <ul style="list-style-type: none"> <li>• For A-tailing reaction</li> </ul>
6	 <b>dATP</b> <ul style="list-style-type: none"> <li>• For A-tailing reaction</li> </ul>
7	 <b>Taq DNA polymerase</b> <ul style="list-style-type: none"> <li>• For A-tailing reaction</li> </ul>

## SMRTbell prep kit 3.0 (102-141-700) + auxiliary components

- Reagents for SMRTbell library construction + cleanup
- Supports 24 samples for manual PureTarget library prep workflows

### SMRTbell prep kit 3.0 + auxiliary components

Component	Description
1	 <b>SMRTbell prep kit 3.0</b> <ul style="list-style-type: none"> <li>• Contains core reagents for SMRTbell template construction</li> </ul>
2	 <b>Low TE buffer</b> <ul style="list-style-type: none"> <li>• For DNA shearing and cleanup</li> </ul>
3	 <b>SMRTbell cleanup beads</b> <ul style="list-style-type: none"> <li>• For DNA cleanup</li> </ul>



## PureTarget cleanup kit (103-708-100)

- Reagents for PureTarget library cleanup
- For use with PureTarget kit 24 (manual) or PureTarget kit 96 (automated) library prep workflows
- Supports cleanup of 24 PureTarget library pools

### PureTarget cleanup kit components

Component	Description
1	 <b>PureTarget cleanup beads kit</b> <ul style="list-style-type: none"> <li>• For performing bead-based cleanup of final PureTarget libraries</li> </ul>
2	 <b>PureTarget cleanup buffer kit</b> <ul style="list-style-type: none"> <li>• Buffers needed for performing cleanup of final PureTarget libraries</li> </ul>



<sup>1</sup> PureTarget kit 2.0 (103-632-900) can only be purchased as part of the PureTarget kit 24 bundle product (103-707-900) or the PureTarget kit 96 bundle product (103-708-000).

<sup>2</sup> PureTarget kit 24 bundle (103-707-900) reagent volumes are **optimized for processing batches of 8 samples and reagent volumes may be insufficient to support batching of fewer than 8 samples.**

<sup>3</sup> Guide RNA panel reagents [ e.g., PureTarget repeat expansion panel 2.0 / PureTarget carrier panel / PureTarget control panel ], PureTarget cleanup kit, and barcoded adapters are sold separately.

# PureTarget kit 96 bundle<sup>1,2</sup> and PureTarget cleanup kit components

PureTarget kit 96 bundle includes PureTarget kit 2.0 and HPPK 96 to support automated PureTarget library prep

## PureTarget kit 2.0 (103-632-900)

- Reagents for PureTarget library preparation
- **Note:** This kit **does not** include guide RNA reagents<sup>3</sup>
- Supports 96 samples for automated PureTarget library prep workflows when bundled with HiFi plex prep kit 96

### PureTarget kit 2.0 components

Component	Description
1	 <b>Cas9 buffer</b> <ul style="list-style-type: none"> <li>• For Cas 9 digestion of gDNA</li> </ul>
2	 <b>Phosphatase</b> <ul style="list-style-type: none"> <li>• For blocking gDNA fragment ends</li> </ul>
3	 <b>Cas9 nuclease</b> <ul style="list-style-type: none"> <li>• For Cas 9 digestion of gDNA</li> </ul>
4	 <b>PureTarget nuclease mix</b> <ul style="list-style-type: none"> <li>• For nuclease Tx of SMRTbell libraries</li> </ul>
5	 <b>dA tail buffer</b> <ul style="list-style-type: none"> <li>• For A-tailing reaction</li> </ul>
6	 <b>dATP</b> <ul style="list-style-type: none"> <li>• For A-tailing reaction</li> </ul>
7	 <b>Taq DNA polymerase</b> <ul style="list-style-type: none"> <li>• For A-tailing reaction</li> </ul>

## HiFi plex prep kit 96 (103-122-800) + auxiliary components

- Reagents for SMRTbell library construction + cleanup
- Supports 96 samples for automated PureTarget library prep workflows

### HiFi plex prep kit 96 + auxiliary components

Component	Description
1	 <b>SMRTbell prep kit 3.0</b> <ul style="list-style-type: none"> <li>• Contains core reagents for SMRTbell template construction</li> </ul>
2	 <b>Low TE buffer</b> <ul style="list-style-type: none"> <li>• For DNA shearing and cleanup</li> </ul>
3	 <b>SMRTbell cleanup beads</b> <ul style="list-style-type: none"> <li>• For DNA cleanup</li> </ul>


## PureTarget cleanup kit (103-708-100)

- Reagents for PureTarget library cleanup
- For use with PureTarget kit 24 (manual) or PureTarget kit 96 (automated) library prep workflows
- Supports cleanup of 24 PureTarget library pools

### PureTarget cleanup kit components

Component	Description
1	 <b>PureTarget cleanup beads kit</b> <ul style="list-style-type: none"> <li>• For performing bead-based cleanup of final PureTarget libraries</li> </ul>
2	 <b>PureTarget cleanup buffer kit</b> <ul style="list-style-type: none"> <li>• Buffers needed for performing cleanup of final PureTarget libraries</li> </ul>

**Note:** PureTarget kit 96 bundle is intended to support **automated** PureTarget library preparation workflows and is **not** intended to support manual workflows



**PureTarget kit 96 bundle**  
103-708-000 (96 rxn)



**PureTarget cleanup kit**  
103-708-100 (24 rxn)

<sup>1</sup> PureTarget kit 2.0 (103-632-900) can only be purchased as part of the PureTarget kit 24 bundle product (103-707-900) or the PureTarget kit 96 bundle product (103-708-000).


<sup>2</sup> PureTarget kit 96 bundle (103-707-900) reagent volumes support processing of **batches of 96 samples in automated library prep workflows (batches <96 are not supported for automated workflows).**

<sup>3</sup> Guide RNA panel reagents [ e.g., PureTarget repeat expansion panel 2.0 / PureTarget carrier panel / PureTarget control panel ], PureTarget cleanup kit, and barcoded adapters are sold separately.

# Other recommended kits & consumables for PureTarget sample DNA extraction, sample multiplexing and HiFi sequencing


Ancillary kits must be purchased separately from PureTarget kit bundle (103-390-400)

## HMW DNA extraction



**Nanobind PanDNA kit**  
(103-260-000)

**Nanobind HT CBB kit**  
(102-762-700)



**Procedure & checklist**  
([Documentation](#))

- **Nanobind PanDNA kit [ 24 rxn ]** is recommended for manual DNA extraction from cultured cells, human blood and saliva for PureTarget panel applications<sup>1</sup>
- **Nanobind HT CBB kit [ 96 rxn ]** is recommended for automated DNA extraction from blood and cells

## Sample multiplexing<sup>2</sup>



**SMRTbell adapter index plate 96A/B/C/D**  
(102-009-200/102-547-800/102-547-900/102-548-000)

	1	2	3	4	5	6	7	8	9	10	11	12
A	bc2001	bc2009	bc2017	bc2025	bc2033	bc2041	bc2049	bc2057	bc2065	bc2073	bc2081	bc2089
B	bc2002	bc2010	bc2018	bc2026	bc2034	bc2042	bc2050	bc2058	bc2066	bc2074	bc2082	bc2090
C	bc2003	bc2011	bc2019	bc2027	bc2035	bc2043	bc2051	bc2059	bc2067	bc2075	bc2083	bc2091
D	bc2004	bc2012	bc2020	bc2028	bc2036	bc2044	bc2052	bc2060	bc2068	bc2076	bc2084	bc2092
E	bc2005	bc2013	bc2021	bc2029	bc2037	bc2045	bc2053	bc2061	bc2069	bc2077	bc2085	bc2093
F	bc2006	bc2014	bc2022	bc2030	bc2038	bc2046	bc2054	bc2062	bc2070	bc2078	bc2086	bc2094
G	bc2007	bc2015	bc2023	bc2031	bc2039	bc2047	bc2055	bc2063	bc2071	bc2079	bc2087	bc2095
H	bc2008	bc2016	bc2024	bc2032	bc2040	bc2048	bc2056	bc2064	bc2072	bc2080	bc2088	bc2096

- Each plate contains **96 indexed SMRTbell adapters** in plate format (1 sample per index)
- E.g., plate 96A includes indexes bc2001–bc2096

## HiFi sequencing



Revio system + SPRQ chemistry supports **8- to 96-plex** PureTarget sample multiplexing



Vega system supports **8- to 48-plex** PureTarget sample multiplexing

- **Revio SPRQ polymerase kit** supports up to **12 Revio SMRT Cells** for PureTarget libraries
- **Vega polymerase kit** supports up to **4 Vega SMRT Cells** for PureTarget libraries

<sup>1</sup> For genomic DNA extraction from blood, we recommend using the **red blood cell lysis protocol** described in *Procedure & checklist – Extracting HMW DNA from human whole blood with RBC lysis using the Nanobind PanDNA kit* (103-377-500) or using the whole blood extraction protocol described in *Procedure & checklist – Extracting HMW DNA from human whole blood using Nanobind kits* (102-573-500).

<sup>2</sup> **Note:** PureTarget automated library preparation procedure (103-740-800) supports up to 96-plex sample multiplexing through use of 96 different SMRTbell indexed adapters.

# PureTarget supported use cases and experimental design considerations

PureTarget repeat expansion panel 2.0 features 38 targets covering a variety of commonly screened tandem repeats associated with neurological disorders

Disease	Targets
Spinocerebellar ataxia (SCA)	<i>ATN1, ATXN1, ATXN2, ATXN3, ATXN7, ATXN8, ATXN10, CACNA1A, PPP2R2B, TBP, BEAN1, DAB1, FGF14, NOP56, ZFHX3</i>
Fragile-X disease (FXS)	<i>FMR1</i>
Fragile X syndrome, FRAXE type	<i>AFF2</i>
Intellectual disability associated with fragile site FRA2A	<i>AFF3</i>
Frontotemporal dementia (FTD), amyotrophic lateral sclerosis (ALS)	<i>C9orf72</i>
Friedreich ataxia (FRDA)	<i>FXN</i>
Cerebellar ataxia, neuropathy, and vestibular areflexia syndrome (CANVAS)	<i>RFC1</i>
Neuronal intranuclear inclusion disease, Alzheimer disease and parkinsonism phenotype (NIID)	<i>NOTCH2NLC</i>
Myotonic dystrophy (DM)	<i>DMPK, CNBP</i>
Huntington disease (HD)	<i>HTT</i>
Huntington's disease-like type2 (HDL2)	<i>JPH3</i>
Fuchs endothelial corneal dystrophy 3 (FECD3)	<i>TCF4</i>
Kennedy Disease, Spinal and bulbar muscular atrophy, (SBMA)	<i>AR</i>
Oculopharyngeal muscular dystrophy (OPMD)	<i>PABPN1</i>
Oculopharyngodistal myopathy (OPDM)	<i>ABCD3, GIPC1, LRP12, RILPL1</i>
Syndactyly (SD5)	<i>HOXD13</i>
Congenital central hypoventilation syndrome (CCHS)	<i>PHOX2B</i>
Creutzfeldt-Jakob disease (CJD)	<i>PRNP</i>
Progressive Myoclonic Epilepsy Type 1 (EPM1) Unverricht-Lundborg Disease (ULD)	<i>CSTB</i>
Familial adult myoclonic epilepsy type 1 (FAME)	<i>SAMD12</i>

18 new targets

## PureTarget repeat expansion panel 2.0 enables genotyping of critical pathogenic repeat expansion loci at scale

- Panel targets 38 repeat expansion loci and captures ~2 kb upstream and downstream of the repeat (total panel size = ~163 kb)
- For normal alleles, resulting sequences are 4–5 kb in length but reads for expanded alleles may be longer
- Sequencing results for samples with large expansions indicate that it is possible to span repeats up to 35 kb in length in a single read

**PacBio**

Application note

**Comprehensive genotyping with the PureTarget repeat expansion panel and HiFi sequencing**

**Introduction**

Tandem repeats are regions of the genome consisting of repetitive units of specific DNA sequences. These regions are hypermutable; they can increase in length across generations and may have variable lengths within somatic tissues of an individual (Paulsen 2018). Many tandem repeats become pathogenic when they exceed a length threshold, which varies from gene to gene, resulting in mutations called repeat expansions (RE). (Bañez et al., 2022). Repeat expansions have been linked to dozens of diseases and cancer, most notably neurodegenerative disorders like Huntington's disease, Fragile-X disorder, spinocerebellar ataxia, and myotonic dystrophy (Bisgaard & Møller, 2021). Disease severity and age of onset of these conditions are often associated with their repeat length (Bañez et al., 2022).

Repeat expansion disease alleles are common, with a recent study estimating more than 1 in 200 individuals harboring disease alleles in their genome (Bañez et al., 2022). Despite this, these regions are challenging to characterize and as such, the majority of patients with rare neurological diseases remain undiagnosed (Bañez et al., 2022). Repeat expansions have historically been

profiled with Southern blotting and PCR-based assays (Tarkenton, 2003) and more recently short-read bioinformatic methods (Dobson et al., 2020; Dabrowski et al., 2022), but these methods cannot accurately size long expansions or simultaneously profile methylation, and are limited in their ability to characterize the base-level sequence of expanded repeats. Given that many repeat expansions are larger than the typical length of short reads, comprehensive genotyping of REs requires a high-throughput, long-read sequencing approach that can reliably manage the high structure variability of these regions.

Powered by the exceptional accuracy of HiFi sequencing and the Tandem Repeat Genotyping Tool (TRGT), the PacBio® PureTarget™ repeat expansion panel offers a more comprehensive genotyping, including methylation, for 38 of the most important tandem repeats impacting human health. This application note demonstrates the performance of PureTarget for repeat expansion sequencing, and presents the PureTarget panel as a scalable, more comprehensive solution for profiling repeat expansions, compared to legacy genotyping and next-generation sequencing methods.

© 2023 PacBio. All rights reserved. Research use only. Not for use in diagnostic procedures. 102-326-614 REV09 SEP2023

PacBio

PureTarget repeat expansion panel 2.0 (103-633-100)

Pooled guide RNAs for PureTarget repeat expansion panel 2.0 (38 genes). Supports 24 samples with manual PureTarget 24 kit and 96 samples with automated PureTarget 96 kit.

**Application note – Comprehensive genotyping with the PureTarget repeat expansion panel and HiFi sequencing (102-326-614)<sup>1</sup>**


# PureTarget supported use cases and experimental design considerations (cont.)

PureTarget carrier panel enables consolidated analysis of 12 challenging genes commonly genotyped with legacy non-NGS assays

Disease	Targets	Variant type
Hemophilia	F8	Inversions
Friedreich ataxia (FRDA)	FXN	Repeat Expansion
Fragile-X disease (FXS)	FMR1	Repeat Expansion
Congenital adrenal hyperplasia	CYP21A2	Small variants and copy number
Classical-like Ehlers-Danlos syndrome	TNXB	Small variants and copy number
Alpha thalassemia	HBA1/2	Deletions
Gaucher disease	GBA	Small variants and copy number
Spinal muscular atrophy	SMN1/2	Small variants and copy number
Early-infantile epileptic encephalopathy (EIEE1) and Partington syndrome (PRTS)	ARX	Repeat Expansion (N=2)
Beta thalassemia	HBB	Small variants
X-linked retinitis pigmentosa	RPGR	Small variants
Fragile X syndrome, FRAAXE type	AFF2	Repeat Expansion

PureTarget carrier panel enables comprehensive, scalable screening of challenging carrier genes that are difficult to resolve with traditional technologies

- Panel targets 12 challenging genes and captures ~2 kb upstream and downstream of the gene (total panel size = ~141 kb)
- PureTarget carrier panel is a highly accurate and streamlined assay that can be run alongside routine short read panels to capture all expanded carrier screening targets.



**PureTarget carrier  
Screening panel  
(103-633-200)**

Pooled guide RNAs for PureTarget carrier panel (12 genes). Supports 24 samples with manual PureTarget 24 kit and 96 samples with automated PureTarget 96 kit.

**Application note – Consolidate challenging genes with PureTarget carrier screen panel (102-326-653)<sup>1</sup>**

# PureTarget supported use cases and experimental design considerations (cont.)

**PureTarget custom panel** option enables users to design target-specific panels to analyze biological targets of interest

**Technical note – A practical guide to amplification-free custom PureTarget panels (102-326-652) provides guidance to customers on creating custom panels**

## Performance specifications

- Max and min panel size that has been tested
- Coverage for different fragment length

## Types of custom designs

- Removing targets (from existing panel designs)
- Adding targets (by spiking in new guide RNA reagents to a control panel or other existing panel designs)
- Tiling and phasing

## Designing guide RNAs

- Recommended design tools, vendors and products
- Success rate, common challenges

## Protocol modifications


- Specific steps to modify in protocol, reagent concentrations etc.
- Spiking in custom guides to control panel
- Protocol modifications for tiled designs

## Analysis recommendations

- SMRT Analysis target enrichment for coverage QC assessment
- Recommend PacBio GitHub tools for genotyping of different types of targets

## Note about using custom panel designs

- PureTarget custom panels should include the PureTarget control panel, which contains a guide RNA mix to target 3 positive control genes (1 X-linked + 2 autosomal)
  - Users may additionally spike in custom gRNA mixes for their own custom panel designs
- In all cases, **we recommend first demonstrating success on the PureTarget repeat expansion panel or PureTarget carrier panel using supported sample types** before designing your own custom panels



## PureTarget carrier Screening panel (103-633-200)

Pooled guide RNAs for PureTarget control panel (3 genes). Supports 24 samples with manual PureTarget 24 kit and 96 samples with automated PureTarget 96 kit.

**Technical note – A practical guide to amplification-free PureTarget custom panels (102-326-652)<sup>1</sup>**

# PureTarget 2.0 performance specifications for supported sample types

PureTarget enables comprehensive characterization of repeat expansions and other difficult genes at scale<sup>1</sup>

Parameter	Specification	Notes
Target gene panel size	<b>PureTarget repeat expansion panel 2.0:</b> 38 genes (163 kb) <b>PureTarget carrier panel:</b> 12 genes (141 kb)	<ul style="list-style-type: none"> <li>See <i>PureTarget Brochure</i> (<a href="#">102-326-609</a>) for list of target genes included in PureTarget repeat expansion panel 2.0 (103-633-100) and PureTarget carrier panel (103-633-200)</li> </ul>
DNA input amount	<b>Automated library prep with PureTarget kit 96:</b> 1–1.5 µg per sample	<ul style="list-style-type: none"> <li>Max. total DNA input = 150 µg per Revio SMRT Cell</li> </ul>
	<b>Manual library prep with PureTarget kit 24:</b> 1–4 µg per sample	<ul style="list-style-type: none"> <li>Max. total DNA input = 50 µg per Revio SMRT Cell or 100 µg per Vega SMRT Cell</li> </ul>
DNA input quality	GQN30kb > 5	<ul style="list-style-type: none"> <li>50% of mass of DNA molecules longer than 30 kb as measured on Femto Pulse (Agilent)</li> </ul>
Target coverage	<b>Mean target coverage:</b> 100-fold or greater per 1 µg DNA per sample	<ul style="list-style-type: none"> <li>Mean coverage per 1 µg of input DNA from supported sample types<sup>1</sup> at max. sample multiplexing level (Revio SPRQ = 96, Vega = 48)</li> <li>Higher coverage is possible with higher DNA input amounts and lower sample multiplexing</li> </ul>
	<b>Minimum target coverage:</b> 20-fold per sample	<ul style="list-style-type: none"> <li>Minimum coverage per 1 µg of input DNA from supported sample types<sup>1</sup></li> </ul>
Sample multiplexing	<b>Vega system:</b> Up to 48-plex	<ul style="list-style-type: none"> <li>For up to 48-plex: PureTarget kit 24 bundle (103-707-900) supports smaller batches in multiples of 8 samples</li> </ul>
	<b>Revio + SPRQ system:</b> Up to 96-plex	<ul style="list-style-type: none"> <li>For 96-plex: PureTarget kit 96 bundle (103-708-000) supports single batches of 96 samples</li> </ul>
Library insert size	4 – 5 kb	<ul style="list-style-type: none"> <li>Inserts with expanded alleles will be longer</li> </ul>
Methylation	5mC in CpG sites detected	<ul style="list-style-type: none"> <li>Methylation probabilities for CpG sites are encoded in BAM output file</li> </ul>
PureTarget library prep protocol (Proc. & Checklist)	<b>Generating PureTarget libraries with PureTarget kit 24 – Manual protocol</b> ( <a href="#">103-740-700</a> )	<ul style="list-style-type: none"> <li>Requires PureTarget kit 24 bundle (103-707-900)</li> </ul>
	<b>Generating PureTarget libraries with PureTarget kit 96 – Automation protocol</b> ( <a href="#">103-740-800</a> )	<ul style="list-style-type: none"> <li>Requires PureTarget kit 96 bundle (103-708-000)</li> </ul>
PureTarget data analysis workflow	<b>PureTarget repeat expansion panel 2.0</b> → SMRT Link PureTarget repeat expansion analysis application	<ul style="list-style-type: none"> <li>Produces target coverage summary QC stats and uses tandem repeat genotyping tool (TRGT) for variant calling and for visualization</li> </ul>
	<b>PureTarget carrier panel</b> → SMRT Link Target enrichment analysis application & PureTarget carrier pipeline (PTCP) analysis	<ul style="list-style-type: none"> <li>SMRT Link Target enrichment analysis produces target coverage QC statistics</li> <li>PTCP analysis workflow available through GitHub or PacBio Compatible partners enables variant calling for tandem repeat regions and hard genes</li> </ul>
	<b>PureTarget custom panel</b> → SMRT Link Target enrichment analysis application & PacBio GitHub software tools for variant calling	<ul style="list-style-type: none"> <li>SMRT Link Target enrichment analysis produces target coverage QC statistics</li> <li>Can use other PacBio GitHub software tools to perform variant calling</li> </ul>

<sup>1</sup> Optimal performance is obtained using PacBio Nanobind extraction kits with **human blood or cell line samples**. Refer to *Application note – Comprehensive genotyping with the PureTarget repeat expansion panel and HiFi sequencing* ([102-326-614](#)) for the most up-to-date list of supported DNA extraction kits and sample types for PureTarget sequencing applications. Individual target coverage is lower for longer target regions or expanded alleles. The same library loaded on Revio SPRQ system will typically give higher coverage than Vega system.



# PureTarget library preparation workflow details

# Procedure & checklist – Generating PureTarget libraries with PureTarget kit 24 – Manual protocol ([103-740-700](#))

Procedure & checklist – Generating PureTarget libraries with PureTarget kit 24 – Manual protocol ([103-740-700](#)) describes the workflow for generating up to 48-plex PureTarget libraries using the **PureTarget kit 24** reagent bundle for sequencing on PacBio long-read systems

## Procedure & checklist contents

1. **DNA input QC recommendations** and **general best practices** for reagent & sample handling.
2. Enzymatic workflow steps for performing **targeted Cas9 digestion** of input genomic DNA samples.
3. Enzymatic workflow steps for **PureTarget SMRTbell library construction** from Cas9-digested gDNA samples.
4. Workflow steps for **final cleanup** of PureTarget SMRTbell libraries using PureTarget cleanup beads
5. Workflow steps for **sample setup ABC<sup>1</sup>** and final dilution to prepare samples for sequencing

**Note:** Procedure 103-740-700 includes instructions for up to 48-plex PureTarget library construction using a **manual workflow**.

For instructions on high-throughput 96-plex PureTarget library construction using an **automated workflow**, refer to **Procedure & checklist – Generating PureTarget libraries with PureTarget kit 96 – Automation protocol** ([103-740-800](#))

Generating PureTarget™ libraries with PureTarget kit 24 – manual protocol **PacBio**

Procedure & checklist

### Before you begin

This procedure describes the workflow for manually generating PureTarget libraries.

Overview	
Samples	8–48 (processed in batches of 8)
Library prep time	8 hours +/- 2 hours for up to 24 samples
DNA input	
Quantity	1 – 4 µg in Elution buffer, TE buffer (pH 8.0, 0.1mM EDTA), or nuclease-free water
DNA size distribution	50% ≥30 kb on Femto Pulse system; DIN >7 on TapeStation system
Sample multiplexing	
Vega™ system	Up to 48 samples
Revio® system with SPRQ	Up to 48 samples

© 2025 PacBio. All rights reserved. Research use only. Not for use in diagnostic procedures.  
PN 103-740-700 REV01 SEP2025 **PacBio**

PacBio [Documentation](#) ([103-740-700](#))

# Key PureTarget library preparation and sequencing workflow steps [ Manual ]

Manual PureTarget library construction steps can be completed within 1 day<sup>1</sup>

Extraction & QC

DNA sample extraction

2 – 2.5 hrs



Nanobind DNA extraction kit



Refer to Nanobind user guide documentation

1 Input gDNA QC

90 min



Qubit 4 fluorometer system (ThermoFisher Sci. Q33238)  
DNA sizing QC tool (Femto pulse system or other tool)

2 gDNA repair<sup>2</sup> & SMRTbell bead cleanup

2 hrs



PureTarget kit 2.0 (103-632-900)

Safe stop

3 Dephosphorylation

30 min



PureTarget cleanup kit (103-708-100)

4 CAS-9 digestion & SMRTbell bead cleanup

1.5 hrs



SMRTbell prep kit 3.0 (102-141-700)

Safe stop

5 dA tailing

10 min



SMRTbell adapter index plate 96A/B/C/D (102-009-200)

6 Adapter ligation & SMRTbell bead cleanup

1 hr

Safe stop

7 Nuclease treatment & SMRTbell bead cleanup

2 hrs



PureTarget repeat expansion panel 2.0 (103-633-100) *or*  
PureTarget carrier panel (103-633-200) *or*  
PureTarget control panel (103-633-300)

8 PureTarget library cleanup

30 min

Safe stop

9 Annealing, binding & cleanup (ABC)

1 hr



Revio SPRQ polymerase kit (103-496-900) *or*  
Vega polymerase kit (103-426-500)

Library construction

1 DAY

ABC

## Manual PureTarget protocol reference

[Procedure & checklist – Generating PureTarget libraries with PureTarget kit 24 – Manual protocol \(103-740-700\)](#)

PacBio

Generating PureTarget™ libraries with PureTarget kit 24 – manual protocol

Procedure & checklist

**Before you begin**

This procedure describes the workflow for manually generating PureTarget libraries.

Overview	
Samples	8–48 (processed in batches of 8)
Library prep time	8 hours +/- 2 hours for up to 24 samples
DNA input	
Quantity	1 – 4 µg in Elution buffer, TE buffer (pH 8.0, 0.1mM EDTA), or nuclease-free water
DNA size distribution	50% ≥30 kb on Femto Pulse system; DIN >7 on TapeStation system
Sample multiplexing	
Vega™ system	Up to 48 samples
Revio® system with SPRQ	Up to 48 samples

© 2025 PacBio. All rights reserved. Research use only. Not for use in diagnostic procedures.  
PN 103-740-700 REV01 SEP2025

PacBio




<sup>1</sup> Example workflow step times shown are for manually processing up to 24 PureTarget samples using PureTarget kit 24 reagent bundle.

<sup>2</sup> This step repairs input genomic DNA, which increases target coverage when using lower quality gDNA. For Nanobind-extracted blood (whole blood or RBC lysis method) gDNA, this step *can* be skipped..

# Key PureTarget library preparation and sequencing workflow steps [ Automated ]

Automated PureTarget library construction steps can be completed in 2 days<sup>1</sup>



Refer to Nanobind user guide documentation

## Automated PureTarget protocol reference

*Procedure & checklist – Generating PureTarget libraries with PureTarget kit 96 – Automation protocol (103-740-800)*

### Generating PureTarget™ libraries with PureTarget kit 96 – automation protocol



Procedure & checklist

#### Before you begin

This procedure describes the workflow for generating PureTarget repeat expansion or carrier libraries.

Overview	
Samples	96
Library prep time	16 hours with the Hamilton NGS STAR MOA system
DNA input	
Quantity	1 - 1.5 µg in Elution buffer, TE buffer (pH 8, 0.1mM EDTA), or nuclease-free water
DNA size distribution	50% ≥30 kb on Femto Pulse system; DIN >7 on TapeStation system
Sample multiplexing	
Revio® system with SPRQ™ chemistry	96 samples

© 2025 PacBio. All rights reserved. Research use only. Not for use in diagnostic procedures.  
PN 103-740-800 REV 01 SEP2025



<sup>1</sup> Example workflow step times shown are for automated processing 96 PureTarget samples using PureTarget kit 96 reagent bundle.  
<sup>2</sup> gDNA repair is required unless using Nanobind-extracted blood genomic DNA using the RBC lysis method or whole blood extraction..

# General best practices recommendations for preparing PureTarget libraries

## Supported input sample types and DNA extraction methods

PureTarget performance is sensitive to DNA quality, which can be impacted by DNA extraction method and sample type used

- It is recommended that users start with high-quality genomic DNA extracted with PacBio Nanobind extraction kits
- When using sample types and extraction methods other than the above, we recommend users:<sup>1</sup>
  - First, demonstrate success using supported sample types, starting with an 8-plex and increasing sample quantity thereafter
  - Introduce new sample types or extraction methods in limited numbers, for example, 3 or fewer new sample types in an 8-plex of otherwise controls

Table 1 (below). Officially supported samples types for PureTarget library preparation workflows.

Supported sample types	Supported DNA extraction methods
Human blood	Extracted by Nanobind RBC lysis PanDNA kit [ <a href="#">103-377-500</a> ], Nanobind Whole blood manual CBB kit/PanDNA kit [ <a href="#">102-573-500</a> ], Nanobind HT 200 µL protocol [ <a href="#">103-028-100</a> ], or Nanobind HT 1 mL protocol [ <a href="#">103-028-100</a> ]
Human B-lymphocyte cell lines	Extracted by Nanobind CBB kit [ <a href="#">103-394-500</a> / <a href="#">102-573-600</a> ], Nanobind PanDNA kit [ <a href="#">103-394-500</a> / <a href="#">102-573-600</a> ], or Nanobind CBB HT kit [ <a href="#">103-028-100</a> ]; Coriell cell line genomic DNA (GQN30kb ≥5)
Human saliva (PureTarget kit 24 only)	<p>Extracted by Nanobind PanDNA kit [ <a href="#">103-544-000</a> ], Nanobind CBB kit [ <a href="#">103-544-000</a> ], or Nanobind CBB HT kit</p> <p><b>NOTE:</b> Saliva DNA samples are only supported for manual library prep workflows using PureTarget kit 24 [ Up to 48-plex]<sup>2</sup></p> <ul style="list-style-type: none"> <li>A diluted AMPure cleanup is <b>required</b> before the start of library prep for this sample type. Please refer to Appendix A1 in <i>Procedure &amp; checklist – Generating PureTarget libraries with PureTarget kit 24 – Manual protocol</i> ( <a href="#">103-740-700</a> ) for details.</li> </ul>

Table 2 (right). Other demonstrated sample types and DNA extraction methods for PureTarget applications. Total sample multiplex level and total DNA input mass into library prep (across all samples) is shown for unsupported sample types and extraction kits listed in table. Demonstrated sample configurations shown in table resulted in successful library prep and sequencing but have not been technically validated. Results may vary. All samples are human-derived.

Other sample types	Demonstrated DNA extraction method	Sample multiplexing level [Total DNA input]
Blood	QIASymphony (QIAGEN)	Up to 96-plex [130 µg]
Blood	Puregene Blood kit (QIAGEN)	Up to 96 plex and [130 µg]
Blood	Biosystems MagMAX DNA Multi-Sample Ultra 2.0 Kit	Up to 32 plex and [40 µg]
Blood	Chemagic DNA blood kit (Revvity)	Up to 24 plex and [48 µg]
Blood	NucleoMag (Macherey-Nagel)	Up to 24-plex and [42 µg]
Blood	Monarch HMW DNA Extraction Kit for Cells & Blood (NEB)	Up to 8 plex and [20 µg]
Blood	QIAGEN Genomic-tips (QIAGEN)	Up to 8 plex and [16 µg]
Blood	FlexiGene DNA Whole Blood Kit (QIAGEN)	Up to 8 plex and [16 µg]
Blood	Bionano SP Blood and Cell Culture DNA Isolation Kit	Up to 8 plex and [16 µg]
Saliva	QIASymphony (QIAGEN)	Up to 48 plex and [75 µg]
Brain tissue	<a href="#">Nanobind DNA extraction from animal tissue protocol</a>	Up to 8-plex and [16 µg]
Skeletal muscle	<a href="#">Nanobind DNA extraction from animal tissue protocol</a>	Less than 8-plex and [<16 µg]
iPSCs	<a href="#">Nanobind DNA extraction from cultured adherent cells protocol</a>	Up to 16-plex and [36 µg]
Fibroblasts	<a href="#">Nanobind DNA extraction from cultured adherent cells protocol</a>	Up to 8-plex and [16 µg]
Myoblasts	<a href="#">Nanobind DNA extraction from cultured adherent cells protocol</a>	Less than 8-plex and [<16 µg]
iPSCs	QIAGEN Genomic-tips (QIAGEN)	Up 8-plex and [21 µg]
Fibroblasts	QIAGEN Genomic-tips (QIAGEN)	Up to 8-plex and [16 µg]
Corneal endothelial (CEC) cell culture	Bionano SP Blood and Cell Culture DNA Isolation Kit	Up to 8-plex and [16 µg]
Fibroblasts	Monarch HMW DNA Extraction Kit for Cells & Blood (NEB)	Up to 8-plex and [20 µg]
Brain tissue	Monarch HMW DNA Extraction Kit for Cells & Blood (NEB)	Less than 8-plex and [<16 µg]

<sup>1</sup> See *Application note – Comprehensive genotyping with PureTarget repeat expansion panel and HiFi sequencing* ( [102-326-614](#) ) for more information about samples that are officially supported or have been demonstrated to result in successful sequencing but have not been technically validated.

<sup>2</sup> Although 48-plex designs are supported for saliva samples, since saliva DNA is typically lower in quality and shorter in length, it is recommended to start at the gDNA repair step with 2 µg of (diluted AMPure-cleaned) saliva DNA per sample, with a maximum of 24-plex on Revio SPRQ or Vega chemistry for the best performance.

# General best practices recommendations for preparing PureTarget libraries (cont.)

## DNA sizing QC

- Agilent **Femto Pulse system**<sup>1</sup> is recommended for the accurate sizing of genomic DNA samples
  - Femto Pulse system enables simple, rapid sizing QC of genomic DNA and SMRTbell libraries, and conserves sample by using femtogram ranges of input DNA
    - Resolves fragments 1,300 bp to 165 kb using gDNA 165 kb Analysis kit (can resolve 100 – 6,000 bp using Ultra Sensitivity NGS kit)
    - Requires <1 ng of sample DNA
    - Can analyze up to 12 samples in <1.5 hrs
    - Outputs quality metrics such as Genomic Quality Number (GQN)<sup>2</sup> to quickly score integrity of HMW gDNA
- Alternative DNA sizing tools (e.g., TapeStation system) may be used if a Femto Pulse system is unavailable
  - However, **caution should be used when interpreting results from other tools that employ constant-field electrophoresis technology**
    - These technologies tend to **inflate the true size of the gDNA (or library)** and should only be used for qualitative assessment of whether an experiment was successful (e.g., intact library) rather than for accurate measurement of fragment size distributions



Femto Pulse system  
(Agilent Technologies)

## DNA quantification QC

- For DNA quantification QC, we recommend using a quantification assay specific for double-stranded DNA (dsDNA) such as the Qubit 1X dsDNA high sensitivity assay kit<sup>3</sup> (Thermo Fisher Scientific)
  - **Note:** We do not recommend quantification with UV-Vis Spectrophotometers (e.g. NanoDrop) that measure all nucleic acids in a sample. For example, measuring all nucleic acid will inflate the true concentration of gDNA in samples.



Qubit 4 fluorometer  
(Thermo Fisher Scientific)

<sup>1</sup> See [Product Note – HiFi WGS sequencing with the Agilent Femto Pulse system \(102-326-561\)](#) for more details.

<sup>2</sup> See [Application Note – Quality Metrics for Nucleic Acids with the Agilent Fragment Analyzer and Femto Pulse Systems \(Agilent 5994-0521EN\)](#)

<sup>3</sup> Alternatively, for high-throughput applications DNA quantification QC may be performed with a microplate reader using the Quant-iT 1X dsDNA high sensitivity assay kit (Thermo Fisher Scientific).

# General best practices recommendations for preparing PureTarget libraries (cont.)

## Recommended genomic DNA input amount and quality

### DNA input quality


- 50% or more of the DNA should be  $\geq 30$  kb for this PureTarget library prep protocol
  - This corresponds to a genome quality number (GQN) of 5.0 or higher at 30 kb based on the Agilent Femto Pulse system

DNA quality	Femto Pulse genome quality score <sup>1</sup>	Notes
DNA size distribution (Femto pulse system)	50% $\geq 30$ kb (GQN30kb $\geq 5.0$ )	<ul style="list-style-type: none"> <li>For low-quality human gDNA, it is still feasible to proceed but lower on-target coverage performance is expected<sup>2</sup></li> </ul>

- Note:** Removal of RNA with RNase is **required** for any DNA extraction method used
  - Failure to remove RNA may result in sequencing inhibition

Any degradation present should be due to shearing from extraction process and **not** from poor sample handling/storage or biochemical processes

(Right figure.) Example DNA sizing QC analysis of a high-quality HG002 human genomic DNA sample using a Femto Pulse system with Genomic DNA 165 kb kit.



<sup>1</sup> If using an Agilent TapeStation for input gDNA sizing QC, then a DIN value  $>7$  is recommended.

<sup>2</sup> **Important:** The HiFi yield and HiFi mean read length of a sequencing run are directly proportional to the quality of the genomic DNA input. To maximize yield and target coverage per SMRT Cell, start with high quality gDNA containing minimal DNA below 10 kb, and with  $>50\%$  mass over 30 kb. High quality gDNA will typically have a higher percent library recovery and HiFi sequencing yield.

# General best practices recommendations for preparing PureTarget libraries (cont.)

## Recommended genomic DNA input amount and quality

### Genomic DNA input amount required for manual library preparation using PureTarget kit 24

- Recommended mass of input gDNA is **2 µg per sample** to ensure there are sufficient gene copies to load and maximize sequencing coverage
  - Manual PureTarget library prep protocol supports input gDNA mass amounts in the range of 1–4 µg per sample.
- Recommend maximum total DNA mass per sequencing platform/chemistry is outlined in the following table across all multiplexed samples
  - For example, a multiplex of 48 samples on Revio (+SPRQ) should not exceed 1 µg per sample on average

Recommended genomic DNA input mass amounts across all multiplexed samples for manual PureTarget library preparation for Revio and Vega system sequencing platforms.

Minimum and maximum total gDNA mass for manual PureTarget library preparation using PureTarget kit 24		
Genomic DNA input for 1 SMRT Cell (8- to 48-plex)	Vega system chemistry	Revio system with SPRQ chemistry
Minimum	16 µg	16 µg
Maximum	100 µg	50 µg

### Genomic DNA input amount required for automated library preparation using PureTarget kit 96

- Automated PureTarget library prep protocol supports **1 – 1.5 µg per sample**
  - 1.3 – 1.5 µg per sample is highly recommended for samples going through the gDNA repair step, to ensure sufficient mass to achieve optimal target coverage.
  - **Maximum mass of 150 µg of gDNA input** (across 96 samples) is recommended for a single Revio SPRQ SMRT cell.

Recommended genomic DNA input mass amounts across all multiplexed samples for automated PureTarget library preparation for Revio system sequencing platform.

Minimum and maximum total gDNA mass for automated PureTarget library preparation using PureTarget kit 96	
Genomic DNA input for 1 SMRT Cell (96-plex)	Revio system with SPRQ chemistry
Minimum	96 µg
Maximum	150 µg

# General best practices recommendations for preparing PureTarget libraries (cont.)

## Recommended genomic DNA input amount and quality

### Expected PureTarget library construction yield for manual and automated workflows

- Overall PureTarget library construction yield is dependent on input gDNA quality and size
  - The recovery from input gDNA to completed SMRTbell library typically ranges between **0.02 – 0.2% for manual workflows** and **between 0.01 – 0.1% for automated workflows** (includes PureTarget library construction using PureTarget kit 2.0 and PureTarget library cleanup using PureTarget cleanup kit)

PureTarget library construction step	Manual PureTarget library prep [ Up to 48-plex ] <sup>1</sup>		Automated PureTarget library prep [ 96-plex ] <sup>2</sup>	
	DNA or SMRTbell step recovery	DNA or SMRTbell overall recovery	DNA or SMRTbell step recovery	DNA or SMRTbell overall recovery
Starting input genomic DNA	100%	<b>100%</b>	100%	<b>100%</b>
Post-gDNA repair & 1.2X SMRTbell bead cleanup	70 – 90%	<b>70 – 90%</b>	70 – 90%	<b>70 – 90%</b>
Post-CAS9 digestion & 1.2X SMRTbell bead cleanup	50 – 80%	<b>35 – 72%</b>	50 – 80%	<b>35 – 72%</b>
Post-adapter ligation & 1.0X SMRTbell bead cleanup	70 – 90%	<b>24 – 65%</b>	70 – 90%	<b>24 – 65%</b>
Post-nuclease treatment & 1.2X SMRTbell bead cleanup <sup>3</sup>	0.2 – 1%	<b>0.05 – 0.5%</b>	0.1 – 0.3%	<b>0.02 – 0.1%</b>
Post-PureTarget library cleanup using PureTarget cleanup kit	40 – 60%	<b>0.02 – 0.2%</b>	40 – 60%	<b>0.005 – 0.05%</b>

<sup>1</sup> Library prep recovery yields are shown for manual PureTarget library prep workflow using PureTarget kit 24 reagent bundle.

<sup>2</sup> Library prep recovery yields are shown for automated PureTarget library prep workflow using PureTarget kit 96 reagent bundle.

<sup>3</sup> **Note:** Automated PureTarget library prep workflow using PureTarget kit 96 employs **two** rounds of nuclease treatment.

# General best practices recommendations for preparing PureTarget libraries (cont.)

## Reagent and sample handling

- Room temperature is defined as any temperature in the range of 18–23°C for this protocol.
- Do not vortex enzymes.
- Bring SMRTbell cleanup beads to room temperature. Always vortex immediately prior to use.
- Bring Qubit reagents to room temperature prior to use.
- **Thaw frozen reagents at room temperature. Place on ice after thawing.**
- **Keep master mixes involving temperature-sensitive reagents on ice until use.**
- Quick-spin all reagents in a microcentrifuge to collect liquid at the bottom prior to use.
- Samples can be stored at the specified temperature at the safe stopping points listed in the protocol.



### PureTarget library prep reagents

[ Store kits at -20°C ]

PureTarget kit 2.0		
Step used	Tube color	Reagent
Dephosphorylation	Blue	Phosphatase
Cas9 digestion	Green	Cas9 Nuclease
Dephosphorylation/Cas9 digestion	Red	Cas9 buffer
dA tailing	Orange	dA tail buffer
dA tailing	Violet	Taq DNA Polymerase
dA tailing	Yellow	dATP
Nuclease treatment	White	Pure Target Nuclease mix
PureTarget panel		
	White	Repeat expansion panel 2.0 or
Cas9 digestion	Purple	Carrier panel or
	Orange	Control panel
SMRTbell prep kit 3.0		
gDNA repair	Purple	Repair buffer
gDNA repair	Blue	End repair mix
gDNA repair	Green	DNA repair mix
Adapter ligation	Yellow	Ligation mix
Adapter ligation	Red	Ligation enhancer
Nuclease treatment	Light green	Nuclease mix
Nuclease treatment	Light purple	Nuclease buffer
SMRTbell adapter index plate		
Adapter ligation	White plate	SMRTbell adapter index plate 96A, 96B, 96C, or 96D

### PureTarget cleanup beads kit<sup>1</sup>

[ Store kit at room temperature ]

Component	Tube color
PureTarget cleanup wash buffer*	Bottle, white
PureTarget cleanup binding buffer	Bottle, white
PureTarget cleanup beads	Clear

### PureTarget cleanup buffer kit

[ Store kit at room temperature ]

Component	Tube color
PureTarget cleanup buffer 1	Red
PureTarget cleanup buffer 2	Blue
PureTarget cleanup buffer 3	Green


### Polymerase kit<sup>2</sup>

[ Store kit at -20°C ]


Component	Tube color
Annealing buffer	Light blue
Standard sequencing primer	Light green
Polymerase buffer	Yellow
Loading buffer	Green
Dilution buffer	Blue
Sequencing polymerase	Purple
Sequencing control	Red

# Multiplexed sample processing workflow for manual library preparation method using PureTarget kit 24 [ For up to 48-plex ]

**Procedure & checklist – Generating PureTarget libraries** ([103-329-400](#))




## To process 8, 16 or 24 samples



**Note:** Each PureTarget kit 24 supports the preparation of **24** samples. Therefore, **two** PureTarget library prep kits are required to generate a Revio/Vega run of **32, 40, or 48** samples.

## To process 32, 40 or 48 samples




- With **manual library prep workflows** using PureTarget kit 24, up to **48** samples can be barcoded and sequenced per Revio or Vega SMRT Cell<sup>1</sup>
- Multiplexed samples should be processed in **batches of 8/16/24/32/40/48** samples
- For runs with 8/16/24 samples (**left figure**), follow protocol steps 1 – 9 to **process all samples in a single workflow**

- For runs with 32/40/48 samples (**right figure**), **process two workflows in parallel:**
  - **Workflow 1:** Process **24 samples** following protocol steps 1 – 7 (PureTarget library construction); and
  - **Workflow 2:** Process **an additional 8, 16 or 24 samples** following protocol steps 1 – 7
- At the end of step 7, **pool** 32/40/48 samples and proceed with PureTarget library cleanup (Step 8) and ABC (Step 9)

# Input DNA quality control

Before you begin, evaluate the size distribution of the input DNA using an Agilent Femto Pulse system or TapeStation system to determine whether it is suitable for the protocol




## 1. Input DNA quality control and dilution

✓ Step	Instructions
1.1	Bring the Qubit 1X dsDNA HS working solution and standards to room temperature.
1.2	Pulse vortex and/or pipette-mix each sample 5 times to homogenize the DNA in solution. Note: For viscous input DNA, it is important to homogenize the extracted DNA prior to start of the protocol. To homogenize the DNA, pulse-vortex and/or pipette-mix full sample volume with standard (not wide bore) tips. These steps will maintain HMW of your DNA but will improve accuracy of quantification and subsequent handling.
1.3	Quick-spin each sample to collect liquid.
1.4	Take a 1 µL aliquot from each sample and dilute with 9 µL of Elution buffer or water.
1.5	Measure DNA concentration with a Qubit fluorometer using the 1X dsDNA HS kit.
1.6	Measure DNA size distribution with a Femto Pulse system using the gDNA 165 kb analysis kit. Or Measure DNA size distribution with a TapeStation system using the Genomic DNA ScreenTape Analysis
1.7	Aliquot or bring 1- 4 µg DNA to a final volume of 43 µL per sample (24- 93 ng/µL) with nuclease-free water or Elution buffer.

SAFE STOPPING POINT - Store at 4°C

- Protocol requires high-quality, high molecular weight (HMW) human gDNA with **≥50% of the mass of DNA in molecules of length ≥30 kb**, or Genome quality number (GQN) at 30 kb of **≥5** based on Agilent Femto Pulse system (or DIN >7 based on Agilent TapeStation).
- For manual workflows, recommended input DNA amount is **2 µg per sample** to ensure sufficient gene copies to load and maximize sequencing coverage. This protocol supports 1-4 µg input DNA per sample.<sup>1</sup>


- For **viscous** input DNA, it is important to homogenize the extracted DNA prior to starting the protocol
- To homogenize the DNA, pulse-vortex and/or pipette-mix full sample volume with standard (not wide bore) tips.
- These steps will maintain HMW of your DNA but will improve accuracy of quantification and subsequent handling



Example Femto Pulse genomic DNA sizing QC analysis results for high-molecular weight genomic DNA extracted from a human whole blood sample using Nanobind PanDNA kit.

# Genomic DNA repair and SMRTbell bead cleanup

This step repairs input genomic DNA, which increases target coverage when using lower quality gDNA.



## 2.1 Genomic DNA repair

✓ Step Instructions

Add the following components to a microcentrifuge tube. Adjust component volumes for the number of samples being prepared, plus overage according to the table below.

gDNA repair master mix				8 reactions*	16 reactions*	24 reactions*
✓ Reagent	Component	Volume per sample				
Purple	Repair buffer	5 µL	48 µL	96 µl	144 µl	
Green	DNA repair mix	1.8 µL	17.3 µL	34.6 µl	51.8 µl	
Blue	End repair mix	0.2 µL	1.9 µL	3.8 µl	5.8 µl	
Total volume		7µL	67.3 µL	134.4 µl	201.6 µl	

\*20% overage included

2.1

2.2 Pipette-mix and quick-spin the gDNA repair master mix.

2.3 Add 7 µL of the gDNA repair master mix to each gDNA sample from step 1.7 for a total reaction volume of 50 µL.

2.4 Pipette-mix and quick-spin.

2.5 Run the **DNA repair** thermocycler program. Set the lid temperature to  $\geq 65^{\circ}\text{C}$ .

Step	Time	Temperature
1	60 min	37°C
2	30 min	52°C
3	Hold	4°C

- **IMPORTANT!** For Nanobind extracted blood (whole blood or RBC lysis method) gDNA, this step can be **skipped**. For **other** sample types or extraction methods, the gDNA repair step is **required** to ensure sufficient coverage.
- **IMPORTANT!** for Nanobind extracted saliva gDNA, a 3.1X 35% AMPure bead cleanup is **required** before the start of library prep.<sup>1</sup>


- Prepare a reaction master mix by adding the required components in the order and volume listed to a new microcentrifuge tube
- Adjust component volumes for the number of samples being prepared (8-plex, 16-plex or 24-plex), plus 15% overage according to table
- **IMPORTANT!** Prepare the master mix immediately before use. Keep the DNA repair mix on ice and immediately return the DNA repair mix back to freezer ( $-20^{\circ}\text{C}$ ) after use.

- Run **DNA repair** thermocycler program
- Set lid temperature to  $\geq 65^{\circ}\text{C}$  (if lid temp. is not programmable, leave at  $95\text{--}105^{\circ}\text{C}$ )

1.2X SMRTbell bead cleanup  
(Steps 2.6 – 2.21)

# Genomic DNA repair and SMRTbell bead cleanup (cont.)

This step repairs input genomic DNA, which increases target coverage when using lower quality gDNA. [ For Nanobind extracted blood (whole blood or RBC lysis method) gDNA, this step can be skipped. ]



## 2.6 1.2X SMRTbell bead cleanup

Step	Instructions
2.6	Add 60 $\mu$ L of resuspended, room-temperature SMRTbell cleanup beads to each sample.
2.7	Pipette-mix the beads until evenly distributed. Bead clumping may occur and is not a concern. Avoid over-pipetting as it may cause DNA/bead mixture to stick to the pipette tip.
2.9	Incubate at room temperature for 10 minutes to allow DNA to bind beads.
2.10	Place on a magnetic separation until beads separate fully from the solution.
2.11	Slowly remove the cleared supernatant without disturbing the beads. Discard the supernatant.
2.12	Slowly dispense 200 $\mu$ L, or enough to cover the beads, of freshly prepared 80% ethanol into each well. After 30 seconds, remove the 80% ethanol and discard.
2.13	Repeat the previous step.
2.14	Remove residual 80% ethanol and discard.
2.15	Remove sample from the magnetic rack. Immediately add 68 $\mu$ L of <b>Elution buffer</b> to each well and resuspend the beads.
2.17	Incubate at room temperature for 5 minutes to elute DNA.
2.18	Place on a magnetic separation rack until beads separate fully from the solution.
2.19	Slowly transfer 68 $\mu$ L cleared eluate without disturbing the beads to a new tube. Discard the old tube strip with beads.
2.20	Optional QC step: Take a 1 $\mu$ L each sample and measure DNA concentration with a Qubit fluorometer using the 1X dsDNA HS kit. <i>Expect a step recovery of 70-90% per sample.</i>
2.21	Proceed to the next step of the protocol. <b>SAFE STOPPING POINT – Store at 4°C overnight or 20°C long term</b>

- **Optional QC step:** Measure eluted DNA concentration using Qubit 1X dsDNA HS kit
- Expected step recovery is **70-90%** per sample

- Add **60  $\mu$ L** of SMRTbell cleanup beads to each sample and pipette-mix DNA/bead mixture until evenly distributed
- **Note: Bead clumping may occur and is not a concern**
- **Avoid over-pipetting** as it may cause DNA/bead mixture to stick to the pipette tip and clog the tip




- Elute sample in **68  $\mu$ L** of EB buffer
- Pipette-mix DNA/bead mixture until beads are fully resuspended and no clumps appear



# Dephosphorylation

This step dephosphorylates the genomic DNA 5' termini, which prevents adapter ligation to non-targeted genomic DNA.



## 3. Dephosphorylation

3.1 Add the following components to a new microcentrifuge tube. Adjust component volumes for the number of samples being prepared, plus overage according to the table below.

Dephosphorylation master mix						
Tube color	Component	Volume per sample	8 reactions*	16 reactions*	24 reactions*	
Red	Cas9 buffer	8 µL	76.8 µL	153.6 µl	230.4 µl	
Blue	Phosphatase	5 µL	48 µL	96 µl	144 µl	
Total volume		13 µL	124.8 µL	249.6 µl	374.4 µl	

\*20% overage included

3.2 Pipette-mix the Dephosphorylation master mix and quick-spin in a microcentrifuge to collect liquid.

3.3 Add 13 µL of the Dephosphorylation master mix to each sample from step 2.19 for total reaction volume of 80 µL.

3.4 Pipette-mix each sample.

3.5 Run the **Dephosphorylation** thermocycler program. Set the lid temperature to  $\geq 95^{\circ}\text{C}$ .

Step	Time	Temperature
1	20 min	37°C
2	5 min	85°C
3	Hold	4°C


3.6 Proceed to the next step of the protocol.

- Prepare a reaction master mix by adding the required components in the order and volume listed to a new microcentrifuge tube
- Adjust component volumes for the number of samples being prepared (8-plex, 16-plex or 24-plex), plus 15% overage according to table

- Run **dephosphorylation** thermocycler program
- Set the lid temperature to  $\geq 95^{\circ}\text{C}$  (if lid temp. is not programmable, leave at 95–105°C)

# Cas9 digestion and SMRTbell bead cleanup

This step enables digestion of double-stranded DNA at targeted regions mediated by Cas9-gRNA ribonucleoprotein complex formation



## 4.1 Cas9 digestion

Step	Instructions																																										
4.1	<p>✓ Add the following components to a new microcentrifuge tube. Adjust component volumes for the number of samples being prepared, plus overage according to the table below.</p> <table border="1"> <thead> <tr> <th colspan="2">Cas9 digestion master mix</th> <th>Volume per sample</th> <th>8 reactions*</th> <th>16 reactions*</th> <th>24 reactions*</th> </tr> </thead> <tbody> <tr> <td>✓</td> <td>Reagent</td> <td>Component</td> <td></td> <td></td> <td></td> </tr> <tr> <td></td> <td></td> <td>Nuclease free water</td> <td>2 µL</td> <td>19.2 µL</td> <td>38.4 µl</td> </tr> <tr> <td></td> <td>Red</td> <td>Cas9 buffer</td> <td>0.5 µL</td> <td>4.8 µL</td> <td>9.6 µl</td> </tr> <tr> <td></td> <td>Green</td> <td>Cas9 Nuclease</td> <td>0.5 µL</td> <td>4.8 µL</td> <td>9.6 µl</td> </tr> <tr> <td></td> <td>variable</td> <td>PureTarget panel**</td> <td>2 µL</td> <td>19.2 µL</td> <td>38.4 µl</td> </tr> <tr> <td></td> <td></td> <td>Total volume</td> <td>5.0 µL</td> <td>45.6 µL</td> <td>91.2 µl</td> </tr> </tbody> </table> <p>*20% overage included **Repeat expansion 2.0, Carrier, or Control panel with custom guide RNA</p>	Cas9 digestion master mix		Volume per sample	8 reactions*	16 reactions*	24 reactions*	✓	Reagent	Component						Nuclease free water	2 µL	19.2 µL	38.4 µl		Red	Cas9 buffer	0.5 µL	4.8 µL	9.6 µl		Green	Cas9 Nuclease	0.5 µL	4.8 µL	9.6 µl		variable	PureTarget panel**	2 µL	19.2 µL	38.4 µl			Total volume	5.0 µL	45.6 µL	91.2 µl
Cas9 digestion master mix		Volume per sample	8 reactions*	16 reactions*	24 reactions*																																						
✓	Reagent	Component																																									
		Nuclease free water	2 µL	19.2 µL	38.4 µl																																						
	Red	Cas9 buffer	0.5 µL	4.8 µL	9.6 µl																																						
	Green	Cas9 Nuclease	0.5 µL	4.8 µL	9.6 µl																																						
	variable	PureTarget panel**	2 µL	19.2 µL	38.4 µl																																						
		Total volume	5.0 µL	45.6 µL	91.2 µl																																						

4.2 Pipette-mix Cas9 digestion master mix and quick-spin in a microcentrifuge to collect liquid.

4.3 Add 5 µL of the Cas9 digestion master mix to each sample for total reaction volume of 85 µL.


4.4 Pipette-mix each sample

Run the **Cas-9 digestion** thermocycler program. Set the lid temperature to  $\geq 65^{\circ}\text{C}$ .

Step	Time	Temperature
4.5	1	60 min
	2	Hold
		4°C

## 4.7 1.2X SMRTbell bead cleanup

Step	Instructions
4.7	Add 102 µL of resuspended, room-temperature SMRTbell cleanup beads to each sample.
4.8	Pipette-mix the beads slowly until evenly distributed. Bead clumping may occur and is not a concern. Avoid over-pipetting as it may cause DNA/bead mixture to stick to the pipette tip.




- Prepare a reaction master mix by adding the required components in the order and volume listed to a new microcentrifuge tube
- Adjust component volumes for the number of samples being prepared (8-plex, 16-plex or 24-plex), plus 15% overage according to table

- Run **Cas-9 digestion** thermocycler program
- Set the lid temperature to  $\geq 65^{\circ}\text{C}$  (if lid temp. is not programmable, leave at  $95\text{--}105^{\circ}\text{C}$ )

- Perform **1.2X** SMRTbell bead cleanup
- **Optional QC step:** Measure eluted DNA concentration using Qubit 1X dsDNA HS kit
- Expected step recovery is **~50-80%** per sample

# dA tail

This step enables A-tailing of DNA 3' ends after Cas9-gRNA digestion at targeted regions



## 5. dA tail

✓ Step Instructions

Add the following components to a new microcentrifuge tube. Adjust component volumes for the number of samples being prepared, plus overage according to the table below.

dA-tailing master mix						
✓	Reagent	Component	Volume per sample	8 reactions*	16 reactions*	24 reactions*
5.1		Nuclease free water	3.5	34 µL	68 µL	102 µL
	Orange	dA tail buffer	5 µL	48 µL	96 µL	144 µL
	Yellow	dATP (100 mM)	0.2 µL	2 µL	4 µL	6 µL
	Violet	Taq DNA polymerase	1.3 µL	12.5 µL	25 µL	37.5 µL
		Total volume	10 µL	96 µL	192 µL	288 µL

\*20% overage included

5.2 Pipette-mix the dA-tailing master mix and quick-spin in a microcentrifuge to collect liquid.

5.3 Add 10 µL of the dA-tailing master mix to each sample for a total reaction volume of 50 µL.

5.4 Pipette-mix each sample.

Run the dA-tail thermocycler program. Set the lid temperature to ≥95°C.

Step	Time	Temperature
1	1 min	37°C
2	5 min	72°C
3	Hold	4°C


5.6 Proceed to the next step of the protocol.

- Prepare a reaction master mix by adding the required components in the order and volume listed to a new microcentrifuge tube
- Adjust component volumes for the number of samples being prepared (8-plex, 16-plex or 24-plex), plus 15% overage according to table

- Run **dA-tail** thermocycler program
- Set the lid temperature to ≥95°C (if lid temp. is not programmable, leave at 95–105°C)

# Adapter ligation and SMRTbell bead cleanup

This step ligates the indexed SMRTbell adapters to the ends of each targeted DNA fragment



## 6.1 Adapter ligation

Step	Instructions																								
6.1	Add 4 $\mu\text{L}$ SMRTbell adapter index to each sample from Section 5 and pipette mix.																								
6.2	Add the following to a new microcentrifuge tube. Adjust component volumes for the number of samples being prepared, plus overage according to the table below. Note: Aspirate and dispense viscous ligation mix slowly to reduce liquid sticking to the exterior and interior wall of the pipette tip.																								
Ligation master mix																									
	<table border="1"> <thead> <tr> <th>Reagent</th> <th>Component</th> <th>Volume per sample</th> <th>8 reactions*</th> <th>16 reactions*</th> <th>24 reactions*</th> </tr> </thead> <tbody> <tr> <td>Yellow</td> <td>Ligation mix</td> <td>30 <math>\mu\text{L}</math></td> <td>264 <math>\mu\text{L}</math></td> <td>528 <math>\mu\text{L}</math></td> <td>792 <math>\mu\text{L}</math></td> </tr> <tr> <td>Red</td> <td>Ligation enhancer</td> <td>1 <math>\mu\text{L}</math></td> <td>8.8 <math>\mu\text{L}</math></td> <td>17.6 <math>\mu\text{L}</math></td> <td>26.4 <math>\mu\text{L}</math></td> </tr> <tr> <td colspan="2">Total volume</td> <td>31 <math>\mu\text{L}</math></td> <td>272.8 <math>\mu\text{L}</math></td> <td>545.6 <math>\mu\text{L}</math></td> <td>818.4 <math>\mu\text{L}</math></td> </tr> </tbody> </table>	Reagent	Component	Volume per sample	8 reactions*	16 reactions*	24 reactions*	Yellow	Ligation mix	30 $\mu\text{L}$	264 $\mu\text{L}$	528 $\mu\text{L}$	792 $\mu\text{L}$	Red	Ligation enhancer	1 $\mu\text{L}$	8.8 $\mu\text{L}$	17.6 $\mu\text{L}$	26.4 $\mu\text{L}$	Total volume		31 $\mu\text{L}$	272.8 $\mu\text{L}$	545.6 $\mu\text{L}$	818.4 $\mu\text{L}$
Reagent	Component	Volume per sample	8 reactions*	16 reactions*	24 reactions*																				
Yellow	Ligation mix	30 $\mu\text{L}$	264 $\mu\text{L}$	528 $\mu\text{L}$	792 $\mu\text{L}$																				
Red	Ligation enhancer	1 $\mu\text{L}$	8.8 $\mu\text{L}$	17.6 $\mu\text{L}$	26.4 $\mu\text{L}$																				
Total volume		31 $\mu\text{L}$	272.8 $\mu\text{L}$	545.6 $\mu\text{L}$	818.4 $\mu\text{L}$																				
*10% overage included																									

- Add 4  $\mu\text{L}$  SMRTbell adapter index plate 96A to each sample and then tap-mix or pipette up and down 10 times (do not vortex)

- Prepare a reaction master mix by adding the required components in the order and volume listed to a new microcentrifuge tube
- Adjust component volumes for the number of samples being prepared (8-plex, 16-plex or 24-plex), plus 15% overage according to table

6.5	Pipette-mix each sample Run the <b>Adapter ligation</b> thermocycler program. Heating of the lid is not necessary.									
6.6	<table border="1"> <thead> <tr> <th>Step</th> <th>Time</th> <th>Temperature</th> </tr> </thead> <tbody> <tr> <td>1</td> <td>30 min</td> <td>20°C or RT</td> </tr> <tr> <td>2</td> <td>Hold</td> <td>4°C</td> </tr> </tbody> </table>	Step	Time	Temperature	1	30 min	20°C or RT	2	Hold	4°C
Step	Time	Temperature								
1	30 min	20°C or RT								
2	Hold	4°C								

- Run **adapter ligation** thermocycler program
- Heating of the lid is not necessary


## 6.7 1X SMRTbell bead cleanup

Step	Instructions
6.7	Add 85 $\mu\text{L}$ (1X) of resuspended, room-temperature SMRTbell cleanup beads to each sample.
6.8	Pipette-mix until evenly distributed.
6.9	Incubate at room temperature for 10 minutes to allow DNA to bind beads.

- Perform 1X SMRTbell bead cleanup<sup>1</sup>
- **Optional QC step:** Measure eluted DNA concentration using Qubit 1X dsDNA HS kit
- Expected step recovery is ~70-90% per sample

# Nuclease treatment and SMRTbell bead cleanup

Nuclease treatment step removes unligated DNA fragments<sup>1</sup> that have not formed complete (intact) SMRTbell templates



## 7.1 Nuclease treatment

✓ Step Instructions

Add the following components to a new microcentrifuge tube. Adjust component volumes for the number of samples being prepared, plus overage according to the table below.

Nuclease master mix					
✓ Reagent	Component	Vol per sample	8 reactions*	16 reactions*	24 reactions*
Light purple	Nuclease buffer	5 µL	44 µL	88 µL	132 µL
Light Green	Nuclease mix	5 µL	44 µL	88 µL	132 µL
White	PureTarget Nuclease	0.5 µL	4.4 µL	8.8 µL	13.2 µL
Total volume		10.5 µL	92.4 µL	184.8 µL	277.2 µL

\*10% overage included

7.1 Pipette-mix the Nuclease master mix and quick-spin in a microcentrifuge to collect liquid.

7.2 Add 10.5 µL of Nuclease master mix to 40 µL of each sample for a total volume of 50.5 µL.

7.3 Pipette-mix each sample.

7.4 Run the **Nuclease treatment** thermocycler program. Set the lid temperature to  $\geq 65^{\circ}\text{C}$ .

Step	Time	Temperature
1	60 min	37° C
2	Hold	4° C

- Prepare a reaction master mix by adding the required components in the order and volume listed to a new microcentrifuge tube
- Adjust component volumes for the number of samples being prepared (8-plex, 16-plex or 24-plex), plus 15% overage according to table

- Run **nuclease treatment** thermocycler program
- Set the lid temperature to  $\geq 65^{\circ}\text{C}$  (if lid temp. is not programmable, leave at 95–105°C)

Sample pooling and 1.2X SMRTbell bead cleanup (Steps 7.6 – 7.19)

Concentrate samples with 1.2X SMRTbell bead cleanup (Steps 7.20 – 7.33)

### Note:


- **Manual** PureTarget library prep workflow using PureTarget kit 24 employs **one round of nuclease treatment**
- **Automated** PureTarget library prep workflow using PureTarget kit 96 employs **two rounds of nuclease treatment**.<sup>2</sup>

<sup>1</sup> Unligated DNA fragments comprise the majority of DNA fragments present prior to this step.

<sup>2</sup> Refer to **Procedure & checklist – Generating PureTarget libraries with PureTarget kit 96 – Automation protocol** (103-740-800) for details regarding nuclease treatment and post-nuclease cleanup. 41

# Nuclease treatment and SMRTbell bead cleanup (cont.)

Pool samples for post-nuclease cleanup and concentration steps using SMRTbell cleanup beads



## 7.6 Sample pooling and 1.2X SMRTbell bead cleanup

Step	Instructions										
7.6	In a 1.5 mL DNA LoBind tube combine nuclease treated libraries from step 7.5 in groups of 8 (8 x 50.5 µL) for a final sample volume of 404 µL.										
7.7	Add 485 µL of SMRTbell cleanup beads to each pooled nuclease-treated sample. Pipette-mix until the beads are evenly distributed.										
7.8	Quick-spin the tube in a microcentrifuge to collect all liquid from the sides of the tubes.										
7.15	Remove from the magnetic rack. Immediately add 55 µL of Elution buffer to each tube and resuspend the beads by pipette mixing until resuspended.										
7.16	Incubate at room temperature for 5 minutes to elute DNA.										
7.17	Place on a magnetic separation rack until beads separate fully from the solution.										
7.18	Slowly transfer 55 µL cleared eluate from each nuclease reaction and pool into a new 1.5 mL DNA LoBind tube. <table border="1" data-bbox="861 753 1243 888"> <thead> <tr> <th>Plex</th> <th>Total pooled volume</th> </tr> </thead> <tbody> <tr> <td>8-plex</td> <td>55 µl</td> </tr> <tr> <td>16-plex</td> <td>110 µl*</td> </tr> <tr> <td>24-plex</td> <td>165 µl*</td> </tr> <tr> <td>48-plex</td> <td>330 µl*</td> </tr> </tbody> </table>	Plex	Total pooled volume	8-plex	55 µl	16-plex	110 µl*	24-plex	165 µl*	48-plex	330 µl*
Plex	Total pooled volume										
8-plex	55 µl										
16-plex	110 µl*										
24-plex	165 µl*										
48-plex	330 µl*										
7.19	Proceed to the step 7.20 if the pooled volume is greater than 55 µL; otherwise proceed to step 7.32 for QC.										

## 7.20 Concentrate samples with 1.2X SMRTbell bead cleanup

Step	Instructions												
7.20	Add 1.2X volume of SMRTbell cleanup beads to the sample from step 7.18. <table border="1" data-bbox="733 1139 1431 1173"> <thead> <tr> <th>Plex</th> <th>Total pooled volume</th> <th>SMRTbell cleanup beads volume</th> </tr> </thead> <tbody> <tr> <td>16-plex</td> <td>110 µl</td> <td>132 µl</td> </tr> <tr> <td>24-plex</td> <td>165 µl</td> <td>198 µl</td> </tr> <tr> <td>48-plex</td> <td>330 µl</td> <td>396 µl</td> </tr> </tbody> </table>	Plex	Total pooled volume	SMRTbell cleanup beads volume	16-plex	110 µl	132 µl	24-plex	165 µl	198 µl	48-plex	330 µl	396 µl
Plex	Total pooled volume	SMRTbell cleanup beads volume											
16-plex	110 µl	132 µl											
24-plex	165 µl	198 µl											
48-plex	330 µl	396 µl											
7.32	QC step: Take a 1 µL aliquot from each sample and measure DNA concentration with a Qubit fluorometer using the 1X dsDNA HS kit. Expect an overall recovery of 0.05-0.5% relative to the starting gDNA total mass.												
7.33	Proceed to the next step of the protocol.												

- Pool nuclease-treated libraries from Step 7.5 in **groups of 8** (i.e., 8 x 50.5 µL) for a final sample vol. of **404 µL<sup>1</sup>**
- Perform **1.2X** SMRTbell bead cleanup by adding **485 µL** of SMRTbell cleanup beads to each 8-plex pool

- Elute each 8-plex pooled sample in **55 µL** of EB buffer

- Slowly transfer **55 µL** of cleared eluate from **all** tubes into a **single** new 1.5 mL DNA LoBind tube
- Proceed to **Step 7.20** if total pooled sample volume **>55 µL**; otherwise proceed directly to **Step 7.32** for QC

- If Step 7.18 total pooled volume >55 µL:** Concentrate samples by adding **1.2X** SMRTbell cleanup beads to pooled samples from Step 7.18


Plex	Total pooled volume	SMRTbell cleanup beads volume
16-plex	110 µl	132 µl
24-plex	165 µl	198 µl
48-plex	330 µl	396 µl

- Elute concentrated samples in **55 µL** of EB buffer


- For manual workflows:** Expect an overall recovery of **0.05-0.5%** relative to the starting input gDNA total mass<sup>1</sup>
- Higher-than-expected yields may indicate inefficient nuclease digestion, or sample incompatibility. Refer to Appendix A2 for recommendations to further digest sample

# PureTarget library cleanup

Perform final library cleanup using PureTarget cleanup kit to prepare the SMRTbell library for sequencing



## PureTarget library cleanup procedure overview




- Perform one cleanup reaction per final pool of 8 – 48 PureTarget samples for manual library prep workflows using PureTarget kit 24<sup>1</sup>
- All reagents in the PureTarget cleanup kit should be stored at **room-temperature**; do not chill or freeze any reagent
- **Vortex** buffers in the PureTarget cleanup kit prior to use
- **Prior to the first use of the PureTarget cleanup beads kit**, add 15 mL of 200 Proof ethanol to PureTarget cleanup wash buffer and mix well

# PureTarget library cleanup (cont.)

Perform final library cleanup using PureTarget cleanup kit to prepare the SMRTbell library for sequencing

## 8. PureTarget Cleanup of the SMRTbell library



Step	Instructions																									
8.1	Add 28 $\mu\text{L}$ PureTarget cleanup buffer 1 (●) to 55 $\mu\text{L}$ SMRTbell library from step 7.31 and pipette-mix.  Add the following components to a new microcentrifuge tube. <table border="1"><thead><tr><th colspan="5">PureTarget cleanup buffer 2+3 mix</th></tr><tr><th>Reagent</th><th>Component</th><th>Volume per sample</th><th colspan="2">Volume (overage included)</th></tr></thead><tbody><tr><td>Blue</td><td>PureTarget cleanup buffer 2</td><td>14 <math>\mu\text{L}</math></td><td>18 <math>\mu\text{L}</math></td><td></td></tr><tr><td>Green</td><td>PureTarget cleanup buffer 3</td><td>14 <math>\mu\text{L}</math></td><td>18 <math>\mu\text{L}</math></td><td></td></tr><tr><td colspan="2">Total volume</td><td>28 <math>\mu\text{L}</math></td><td>36 <math>\mu\text{L}</math></td><td></td></tr></tbody></table>	PureTarget cleanup buffer 2+3 mix					Reagent	Component	Volume per sample	Volume (overage included)		Blue	PureTarget cleanup buffer 2	14 $\mu\text{L}$	18 $\mu\text{L}$		Green	PureTarget cleanup buffer 3	14 $\mu\text{L}$	18 $\mu\text{L}$		Total volume		28 $\mu\text{L}$	36 $\mu\text{L}$	
PureTarget cleanup buffer 2+3 mix																										
Reagent	Component	Volume per sample	Volume (overage included)																							
Blue	PureTarget cleanup buffer 2	14 $\mu\text{L}$	18 $\mu\text{L}$																							
Green	PureTarget cleanup buffer 3	14 $\mu\text{L}$	18 $\mu\text{L}$																							
Total volume		28 $\mu\text{L}$	36 $\mu\text{L}$																							
8.2	Pulse vortex to mix and quick-spin. <b>Important: PureTarget cleanup buffer 2 and buffer 3 must be pre-mixed prior to adding to sample. If added individually to sample, sample recovery may be low.</b>																									
8.3	Add 28 $\mu\text{L}$ PureTarget cleanup buffer 2+3 mix to the mixture from step 8.1 and pipette-mix. The solution will turn opaque and cloudy after homogenous mixing. Incubate at room temperature for 1 minute.																									
8.4	Centrifuge at 13,000 $\times g$ for 2 minutes at room temperature.																									
8.5	Taking care to avoid the pellet, transfer 95 $\mu\text{L}$ of the clear supernatant to a clean 1.5 mL LoBind tube. Discard the tube with pellet.																									
8.6	Thoroughly resuspend PureTarget cleanup beads by pulse-vortexing and quick spin. Add 5 $\mu\text{L}$ of fully resuspended PureTarget cleanup beads to 95 $\mu\text{L}$ of the sample from step 8.5.																									
8.7	Add 100 $\mu\text{L}$ of PureTarget cleanup binding buffer to the sample from step 8.6.																									
8.8	Pipette-mix until evenly distributed. Incubate at room temperature for 15 mins, with <b>periodic pipette-mixing to ensure that beads are resuspended during the entire incubation.</b>																									
8.9	Quick-spin the tube in a microcentrifuge to collect all liquid from the sides of the tubes.																									
8.10	Place the tube in a magnetic separation rack until beads separate fully from the solution.																									
8.11	Slowly remove the cleared supernatant without disturbing the beads. Discard the supernatant.																									
8.12	Slowly dispense 500 $\mu\text{L}$ of PureTarget cleanup wash buffer (with ethanol added) into the tube. Remove the tube from the magnetic separation rack and gently resuspend the beads. Close the tube cap and invert the tube 5 times to wash off any residual PureTarget Cleanup Binding Buffer.																									
8.13	Place the tube in a magnetic separation rack until beads separate fully from the solution.																									
8.14	Slowly remove the cleared supernatant without disturbing the beads. Discard the supernatant.																									

- Add **PureTarget cleanup buffers 1, 2 and 3** to SMRTbell library sample (**55  $\mu\text{L}$** ) in the order shown
  - Prepare a master mix of PureTarget cleanup Buffer 2 and Buffer 3 by combining equal volumes of each (16  $\mu\text{L}$  per buffer)
  - Add PureTarget cleanup buffer 1 to SMRTbell library **before** adding the master mix containing PureTarget cleanup buffers 2 & 3
- **Note:** After adding the PureTarget cleanup buffer 2/3 mixture to the sample, the solution will turn **opaque and cloudy** after homogenous mixing; a white jelly-like pellet will form after centrifugation




- Prior to the first use of the PureTarget cleanup beads kit, add 15 mL of 200 Proof ethanol to PureTarget cleanup wash buffer and mix well

# PureTarget library cleanup (cont.)

Perform final library cleanup using PureTarget cleanup kit to prepare the SMRTbell library for sequencing

## 8. PureTarget Cleanup of the SMRTbell library (cont.)



Step	Instructions
8.15	Slowly dispense 500 $\mu$ L of fresh 80% ethanol into the tube without disturbing the beads.
8.16	After 30 seconds, remove the 80% ethanol and discard.
8.17	<b>It is critical to remove residual 80% ethanol for efficient elution:</b> <ul style="list-style-type: none"><li>Remove the tube from the magnetic separation rack.</li><li>Quick-spin the tube in a microcentrifuge.</li><li>Place the tube back in a magnetic separation rack until beads separate fully from the solution.</li><li>Remove residual 80% ethanol and discard.</li><li>Air dry the bead pellet for 2 mins.</li></ul>
8.18	Remove the tube from the magnetic rack. Add 26 $\mu$ L of Elution buffer and resuspend the beads by gentle vortexing.
8.19	Incubate at room temperature for 15 minutes to elute DNA, with gentle vortexing every 5 mins to fully resuspend the beads.
8.20	Place the tube in a magnetic separation rack until beads separate fully from the solution.
8.21	Slowly transfer 26 $\mu$ L cleared supernatant without disturbing the beads to a new 1.5 mL DNA LoBind tube.
8.22	QC step: Take a 1 $\mu$ L aliquot from the sample and measure DNA concentration with a Qubit fluorometer using the 1X dsDNA HS kit.  Expect a step recovery of 40-60%, and an overall recovery of 0.02-0.2% relative to starting gDNA total mass. For example, starting with 1.0 $\mu$ g gDNA per sample for a 48-plex prep to load on Revio with SPRQ, the final mass of the pooled 48-plex library recovered at this step is expected to be in the range of 10 ng– 100 ng.

SAFE STOPPING POINT – Store at 4°C overnight or -20°C long term

Proceed to **Step 9** to perform sample setup (ABC) and final dilution (Do not use SMRT Link Sample Setup)

Elute cleaned library (containing up to 48-plex) into 26  $\mu$ L of EB

- Take 1  $\mu$ L to perform DNA concentration QC using Qubit dsDNA HS kit
- Use remaining 25  $\mu$ L to proceed with sample setup (ABC<sup>1</sup>)

- **QC step:** Measure DNA concentration using Qubit 1X dsDNA HS kit
- **For manual workflows:** Expected recovery is ~0.02 – 0.2% relative to starting input gDNA mass (e.g., ~10 ng – 100 ng of final library per 48-plex when starting with 1.0  $\mu$ g per sample)<sup>1</sup>
- **Note:** Significantly higher recoveries could indicate a sample type incompatibility or an issue with the nuclease treatment step

For primer annealing, polymerase binding & complex cleanup (ABC) and final dilution steps, follow sample setup instructions provided in [Procedure & checklist – Generating PureTarget libraries with PureTarget kit 24 – Manual protocol](#) (103-740-700) [ or [Procedure & checklist – Generating PureTarget libraries with PureTarget kit 96 – Automation protocol](#) (103-740-800) if using automation ] – Do not use SMRT Link Sample Setup



# PureTarget sequencing preparation workflow details

# Sample setup and run design recommendations for PureTarget libraries

## Key sample setup and run design setup parameters for Revo and Vega systems

Workflow	Key setup parameters	Vega system recommended settings	Revo system recommended settings
		for PureTarget libraries	
Sample setup	Library type	Standard	
	Primer	Standard sequencing primer	
	Polymerase kit	Vega polymerase kit	Revo SPRQ polymerase kit
Run design	Application type	PureTarget repeat expansion / PureTarget carrier / PureTarget custom	
	Library type	Standard	
	Insert size (bp)	5000 <sup>1</sup>	
	Library concentration (pM)	<user-specified> <sup>2</sup>	
	Movie acquisition time	24 hours	
	Use adaptive loading	N/A	NO <sup>3</sup>
	Data options	Include base kinetics = NO Consensus Mode = MOLECULE	
	Analysis options <sup>4</sup>	<b>PureTarget repeat expansion panel 2.0:</b> Add Analysis = YES Analysis Workflow = PureTarget repeat expansion	
		<b>PureTarget carrier panel:</b> Add Analysis = YES Analysis Workflow = Target Enrichment <sup>5</sup>	
		<b>PureTarget custom panel:</b> Add Analysis = YES Analysis Workflow = Target Enrichment <sup>5</sup>	

For PureTarget samples, specify **Insert Size = 5000 bp**<sup>1</sup>

<sup>1</sup> **Note:** This '5000 bp' value only serves as a placeholder since the actual average insert size distribution of PureTarget repeat expansion samples may vary.

If you calculate your PureTarget library concentration, enter it – but if not, enter "0"<sup>2</sup>

<sup>2</sup> **Note:** For PureTarget samples, run conditions are not impacted by insert size or library concentration.

**IMPORTANT!** For Revo system, specify **Use Adaptive Loading = NO**<sup>3</sup>


<sup>3</sup> **Note:** In SMRT Link v25.3, Adaptive Loading is ON by default for all Revo system run designs. For PureTarget libraries, specify **Adaptive Loading = NO** to enable correct sample immobilization conditions to be used on the Revo system. PureTarget samples should not be included in the same run design as other sample types that require Adaptive Loading to be enabled (If an attempt is made to save a Revo run design with Adaptive Loading = YES and the run includes one or more PureTarget samples, then an error window will appear prompting the user to specify NO for Use Adaptive Loading.)

<sup>4</sup> Users who prefer to use command line tools may configure SMRT Link to perform automatic demultiplexing only. Demultiplexed BAM files may then be transferred for downstream command line analysis.

<sup>5</sup> SMRT Link Target Enrichment analysis application features a generalizable workflow that can be used to evaluate PureTarget panel performance QC (see SMRT Link user guide for further details.)

# Sample setup workflow overview for PureTarget libraries

Follow sample setup instructions for PureTarget libraries in **Procedure & checklist** library prep protocol – Do not use SMRT Link Sample Setup



- Use **entire volume (25 µL)** of final PureTarget library per ABC reaction
- After complex cleanup, elute samples in Loading buffer and add specified volume of diluted sequencing control  
→ Do not use SMRT Link Loading Calculator

## ABC procedure for PureTarget kit 24 manual workflows (8-, 16-, 24- or 48-plex)

**P & C – Generating PureTarget libraries with PureTarget kit 24 – Manual protocol (103-740-700)**



Perform **one** ABC reaction per final pool of 8 – 48 samples in 25 µL)

Step	Instructions																				
9.1	Prepare the appropriate volume of master mix with 10% overage using the per reaction volumes listed below.																				
<table border="1"> <thead> <tr> <th colspan="4">Annealing mix</th> </tr> <tr> <th>Tube</th> <th>Component</th> <th>Volume</th> <th></th> </tr> </thead> <tbody> <tr> <td>Light blue</td> <td>Annealing buffer</td> <td>12.5 µL</td> <td></td> </tr> <tr> <td>Light green</td> <td>Standard sequencing primer</td> <td>12.5 µL</td> <td></td> </tr> <tr> <td colspan="2">Total volume</td> <td>25 µL</td> <td></td> </tr> </tbody> </table>		Annealing mix				Tube	Component	Volume		Light blue	Annealing buffer	12.5 µL		Light green	Standard sequencing primer	12.5 µL		Total volume		25 µL	
Annealing mix																					
Tube	Component	Volume																			
Light blue	Annealing buffer	12.5 µL																			
Light green	Standard sequencing primer	12.5 µL																			
Total volume		25 µL																			

9.2 Pipette-mix the Annealing mix and quick spin to collect liquid.

Prepare sample for sequencing. Add the following components for the appropriate sequencing platform.

Step	Sequencing sample and control		
	Component	Revio SPRQ	Vega
10.4	Sample from step 9.2.2	23 µL	98 µL
	Diluted internal sequencing control (Dilution 3)	1 µL	3 µL
	Total volume	24 µL	101 µL

10.5 Load the Sequencing sample and control prepared in step 10.4. Load the exact volume indicated for the appropriate sequencing platform.

Revio SPRQ: 23 µL  
Vega: 100 µL

Sequence pooled samples on one Revio or Vega SMRT Cell



## ABC procedure for PureTarget kit 96 automated workflows (96-plex)

**P & C – Generating PureTarget libraries with PureTarget kit 96 – Automation protocol (103-740-800)**



Perform **one** ABC reaction per final pool of 96 samples in 25 µL)

Step	Instructions																									
9.1	Prepare the appropriate volume of master mix with 10% overage using the per reaction volumes listed below.																									
<table border="1"> <thead> <tr> <th colspan="5">Annealing mix</th> </tr> <tr> <th>Tube</th> <th>Component</th> <th>Volume</th> <th>Volume (overage included)</th> <th></th> </tr> </thead> <tbody> <tr> <td>Light blue</td> <td>Annealing buffer</td> <td>12.5 µL</td> <td>15 µL</td> <td></td> </tr> <tr> <td>Light green</td> <td>Standard sequencing primer</td> <td>12.5 µL</td> <td>15 µL</td> <td></td> </tr> <tr> <td colspan="2">Total volume</td> <td>25 µL</td> <td>30 µL</td> <td></td> </tr> </tbody> </table>		Annealing mix					Tube	Component	Volume	Volume (overage included)		Light blue	Annealing buffer	12.5 µL	15 µL		Light green	Standard sequencing primer	12.5 µL	15 µL		Total volume		25 µL	30 µL	
Annealing mix																										
Tube	Component	Volume	Volume (overage included)																							
Light blue	Annealing buffer	12.5 µL	15 µL																							
Light green	Standard sequencing primer	12.5 µL	15 µL																							
Total volume		25 µL	30 µL																							

9.2 Pipette-mix the Annealing mix and quick spin to collect liquid.

Prepare sample for sequencing. Add the following components to the Revio SPRQ sequencing reagent plate.

Step	Sequencing reagent plate sample	
	Component	Volume
10.4	Sample from step 10.2.2	23 µL
	Diluted internal sequencing control (Dilution 3)	1 µL
	Total volume	24 µL

10.5 Load **exactly 23 µL** of sample prepared in 10.4 per sequencing well or store at 4°C for up to 24 hours before use. For long-term storage, store the sample at -20°C.

Sequence pooled samples on one Revio SMRT Cell



# SMRT Link Run Design options for PureTarget libraries – Sample indexing

Specifying sample indexing (barcoding) information for Revio and Vega systems

**Samples**

Sample is indexed  YES  NO


Indexes Required: SMRTbell adapter indexes

Biosample names Required: Interactively From a File

Default = YES for Sample is indexed

Specify Indexes FASTA = SMRTbell adapter indexes

Example PureTarget library molecule containing SMRTbell indexed adapters<sup>1</sup> at both ends



Both forward & reverse terminal SMRTbell adapters contain the same barcode sequence

Example interactive biosample name specification for a multiplexed PureTarget library sample

Barcode Selector and Sample Name Editor

Barcode ID	Bio Sample ID
------------	---------------

Barcode Selector and Sample Name Editor

Barcode ID	Bio Sample ID
bc2001--bc2001	PureTarget_Sample_1
bc2002--bc2002	PureTarget_Sample_2

SMRT Link



SMRTbell adapter indexes

```
>bc2001
ATCGTGCGACGAGTAT
>bc2002
TGCATGTCATGAGTAT
>bc2003
ACGAGTGCTCGAGTAT
>bc2004
TGCAGTGCTCGAGTAT
```

# SMRT Link Run Design options for PureTarget libraries – Data options

Specifying data options for Revio and Vega systems

**Data Options**

Include Base Kinetics   YES  NO

Consensus Mode  MOLECULE  STRAND

Assign Data To Project  General Project

Specify Include Base Kinetics = NO

Specify Consensus Mode = MOLECULE<sup>1</sup>

Can leave **Include Base Kinetics** and **Consensus Mode** fields at their **default settings** for PureTarget library samples

# SMRT Link Run Design options for PureTarget libraries – Analysis options

## 1. Specifying PureTarget analysis options for Revio and Vega systems – Add analysis

### PureTarget repeat expansion panel

Application Required PureTarget repeat expansion

Analysis Options

Add Analysis  YES  NO

Analysis Name Required

Select Analysis Workflow Required PureTarget repeat expansion

Reference Genome Required Human Genome hg38, with Gencode v39 annotations

Target and repeat definitions Required PureTarget repeat expansion panel 2.0

Advanced Parameters

Add Analysis  YES  NO

Default = YES for Add Analysis<sup>1</sup>

### PureTarget carrier panel

Application Required PureTarget carrier

Analysis Options

Add Analysis  YES  NO

Analysis Name Required

Select Analysis Workflow Required Target Enrichment

Reference Set Required Human Genome hg38, with Gencode v39 annotations

Target Set Required PureTarget carrier panel 1.0

Mark PCR d

Advanced Parameters

Add Analysis  YES  NO

Default = YES for Add Analysis<sup>1</sup>

### PureTarget custom panel

Application Required PureTarget custom

Analysis Options

Add Analysis  YES  NO

Analysis Name Required

Select Analysis Workflow Required Target Enrichment

Reference Set Required Human Genome hg38, with Gencode v39 annotations

Target Set Required

Mark PCR d

Advanced Parameters

Add Analysis  YES  NO

Default = YES for Add Analysis<sup>1</sup>

<sup>1</sup> Users have two options for analysis when setting up sequencing runs in SMRT Link. For the fastest turnaround time, users can specify to add the desired PureTarget analysis in their run design and analysis will be automatically performed when sequencing is complete. Alternatively, users who prefer command line analysis may configure SMRT Link to perform automatic demultiplexing only. Demultiplexed BAM files may then be transferred for downstream command line analysis using appropriate PacBio or PacBio Compatible software tools.

# SMRT Link Run Design options for PureTarget libraries – Analysis options

## 2. Specifying PureTarget analysis options for Revio and Vega systems – Analysis workflow selection<sup>1</sup>

### PureTarget repeat expansion panel

Application Required PureTarget repeat expansion

Analysis Options

Add Analysis  YES  NO

Analysis Name Required

Select Analysis Workflow Required PureTarget repeat expansion

Reference Genome Required Human Genome hg38, with Gencode v39 annotations

Target and repeat definitions Required PureTarget repeat expansion panel 2.0

Advanced Parameters

Select Analysis Workflow Required PureTarget repeat expansion

- SMRT Link **PureTarget repeat expansion analysis application** produces **target enrichment summary QC statistics** and uses tandem repeat genotyping tool (TRGT) for variant calling and for visualization

### PureTarget carrier panel

Application Required PureTarget carrier

Analysis Options

Add Analysis  YES  NO

Analysis Name Required

Select Analysis Workflow Required Target Enrichment

Reference Set Required Human Genome hg38, with Gencode v39 annotations

Target Set Required PureTarget carrier panel 1.0

Mark PCR duplicates  ON  OFF

Advanced Parameters

Select Analysis Workflow Required Target Enrichment

- SMRT Link **Target Enrichment analysis application** produces **target coverage summary QC statistics** for PureTarget libraries

### PureTarget custom panel

Application Required PureTarget custom

Analysis Options

Add Analysis  YES  NO

Analysis Name Required

Select Analysis Workflow Required Target Enrichment

Reference Set Required Human Genome hg38, with Gencode v39 annotations

Target Set Required

Mark PCR duplicates  ON  OFF

Advanced Parameters

Select Analysis Workflow Required Target Enrichment

- SMRT Link **Target Enrichment analysis application** produces **target coverage summary QC statistics** for PureTarget libraries

# SMRT Link Run Design options for PureTarget libraries – Analysis options

## 3. Specifying PureTarget analysis options for Revio and Vega systems – Reference set selection<sup>1</sup>

### PureTarget repeat expansion panel

Application Required PureTarget repeat expansion

Analysis Options

Add Analysis  YES  NO

Analysis Name Required

Select Analysis Workflow Required PureTarget repeat expansion

Reference Genome Required Human Genome hg38, with Gencode v39 annotations

Target and repeat definitions Required PureTarget repeat expansion panel 2.0

Advanced Parameters

Reference Genome Required Human Genome hg38, with Gencode v39 annotations

- Specify a **reference genome** against which to align the reads
- Default set = **Human Genome hg38, with Gencode v39 annotations**

### PureTarget carrier panel

Application Required PureTarget carrier

Analysis Options

Add Analysis  YES  NO

Analysis Name Required

Select Analysis Workflow Required Target Enrichment

Reference Set Required Human Genome hg38, with Gencode v39 annotations

Target Set Required PureTarget carrier panel 1.0

Mark PCR duplicates  ON  OFF

Reference Set Required Human Genome hg38, with Gencode v39 annotations

- Specify a **reference set** against which to align the reads
- Default set = **Human Genome hg38, with Gencode v39 annotations**

### PureTarget custom panel

Application Required PureTarget custom

Analysis Options

Add Analysis  YES  NO

Analysis Name Required

Select Analysis Workflow Required Target Enrichment

Reference Set Required Human Genome hg38, with Gencode v39 annotations

Target Set Required

Mark PCR duplicates  ON  OFF

Reference Set Required Human Genome hg38, with Gencode v39 annotations

- Specify a **reference set** against which to align the reads
- Default set = **Human Genome hg38, with Gencode v39 annotations**

# SMRT Link Run Design options for PureTarget libraries – Analysis options

## 4. Specifying PureTarget analysis options for Revio and Vega systems – Target set selection<sup>1</sup>

### PureTarget repeat expansion panel

Application Required PureTarget repeat expansion

Analysis Options

Add Analysis  YES  NO

Analysis Name Required

Select Analysis Workflow Required PureTarget repeat expansion

Reference Genome Required Human Genome hg38, with Gencode v39 annotations

Target and repeat definitions Required PureTarget repeat expansion panel 2.0

Advanced Parameters

Target and repeat definitions Required PureTarget repeat expansion panel 2.0

- Specify a target and repeat definition browser extensible data (BED<sup>2</sup>) dataset
  - Default set = PureTarget repeat expansion panel 2.0

### PureTarget carrier panel

Application Required PureTarget carrier

Analysis Options

Add Analysis  YES  NO

Analysis Name Required

Select Analysis Workflow Required Target Enrichment

Reference Set Required Human Genome hg38, with Gencode v39 annotations

Target Set Required PureTarget carrier panel 1.0

Mark PCR duplicates  ON  OFF

Advanced Parameters

Target Set Required PureTarget carrier panel 1.0

- Specify a target definition BED<sup>2</sup> dataset
  - Default set = PureTarget carrier panel 1.0

### PureTarget custom panel

Application Required PureTarget custom

Analysis Options

Add Analysis  YES  NO

Analysis Name Required

Select Analysis Workflow Required Target Enrichment

Reference Set Required Human Genome hg38, with Gencode v39 annotations

Target Set Required < User-specified >

Mark PCR duplicates  ON  OFF

Advanced Parameters

Target Set Required < User-specified >

- Specify a target definition BED<sup>2</sup> dataset
  - User must manually specify the desired BED set (no default set is specified)
  - Recommend to include positive control target sequences corresponding to PureTarget control panel (103-633-300)

# SMRT Link Run Design options for PureTarget libraries – Analysis options

## 5. Specifying PureTarget analysis options for Revio and Vega systems – Mark PCR duplicates<sup>1</sup>

### PureTarget repeat expansion panel

Application Required PureTarget repeat expansion

Analysis Options

Add Analysis  YES  NO

Analysis Name Required

Select Analysis Workflow Required PureTarget repeat expansion

Reference Genome Required Human Genome hg38, with Gencode v39 annotations

Target and repeat definitions Required PureTarget repeat expansion panel 2.0

Advanced Parameters

No Mark PCR duplicates option

- Mark PCR duplicates option is not displayed
  - PureTarget repeat expansion analysis application workflow does not include mark PCR duplicates

### PureTarget carrier panel

Application Required PureTarget carrier

Analysis Options

Add Analysis  YES  NO

Analysis Name Required

Select Analysis Workflow Required Target Enrichment

Reference Set Required Human Genome hg38, with Gencode v39 annotations

Target Set Required PureTarget carrier panel 1.0

Mark PCR duplicates  ON  OFF

Advanced Parameters

Mark PCR duplicates  ON  OFF

- Mark PCR duplicates option is displayed
  - For PureTarget carrier analysis application workflow, Mark PCR duplicates = OFF (default)

### PureTarget custom panel

Application Required PureTarget custom

Analysis Options

Add Analysis  YES  NO

Analysis Name Required

Select Analysis Workflow Required Target Enrichment

Reference Set Required Human Genome hg38, with Gencode v39 annotations

Target Set Required

Mark PCR duplicates  ON  OFF

Advanced Parameters

Mark PCR duplicates  ON  OFF

- Mark PCR duplicates option is displayed
  - For PureTarget custom analysis application workflow, Mark PCR duplicates = OFF (default)

# SMRT Link Run Design options for PureTarget libraries – Analysis options

## 6. Specifying PureTarget analysis options for Revio and Vega systems – Advanced Parameters<sup>1</sup>

### PureTarget repeat expansion panel

Application Required PureTarget repeat expansion

Analysis Options

Add Analysis  YES  NO

Analysis Name Required

Select Analysis Workflow Required PureTarget repeat expansion

Reference Genome Required Human Genome hg38, with Gencode v39 annotations

Target and repeat definitions Required PureTarget repeat expansion panel 2.0

Advanced Parameters

Optionally specify Advanced Parameters

- Example Advanced Parameters:
  - Optionally specify sample karyotypes (XX or XY; default = XX for unlisted samples)<sup>2</sup>

### PureTarget carrier panel

Application Required PureTarget carrier

Analysis Options

Add Analysis  YES  NO

Analysis Name Required

Select Analysis Workflow Required Target Enrichment

Reference Set Required Human Genome hg38, with Gencode v39 annotations

Target Set Required PureTarget carrier panel 1.0

Mark PCR duplicates  ON  OFF

Advanced Parameters

Optionally specify Advanced Parameters

- Example Advanced Parameters:
  - Include Fail Reads (Leave default = ON)<sup>3</sup>
  - Perform variant calling (Leave default = OFF)<sup>4</sup>
  - Padding around regions (bp) (Leave default = 0)<sup>5</sup>

### PureTarget custom panel

Application Required PureTarget custom

Analysis Options

Add Analysis  YES  NO

Analysis Name Required

Select Analysis Workflow Required Target Enrichment

Reference Set Required Human Genome hg38, with Gencode v39 annotations

Target Set Required

Mark PCR duplicates  ON  OFF

Advanced Parameters

Optionally specify Advanced Parameters

- Example Advanced Parameters:
  - Include Fail Reads (Set to = ON)<sup>3</sup>
  - Perform variant calling (Leave default = OFF)<sup>4</sup>
  - Padding around regions (bp) (Set to = 0)<sup>5</sup>

<sup>1</sup> See *SMRT Link User Guide* ([Documentation](#)) for detailed descriptions of parameter settings for PureTarget analysis applications.

<sup>2</sup> Optionally specify sample karyotypes since ploidy is considered when genotyping X-chromosome repeats. CSV header should be 'biosample, karyotype', followed by one sample and karyotype per line.

<sup>3</sup> Include Fail Reads should be set to ON for PureTarget applications. Use of fail reads increases sensitivity to detect long repeat expansions in PureTarget libraries.

<sup>4</sup> Perform variant calling should be set to OFF for PureTarget applications. For PureTarget carrier panel analysis, optionally call variants using the PureTarget Carrier Pipeline analysis workflow.

<sup>5</sup> Padding around regions should be set to 0 for PureTarget applications. This parameter adds padding on each side of regions specified in the target BED file.

# SMRT Link Run Design options for PureTarget libraries – Analysis options

## 6. Specifying PureTarget analysis options for Revio and Vega systems – Advanced Parameters<sup>1</sup>

### PureTarget repeat expansion panel

Application Required PureTarget repeat expansion

Analysis Options

Add Analysis  YES  NO

Analysis Name Required

Select Analysis Workflow Required PureTarget repeat expansion

Reference Genome Required Human Genome hg38, with Gencode v39 annotations

Target and repeat definitions Required PureTarget repeat expansion panel 2.0

Advanced Parameters

Optionally specify **Advanced Parameters**

- Example Advanced Parameters:
  - Optionally specify sample karyotypes (XX or XY; default = XX for unlisted samples)<sup>2</sup>

Advanced Analysis Parameters

Sample karyotypes <sup>i</sup>  
Choose file Browse

Target order in boxplots <sup>i</sup>  
BED order

Add Task Memory (MB) <sup>i</sup>  
0

Compute Settings <sup>i</sup>  
-- select --

Ok Cancel

# SMRT Link Run Design options for PureTarget libraries – Analysis options

## 6. Specifying PureTarget analysis options for Revio and Vega systems – Advanced Parameters<sup>1</sup>

**PureTarget carrier panel**

Application Required: PureTarget carrier

**Analysis Options**

Add Analysis:  YES  NO

Analysis Name Required:

Select Analysis Workflow Required: Target Enrichment

Reference Set Required: Human Genome hg38, with Gencode v39 annotations

Target Set Required: PureTarget carrier panel 1.0

Mark PCR duplicates:  ON  OFF

**Advanced Parameters**

Optionally specify **Advanced Parameters**

### • Example Advanced Parameters:

- Include Fail Reads (Leave default = ON)<sup>2</sup>
- Perform variant calling (Leave default = OFF)<sup>3</sup>
- Padding around regions (bp) (Leave default = 0)<sup>4</sup>

**Advanced Analysis Parameters**

**Include Fail Reads** ⓘ  
 ON  OFF  
**Leave at default setting (ON)**

**Perform variant calling** ⓘ  
 ON  OFF  
**Leave at default setting (OFF)**

**Minimum Compressed Identity (%)** ⓘ

**Minimum Mapping Length (bp)** ⓘ

**Target order in boxplots** ⓘ  
BED order

**Advanced pbsv Options** ⓘ

**Padding around regions (bp)** ⓘ  
  
**Leave at default setting (0)**

**Minimum % of Reads that Support Variant (any one sample)** ⓘ

**Minimum Reads that Support Variant (total over all samples)** ⓘ

**Maximum number of parallel tasks** ⓘ

**Use GPU if available** ⓘ  
 ON  OFF

**Compute Settings** ⓘ  
-- select --

**Min. CCS Predicted Accuracy (Phred Scale)** ⓘ

**Identify Duplicates Across Sequencing Libraries** ⓘ  
 ON  OFF

**Minimum Length of Structural Variant (bp)** ⓘ

**Minimum Reads that Support Variant (any one sample)** ⓘ

**Add task memory (MB)** ⓘ

Ok Cancel

<sup>1</sup> See *SMRT Link User Guide* ([Documentation](#)) for detailed descriptions of parameter settings for PureTarget analysis applications.

<sup>2</sup> Include Fail Reads should be set to ON for PureTarget applications. Use of fail reads increases sensitivity to detect long repeat expansions in PureTarget libraries.

<sup>3</sup> Perform variant calling should be set to OFF for PureTarget applications. For PureTarget carrier panel analysis, optionally call variants using the PureTarget Carrier Pipeline analysis workflow.

<sup>4</sup> Padding around regions should be set to 0 for PureTarget applications. This parameter adds padding on each side of regions specified in the target BED file.

# SMRT Link Run Design options for PureTarget libraries – Analysis options

## 6. Specifying PureTarget analysis options for Revio and Vega systems – Advanced Parameters<sup>1</sup>

### PureTarget custom panel

Application  
Required PureTarget custom

Analysis Options

Add Analysis  YES  NO

Analysis Name  
Required

Select Analysis Workflow  
Required Target Enrichment

Reference Set Required  
Human Genome hg38, with Gencode v39 annotations

Target Set Required

Mark PCR duplicates  ON  OFF

Advanced Parameters

Optionally specify Advanced Parameters

- Example Advanced Parameters:
  - Include Fail Reads (Set to = ON)<sup>2</sup>
  - Perform variant calling (Leave default = OFF)<sup>3</sup>
  - Padding around regions (bp) (Set to = 0)<sup>4</sup>

Advanced Analysis Parameters

Include Fail Reads  ON  OFF

Perform variant calling  ON  OFF

Minimum Compressed Identity (%)

Minimum Length (bp)

Target order in boxplots  
BED order

Advanced pbsv Options

Minimum Reads that Support Variant (total over all samples)  
3

Maximum number of parallel tasks  
5120

Use GPU if available  ON  OFF

Compute Settings  
-- select --

Min. CCS Predicted Accuracy (Phred Scale)  
20

Identify Duplicates Across Sequencing Libraries  
 ON  OFF

Minimum Length of Structural Variant (bp)  
20

Minimum Reads that Support Variant (any one sample)  
3

Add task memory (MB)  
0

Ok Cancel

<sup>1</sup> See *SMRT Link User Guide* ([Documentation](#)) for detailed descriptions of parameter settings for PureTarget analysis applications.

<sup>2</sup> Include Fail Reads should be set to ON for PureTarget applications. Use of fail reads increases sensitivity to detect long repeat expansions in PureTarget libraries.

<sup>3</sup> Perform variant calling should be set to OFF for PureTarget applications. For PureTarget custom analysis, optionally perform call variants using other PacBio variant calling tools.

<sup>4</sup> Padding around regions should be set to 0 for PureTarget applications. This parameter adds padding on each side of regions specified in the target BED file.



# PureTarget data analysis workflow overview

# PureTarget data analysis software and workflow options



PureTarget applications leverage PacBio and PacBio Compatible Partner software for flexible secondary and tertiary analysis options

Run Design & Monitoring


Panel performance

Variant calling

Variant interpretation



```
##fileformat=VCFv4.2
##FILTER=<ID=PASS,Description="PASS">
##INFO=<ID=TRID,Number=1,Type="Text">
##INFO=<ID=END,Number=1,Type="Text">
##INFO=<ID=MOTIFS,Number=1,Type="Text">
##INFO=<ID=STRUC,Number=1,Type="Text">
##FORMAT=<ID=GT,Number=1,Type="Text">
```



Repeat expansion panel<sup>1</sup>

SMRT Link run design<sup>4</sup>

SMRT Link PureTarget repeat expansion

Carrier panel<sup>2</sup>

SMRT Link run design<sup>4</sup>

PureTarget carrier pipeline (PTCP)

Tertiary reporting



SMRT Link run design<sup>4</sup>

SMRT Link Target Enrichment

For lab techs to quickly assess sequencing performance

Custom panel<sup>3</sup>

SMRT Link run design<sup>4</sup>


SMRT Link Target Enrichment

<sup>1</sup> See *Application note – Comprehensive genotyping with the PureTarget repeat expansion panel and HiFi sequencing* ([102-326-614](#)) for descriptions of analysis software and workflow options.  
<sup>2</sup> See *Application note – Consolidate challenging genes with PureTarget carrier screen panel* ([102-326-653](#)) for descriptions of analysis software and workflow options.  
<sup>3</sup> See *Technical Note – A practical guide to amplification-free PureTarget custom panels* ([102-326-652](#)) for descriptions of analysis software and workflow options.  
<sup>4</sup> PureTarget run designs may be created using either SMRT Link (on-premise) or SMRT Link Cloud software.

# PureTarget repeat expansion panel 2.0 data analysis workflow overview

Analysis of PureTarget repeat expansion panel libraries can be performed in SMRT Link using the PureTarget repeat expansion analysis workflow or at the command line<sup>1</sup>

## PureTarget repeat expansion panel SMRT Link analysis workflow<sup>2</sup>



# PureTarget repeat expansion panel 2.0 data analysis workflow overview (cont.)

Target and repeat definitions specification for SMRT Link PureTarget repeat expansion panel QC evaluation and variant calling analysis<sup>1</sup>

## BED file target and repeat specification



PureTarget\_repeat\_expansion\_panel\_2.0.repeat\_definition.GRCh38.bed

### Example tandem repeat definition

```
chr4 3074876 3074966 ID=HTT,MOTIFS=CAG,CCG;STRUC=<TR>
```

- Repeat region has coordinates chr4:3074876-3074966
- Identifier is *HTT*
- **Note:** SMRT Link BED files must contain a fourth column with ID=NAME
- MOTIFS are the expected motifs the regions contain; CAG and CCG
- STRUC field can be always set to value <TR>






SMRT Link v25.3

SMRT Analysis

Default target BED set is PureTarget repeat expansion panel 2.0

- Set to run as “One analysis for all Data Sets” to summarize across indexed samples
- Only reads that map within the target regions in the BED file are included in the analysis<sup>2</sup>




<sup>1</sup> See *SMRT Link User Guide* ([Documentation](#)) for detailed descriptions of parameter settings for PureTarget repeat expansion analysis application.

<sup>2</sup> To “in-silico” mask data from targets included in the 38-gene PureTarget repeat expansion panel 2.0, create a new BED dataset without that target.


<sup>3</sup> Optionally specify sample karyotypes since ploidy is considered when genotyping X-chromosome repeats. CSV header should be ‘biosample, karyotype’, followed by one sample and karyotype per line.

# PureTarget repeat expansion panel 2.0 data analysis workflow overview (cont.)

## SMRT Link PureTarget repeat expansion analysis outputs – Target coverage QC plots




### Read Categories<sup>1</sup>




- Histogram mapping the length of On target, Off target, and Unmapped reads in the sample

### Target Coverage<sup>2</sup>



- Box plot for each target regions of mean coverage across all samples analyzed

### Sample Coverage<sup>2</sup>



- Box plot for each sample of mean coverage across all target regions

<sup>1</sup> **Note:** A read is defined as being on-target if its alignment region in the reference genome has a non-empty overlap with any defined target in the input BED file. Histogram bars for the different categories are stacked.

<sup>2</sup> **Note:** PureTarget repeat expansion panel application reports **mean read coverage** for target coverage and sample coverage QC plots.

# PureTarget repeat expansion panel 2.0 data analysis workflow overview (cont.)

## SMRT Link PureTarget repeat expansion analysis outputs – File downloads

SMRT Link v25.3

SMRT Analysis



Select Analysis Workflow

Required

PureTarget repeat expansion

▼Data

File Downloads



File downloads

Edit Output File Name Prefix

Example:analysis-GV-18-12105

File ↑	Type ⌵	Size
Mapped BAM	zip	462 MB
Repeat expansion genotypes	csv	185 KB
SMRT Link Log	log	40 KB
Sample summary	csv	4 KB
TRGT plots meth-allele	zip	46 MB
TRGT plots meth-waterfall	zip	49 MB
TRGT plots motifs-allele	zip	24 MB
TRGT plots motifs-waterfall	zip	26 MB
TRGT spanning BAM	zip	22 MB

### • Repeat expansion genotypes

- CSV file containing repeat unit sequence, allele count, min/max/consensus repeat array length, motif count, and motif span information.

### • TRGT plots meth-allele

- Depicts consensus repeat alleles and reads aligning to them. Bases in repeats are colored by methylation levels.

### • TRGT plots meth-waterfall

- Depicts portions of reads spanning the repeat without aligning them, which is convenient for showing mosaicism. Bases in repeats are colored by methylation levels.

### • TRGT plots motif-allele

- Depicts consensus repeat alleles and reads aligning to them. Bases in repeats are colored by repeat motif.

### • TRGT plots motif-waterfall

- Depicts portions of reads spanning the repeat without aligning them, which is convenient for showing mosaicism. Bases in repeats are colored by repeat motif.

# PureTarget repeat expansion panel 2.0 data analysis workflow overview (cont.)

## SMRT Link PureTarget repeat expansion analysis outputs – File downloads (cont.)

### trgt plot enables visualizations with TRGT

- Create **allele plots** that depict repeat alleles and reads aligning to them
- Create **waterfall plots** that depict portions of HiFi reads spanning the repeat without aligning them; waterfall plots are especially convenient for visualizing mosaicism.

Edit Output File Name Prefix

Example: analysis-GV-18-12105

File ↑	Type IT	Size
Mapped BAM	zip	462 MB
Repeat expansion genotypes	csv	185 KB
SMRT Link Log	log	40 KB
Sample summary	csv	4 KB
TRGT plots meth-allele	zip	46 MB
TRGT plots meth-waterfall	zip	49 MB
TRGT plots motifs-allele	zip	24 MB
TRGT plots motifs-waterfall	zip	26 MB
TRGT spanning BAM	zip	22 MB


### Repeat expansion genotypes

Sample	Bio_Sample_1
CNBP repeat unit	CAGG:CAGA:CA
CNBP read count allele 0	146
CNBP consensus size allele 0	139
CNBP min size allele 0	129
CNBP max size allele 0	317
CNBP motif counts allele 0	15_8_23
CNBP motif spans allele 0	0(0-61)_1(61-93)_2(93-139)


### PureTarget motifs analysis

■ CAG ■ CCG

#### Allele plot




#### Waterfall plot




### PureTarget methylation analysis

■ Methylated ■ Unmethylated




Profile at *FMR1* shows consistent methylation of expanded allele in a female carrier (NA07537) sample

#### Allele plot



#### Waterfall plot




For this sample, the majority of reads in both *ATXN7* alleles contain mostly unmethylated cytosines at CpG sites

Repeat structures are color-coded in allele and waterfall plots (*ATXN7* example shown)

# PureTarget carrier panel data analysis workflow overview

Evaluation of PureTarget carrier panel performance QC can be performed in SMRT Link using the Target Enrichment analysis workflow or at the command line<sup>1</sup>

## PureTarget carrier panel SMRT Link QC analysis workflow<sup>2</sup>



<sup>1</sup> See *Application note – Consolidate challenging genes with PureTarget carrier screening panel* (102-326-653) for detailed descriptions of analysis workflow options for PureTarget carrier applications.

<sup>2</sup> See *SMRT Link User Guide* ([Documentation](#)) for detailed descriptions of parameter settings for Target Enrichment analysis application.

# PureTarget carrier panel data analysis workflow overview (cont.)

## Target definitions specification for SMRT Link PureTarget carrier panel QC evaluation<sup>1</sup>

### BED file target specification



PureTarget\_carrier\_panel\_1.0.cut\_site.GRCh38.bed

### Example target cut site definition

```
chr4 3074876 3074966 ID=GBA
```

- Target cut site has coordinates chr4:3074876-3074966
- Identifier is *GBA*
- **Note:** SMRT Link BED files must contain a fourth column with ID=NAME
- See *SMRT Link user guide* for more details<sup>1</sup>

Set to run as "One analysis for all Data Sets"

PacBio Select a Module

SMRT Analysis / Create New Analysis

1. Select Data 2. Select Analysis

Analysis Application Required

Target Enrichment

Import Analysis Settings Export

Associated Inputs

Reference Set Required

Human Genome hg38, with Gencode v39 annotatic

Target Set Required

PureTarget carrier panel 1.0

Mark PCR duplicates

ON OFF

Advanced Parameters






SMRT Link v25.3

SMRT Analysis

Select PureTarget carrier panel 1.0 target BED set

Set Include Fail Reads = ON<sup>2</sup>  
Leave Perform variant calling = OFF<sup>3</sup>  
Set Padding around regions (bp) = 0<sup>4</sup>

- Set to run as "One analysis for all Data Sets" to summarize across indexed samples



<sup>1</sup> See *SMRT Link User Guide (Documentation)* for detailed descriptions of parameter settings for PureTarget Target Enrichment analysis application.

<sup>2</sup> Include Fail Reads should be set to ON for PureTarget applications. Use of fail reads increases sensitivity to detect long repeat expansions in PureTarget libraries.

<sup>3</sup> Perform variant calling should be set to OFF for PureTarget applications. For PureTarget carrier panel analysis, optionally perform genotyping using the PureTarget carrier panel analysis workflow.

<sup>4</sup> Padding around regions should be set to 0 for PureTarget applications. This parameter adds padding on each side of regions specified in the target BED file.

# PureTarget carrier panel data analysis workflow overview (cont.)

## SMRT Link PureTarget carrier panel QC analysis outputs – Target coverage QC plots

SMRT Link v25.3



Select Analysis Workflow  
Required

Target Enrichment

Target Enrichment

Summary Metrics

Read Categories


Target Coverage

Sample Coverage




QC plots & metrics

### Read Categories<sup>1</sup>




- Histogram mapping the length of On target, Off target, and Unmapped reads in the sample

### Target Coverage<sup>2</sup>



- Box plot for each target regions of mean coverage across all samples analyzed

### Sample Coverage<sup>2</sup>



- Box plot for each sample of mean coverage across all target regions

<sup>1</sup> **Note:** A read is defined as being on-target if its alignment region in the reference genome has a non-empty overlap with any defined target in the input BED file. Histogram bars for the different categories are stacked. For PureTarget samples, the number reads in the 'Duplicates' category is typically zero since PCR duplicates are not expected to be present in PureTarget libraries and Mark PCR Duplicates is therefore not performed..

<sup>2</sup> **Note:** For PureTarget carrier panels and PureTarget custom panels, SMRT Link Target Enrichment analysis workflow reports **mean base coverage** for target coverage and sample coverage QC plots.

# PureTarget carrier panel data analysis workflow overview (cont.)

## SMRT Link PureTarget carrier panel QC analysis outputs – File downloads

SMRT Link v25.3



Select Analysis Workflow

Required

Target Enrichment

▼Data

File Downloads



File downloads

Edit Output File Name Prefix

Example:analysis-[multiple]-123320

File	Size	Type
Mapped+Haplotagged BAM - BioSample146	26 MB	bam
Mapped+Haplotagged BAM - BioSample108	25 MB	bam
Mapped+Haplotagged BAM - BioSample162	36 MB	bam
Mapped+Haplotagged BAM Index - BioSample185	2 MB	bam_bai
Target Coverage By Sample CSV	22 KB	csv
Sample Summary CSV	8 KB	csv

### • Target Coverage By Sample CSV

- Comma-delimited text file of the matrix of all mean coverage values for each target (rows) and sample (columns). The data in this file is used to generate the Target Coverage and Sample Coverage box plot.

	A	B	C	D	E	F	G	H	I	J	K
1	targetName	BioSample100	BioSample101	BioSample102	BioSample103	BioSample104	BioSample105	BioSample106	BioSample107	BioSample108	BioSample109
2	GBAP	143.039226	235.618455	159.680491	205.643511	164.104372	187.5212	172.4631	201.5143	157.0218	179.8754
3	GBA	109.468538	183.412531	176.706672	152.876051	131.967724	185.0193	177.8427	140.387	122.852	102.0357
4	HBB	196.493551	295.325873	201.311077	282.075683	239.581563	258.0986	263.5472	273.1737	182.9837	230.0867
5	HBA1/hs-40	187.466172	257.382333	253.085188	192.052596	201.34637	248.0074	207.1164	244.3817	272.1106	253.6923
6	HBA1/2-smalldel	42.775765	50.690922	52.119821	86.559026	43.743181	37.81439	83.11514	80.20248	41.1167	78.1674
7	HBA1/2-mediumdel	14.634827	23.264869	35.464813	25.631649	18.096078	22.81532	20.14724	25.31156	31.40876	14.31733
8	HBA1/2-largedel	2.558176	2.955431	4.608815	3.48891	2.294498	5.938984	2.33284	4.378707	2.198251	3.434638
9	SMN2	46.470405	51.342962	29.593529	76.550765	0.248159	36.29135	50.87199	81.42325	26.48485	33.91766

### • Sample Summary CSV

- CSV version of the data displayed in the Target Enrichment Sample Summary table.

	A	B	C	D	E	F
1	Sample	Number of Reads	Median Read Length (bp)	Median Read Quality	Mean Target Coverage	Percent of Targets with ≥10-fold Coverage
2	BioSample97	11,241	4,951	Q43	211	100.00%
3	BioSample98	107,054	4,662	Q43	220.95	100.00%
4	BioSample99	1,093,548	3,542	Q47	197.15	100.00%
5	BioSample100	8,570	5,009	Q42	176.25	100.00%
6	BioSample101	12,882	4,946	Q42	221.05	100.00%
7	BioSample102	151,036	3,832	Q46	211.55	100.00%

# PureTarget carrier panel data analysis workflow overview (cont.)

PureTarget carrier panel variant calling analysis can be performed using PureTarget carrier pipeline (PTCP) analysis software available from GitHub and PacBio Compatible partners<sup>1</sup>

## PureTarget carrier panel variant calling analysis workflow<sup>2</sup>






# PureTarget carrier panel data analysis workflow overview (cont.)

## PureTarget carrier panel (PTCP) variant calling analysis setup

GitHub

Customer  
HPC



Deployment



DNAnexus  
platform



DNAnexus

### Configure workflow

```
git clone \
https://github.com/PacificBiosciences
/ptcp.git
cd ptcp
conda create -n ptcp \
python=3.12.11
conda activate ptcp python -m pip
install \
-r requirements.txt
```

```
docker pull quay.io/pacbio/ptcp:latest
```

### Create Sample Sheet

```
bam_name,bam_id,sex
movie.bc2042.bam,movie.bc2042,F
movie.bc2043.bam,movie.bc2043,F
movie.bc2044.bam,movie.bc2044,F
movie.bc2045.bam,movie.bc2045,F
movie.bc2046.bam,movie.bc2046,F
movie.bc2047.bam,movie.bc2047,F
movie.bc2048.bam,movie.bc2048,F
movie.bc2049.bam,movie.bc2049,F
```


### Create inputs JSON

```
{
  "ptcp.sample_sheet": "inputs/sample_sheet.csv",
  "ptcp.ref_fasta": "data/reference/GRCh38.fa",
  "ptcp.ref_index": "data/reference/GRCh38.fa.fai",
  "ptcp.trgt_bed": "meta/repeat_definition.GRCh38.bed",
  "ptcp.paraphase_config_yaml": "meta/config.GRCh38.yaml",
  "ptcp.genome_version": "38",
  "ptcp.ptcp_qc_bed": "meta/ptcp-qc.GRCh38.bed",
  "ptcp.hifi_reads": [
    "/data/reads/SAMPLE1.hifi_reads.bc2042.bam"
  ],
  "ptcp.fail_reads": [
    "/data/reads/SAMPLE1.fail_reads.bc2042.bam"
  ]
}
```

### Upload data

- Import from AWS S3
- Local import using Upload Agent


### Select workflow



### Create Sample Sheet

```
bam_name,bam_id,sex
movie.bc2042.bam,movie.bc2042,F
movie.bc2043.bam,movie.bc2043,F
movie.bc2044.bam,movie.bc2044,F
movie.bc2045.bam,movie.bc2045,F
movie.bc2046.bam,movie.bc2046,F
movie.bc2047.bam,movie.bc2047,F
movie.bc2048.bam,movie.bc2048,F
movie.bc2049.bam,movie.bc2049,F
```

### Select inputs




# PureTarget carrier panel data analysis workflow overview (cont.)

## PureTarget carrier panel (PTCP) variant calling analysis execution and outputs

GitHub

Customer  
HPC



Deployment



DNAnexus  
platform




### Run workflow

```
conda activate ptcp
cd /Path/to/ptcp

miniwdl run \
  --verbose \
  --dir /Path/to/output_dir \
  --cfg /Path/to/miniwdl.cfg \
  --input /Path/to/ptcp_inputs.json \
  main.wdl
```

### Start analysis



### Per sample outputs for tandem repeat analysis, gene phasing, and structural variant detection

```
bc2048.mapped.bam      <- pbmm2 mapped reads
bc2048.mapped.bam.bai
bc2048.repeats.bam     <- pbmm2 mapped reads intersecting the repeats
bc2048.repeats.bam.bai
bc2048.trgt.sorted.spanning.bam  <- TRGT BAM with spanning reads
bc2048.trgt.sorted.spanning.bam.bai
bc2048.trgt.vcf        <- Tandem repeat genotypes generated by TRGT
bc2048.meth_allele.trgt_plots.zip <- Tandem repeat plots generated by TRGT
bc2048.meth_waterfall.trgt_plots.zip <- ...
bc2048.motifs_allele.trgt_plots.zip <- ...
bc2048.motifs_waterfall.trgt_plots.zip <- ...
bc2048.f8inversion.json <- F8 inversion calls with extended information
bc2048.f8inversion.vcf <- F8 inversion calls
bc2048.paraphase.bam   <- Paraphase re-aligned BAM with annotations
bc2048.paraphase.bam.bai
bc2048.paraphase.json  <- Paraphase results per target region
bc2048_cyp21.vcf      <- Paraphase VCF with small variant calls per target per haplotype
bc2048.havanno.json   <- Annotated small variants called by Paraphase
```

### Coverage and quality metrics per sample and aggregated report across all sample

```
qc.aggregate.json <- per target coverage
qc.bc2045.json    <- per sample per target coverage
qc.bc2046.json    <- ...
qc.bc2047.json    <- ...
qc.bc2048.json    <- ...
qc.bc2049.json    <- ...
```

# PureTarget carrier panel data analysis workflow overview (cont.)

PureTarget carrier panel (PTCP) variant calling analysis outputs **read count** as a measure of depth

## Coverage and quality metrics per sample and aggregated report across all sample

```
qc.aggregate.json <- per target
qc.bc2045.json <- per sample+target
qc.bc2046.json <- ...
qc.bc2047.json <- ...
qc.bc2048.json <- ...
qc.bc2049.json <- ...
```

```
... qc.aggregate.json
"hbh": {
  "samples_with_coverage": 16,
  "total_reads": 4388,
  "coverage_stats": {
    "min": 130.0,
    "q1": 234.5,
    "median": 271.0,
    "mean": 274.25,
    "q3": 318.0,
    "max": 421.0
  }
},
...
```


```
... qc.bc2045.json
{
  "sample_name": "bc2045",
  "genome_version": "38",
  "targets_bed": "ptcp-qc.GRCh38.bed",
  "ptcp-qc_version": "0.8.0",
  "timestamp": "2025-09-16T03:09:18.573904366+00:00",
  "stats": {
    "total_reads": 23230,
    "on_target": {
      "total": 2773,
      "hifi": 2445,
      "fail": 328
    },
    "off_target": {
      "total": 20402,
      "hifi": 19195,
      "fail": 1207
    },
    "unmapped": {
      "total": 55,
      "hifi": 5,
      "fail": 50
    },
    "loci": {
      "count": 24,
      "coverages": {
        "stats": {
          "min": 1.0,
          "q1": 2.0,
          "median": 47.0,
          "mean": 115.541664,
          "q3": 215.75,
          "max": 400.0
        }
      }
    }
  }
},
...
```

```
... qc.bc2045.json
...
  },
  "paraphase_results": {
    "gba": {
      "total_reads": 334,
      "unique_reads": 334,
      "total_cn": 4,
      "haplotypes": {
        "gba_hap1": {
          "total_count": 58,
          "unique_count": 58,
          "fractional_count": 58.0
        },
        "gba_hap2": {
          "total_count": 201,
          "unique_count": 201,
          "fractional_count": 201.0,
          "n_copy": 2
        },
        "gba_hap3": {
          "total_count": 75,
          "unique_count": 75,
          "fractional_count": 75.0
        }
      }
    }
  }
},
...
}
```

# PureTarget custom panel data analysis workflow overview

Evaluation of PureTarget custom panel performance QC can be performed in SMRT Link using the Target enrichment analysis workflow or at the command line<sup>1</sup>

## PureTarget custom panel SMRT Link QC analysis workflow<sup>2</sup>



# PureTarget custom panel data analysis workflow overview (cont.)

## Target definitions specification for SMRT Link PureTarget custom panel QC evaluation<sup>1</sup>

### BED file target specification



\*. PureTarget\_control\_panel.cut\_site.GRCh38.bed

### Example target cut site definitions for PureTarget control panel

chrX	154482054	154488569	ID=UBL4A
chr7	5525920	5532077	ID=ACTB
chr12	6534182	6539602	ID=GAPDH

- Target cut site are located on Chr X, Chr 7 and Chr 12
- Identifiers are *UBL4A*, *ACTB* and *GAPDH*
- Note: SMRT Link BED files must contain a fourth column with ID=NAME. See the user guide for more details<sup>1</sup>

Set to run as “One analysis for all Data Sets”



SMRT Link v25.3

SMRT Analysis

Specify to use the desired user-created target BED set

PacBio Select a Module

SMRT Analysis / Create New Analysis

1. Select Data 2. Select Analysis

Analysis Application Required

Target Enrichment

Import Analysis Settings Export

Associated Inputs

Reference Set Required

Human Genome hg38, with Gencode v39 annotatic

Target Set Required




< User-specified >

Mark PCR duplicates

ON OFF

Advanced Parameters

- Set to run as “One analysis for all Data Sets” to summarize across indexed samples
- To analyze data from a custom panel, create a new BED dataset containing positive control target sequences corresponding to **PureTarget control panel** (103-633-300) along with additional desired custom targets<sup>2</sup>



# PureTarget custom panel data analysis workflow overview (cont.)

## SMRT Link PureTarget custom panel QC analysis outputs

SMRT Link v25.3

SMRT Analysis



Select Analysis Workflow

Required

Target Enrichment

Target Enrichment

Summary Metrics

Read Categories


Target Coverage

Sample Coverage




QC plots & metrics

### Read Categories<sup>1</sup>




- Histogram mapping the length of On target, Off target, and Unmapped reads in the sample

### Target Coverage<sup>2</sup>



- Box plot for each target regions of mean coverage across all samples analyzed

### Sample Coverage<sup>2</sup>



- Box plot for each sample of mean coverage across all target regions

<sup>1</sup> **Note:** A read is defined as being on-target if its alignment region in the reference genome has a non-empty overlap with any defined target in the input BED file. Histogram bars for the different categories are stacked. For PureTarget samples, the number reads in the 'Duplicates' category is typically zero since PCR duplicates are not expected to be present in PureTarget libraries and Mark PCR Duplicates is therefore not performed.

<sup>2</sup> **Note:** For PureTarget carrier panels and PureTarget custom panels, SMRT Link Target Enrichment analysis workflow reports **mean base coverage** for target coverage and sample coverage QC plots.

# PureTarget custom panel data analysis workflow overview (cont.)

## SMRT Link PureTarget custom panel QC analysis outputs – File downloads



Edit Output File Name Prefix **Example:analysis-[multiple]-123320**

File	Size	Type
Mapped+Haplotagged BAM - BioSample146	26 MB	bam
Mapped+Haplotagged BAM - BioSample108	25 MB	bam
Mapped+Haplotagged BAM - BioSample162	36 MB	bam
Mapped+Haplotagged BAM Index - BioSample185	2 MB	bam_bai
Target Coverage By Sample CSV	22 KB	csv
Sample Summary CSV	8 KB	csv
Sample Mapped BAMs	5 GB	zip
Analysis Log	441 KB	log
SMRT Link Log	28 KB	log

- **Target Coverage By Sample CSV**
  - Comma-delimited text file of the matrix of all mean coverage values for each target (rows) and sample (columns). The data in this file is used to generate the Target Coverage and Sample Coverage box plot.

	A	B	C	D	E	F	G	H	I
1	targetName	BioSample1	BioSample2	BioSample3	BioSample4	BioSample5	BioSample6	BioSample7	BioSample8
2	UBL4A	690.897467	679.900844	774.944129	765.051266	698.799233	797.569916	753.020107	741.452341
3	ACTB	1290.993503	1394.134319	1493.070814	1414.276108	1371.073412	1434.178983	1348.245737	1484.30242
4	GAPDH	1777.403321	1812.297048	2021.746494	2020.251476	2075.365129	1974.968635	1951.868081	1901.698708
5									
6									
7									
8									

- **Sample Summary CSV**
  - CSV version of the data displayed in the Target Enrichment Sample Summary table.

	A	B	C	D	E	F
1	Sample	Number of Reads	Median Read Length (bp)	Median Read Quality	Mean Target Coverage	Percent of Targets with a%¥10-fold Coverage
2	BioSample1	12,576	5,419	Q41	1236.67	100.00%
3	BioSample2	12,188	5,419	Q41	1275.67	100.00%
4	BioSample3	12,580	5,419	Q41	1408	100.00%
5	BioSample4	12,592	5,419	Q40	1373	100.00%
6	BioSample5	12,428	5,419	Q41	1354.67	100.00%
7	BioSample6	13,134	5,419	Q40	1354	100.00%

# SMRT Link Target Enrichment analysis application is a generalizable workflow that can be used to evaluate targeted sequencing performance QC

Recommended Target Enrichment analysis settings for amplicon, hybrid capture, and PureTarget sequencing applications

Targeted sequencing application	Target Enrichment analysis options			
	Include Fail Reads	Padding around regions (bp)	Mark PCR Duplicates	Variant Calling <sup>1</sup>
Amplicon sequencing	OFF (default)	0	OFF (default)	Optional
Hybrid capture sequencing <sup>2</sup>	OFF (default)	3000 (default)	ON	Optional
PureTarget sequencing (for PureTarget carrier screening & custom panels)	ON	0	OFF (default)	OFF (default)

<sup>1</sup> Singularity is only required if Variant Calling is performed.

<sup>2</sup> **Note:** For QIAGEN hybrid capture long read panels, HiFi data analysis should be performed using QIAGEN software tools (e.g., [QIAGEN GeneGlobe Data Analysis portal](#) offers simplified downstream data analysis for HLA haplotyping and [QIAGEN CLC Genomics Workbench](#) is used to detect large structural variants in hereditary cancers).




# PureTarget example sequencing performance data

# Example PureTarget library preparation QC results

Example PureTarget repeat expansion panel library QC results for human blood gDNA samples

## Input genomic DNA sizing QC



Example Femto Pulse genomic DNA sizing QC analysis results for high-molecular weight genomic DNA extracted from a human whole blood sample using Nanobind PanDNA kit.

## PureTarget library preparation step yields

Library preparation step	DNA or SMRTbell overall recovery	
	Manual	Automated
Starting input genomic DNA	2,000 ng per sample	1,300 ng per sample
Post-gDNA repair & cleanup	1600 ng (80%) per sample	1040 ng (80%) per sample
Post-Cas9 digestion & cleanup	1200 ng (60%) per sample	780 ng (60%) per sample
Post-Adapter ligation & cleanup	1000 ng (50%) per sample	650 ng (50%) per sample
Post-nuclease treatment & cleanup <sup>1</sup>	288 ng (0.3%) per 48-plex	124.8 ng (0.1%) per 96-plex
Post-PureTarget library cleanup	96 ng (0.1%) per 48-plex	62.4 ng (0.05%) per 96-plex

Example library preparation yield results for PureTarget repeat expansion library prepared from a human whole blood sample using manual or automated workflows. For manual workflows, expected overall recovery is typically ~0.02-0.2% relative to starting gDNA total mass. For automated workflows, expected overall recovery is typically ~0.005-0.05%.<sup>1,2</sup>

Final PureTarget library yield is typically sufficient to load 1 SMRT Cell on Revo system + SPRQ chemistry or 1 SMRT Cell on Vega system


<sup>1</sup> **Note:** Manual PureTarget library prep workflow using PureTarget kit 24 employs one round of nuclease treatment. Automated PureTarget library prep workflow using PureTarget kit 96 employs two rounds of nuclease treatment.

<sup>2</sup> **Note:** It possible to observe PureTarget library prep recoveries outside of these expected ranges and still obtain acceptable sequencing yields and good target coverage performance.


# Example 96-plex PureTarget repeat expansion panel library sequencing performance (Revio system + SPRQ chemistry)

96-plex PureTarget repeat expansion panel example data for human whole blood DNA samples<sup>1</sup>


## Raw Data Report



## HiFi Read Length



## Target Enrichment Summary Metrics



Raw Base Yield	395 Gb
Mean Polymerase Read Length	81.3 kb
P0	80%
P1	19%
P2	0%

Example metrics for a 96-plex PureTarget repeat expansion panel library run on a Revio system with SPRQ chemistry using a 24-hrs movie time. Revio system *P1* range for 96-plex PureTarget libraries was typically ~15%–40%.

HiFi Reads	2.5 M
HiFi Base Yield	11.4 Gb
Mean HiFi Read Length	4.5 kb
Median HiFi Read Quality	Q49
HiFi Read Mean # of Passes	33

For 96-plex PureTarget libraries, per-Revio SMRT Cell HiFi read counts were typically ~1 M–2 M depending on the final library insert size and *P1* loading performance.


Total Bases*	11,812,292,440
Total Reads*	2,566,000
Median Read Length*	4,152 bp
Median Read Quality*	Q48
Sample Count	96
Target Regions	38

\* Includes HiFi + non-HiFi data. For 96-plex PureTarget libraries, per-Revio SMRT Cell total read counts were typically ~1 M–2 M.


# Example 96-plex PureTarget carrier panel library sequencing performance (Revio system + SPRQ chemistry)

96-plex PureTarget carrier panel example data for human whole blood DNA samples<sup>1</sup>


## Raw Data Report



## HiFi Read Length



## Target Enrichment Summary Metrics



Raw Base Yield	362 Gb
Mean Polymerase Read Length	53.7 kb
P0	73%
P1	27%
P2	1%

Example metrics for a 96-plex PureTarget carrier panel library run on a Revio system with SPRQ chemistry using a 24-hrs movie time. Revio system P1 range for 96-plex PureTarget libraries was typically ~15%–40%.

HiFi Reads	2.2 M
HiFi Base Yield	12.2 Gb
Mean HiFi Read Length	5.5 kb
Median HiFi Read Quality	Q46
HiFi Read Mean # of Passes	31

For 96-plex PureTarget libraries, per-Revio SMRT Cell HiFi read counts were typically ~1 M–2 M depending on the final library insert size and P1 loading performance.


Total Bases*	13,043,277,783
Total Reads*	2,337,184
Median Read Length*	4,281 bp
Median Read Quality*	Q45
Sample Count	96
Target Regions	20

\* Includes HiFi + non-HiFi data. For 96-plex PureTarget libraries, per-Revio SMRT Cell total read counts were typically ~1 M–2 M.


# Example 48-plex PureTarget repeat expansion panel library sequencing performance (Revio system + SPRQ chemistry)

48-plex PureTarget repeat expansion panel example data for human whole blood DNA samples<sup>1</sup>


## Raw Data Report



## HiFi Read Length



## Target Enrichment Summary Metrics



Raw Base Yield	202 Gb
Mean Polymerase Read Length	56.1 kb
P0	85%
P1	14%
P2	0%

Example metrics for a 48-plex PureTarget repeat expansion panel library run on a Revio system with SPRQ chemistry using a 24-hrs movie time. Revio system *P1* range for 48-plex PureTarget libraries was typically ~15%–50%.

HiFi Reads	1.1 M
HiFi Base Yield	7.5 Gb
Mean HiFi Read Length	6.6 kb
Median HiFi Read Quality	Q43
HiFi Read Mean # of Passes	26

For 48-plex PureTarget libraries, per-Revio SMRT Cell HiFi read counts were typically ~1 M–2 M depending on the final library insert size and *P1* loading performance.


Total Bases*	7,692,021,081
Total Reads*	1,160,267
Median Read Length*	4,783 bp
Median Read Quality*	Q42
Sample Count	48
Target Regions	38

\* Includes HiFi + non-HiFi data. For 48-plex PureTarget libraries, per-Revio SMRT Cell total read counts were typically ~1 M–2 M.


# Example 48-plex PureTarget carrier panel library sequencing performance (Revio system + SPRQ chemistry)

48-plex PureTarget carrier panel example data for human whole blood DNA samples<sup>1</sup>


## Raw Data Report



## HiFi Read Length



## Target Enrichment Summary Metrics



Raw Base Yield	194 Gb
Mean Polymerase Read Length	55.2 kb
P0	86%
P1	14%
P2	0%

Example metrics for a 48-plex PureTarget carrier panel library run on a Revio system with SPRQ chemistry using a 24-hrs movie time. Revio system P1 range for 48-plex PureTarget libraries was typically ~15%–50%.

HiFi Reads	1.1 M
HiFi Base Yield	8.8 Gb
Mean HiFi Read Length	8.3 kb
Median HiFi Read Quality	Q39
HiFi Read Mean # of Passes	22

For 48-plex PureTarget libraries, per-Revio SMRT Cell HiFi read counts were typically ~1 M–2 M depending on the final library insert size and P1 loading performance.


Total Bases*	9,099,180,991
Total Reads*	1,076,225
Median Read Length*	6,825 bp
Median Read Quality*	Q37
Sample Count	48
Target Regions	20

\* Includes HiFi + non-HiFi data. For 48-plex PureTarget libraries, per-Revio SMRT Cell total read counts were typically ~1 M–2 M.


# Example 48-plex PureTarget repeat expansion panel library sequencing performance (Vega system)

48-plex PureTarget repeat expansion panel example data for human blood & cell line DNA samples<sup>1</sup>


## Raw Data Report



## HiFi Read Length



## Target Enrichment Summary Metrics



Raw Base Yield	167 Gb
Mean Polymerase Read Length	47.2 kb
Loading level	22%

Example metrics for a 16-plex PureTarget repeat expansion panel library run on a Vega system using a 24-hrs movie time. Vega system loading level range for 16-plex PureTarget libraries was typically ~20–40%.

HiFi Reads	1.0 M
HiFi Base Yield	6.4 Gb
Mean HiFi Read Length	6.2 kb
Median HiFi Read Quality	Q41
HiFi Read Mean # of Passes	24

For 16-plex PureTarget libraries, per-Vega SMRT Cell HiFi read counts were typically ~1 M–2 M depending on the final library insert size and loading level performance.

Total Bases*	6,548,671,466
Total Reads*	1,042,657
Median Read Length*	4,799 bp
Median Read Quality*	Q40
Sample Count	48
Target Regions	38


\* Includes HiFi + non-HiFi data. For 16-plex PureTarget libraries, per-Vega SMRT Cell total read counts were typically ~1 M–2 M.

<sup>1</sup> Human blood DNA samples were extracted using PacBio Nanobind kits with or without RBC lysis. Cell line DNA samples were Coriell lymphoblastoid cell DNA with GQN30kb > 5. 2 µg of input DNA per sample was used for PureTarget library preparation. **Note:** Each sample was processed with both the PureTarget repeat expansion panel 2.0 (103-633-100) and the PureTarget control panel (103-633-300). Read lengths, reads/data per SMRT Cell and other sequencing performance results can vary depending on DNA sample quality, insert size, sample loading performance & movie time. Note: Shorter library insert sizes (<5 kb), lower DNA quality samples, and suboptimal sample loading performance may result in lower data yields per SMRT Cell.


# Example 16-plex PureTarget repeat expansion panel library sequencing performance (Vega system)

16-plex PureTarget repeat expansion panel example data for human cell line DNA samples<sup>1</sup>


## Raw Data Report



## HiFi Read Length



## Target Enrichment Summary Metrics



Raw Base Yield	93 Gb
Mean Polymerase Read Length	54.0 kb
Loading level	9%

Example metrics for a 16-plex PureTarget repeat expansion panel library run on a Vega system using a 24-hrs movie time. Vega system loading level range for 16-plex PureTarget libraries was typically ~10–15%.

HiFi Reads	0.5 M
HiFi Base Yield	3.4 Gb
Mean HiFi Read Length	6.8 kb
Median HiFi Read Quality	Q39
HiFi Read Mean # of Passes	27

For 16-plex PureTarget libraries, per-Vega SMRT Cell HiFi read counts were typically ~0.5 M–1 M depending on the final library insert size and loading level performance.


Total Bases*	3,198,060,800
Total Reads*	469,669
Median Read Length*	4,880 bp
Median Read Quality*	Q39
Sample Count	16
Target Regions	38

\* Includes HiFi + non-HiFi data. For 16-plex PureTarget libraries, per-Vega SMRT Cell total read counts were typically ~0.5 M–1 M.


# Example 16-plex PureTarget carrier panel library sequencing performance (Vega system)

16-plex PureTarget carrier panel example data for human cell line DNA samples<sup>1</sup>


## Raw Data Report



## HiFi Read Length



## Target Enrichment Summary Metrics



Raw Base Yield	126 Gb
Mean Polymerase Read Length	58.2 kb
Loading level	11%

Example metrics for a 16-plex PureTarget carrier panel library run on a Vega system using a 24-hrs movie time. Vega system loading level range for 16-plex PureTarget libraries was typically ~10–15%.

HiFi Reads	0.7 M
HiFi Base Yield	5.2 Gb
Mean HiFi Read Length	7.5 kb
Median HiFi Read Quality	Q37
HiFi Read Mean # of Passes	23

For 16-plex PureTarget libraries, per-Vega SMRT Cell HiFi read counts were typically ~0.5 M–1 M depending on the final library insert size and loading level performance.


Total Bases*	5,132,040,788
Total Reads*	666,609
Median Read Length*	6,000 bp
Median Read Quality*	Q36
Sample Count	16
Target Regions	20

\* Includes HiFi + non-HiFi data. For 16-plex PureTarget libraries, per-Vega SMRT Cell total read counts were typically ~0.5 M–1 M.

# Recommended guidance for evaluating PureTarget repeat expansion panel sequencing run performance


When evaluating PureTarget runs, it is generally more useful to examine the **secondary analysis results** (e.g., on-target coverage) since primary sequencing metrics like Productivity ( $P0$ ,  $P1$ ,  $P2$ ) or Loading level are mostly dominated by 'background' non-targeted reads

Read Categories<sup>1</sup>




- Histogram mapping the length of On target, Off target, and Unmapped reads in the sample

Target Coverage<sup>2</sup>



- Box plot for each target regions of mean coverage across all samples analyzed

Sample Coverage<sup>2</sup>



- Box plot for each sample of mean coverage across all target regions

# Recommended guidance for evaluating PureTarget repeat expansion panel sequencing run performance (cont.)

For Revio system, Barcode Counts preview metrics<sup>1</sup> in SMRT Link Run Details report are useful for early evaluation of PureTarget sample demultiplexing performance and per-sample mean HiFi read length

The screenshot shows the PacBio Run Details interface. The 'Overview' section displays run metadata such as 'Run Created', 'Run Start', 'Run Complete', 'Created By', 'Instrument Name', 'Completed Cells', 'Failed Cells', 'Time remaining for PostProcessing', 'Transfer Status', 'Run ID', 'Instrument SN', 'Instrument Control SW Version', 'Instrument Chemistry Build', and 'Primary SW Version'. Below this is a 'Consumables' table. The 'Barcode Counts' section is highlighted with a pink box and shows a table with columns for Barcode ID, HiFi reads, and HiFi read length (mean, bp).

This section provides a detailed view of the Barcode Counts table. It includes filters for 'Well name' (PHX\_SEQ-12) and 'Time point' (23 hr). A search bar is present. The table lists various barcode IDs, their corresponding HiFi reads percentage, and their mean HiFi read length in base pairs. Three callouts provide definitions for the terms used in the table:

- Barcode:** An individual barcode detected in the sample, as well as unbarcoded reads
- HiFi reads:** An estimate of the % of reads with each barcode, as well as the % of unbarcoded reads
- HiFi read length, mean:** An estimate of the average HiFi read length for each barcode or for unbarcoded reads

Barcode ID	HiFi reads	HiFi read length (mean, bp)
bc2002--bc2002	7.1%	5747
bc2003--bc2003	3.5%	9374
bc2004--bc2004	4.7%	8112
bc2005--bc2005	3.5%	5707
bc2007--bc2007	2.4%	6627
bc2008--bc2008	2.4%	6063
bc2009--bc2009	1.2%	8257
bc2011--bc2011	3.5%	6463
bc2012--bc2012	4.7%	10685
bc2013--bc2013	1.2%	3381
bc2014--bc2014	2.4%	4243
bc2047--bc2047	3.5%	4884
bc2048--bc2048	1.2%	4294
Not Barcoded	5.9%	6508

# Recommended guidance for evaluating PureTarget repeat expansion panel sequencing run performance (cont.)

View Barcode demultiplexing results in SMRT Link SMRT Analysis to perform more detailed evaluation of PureTarget sample demultiplexing performance

**PacBio Data Management**  
 PHX\_SEQ-12-Cell1 (all samples)  
 Copy Analyze Export Delete

**Barcode Data**

Sample Name	Barcode	Barcode Qua...	HiFi Reads	HiFi Read Le...	HiFi Read Q...	HiFi Yield (bp)	Polymerase
Bio Sample 1	bc2001--bc20...	97.9	23,838	6,498	Q39	154,912,577	121,916
Bio Sample 2	bc2002--bc20...	96.3	23,829	6,606	Q38	157,423,933	121,269
Bio Sample 3	bc2003--bc20...	97.7	25,236	6,484	Q38	163,630,897	120,873
Bio Sample 4	bc2004--bc20...	96.8	23,954	6,630	Q38	158,825,582	121,598
Bio Sample 5	bc2005--bc20...	97.0	24,809	6,538	Q38	162,760,987	121,678
Bio Sample 6	bc2006--bc20...	97.3	22,362	6,738	Q38	152,786,566	121,731
Bio Sample 7	bc2007--bc20...	97.7	23,388	6,655	Q38	141,695,083	120,761
Bio Sample 8	bc2008--bc20...	97.5	22,407	6,662	Q38	150,100,389	122,274
Bio Sample 9	bc2009--bc20...	97.4	27,277	6,828	Q38	186,250,674	121,427
Bio Sample 10	bc2010--bc20...	97.3	22,084	6,641	Q38	146,668,330	120,977
Bio Sample 11	bc2011--bc20...	97.1	22,351	6,539	Q38	146,158,985	120,805
Bio Sample 12	bc2012--bc20...	97.6	21,407	6,594	Q39	141,171,468	121,334
Bio Sample 13	bc2013--bc20...	97.4	22,347	6,675	Q39	147,845,367	121,793
Bio Sample 14	bc2014--bc20...	97.6	24,733	6,561	Q39	162,281,425	121,368
Bio Sample 15	bc2015--bc20...	97.6	24,921	6,664	Q38	166,086,764	122,085
Bio Sample 16	bc2016--bc20...	97.5	19,922	6,528	Q38	130,058,233	121,249
Bio Sample 17	bc2017--bc20...	97.5	22,234	7,150	Q38	158,979,009	122,790
Bio Sample 18	bc2018--bc20...	97.9	25,923	6,775	Q38	175,644,772	121,687
Bio Sample 19	bc2019--bc20...	97.5	23,040	6,747	Q38	155,467,070	122,291
Bio Sample 20	bc2020--bc20...	97.9	26,138	6,646	Q38	173,730,992	121,650
Bio Sample 21	bc2021--bc20...	97.4	24,954	6,750	Q38	168,449,825	122,122
Bio Sample 22	bc2022--bc20...	97.9	22,402	6,627	Q38	148,478,815	121,233
Bio Sample 23	bc2023--bc20...	97.4	22,389	6,692	Q38	149,847,373	121,316
Bio Sample 24	bc2024--bc20...	97.8	24,244	6,677	Q38	161,888,862	121,594
Bio Sample 25	bc2025--bc20...	97.5	14,985	6,874	Q39	100,264,900	120,846
Bio Sample 26	bc2026--bc20...	97.4	15,678	7,042	Q39	110,418,581	121,313
Bio Sample 27	bc2027--bc20...	96.3	14,540	7,040	Q39	102,663,839	122,298
Bio Sample 28	bc2028--bc20...	97.8	14,674	6,879	Q39	100,948,914	121,355
Bio Sample 29	bc2029--bc20...	96.0	16,160	7,189	Q39	116,341,083	121,688
Bio Sample 30	bc2030--bc20...	97.5	15,112	7,117	Q39	122,556,883	121,239
Bio Sample 31	bc2031--bc20...	97.4	15,283	7,312	Q39	111,762,309	122,443
Bio Sample 32	bc2032--bc20...	97.3	14,740	7,177	Q39	105,480,469	121,769
Bio Sample 33	bc2033--bc20...	97.3	15,329	6,873	Q39	106,839,195	121,794
Bio Sample 34	bc2034--bc20...	97.4	15,866	7,034	Q39	111,621,343	121,645
Bio Sample 35	bc2035--bc20...	97.3	15,809	7,117	Q39	112,514,643	122,382
Bio Sample 36	bc2036--bc20...	97.5	14,012	7,077	Q39	99,166,762	121,531
Bio Sample 37	bc2037--bc20...	97.4	15,376	7,244	Q39	111,695,500	122,290
Bio Sample 38	bc2038--bc20...	97.6	14,470	7,544	Q39	103,382,534	122,633
Bio Sample 39	bc2039--bc20...	97.5	15,277	7,286	Q39	111,376,012	121,992
Bio Sample 40	bc2040--bc20...	96.7	15,195	7,323	Q39	115,281,396	122,716
Bio Sample 41	bc2041--bc20...	97.2	13,970	6,889	Q39	93,224,038	120,817
Bio Sample 42	bc2042--bc20...	97.5	12,714	6,726	Q39	85,516,208	120,420
Bio Sample 43	bc2043--bc20...	97.5	14,753	7,045	Q39	103,035,396	121,444
Bio Sample 44	bc2044--bc20...	97.6	15,508	7,169	Q39	111,183,795	121,566
Bio Sample 45	bc2045--bc20...	97.6	13,848	7,076	Q39	98,001,581	120,954
Bio Sample 46	bc2046--bc20...	97.3	14,362	7,107	Q39	102,074,985	121,586
Bio Sample 47	bc2047--bc20...	97.3	15,416	7,168	Q39	110,509,529	121,808
Bio Sample 48	bc2048--bc20...	96.9	13,819	7,156	Q39	98,899,414	121,848
No Name	Not Barcoded	0.0	36,450	8,306	Q32	302,770,918	119,661

**PacBio Data Management**  
 PHX\_SEQ-12-Cell1 (all samples)  
 Copy Analyze... Export Delete


**Barcode Data**

Sample Name	Barcode	Barcode Qua...	HiFi Reads	HiFi Read Le...	HiFi Read Q...	HiFi Yield (bp)	Polymerase
Bio Sample 1	bc2001--bc20...	97.9	23,838	6,498	Q39	154,912,577	121,916
Bio Sample 2	bc2002--bc20...	96.3	23,829	6,606	Q38	157,423,933	121,269
Bio Sample 14	bc2014--bc20...	97.6	24,733	6,561	Q39	162,281,425	121,368
Bio Sample 15	bc2015--bc20...	97.6	24,921	6,664	Q38	166,086,764	122,085
Bio Sample 16	bc2016--bc20...	97.5	19,922	6,528	Q38	130,058,233	121,249
Bio Sample 17	bc2017--bc20...	97.5	22,234	7,150	Q38	158,979,009	122,790
Bio Sample 18	bc2018--bc20...	97.9	25,923	6,775	Q38	175,644,772	121,687
Bio Sample 19	bc2019--bc20...	97.5	23,040	6,747	Q38	155,467,070	122,291
Bio Sample 20	bc2020--bc20...	97.9	26,138	6,646	Q38	173,730,992	121,650
Bio Sample 21	bc2021--bc20...	97.4	24,954	6,750	Q38	168,449,825	122,122
Bio Sample 22	bc2022--bc20...	97.9	22,402	6,627	Q38	148,478,815	121,233
Bio Sample 23	bc2023--bc20...	97.4	22,389	6,692	Q38	149,847,373	121,316
Bio Sample 24	bc2024--bc20...	97.8	24,244	6,677	Q38	161,888,862	121,594
Bio Sample 48	bc2048--bc20...	96.9	13,819	7,156	Q39	98,899,414	121,848
No Name	Not Barcoded	0.0	36,450	8,306	Q32	302,770,918	119,661

# Recommended guidance for evaluating PureTarget repeat expansion panel sequencing run performance (cont.)

On-target coverage is affected by input gDNA quantity, input gDNA quality and multiplex level<sup>1</sup>


## Input DNA quantity and read coverage<sup>1</sup>



**DNA input quantity versus coverage.** A total mass of 58 µg was prepared from 29 samples and sequenced on a Revo system with SPRQ chemistry.<sup>2</sup>


Sample coverage for target regions increases with quantity of DNA used in PureTarget library prep for a given multiplex level

## Input DNA quality and read coverage<sup>2</sup>



**Importance of DNA quality for target coverage.** PureTarget repeat expansion panel 1.0 libraries were prepared with 2 µg DNA input each and sequenced in a 16-plex on a Revo system with v1 chemistry.<sup>3</sup>

Higher on-target coverage is obtained for samples of higher quality, measured as GQN at 30 kb with the Femto Pulse system



Comparison of a) high- and b) low-quality gDNA samples illustrates that more on-target reads are observed at an *FMR1* repeat expansion for the sample with higher molecular weight DNA

<sup>1</sup> See *Application note — Comprehensive genotyping with the PureTarget repeat expansion panel and HiFi sequencing* ([102-326-614](#)).

<sup>2</sup> DNA was extracted from iPSCs using Nanobind PanDNA kit and QIAGEN Genomic-Tip.

<sup>3</sup> DNA samples were extracted with the MadMAX kit from Applied Biosystems. Mean target coverage is calculated using SMRT Link v13.1 PureTarget repeat expansion analysis application.

# Example coverage performance for PureTarget repeat expansion panel library (96-plex, Revio system + SPRQ chemistry, Nanobind whole blood DNA)

## DNA sample preparation

- 1.2 µg DNA per sample
- DNA samples were extracted from human whole blood with Nanobind kits

## PureTarget library preparation

- Procedure & checklist – Generating PureTarget libraries with PureTarget kit 96 – Automation protocol (103-740-800)
- PureTarget repeat expansion panel 2.0 (38 targets) (103-633-100)
- 96-plex

## Sequencing run design

- Revio system with SPRQ chemistry
- 24 hrs movie time
- No adaptive loading

## Genotyping analysis

- SMRT Link PureTarget repeat expansion analysis application (default settings)

## Coverage (per target per sample)

- Mean = 457 (N = 96)

On-target coverage results for one representative human blood sample from a 96-plex PureTarget repeat expansion panel library.

Gene	Motif	Motif repeat number and PureTarget sequencing coverage <sup>1</sup>			
		Allele 1 motif count	Allele 1 read count	Allele 2 motif count	Allele 2 read count
ABCD3	GCC	7	224	7	228
AFF2	GCC	37	320	41	339
AFF3	GCC	8	290	22	302
AR	GCA	23	457	28	422
ATN1	CAG	17	328	20	299
ATXN1	CTG	27	333	27	286
ATXN10	ATTCT	14	328	13	335
ATXN2	CTG	20	200	21	218
ATXN3	CTG	20	306	23	256
ATXN7	CAG:CCG	10_3	350	10_3	346
ATXN8OS	CTA:CTG	8_10	314	10_18	303
BEAN1	TGGAA:TAGAA:AATAA	0_0_10	276	0_0_17	282
C9orf72	GGCCCC	2	312	2	300
CACNA1A	CTG	13	227	12	236
CNBP	CAGG:CAGA:CA	13_10_15	168	13_10_21	136
CSTB	CGCGGGGCGGGG	3	241	3	241
DAB1	AAAAT:GAAAT	12_0	301	22_0	290
DMPK	CAG	12	360	14	388
FGF14	GAA:GAAGGA:GAAGAAGAAGAAGCA:AAGGAG	17_0_0_0	252	105_1_2_0	205
FMR1	CGG:AGG	27_2	272	49_2	262
FXN	A:GAA	18_9	332	16_16	325
GIPC1	CCG	11	267	10	261
HOXD13	GCN	15	293	15	333
HTT	CAG:CCG	18_8	310	16_11	285
JPH3	CTG	15	279	15	270
LRP12	CGC	8	314	12	285
NOP56	GGCCTG:CGCCTG	6_2	273	8_2	295
NOTCH2NLC	GGC	10	262	15	262
PABPN1	GCN	10	358	10	327
PHOX2B	GCN	20	150	20	110
PPP2R2B	GCT	10	312	10	311
PRNP	CCTCAGGGCGGTGGTGGCTGGGGGCAG:CCTCATGGTGGTGGCTGGGGGCAG:GGTGGTGGCTGGGGGCAGCCTCAT	1_4_0	179	1_4_0	178
RFC1	AAGGG:ACAGG:AGGGC:AAGGC:AGAGG:AAAAG:AAAGG:AAGAG:AAAGGG	0_0_0_0_0_16_0_0_0	270	0_0_0_0_0_92_0_0_0	231
RILPL1	GGC	10	269	10	272
SAMD12	TAAAA:TGAAA	20_0	200	15_0	189
TBP	GCA	32	120	33	109
TCF4	CAG	11	245	23	254
ZFX3	GCC	18	349	18	348

# Example coverage performance for PureTarget repeat expansion panel library (48-plex, Revio system + SPRQ chemistry, Nanobind whole blood DNA)

## DNA sample preparation

- 1.0 µg DNA per sample
- Blood DNA samples were extracted from human whole blood (with or without RBC lysis) using Nanobind kits
- Cell line DNA samples were Coriell lymphoblastoid cell DNA with GQN30kb > 5

## PureTarget library preparation

- Procedure & checklist – Generating PureTarget libraries with PureTarget kit 24 – Manual protocol (103-740-700)
- PureTarget repeat expansion panel 2.0 (38 targets) (103-633-100)
- 48-plex

## Sequencing run design

- Revio system with SPRQ chemistry
- 24 hrs movie time
- No adaptive loading

## Genotyping analysis

- SMRT Link PureTarget repeat expansion analysis application (default settings)

## Coverage (per target per sample)

- Mean = 267 (N = 48)

On-target coverage results for one representative human blood sample from a 48-plex PureTarget repeat expansion panel library.

Gene	Motif	Motif repeat number and PureTarget sequencing coverage <sup>1</sup>			
		Allele 1 motif count	Allele 1 read count	Allele 2 motif count	Allele 2 read count
ABCD3	GCC	7	241	7	240
AFF2	GCC	37	237	37	233
AFF3	GCC	8	196	10	164
AR	GCA	24	210	30	240
ATN1	CAG	17	207	17	207
ATXN1	CTG	25	201	28	182
ATXN10	ATTCT	14	214	16	198
ATXN2	CTG	20	197	21	225
ATXN3	CTG	21	200	24	209
ATXN7	CAG:CCG	10_3	225	10_3	193
ATXN8OS	CTA:CTG	8_11	239	9_16	229
BEAN1	TGGAA:TAGAA:AATAA	0_0_9	210	0_0_10	193
C9orf72	GGCCCC	4	203	5	200
CACNA1A	CTG	11	172	12	166
CNBP	CAGG:CAGA:CA	15_10_18	156	15_9_21	160
CSTB	CGCGGGCGGGG	3	189	2	183
DAB1	AAAAT:GAAAT	15_0	215	13_0	241
DMPK	CAG	11	231	12	211
FGF14	GAA:GAAGGA:GAAGAAGAAGAAGCA:AAGGAG	9_0_0_0	206	23_0_1_0	228
FMR1	CGG:AGG	19_1	198	28_2	212
FXN	A:GAA	16_8	214	18_9	226
GIPC1	CCG	11	223	41	208
HOXD13	GCN	15	242	15	235
HTT	CAG:CCG	18_8	180	17_11	178
JPH3	CTG	15	217	14	190
LRP12	CGC	11	266	12	236
NOP56	GGCCTG:CGCCTG	4_2	217	6_2	250
NOTCH2NLC	GGC	11	220	16	214
PABPN1	GCN	10	203	10	206
PHOX2B	GCN	20	151	20	145
PPP2R2B	GCT	10	209	13	207
PRNP	CCTCAGGGCGGTGGTGGCTGGGGGCAG:CCTCATGGTGGTGGCTGGGGGCAG:GGTGGTGGCTGGGGGCAGCCTCAT	1_4_0	145	1_3_0	171
RFC1	AAGGG:ACAGG:AGGGC:AAGGC:AGAGG:AAAAG:AAAGG:AAGAG:AAAGGG	0_0_0_0_0_11_0_0_0	170	0_0_0_0_0_90_0_0_0	177
RILPL1	GGC	10	236	10	206
SAMD12	TAAAA:TGAAA	17_0	169	21_0	163
TBP	GCA	29	161	29	177
TCF4	CAG	16	179	26	188
ZFX3	GCC	18	220	18	219

# Example coverage performance for PureTarget repeat expansion panel library (48-plex, Vega system, Nanobind whole blood DNA)

## DNA sample preparation

- 2.0 µg DNA per sample
- Blood DNA samples were extracted from human whole blood (with or without RBC lysis) using Nanobind kits
- Cell line DNA samples were Coriell lymphoblastoid cell DNA with GQN30kb > 5

## PureTarget library preparation

- Procedure & checklist – Generating PureTarget libraries with PureTarget kit 24 – Manual protocol (103-740-700)
- PureTarget repeat expansion panel 2.0 (38 targets) (103-633-100)
- 48-plex

## Sequencing run design

- Vega system
- 24 hrs movie time

## Genotyping analysis

- SMRT Link PureTarget repeat expansion analysis application (default settings)

## Coverage (per target per sample)

- Mean = 223 (N = 48)

On-target coverage results for one representative human blood sample from a 48-plex PureTarget repeat expansion panel library.

Gene	Motif	Motif repeat number and PureTarget sequencing coverage <sup>1</sup>			
		Allele 1 motif count	Allele 1 read count	Allele 2 motif count	Allele 2 read count
ABCD3	GCC	7	201	7	200
AFF2	GCC	37	170	37	191
AFF3	GCC	8	167	10	173
AR	GCA	24	168	30	183
ATN1	CAG	17	188	17	187
ATXN1	CTG	25	187	28	191
ATXN10	ATTCT	14	144	16	130
ATXN2	CTG	20	188	21	170
ATXN3	CTG	21	155	24	185
ATXN7	CAG:CCG	10_3	155	10_3	10_3
ATXN8OS	CTA:CTG	8_11	157	9_16	126
BEAN1	TGGAA:TAGAA:AATAA	0_0_9	177	0_0_10	184
C9orf72	GGCCCC	4	149	5	173
CACNA1A	CTG	11	186	12	171
CNBP	CAGG:CAGA:CA	15_10_18	141	15_9_21	119
CSTB	CGCGGGCGGGG	3	158	2	153
DAB1	AAAAT:GAAAT	15_0	151	13_0	156
DMPK	CAG	11	181	12	162
FGF14	GAA:GAAGGA:GAAGAAGAAGAAGCA:AAGGAG	9_0_0_0	160	23_0_1_0	187
FMR1	CGG:AGG	19_1	169	28_2	186
FXN	A:GAA	16_8	195	18_9	191
GIPC1	CCG	11	185	41	191
HOXD13	GCN	15	185	15	172
HTT	CAG:CCG	18_8	140	17_11	142
JPH3	CTG	15	212	14	178
LRP12	CGC	11	198	12	216
NOP56	GGCCTG:CGCCTG	4_2	176	6_2	177
NOTCH2NLC	GGC	11	209	16	180
PABPN1	GCN	10	161	10	165
PHOX2B	GCN	20	127	20	138
PPP2R2B	GCT	10	151	13	168
PRNP	CCTCAGGGCGGTTGGTGGCTGGGGGCAG:CCTCATGGTGGTGGCTGGGGGCAG:GGTGGTGGCTGGGGGCAGCCTCAT	1_4_0	156	1_3_0	156
RFC1	AAGGG:ACAGG:AGGGC:AAGGC:AGAGG:AAAAG:AAAGG:AAGAG:AAAGGG	0_0_0_0_0_11_0_0_0	157	0_0_0_0_0_83_0_0_0	138
RILPL1	GGC	10	242	10	207
SAMD12	TAAAA:TGAAA	17_0	144	21_0	132
TBP	GCA	29	123	29	128
TCF4	CAG	16	161	26	171
ZFX3	GCC	18	185	18	184

# Example coverage performance for PureTarget carrier panel library (96-plex, Revio system + SPRQ chemistry, Non-Nanobind whole blood DNA)

## DNA sample preparation

- 1.2 µg DNA per sample
- DNA samples were extracted from human whole blood with a custom (non-Nanobind) third-party DNA extraction method

## PureTarget library preparation

- *Procedure & checklist – Generating PureTarget libraries with PureTarget kit 96 – Automation protocol (103-740-800)*
- PureTarget carrier panel (12 genes) (103-633-200)
- 96-plex

## Sequencing run design

- Revio system with SPRQ chemistry
- 24 hrs movie time
- No adaptive loading

## Panel QC analysis

- SMRT Link Target Enrichment analysis application [ Include Fail Reads = On; Padding around regions (bp) = 0 ]

## Coverage (per target per sample)

- Mean = 294 (N = 96)

On-target coverage results for one representative human blood sample from a 96-plex PureTarget carrier panel library.

Target	Variant type(s) commonly detected <sup>1,2</sup>	PureTarget sequencing coverage
AFF2	Repeat expansion	627
ARX	Repeat expansion	130
CYP21A2	Small variants and copy number variants	483
CYP21A2P	Small variants and copy number variants	761
F8_intron1_1	Inversions	589
F8_intron1_2	Inversions	704
F8_intron22_1	Inversions	123
F8_intron22_3	Inversions	217
FMR1	Repeat expansion	629
FXN	Repeat expansion	641
GBA	Small variants and copy number variants	282
GBAP	Small variants and copy number variants	322
HBA1/2-largedel	Large deletions	5
HBA1/2-mediumdel	Medium deletions	29
HBA1/2-smalldel	Small deletions	101
HBA1/hs-40	Deletions	458
HBB	Small variants	460
RPGR	Small variants	571
SMN1	Small variants and copy number variants	83
SMN2	Small variants and copy number variants	91

<sup>1</sup> Example variant types listed in the table can be analyzed using the PureTarget Carrier Pipeline analysis workflow available in GitHub and DNAnexus.

<sup>2</sup> **Note:** *TNXB* and *CYP21A2* genes overlap and so sequencing coverage for this region is reported under the *CYP21A2* target in SMRT Link PureTarget carrier panel coverage reports. PureTarget carrier panel includes multiple guide RNA pairs to target *F8* introns 1 and 22 regions to detect inversion variant types (Inv1 and Inv22) resulting from intra-chromosomal translocation events occurring within *F8* intron 1 and intron 22 hotspots. Similarly, multiple guide RNAs are used (in a nested design spanning ~40 kb) to target the *HBA 1/2* region in order to enable detection of small, medium and large deletion events.

# Example coverage performance for PureTarget carrier panel library (48-plex, Revio system + SPRQ chemistry, Nanobind saliva DNA)

## DNA sample preparation

- 1 µg DNA per sample
- DNA samples were extracted from human saliva using Nanobind kits

## PureTarget library preparation

- *Procedure & checklist – Generating PureTarget libraries with PureTarget kit 24 – Manual protocol* (103-329-400)
- PureTarget carrier panel (12 genes) (103-633-200)
- 48-plex

## Sequencing run design

- Revio system with SPRQ chemistry
- 24 hrs movie time
- No adaptive loading

## Panel QC analysis

- SMRT Link Target Enrichment analysis application [ Include Fail Reads = On; Padding around regions (bp) = 0 ]

## Coverage (per target per sample)

- Mean = 250 (N = 48)

On-target coverage results for one representative human blood sample from a 24-plex PureTarget carrier panel library.

Target	Variant type(s) commonly detected <sup>1,2</sup>	PureTarget sequencing coverage
AFF2	Repeat expansion	353
ARX	Repeat expansion	73
CYP21A2	Small variants and copy number variants	673
CYP21A2P	Small variants and copy number variants	483
F8_intron1_1	Inversions	252
F8_intron1_2	Inversions	318
F8_intron22_1	Inversions	36
F8_intron22_3	Inversions	104
FMR1	Repeat expansion	334
FXN	Repeat expansion	590
GBA	Small variants and copy number variants	262
GBAP	Small variants and copy number variants	313
HBA1/2-largedel	Large deletions	5
HBA1/2-mediumdel	Medium deletions	37
HBA1/2-smalldel	Small deletions	230
HBA1/hs-40	Deletions	579
HBB	Small variants	560
RPGR	Small variants	237
SMN1	Small variants and copy number variants	128
SMN2	Small variants and copy number variants	114

<sup>1</sup> Example variant types listed in the table can be analyzed using the PureTarget Carrier Pipeline analysis workflow available in GitHub and DNAnexus.

<sup>2</sup> **Note:** *TNXB* and *CYP21A2* genes overlap and so sequencing coverage for this region is reported under the *CYP21A2* target in SMRT Link PureTarget carrier panel coverage reports. PureTarget carrier panel includes multiple guide RNA pairs to target *F8* introns 1 and 22 regions to detect inversion variant types (Inv1 and Inv22) resulting from intra-chromosomal translocation events occurring within *F8* intron 1 and intron 22 hotspots. Similarly, multiple guide RNAs are used (in a nested design spanning ~40 kb) to target the *HBA 1/2* region in order to enable detection of small, medium and large deletion events.

# Example coverage performance for PureTarget carrier panel library (24-plex, Revio system + SPRQ chemistry, Nanobind saliva DNA)

## DNA sample preparation

- 2 µg DNA per sample
- DNA samples were extracted from human saliva using Nanobind kits

## PureTarget library preparation

- *Procedure & checklist – Generating PureTarget libraries with PureTarget kit 24 – Manual protocol* (103-329-400)
- PureTarget carrier panel (12 genes) (103-633-200)
- 24-plex

## Sequencing run design

- Revio system with SPRQ chemistry
- 24 hrs movie time
- No adaptive loading

## Panel QC analysis

- SMRT Link Target Enrichment analysis application [ Include Fail Reads = On; Padding around regions (bp) = 0 ]

## Coverage (per target per sample)

- Mean = 281 (N = 24)

On-target coverage results for one representative human blood sample from a 24-plex PureTarget carrier panel library.

Target	Variant type(s) commonly detected <sup>1,2</sup>	PureTarget sequencing coverage
AFF2	Repeat expansion	541
ARX	Repeat expansion	129
CYP21A2	Small variants and copy number variants	321
CYP21A2P	Small variants and copy number variants	355
F8_intron1_1	Inversions	388
F8_intron1_2	Inversions	561
F8_intron22_1	Inversions	128
F8_intron22_3	Inversions	163
FMR1	Repeat expansion	502
FXN	Repeat expansion	539
GBA	Small variants and copy number variants	216
GBAP	Small variants and copy number variants	232
HBA1/2-largedel	Large deletions	7
HBA1/2-mediumdel	Medium deletions	40
HBA1/2-smalldel	Small deletions	193
HBA1/hs-40	Deletions	477
HBB	Small variants	449
RPGR	Small variants	431
SMN1	Small variants and copy number variants	122
SMN2	Small variants and copy number variants	61

<sup>1</sup> Example variant types listed in the table can be analyzed using the PureTarget Carrier Pipeline analysis workflow available in GitHub and DNAnexus.

<sup>2</sup> **Note:** *TNXB* and *CYP21A2* genes overlap and so sequencing coverage for this region is reported under the *CYP21A2* target in SMRT Link PureTarget carrier panel coverage reports. PureTarget carrier panel includes multiple guide RNA pairs to target *F8* introns 1 and 22 regions to detect inversion variant types (Inv1 and Inv22) resulting from intra-chromosomal translocation events occurring within *F8* intron 1 and intron 22 hotspots. Similarly, multiple guide RNAs are used (in a nested design spanning ~40 kb) to target the *HBA 1/2* region in order to enable detection of small, medium and large deletion events.

# Example coverage performance for PureTarget carrier panel library (24-plex, Revio system + SPRQ chemistry, Nanobind whole blood DNA)

## DNA sample preparation

- 1 µg DNA per sample
- DNA samples were extracted from human whole blood with RBC lysis using Nanobind kits

## PureTarget library preparation

- *Procedure & checklist – Generating PureTarget libraries with PureTarget kit 24 – Manual protocol* (103-329-400)
- PureTarget carrier panel (12 genes) (103-633-200)
- 24-plex

## Sequencing run design

- Revio system with SPRQ chemistry
- 24 hrs movie time
- No adaptive loading

## Panel QC analysis

- SMRT Link Target Enrichment analysis application [ Include Fail Reads = On; Padding around regions (bp) = 0 ]

## Coverage (per target per sample)

- Mean = 747 (N = 24)

On-target coverage results for one representative human blood sample from a 24-plex PureTarget carrier panel library.

Target	Variant type(s) commonly detected <sup>1,2</sup>	PureTarget sequencing coverage
AFF2	Repeat expansion	1427
ARX	Repeat expansion	447
CYP21A2	Small variants and copy number variants	956
CYP21A2P	Small variants and copy number variants	1900
F8_intron1_1	Inversions	1101
F8_intron1_2	Inversions	1334
F8_intron22_1	Inversions	392
F8_intron22_3	Inversions	602
FMR1	Repeat expansion	1155
FXN	Repeat expansion	1230
GBA	Small variants and copy number variants	405
GBAP	Small variants and copy number variants	539
HBA1/2-largedel	Large deletions	5
HBA1/2-mediumdel	Medium deletions	208
HBA1/2-smalldel	Small deletions	240
HBA1/hs-40	Deletions	1281
HBB	Small variants	958
RPGR	Small variants	942
SMN1	Small variants and copy number variants	363
SMN2	Small variants and copy number variants	353

<sup>1</sup> Example variant types listed in the table can be analyzed using the PureTarget Carrier Pipeline analysis workflow available in GitHub and DNAnexus.

<sup>2</sup> **Note:** *TNXB* and *CYP21A2* genes overlap and so sequencing coverage for this region is reported under the *CYP21A2* target in SMRT Link PureTarget carrier panel coverage reports. PureTarget carrier panel includes multiple guide RNA pairs to target *F8* introns 1 and 22 regions to detect inversion variant types (Inv1 and Inv22) resulting from intra-chromosomal translocation events occurring within *F8* intron 1 and intron 22 hotspots. Similarly, multiple guide RNAs are used (in a nested design spanning ~40 kb) to target the *HBA 1/2* region in order to enable detection of small, medium and large deletion events.

# Example coverage performance for PureTarget carrier panel library (24-plex, Vega system , Nanobind whole blood DNA)

## DNA sample preparation

- 2 µg DNA per sample
- DNA samples were extracted from human whole blood with RBC lysis using Nanobind kits

## PureTarget library preparation

- *Procedure & checklist – Generating PureTarget libraries with PureTarget kit 24 – Manual protocol* (103-329-400)
- PureTarget carrier panel (12 genes) (103-633-200)
- 24-plex

## Sequencing run design

- Vega system
- 24 hrs movie time

## Panel QC Analysis

- SMRT Link Target Enrichment analysis application [ Include Fail Reads = On; Padding around regions (bp) = 0 ]

## Coverage (per target per sample)

- Mean = 592 (N = 24)

On-target coverage results for one representative human blood sample from a 24-plex PureTarget carrier panel library.

Target	Variant type(s) commonly detected <sup>1,2</sup>	PureTarget sequencing coverage
AFF2	Repeat expansion	1064
ARX	Repeat expansion	329
CYP21A2	Small variants and copy number variants	919
CYP21A2P	Small variants and copy number variants	1662
F8_intron1_1	Inversions	860
F8_intron1_2	Inversions	977
F8_intron22_1	Inversions	288
F8_intron22_3	Inversions	426
FMR1	Repeat expansion	847
FXN	Repeat expansion	863
GBA	Small variants and copy number variants	365
GBAP	Small variants and copy number variants	496
HBA1/2-largedel	Large deletions	8
HBA1/2-mediumdel	Medium deletions	169
HBA1/2-smalldel	Small deletions	218
HBA1/hs-40	Deletions	925
HBB	Small variants	765
RPGR	Small variants	856
SMN1	Small variants and copy number variants	258
SMN2	Small variants and copy number variants	244

<sup>1</sup> Example variant types listed in the table can be analyzed using the PureTarget Carrier Pipeline analysis workflow available in GitHub and DNAnexus.

<sup>2</sup> **Note:** *TNXB* and *CYP21A2* genes overlap and so sequencing coverage for this region is reported under the *CYP21A2* target in SMRT Link PureTarget carrier panel coverage reports. PureTarget carrier panel includes multiple guide RNA pairs to target *F8* introns 1 and 22 regions to detect inversion variant types (Inv1 and Inv22) resulting from intra-chromosomal translocation events occurring within *F8* intron 1 and intron 22 hotspots. Similarly, multiple guide RNAs are used (in a nested design spanning ~40 kb) to target the *HBA 1/2* region in order to enable detection of small, medium and large deletion events.



# Technical documentation & applications support resources


# Technical resources for PureTarget library preparation, sequencing & data analysis

## DNA sample preparation literature & other resources

- Brochure – Nanobind PanDNA kit ([102-326-604](#))
- Procedure & checklist – Extracting HMW DNA from human whole blood using Nanobind kits ([102-573-500](#))
- Procedure & checklist – Extracting HMW DNA from human whole blood with RBC lysis using Nanobind kits ([103-377-500](#))
- Procedure & checklist – Extracting HMW DNA from cultured suspension cells using Nanobind kits ([103-394-500](#))
- Procedure & checklist – Extracting HMW DNA from cultured adherent cells using Nanobind kits ([102-573-600](#))

## PureTarget library preparation literature & other resources

- Application note – Comprehensive genotyping with PureTarget repeat expansion panel and HiFi sequencing ([102-326-614](#))
- Application note – Consolidate challenging genes with PureTarget carrier screen panel ([102-326-653](#))
- Brochure – Comprehensive genotyping with PureTarget repeat expansion panel ([102-326-609](#))
- Guide & overview – Automated PureTarget kit 96 for the Hamilton NGS STAR MOA system ([103-740-900](#))
- Procedure & checklist – Generating PureTarget libraries with PureTarget kit 24 – Manual protocol ([103-740-700](#))
- Procedure & checklist – Generating PureTarget libraries with PureTarget kit 96 – Automation protocol ([103-740-800](#))
- Technical note – A practical guide to amplification-free PureTarget custom panels ([102-326-652](#))
- Technical overview – PureTarget library preparation using PureTarget kit 2.0 ([103-742-600](#))



# Technical resources for PureTarget library preparation, sequencing & data analysis (cont.)

## Data analysis resources

- Application note – Comprehensive genotyping with PureTarget repeat expansion panel and HiFi sequencing ([102-326-614](#))
- Application note – Consolidate challenging genes with PureTarget carrier screen panel ([102-326-653](#))
- PacBio GitHub Wiki – PureTarget Carrier Pipeline (PTCP) analysis workflow [ [Link](#) ]
- SMRT Link software installation guide [ [Link](#) ]
- SMRT Link user guide [ [Link](#) ]
- SMRT Tools reference guide [ [Link](#) ]

## Publications and posters

- PacBio ESHG poster (2025) - Targeted long-read sequencing of native DNA for genetic disease diagnostic and screening research [ [Link](#) ]
- PacBio ACMG poster (2025) – Targeted long-read sequencing of native DNA for genetic disease diagnostic and screening research [ [Link](#) ]
- PacBio AMP poster (2024) – Targeted long-read sequencing of native DNA for comprehensive characterization of repeat expansions [ [Link](#) ]

## Videos

- PacBio YouTube video (2025) – PureTarget in action: Scalable and accurate repeat expansion sequencing [ [Link](#) ]
- PacBio YouTube video (2025) – PureTarget: PCR-Free Enrichment with Methylation for Challenging Genes [ [Link](#) ]
- PacBio YouTube video (2025) – TRGT: profiling variation in and around tandem repeats [ [Link](#) ]



# Technical resources for PureTarget library preparation, sequencing & data analysis (cont.)

## Example PacBio data sets<sup>1</sup>

Application	Dataset	Data type	PacBio system
PureTarget 2.0	Carrier panel Nanobind 96plex [ <a href="#">Link</a> ]	HiFi long read	Revio + SPRQ
	Carrier panel Coriell 16plex [ <a href="#">Link</a> ]	HiFi long read	Revio + SPRQ
	Carrier panel Coriell 16plex [ <a href="#">Link</a> ]	HiFi long read	Vega system
	Repeat Expansion 2.0 Nanobind Coriell 48plex [ <a href="#">Link</a> ]	HiFi long read	Revio + SPRQ
	Repeat Expansion 2.0 Nanobind Coriell 48plex [ <a href="#">Link</a> ]	HiFi long read	Vega system



[www.pacb.com](http://www.pacb.com)

Revision history (description)	Version	Date
Initial release	01	September 2025

Research use only. Not for use in diagnostic procedures. © 2025 Pacific Biosciences of California, Inc. (“PacBio”). All rights reserved. Information in this document is subject to change without notice. PacBio assumes no responsibility for any errors or omissions in this document. Certain notices, terms, conditions and/or use restrictions may pertain to your use of PacBio products and/or third-party products. Refer to the applicable PacBio terms and conditions of sale and to the applicable license terms at [pacb.com/license](http://pacb.com/license). Pacific Biosciences, the PacBio logo, PacBio, Circulomics, Omniome, SMRT, SMRTbell, Iso-Seq, Sequel, Nanobind, SBB, Revio, Onso, Apton, Kinnex, PureTarget, SPRQ, and Vega are trademarks of PacBio.