




Abstract

Microbial genome sequencing can be done quickly, easily, and efficiently with PacBio sequencing instruments, resulting in complete *de novo* assemblies. Alternative protocols have been developed to reduce the amount of purified DNA required for SMRT Sequencing to broaden applicability to lower-abundance samples. If 50-100 ng of microbial DNA is available, a 10-20 kb SMRTbell library can be made. The resulting library can be loaded onto one or more SMRT Cells, yielding more than enough data for complete assembly of microbial genomes using the SMRT Portal assembly pipeline HGAP, plus epigenetic analysis. The entire process can be done in 2 to 3 days by standard laboratory personnel.

This approach is particularly important for analysis of microbial communities, in which genomic DNA is often limited. For simple microbiomes, it may be possible to obtain complete or near complete genome sequences of abundant members with low-input 10 kb libraries. For more complex communities, or when only a few ng of input DNA are available, a 2 kb sheared library can be prepared to generate 1-2 kb reads with >99.9% accuracy with Circular Consensus Sequencing. These very accurate, single-molecule reads can be used to provide information about the microbial composition of a community, as well as biochemical pathways present in the sample.


SMRTbell Library Prep Workflow



Highly Accurate Single-Molecule Sequencing

MULTIPLE READS FROM A SINGLE MOLECULE

As a function of the SMRTbell adapters, multiple single-pass reads are generated from an individual molecule. Combining these subreads corrects for random errors and results in a highly accurate single-molecule consensus sequence. Data can be filtered to an accuracy of 99.99% with CCS2 analysis.



Library Prep Options for Low-Input Sequencing


Library Size	Protocol	Input	SMRT Cells	Total Bases
10 kb	Standard	1 µg (minimum)	≥150	≥90 Gb
10 kb	Low Input	100 ng	≥4	≥2.4 Gb
2 kb	Standard	500 ng (minimum)	≥400	≥240 Gb
2 kb	Very Low Input	10 ng	≥4	≥2.4 Gb

Results for the PacBio RS II, with conservative estimates of library and SMRT Cell yield.

2 kb SMRTbell Libraries from 10 ng Input DNA


MICROBIOME PROFILING

DNA was purified from an environmental (lake) sample and prepared for sequencing using the 2 kb, very low input SMRTbell Library Prep Protocol. Data was used to determine genes, which were aligned to a reference database.



Lonza FlashGel and Bioanalyzer QC of sheared microbiome DNA prior to SMRTbell library prep


B Primary analysis of PacBio data, 2 kb sheared lake microbiome samples




C Circular consensus Reads of Insert yield, SMRT Analysis

Sample	Primary		90% Accuracy		99% Accuracy	
	P1 Reads	Total Bases	# of Reads	Total Bases	# of Reads	
Lake Microbiome	90 K	74 Mb	64 K	56 Mb	48 K	
Mock Community	114 K	90 Mb	82 K	66 Mb	60 K	

D Taxonomic groups found in the lake microbiome sample



Bioinformatic Workflow for Profiling Microbial Communities




MULTIPLEXING EXTREMELY LOW INPUT SAMPLES

For sub-ng inputs, a small number of highly accurate 1-2 kb sequences may be obtained by pooling libraries with barcoded adapters. Here is an example of where an undetectable amount of DNA was pooled with a sample with ~4 ng input following ligation of barcoded adapters

Sample	Input DNA (ng)	Bar-code	# Cells	Total Bases	P1 Reads	CCS Reads (≥99% accuracy)
Sample 1	3.9	34	3	2.8 Gb	195,876	104,631
Sample 2	Un-detectable*	51				71

*With Qubit dsDNA High Sensitivity Assay



Profiling results

- Sample 1 and Sample 2 communities were quite different
- Sample 2 reads included perfect to near-perfect kb-long alignments

Primary analysis results for one PacBio RS II SMRT Cell; all 3 SMRT Cells had similar results.

2 - 20 kb Low-Input SMRTbell Library Protocols

Procedure & Checklist - 10 kb to 20 kb Template Preparation and Sequencing with Low (100 ng) Input DNA

Before You Begin
To perform this procedure, you must have the PacBio®:

- Template Prep Kit
- Akit® P1 Beads
- DNA Polymerase Binding Kit (P1 v2 or later)
- MagBead Kit
- DNA Sequencing Kit
- SMRT® Cells

This procedure can be used to prepare 10-20 kb libraries from 100 ng to 200 ng of sheared and concentrated DNA, minimally 100 ng into shearing. Note: for input amounts between 200 ng and 1 µg, the standard 10 kb library prep protocol may be used.

Insert Size Target	Insert Size Range	Input DNA Amount	Ligation	DNA Damage Repair
10 kb to 20 kb	8 kb to 22 kb	100 to 200 ng	Blank	Required

For the full protocols, visit www.pacb.com/support/documentation

Procedure & Checklist - Very Low (10 ng) Input 2 kb Template Preparation and Sequencing with Carrier DNA

Before You Begin
To perform this procedure, you must have the PacBio®:

- Template Prep Kit
- Akit® P1 Beads
- DNA Polymerase Binding Kit (P1 v2 or later)
- MagBead Kit
- DNA Sequencing Kit
- SMRT® Cells

Additional required materials:

- Carrier plasmid DNA, 10 µg (pBR322 or pUC, NEB or Thermo Scientific)
- Exo III, Exo VI, and Template prep buffer from the DNA Template Prep Kit may be used to prepare carrier DNA (Bioss)


This procedure can be used to prepare 1 kb to 3 kb libraries from 5 ng to 50 ng of sheared and concentrated DNA, or from 100 to 200 ng of input DNA. Note that when preparing libraries with a low DNA amount you must use MagBead loading for sequencing.

10 kb SMRTbell Libraries from 100 ng Input DNA


10-20 KB LOW-INPUT LIBRARY FROM A SIMPLE COMMUNITY: BIOCATHODE DNA

DNA was isolated from biofilm collected on an electrode. SMRTbell libraries were prepared using the 10-20 kb Low-Input Shared Protocol (125 ng input) or the standard 10 kb Template Preparation and Sequencing protocol (1.25 µg input).

A. QC of sheared biofilm DNA prior to SMRTbell library prep



B. Primary analysis of PacBio data, 10 kb sheared biofilm samples




C. PacBio RS II sequencing yield

Sample	Number of SMRT Cells ¹	Number of P1 ² Reads	Total Bases	Avg Insert Length	Avg Polymerase Read Length
Low (125 ng) Input	2	97 K	1.3 Gb	4.6 kb	13.7 kb
Std (1.25 µg) Input	7	579 K	9.2 Gb	5.2 kb	15.9 kb

¹Did not sequence entire library

²P1 reads contain usable sequence information

D. Genome assembly of the two most prevalent species



Results

- Microbial genomes have been assembled from low input amounts (100 ng) of DNA using the 10 kb – 20 kb SMRTbell library prep protocol, sequenced on the PacBio RS II.
- PacBio community profiling can be done with very low amounts (10 ng or less) of DNA from microbiome samples prepared with the 2 kb, very low input SMRTbell library protocol. Additionally, samples with sub-ng amounts of DNA may yield a small number of very accurate reads when multiplexed with samples with a few ng.

Acknowledgements

The authors would like to thank:

- Dr. Sarah Strycharz-Glaven, Research Biologist at the Naval Research Laboratory, for the biocathode samples.
- Dr. Tanja Woyke, the Microbial Genomics Program Lead at the DOE Joint Genome Institute, for the lake, plant and mock metagenomic samples.