

---

**Fixed Issues**    **RS System Software**

1. Fixed a failure involving SMRT Cells with more than 2 Gigabases of data. (26131)
2. Fixed out-of-memory issues related to Primary Analysis. (26091, 26332, 26628)
3. Fixed timeout issues with excess control reads. (26062)
4. Increased the maximum movie time to 6 hours. (26726)

**RS Remote**

1. Users can now resubmit failed primary jobs from RS Touch. (23474)

For Research Use Only. Not for use in diagnostic procedures. © Copyright 2015, Pacific Biosciences of California, Inc. All rights reserved. Information in this document is subject to change without notice. Pacific Biosciences assumes no responsibility for any errors or omissions in this document. Certain notices, terms, conditions and/or use restrictions may pertain to your use of Pacific Biosciences products and/or third party products. Please refer to the applicable Pacific Biosciences Terms and Conditions of Sale and the applicable license terms at <http://www.pacificbiosciences.com/licenses.html>.

Pacific Biosciences, the Pacific Biosciences logo, PacBio, SMRT, SMRTbell and Iso-Seq are trademarks of Pacific Biosciences. BluePippin and SageELF are trademarks of Sage Science, Inc. NGS-go and NGSengine are trademarks of GenDx. All other trademarks are the sole property of their respective owners.