

Novel Insights into Microbial and Viral Complexity




Fully characterize complete microbial genomes and populations affordably with SMRT® Sequencing, simultaneously detailing base modifications and mobile elements, quantifying low-level variants and achieving strain-level resolution within communities.

Complete Microbial Characterization in a Single Experiment

“We can now do it [investigate plasmid transfer] in the context of fully sequenced genomes because we now know every plasmid and what it carries. We’re now doing a study that I think is fully informed.”


- Julie Segre, NHGRI¹

- Generate gold-standard reference genomes
- Assess contribution of base modification on pathogenesis
- Understand the role of mobile elements in drug resistance and transmission



Long-read sequencing can resolve complete plasmids in order to track the dissemination of antimicrobial-resistance genes among pathogenic bacteria²

Investigate Viral Population Dynamics



- De-convolute complex mixtures of unique haplotypes
- Track evolution and phylogeny of viral populations
- Identify and quantify minor variants

“To understand within-host viral evolution, long high-quality sequences are absolutely necessary. PacBio is the only NGS technology that hits the sweet spot for length, accuracy and depth.”

- Ben Murrell, UCSD CFAR³

HIV *env* sequences expand and diversify throughout infection³


Detailed Examination of Microbial Communities

“SMRT Sequencing is a single molecule technique that can generate long reads (10-15 kb), is highly accurate and can distinguish methylated bases from the normal A,C,G,T. This latter property is unique as no other method can do that for N⁶-methyladenine or N⁴-methylcytosine without additional chemistry being involved.”

- Rich Roberts, New England Biolabs⁴

- Discover novel genes with longer reads and contigs
- Identify and assign functions to community members with better genome assemblies
- Leverage epigenomic data to cluster contigs, including those from closely related strains

Full-length 16S



Species-level resolution of soil sample⁵

Whole-Genome Shotgun



De novo assembly string graph of Human Microbiome Mock Community B (HM-276D), from the Human Microbiome Project

Key References

1. GenomeWeb. (2014) PacBio Long Reads Assess Antibiotic Resistance Plasmids in Wake of Enterobacteriaceae Outbreak. Retrieved February 3, 2015 from <https://www.genomeweb.com/sequencing/pacbio-long-reads-assess-antibiotic-resistance-plasmids-wake-enterobacteriaceae>
2. Beatson and Walker. (2014) [Tracking Antibiotic Resistance](#). *Science*. **345** (6203):1454-1455.
3. Murrell, et al. (2015) Full-length *env* Deep Sequencing in a Donor with Broadly Neutralizing V1/V2 Antibodies. Conference on Retroviruses and Opportunistic Infections (2015 CROI) Poster Presentation.
4. Naomi Attar. (2013) Rich Roberts discusses single-molecule sequencing technology. Retrieved February 3, 2015 from <http://www.biomedcentral.com/biome/rich-roberts-discusses-single-molecule-sequencing-technology/>
5. Bowman, et al. (2013) Analysis of Full-Length Metagenomic 16S Genes by SMRT[®] Sequencing. American Society for Microbiology (2013 ASM) Poster Presentation. Retrieved February 5, 2015 from http://files.pacb.com/pdf/Analysis_of_Full_Length_Metagenomic_16S_Genes_by_SMRT_Sequencing.pdf



For Research Use Only. Not for use in diagnostic procedures. © Copyright 2015, Pacific Biosciences of California, Inc. All rights reserved. Information in this document is subject to change without notice. Pacific Biosciences assumes no responsibility for any errors or omissions in this document. Certain notices, terms, conditions and/or use restrictions may pertain to your use of Pacific Biosciences products and/or third party products. Please refer to the applicable Pacific Biosciences Terms and Conditions of Sale and to the applicable license terms at <http://www.pacificbiosciences.com/licenses.html>.

Pacific Biosciences, the Pacific Biosciences logo, PacBio, SMRT, SMRTbell and Iso-Seq are trademarks of Pacific Biosciences in the United States and/or certain other countries. All other trademarks are the sole property of their respective owners.

PN: VM102-020915