

**Description**

The DNA Internal Control Complex P6 contains a fixed template of 2 kb length, bound to DNA Polymerase P6 for use as an internal sequencing control.

Reagent quantities provide spike-in controls for a minimum of 24 samples.

**Part Number**

100-356-500

**Quantities**

Minimum of 24 reactions

**Intended Use**

For Research Use Only. Not for use in diagnostic procedures.

Certain notices, terms, conditions and/or use restrictions may pertain to your use of Pacific Biosciences® products and/or third party products. Please refer to the applicable Pacific Biosciences Terms and Conditions of Sale and to the applicable license terms at <http://www.pacificbiosciences.com/licenses.html>.

**Storage and Handling Requirements**

Product should be stored at -20°C.  
Refer to Sample Preparation and Sequencing Guide for additional recommendations for reagent handling.

**Shelf Life**

Refer to expiration date on kit label.

**Support**

For additional product and technical information, such as Safety Data Sheets (SDS) or procedures, please visit our website at [www.pacificbiosciences.com](http://www.pacificbiosciences.com) and on our customer portal at [www.pacbioportal.com](http://www.pacbioportal.com).

**Legal**

© Copyright 2010-2014, Pacific Biosciences of California, Inc. All rights reserved.

Information in this document is subject to change without notice. Pacific Biosciences assumes no responsibility for any errors or omissions in this document.

PACIFIC BIOSCIENCES DISCLAIMS ALL WARRANTIES WITH RESPECT TO THIS DOCUMENT, EXPRESS, STATUTORY, IMPLIED OR OTHERWISE, INCLUDING, BUT NOT LIMITED TO, ANY WARRANTIES OF MERCHANTABILITY, SATISFACTORY QUALITY, NONINFRINGEMENT OR FITNESS FOR A PARTICULAR PURPOSE. IN NO EVENT SHALL PACIFIC BIOSCIENCES BE LIABLE, WHETHER IN CONTRACT, TORT, WARRANTY, PURSUANT TO ANY STATUTE, OR ON ANY OTHER BASIS FOR SPECIAL, CONSEQUENTIAL, INCIDENTAL, EXEMPLARY OR INDIRECT DAMAGES IN CONNECTION WITH (OR ARISING FROM) THIS DOCUMENT, WHETHER OR NOT FORESEEABLE AND WHETHER OR NOT PACIFIC BIOSCIENCES IS ADVISED OF THE POSSIBILITY OF SUCH DAMAGES.

Trademarks:

Pacific Biosciences, the Pacific Biosciences logo, PacBio, SMRT, SMRTbell, and Iso-Seq are trademarks of Pacific Biosciences in the United States and/or certain other countries.

All other trademarks are the sole property of their respective owners.