

SMRT Link Cloud

Skip the setup. Start sequencing.

Stop waiting weeks for IT setup and server approvals.

With **SMRT® Link Cloud**, you can get your Revio® or Vega™ system running in hours.






SMRT Link Cloud is a **free, PacBio®-hosted version of SMRT Link** that removes the need for IT maintenance. Your system stays secure, current, and ready to run, fully managed by PacBio.

Access your familiar SMRT Link dashboard from any browser or any operating system to set up and monitor runs. PacBio handles updates and security so you can focus on results, not infrastructure.

Why SMRT Link Cloud?

- **Start fast.** Software is setup in hours.
- **Save money.** Skip \$1K–\$10K in hardware costs.
- **Stay current.** Always have access to the latest features.
- **Run reliably.** PacBio-hosted infrastructure ensures uptime and performance.
- **Keep control.** Your data stays yours. Store it locally or in your own cloud.
- **Stay flexible.** Use as your main, backup, or test environment.

How it works

Function	Who manages it?	What you get
 Hosting and maintenance	PacBio	No setup or upkeep. PacBio keeps everything running smoothly.
 Run setup and monitoring	You	Control your runs from any browser with live performance data.
 Instrument connection (Revio/Vega)	You	Connect directly to start runs instantly.
 Data storage	You	Keep full ownership of your data. Store it locally or in your own cloud (S3, Google Cloud, or Azure).
 Analysis	You	Run analysis locally or on your compute system. 'Bring Your Own Compute' support planned for 2026.
 Cost and maintenance	PacBio	Free for all PacBio customers. Setup takes about one business day after site prep.

Built to help every lab get sequencing started faster without IT or delays.



Ready to start sequencing?

Contact your FABS or [PacBio Tech Support](#) to activate your free SMRT Link Cloud workspace and begin sequencing in hours.

Research use only. Not for use in diagnostic procedures. © 2025 Pacific Biosciences of California, Inc. ("PacBio"). All rights reserved. Information in this document is subject to change without notice. PacBio assumes no responsibility for any errors or omissions in this document. Certain notices, terms, conditions and/or use restrictions may pertain to your use of PacBio products and/or third-party products. Refer to the applicable PacBio terms and conditions of sale and to the applicable license terms at [pacb.com/license](#). Pacific Biosciences, the PacBio logo, PacBio, Circulomics, Omniome, SMRT, SMRTbell, Iso-Seq, Sequel, Nanobind, SBB, Revio, Onso, Apton, Kinnex, PureTarget, SPRQ, and Vega are trademarks of PacBio.