

Description

Barcoded Universal F/R Primers Plate – 96 v2 is a 96-well plate containing barcoded primers. Each well on the plate contains a primer pair consisting of forward and reverse primers having the same barcode. Each barcode is present in only one well for a total of 96 different barcodes on the plate. Each primer has a universal sequence on the 3' end to allow for incorporation of barcodes into templates through a two-step PCR strategy. Reagent quantities support a single use of each of the 96 primer pairs in library preparations. The barcode sequences have been selected for optimized demultiplexing on both the Sequel® System and PacBio® RS II System.

Part Number

101-629-100

Quantity

1 reaction for each barcoded primer pair

Intended Use

For Research Use Only. Not for use in diagnostic procedures.

Certain notices, terms, conditions and/or use restrictions may pertain to your use of Pacific Biosciences® products and/or third-party products. Please refer to the applicable Pacific Biosciences Terms and Conditions of Sale and to the applicable license terms at

<https://www.pacb.com/legal-and-trademarks/terms-and-conditions-of-sale/>

Storage and Handling Requirements

Product should be stored at -20°C.

Refer to Sample Preparation and Sequencing Guide for additional recommendations for reagent handling.

Shelf Life

Refer to expiration date on kit label.

Support

For additional product and technical information, such as Safety Data Sheets (SDS) or procedures, please visit our website at www.pacificbiosciences.com and our customer portal at www.pacbioportal.com.

Legal

© Copyright 2010-2018, Pacific Biosciences of California, Inc. All rights reserved.

Information in this document is subject to change without notice. Pacific Biosciences assumes no responsibility for any errors or omissions in this document.

PACIFIC BIOSCIENCES DISCLAIMS ALL WARRANTIES WITH RESPECT TO THIS DOCUMENT, EXPRESS, STATUTORY, IMPLIED OR OTHERWISE, INCLUDING, BUT NOT LIMITED TO, ANY WARRANTIES OF MERCHANTABILITY, SATISFACTORY QUALITY, NONINFRINGEMENT OR FITNESS FOR A PARTICULAR PURPOSE. IN NO EVENT SHALL PACIFIC BIOSCIENCES BE LIABLE, WHETHER IN CONTRACT, TORT, WARRANTY, PURSUANT TO ANY STATUTE, OR ON ANY OTHER BASIS FOR SPECIAL, CONSEQUENTIAL, INCIDENTAL, EXEMPLARY OR INDIRECT DAMAGES IN CONNECTION WITH (OR ARISING FROM) THIS DOCUMENT, WHETHER OR NOT FORESEEABLE AND WHETHER OR NOT PACIFIC BIOSCIENCES IS ADVISED OF THE POSSIBILITY OF SUCH DAMAGES.

Trademarks:

Pacific Biosciences, the Pacific Biosciences logo, PacBio, SMRT, SMRTbell, Iso-Seq, and Sequel are trademarks of Pacific Biosciences in the United States and/or certain other countries.

All other trademarks are the sole property of their respective owners.