# Iso-Seq® express oligo kit

## **Description**

The Iso-Seq express oligo kit contains template switching oligo (TSO) and cDNA PCR primer for Iso-Seq library preparation.

## Kit configuration & layout

## Configuration

#	Component	Part number	Qty	Color	Volume
1	Iso-Seq express template switching oligo	101-738-400	1	blue	40 μL
2	Iso-Seq express cDNA PCR primer	101-738-500	1	green	80 μL

#### Layout









24 reactions

## Storage and handling requirements

Product should be stored frozen between -70°C to -80°C.

Refer to protocols for additional recommendations for reagent handling.

#### Shelf life

Refer to expiration date on kit label.

#### Intended use

Research use only. Not for use in diagnostic procedures.

## **Support**

For additional product and technical information, such as safety data sheets (SDS) or procedures, visit our website.

## Legal

Information in this document is subject to change without notice. PacBio assumes no responsibility for any errors or omissions in this document. Certain notices, terms, conditions and/or use restrictions may pertain to your use of PacBio products and/or third-party products. Refer to the applicable PacBio terms and conditions of sale and to the applicable license terms at pacb.com/license. Pacific Biosciences, the PacBio logo, PacBio, Circulomics, Omniome, SMRT, SMRTbell, Iso-Seq, Sequel, Nanobind, and SBB are trademarks of PacBio. All other trademarks are the sole property of their respective owners.

