

Description

The Iso-Seq™ Express Oligo Kit contains two oligonucleotides to support the Iso-Seq Express workflow. The oligonucleotides are the PacBio Template Switching Oligo (TSO) and PacBio cDNA PCR primer. The oligos are designed to replace the equivalent oligos in the NEBNext® Single Cell/Low Input cDNA Synthesis & Amplification Module. The reagents are sufficient for approximately 24 reactions.

Kit components include:

Subcomponent Name	Part Number	Number	Color	Quantity (µl)
Iso-Seq Express Template Switching Oligo	101-738-400	1 tube	Blue	40 µl
Iso-Seq Express cDNA PCR Primer	101-738-500	1 tube	Green	80 µl

Part Number

101-737-500

Quantity

24 reactions

Intended Use

For Research Use Only. Not for use in diagnostic procedures.

Certain notices, terms, conditions and/or use restrictions may pertain to your use of Pacific Biosciences® products and/or third-party products. Please refer to the applicable Pacific Biosciences Terms and Conditions of Sale and to the applicable license terms at <https://www.pacb.com/legal-and-trademarks/terms-and-conditions-of-sale/>

Storage and Handling Requirements

Product should be stored frozen between -70°C to -80°C.

Refer to Sample Preparation and Sequencing Guide for additional recommendations for reagent handling.

Shelf Life

Refer to expiration date on kit label.

Support

For additional product and technical information, such as Safety Data Sheets (SDS) or procedures, please visit our website at <http://www.pacb.com/support/documentation/>

Legal

© Copyright 2019, Pacific Biosciences of California, Inc. All rights reserved.

Information in this document is subject to change without notice. Pacific Biosciences assumes no responsibility for any errors or omissions in this document.

PACIFIC BIOSCIENCES DISCLAIMS ALL WARRANTIES WITH RESPECT TO THIS DOCUMENT, EXPRESS, STATUTORY, IMPLIED OR OTHERWISE, INCLUDING, BUT NOT LIMITED TO, ANY WARRANTIES OF MERCHANTABILITY, SATISFACTORY QUALITY, NONINFRINGEMENT OR FITNESS FOR A PARTICULAR PURPOSE. IN NO EVENT SHALL PACIFIC BIOSCIENCES BE LIABLE, WHETHER IN CONTRACT, TORT, WARRANTY, PURSUANT TO ANY STATUTE, OR ON ANY OTHER BASIS FOR SPECIAL, CONSEQUENTIAL, INCIDENTAL, EXEMPLARY OR INDIRECT DAMAGES IN CONNECTION WITH (OR ARISING FROM) THIS DOCUMENT, WHETHER OR NOT FORESEEABLE AND WHETHER OR NOT PACIFIC BIOSCIENCES IS ADVISED OF THE POSSIBILITY OF SUCH DAMAGES.

Trademarks:

Pacific Biosciences, the Pacific Biosciences logo, PacBio, SMRT, SMRTbell, Iso-Seq, and Sequel are trademarks of Pacific Biosciences in the United States and/or certain other countries.

All other trademarks are the sole property of their respective owners.