

Description

The Kinnex capture beads kit contains beads and buffers to capture the 3' or 5' of 10x single-cell cDNA molecules.

Kit configuration & layout

Configuration

#	Component	Part number	Qty	Color	Volume
1	Kinnex capture beads	103-144-900	1	Clear	160 µL
2	Kinnex bead binding buffer	103-145-400	1	Clear	1250 µL
3-4	Kinnex bead washing buffer	103-145-500	2	Clear	1100 µL

Layout



Quantities

12 reactions

Refer to protocols for additional recommendations for reagent handling.

Shelf life

Refer to expiration date on kit label.

Intended use

Research use only. Not for use in diagnostic procedures.

Support

For additional product and technical information, such as safety data sheets (SDS) or procedures, visit our [website](#).

Legal

Information in this document is subject to change without notice. PacBio assumes no responsibility for any errors or omissions in this document. Certain notices, terms, conditions and/or use restrictions may pertain to your use of PacBio products and/or third-party products. Refer to the applicable PacBio terms and conditions of sale and to the applicable license terms at pacb.com/license. Pacific Biosciences, the PacBio logo, PacBio, Circulomics, Omniome, SMRT, SMRTbell, Iso-Seq, Sequel, Nanobind, SBB, Revio, and Onso are trademarks of PacBio.

Storage and handling requirements

Product should be stored between 2 °C to 8 °C.