# **Description**

The Kinnex concatenation kit contains reagents and oligos to be used with other Kinnex PCR kits to create concatenated barcoded Kinnex libraries.

## Kit configuration & layout

#### Configuration

#	Component	Part number	Qty	Color	Volume
1	Kinnex ligase	103-110-500	1	yellow	80 µL
2	Kinnex enzyme	103-110-400	1	red	53 µL
3	Kinnex array and repair buffer	103-110-300	1	white	172 µL
4	Kinnex adapter bc01 mix	103-109-600	1	blue	10 μL
5	Kinnex adapter bc02 mix	103-109-700	1	blue	10 μL
6	Kinnex adapter bc03 mix	103-109-800	1	blue	10 μL
7	Kinnex adapter bc04 mix	103-109-900	1	blue	10 μL
8	DNA repair mix	103-110-000	1	green	27 μL
9	Nuclease buffer	103-110-200	1	light purple	66 µL
10	Nuclease mix	103-110-100	1	light green	66 µL

#### Layout



# **Quantities**

12 reactions

## Storage and handling requirements

Product should be stored frozen between -15 to -25 °C.

Refer to protocols for additional recommendations for reagent handling.

#### Shelf life

Refer to expiration date on kit label.

#### Intended use

Research use only. Not for use in diagnostic procedures.

## **Support**

For additional product and technical information, such as safety data sheets (SDS) or procedures, visit our website.

## Legal

Research use only. Not for use in diagnostic procedures. © 2023 Pacific Biosciences of California, Inc. ("PacBio"). Information in this document is subject to change without notice. PacBio assumes no responsibility for any errors or omissions in this document. Certain notices, terms, conditions and/or use restrictions may pertain to your use of PacBio products and/or third-party products. Refer to the applicable PacBio terms and conditions of sale and to the applicable license terms at <a href="macb.com/license">pacb.com/license</a>. Pacific Biosciences, the PacBio logo, PacBio, Circulomics, Omniome, SMRT, SMRTbell, Iso-Seq, Sequel, Nanobind, SBB, Revio, Onso, Apton, and Kinnex are trademarks of PacBio.



Kit insert PN: 103-293-300 Rev: 02