

## Description

The LongPlex multiplexing kit box 1 contains reagents to be used in conjunction with the LongPlex multiplexing kit box 2 (301304) to simultaneously fragment and index DNA to generate 8-10 kb fragments.

## Kit configuration

### Configuration

#	Component	Part number	Qty	Color	Volume
1	UDI Tagging Reagent Plate	301308	1	Green	12 µL per well
2	Library Primer Mix (LPM)	300134	1	Clear	700 µL

## Quantities

96 reactions

## Storage and handling requirements

Product should be stored frozen between -15 to -25 °C.

Stable for up to 12 freeze/thaw cycles.

Packaging information:

- Type: Box
- Dimensions: 5.3"L X 3.0"H X 3.5"W

Refer to protocols for additional recommendations for product handling.

## Shelf life

Refer to expiration date on product label.

## Intended use

Research use only. Not for use in diagnostic procedures.

## Support

For additional product and technical information, such as safety data sheets (SDS) or procedures, visit our [website](#). For region-specific SDS, search the [SDS library](#).

## Legal

Information in this document is subject to change without notice. PacBio assumes no responsibility for any errors or omissions in this document. Certain notices, terms, conditions and/or use restrictions may pertain to your use of PacBio products and/or third-party products. Refer to the applicable PacBio terms and conditions of sale and to the applicable license terms at [pacb.com/license](http://pacb.com/license). Pacific Biosciences, the PacBio logo, PacBio, Circulomics, Omniome, SMRT, SMRTbell, Iso-Seq, Sequel, Nanobind, SBB, Revio, Onso, Apton, Kinnex, PureTarget, SPRQ, and Vega are trademarks of PacBio. LongPlex is a trademark of seqWell.