

### Description

The MagBead Kit v2 contains reagent supplies to capture DNA polymerase/template complexes for loading onto the PacBio® System.

The MagBead Kit v2 (P/N 100-676-500) includes:

- MagBeads (P/N 100-125-900)

Subcomponent Name	Part Number	Number	Color	Quantity
MagBeads	100-125-900	1 bottle	N/A	2.88 mL

- MagBead Binding Buffer Kit v2 (P/N 101-046-400)

Subcomponent Name	Part Number	Number	Color	Quantity
MagBead Binding Buffer v2	100-672-100	9	Red	1.7 mL

- MagBead Wash Buffer Kit v2 (P/N 101-046-500)

Subcomponent Name	Part Number	Number	Color	Quantity
MagBead Wash Buffer v2	100-672-200	9	Yellow	1.7 mL

### Part Number

The MagBead Kit v2 (P/N 100-676-500) is shipped as three separate parts: MagBeads (P/N 100-125-900), MagBead Binding Buffer Kit v2 (P/N 101-046-400), and MagBead Wash Buffer Kit v2 (P/N: 101-046-500).

### Quantity

Refer to the Pacific Biosciences Binding Calculator for reaction numbers per kit.

### Intended Use

For Research Use Only. Not for use in diagnostic procedures.

Certain notices, terms, conditions and/or use restrictions may pertain to your use of Pacific Biosciences® products and/or third party products. Please refer to the applicable Pacific Biosciences Terms and Conditions of Sale and to the applicable license terms at <http://www.pacb.com/licenses.html>

### Storage and Handling Requirements

MagBeads (P/N 100-125-900) should be stored refrigerated between 2-8 °C upon arrival. DO NOT FREEZE.

MagBead Binding Buffer Kit v2 (P/N: 101-046-400) and MagBead Wash Buffer Kit v2 (P/N: 101-046-500) should be stored frozen between -15 to -25 °C.

### Shelf Life

Refer to expiration dates on product labels.

### Support

For additional product and technical information, such as Safety Data Sheets (SDS) or procedures, please visit our website at [www.pacificbiosciences.com](http://www.pacificbiosciences.com) or our customer portal at [www.pacbioportal.com](http://www.pacbioportal.com).

### Legal

© Copyright 2010-2017, Pacific Biosciences of California, Inc. All rights reserved.

Information in this document is subject to change without notice. Pacific Biosciences assumes no responsibility for any errors or omissions in this document.

PACIFIC BIOSCIENCES DISCLAIMS ALL WARRANTIES WITH RESPECT TO THIS DOCUMENT, EXPRESS, STATUTORY, IMPLIED OR OTHERWISE, INCLUDING, BUT NOT LIMITED TO, ANY WARRANTIES OF MERCHANTABILITY, SATISFACTORY QUALITY, NONINFRINGEMENT OR FITNESS FOR A PARTICULAR PURPOSE. IN NO EVENT SHALL PACIFIC BIOSCIENCES BE LIABLE, WHETHER IN CONTRACT, TORT, WARRANTY, PURSUANT TO ANY STATUTE, OR ON ANY OTHER BASIS FOR SPECIAL, CONSEQUENTIAL, INCIDENTAL, EXEMPLARY OR INDIRECT DAMAGES IN CONNECTION WITH (OR ARISING FROM) THIS DOCUMENT, WHETHER OR NOT FORESEEABLE AND WHETHER OR NOT PACIFIC BIOSCIENCES IS ADVISED OF THE POSSIBILITY OF SUCH DAMAGES.

Trademarks:

Pacific Biosciences, the Pacific Biosciences logo, PacBio, SMRT, SMRTbell, Iso-Seq, and Sequel are trademarks of Pacific Biosciences in the United States and/or certain other countries.

All other trademarks are the sole property of their respective owners.