

Description

The Nanobind PanDNA kit 4C contains reagents to be used in conjunction with the Nanobind PanDNA kit RT (103-260-300) to complete HMW DNA extraction from cells, blood, saliva, bacteria, tissues, insects, and plants.

Kit configuration

Configuration

#	Component	Part number	Qty	Color	Volume
1	RNase A	102-278-800	1	blue	0.5 mL
2	Buffer CT	102-280-300	1	clear	42 mL
3	Buffer RBC 10X	103-232-000	1	clear	5 mL

Quantities

24 reactions

Storage and handling requirements

Product should be stored refrigerated between 2 to 8 °C.

Refer to protocols for additional recommendations for product handling.

Shelf life

Refer to expiration date on product label.

Intended use

Research use only. Not for use in diagnostic procedures.

Support

For additional product and technical information, such as safety data sheets (SDS) or procedures, visit our [website](#).

Legal

Information in this document is subject to change without notice. PacBio assumes no responsibility for any errors or omissions in this document. Certain notices, terms, conditions and/or use restrictions may pertain to your use of PacBio products and/or third-party products. Refer to the applicable PacBio terms and conditions of sale and to the applicable license terms at [pacb.com/license](https://www.pacb.com/license). Pacific Biosciences, the PacBio logo, PacBio, Circulomics, Omniome, SMRT, SMRTbell, Iso-Seq, Sequel, Nanobind, SBB, Revio, Onso, Apton, Kinnex, PureTarget, SPRQ, and Vega are trademarks of PacBio.