

Description

No-Amp Accessory Kit contains Annealing Buffer, Sequencing Primer v4, and 10X Primer Buffer v2. Annealing buffer is to be used during SMRTbell® Adapter hairpin annealing. Sequencing Primer v4 and 10X Primer Buffer v2 is used to anneal Sequencing Primer v4 to final SMRTbell.

Kit components include:

Subcomponent Name	Part Number	Number	Color	Quantity (µl)
10X Annealing Buffer	101-789-000	2 tubes	Light Green	1000 µl
Sequencing Primer v4	101-654-600	1 tube	Light Blue	18 µl
10X Primer Buffer v2	001-560-849	1 tube	Blue	1000 µl

Part Number

101-788-900

Quantity

1 Kit

Intended Use

For Research Use Only. Not for use in diagnostic procedures.

Certain notices, terms, conditions and/or use restrictions may pertain to your use of Pacific Biosciences® products and/or third-party products. Please refer to the applicable Pacific Biosciences Terms and Conditions of Sale and to the applicable license terms at <https://www.pacb.com/legal-and-trademarks/terms-and-conditions-of-sale>.

Storage and Handling Requirements

Product should be stored frozen between -15 to -25 °C.

Refer to Sample Preparation and Sequencing Guide for additional recommendations for reagent handling.

Shelf Life

Refer to expiration date on kit label.

Support

For additional product and technical information, such as Safety Data Sheets (SDS) or procedures, please visit our website at <http://www.pacb.com/support/documentation/>

Legal

© Copyright 2010-2019, Pacific Biosciences of California, Inc. All rights reserved.

Information in this document is subject to change without notice. Pacific Biosciences assumes no responsibility for any errors or omissions in this document.

PACIFIC BIOSCIENCES DISCLAIMS ALL WARRANTIES WITH RESPECT TO THIS DOCUMENT, EXPRESS, STATUTORY, IMPLIED OR OTHERWISE, INCLUDING, BUT NOT LIMITED TO, ANY WARRANTIES OF MERCHANTABILITY, SATISFACTORY QUALITY, NONINFRINGEMENT OR FITNESS FOR A PARTICULAR PURPOSE. IN NO EVENT SHALL PACIFIC BIOSCIENCES BE LIABLE, WHETHER IN CONTRACT, TORT, WARRANTY, PURSUANT TO ANY STATUTE, OR ON ANY OTHER BASIS FOR SPECIAL, CONSEQUENTIAL, INCIDENTAL, EXEMPLARY OR INDIRECT DAMAGES IN CONNECTION WITH (OR ARISING FROM) THIS DOCUMENT, WHETHER OR NOT FORESEEABLE AND WHETHER OR NOT PACIFIC BIOSCIENCES IS ADVISED OF THE POSSIBILITY OF SUCH DAMAGES.

Trademarks:

Pacific Biosciences, the Pacific Biosciences logo, PacBio, SMRT, SMRTbell, Iso-Seq, and Sequel are trademarks of Pacific Biosciences in the United States and/or certain other countries.

All other trademarks are the sole property of their respective owners.