Onso[™] DNA library prep kit

Description

Library prep kit used for creation of Onso-compatible libraries from fragmented DNA. Contains reagents for end-repair, A-tailing, and adapter ligation of DNA fragments (96 reactions).

Kit configuration & layout

Configuration

#	Component	Part number	Qty	Color	Volume
1	End prep mix	102-554-200	1	Orange	288 µL
2	End prep buffer	102-554-300	1	green	768 µL
3	Ligation mix	102-554-400	1	Yellow	960 µL
4	Ligation mix	102-554-400	1	Yellow	960 µL
5	Ligation mix	102-554-400	1	Yellow	960µL
6	Ligation enhancer	102-554-500	1	Red	96 µL
7	EDTA	102-554-600	1	White	1150µL
8	EDTA	102-554-600	1	White	1150µL

Layout



Quantities

96 reactions

Storage and handling requirements

Product should be stored frozen between -15 to -25 °C.

Refer to protocols for additional recommendations for reagent handling.

Shelf life

Refer to expiration date on kit label.

Intended use

Research use only. Not for use in diagnostic procedures.

Support

For additional product and technical information, such as safety data sheets (SDS) or procedures, visit our <u>website</u>.

Legal

Information in this document is subject to change without notice. PacBio assumes no responsibility for any errors or omissions in this document. Certain notices, terms, conditions and/or use restrictions may pertain to your use of PacBio products and/or third-party products. Refer to the applicable PacBio terms and conditions of sale and to the applicable license terms at <u>pacb.com/license</u>. Pacific Biosciences, the PacBio logo, PacBio, Circulomics, Omniome, SMRT, SMRTbell, Iso-Seq, Sequel, Nanobind, SBB, Revio, and Onso are trademarks of PacBio.

