

## Description

Onso cluster reagent plate contains all necessary reagents to dilute and denature DNA libraries, and enables on flow cell cluster generation.

## Kit configuration

#	Component	Part number	Qty	Color	Volume
1	Tube, Dil buffer B	102-903-700	1	white	1000 µL
2	Onso cluster reagent plate	102-752-200	1	N/A	1 Qty

## Quantities

1 Qty

## Storage and handling requirements

Product should be stored frozen between -15 to -25 °C.

Refer to protocols for additional recommendations for reagent handling.

## Shelf life

Refer to expiration date on kit label.

## Intended use

Research use only. Not for use in diagnostic procedures.

## Support

For additional product and technical information, such as safety data sheets (SDS) or procedures, visit our [website](#).

## Legal

Research use only. Not for use in diagnostic procedures. © 2023 Pacific Biosciences of California, Inc. ("PacBio"). Information in this document is subject to change without notice. PacBio assumes no responsibility for any errors or omissions in this document. Certain notices, terms, conditions and/or use restrictions may pertain to your use of PacBio products and/or third-party products. Refer to the applicable PacBio terms and conditions of sale and to the applicable license terms at [pacb.com/license](https://pacb.com/license). Pacific Biosciences, the PacBio logo, PacBio, Circulomics, Omniome, SMRT, SMRTbell, Iso-Seq, Sequel, Nanobind, SBB, Revio and Onso are trademarks of PacBio.