

Onso™ sequencing reagent pack, 200 cycles

Product PN: 102-656-400

Description

The Onso sequencing reagent pack, 200 cycles, contains sequencing reagents to enable up to 200 sequencing-by-binding (SBB) cycles on the Onso system. Includes support for single- and dual-index reads.

Kit configuration & layout

Configuration

#	Component	Part number	Qty	Color	Volume
1	Buffer C	102-379-600	1	N/A	300 mL
2	Incorporation Buffer S	102-807-400	1	N/A	135 mL
3	Buffer A	102-385-100	1	N/A	155 mL
4	Imaging Buffer	102-384-600	1	N/A	300 mL
5	Stripping Buffer	102-384-700	1	N/A	300 mL
6	Buffer C	102-379-600	1	N/A	300 mL
7	Incorporation Buffer S	102-807-400	1	N/A	135 mL
8	Buffer B	102-385-200	1	N/A	155 mL
9	QC Buffer	102-384-800	1	N/A	250 mL
10	Stripping Buffer	102-384-700	1	N/A	300 mL
11	Change Strand Reagent	102-385-900	1	N/A	4 mL
12	Blocking Buffer	102-385-300	1	N/A	6 mL
13	Scale Insert Primer	102-385-800	1	N/A	4 mL
14	Index Primer 2	102-386-000	1	N/A	4 mL
15	Insert Primer 1	102-385-500	1	N/A	4 mL
16	Index Primer 1	102-385-600	1	N/A	4 mL
17-20	Primer Buffer	102-385-400	4	N/A	4 mL

Layout



Quantities

200 cycles

Storage and handling requirements

Product should be stored frozen between -15 to -25 °C.

Refer to protocols for additional recommendations for reagent handling.

Shelf life

Refer to expiration date on kit label.

Intended use

Research use only. Not for use in diagnostic procedures.

Support

For additional product and technical information, such as safety data sheets (SDS) or procedures, visit our [website](#).

Legal

Information in this document is subject to change without notice. PacBio assumes no responsibility for any errors or omissions in this document. Certain notices, terms, conditions and/or use restrictions may pertain to your use of PacBio products and/or third-party products. Refer to the applicable PacBio terms and conditions of sale and to the applicable license terms at pacb.com/license. Pacific Biosciences, the PacBio logo, PacBio, Circulomics, Omniome, SMRT, SMRTbell, Iso-Seq, Sequel, Nanobind, SBB, Revio, and Onso are trademarks of PacBio.