

## Description

The Revio SMRT Cell tray - Nx contains 4 Revio SMRT Cells - Nx for sequencing on the Revio system.

## Kit configuration

### Configuration

#	Component	Part number	Qty	Color	Volume
1	Revio SMRT Cell - Nx	N/A	1	N/A	N/A
2	Revio SMRT Cell - Nx	N/A	1	N/A	N/A
3	Revio SMRT Cell - Nx	N/A	1	N/A	N/A
4	Revio SMRT Cell - Nx	N/A	1	N/A	N/A

## Quantities

1 SMRT Cell tray - Nx containing 4 SMRT Cells - Nx

## Storage and handling requirements

During shipment, product should be stored between -5 to 50°C.

Once delivered, product should be stored at ambient temperature between 8 to 40°C.

Packaging information:

- Type: Box
- Dimensions (inches): 8.1" x 1.5" x 2.3"

Refer to protocols for additional recommendations for product handling.

## Shelf life

Refer to expiration date on product label.

## Intended use

Research use only. Not for use in diagnostic procedures.

## Support

For additional product and technical information, such as safety data sheets (SDS) or procedures, visit our [website](#). For region-specific SDS, search the [SDS library](#).

## Legal

Information in this document is subject to change without notice. PacBio assumes no responsibility for any errors or omissions in this document. Certain notices, terms, conditions and/or use restrictions may pertain to your use of PacBio products and/or third-party products. Refer to the applicable PacBio terms and conditions of sale and to the applicable license terms at [pacb.com/license](http://pacb.com/license). Pacific Biosciences, the PacBio logo, PacBio, Circulomics, Omniome, SMRT, SMRTbell, Iso-Seq, Sequel, Nanobind, SBB, Revio, Onso, Apton, Kinnex, PureTarget, SPRQ, and Vega are trademarks of PacBio.