

## Description

Reagents for binding sequencing polymerase to SMRTbell® library. Recommended for SMRTbell libraries ≥3 kb. Supports 24 Revio SMRT® Cells.

## Kit configuration & layout

### Configuration

#	Component	Part number	Qty	Color	Volume
1	Annealing buffer	102-797-600	1	light blue	400 µL
2	Sequencing primer	102-797-700	1	light green	400 µL
3	Polymerase buffer	102-797-800	1	yellow	1.8 mL
4	Sequencing polymerase	102-797-300	1	purple	35 µL
5-7	Dilution buffer	102-797-900	3	blue	1.8 mL
8	Sequencing control	102-798-000	1	red	35 µL
9-10	Loading buffer	102-797-500	2	green	1.6 mL

### Layout



## Quantities

24 reactions

## Storage and handling requirements

Product should be stored frozen between -15°C to -25°C.

Refer to protocols for additional recommendations for reagent handling.

## Shelf life

Refer to expiration date on kit label.

## Intended use

Research use only. Not for use in diagnostic procedures.

## Support

For additional product and technical information, such as safety data sheets (SDS) or procedures, visit our [website](#).

## Legal

Research use only. Not for use in diagnostic procedures. © 2023 Pacific Biosciences of California, Inc. ("PacBio"). Information in this document is subject to change without notice. PacBio assumes no responsibility for any errors or omissions in this document. Certain notices, terms, conditions and/or use restrictions may pertain to your use of PacBio products and/or third-party products. Refer to the applicable PacBio terms and conditions of sale and to the applicable license terms at [pacb.com/license](https://pacb.com/license). Pacific Biosciences, the PacBio logo, PacBio, Circulomics, Omniome, SMRT, SMRTbell, Iso-Seq, Sequel, Nanobind, SBB, Revio, and Onso are trademarks of PacBio.