

Description

The SMRTbell™ Clean Up Column Kit (PN 100-929-600) contains reagent supplies to assist in the removal of free sequencing polymerase for up to 5 samples. The SMRTbell™ Clean Up Column Kit consists of 1 of each of the following parts: SMRTbell™ Clean Up Columns (PN: 100-912-000) and SMRTbell™ Conditioning Buffer Kit (PN: 100-912-100). The SMRTbell Conditioning Buffer Kit is used to perform a buffer exchange of the native buffer in the Clean Up Column for a polymerase-compatible buffer.

SMRTbell™ Clean Up Kit (100-929-600) includes:

- SMRTbell™ Clean Up Columns (100-912-000)

Subcomponent Name	Number	Container
SMRTbell™ Clean Up Columns	5	Bag
Collection Tubes	10	Bag
Vial Stoppers	5	Bag

- SMRTbell™ Conditioning Buffer Kit (100-912-100)

Subcomponent Name	Part Number	Number	Color
MagBead Wash Buffer v2	100-672-200	9	Yellow

Part Number

The SMRTbell™ Clean Up Column Kit (PN: 100-929-600) is shipped as two separate parts: SMRTbell™ Clean Up Columns (PN: 100-912-000) and SMRTbell™ Conditioning Buffer Kit (PN: 100-912-100).

Quantities

5 reactions

Intended Use

For Research Use Only. Not for use in diagnostic procedures.

Certain notices, terms, conditions and/or use restrictions may pertain to your use of Pacific Biosciences® products and/or third party products. Please refer to the applicable Pacific Biosciences Terms and Conditions of Sale and to the applicable license terms at <http://www.pacb.com/licenses.html>

Storage and Handling Requirements

SMRTbell™ Clean Up Columns (PN: 100-912-000) should be stored between 15 to 30 °C. DO NOT FREEZE.
SMRTbell™ Conditioning Buffer Kit (PN: 100-912-100) should be stored frozen between -15 to -25 °C

Shelf Life

Refer to expiration dates on product labels.

Support

For additional product and technical information, such as Safety Data Sheets (SDS) or procedures, please visit our website at <http://www.pacb.com/support/documentation/>.

Legal

© Copyright 2010-2016, Pacific Biosciences of California, Inc. All rights reserved.

Information in this document is subject to change without notice. Pacific Biosciences assumes no responsibility for any errors or omissions in this document.

PACIFIC BIOSCIENCES DISCLAIMS ALL WARRANTIES WITH RESPECT TO THIS DOCUMENT, EXPRESS, STATUTORY, IMPLIED OR OTHERWISE, INCLUDING, BUT NOT LIMITED TO, ANY WARRANTIES OF MERCHANTABILITY, SATISFACTORY QUALITY, NONINFRINGEMENT OR FITNESS FOR A PARTICULAR PURPOSE. IN NO EVENT SHALL PACIFIC BIOSCIENCES BE LIABLE, WHETHER IN CONTRACT, TORT, WARRANTY, PURSUANT TO ANY STATUTE, OR ON ANY OTHER BASIS FOR SPECIAL, CONSEQUENTIAL, INCIDENTAL, EXEMPLARY OR INDIRECT DAMAGES IN CONNECTION WITH (OR ARISING FROM) THIS DOCUMENT, WHETHER OR NOT FORESEEABLE AND WHETHER OR NOT PACIFIC BIOSCIENCES IS ADVISED OF THE POSSIBILITY OF SUCH DAMAGES.

Trademarks:

Pacific Biosciences, the Pacific Biosciences logo, PacBio, SMRT, SMRTbell, Iso-Seq, and Sequel are trademarks of Pacific Biosciences in the United States and/or certain other countries.

All other trademarks are the sole property of their respective owners.