

Description

The Sequel II SMRT Cell ctrl complex 1.0 contains Lambda library DNA template sheared to a peak size of 2 kb and bound to Sequel II polymerase 1.0. It acts as a complete Sequel II SMRT Cell control.

Kit configuration & layout

Configuration

#	Component	Part number	Qty	Color	Volume
1	Sequel II SMRT Cell control complex 1.0	101-756-900	1	purple	30 µL
2-4	Sequel complex dilution buffer	101-623-600	3	blue	1600 µL

Layout



Quantities

30 µL

Storage and handling requirements

Product should be stored frozen between -15 to -25 °C.

Refer to protocols for additional recommendations for reagent handling.

Shelf life

Refer to expiration date on kit label.

Intended use

Research use only. Not for use in diagnostic procedures.

Support

For additional product and technical information, such as safety data sheets (SDS) or procedures, visit our [website](#).

Legal

Information in this document is subject to change without notice. PacBio assumes no responsibility for any errors or omissions in this document. Certain notices, terms, conditions and/or use restrictions may pertain to your use of PacBio products and/or third-party products. Refer to the applicable PacBio terms and conditions of sale and to the applicable license terms at pacb.com/license. Pacific Biosciences, the PacBio logo, PacBio, Circulomics, Omniome, SMRT, SMRTbell, Iso-Seq, Sequel, Nanobind, and SBB are trademarks of PacBio. All other trademarks are the sole property of their respective owners.