

Description

This buffer consists of 10 mM Tris-HCl (pH 8.0) and 0.1 mM EDTA.

Subcomponent Name	Part Number	Number	Fill Volume (mL)
Low TE Buffer	102-178-400	1 bottle	10 mL

Part Number

102-178-400

Quantity

10 mL

Intended Use

For Research Use Only. Not for use in diagnostic procedures.

Certain notices, terms, conditions and/or use restrictions may pertain to your use of Pacific Biosciences® products and/or third-party products. Please refer to the applicable Pacific Biosciences Terms and Conditions of Sale and to the applicable license terms at <https://www.pacb.com/legal-and-trademarks/terms-and-conditions-of-sale>.

Storage and Handling Requirements

Product should be stored at room temperature (15-30°C) upon receipt.

Shelf Life

Refer to expiration date on kit label.

Support

For additional product and technical information, such as Safety Data Sheets (SDS) or procedures, please visit our website at <http://www.pacb.com/support/documentation/>

Legal

© Copyright 2021, Pacific Biosciences of California, Inc. All rights reserved.

Information in this document is subject to change without notice. Pacific Biosciences assumes no responsibility for any errors or omissions in this document.

PACIFIC BIOSCIENCES DISCLAIMS ALL WARRANTIES WITH RESPECT TO THIS DOCUMENT, EXPRESS, STATUTORY, IMPLIED OR OTHERWISE, INCLUDING, BUT NOT LIMITED TO, ANY WARRANTIES OF MERCHANTABILITY, SATISFACTORY QUALITY, NONINFRINGEMENT OR FITNESS FOR A PARTICULAR PURPOSE. IN NO EVENT SHALL PACIFIC BIOSCIENCES BE LIABLE, WHETHER IN CONTRACT, TORT, WARRANTY, PURSUANT TO ANY STATUTE, OR ON ANY OTHER BASIS FOR SPECIAL, CONSEQUENTIAL, INCIDENTAL, EXEMPLARY OR INDIRECT DAMAGES IN CONNECTION WITH (OR ARISING FROM) THIS DOCUMENT, WHETHER OR NOT FORESEEABLE AND WHETHER OR NOT PACIFIC BIOSCIENCES IS ADVISED OF THE POSSIBILITY OF SUCH DAMAGES.

Trademarks:

Pacific Biosciences, the Pacific Biosciences logo, PacBio, SMRT, SMRTbell, Iso-Seq, and Sequel are trademarks of Pacific Biosciences in the United States and/or certain other countries.

All other trademarks are the sole property of their respective owners.