
This system version supports Sequel II Systems **only**.

- | | |
|---------------------|--|
| Improvements | Sequel® Instrument Control Software (ICS) <ul style="list-style-type: none">• Various automation and reliability improvements. |
| Fixed Issues | Sequel® Instrument Control Software (ICS) <ul style="list-style-type: none">• The context folders are no longer missing the <code>LaserAlign_Details.csv</code> and <code>LaserAlign_Summary.csv</code> report files.• Fixed an issue where runs yielding very large primary data volumes (600 Gb and above) would occasionally fail Post-Primary Analysis due to system resource limits. However, Post-Primary Analysis may take longer than usual for such runs. |
| Known Issues | <ul style="list-style-type: none">• When performing a run using two Sequencing Plates, each plate must have a separate Sequel SMRT Cell Oil tube. |

For Research Use Only. Not for use in diagnostic procedures. © Copyright 2019-2020, Pacific Biosciences of California, Inc. All rights reserved. Information in this document is subject to change without notice. Pacific Biosciences assumes no responsibility for any errors or omissions in this document. Certain notices, terms, conditions and/or use restrictions may pertain to your use of Pacific Biosciences products and/or third party products. Please refer to the applicable Pacific Biosciences Terms and Conditions of Sale and the applicable license terms at <http://www.pacb.com/legal-and-trademarks/product-license-and-use-restrictions/>.

Pacific Biosciences, the Pacific Biosciences logo, PacBio, SMRT, SMRTbell, Iso-Seq and Sequel are trademarks of Pacific Biosciences. BluePippin and SageELF are trademarks of Sage Science, Inc. NGS-go and NGSengine are trademarks of GenDx. FEMTO Pulse and Fragment Analyzer are trademarks of Agilent Technologies Inc. All other trademarks are the sole property of their respective owners.

P/N 101-913-800 Version 01 (April 2020)