

## New features

- Output FASTQ files now use the following file naming convention to improve compatibility with downstream secondary analysis tools:

```
<Sample Name>_S<Sample Number>_L00<Lane Number>_R<Read Number>_001.fastq.gz
```

**Note:** Contact your FSE (Field Service Engineer) or PacBio Technical Support to perform the Instrument Control Software installation.

Research use only. Not for use in diagnostic procedures. © 2025 Pacific Biosciences of California, Inc. (“PacBio”). All rights reserved. Information in this document is subject to change without notice. PacBio assumes no responsibility for any errors or omissions in this document. Certain notices, terms, conditions and/or use restrictions may pertain to your use of PacBio products and/or third-party products. Refer to the applicable PacBio terms and conditions of sale and to the applicable license terms at [pacb.com/license](https://pacb.com/license).

Pacific Biosciences, the PacBio logo, PacBio, Circulomics, Omniome, SMRT, SMRTbell, Iso-Seq, Sequel, Nanobind, SBB, Revio, Onso, Apton, Kinnex, PureTarget, SPRQ, and Vega are trademarks of PacBio.

© 2025 PacBio. All rights reserved. Research use only. Not for use in diagnostic procedures.

PN 103-617-200 REV01 (January 2025)